

(12) INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY

(19) World Intellectual Property Organization  
International Bureau



(43) International Publication Date  
15 March 2001 (15.03.2001)

PCT

(10) International Publication Number  
**WO 01/18542 A2**

(51) International Patent Classification<sup>7</sup>: **G01N 33/574**

One Kendall Square Bldg. 700, Cambridge, MA 02139 (US).

(21) International Application Number: **PCT/US00/24199**

(72) Inventors; and

(22) International Filing Date:

1 September 2000 (01.09.2000)

(75) Inventors/Applicants (for US only): **LEE, John** [US/US]; 119 Walnut Street, Somerville, MA 02145 (US). **THOMPSON, Pamela** [US/US]; 83 Beech Street #2, Belmont, MA 02478 (US). **LILLIE, James** [US/US]; 119 Walnut Street, Somerville, MA 02145 (US).

(25) Filing Language:

English

(26) Publication Language:

English

(74) Agent: **SMITH, DeAnn, F.**; Lahive & Cockfield, LLP, 28 State Street, Boston, MA 02109 (US).

(30) Priority Data:

60/152,547	3 September 1999 (03.09.1999)	US
60/190,347	16 March 2000 (16.03.2000)	US
60/191,321	21 March 2000 (21.03.2000)	US
60/208,382	31 May 2000 (31.05.2000)	US
60/220,467	20 July 2000 (20.07.2000)	US

(81) Designated States (*national*): AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CR, CU, CZ, DE, DK, DM, DZ, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, TZ, UA, UG, US, UZ, VN, YU, ZA, ZW.

(63) Related by continuation (CON) or continuation-in-part (CIP) to earlier applications:

US	60/152,547 (CON)
Filed on	3 September 1999 (03.09.1999)
US	60/190,347 (CON)
Filed on	16 March 2000 (16.03.2000)
US	60/191,321 (CON)
Filed on	21 March 2000 (21.03.2000)
US	60/208,382 (CON)
Filed on	31 May 2000 (31.05.2000)
US	60/220,467 (CON)
Filed on	20 July 2000 (20.07.2000)

(84) Designated States (*regional*): ARIPO patent (GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG).

Published:

— Without international search report and to be republished upon receipt of that report.

(71) Applicant (for all designated States except US): **MILLENNIUM PREDICTIVE MEDICINE, INC.** [US/US];

For two-letter codes and other abbreviations, refer to the "Guidance Notes on Codes and Abbreviations" appearing at the beginning of each regular issue of the PCT Gazette.

(54) Title: COMPOSITIONS, KITS, AND METHODS FOR IDENTIFICATION, ASSESSMENT, PREVENTION, AND THERAPY OF OVARIAN CANCER

(57) Abstract: The invention relates to compositions, kits, and methods for detecting, characterizing, preventing, and treating human ovarian cancers. A variety of markers are provided, wherein changes in the levels of expression of one or more of the markers is correlated with the presence of ovarian cancer.



WO 01/18542 A2

**COMPOSITIONS, KITS, AND METHODS FOR  
IDENTIFICATION, ASSESSMENT, PREVENTION, AND THERAPY OF  
OVARIAN CANCER**

5

**RELATED APPLICATIONS**

The present application claims priority to U.S. provisional patent application serial no. 60/152,547, filed on September 3, 1999, U.S. provisional patent application serial no. 60/190,347, filed on March 16, 2000, U.S. provisional patent application serial no. 60/191,321, filed on March 21, 2000, U.S. provisional patent application serial no. 60/208,382, filed on May 31, 2000 and U.S. provisional patent application serial no. 60/220,467, filed on July 20, 2000, all of which are expressly incorporated by reference.

**FIELD OF THE INVENTION**

The field of the invention is ovarian cancer, including diagnosis,  
15 characterization, management, and therapy of ovarian cancer.

**BACKGROUND OF THE INVENTION**

Ovarian cancer is responsible for significant morbidity and mortality in populations around the world. Ovarian cancer is classified, on the basis of clinical and pathological features, in three groups, namely epithelial ovarian cancer (EOC; >90% of  
20 ovarian cancer in Western countries), germ cell tumors (*circa* 2-3% of ovarian cancer), and stromal ovarian cancer (*circa* 5% of ovarian cancer; Ozols *et al.*, 1997, *Cancer Principles and Practice of Oncology*, 5th ed., DeVita *et al.*, Eds. pp. 1502). Relative to EOC, germ cell tumors and stromal ovarian cancers are more easily detected and treated  
25 at an early stage, translating into higher/better survival rates for patients afflicted with these two types of ovarian cancer.

There are numerous types of ovarian tumors, some of which are benign, and others of which are malignant. Treatment (including non-treatment) options and predictions of patient outcome depend on accurate classification of the ovarian cancer.  
30 Ovarian cancers are named according to the type of cells from which the cancer is derived and whether the ovarian cancer is benign or malignant. Recognized histological tumor types include, for example, serous, mucinous, endometrioid, and clear cell

tumors. In addition, ovarian cancers are classified according to recognized grade and stage scales.

In grade I, the tumor tissue is well differentiated from normal ovarian tissue. In grade II, tumor tissue is moderately well differentiated. In grade III, the tumor tissue is poorly differentiated from normal tissue, and this grade correlates with a less favorable prognosis than grades I and II. Stage I is generally confined within the capsule surrounding one (stage IA) or both (stage IB) ovaries, although in some stage I (*i.e.* stage IC) cancers, malignant cells may be detected in ascites, in peritoneal rinse fluid, or on the surface of the ovaries. Stage II involves extension or metastasis of the tumor from one or both ovaries to other pelvic structures. In stage IIA, the tumor extends or has metastasized to the uterus, the fallopian tubes, or both. Stage IIB involves extension of the tumor to the pelvis. Stage IIC is stage IIA or IIB in which malignant cells may be detected in ascites, in peritoneal rinse fluid, or on the surface of the ovaries. In stage III, the tumor comprises at least one malignant extension to the small bowel or the omentum, has formed extrapelvic peritoneal implants of microscopic (stage IIIA) or macroscopic (< 2 centimeter diameter, stage IIIB; > 2 centimeter diameter, stage IIIC) size, or has metastasized to a retroperitoneal or inguinal lymph node (an alternate indicator of stage IIIC). In stage IV, distant (*i.e.* non-peritoneal) metastases of the tumor can be detected.

The durations of the various stages of ovarian cancer are not presently known, but are believed to be at least about a year each (Richart *et al.*, 1969, *Am. J. Obstet. Gynecol.* 105:386). Prognosis declines with increasing stage designation. For example, 5-year survival rates for patients diagnosed with stage I, II, III, and IV ovarian cancer are 80%, 57%, 25%, and 8%, respectively.

Despite being the third most prevalent gynecological cancer, ovarian cancer is the leading cause of death among those afflicted with gynecological cancers. The disproportionate mortality of ovarian cancer is attributable to a substantial absence of symptoms among those afflicted with early-stage ovarian cancer and to difficulty diagnosing ovarian cancer at an early stage. Patients afflicted with ovarian cancer most often present with non-specific complaints, such as abnormal vaginal bleeding, gastrointestinal symptoms, urinary tract symptoms, lower abdominal pain, and generalized abdominal distension. These patients rarely present with paraneoplastic

symptoms or with symptoms which clearly indicate their affliction. Presently, less than about 40% of patients afflicted with ovarian cancer present with stage I or stage II.

Management of ovarian cancer would be significantly enhanced if the disease could be detected at an earlier stage, when treatments are much more generally efficacious.

5           Ovarian cancer may be diagnosed, in part, by collecting a routine medical history from a patient and by performing physical examination, x-ray examination, and chemical and hematological studies on the patient. Hematological tests which may be indicative of ovarian cancer in a patient include analyses of serum levels of proteins designated CA125 and DF3 and plasma levels of lysophosphatidic acid (LPA).

10          Palpation of the ovaries and ultrasound techniques (particularly including endovaginal ultrasound and color Doppler flow ultrasound techniques) can aid detection of ovarian tumors and differentiation of ovarian cancer from benign ovarian cysts. However, a definitive diagnosis of ovarian cancer typically requires performing exploratory laparotomy of the patient.

15           Potential tests for the detection of ovarian cancer (*e.g.*, screening, reflex or monitoring) may be characterized by a number of factors. The "sensitivity" of an assay refers to the probability that the test will yield a positive result in an individual afflicted with ovarian cancer. The "specificity" of an assay refers to the probability that the test will yield a negative result in an individual not afflicted with ovarian cancer. The  
20          "positive predictive value" (PPV) of an assay is the ratio of true positive results (*i.e.* positive assay results for patients afflicted with ovarian cancer) to all positive results (*i.e.* positive assay results for patients afflicted with ovarian cancer + positive assay results for patients not afflicted with ovarian cancer). It has been estimated that in order for an assay to be an appropriate population-wide screening tool for ovarian cancer the  
25          assay must have a PPV of at least about 10% (Rosenthal *et al.*, 1998, *Sem. Oncol.* 25:315-325). It would thus be desirable for a screening assay for detecting ovarian cancer in patients to have a high sensitivity and a high PPV. Monitoring and reflex tests would also require appropriate specifications.

            Owing to the cost, limited sensitivity, and limited specificity of known methods  
30          of detecting ovarian cancer, screening is not presently performed for the general population. In addition, the need to perform laparotomy in order to diagnose ovarian cancer in patients who screen positive for indications of ovarian cancer limits the



desirability of population-wide screening, such that a PPV even greater than 10% would be desirable.

Prior use of serum CA125 level as a diagnostic marker for ovarian cancer indicated that this method exhibited insufficient specificity for use as a general  
5 screening method. Use of a refined algorithm for interpreting CA125 levels in serial retrospective samples obtained from patients improved the specificity of the method without shifting detection of ovarian cancer to an earlier stage (Skakes, 1995, *Cancer* 76:2004). Screening for LPA to detect gynecological cancers including ovarian cancer exhibited a sensitivity of about 96% and a specificity of about 89%. However, CA125-  
10 based screening methods and LPA-based screening methods are hampered by the presence of CA125 and LPA, respectively, in the serum of patients afflicted with conditions other than ovarian cancer. For example, serum CA125 levels are known to be associated with menstruation, pregnancy, gastrointestinal and hepatic conditions such as colitis and cirrhosis, pericarditis, renal disease, and various non-ovarian malignancies.  
15 Serum LPA is known, for example, to be affected by the presence of non-ovarian gynecological malignancies. A screening method having a greater specificity for ovarian cancer than the current screening methods for CA125 and LPA could provide a population-wide screening for early stage ovarian cancer.

Presently greater than about 60% of ovarian cancers diagnosed in patients are  
20 stage III or stage IV cancers. Treatment at these stages is largely limited to cytoreductive surgery (when feasible) and chemotherapy, both of which aim to slow the spread and development of metastasized tumor. Substantially all late stage ovarian cancer patients currently undergo combination chemotherapy as primary treatment, usually a combination of a platinum compound and a taxane. Median survival for  
25 responding patients is about one year. Combination chemotherapy involving agents such as doxorubicin, cyclophosphamide, cisplatin, hexamethylmelamine, paclitaxel, and methotrexate may improve survival rates in these groups, relative to single-agent therapies. Various recently-developed chemotherapeutic agents and treatment regimens have also demonstrated usefulness for treatment of advanced ovarian cancer. For  
30 example, use of the topoisomerase I inhibitor topotecan, use of amifostine to minimize chemotherapeutic side effects, and use of intraperitoneal chemotherapy for patients having peritoneally implanted tumors have demonstrated at least limited utility.

Presently, however, the 5-year survival rate for patients afflicted with stage III ovarian cancer is 25%, and the survival rate for patients afflicted with stage IV ovarian cancer is 8%.

In summary, the earlier ovarian cancer is detected, the aggressiveness of  
5 therapeutic intervention and the side effects associated with therapeutic intervention are minimized. More importantly, the earlier the cancer is detected, the survival rate and quality of life of ovarian cancer patients is enhanced. Thus, a pressing need exists for methods of detecting ovarian cancer as early as possible. There also exists a need for  
10 methods of detecting recurrence of ovarian cancer as well as methods for predicting and monitoring the efficacy of treatment. The present invention satisfies these needs.

### **SUMMARY OF THE INVENTION**

The invention relates to a method of assessing whether a patient is afflicted with ovarian cancer. This method comprises the step of comparing the level of expression of  
15 a marker in a patient sample, wherein the marker is listed in Tables 1-11, and the normal level of expression of the marker in a control, *e.g.*, a sample from a patient without ovarian cancer. A significant difference between the level of expression of the marker in the patient sample and the normal level is an indication that the patient is afflicted with ovarian cancer. In a preferred embodiment, the marker is listed in Tables 2B or 2C  
20 (which are subsets of the markers listed in Table 2A), in Tables 3B or 3C (which are subsets of the markers listed in Table 3A), in Tables 4A or 5A (which are subsets of the markers listed in Tables 4 and 5, respectively), in Table 6A, in Tables 7A-7E or in Table 8. Preferably, a protein corresponding to the marker is a secreted protein or is predicted to correspond to a secreted protein (see, *e.g.* Tables 2D, 4A, 7A-7E). Alternatively, the  
25 marker can correspond to a protein which is normally expressed in ovarian tissue at a detectable level, to one having an extracellular portion, or both (see *e.g.*, Table 8).

In one method, the marker(s) are preferably selected such that the positive predictive value of the method is at least about 10%. Also preferred are embodiments of the method wherein the marker is over- or under-expressed by at least two-fold in at  
30 least about 20% of stage I ovarian cancer patients, stage II ovarian cancer patients, stage III ovarian cancer patients, stage IV ovarian cancer patients, grade I ovarian cancer patients, grade II ovarian cancer patients, grade III ovarian cancer patients, epithelial

ovarian cancer patients, stromal ovarian cancer patients, germ cell ovarian cancer patients, malignant ovarian cancer patients, benign ovarian patients, serous neoplasm ovarian cancer patients, mucinous neoplasm ovarian cancer patients, endometrioid neoplasm ovarian cancer patients and/or clear cell neoplasm ovarian cancer patients.

5           In one embodiment of the methods of the present invention, the patient sample is an ovary-associated body fluid. Such fluids include, for example, blood fluids, lymph, ascitic fluids, gynecological fluids, cystic fluids, urine, and fluids collected by peritoneal rinsing. In another embodiment, the sample comprises cells obtained from the patient. In this embodiment, the cells may be found in a fluid selected from the group consisting  
10 of a fluid collected by peritoneal rinsing, a fluid collected by uterine rinsing, a uterine fluid, a uterine exudate, a pleural fluid, and an ovarian exudate. In another embodiment, the patient sample is *in vivo*.

          In accordance with the methods of the present invention, the level of expression of the marker in a sample can be assessed, for example, by detecting the presence in the  
15 sample of :

- a protein corresponding to the marker or a fragment of the protein (*e.g.* using a reagent, such as an antibody, an antibody derivative, or an antibody fragment, which binds specifically with the protein)
- a metabolite which is produced directly (*i.e.*, catalyzed) or indirectly by a  
20 protein corresponding to the marker
- a transcribed polynucleotide (*e.g.* an mRNA or a cDNA), or fragment thereof, having at least a portion with which the marker is substantially homologous (*e.g.* by contacting a mixture of transcribed polynucleotides obtained from the sample with a substrate having one or more of the markers  
25 listed in Tables 1-11 fixed thereto at selected positions)
- a transcribed polynucleotide or fragment thereof, wherein the polynucleotide anneals with the marker under stringent hybridization conditions.

          The methods of the present invention are particularly useful for patients with an  
30 identified pelvic mass or symptoms associated with ovarian cancer. The methods of the present invention can also be of particular use with patients having an enhanced risk of developing ovarian cancer (*e.g.*, patients having a familial history of ovarian cancer,

- 7 -

patients identified as having a mutant oncogene, and patients at least about 50 years of age). The methods of the present invention may further be of particular use in monitoring the efficacy of treatment of an ovarian cancer patient (e.g. the efficacy of chemotherapy).

5           The methods of the present invention may be performed using a plurality (e.g. 2, 3, 5, or 10 or more) of markers. According to a method involving a plurality of markers, the level of expression in the sample of each of a plurality of markers independently selected from the markers listed in Tables 1-11 is compared with the normal level of expression of each of the plurality of markers in samples of the same type obtained from  
10 control humans not afflicted with ovarian cancer. A significantly enhanced level of expression of one or more of the markers listed in Tables 1, 1A, 2A, 4 and 6, 6A, 7A, 7B, 7D and 8, a significantly reduced level of expression of one or more of the markers listed in Tables 3A, 5, 7C and 7E, or some combination thereof, in the sample, relative to the corresponding normal levels, is an indication that the patient is afflicted with  
15 ovarian cancer. The markers of Tables 1-11 may also be used in combination with known ovarian cancer markers in the methods of the present invention.

In a preferred method of assessing whether a patient is afflicted with ovarian cancer (e.g., new detection ("screening"), detection of recurrence, reflex testing), the method comprises comparing:

- 20           a) the level of expression of a marker in a patient sample, wherein at least one marker is selected from the markers of Tables 1-11 and,  
            b) the normal level of expression of the marker in a control non-ovarian cancer sample.

A significant difference between the level of expression of the marker in the patient  
25 sample and the normal level is an indication that the patient is afflicted with ovarian cancer.

The methods of the present invention further include a method of assessing the efficacy of a test compound for inhibiting ovarian cancer in a patient. This method comprises comparing:

- 8 -

a) expression of a marker in a first sample obtained from the patient and maintained in the presence of the test compound, wherein the marker is selected from the group consisting of the markers listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8, and

5           b) expression of the marker in a second sample obtained from the patient and maintained in the absence of the test compound.

A significantly lower level of expression of the marker in the first sample, relative to the second sample, is an indication that the test compound is efficacious for inhibiting ovarian cancer in the patient. For example, the first and second samples can be portions  
10 of a single sample obtained from the patient or portions of pooled samples obtained from the patient.

The invention still further includes a method of assessing the efficacy of a test compound for inhibiting ovarian cancer in a patient. This method comprises comparing:

15           a) expression of a marker in a first sample obtained from the patient and maintained in the presence of the test compound, wherein the marker is selected from the group consisting of the markers listed in Tables 3A, 5, 7C and 7E, and

b) expression of the marker in a second sample obtained from the patient and maintained in the absence of the test compound.

20 A significantly enhanced level of expression of the marker in the first sample, relative to the second sample, is an indication that the test compound is efficacious for inhibiting the ovarian cancer in the patient.

The invention further relates to a method of assessing the efficacy of a therapy for inhibiting ovarian cancer in a patient. This method comprises comparing:

25           a) expression of a marker in a first sample obtained from the patient prior to providing at least a portion of the therapy to the patient, wherein the marker is selected from the group consisting of the markers listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8, and

30           b) expression of the marker in a second sample obtained from the patient following provision of the portion of the therapy.

A significantly lower level of expression of the marker in the second sample, relative to the first sample, is an indication that the therapy is efficacious for inhibiting ovarian cancer in the patient.

The invention further includes a method of assessing the efficacy of a therapy for inhibiting ovarian cancer in a patient, comprising comparing:

- a) expression of a marker in a first sample obtained from the patient prior to providing at least a portion of the therapy to the patient, wherein the marker is selected from the group consisting of the markers listed in Tables 3A, 5, 7C and 7E, and
- b) expression of the marker in a second sample obtained from the patient following provision of the portion of the therapy.

A significantly enhanced level of expression of the marker in the second sample, relative to the first sample, is an indication that the therapy is efficacious for inhibiting ovarian cancer in the patient.

It will be appreciated that in these methods the "therapy" may be any traditional therapy for treating ovarian cancer including, but not limited to, chemotherapy, radiation therapy and surgical removal of tissue, e.g., an ovarian tumor. Thus, the methods of the invention may be used to evaluate a patient before, during and after therapy, for example, to evaluate the reduction in tumor burden.

The present invention therefore further comprises a method for monitoring the progression of ovarian cancer in a patient, the method comprising:

- a) detecting in a patient sample at a first time point, the expression of a marker, wherein the marker is selected from the group consisting of the markers listed in Tables 1-11;
- b) repeating step a) at a subsequent time point in time; and
- c) comparing the level of expression detected in steps a) and b), and therefrom monitoring the progression of ovarian cancer in the patient.

The invention also includes a method of selecting a composition for inhibiting ovarian cancer in a patient. This method comprises the steps of:

- a) obtaining a sample comprising cancer cells from the patient;
- b) separately maintaining aliquots of the sample in the presence of a plurality of test compositions;

- 10 -

c) comparing expression of a marker listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8 in each of the aliquots; and

d) selecting one of the test compositions which induces a lower level of expression of the marker in the aliquot containing that test composition, relative to other test compositions.

5

The invention further includes a method of selecting a composition for inhibiting ovarian cancer in a patient. This method comprises the steps of:

a) obtaining a sample comprising cancer cells from the patient;

b) separately maintaining aliquots of the sample in the presence of a plurality of test compositions;

10

c) comparing expression of a marker listed in Tables 3A, 5, 7C and 7E in each of the aliquots; and

d) selecting one of the test compositions which induces an enhanced level of expression of the marker in the aliquot containing that test

15

composition, relative to other test compositions.

In addition, the invention includes a method of inhibiting ovarian cancer in a patient. This method comprises the steps of:

a) obtaining a sample comprising cancer cells from the patient;

b) separately maintaining aliquots of the sample in the presence of a plurality of test compositions;

20

c) comparing expression of a marker listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8 in each of the aliquots; and

d) administering to the patient at least one of the test compositions which induces a lower level of expression of the marker in the aliquot

25

containing that test composition, relative to other test compositions.

The invention also includes a method of inhibiting ovarian cancer in a patient. This method comprises the steps of:

a) obtaining a sample comprising cancer cells from the patient;

b) separately maintaining aliquots of the sample in the presence of a plurality of test compositions;

30

c) comparing expression of a marker listed in Tables 3A, 5, 7C and 7E, in each of the aliquots; and

- 11 -

d) administering to the patient at least one of the test compositions which induces an enhanced expression of the marker in the aliquot containing that test composition, relative to other test compositions.

The invention also includes a kit for assessing whether a patient is afflicted with  
5 ovarian cancer. This kit comprises reagents for assessing expression of a marker listed in Tables 1-11.

In another aspect, the invention relates to a kit for assessing the suitability of each of a plurality of compounds for inhibiting an ovarian cancer in a patient. The kit comprises a reagent for assessing expression of a marker listed in Tables 1-11, and may  
10 also comprise a plurality of compounds.

In another aspect, the invention relates to a kit for assessing the presence of ovarian cancer cells. This kit comprises an antibody, wherein the antibody binds specifically with a protein corresponding to a marker listed in Tables 1-11. The kit may also comprise a plurality of antibodies, wherein the plurality binds specifically with a  
15 protein corresponding to a different marker listed in Tables 1-11.

The invention also includes a kit for assessing the presence of ovarian cancer cells, wherein the kit comprises a nucleic acid probe. The probe binds specifically with a transcribed polynucleotide corresponding to a marker listed in Tables 1-11. The kit may also comprise a plurality of probes, wherein each of the probes binds specifically  
20 with a transcribed polynucleotide corresponding to a different marker listed in Tables 1-11.

The invention further relates to a method of making an isolated hybridoma which produces an antibody useful for assessing whether a patient is afflicted with ovarian cancer. The method comprises isolating a protein corresponding to a marker listed in  
25 Tables 1-11, immunizing a mammal using the isolated protein, isolating splenocytes from the immunized mammal, fusing the isolated splenocytes with an immortalized cell line to form hybridomas, and screening individual hybridomas for production of an antibody which specifically binds with the protein to isolate the hybridoma. The invention also includes an antibody produced by this method.



- 12 -

The invention further includes a method of assessing the ovarian carcinogenic potential of a test compound. This method comprises the steps of:

- a) maintaining separate aliquots of ovarian cells in the presence and absence of the test compound; and
- 5        b) comparing expression of a marker in each of the aliquots.

The marker is selected from those listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8. A significantly enhanced level of expression of the marker in the aliquot maintained in the presence of (or exposed to) the test compound, relative to the aliquot maintained in the absence of the test compound, is an indication that the test compound  
10       possesses ovarian carcinogenic potential.

The invention includes another method of assessing the ovarian carcinogenic potential of a test compound. This method comprises the steps of:

- a) maintaining separate aliquots of ovarian cells in the presence and absence of the test compound; and
- 15       b) comparing expression of a marker in each of the aliquots.

In this method, the marker is selected from those listed in Tables 3A, 5, 7C and 7E. A significantly lower level of expression of the marker in the aliquot maintained in the presence of the test compound, relative to the aliquot maintained in the absence of the test compound, is an indication that the test compound possesses ovarian  
20       carcinogenic potential.

Additionally, the invention includes a kit for assessing the ovarian carcinogenic potential of a test compound. The kit comprises ovarian cells and a reagent for assessing expression of a marker in each of the aliquots. The marker is selected from those listed in Tables 1-11.

25       The invention further relates to a method of treating a patient afflicted with ovarian cancer. This method comprises providing to cells of the patient a protein corresponding to a marker listed in Tables 3A, 5, 7C and 7E. The protein can be provided to the cells, for example, by providing a vector comprising a polynucleotide encoding the protein to the cells.

The invention includes another method of treating a patient afflicted with ovarian cancer. This method comprises providing to cells of the patient an antisense oligonucleotide complementary to a polynucleotide corresponding to a marker listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8.

5           The invention includes a method of inhibiting ovarian cancer in a patient at risk for developing ovarian cancer. This method comprises inhibiting expression or overexpression of a gene corresponding to a marker listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8.

10           The invention includes another method of inhibiting ovarian cancer in a patient at risk for developing ovarian cancer. This method comprises enhancing expression of a gene corresponding to a marker listed in Tables 3A, 5, 7C and 7E.

It will be appreciated that the methods and kits of the present invention may also include known cancer markers including known ovarian cancer markers. It will further be appreciated that the methods and kits may be used to identify cancers other than  
15   ovarian cancer.

### **DETAILED DESCRIPTION OF THE INVENTION**

The invention relates to newly discovered correlations between expression of certain markers and the cancerous state of ovarian cells. It has been discovered that the  
20   level of expression of individual markers and combinations of markers described herein correlates with the presence of ovarian cancer in a patient. Methods are provided for detecting the presence of ovarian cancer in a sample, the absence of ovarian cancer in a sample, the stage of an ovarian cancer, and with other characteristics of ovarian cancer that are relevant to prevention, diagnosis, characterization, and therapy of ovarian cancer  
25   in a patient.

#### **Definitions**

As used herein, each of the following terms has the meaning associated with it in this section.

30           The articles "a" and "an" are used herein to refer to one or to more than one (*i.e.* to at least one) of the grammatical object of the article. By way of example, "an element" means one element or more than one element.

A "marker" of the invention is a naturally-occurring polymer corresponding to at least one of the nucleic acids listed in Tables 1-11. In particular, a marker of the invention may be a nucleic acid molecule comprising a sequence listed in Tables 1-11 or a sequence which hybridizes under high stringency conditions with a polynucleotide  
5 sequence listed in Tables 1-11 ("nucleic acid marker"). Nucleic acid markers include, without limitation, sense and anti-sense strands of genomic DNA (*i.e.* including any introns occurring therein), RNA generated by transcription of genomic DNA (*i.e.* prior to splicing), RNA generated by splicing of RNA transcribed from genomic DNA, and proteins generated by translation of spliced RNA (*i.e.* including proteins both before and  
10 after cleavage of normally cleaved regions such as transmembrane signal sequences). As used herein, "marker" may also include a cDNA made by reverse transcription of an RNA generated by transcription of genomic DNA (including spliced RNA). A marker of the invention also may be a protein encoded by, for example, a nucleic acid marker.

The term "probe" refers to any molecule which is capable of selectively binding  
15 to a specifically intended target molecule, for example a marker of the invention. Probes can be either synthesized by one skilled in the art, or derived from appropriate biological preparations. For purposes of detection of the target molecule, probes may be specifically designed to be labeled, as described herein. Examples of molecules that can be utilized as probes include, but are not limited to, RNA, DNA, proteins, antibodies,  
20 and organic monomers.

An "ovary-associated" body fluid is a fluid which, when in the body of a patient, contacts or passes through ovarian cells or into which cells or proteins shed from ovarian cells *e.g.* ovarian epithelium, are capable of passing. Exemplary ovary-associated body fluids include blood fluids, lymph, ascites, gynecological fluids, cystic  
25 fluid, urine, and fluids collected by peritoneal rinsing.

The "normal" level of expression of a marker is the level of expression of the marker in ovarian cells of a patient, *e.g.* a human, not afflicted with ovarian cancer.

"Over-expression" and "under-expression" of a marker refer to expression of the marker of a patient at a greater or lesser level, respectively, than normal level of  
30 expression of the marker (*e.g.* at least two-fold greater or lesser level).

As used herein, the term "promoter/regulatory sequence" means a nucleic acid sequence which is required for expression of a gene product operably linked to the promoter/regulatory sequence. In some instances, this sequence may be the core promoter sequence and in other instances, this sequence may also include an enhancer sequence and other regulatory elements which are required for expression of the gene product. The promoter/regulatory sequence may, for example, be one which expresses the gene product in a tissue-specific manner.

A "constitutive" promoter is a nucleotide sequence which, when operably linked with a polynucleotide which encodes or specifies a gene product, causes the gene product to be produced in a living human cell under most or all physiological conditions of the cell.

An "inducible" promoter is a nucleotide sequence which, when operably linked with a polynucleotide which encodes or specifies a gene product, causes the gene product to be produced in a living human cell substantially only when an inducer which corresponds to the promoter is present in the cell.

A "tissue-specific" promoter is a nucleotide sequence which, when operably linked with a polynucleotide which encodes or specifies a gene product, causes the gene product to be produced in a living human cell substantially only if the cell is a cell of the tissue type corresponding to the promoter.

A "transcribed polynucleotide" is a polynucleotide (*e.g.* an RNA, a cDNA, or an analog of one of an RNA or cDNA) which is complementary to or homologous with all or a portion of a mature RNA made by transcription of a genomic DNA corresponding to a marker of the invention and normal post-transcriptional processing (*e.g.* splicing), if any, of the transcript.

"Complementary" refers to the broad concept of sequence complementarity between regions of two nucleic acid strands or between two regions of the same nucleic acid strand. It is known that an adenine residue of a first nucleic acid region is capable of forming specific hydrogen bonds ("base pairing") with a residue of a second nucleic acid region which is antiparallel to the first region if the residue is thymine or uracil. Similarly, it is known that a cytosine residue of a first nucleic acid strand is capable of base pairing with a residue of a second nucleic acid strand which is antiparallel to the first strand if the residue is guanine. A first region of a nucleic acid is complementary to

a second region of the same or a different nucleic acid if, when the two regions are arranged in an antiparallel fashion, at least one nucleotide residue of the first region is capable of base pairing with a residue of the second region. Preferably, the first region comprises a first portion and the second region comprises a second portion, whereby,  
5 when the first and second portions are arranged in an antiparallel fashion, at least about 50%, and preferably at least about 75%, at least about 90%, or at least about 95% of the nucleotide residues of the first portion are capable of base pairing with nucleotide residues in the second portion. More preferably, all nucleotide residues of the first portion are capable of base pairing with nucleotide residues in the second portion.

10 "Homologous" as used herein, refers to nucleotide sequence similarity between two regions of the same nucleic acid strand or between regions of two different nucleic acid strands. When a nucleotide residue position in both regions is occupied by the same nucleotide residue, then the regions are homologous at that position. A first region is homologous to a second region if at least one nucleotide residue position of each  
15 region is occupied by the same residue. Homology between two regions is expressed in terms of the proportion of nucleotide residue positions of the two regions that are occupied by the same nucleotide residue. By way of example, a region having the nucleotide sequence 5'-ATTGCC-3' and a region having the nucleotide sequence 5'-TATGGC-3' share 50% homology. Preferably, the first region comprises a first portion and the second region comprises a second portion, whereby, at least about 50%, and  
20 preferably at least about 75%, at least about 90%, or at least about 95% of the nucleotide residue positions of each of the portions are occupied by the same nucleotide residue. More preferably, all nucleotide residue positions of each of the portions are occupied by the same nucleotide residue.

25 A marker is "fixed" to a substrate if it is covalently or non-covalently associated with the substrate such the substrate can be rinsed with a fluid (*e.g.* standard saline citrate, pH 7.4) without a substantial fraction of the marker dissociating from the substrate.

As used herein, a "naturally-occurring" nucleic acid molecule refers to an RNA  
30 or DNA molecule having a nucleotide sequence that occurs in nature (*e.g.* encodes a natural protein).

Expression of a marker in a patient is "significantly" higher or lower than the normal level of expression of a marker if the level of expression of the marker is greater or less, respectively, than the normal level by an amount greater than the standard error of the assay employed to assess expression, and preferably at least twice, and more preferably three, four, five or ten times that amount. Alternately, expression of the marker in the patient can be considered "significantly" higher or lower than the normal level of expression if the level of expression is at least about two, and preferably at least about three, four, or five times, higher or lower, respectively, than the normal level of expression of the marker.

10 Ovarian cancer is "inhibited" if at least one symptom of the cancer is alleviated, terminated, slowed, or prevented. As used herein, ovarian cancer is also "inhibited" if recurrence or metastasis of the cancer is reduced, slowed, delayed, or prevented.

A kit is any manufacture (*e.g.* a package or container) comprising at least one reagent, *e.g.* a probe, for specifically detecting a marker of the invention, the manufacture being promoted, distributed, or sold as a unit for performing the methods of the present invention.

### Description

The present invention is based, in part, on identification of markers which are expressed at a different level in ovarian cancer cells than they are in normal (*i.e.* non-cancerous) ovarian cells. The markers of the invention correspond to nucleic acid and polypeptide molecules which can be detected in one or both of normal and cancerous ovarian cells. The presence, absence, or level of expression of one or more of these markers in ovarian cells is herein correlated with the cancerous state of the tissue. The invention thus includes compositions, kits, and methods for assessing the cancerous state of ovarian cells (*e.g.* cells obtained from a human, cultured human cells, archived or preserved human cells and *in vivo* cells).

The compositions, kits, and methods of the invention have the following uses, among others:

- 30 1) assessing whether a patient is afflicted with ovarian cancer;
- 2) assessing the stage of ovarian cancer in a human patient;
- 3) assessing the grade of ovarian cancer in a patient;

- 18 -

- 4) assessing the benign or malignant nature of ovarian cancer in a patient;
- 5) assessing the histological type of neoplasm (*e.g.* serous, mucinous, endometrioid, or clear cell neoplasm) associated with ovarian cancer in a patient;
- 6) making an isolated hybridoma which produces an antibody useful for assessing whether a patient is afflicted with ovarian cancer;
- 7) assessing the presence of ovarian cancer cells;
- 8) assessing the efficacy of one or more test compounds for inhibiting ovarian cancer in a patient;
- 9) assessing the efficacy of a therapy for inhibiting ovarian cancer in a patient;
- 10) monitoring the progression of ovarian cancer in a patient;
- 11) selecting a composition or therapy for inhibiting ovarian cancer in a patient;
- 12) treating a patient afflicted with ovarian cancer;
- 13) inhibiting ovarian cancer in a patient;
- 14) assessing the ovarian carcinogenic potential of a test compound; and
- 15) inhibiting an ovarian cancer in a patient at risk for developing ovarian cancer.

The invention thus includes a method of assessing whether a patient is afflicted with ovarian cancer. This method comprises comparing the level of expression of a marker in a patient sample and the normal level of expression of the marker in a control, *e.g.*, a non-ovarian cancer sample. A significant difference between the level of expression of the marker in the patient sample and the normal level is an indication that the patient is afflicted with ovarian cancer. The marker is selected from the group consisting of the markers listed in Tables 1-11. The markers listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8 are expressed at a greater level in ovarian cancer cells than in normal ovarian cells. The markers listed in Tables 3A, 5, 7C and 7E are expressed at a lower level in ovarian cancer cells than in normal ovarian cells. Although one or more

molecules corresponding to the markers listed in Tables 1-11 may have been described by others, the significance of the level of expression of these markers with regard to the cancerous state of ovarian cells has not previously been recognized.

Tables 1 and 1A list markers that were identified in subtractive libraries and  
5 which are preferentially expressed in ovarian cancer cells over normal (*i.e.*, non-cancerous) ovarian cells.

Table 2A lists markers, expression of which was increased by at least 5-fold in at least one of twenty-three ovarian cancer samples tested, relative to its expression in normal (*i.e.* non-cancerous) ovarian samples. Table 2B lists markers, expression of  
10 which was increased by at least 2-fold in all twenty-three ovarian cancer samples tested, relative to its expression in normal ovarian samples. Table 2C lists markers, expression of which was increased by at least 5-fold in at least 6 of the 23 ovarian cancer samples tested, relative to its expression in normal ovarian cells. Table 2D lists markers,  
15 expression of which was increased by at least 5-fold in at least 6 of the 23 ovarian cancer samples, relative to expression in normal ovarian samples. In a preferred embodiment, proteins corresponding to the markers of Table 2D as well as fragments of the proteins, serve as antigens for antibody production, based upon proteomic studies, sequence analysis and/or literature references

Table 3A lists markers, expression of which was decreased by at least 5-fold in  
20 at least one of twenty-three ovarian cancer samples tested, relative to its expression in normal (*i.e.*, non-cancerous) ovarian cells. Table 3B lists markers, expression of which was decreased by at least 2-fold in all twenty-three ovarian cancer samples tested, relative to its expression in normal ovarian cells. Table 3C lists markers, expression of which was decreased by at least 5-fold in at least 6 of the 23 ovarian cancer samples  
25 tested, relative to its expression in normal ovarian cells.

Tables 4 and 5 list markers, expression of which was either increased (Table 4) or decreased (Table 5) in ovarian cancer samples, relative to expression in normal (*i.e.*, non-cancerous) ovarian samples. In particular, expression of the markers in 37 tumors (7 endometrioid tumors, 5 clear cell tumors and 25 serous tumors) was evaluated. A  
30 ranking system based on the sum of the number of tumors multiplied by the fold regulation (for 2-fold, 3-fold, 5-fold and 10-fold regulation), divided by the total number



of tumors, was employed. A rank score was generated for four categories, endometroid tumors, clear cell tumors, serous tumors and overall.

For example, for # 19109 in Table 4A (first marker listed):

# of tumors > 2-fold:  $36 = (2 \times 0) = 0$

5 # of tumors > 3-fold:  $36 = (3 \times 1) = 3$

# of tumors > 5-fold:  $35 = (5 \times 3) = 15$

# of tumors > 10-fold:  $32 = (10 \times 32) = 320$

The sum is 3 plus 15 plus 320, which equals 338. The score is therefore 338 divided by 37, which equals 9.1.

10 The markers of Table 4 had a score of greater than 1.5 for endometroid tumors, greater than 1.5 for clear cell tumors, greater than 1 for serous tumors, or greater than 0.8 overall. Table 4A shows the markers of Table 4 with a score of greater than 3 in any of the four categories.

The markers of Table 5 had a score of greater than 2.5 for endometroid tumors, 15 greater than 2.5 for clear cell tumors, greater than 2 for serous tumors, or greater than 1.75 overall. Table 5A shows the markers of Table 5 with a score of greater than 3 in any of the four categories.

Tables 6 and 6A list markers that were identified in subtractive libraries and which are preferentially expressed in ovarian cancer cells over normal (*i.e.*, non- 20 cancerous) ovarian cells.

Tables 7A-7E list markers that were identified in proteomic studies. The markers of Table 7A are secreted or membrane proteins, expression of which was increased at least 5-fold in two or more ovarian cancer samples or cell lines, relative to at least a 10-fold decrease in expression in normal ovarian samples.

25 The markers of Table 7B are secreted or membrane proteins, expression of which was increased in one ovarian cancer sample cell line, relative to expression in normal ovarian samples, where the medium expression of normals equaled 0 (the expression level of the ovarian cancer sample and cell lines was divided by 0.001, rather than 0).

30 The markers of Table 7C are preferred secreted or membrane proteins, expression of which was decreased in ovarian cancer samples and cell lines, relative to expression in normal ovarian samples.

The markers of Table 7D are secreted or membrane proteins present in ovarian cancer cell supernatants.

The markers of Table 7E are secreted or membrane proteins present in normal cell supernatants.

5 Table 8 lists novel genes that are overexpressed in ovarian cancer samples, relative to expression in normal ovarian samples.

Table 9 summarizes TaqMan® expression data for the novel genes of Table 8.

Tables 10A-10N summarize Northern Blot analysis of the novel genes of Table 8.

10 Table 11 summarizes LightCycler data and RT-PCR data for various markers of the present invention.

Any marker or combination of markers listed in Tables 1-11, as well as any known markers in combination with the markers set forth in Tables 1-11, may be used in the compositions, kits, and methods of the present invention. Use of markers listed in  
15 Tables 2B, 2C, 2D, 3B, 3C, 4A, 5A, 6A, 7A-7E and 8 are preferred, wherein use of markers listed in Tables 2C, 2D, 3C, 6A, 7A-7C and 8 are more preferred. In general, it is preferable to use markers for which the difference between the level of expression of the marker in ovarian cancer cells and the level of expression of the same marker in normal ovarian cells is as great as possible. Although this difference can be as small as  
20 the limit of detection of the method for assessing expression of the marker, it is preferred that the difference be at least greater than the standard error of the assessment method, and preferably a difference of at least 2-, 3-, 4-, 5-, 6-, 7-, 8-, 9-, 10-, 15-, 20-, 25-, 100-, 500-, 1000-fold or greater.

It is recognized that certain markers correspond to proteins which are secreted  
25 from ovarian cells (*i.e.* one or both of normal and cancerous cells) to the extracellular space surrounding the cells (see, *e.g.* Tables 2D, 7A-7E and 8). These markers are preferably used in certain embodiments of the compositions, kits, and methods of the invention, owing to the fact that the protein corresponding to each of these markers can be detected in an ovary-associated body fluid sample, which may be more easily  
30 collected from a human patient than a tissue biopsy sample. In addition, preferred *in vivo* techniques for detection of a protein corresponding to a marker of the invention include introducing into a subject a labeled antibody directed against the protein. For

example, the antibody can be labeled with a radioactive marker whose presence and location in a subject can be detected by standard imaging techniques.

Although not every marker corresponding to a secreted protein is indicated as such in the Tables herein, it is a simple matter for the skilled artisan to determine

5 whether any particular marker corresponds to a secreted protein. In order to make this determination, the protein corresponding to a marker is expressed in a test cell (*e.g.* a cell of an ovarian cell line), extracellular fluid is collected, and the presence or absence of the protein in the extracellular fluid is assessed (*e.g.* using a labeled antibody which binds specifically with the protein).

10 The following is an example of a method which can be used to detect secretion of a protein corresponding to a marker of the invention. About  $8 \times 10^5$  293T cells are incubated at 37°C in wells containing growth medium (Dulbecco's modified Eagle's medium {DMEM} supplemented with 10% fetal bovine serum) under a 5% (v/v) CO<sub>2</sub>, 95% air atmosphere to about 60-70% confluence. The cells are then transfected using a  
15 standard transfection mixture comprising 2 micrograms of DNA comprising an expression vector encoding the protein and 10 microliters of LipofectAMINE™ (GIBCO/BRL Catalog no. 18342-012) per well. The transfection mixture is maintained for about 5 hours, and then replaced with fresh growth medium and maintained in an air atmosphere. Each well is gently rinsed twice with DMEM which does not contain  
20 methionine or cysteine (DMEM-MC; ICN Catalog no. 16-424- 54). About 1 milliliter of DMEM-MC and about 50 microcuries of Trans-<sup>35</sup>S™ reagent (ICN Catalog no. 51006) are added to each well. The wells are maintained under the 5% CO<sub>2</sub> atmosphere described above and incubated at 37°C for a selected period. Following incubation, 150 microliters of conditioned medium is removed and centrifuged to remove floating cells  
25 and debris. The presence of the protein in the supernatant is an indication that the protein is secreted.

Examples of ovary-associated body fluids include blood fluids (*e.g.* whole blood, blood serum, blood having platelets removed therefrom, etc.), lymph, ascitic fluids, gynecological fluids (*e.g.* ovarian, fallopian, and uterine secretions, menses, vaginal  
30 douching fluids, fluids used to rinse cervical cell samples, etc.), cystic fluid, urine, and fluids collected by peritoneal rinsing (*e.g.* fluids applied and collected during laparoscopy or fluids instilled into and withdrawn from the peritoneal cavity of a human

patient). In these embodiments, the level of expression of the marker can be assessed by assessing the amount (*e.g.* absolute amount or concentration) of the marker in an ovary-associated body fluid obtained from a patient. The fluid can, of course, be subjected to a variety of well-known post-collection preparative and storage techniques (*e.g.* storage, freezing, ultrafiltration, concentration, evaporation, centrifugation, etc.) prior to assessing the amount of the marker in the fluid.

Many ovary-associated body fluids (*i.e.* usually excluding urine) can have ovarian cells, *e.g.* ovarian epithelium, therein, particularly when the ovarian cells are cancerous, and, more particularly, when the ovarian cancer is metastasizing. Cell-containing fluids which can contain ovarian cancer cells include, but are not limited to, peritoneal ascites, fluids collected by peritoneal rinsing, fluids collected by uterine rinsing, uterine fluids such as uterine exudate and menses, pleural fluid, and ovarian exudates. Thus, the compositions, kits, and methods of the invention can be used to detect expression of markers corresponding to proteins having at least one portion which is displayed on the surface of cells which express it. Examples of such proteins are indicated in the Tables herein. Although not every protein having at least one cell-surface portion is indicated in the Tables, it is a simple matter for the skilled artisan to determine whether the protein corresponding to any particular marker comprises a cell-surface protein. For example, immunological methods may be used to detect such proteins on whole cells, or well known computer-based sequence analysis methods (*e.g.* the SIGNALP program; Nielsen *et al.*, 1997, *Protein Engineering* 10:1-6) may be used to predict the presence of at least one extracellular domain (*i.e.* including both secreted proteins and proteins having at least one cell-surface domain). Expression of a marker corresponding to a protein having at least one portion which is displayed on the surface of a cell which expresses it may be detected without necessarily lysing the cell (*e.g.* using a labeled antibody which binds specifically with a cell-surface domain of the protein).

Expression of a marker of the invention may be assessed by any of a wide variety of well known methods for detecting expression of a transcribed molecule or its corresponding protein. Non-limiting examples of such methods include immunological methods for detection of secreted, cell-surface, cytoplasmic, or nuclear proteins, protein purification methods, protein function or activity assays, nucleic acid hybridization

methods, nucleic acid reverse transcription methods, and nucleic acid amplification methods.

In a preferred embodiment, expression of a marker is assessed using an antibody (e.g. a radio-labeled, chromophore-labeled, fluorophore-labeled, or enzyme-labeled antibody), an antibody derivative (e.g. an antibody conjugated with a substrate or with the protein or ligand of a protein-ligand pair {e.g. biotin-streptavidin} ), or an antibody fragment (e.g. a single-chain antibody, an isolated antibody hypervariable domain, etc.) which binds specifically with a protein corresponding to the marker, such as the protein encoded by the open reading frame corresponding to the marker or such a protein which has undergone all or a portion of its normal post-translational modification.

In another preferred embodiment, expression of a marker is assessed by preparing mRNA/cDNA (i.e. a transcribed polynucleotide) from cells in a patient sample, and by hybridizing the mRNA/cDNA with a reference polynucleotide which is a complement of a polynucleotide comprising the marker, and fragments thereof. cDNA can, optionally, be amplified using any of a variety of polymerase chain reaction methods prior to hybridization with the reference polynucleotide; preferably, it is not amplified. Expression of one or more markers can likewise be detected using quantitative PCR to assess the level of expression of the marker(s). Alternatively, any of the many known methods of detecting mutations or variants (e.g. single nucleotide polymorphisms, deletions, etc.) of a marker of the invention may be used to detect occurrence of a marker in a patient.

In a related embodiment, a mixture of transcribed polynucleotides obtained from the sample is contacted with a substrate having fixed thereto a polynucleotide complementary to or homologous with at least a portion (e.g. at least 7, 10, 15, 20, 25, 30, 40, 50, 100, 500, or more nucleotide residues) of a marker of the invention. If polynucleotides complementary to or homologous with are differentially detectable on the substrate (e.g. detectable using different chromophores or fluorophores, or fixed to different selected positions), then the levels of expression of a plurality of markers can be assessed simultaneously using a single substrate (e.g. a "gene chip" microarray of polynucleotides fixed at selected positions). When a method of assessing marker expression is used which involves hybridization of one nucleic acid with another, it is preferred that the hybridization be performed under stringent hybridization conditions.

Because the compositions, kits, and methods of the invention rely on detection of a difference in expression levels of one or more markers of the invention, it is preferable that the level of expression of the marker is significantly greater than the minimum detection limit of the method used to assess expression in at least one of normal ovarian cells and cancerous ovarian cells.

Preferably, at least one of the marker(s) used in the compositions, kits, and methods of the invention is a marker for which the "Tissue Prominence," as indicated in the Tables herein, includes, without limitation, an epithelial tissue such as ovarian, stomach, foreskin, colon, uterus, esophagus, synovial membrane, small intestine, breast, skin, cervix, adrenal gland, eye, gall bladder, lung, placenta, prostate and retina tissues. Preferably, the marker is one for which ovary is listed among the Tissue Prominence tissues in one or more of the Tables.

The chromosomal location corresponding to each of a number of the markers listed in the Tables herein is known and is also listed in the Tables. In addition, the chromosomal locations of a number of loci and chromosomal regions associated with ovarian cancers are known (Lynch *et al.*, 1998, *Sem. Oncol.* 25:265-280). For example, *AKT2* is located on chromosome 19 at q13.1-13.2, copy number increases have been observed at 8q24, 20q13.2-qter, 3q26.3, 1q32, 20p, 9p21-pter, 12p, and 5p14-pter, DNA amplifications have been observed at 8q24, 3q26.3, and 20q13.3, *c-MYC* is located at 8q24, *MYBL2* is located at 20q13.1, *EVII* is located at 3q26, loss of heterozygosity has been observed on chromosomes 6, 9, 13q, 17, 18q, 19p, 22q and Xp, including at locations 17p(p13.3, 13.1), 17q(q21, q22-q23), 18q (q21.3-qter), 6q(q26-q27), 11q(q23.3-qter), and 11p(p13-p15.5), *TP53* is located at 17p13.1, *BRCA1* is located at 17q21, the prohibitin gene and *NM23* are both located at 17q23-24, *NF1* is located at 17q11, and *ERBB2* is located at 17q21. At least one previously unidentified gene which contributes to development of ovarian cancer has been suggested to reside on chromosome 17 (Lynch *et al.*, *supra*), particularly on 17p, and more particularly in the vicinity of 17p13.3. Thus, markers which map to one or more of these chromosomal locations, or to a location relatively near one of these locations are preferred for use in the compositions, kits, and methods of the invention.

- 26 -

It is understood that by routine screening of additional patient samples using one or more of the markers of the invention, it will be realized that certain of the markers are over- or under-expressed in cancers of various types, including specific ovarian cancers, as well as other cancers such as breast cancer, cervical cancer, etc. For example, it will

5 be confirmed that some of the markers of the invention are over- or under-expressed in most (*i.e.* 50% or more) or substantially all (*i.e.* 80% or more) of ovarian cancer. Furthermore, it will be confirmed that certain of the markers of the invention are associated with ovarian cancer of various stages (*i.e.* stage I, II, III, and IV ovarian cancers, as well as subclassifications IA, IB, IC, IIA, IIB, IIC, IIIA, IIIB, and IIIC, using

10 the FIGO Stage Grouping system for primary carcinoma of the ovary; 1987, *Am. J. Obstet. Gynecol.* 156:236), of various histologic subtypes (*e.g.* serous, mucinous, endometrioid, and clear cell subtypes, as well as subclassifications and alternate classifications adenocarcinoma, papillary adenocarcinoma, papillary

15 cystadenocarcinoma, surface papillary carcinoma, malignant adenofibroma, cystadenofibroma, adenocarcinoma, cystadenocarcinoma, adenoacanthoma, endometrioid stromal sarcoma, mesodermal (Müllerian) mixed tumor, mesonephroid tumor, malignant carcinoma, Brenner tumor, mixed epithelial tumor, and undifferentiated carcinoma, using the WHO/FIGO system for classification of malignant ovarian tumors; Scully, *Atlas of Tumor Pathology*, 3d series, Washington DC), and

20 various grades (*i.e.* grade I {well differentiated} , grade II {moderately well differentiated}, and grade III {poorly differentiated from surrounding normal tissue} ). In addition, as a greater number of patient samples are assessed for expression of the markers of the invention and the outcomes of the individual patients from whom the samples were obtained are correlated, it will also be confirmed that altered expression of

25 certain of the markers of the invention are strongly correlated with malignant cancers and that altered expression of other markers of the invention are strongly correlated with benign tumors. The compositions, kits, and methods of the invention are thus useful for characterizing one or more of the stage, grade, histological type, and benign/malignant nature of ovarian cancer in patients. In addition, these compositions, kits, and methods

30 can be used to detect and differentiate epithelial, stromal, and germ cell ovarian cancers.

When the compositions, kits, and methods of the invention are used for characterizing one or more of the stage, grade, histological type, and benign/malignant nature of ovarian cancer in a patient, it is preferred that the marker or panel of markers of the invention is selected such that a positive result is obtained in at least about 20%,  
5 and preferably at least about 40%, 60%, or 80%, and more preferably in substantially all patients afflicted with an ovarian cancer of the corresponding stage, grade, histological type, or benign/malignant nature. Preferably, the marker or panel of markers of the invention is selected such that a PPV of greater than about 10% is obtained for the general population (more preferably coupled with an assay specificity greater than  
10 99.5%).

When a plurality of markers of the invention are used in the compositions, kits, and methods of the invention, the level of expression of each marker in a patient sample can be compared with the normal level of expression of each of the plurality of markers in non-cancerous samples of the same type, either in a single reaction mixture (*i.e.* using  
15 reagents, such as different fluorescent probes, for each marker) or in individual reaction mixtures corresponding to one or more of the markers. In one embodiment, a significantly enhanced level of expression of more than one of the plurality of markers in the sample, relative to the corresponding normal levels, is an indication that the patient is afflicted with ovarian cancer. In another embodiment, a significantly lower  
20 level of expression in the sample of each of the plurality of markers, relative to the corresponding normal levels, is an indication that the patient is afflicted with ovarian cancer. In yet another embodiment, a significantly enhanced level of expression of one or more marks and a significantly lower level of expression of one or more markers in a sample relative to the corresponding normal levels, is an indication that the patient is  
25 afflicted with ovarian cancer. When a plurality of markers is used, it is preferred that 2, 3, 4, 5, 8, 10, 12, 15, 20, 30, or 50 or more individual markers be used, wherein fewer markers are preferred.

In order to maximize the sensitivity of the compositions, kits, and methods of the invention (*i.e.* by interference attributable to cells of non-ovarian origin in a patient  
30 sample), it is preferable that the marker of the invention used therein be a marker which has a restricted tissue distribution, *e.g.*, normally not expressed in a non-epithelial tissue, and more preferably a marker which is normally not expressed in a non-ovarian tissue.



Only a small number of markers are known to be associated with ovarian cancers (*e.g.* *AKT2*, *Ki-RAS*, *ERBB2*, *c-MYC*, *RBI*, and *TP53*; Lynch, *supra*). These markers are not, of course, included among the markers of the invention, although they may be used together with one or more markers of the invention in a panel of markers, for example.

- 5 It is well known that certain types of genes, such as oncogenes, tumor suppressor genes, growth factor-like genes, protease-like genes, and protein kinase-like genes are often involved with development of cancers of various types. Thus, among the markers of the invention, use of those which correspond to proteins which resemble known proteins encoded by known oncogenes and tumor suppressor genes, and those which correspond  
10 to proteins which resemble growth factors, proteases, and protein kinases are preferred.

- Known oncogenes and tumor suppressor genes include, for example, *abl*, *abr*, *akt2*, *apc*, *bcl2 $\alpha$* , *bcl2 $\beta$* , *bcl3*, *bcr*, *brca1*, *brca2*, *cbl*, *ccnd1*, *cdc42*, *cdk4*, *crk- II*, *csflr/fms*, *dbl*, *dcc*, *dpc4/smad4*, *e-cad*, *e2f1/rbap*, *egfr/erbb-1*, *elk1*, *elk3*, *eph*, *erg*, *ets1*, *ets2*, *fer*, *fgr/src2*, *flil/erbb2*, *fos*, *fps/fes*, *fra1*, *fra2*, *fyn*, *hck*, *hek*, *her2/erbb- 2/neu*,  
15 *her3/erbb-3*, *her4/erbb-4*, *hras1*, *hst2*, *hstf1*, *igfbp2*, *ink4a*, *ink4b*, *int2/fgf3*, *jun*, *junb*, *jund*, *kip2*, *kit*, *kras2a*, *kras2b*, *lck*, *lyn*, *mas*, *max*, *mcc*, *mdm2*, *met*, *mlh1*, *mmp10*, *mos*, *msh2*, *msh3*, *msh6*, *myb*, *myba*, *mybb*, *myc*, *mycl1*, *mycn*, *nfl*, *nf2*, *nme2*, *nras*, *p53*, *pdgfb*, *phb*, *pim1*, *pms1*, *pms2*, *ptc*, *pten*, *raf1*, *rap1a*, *rb1*, *rel*, *ret*, *ros1*, *ski*, *src1*, *tall*, *tgfb2*, *tgfb3*, *tgfb3*, *thra1*, *thrb*, *tiam1*, *timp3*, *tjp1*, *tp53*, *trk*, *vav*, *vhl*, *vil2*, *waf1*, *wnt1*,  
20 *wnt2*, *wl1*, and *yes1* (Hesketh, 1997, In: *The Oncogene and Tumour Suppressor Gene Facts Book*, 2nd Ed., Academic Press; Fishel *et al.*, 1994, *Science* 266:1403-1405).

- Known growth factors include platelet-derived growth factor alpha, platelet-derived growth factor beta (simian sarcoma viral {v-sis} oncogene homolog), thrombopoietin (myeloproliferative leukemia virus oncogene ligand, megakaryocyte  
25 growth and development factor), erythropoietin, B cell growth factor, macrophage stimulating factor 1 (hepatocyte growth factor-like protein), hepatocyte growth factor (hepapoietin A), insulin-like growth factor 1 (somatomedia C), hepatoma-derived growth factor, amphiregulin (schwannoma-derived growth factor), bone morphogenetic proteins 1, 2, 3, 3 beta, and 4, bone morphogenetic protein 7 (osteogenic protein 1), bone  
30 morphogenetic protein 8 (osteogenic protein 2), connective tissue growth factor, connective tissue activation peptide 3, epidermal growth factor (EGF), teratocarcinoma-derived growth factor 1, endothelin, endothelin 2, endothelin 3, stromal cell-derived

factor 1, vascular endothelial growth factor (VEGF), VEGF-B, VEGF-C, placental growth factor (vascular endothelial growth factor-related protein), transforming growth factor alpha, transforming growth factor beta 1 and its precursors, transforming growth factor beta 2 and its precursors, fibroblast growth factor 1 (acidic), fibroblast growth factor 2 (basic), fibroblast growth factor 5 and its precursors, fibroblast growth factor 6 and its precursors, fibroblast growth factor 7 (keratinocyte growth factor), fibroblast growth factor 8 (androgen-induced), fibroblast growth factor 9 (glia-activating factor), pleiotrophin (heparin binding growth factor 8, neurite growth-promoting factor 1), brain-derived neurotrophic factor, and recombinant glial growth factor 2.

10 Known proteases include interleukin-1 beta convertase and its precursors, Mch6 and its precursors, Mch2 isoform alpha, Mch4, Cpp32 isoform alpha, Lice2 gamma cysteine protease, Ich-1S, Ich-1L, Ich-2 and its precursors, TY protease, matrix metalloproteinase 1 (interstitial collagenase), matrix metalloproteinase 2 (gelatinase A, 72kD gelatinase, 72kD type IV collagenase), matrix metalloproteinase 7 (matrilysin),  
15 matrix metalloproteinase 8 (neutrophil collagenase), matrix metalloproteinase 12 (macrophage elastase), matrix metalloproteinase 13 (collagenase 3), metalloproteinase 1, cysteine-rich metalloproteinase (disintegrin) and its precursors, subtilisin-like protease Pc8 and its precursors, chymotrypsin, snake venom-like protease, cathepsin L, cathepsin D (lysosomal aspartyl protease), stromelysin, aminopeptidase N, plasminogen, tissue  
20 plasminogen activator, plasminogen activator inhibitor type II, and urokinase-type plasminogen activator.

Known protein kinases include DAP kinase, serine/threonine protein kinases NIK, PK428, Krs-2, SAK, and EMK, interferon-inducible double stranded RNA dependent protein kinase, FAST kinase, AIM1, IPL1-like midbody-associated protein  
25 kinase-1, NIMA-like protein kinase 1 (NLK1), the cyclin-dependent kinases (cdk1-10), checkpoint kinase Chk1, Nek3 protein kinase, BMK1 beta kinase, Clk1, Clk2, Clk3, extracellular signal-regulated kinases 1, 3, and 6, cdc28 protein kinase 1, cdc28 protein kinase 2, pLK, Myt1, c-Jun N-terminal kinase 2, Cam kinase 1, the MAP kinases, insulin-stimulated protein kinase 1, beta-adrenergic receptor kinase 2, ribosomal protein  
30 S6 kinase, kinase suppressor of ras-1 (KSR1), putative serine/threonine protein kinase Prk, PkB kinase, cAMP-dependent protein kinase, cGMP-dependent protein kinase, type II cGMP-dependent protein kinase, protein kinases Dyrk2, Dyrk3, and Dyrk4, Rho-

associated coiled-coil containing protein kinase p160ROCK, protein tyrosine kinase t-Ror1, Ste20-related kinases, cell adhesion kinase beta, protein kinase 3, stress-activated protein kinase 4, protein kinase Zpk, serine kinase hPAK65, dual specificity mitogen-activated protein kinases 1 and 2, casein kinase I gamma 2, p21-activated protein kinase

5 Pak1, lipid-activated protein kinase PRK2, focal adhesion kinase, dual-specificity tyrosine-phosphorylation regulated kinase, myosin light chain kinase, serine kinases SRPK2, TESK1, and VRK2, B lymphocyte serine/threonine protein kinase, stress-activated protein kinases JNK1 and JNK2, phosphorylase kinase, protein tyrosine kinase Tec, Jak2 kinase, protein kinase Ndr, MEK kinase 3, SHB adaptor protein (a Src

10 homology 2 protein), agammaglobulinaemia protein-tyrosine kinase (Atk), protein kinase ATR, guanylate kinase 1, thrombopoietin receptor and its precursors, DAG kinase epsilon, and kinases encoded by oncogenes or viral oncogenes such as v-fgr (Gardner-Rasheed), v-abl (Abelson murine leukemia viral oncogene homolog 1), v-arg (Abelson murine leukemia viral oncogene homolog, Abelson-related gene), v-fes and v-

15 fps (feline sarcoma viral oncogene and Fujinami avian sarcoma viral oncogene homologs), proto-oncogene *c-cot*, oncogene *pim-1*, and oncogene *mas1*.

Previously known proteins (and, of course, the genes, transcripts, mRNAs, etc. corresponding to those proteins) designated NES1, HE4, and neurosin, are included as markers. NES1 protein is also known as protease serine-like 1 and normal epithelial

20 cell-specific protein, and has been assigned Swiss-Prot accession number O43240 and GenBank accession number AF024605. The amino acid sequence of NES1 protein and the nucleotide sequence of a cDNA encoding it have also been described in U.S. Patent 5,736,377. Association of NES1 protein expression and occurrence of cancer has been described, for example, in U.S. Patent 5,843,694. However, these references (and

25 others, *e.g.* Liu *et al.*, 1996, *Cancer Res.* 56:3371-3379; Luo *et al.*, 1998, *Biochem. Biophys. Res. Comm.* 247:580-586; Goyal *et al.*, 1998, *Cancer Res.* 58:4782-4786) indicate that NES1 expression is down-regulated in cancer patients. In contrast, the present inventors have discovered that NES1 expression is up-regulated in ovarian cancer samples (*e.g.* in later stage {*i.e.* stage 3 or 4} ovarian cancer cell lines).

30 HE4 protein is also known as major epididymis-specific protein E4 and epididymal secretory protein E4, and has been assigned Swiss-Prot accession number Q14508 and GenBank accession number X63187. The amino acid sequence and the

corresponding cDNA nucleotide sequence were also disclosed in Kirchhoff *et al.* (1991) *Biol. Reprod.* 45:350-357. A possible association between expression of HE4 and occurrence of ovarian cancer was disclosed, for example in Wang *et al.* (1999) *Gene* 229:101-108.

5       Neurosin is also known as protease M, zyme, and SP59, and has been assigned Swiss-Prot accession number Q92876 and GenBank accession number U62801. The amino acid sequence of neurosin and the corresponding cDNA nucleotide sequence were also disclosed in Anisowicz *et al.* (1996) *Mol. Med.* 2:624-636. The same reference discloses a possible association between expression of neurosin and  
10       occurrence of ovarian cancer.

It is recognized that the compositions, kits, and methods of the invention will be of particular utility to patients having an enhanced risk of developing ovarian cancer and their medical advisors. Patients recognized as having an enhanced risk of developing ovarian cancer include, for example, patients having a familial history of ovarian cancer,  
15       patients identified as having a mutant oncogene (*i.e.* at least one allele), and patients of advancing age (*i.e.* women older than about 50 or 60 years).

The level of expression of a marker in normal (*i.e.* non-cancerous) human ovarian tissue can be assessed in a variety of ways. In one embodiment, this normal level of expression is assessed by assessing the level of expression of the marker in a  
20       portion of ovarian cells which appears to be non-cancerous and by comparing this normal level of expression with the level of expression in a portion of the ovarian cells which is suspected of being cancerous. For example, when laparoscopy or other medical procedure, reveals the presence of a lump on one portion of a patient's ovary, but not on another portion of the same ovary or on the other ovary, the normal level of  
25       expression of a marker may be assessed using one or both of the non-affected ovary and a non-affected portion of the affected ovary, and this normal level of expression may be compared with the level of expression of the same marker in an affected portion (*i.e.* the lump) of the affected ovary. Alternately, and particularly as further information becomes available as a result of routine performance of the methods described herein,  
30       population-average values for normal expression of the markers of the invention may be used. In other embodiments, the 'normal' level of expression of a marker may be determined by assessing expression of the marker in a patient sample obtained from a

non-cancer-afflicted patient, from a patient sample obtained from a patient before the suspected onset of ovarian cancer in the patient, from archived patient samples, and the like.

The invention includes compositions, kits, and methods for assessing the  
5 presence of ovarian cancer cells in a sample (*e.g.* an archived tissue sample or a sample obtained from a patient). These compositions, kits, and methods are substantially the same as those described above, except that, where necessary, the compositions, kits, and methods are adapted for use with samples other than patient samples. For example, when the sample to be used is a paraffinized, archived human tissue sample, it can be  
10 necessary to adjust the ratio of compounds in the compositions of the invention, in the kits of the invention, or the methods used to assess levels of marker expression in the sample. Such methods are well known in the art and within the skill of the ordinary artisan.

The invention includes a kit for assessing the presence of ovarian cancer cells  
15 (*e.g.* in a sample such as a patient sample). The kit comprises a plurality of reagents, each of which is capable of binding specifically with a nucleic acid or polypeptide corresponding to a marker of the invention. Suitable reagents for binding with a polypeptide corresponding to a marker of the invention include antibodies, antibody derivatives, antibody fragments, and the like. Suitable reagents for binding with a  
20 nucleic acid (*e.g.* a genomic DNA, an mRNA, a spliced mRNA, a cDNA, or the like) include complementary nucleic acids. For example, the nucleic acid reagents may include oligonucleotides (labeled or non-labeled) fixed to a substrate, labeled oligonucleotides not bound with a substrate, pairs of PCR primers, molecular beacon probes, and the like.

25 The kit of the invention may optionally comprise additional components useful for performing the methods of the invention. By way of example, the kit may comprise fluids (*e.g.* SSC buffer) suitable for annealing complementary nucleic acids or for binding an antibody with a protein with which it specifically binds, one or more sample compartments, an instructional material which describes performance of a method of the  
30 invention, a sample of normal ovarian cells, a sample of ovarian cancer cells, and the like.

The invention also includes a method of making an isolated hybridoma which produces an antibody useful for assessing whether patient is afflicted with an ovarian cancer. In this method, a protein corresponding to a marker of the invention or a fragment of the protein is isolated (*e.g.* by purification from a cell in which it is expressed or by transcription and translation of a nucleic acid encoding the protein *in vivo* or *in vitro* using known methods). A vertebrate, preferably a mammal such as a mouse, rat, rabbit, or sheep, is immunized using the isolated protein or fragment thereof. The vertebrate may optionally (and preferably) be immunized at least one additional time with the isolated protein or fragment, so that the vertebrate exhibits a robust immune response to the protein. Splenocytes are isolated from the immunized vertebrate and fused with an immortalized cell line to form hybridomas, using any of a variety of methods well known in the art. Hybridomas formed in this manner are then screened using standard methods to identify one or more hybridomas which produce an antibody which specifically binds with the protein. The invention also includes hybridomas made by this method and antibodies made using such hybridomas. An antibody of the invention may also be used as a therapeutic agent for treating cancers, particularly ovarian cancers (see *e.g.*, Table 8).

The invention also includes a method of assessing the efficacy of a test compound for inhibiting ovarian cancer cells. As described above, differences in the level of expression of the markers of the invention correlate with the cancerous state of ovarian cells. Although it is recognized that changes in the levels of expression of certain of the markers of the invention likely result from the cancerous state of ovarian cells, it is likewise recognized that changes in the levels of expression of other of the markers of the invention induce, maintain, and promote the cancerous state of those cells. Thus, compounds which inhibit an ovarian cancer in a patient will cause the level of expression of one or more of the markers of the invention to change to a level nearer the normal level of expression for that marker (*i.e.* the level of expression for the marker in non-cancerous ovarian cells).

This method thus comprises comparing expression of a marker in a first ovarian cell sample and maintained in the presence of the test compound and expression of the marker in a second ovarian cell sample and maintained in the absence of the test compound. A significant increase in the level of expression of a marker listed in Table

3A, 5, 7C and/or 7E, or a significant decrease in the level of expression of a marker listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and/or 8, is an indication that the test compound inhibits ovarian cancer. The ovarian cell samples may, for example, be aliquots of a single sample of normal ovarian cells obtained from a patient, pooled  
5 samples of normal ovarian cells obtained from a patient, cells of a normal ovarian cell line, aliquots of a single sample of ovarian cancer cells obtained from a patient, pooled samples of ovarian cancer cells obtained from a patient, cells of an ovarian cancer cell line, or the like. In one embodiment, the samples are ovarian cancer cells obtained from a patient and a plurality of compounds known to be effective for inhibiting various  
10 ovarian cancers are tested in order to identify the compound which is likely to best inhibit the ovarian cancer in the patient.

This method may likewise be used to assess the efficacy of a therapy for inhibiting ovarian cancer in a patient. In this method, the level of expression of one or more markers of the invention in a pair of samples (one subjected to the therapy, the  
15 other not subjected to the therapy) is assessed. As with the method of assessing the efficacy of test compounds, if the therapy induces a significant decrease in the level of expression of a marker listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and/or 8, or blocks induction of a marker listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and/or 8, or if the therapy induces a significant enhancement of the level of expression of a  
20 marker listed in Tables 3A, 5, 7C and 7E, then the therapy is efficacious for inhibiting ovarian cancer. As above, if samples from a selected patient are used in this method, then alternative therapies can be assessed *in vitro* in order to select a therapy most likely to be efficacious for inhibiting ovarian cancer in the patient.

As described herein, ovarian cancer in patients is associated with an increase in  
25 the level of expression of one or more markers listed in either or both of Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and/or 8, with a decrease in the level of expression of one or more markers listed in Table 3A, 5, 7C and 7E, or with both. While, as discussed above, some of these changes in expression level result from occurrence of the ovarian cancer, others of these changes induce, maintain, and promote the cancerous state of ovarian  
30 cancer cells. Thus, ovarian cancer characterized by an increase in the level of expression of one or more markers listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and/or 8 can be inhibited by inhibiting expression of those markers. Likewise, ovarian

cancer characterized by a decrease in the level of expression of one or more markers listed in Table 3A, 5, 7C and 7E can be inhibited by enhancing expression of those markers.

Expression of a marker listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8  
5 can be inhibited in a number of ways generally known in the art. For example, an antisense oligonucleotide can be provided to the ovarian cancer cells in order to inhibit transcription, translation, or both, of the marker(s). Alternately, a polynucleotide encoding an antibody, an antibody derivative, or an antibody fragment which specifically binds the protein corresponding to the marker, and operably linked with an  
10 appropriate promoter/regulator region, can be provided to the cell in order to generate intracellular antibodies which will inhibit the function or activity of the protein. The expression and/or function of a marker may also be inhibited by treating the ovarian cancer cell with a heterologous antibody or antibody derivative that specifically binds the protein corresponding to the marker. Using the methods described herein, a variety  
15 of molecules, particularly including molecules sufficiently small that they are able to cross the cell membrane, can be screened in order to identify molecules which inhibit expression of the marker(s). The compound so identified can be provided to the patient in order to inhibit expression of the marker(s) in the ovarian cancer cells of the patient.

Expression of a marker listed in Tables 3A, 5, 7C and 7E can be enhanced in a  
20 number of ways generally known in the art. For example, a polynucleotide encoding the marker and operably linked with an appropriate promoter/regulator region can be provided to ovarian cancer cells of the patient in order to induce enhanced expression of the protein (and mRNA) corresponding to the marker therein. Alternatively, if the protein is capable of crossing the cell membrane, inserting itself in the cell membrane,  
25 or is normally a secreted protein, then expression of the protein can be enhanced by providing the protein (*e.g.* directly or by way of the bloodstream or another ovary-associated fluid) to ovarian cancer cells in the patient.

As described above, the cancerous state of human ovarian cells is correlated with changes in the levels of expression of the markers of the invention. Thus, compounds  
30 which induce increased expression of one or more of the markers listed in either or both of Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8, decreased expression of one or more of the markers listed in either or both of Tables 3A, 5, 7C and 7E, or both, can induce



- 36 -

ovarian cell carcinogenesis. The invention includes a method for assessing the human ovarian cell carcinogenic potential of a test compound. This method comprises maintaining separate aliquots of human ovarian cells in the presence and absence of the test compound. Expression of a marker of the invention in each of the aliquots is compared. A significant increase in the level of expression of a marker listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8, or a significant decrease in the level of expression of a marker listed in Tables 3A, 5, 7C and 7E in the aliquot maintained in the presence of the test compound (relative to the aliquot maintained in the absence of the test compound) is an indication that the test compound possesses human ovarian cell carcinogenic potential. The relative carcinogenic potentials of various test compounds can be assessed by comparing the degree of enhancement or inhibition of the level of expression of the relevant markers, by comparing the number of markers for which the level of expression is enhanced or inhibited, or by comparing both.

Various aspects of the invention are described in further detail in the following subsections.

#### I. Isolated Nucleic Acid Molecules

One aspect of the invention pertains to isolated nucleic acid molecules that correspond to a marker of the invention, including nucleic acids which encode a polypeptide corresponding to a marker of the invention or a portion of such a polypeptide. Isolated nucleic acids of the invention also include nucleic acid molecules sufficient for use as hybridization probes to identify nucleic acid molecules that correspond to a marker of the invention, including nucleic acids which encode a polypeptide corresponding to a marker of the invention, and fragments of such nucleic acid molecules, *e.g.*, those suitable for use as PCR primers for the amplification or mutation of nucleic acid molecules. As used herein, the term "nucleic acid molecule" is intended to include DNA molecules (*e.g.*, cDNA or genomic DNA) and RNA molecules (*e.g.*, mRNA) and analogs of the DNA or RNA generated using nucleotide analogs. The nucleic acid molecule can be single-stranded or double-stranded, but preferably is double-stranded DNA.

An "isolated" nucleic acid molecule is one which is separated from other nucleic acid molecules which are present in the natural source of the nucleic acid molecule. Preferably, an "isolated" nucleic acid molecule comprises a protein-coding sequence and is free of sequences which naturally flank the coding sequence in the genomic DNA of the organism from which the nucleic acid is derived. For example, in various  
5       embodiments, the isolated nucleic acid molecule can contain less than about 5 kB, 4 kB, 3 kB, 2 kB, 1 kB, 0.5 kB or 0.1 kB of nucleotide sequences which naturally flank the nucleic acid molecule in genomic DNA of the cell from which the nucleic acid is derived. Moreover, an "isolated" nucleic acid molecule, such as a cDNA molecule, can  
10       be substantially free of other cellular material, or culture medium when produced by recombinant techniques, or substantially free of chemical precursors or other chemicals when chemically synthesized.

A nucleic acid molecule of the present invention, *e.g.*, a nucleic acid encoding a protein corresponding to a marker listed in one or more of Tables 1-11, can be isolated  
15       using standard molecular biology techniques and the sequence information in the database records described herein. Using all or a portion of such nucleic acid sequences, nucleic acid molecules of the invention can be isolated using standard hybridization and cloning techniques (*e.g.*, as described in Sambrook *et al.*, ed., *Molecular Cloning: A Laboratory Manual*, 2nd ed., Cold Spring Harbor Laboratory Press, Cold Spring  
20       Harbor, NY, 1989).

A nucleic acid molecule of the invention can be amplified using cDNA, mRNA, or genomic DNA as a template and appropriate oligonucleotide primers according to standard PCR amplification techniques. The nucleic acid so amplified can be cloned into an appropriate vector and characterized by DNA sequence analysis. Furthermore,  
25       oligonucleotides corresponding to all or a portion of a nucleic acid molecule of the invention can be prepared by standard synthetic techniques, *e.g.*, using an automated DNA synthesizer.

In another preferred embodiment, an isolated nucleic acid molecule of the invention comprises a nucleic acid molecule which has a nucleotide sequence  
30       complementary to the nucleotide sequence of a nucleic acid corresponding to a marker of the invention or to the nucleotide sequence of a nucleic acid encoding a protein which corresponds to a marker of the invention. A nucleic acid molecule which is

complementary to a given nucleotide sequence is one which is sufficiently complementary to the given nucleotide sequence that it can hybridize to the given nucleotide sequence thereby forming a stable duplex.

Moreover, a nucleic acid molecule of the invention can comprise only a portion  
5 of a nucleic acid sequence, wherein the full length nucleic acid sequence comprises a marker of the invention or which encodes a polypeptide corresponding to a marker of the invention. Such nucleic acids can be used, for example, as a probe or primer. The probe/primer typically is used as one or more substantially purified oligonucleotides. The oligonucleotide typically comprises a region of nucleotide sequence that hybridizes  
10 under stringent conditions to at least about 7, preferably about 15, more preferably about 25, 50, 75, 100, 125, 150, 175, 200, 250, 300, 350, or 400 or more consecutive nucleotides of a nucleic acid of the invention.

Probes based on the sequence of a nucleic acid molecule of the invention can be used to detect transcripts or genomic sequences corresponding to one or more markers  
15 of the invention. The probe comprises a label group attached thereto, *e.g.*, a radioisotope, a fluorescent compound, an enzyme, or an enzyme co-factor. Such probes can be used as part of a diagnostic test kit for identifying cells or tissues which mis-express the protein, such as by measuring levels of a nucleic acid molecule encoding the protein in a sample of cells from a subject, *e.g.*, detecting mRNA levels or determining  
20 whether a gene encoding the protein has been mutated or deleted.

The invention further encompasses nucleic acid molecules that differ, due to degeneracy of the genetic code, from the nucleotide sequence of nucleic acids encoding a protein which corresponds to a marker of the invention, and thus encode the same protein.

25 In addition to the nucleotide sequences described in the GenBank and IMAGE Consortium database records described herein, it will be appreciated by those skilled in the art that DNA sequence polymorphisms that lead to changes in the amino acid sequence can exist within a population (*e.g.*, the human population). Such genetic polymorphisms can exist among individuals within a population due to natural allelic  
30 variation. An allele is one of a group of genes which occur alternatively at a given genetic locus. In addition, it will be appreciated that DNA polymorphisms that affect

RNA expression levels can also exist that may affect the overall expression level of that gene (e.g., by affecting regulation or degradation).

As used herein, the phrase "allelic variant" refers to a nucleotide sequence which occurs at a given locus or to a polypeptide encoded by the nucleotide sequence.

- 5       As used herein, the terms "gene" and "recombinant gene" refer to nucleic acid molecules comprising an open reading frame encoding a polypeptide corresponding to a marker of the invention. Such natural allelic variations can typically result in 0.1 –0.5 % variance in the nucleotide sequence of a given gene. Alternative alleles can be identified by sequencing the gene of interest in a number of different individuals.
- 10       This can be readily carried out by using hybridization probes to identify the same genetic locus in a variety of individuals. Any and all such nucleotide variations and resulting amino acid polymorphisms or variations that are the result of natural allelic variation and that do not alter the functional activity are intended to be within the scope of the invention.
- 15       In another embodiment, an isolated nucleic acid molecule of the invention is at least 7, 15, 20, 25, 30, 40, 60, 80, 100, 150, 200, 250, 300, 350, 400, 450, 550, 650, 700, 800, 900, 1000, 1200, 1400, 1600, 1800, 2000, 2200, 2400, 2600, 2800, 3000, 3500, 4000, 4500, or more nucleotides in length and hybridizes under stringent conditions to a nucleic acid corresponding to a marker of the invention or to a nucleic acid encoding a
- 20       protein corresponding to a marker of the invention. As used herein, the term "hybridizes under stringent conditions" is intended to describe conditions for hybridization and washing under which nucleotide sequences at least 75% (80%, 85%, preferably 90%) identical to each other typically remain hybridized to each other. Such stringent conditions are known to those skilled in the art and can be found in sections 6.3.1-6.3.6
- 25       of *Current Protocols in Molecular Biology*, John Wiley & Sons, N.Y. (1989). A preferred, non-limiting example of stringent hybridization conditions for annealing two single-stranded DNA each of which is at least about 100 bases in length and/or for annealing a single-stranded DNA and a single-stranded RNA each of which is at least about 100 bases in length, are hybridization in 6X sodium chloride/sodium citrate (SSC)
- 30       at about 45°C, followed by one or more washes in 0.2X SSC, 0.1% SDS at 50-65°C. Further preferred hybridization conditions are taught in Lockhart, *et al.*, *Nature Biotechnology*, Volume 14, 1996 August:1675-1680; Breslauer, *et al.*, *Proc. Natl. Acad.*

Sci. USA, Volume 83, 1986 June: 3746-3750; Van Ness, *et al.*, Nucleic Acids Research, Volume 19, No. 19, 1991 September: 5143-5151; McGraw, *et al.*, BioTechniques, Volume 8, No. 6 1990: 674-678; and Milner, *et al.*, Nature Biotechnology, Volume 15, 1997 June: 537-541, all expressly incorporated by reference.

5           In addition to naturally-occurring allelic variants of a nucleic acid molecule of the invention that can exist in the population, the skilled artisan will further appreciate that sequence changes can be introduced by mutation thereby leading to changes in the amino acid sequence of the encoded protein, without altering the biological activity of the protein encoded thereby. For example, one can make nucleotide substitutions  
10   leading to amino acid substitutions at "non-essential" amino acid residues. A "non-essential" amino acid residue is a residue that can be altered from the wild-type sequence without altering the biological activity, whereas an "essential" amino acid residue is required for biological activity. For example, amino acid residues that are not conserved or only semi-conserved among homologs of various species may be non-  
15   essential for activity and thus would be likely targets for alteration. Alternatively, amino acid residues that are conserved among the homologs of various species (*e.g.*, murine and human) may be essential for activity and thus would not be likely targets for alteration.

          Accordingly, another aspect of the invention pertains to nucleic acid molecules  
20   encoding a polypeptide of the invention that contain changes in amino acid residues that are not essential for activity. Such polypeptides differ in amino acid sequence from the naturally-occurring proteins which correspond to the markers of the invention, yet retain biological activity. In one embodiment, such a protein has an amino acid sequence that is at least about 40% identical, 50%, 60%, 70%, 80%, 90%, 95%, or 98% identical to the  
25   amino acid sequence of one of the proteins which correspond to the markers of the invention.

          An isolated nucleic acid molecule encoding a variant protein can be created by introducing one or more nucleotide substitutions, additions or deletions into the nucleotide sequence of nucleic acids of the invention, such that one or more amino acid  
30   residue substitutions, additions, or deletions are introduced into the encoded protein. Mutations can be introduced by standard techniques, such as site-directed mutagenesis and PCR-mediated mutagenesis. Preferably, conservative amino acid substitutions are

- made at one or more predicted non-essential amino acid residues. A "conservative amino acid substitution" is one in which the amino acid residue is replaced with an amino acid residue having a similar side chain. Families of amino acid residues having similar side chains have been defined in the art. These families include amino acids
- 5 with basic side chains (*e.g.*, lysine, arginine, histidine), acidic side chains (*e.g.*, aspartic acid, glutamic acid), uncharged polar side chains (*e.g.*, glycine, asparagine, glutamine, serine, threonine, tyrosine, cysteine), non-polar side chains (*e.g.*, alanine, valine, leucine, isoleucine, proline, phenylalanine, methionine, tryptophan), beta-branched side chains (*e.g.*, threonine, valine, isoleucine) and aromatic side chains (*e.g.*, tyrosine,
- 10 phenylalanine, tryptophan, histidine). Alternatively, mutations can be introduced randomly along all or part of the coding sequence, such as by saturation mutagenesis, and the resultant mutants can be screened for biological activity to identify mutants that retain activity. Following mutagenesis, the encoded protein can be expressed recombinantly and the activity of the protein can be determined.
- 15 The present invention encompasses antisense nucleic acid molecules, *i.e.*, molecules which are complementary to a sense nucleic acid of the invention, *e.g.*, complementary to the coding strand of a double-stranded cDNA molecule corresponding to a marker of the invention or complementary to an mRNA sequence corresponding to a marker of the invention. Accordingly, an antisense nucleic acid of
- 20 the invention can hydrogen bond to (*i.e.* anneal with) a sense nucleic acid of the invention. The antisense nucleic acid can be complementary to an entire coding strand, or to only a portion thereof, *e.g.*, all or part of the protein coding region (or open reading frame). An antisense nucleic acid molecule can also be antisense to all or part of a non-coding region of the coding strand of a nucleotide sequence encoding a polypeptide of
- 25 the invention. The non-coding regions ("5' and 3' untranslated regions") are the 5' and 3' sequences which flank the coding region and are not translated into amino acids.
- An antisense oligonucleotide can be, for example, about 5, 10, 15, 20, 25, 30, 35, 40, 45, or 50 or more nucleotides in length. An antisense nucleic acid of the invention can be constructed using chemical synthesis and enzymatic ligation reactions using
- 30 procedures known in the art. For example, an antisense nucleic acid (*e.g.*, an antisense oligonucleotide) can be chemically synthesized using naturally occurring nucleotides or variously modified nucleotides designed to increase the biological stability of the

molecules or to increase the physical stability of the duplex formed between the antisense and sense nucleic acids, *e.g.*, phosphorothioate derivatives and acridine substituted nucleotides can be used. Examples of modified nucleotides which can be used to generate the antisense nucleic acid include 5-fluorouracil, 5-bromouracil, 5-chlorouracil, 5-iodouracil, hypoxanthine, xanthine, 4-acetylcytosine, 5-  
5 (carboxyhydroxymethyl) uracil, 5-carboxymethylaminomethyl-2-thiouridine, 5-carboxymethylaminomethyluracil, dihydrouracil, beta-D-galactosylqueosine, inosine, N6-isopentenyladenine, 1-methylguanine, 1-methylinosine, 2,2-dimethylguanine, 2-methyladenine, 2-methylguanine, 3-methylcytosine, 5-methylcytosine, N6-adenine, 7-  
10 methylguanine, 5-methylaminomethyluracil, 5-methoxyaminomethyl-2-thiouracil, beta-D-mannosylqueosine, 5'-methoxycarboxymethyluracil, 5-methoxyuracil, 2-methylthio-N6-isopentenyladenine, uracil-5-oxyacetic acid (v), wybutoxosine, pseudouracil, queosine, 2-thiocytosine, 5-methyl-2-thiouracil, 2-thiouracil, 4-thiouracil, 5-methyluracil, uracil-5-oxyacetic acid methylester, uracil-5-oxyacetic acid (v), 5-methyl-  
15 2-thiouracil, 3-(3-amino-3-N-2-carboxypropyl) uracil, (acp3)w, and 2,6-diaminopurine. Alternatively, the antisense nucleic acid can be produced biologically using an expression vector into which a nucleic acid has been sub-cloned in an antisense orientation (*i.e.*, RNA transcribed from the inserted nucleic acid will be of an antisense orientation to a target nucleic acid of interest, described further in the following  
20 subsection).

The antisense nucleic acid molecules of the invention are typically administered to a subject or generated *in situ* such that they hybridize with or bind to cellular mRNA and/or genomic DNA encoding a polypeptide corresponding to a selected marker of the invention to thereby inhibit expression of the marker, *e.g.*, by inhibiting transcription  
25 and/or translation. The hybridization can be by conventional nucleotide complementarity to form a stable duplex, or, for example, in the case of an antisense nucleic acid molecule which binds to DNA duplexes, through specific interactions in the major groove of the double helix. Examples of a route of administration of antisense nucleic acid molecules of the invention includes direct injection at a tissue site or  
30 infusion of the antisense nucleic acid into an ovary-associated body fluid. Alternatively, antisense nucleic acid molecules can be modified to target selected cells and then administered systemically. For example, for systemic administration, antisense

molecules can be modified such that they specifically bind to receptors or antigens expressed on a selected cell surface, *e.g.*, by linking the antisense nucleic acid molecules to peptides or antibodies which bind to cell surface receptors or antigens. The antisense nucleic acid molecules can also be delivered to cells using the vectors described herein.

- 5 To achieve sufficient intracellular concentrations of the antisense molecules, vector constructs in which the antisense nucleic acid molecule is placed under the control of a strong pol II or pol III promoter are preferred.

An antisense nucleic acid molecule of the invention can be an  $\alpha$ -anomeric nucleic acid molecule. An  $\alpha$ -anomeric nucleic acid molecule forms specific double-  
10 stranded hybrids with complementary RNA in which, contrary to the usual  $\alpha$ -units, the strands run parallel to each other (Gaultier *et al.*, 1987, *Nucleic Acids Res.* 15:6625-6641). The antisense nucleic acid molecule can also comprise a 2'-o-methylribonucleotide (Inoue *et al.*, 1987, *Nucleic Acids Res.* 15:6131-6148) or a chimeric RNA-DNA analogue (Inoue *et al.*, 1987, *FEBS Lett.* 215:327-330).

- 15 The invention also encompasses ribozymes. Ribozymes are catalytic RNA molecules with ribonuclease activity which are capable of cleaving a single-stranded nucleic acid, such as an mRNA, to which they have a complementary region. Thus, ribozymes (*e.g.*, hammerhead ribozymes as described in Haselhoff and Gerlach, 1988, *Nature* 334:585-591) can be used to catalytically cleave mRNA transcripts to thereby  
20 inhibit translation of the protein encoded by the mRNA. A ribozyme having specificity for a nucleic acid molecule encoding a polypeptide corresponding to a marker of the invention can be designed based upon the nucleotide sequence of a cDNA corresponding to the marker. For example, a derivative of a *Tetrahymena* L-19 IVS RNA can be constructed in which the nucleotide sequence of the active site is  
25 complementary to the nucleotide sequence to be cleaved (see Cech *et al.* U.S. Patent No. 4,987,071; and Cech *et al.* U.S. Patent No. 5,116,742). Alternatively, an mRNA encoding a polypeptide of the invention can be used to select a catalytic RNA having a specific ribonuclease activity from a pool of RNA molecules (see, *e.g.*, Bartel and Szostak, 1993, *Science* 261:1411-1418).

- 30 The invention also encompasses nucleic acid molecules which form triple helical structures. For example, expression of a polypeptide of the invention can be inhibited by targeting nucleotide sequences complementary to the regulatory region of the gene



encoding the polypeptide (*e.g.*, the promoter and/or enhancer) to form triple helical structures that prevent transcription of the gene in target cells. See generally Helene (1991) *Anticancer Drug Des.* 6(6):569-84; Helene (1992) *Ann. N.Y. Acad. Sci.* 660:27-36; and Maher (1992) *Bioassays* 14(12):807-15.

5           In various embodiments, the nucleic acid molecules of the invention can be modified at the base moiety, sugar moiety or phosphate backbone to improve, *e.g.*, the stability, hybridization, or solubility of the molecule. For example, the deoxyribose phosphate backbone of the nucleic acids can be modified to generate peptide nucleic acids (see Hyrup *et al.*, 1996, *Bioorganic & Medicinal Chemistry* 4(1): 5-23). As used  
10   herein, the terms "peptide nucleic acids" or "PNAs" refer to nucleic acid mimics, *e.g.*, DNA mimics, in which the deoxyribose phosphate backbone is replaced by a pseudopeptide backbone and only the four natural nucleobases are retained. The neutral backbone of PNAs has been shown to allow for specific hybridization to DNA and RNA under conditions of low ionic strength. The synthesis of PNA oligomers can be  
15   performed using standard solid phase peptide synthesis protocols as described in Hyrup *et al.* (1996), *supra*; Perry-O'Keefe *et al.* (1996) *Proc. Natl. Acad. Sci. USA* 93:14670-675.

          PNAs can be used in therapeutic and diagnostic applications. For example, PNAs can be used as antisense or antigene agents for sequence-specific modulation of  
20   gene expression by, *e.g.*, inducing transcription or translation arrest or inhibiting replication. PNAs can also be used, *e.g.*, in the analysis of single base pair mutations in a gene by, *e.g.*, PNA directed PCR clamping; as artificial restriction enzymes when used in combination with other enzymes, *e.g.*, S1 nucleases (Hyrup (1996), *supra*; or as probes or primers for DNA sequence and hybridization (Hyrup, 1996, *supra*; Perry-  
25   O'Keefe *et al.*, 1996, *Proc. Natl. Acad. Sci. USA* 93:14670-675).

          In another embodiment, PNAs can be modified, *e.g.*, to enhance their stability or cellular uptake, by attaching lipophilic or other helper groups to PNA, by the formation of PNA-DNA chimeras, or by the use of liposomes or other techniques of drug delivery known in the art. For example, PNA-DNA chimeras can be generated which can  
30   combine the advantageous properties of PNA and DNA. Such chimeras allow DNA recognition enzymes, *e.g.*, RNASE H and DNA polymerases, to interact with the DNA portion while the PNA portion would provide high binding affinity and specificity.

PNA-DNA chimeras can be linked using linkers of appropriate lengths selected in terms of base stacking, number of bonds between the nucleobases, and orientation (Hyrup, 1996, *supra*). The synthesis of PNA-DNA chimeras can be performed as described in Hyrup (1996), *supra*, and Finn *et al.* (1996) *Nucleic Acids Res.* 24(17):3357-63. For  
5 example, a DNA chain can be synthesized on a solid support using standard phosphoramidite coupling chemistry and modified nucleoside analogs. Compounds such as 5'-(4-methoxytrityl)amino-5'-deoxy-thymidine phosphoramidite can be used as a link between the PNA and the 5' end of DNA (Mag *et al.*, 1989, *Nucleic Acids Res.* 17:5973-88). PNA monomers are then coupled in a step-wise manner to produce a  
10 chimeric molecule with a 5' PNA segment and a 3' DNA segment (Finn *et al.*, 1996, *Nucleic Acids Res.* 24(17):3357-63). Alternatively, chimeric molecules can be synthesized with a 5' DNA segment and a 3' PNA segment (Peterser *et al.*, 1975, *Bioorganic Med. Chem. Lett.* 5:1119-11124).

In other embodiments, the oligonucleotide can include other appended groups  
15 such as peptides (*e.g.*, for targeting host cell receptors *in vivo*), or agents facilitating transport across the cell membrane (see, *e.g.*, Letsinger *et al.*, 1989, *Proc. Natl. Acad. Sci. USA* 86:6553-6556; Lemaitre *et al.*, 1987, *Proc. Natl. Acad. Sci. USA* 84:648-652; PCT Publication No. WO 88/09810) or the blood-brain barrier (see, *e.g.*, PCT Publication No. WO 89/10134). In addition, oligonucleotides can be modified with  
20 hybridization-triggered cleavage agents (see, *e.g.*, Krol *et al.*, 1988, *Bio/Techniques* 6:958-976) or intercalating agents (see, *e.g.*, Zon, 1988, *Pharm. Res.* 5:539-549). To this end, the oligonucleotide can be conjugated to another molecule, *e.g.*, a peptide, hybridization triggered cross-linking agent, transport agent, hybridization-triggered cleavage agent, etc.

25 The invention also includes molecular beacon nucleic acids having at least one region which is complementary to a nucleic acid of the invention, such that the molecular beacon is useful for quantitating the presence of the nucleic acid of the invention in a sample. A "molecular beacon" nucleic acid is a nucleic acid comprising a pair of complementary regions and having a fluorophore and a fluorescent quencher  
30 associated therewith. The fluorophore and quencher are associated with different portions of the nucleic acid in such an orientation that when the complementary regions are annealed with one another, fluorescence of the fluorophore is quenched by the

quencher. When the complementary regions of the nucleic acid are not annealed with one another, fluorescence of the fluorophore is quenched to a lesser degree. Molecular beacon nucleic acids are described, for example, in U.S. Patent 5,876,930.

## 5 II. Isolated Proteins and Antibodies

One aspect of the invention pertains to isolated proteins which correspond to individual markers of the invention, and biologically active portions thereof, as well as polypeptide fragments suitable for use as immunogens to raise antibodies directed against a polypeptide corresponding to a marker of the invention. In one embodiment,  
10 the native polypeptide corresponding to a marker can be isolated from cells or tissue sources by an appropriate purification scheme using standard protein purification techniques. In another embodiment, polypeptides corresponding to a marker of the invention are produced by recombinant DNA techniques. Alternative to recombinant expression, a polypeptide corresponding to a marker of the invention can be synthesized  
15 chemically using standard peptide synthesis techniques.

An "isolated" or "purified" protein or biologically active portion thereof is substantially free of cellular material or other contaminating proteins from the cell or tissue source from which the protein is derived, or substantially free of chemical precursors or other chemicals when chemically synthesized. The language  
20 "substantially free of cellular material" includes preparations of protein in which the protein is separated from cellular components of the cells from which it is isolated or recombinantly produced. Thus, protein that is substantially free of cellular material includes preparations of protein having less than about 30%, 20%, 10%, or 5% (by dry weight) of heterologous protein (also referred to herein as a "contaminating protein").  
25 When the protein or biologically active portion thereof is recombinantly produced, it is also preferably substantially free of culture medium, *i.e.*, culture medium represents less than about 20%, 10%, or 5% of the volume of the protein preparation. When the protein is produced by chemical synthesis, it is preferably substantially free of chemical precursors or other chemicals, *i.e.*, it is separated from chemical precursors or other  
30 chemicals which are involved in the synthesis of the protein. Accordingly such preparations of the protein have less than about 30%, 20%, 10%, 5% (by dry weight) of chemical precursors or compounds other than the polypeptide of interest.

Biologically active portions of a polypeptide corresponding to a marker of the invention include polypeptides comprising amino acid sequences sufficiently identical to or derived from the amino acid sequence of the protein corresponding to the marker (e.g., the amino acid sequence listed in the GenBank and IMAGE Consortium database records described herein), which include fewer amino acids than the full length protein, and exhibit at least one activity of the corresponding full-length protein. Typically, biologically active portions comprise a domain or motif with at least one activity of the corresponding protein. A biologically active portion of a protein of the invention can be a polypeptide which is, for example, 10, 25, 50, 100 or more amino acids in length.

Moreover, other biologically active portions, in which other regions of the protein are deleted, can be prepared by recombinant techniques and evaluated for one or more of the functional activities of the native form of a polypeptide of the invention.

Preferred polypeptides have the amino acid sequence listed in the one of the GenBank and IMAGE Consortium database records described herein. Other useful proteins are substantially identical (e.g., at least about 40%, preferably 50%, 60%, 70%, 80%, 90%, 95%, or 99%) to one of these sequences and retain the functional activity of the protein of the corresponding naturally-occurring protein yet differ in amino acid sequence due to natural allelic variation or mutagenesis.

To determine the percent identity of two amino acid sequences or of two nucleic acids, the sequences are aligned for optimal comparison purposes (e.g., gaps can be introduced in the sequence of a first amino acid or nucleic acid sequence for optimal alignment with a second amino or nucleic acid sequence). The amino acid residues or nucleotides at corresponding amino acid positions or nucleotide positions are then compared. When a position in the first sequence is occupied by the same amino acid residue or nucleotide as the corresponding position in the second sequence, then the molecules are identical at that position. The percent identity between the two sequences is a function of the number of identical positions shared by the sequences (i.e., % identity = # of identical positions/total # of positions (e.g., overlapping positions) x100). In one embodiment the two sequences are the same length.

The determination of percent identity between two sequences can be accomplished using a mathematical algorithm. A preferred, non-limiting example of a mathematical algorithm utilized for the comparison of two sequences is the algorithm of

Karlin and Altschul (1990) *Proc. Natl. Acad. Sci. USA* 87:2264-2268, modified as in Karlin and Altschul (1993) *Proc. Natl. Acad. Sci. USA* 90:5873-5877. Such an algorithm is incorporated into the NBLAST and XBLAST programs of Altschul, *et al.* (1990) *J. Mol. Biol.* 215:403-410. BLAST nucleotide searches can be performed with the NBLAST program, score = 100, wordlength = 12 to obtain nucleotide sequences homologous to a nucleic acid molecules of the invention. BLAST protein searches can be performed with the XBLAST program, score = 50, wordlength = 3 to obtain amino acid sequences homologous to a protein molecules of the invention. To obtain gapped alignments for comparison purposes, Gapped BLAST can be utilized as described in Altschul *et al.* (1997) *Nucleic Acids Res.* 25:3389-3402. Alternatively, PSI-Blast can be used to perform an iterated search which detects distant relationships between molecules. When utilizing BLAST, Gapped BLAST, and PSI-Blast programs, the default parameters of the respective programs (*e.g.*, XBLAST and NBLAST) can be used. See <http://www.ncbi.nlm.nih.gov>. Another preferred, non-limiting example of a mathematical algorithm utilized for the comparison of sequences is the algorithm of Myers and Miller, (1988) *Comput Appl Biosci*, 4:11-7. Such an algorithm is incorporated into the ALIGN program (version 2.0) which is part of the GCG sequence alignment software package. When utilizing the ALIGN program for comparing amino acid sequences, a PAM120 weight residue table, a gap length penalty of 12, and a gap penalty of 4 can be used. Yet another useful algorithm for identifying regions of local sequence similarity and alignment is the FASTA algorithm as described in Pearson and Lipman (1988) *Proc. Natl. Acad. Sci. USA* 85:2444-2448. When using the FASTA algorithm for comparing nucleotide or amino acid sequences, a PAM120 weight residue table can, for example, be used with a *k*-tuple value of 2.

25       The percent identity between two sequences can be determined using techniques similar to those described above, with or without allowing gaps. In calculating percent identity, only exact matches are counted.

30       The invention also provides chimeric or fusion proteins corresponding to a marker of the invention. As used herein, a "chimeric protein" or "fusion protein" comprises all or part (preferably a biologically active part) of a polypeptide corresponding to a marker of the invention operably linked to a heterologous polypeptide (*i.e.*, a polypeptide other than the polypeptide corresponding to the marker).

Within the fusion protein, the term "operably linked" is intended to indicate that the polypeptide of the invention and the heterologous polypeptide are fused in-frame to each other. The heterologous polypeptide can be fused to the amino-terminus or the carboxyl-terminus of the polypeptide of the invention.

5           One useful fusion protein is a GST fusion protein in which a polypeptide corresponding to a marker of the invention is fused to the carboxyl terminus of GST sequences. Such fusion proteins can facilitate the purification of a recombinant polypeptide of the invention.

10           In another embodiment, the fusion protein contains a heterologous signal sequence at its amino terminus. For example, the native signal sequence of a polypeptide corresponding to a marker of the invention can be removed and replaced with a signal sequence from another protein. For example, the gp67 secretory sequence of the baculovirus envelope protein can be used as a heterologous signal sequence (Ausubel *et al.*, ed., *Current Protocols in Molecular Biology*, John Wiley & Sons, NY, 15 1992). Other examples of eukaryotic heterologous signal sequences include the secretory sequences of melittin and human placental alkaline phosphatase (Stratagene; La Jolla, California). In yet another example, useful prokaryotic heterologous signal sequences include the phoA secretory signal (Sambrook *et al.*, *supra*) and the protein A secretory signal (Pharmacia Biotech; Piscataway, New Jersey).

20           In yet another embodiment, the fusion protein is an immunoglobulin fusion protein in which all or part of a polypeptide corresponding to a marker of the invention is fused to sequences derived from a member of the immunoglobulin protein family. The immunoglobulin fusion proteins of the invention can be incorporated into pharmaceutical compositions and administered to a subject to inhibit an interaction 25 between a ligand (soluble or membrane-bound) and a protein on the surface of a cell (receptor), to thereby suppress signal transduction *in vivo*. The immunoglobulin fusion protein can be used to affect the bioavailability of a cognate ligand of a polypeptide of the invention. Inhibition of ligand/receptor interaction can be useful therapeutically, both for treating proliferative and differentiative disorders and for modulating (*e.g.* 30 promoting or inhibiting) cell survival. Moreover, the immunoglobulin fusion proteins of the invention can be used as immunogens to produce antibodies directed against a

- 50 -

polypeptide of the invention in a subject, to purify ligands and in screening assays to identify molecules which inhibit the interaction of receptors with ligands.

Chimeric and fusion proteins of the invention can be produced by standard recombinant DNA techniques. In another embodiment, the fusion gene can be  
5 synthesized by conventional techniques including automated DNA synthesizers. Alternatively, PCR amplification of gene fragments can be carried out using anchor primers which give rise to complementary overhangs between two consecutive gene fragments which can subsequently be annealed and re-amplified to generate a chimeric gene sequence (see, *e.g.*, Ausubel *et al.*, *supra*). Moreover, many expression vectors are  
10 commercially available that already encode a fusion moiety (*e.g.*, a GST polypeptide). A nucleic acid encoding a polypeptide of the invention can be cloned into such an expression vector such that the fusion moiety is linked in-frame to the polypeptide of the invention.

A signal sequence can be used to facilitate secretion and isolation of the secreted  
15 protein or other proteins of interest. Signal sequences are typically characterized by a core of hydrophobic amino acids which are generally cleaved from the mature protein during secretion in one or more cleavage events. Such signal peptides contain processing sites that allow cleavage of the signal sequence from the mature proteins as they pass through the secretory pathway. Thus, the invention pertains to the described  
20 polypeptides having a signal sequence, as well as to polypeptides from which the signal sequence has been proteolytically cleaved (*i.e.*, the cleavage products). In one embodiment, a nucleic acid sequence encoding a signal sequence can be operably linked in an expression vector to a protein of interest, such as a protein which is ordinarily not secreted or is otherwise difficult to isolate. The signal sequence directs secretion of the  
25 protein, such as from a eukaryotic host into which the expression vector is transformed, and the signal sequence is subsequently or concurrently cleaved. The protein can then be readily purified from the extracellular medium by art recognized methods.

Alternatively, the signal sequence can be linked to the protein of interest using a sequence which facilitates purification, such as with a GST domain.

30 The present invention also pertains to variants of the polypeptides corresponding to individual markers of the invention. Such variants have an altered amino acid sequence which can function as either agonists (mimetics) or as antagonists. Variants

can be generated by mutagenesis, *e.g.*, discrete point mutation or truncation. An agonist can retain substantially the same, or a subset, of the biological activities of the naturally occurring form of the protein. An antagonist of a protein can inhibit one or more of the activities of the naturally occurring form of the protein by, for example, competitively  
5 binding to a downstream or upstream member of a cellular signaling cascade which includes the protein of interest. Thus, specific biological effects can be elicited by treatment with a variant of limited function. Treatment of a subject with a variant having a subset of the biological activities of the naturally occurring form of the protein can have fewer side effects in a subject relative to treatment with the naturally occurring  
10 form of the protein.

Variants of a protein of the invention which function as either agonists (mimetics) or as antagonists can be identified by screening combinatorial libraries of mutants, *e.g.*, truncation mutants, of the protein of the invention for agonist or antagonist activity. In one embodiment, a variegated library of variants is generated by  
15 combinatorial mutagenesis at the nucleic acid level and is encoded by a variegated gene library. A variegated library of variants can be produced by, for example, enzymatically ligating a mixture of synthetic oligonucleotides into gene sequences such that a degenerate set of potential protein sequences is expressible as individual polypeptides, or alternatively, as a set of larger fusion proteins (*e.g.*, for phage display). There are a  
20 variety of methods which can be used to produce libraries of potential variants of the polypeptides of the invention from a degenerate oligonucleotide sequence. Methods for synthesizing degenerate oligonucleotides are known in the art (see, *e.g.*, Narang, 1983, *Tetrahedron* 39:3; Itakura *et al.*, 1984, *Annu. Rev. Biochem.* 53:323; Itakura *et al.*, 1984, *Science* 198:1056; Ike *et al.*, 1983 *Nucleic Acid Res.* 11:477).

25 In addition, libraries of fragments of the coding sequence of a polypeptide corresponding to a marker of the invention can be used to generate a variegated population of polypeptides for screening and subsequent selection of variants. For example, a library of coding sequence fragments can be generated by treating a double stranded PCR fragment of the coding sequence of interest with a nuclease under  
30 conditions wherein nicking occurs only about once per molecule, denaturing the double stranded DNA, renaturing the DNA to form double stranded DNA which can include sense/antisense pairs from different nicked products, removing single stranded portions



from reformed duplexes by treatment with S1 nuclease, and ligating the resulting fragment library into an expression vector. By this method, an expression library can be derived which encodes amino terminal and internal fragments of various sizes of the protein of interest.

5           Several techniques are known in the art for screening gene products of combinatorial libraries made by point mutations or truncation, and for screening cDNA libraries for gene products having a selected property. The most widely used techniques, which are amenable to high through-put analysis, for screening large gene libraries typically include cloning the gene library into replicable expression vectors,  
10 transforming appropriate cells with the resulting library of vectors, and expressing the combinatorial genes under conditions in which detection of a desired activity facilitates isolation of the vector encoding the gene whose product was detected. Recursive ensemble mutagenesis (REM), a technique which enhances the frequency of functional mutants in the libraries, can be used in combination with the screening assays to identify  
15 variants of a protein of the invention (Arkin and Yourvan, 1992, *Proc. Natl. Acad. Sci. USA* 89:7811-7815; Delgrave *et al.*, 1993, *Protein Engineering* 6(3):327- 331).

          An isolated polypeptide corresponding to a marker of the invention, or a fragment thereof, can be used as an immunogen to generate antibodies using standard techniques for polyclonal and monoclonal antibody preparation. The full-length  
20 polypeptide or protein can be used or, alternatively, the invention provides antigenic peptide fragments for use as immunogens. The antigenic peptide of a protein of the invention comprises at least 8 (preferably 10, 15, 20, or 30 or more) amino acid residues of the amino acid sequence of one of the polypeptides of the invention, and encompasses an epitope of the protein such that an antibody raised against the peptide forms a specific  
25 immune complex with a marker of the invention to which the protein corresponds. Preferred epitopes encompassed by the antigenic peptide are regions that are located on the surface of the protein, *e.g.*, hydrophilic regions. Hydrophobicity sequence analysis, hydrophilicity sequence analysis, or similar analyses can be used to identify hydrophilic regions.

30           An immunogen typically is used to prepare antibodies by immunizing a suitable (*i.e.* immunocompetent) subject such as a rabbit, goat, mouse, or other mammal or vertebrate. An appropriate immunogenic preparation can contain, for example,

recombinantly-expressed or chemically-synthesized polypeptide. The preparation can further include an adjuvant, such as Freund's complete or incomplete adjuvant, or a similar immunostimulatory agent.

Accordingly, another aspect of the invention pertains to antibodies directed  
5 against a polypeptide of the invention. The terms "antibody" and "antibody substance" as used interchangeably herein refer to immunoglobulin molecules and immunologically active portions of immunoglobulin molecules, *i.e.*, molecules that contain an antigen binding site which specifically binds an antigen, such as a polypeptide of the invention, *e.g.*, an epitope of a polypeptide of the invention. A molecule which specifically binds  
10 to a given polypeptide of the invention is a molecule which binds the polypeptide, but does not substantially bind other molecules in a sample, *e.g.*, a biological sample, which naturally contains the polypeptide. Examples of immunologically active portions of immunoglobulin molecules include F(ab) and F(ab')<sub>2</sub> fragments which can be generated by treating the antibody with an enzyme such as pepsin. The invention provides  
15 polyclonal and monoclonal antibodies. The term "monoclonal antibody" or "monoclonal antibody composition", as used herein, refers to a population of antibody molecules that contain only one species of an antigen binding site capable of immunoreacting with a particular epitope.

Polyclonal antibodies can be prepared as described above by immunizing a  
20 suitable subject with a polypeptide of the invention as an immunogen. Preferred polyclonal antibody compositions are ones that have been selected for antibodies directed against a polypeptide or polypeptides of the invention. Particularly preferred polyclonal antibody preparations are ones that contain only antibodies directed against a polypeptide or polypeptides of the invention. Particularly preferred immunogen  
25 compositions are those that contain no other human proteins such as, for example, immunogen compositions made using a non-human host cell for recombinant expression of a polypeptide of the invention. In such a manner, the only human epitope or epitopes recognized by the resulting antibody compositions raised against this immunogen will be present as part of a polypeptide or polypeptides of the invention.

30 The antibody titer in the immunized subject can be monitored over time by standard techniques, such as with an enzyme linked immunosorbent assay (ELISA) using immobilized polypeptide. If desired, the antibody molecules can be harvested or

isolated from the subject (*e.g.*, from the blood or serum of the subject) and further purified by well-known techniques, such as protein A chromatography to obtain the IgG fraction. Alternatively, antibodies specific for a protein or polypeptide of the invention can be selected or (*e.g.*, partially purified) or purified by, *e.g.*, affinity chromatography.

5 For example, a recombinantly expressed and purified (or partially purified) protein of the invention is produced as described herein, and covalently or non-covalently coupled to a solid support such as, for example, a chromatography column. The column can then be used to affinity purify antibodies specific for the proteins of the invention from a sample containing antibodies directed against a large number of different epitopes,

10 thereby generating a substantially purified antibody composition, *i.e.*, one that is substantially free of contaminating antibodies. By a substantially purified antibody composition is meant, in this context, that the antibody sample contains at most only 30% (by dry weight) of contaminating antibodies directed against epitopes other than those of the desired protein or polypeptide of the invention, and preferably at most 20%,

15 yet more preferably at most 10%, and most preferably at most 5% (by dry weight) of the sample is contaminating antibodies. A purified antibody composition means that at least 99% of the antibodies in the composition are directed against the desired protein or polypeptide of the invention.

At an appropriate time after immunization, *e.g.*, when the specific antibody titers

20 are highest, antibody-producing cells can be obtained from the subject and used to prepare monoclonal antibodies by standard techniques, such as the hybridoma technique originally described by Kohler and Milstein (1975) *Nature* 256:495-497, the human B cell hybridoma technique (see Kozbor *et al.*, 1983, *Immunol. Today* 4:72), the EBV-hybridoma technique (see Cole *et al.*, pp. 77-96 In *Monoclonal Antibodies and Cancer*

25 *Therapy*, Alan R. Liss, Inc., 1985) or trioma techniques. The technology for producing hybridomas is well known (see generally *Current Protocols in Immunology*, Coligan *et al.* ed., John Wiley & Sons, New York, 1994). Hybridoma cells producing a monoclonal antibody of the invention are detected by screening the hybridoma culture supernatants for antibodies that bind the polypeptide of interest, *e.g.*, using a standard

30 ELISA assay.

Alternative to preparing monoclonal antibody-secreting hybridomas, a monoclonal antibody directed against a polypeptide of the invention can be identified and isolated by screening a recombinant combinatorial immunoglobulin library (*e.g.*, an antibody phage display library) with the polypeptide of interest. Kits for generating and screening phage display libraries are commercially available (*e.g.*, the Pharmacia *Recombinant Phage Antibody System*, Catalog No. 27-9400-01; and the Stratagene *SurfZAP Phage Display Kit*, Catalog No. 240612). Additionally, examples of methods and reagents particularly amenable for use in generating and screening antibody display library can be found in, for example, U.S. Patent No. 5,223,409; PCT Publication No. WO 92/18619; PCT Publication No. WO 91/17271; PCT Publication No. WO 92/20791; PCT Publication No. WO 92/15679; PCT Publication No. WO 93/01288; PCT Publication No. WO 92/01047; PCT Publication No. WO 92/09690; PCT Publication No. WO 90/02809; Fuchs *et al.* (1991) *Bio/Technology* 9:1370-1372; Hay *et al.* (1992) *Hum. Antibod. Hybridomas* 3:81-85; Huse *et al.* (1989) *Science* 246:1275-1281; Griffiths *et al.* (1993) *EMBO J.* 12:725-734.

Additionally, recombinant antibodies, such as chimeric and humanized monoclonal antibodies, comprising both human and non-human portions, which can be made using standard recombinant DNA techniques, are within the scope of the invention. A chimeric antibody is a molecule in which different portions are derived from different animal species, such as those having a variable region derived from a murine mAb and a human immunoglobulin constant region. (See, *e.g.*, Cabilly *et al.*, U.S. Patent No. 4,816,567; and Boss *et al.*, U.S. Patent No. 4,816,397, which are incorporated herein by reference in their entirety.) Humanized antibodies are antibody molecules from non-human species having one or more complementarily determining regions (CDRs) from the non-human species and a framework region from a human immunoglobulin molecule. (See, *e.g.*, Queen, U.S. Patent No. 5,585,089, which is incorporated herein by reference in its entirety.) Such chimeric and humanized monoclonal antibodies can be produced by recombinant DNA techniques known in the art, for example using methods described in PCT Publication No. WO 87/02671; European Patent Application 184,187; European Patent Application 171,496; European Patent Application 173,494; PCT Publication No. WO 86/01533; U.S. Patent No. 4,816,567; European Patent Application 125,023; Better *et al.* (1988) *Science* 240:1041-

1043; Liu *et al.* (1987) *Proc. Natl. Acad. Sci. USA* 84:3439-3443; Liu *et al.* (1987) *J. Immunol.* 139:3521-3526; Sun *et al.* (1987) *Proc. Natl. Acad. Sci. USA* 84:214-218; Nishimura *et al.* (1987) *Cancer Res.* 47:999-1005; Wood *et al.* (1985) *Nature* 314:446-449; and Shaw *et al.* (1988) *J. Natl. Cancer Inst.* 80:1553-1559; Morrison (1985) *Science* 229:1202-1207; Oi *et al.* (1986) *Bio/Techniques* 4:214; U.S. Patent 5,225,539; Jones *et al.* (1986) *Nature* 321:552-525; Verhoeyan *et al.* (1988) *Science* 239:1534; and Beidler *et al.* (1988) *J. Immunol.* 141:4053-4060.

Antibodies of the invention may be used as therapeutic agents in treating cancers. In a preferred embodiment, completely human antibodies of the invention are used for therapeutic treatment of human cancer patients, particularly those having an ovarian cancer. Such antibodies can be produced, for example, using transgenic mice which are incapable of expressing endogenous immunoglobulin heavy and light chains genes, but which can express human heavy and light chain genes. The transgenic mice are immunized in the normal fashion with a selected antigen, *e.g.*, all or a portion of a polypeptide corresponding to a marker of the invention. Monoclonal antibodies directed against the antigen can be obtained using conventional hybridoma technology. The human immunoglobulin transgenes harbored by the transgenic mice rearrange during B cell differentiation, and subsequently undergo class switching and somatic mutation. Thus, using such a technique, it is possible to produce therapeutically useful IgG, IgA and IgE antibodies. For an overview of this technology for producing human antibodies, see Lonberg and Huszar (1995) *Int. Rev. Immunol.* 13:65-93). For a detailed discussion of this technology for producing human antibodies and human monoclonal antibodies and protocols for producing such antibodies, see, *e.g.*, U.S. Patent 5,625,126; U.S. Patent 5,633,425; U.S. Patent 5,569,825; U.S. Patent 5,661,016; and U.S. Patent 5,545,806. In addition, companies such as Abgenix, Inc. (Freemont, CA), can be engaged to provide human antibodies directed against a selected antigen using technology similar to that described above.

Completely human antibodies which recognize a selected epitope can be generated using a technique referred to as "guided selection." In this approach a selected non-human monoclonal antibody, *e.g.*, a murine antibody, is used to guide the selection of a completely human antibody recognizing the same epitope (Jespers *et al.*, 1994, *Bio/technology* 12:899-903).

- An antibody directed against a polypeptide corresponding to a marker of the invention (*e.g.*, a monoclonal antibody) can be used to isolate the polypeptide by standard techniques, such as affinity chromatography or immunoprecipitation. Moreover, such an antibody can be used to detect the marker (*e.g.*, in a cellular lysate or
- 5 cell supernatant) in order to evaluate the level and pattern of expression of the marker. The antibodies can also be used diagnostically to monitor protein levels in tissues or body fluids (*e.g.* in an ovary-associated body fluid) as part of a clinical testing procedure, *e.g.*, to, for example, determine the efficacy of a given treatment regimen. Detection can be facilitated by coupling the antibody to a detectable substance.
- 10 Examples of detectable substances include various enzymes, prosthetic groups, fluorescent materials, luminescent materials, bioluminescent materials, and radioactive materials. Examples of suitable enzymes include horseradish peroxidase, alkaline phosphatase,  $\beta$ -galactosidase, or acetylcholinesterase; examples of suitable prosthetic group complexes include streptavidin/biotin and avidin/biotin; examples of suitable
- 15 fluorescent materials include umbelliferone, fluorescein, fluorescein isothiocyanate, rhodamine, dichlorotriazinylamine fluorescein, dansyl chloride or phycoerythrin; an example of a luminescent material includes luminol; examples of bioluminescent materials include luciferase, luciferin, and aequorin, and examples of suitable radioactive material include  $^{125}\text{I}$ ,  $^{131}\text{I}$ ,  $^{35}\text{S}$  or  $^3\text{H}$ .
- 20 Further, an antibody (or fragment thereof) can be conjugated to a therapeutic moiety such as a cytotoxin, a therapeutic agent or a radioactive metal ion. A cytotoxin or cytotoxic agent includes any agent that is detrimental to cells. Examples include taxol, cytochalasin B, gramicidin D, ethidium bromide, emetine, mitomycin, etoposide, tenoposide, vincristine, vinblastine, colchicin, doxorubicin, daunorubicin, dihydroxy
- 25 anthracin dione, mitoxantrone, mithramycin, actinomycin D, 1-dehydrotestosterone, glucocorticoids, procaine, tetracaine, lidocaine, propranolol, and puromycin and analogs or homologs thereof. Therapeutic agents include, but are not limited to, antimetabolites (*e.g.*, methotrexate, 6-mercaptopurine, 6-thioguanine, cytarabine, 5-fluorouracil decarbazine), alkylating agents (*e.g.*, mechlorethamine, thioepa chlorambucil,
- 30 melphalan, carmustine (BSNU) and lomustine (CCNU), cyclophosphamide, busulfan, dibromomannitol, streptozotocin, mitomycin C, and cis-dichlorodiamine platinum (II) (DDP) cisplatin), anthracyclines (*e.g.*, daunorubicin (formerly daunomycin) and

doxorubicin), antibiotics (*e.g.*, dactinomycin (formerly actinomycin), bleomycin, mithramycin, and anthramycin (AMC)), and anti-mitotic agents (*e.g.*, vincristine and vinblastine).

The conjugates of the invention can be used for modifying a given biological response, the drug moiety is not to be construed as limited to classical chemical therapeutic agents. For example, the drug moiety may be a protein or polypeptide possessing a desired biological activity. Such proteins may include, for example, a toxin such as abrin, ricin A, pseudomonas exotoxin, or diphtheria toxin; a protein such as tumor necrosis factor, .alpha.-interferon, .beta.-interferon, nerve growth factor, platelet derived growth factor, tissue plasminogen activator; or, biological response modifiers such as, for example, lymphokines, interleukin-1 ("IL-1"), interleukin-2 ("IL-2"), interleukin-6 ("IL-6"), granulocyte macrophase colony stimulating factor ("GM-CSF"), granulocyte colony stimulating factor ("G-CSF"), or other growth factors.

Techniques for conjugating such therapeutic moiety to antibodies are well known, see, *e.g.*, Arnon et al., "Monoclonal Antibodies For Immunotargeting Of Drugs In Cancer Therapy", in *Monoclonal Antibodies And Cancer Therapy*, Reisfeld et al. (eds.), pp. 243-56 (Alan R. Liss, Inc. 1985); Hellstrom et al., "Antibodies For Drug Delivery", in *Controlled Drug Delivery* (2nd Ed.), Robinson et al. (eds.), pp. 623-53 (Marcel Dekker, Inc. 1987); Thorpe, "Antibody Carriers Of Cytotoxic Agents In Cancer Therapy: A Review", in *Monoclonal Antibodies '84: Biological And Clinical Applications*, Pinchera et al. (eds.), pp. 475-506 (1985); "Analysis, Results, And Future Prospective Of The Therapeutic Use Of Radiolabeled Antibody In Cancer Therapy", in *Monoclonal Antibodies For Cancer Detection And Therapy*, Baldwin et al. (eds.), pp. 303-16 (Academic Press 1985), and Thorpe et al., "The Preparation And Cytotoxic Properties Of Antibody-Toxin Conjugates", *Immunol. Rev.*, 62:119-58 (1982).

Alternatively, an antibody can be conjugated to a second antibody to form an antibody heteroconjugate as described by Segal in U.S. Patent No. 4,676,980.

Accordingly, in one aspect, the invention provides substantially purified antibodies or fragments thereof, and non-human antibodies or fragments thereof, which antibodies or fragments specifically bind to a polypeptide comprising an amino acid sequence selected from the group consisting of the amino acid sequences of the present invention, an amino acid sequence encoded by the cDNA of the present invention, a

fragment of at least 15 amino acid residues of an amino acid sequence of the present invention, an amino acid sequence which is at least 95% identical to the amino acid sequence of the present invention (wherein the percent identity is determined using the ALIGN program of the GCG software package with a PAM120 weight residue table, a gap length penalty of 12, and a gap penalty of 4) and an amino acid sequence which is encoded by a nucleic acid molecule which hybridizes to a nucleic acid molecule consisting of the nucleic acid molecules of the present invention, or a complement thereof, under conditions of hybridization of 6X SSC at 45°C and washing in 0.2 X SSC, 0.1% SDS at 65°C. In various embodiments, the substantially purified antibodies of the invention, or fragments thereof, can be human, non-human, chimeric and/or humanized antibodies.

In another aspect, the invention provides non-human antibodies or fragments thereof, which antibodies or fragments specifically bind to a polypeptide comprising an amino acid sequence selected from the group consisting of: the amino acid sequence of the present invention, an amino acid sequence encoded by the cDNA of the present invention, a fragment of at least 15 amino acid residues of the amino acid sequence of the present invention, an amino acid sequence which is at least 95% identical to the amino acid sequence of the present invention (wherein the percent identity is determined using the ALIGN program of the GCG software package with a PAM120 weight residue table, a gap length penalty of 12, and a gap penalty of 4) and an amino acid sequence which is encoded by a nucleic acid molecule which hybridizes to a nucleic acid molecule consisting of the nucleic acid molecules of the present invention, or a complement thereof, under conditions of hybridization of 6X SSC at 45°C and washing in 0.2 X SSC, 0.1% SDS at 65°C. Such non-human antibodies can be goat, mouse, sheep, horse, chicken, rabbit, or rat antibodies. Alternatively, the non-human antibodies of the invention can be chimeric and/or humanized antibodies. In addition, the non-human antibodies of the invention can be polyclonal antibodies or monoclonal antibodies.

In still a further aspect, the invention provides monoclonal antibodies or fragments thereof, which antibodies or fragments specifically bind to a polypeptide comprising an amino acid sequence selected from the group consisting of the amino acid sequences of the present invention, an amino acid sequence encoded by the cDNA of the



present invention, a fragment of at least 15 amino acid residues of an amino acid sequence of the present invention, an amino acid sequence which is at least 95% identical to an amino acid sequence of the present invention (wherein the percent identity is determined using the ALIGN program of the GCG software package with a  
5 PAM120 weight residue table, a gap length penalty of 12, and a gap penalty of 4) and an amino acid sequence which is encoded by a nucleic acid molecule which hybridizes to a nucleic acid molecule consisting of the nucleic acid molecules of the present invention, or a complement thereof, under conditions of hybridization of 6X SSC at 45°C and washing in 0.2 X SSC, 0.1% SDS at 65°C. The monoclonal antibodies can be human,  
10 humanized, chimeric and/or non-human antibodies.

The substantially purified antibodies or fragments thereof may specifically bind to a signal peptide, a secreted sequence, an extracellular domain, a transmembrane or a cytoplasmic domain or cytoplasmic membrane of a polypeptide of the invention. In a particularly preferred embodiment, the substantially purified antibodies or fragments  
15 thereof, the non-human antibodies or fragments thereof, and/or the monoclonal antibodies or fragments thereof, of the invention specifically bind to a secreted sequence or an extracellular domain of the amino acid sequences of the present invention.

Any of the antibodies of the invention can be conjugated to a therapeutic moiety or to a detectable substance. Non-limiting examples of detectable substances that can be  
20 conjugated to the antibodies of the invention are an enzyme, a prosthetic group, a fluorescent material, a luminescent material, a bioluminescent material, and a radioactive material.

The invention also provides a kit containing an antibody of the invention conjugated to a detectable substance, and instructions for use. Still another aspect of the  
25 invention is a pharmaceutical composition comprising an antibody of the invention and a pharmaceutically acceptable carrier. In preferred embodiments, the pharmaceutical composition contains an antibody of the invention, a therapeutic moiety, and a pharmaceutically acceptable carrier.

Still another aspect of the invention is a method of making an antibody that  
30 specifically recognizes a polypeptide of the present invention, the method comprising immunizing a mammal with a polypeptide. The polypeptide used as an immunogen comprises an amino acid sequence selected from the group consisting of the amino acid

sequence of the present invention, an amino acid sequence encoded by the cDNA of the nucleic acid molecules of the present invention, a fragment of at least 15 amino acid residues of the amino acid sequence of the present invention, an amino acid sequence which is at least 95% identical to the amino acid sequence of the present invention  
5 (wherein the percent identity is determined using the ALIGN program of the GCG software package with a PAM120 weight residue table, a gap length penalty of 12, and a gap penalty of 4) and an amino acid sequence which is encoded by a nucleic acid molecule which hybridizes to a nucleic acid molecule consisting of the nucleic acid molecules of the present invention, or a complement thereof, under conditions of  
10 hybridization of 6X SSC at 45°C and washing in 0.2 X SSC, 0.1% SDS at 65°C.

After immunization, a sample is collected from the mammal that contains an antibody that specifically recognizes the polypeptide. Preferably, the polypeptide is recombinantly produced using a non-human host cell. Optionally, the antibodies can be further purified from the sample using techniques well known to those of skill in the art.  
15 The method can further comprise producing a monoclonal antibody- producing cell from the cells of the mammal. Optionally, antibodies are collected from the antibody-producing cell.

### III. Recombinant Expression Vectors and Host Cells

20 Another aspect of the invention pertains to vectors, preferably expression vectors, containing a nucleic acid encoding a polypeptide corresponding to a marker of the invention (or a portion of such a polypeptide). As used herein, the term "vector" refers to a nucleic acid molecule capable of transporting another nucleic acid to which it has been linked. One type of vector is a "plasmid", which refers to a circular double  
25 stranded DNA loop into which additional DNA segments can be ligated. Another type of vector is a viral vector, wherein additional DNA segments can be ligated into the viral genome. Certain vectors are capable of autonomous replication in a host cell into which they are introduced (*e.g.*, bacterial vectors having a bacterial origin of replication and episomal mammalian vectors). Other vectors (*e.g.*, non-episomal mammalian vectors)  
30 are integrated into the genome of a host cell upon introduction into the host cell, and thereby are replicated along with the host genome. Moreover, certain vectors, namely expression vectors, are capable of directing the expression of genes to which they are

operably linked. In general, expression vectors of utility in recombinant DNA techniques are often in the form of plasmids (vectors). However, the invention is intended to include such other forms of expression vectors, such as viral vectors (*e.g.*, replication defective retroviruses, adenoviruses and adeno-associated viruses), which  
5 serve equivalent functions.

The recombinant expression vectors of the invention comprise a nucleic acid of the invention in a form suitable for expression of the nucleic acid in a host cell. This means that the recombinant expression vectors include one or more regulatory sequences, selected on the basis of the host cells to be used for expression, which is  
10 operably linked to the nucleic acid sequence to be expressed. Within a recombinant expression vector, "operably linked" is intended to mean that the nucleotide sequence of interest is linked to the regulatory sequence(s) in a manner which allows for expression of the nucleotide sequence (*e.g.*, in an *in vitro* transcription/translation system or in a host cell when the vector is introduced into the host cell). The term "regulatory  
15 sequence" is intended to include promoters, enhancers and other expression control elements (*e.g.*, polyadenylation signals). Such regulatory sequences are described, for example, in Goeddel, *Methods in Enzymology: Gene Expression Technology* vol.185, Academic Press, San Diego, CA (1991). Regulatory sequences include those which direct constitutive expression of a nucleotide sequence in many types of host cell and  
20 those which direct expression of the nucleotide sequence only in certain host cells (*e.g.*, tissue-specific regulatory sequences). It will be appreciated by those skilled in the art that the design of the expression vector can depend on such factors as the choice of the host cell to be transformed, the level of expression of protein desired, and the like. The expression vectors of the invention can be introduced into host cells to thereby produce  
25 proteins or peptides, including fusion proteins or peptides, encoded by nucleic acids as described herein.

The recombinant expression vectors of the invention can be designed for expression of a polypeptide corresponding to a marker of the invention in prokaryotic (*e.g.*, *E. coli*) or eukaryotic cells (*e.g.*, insect cells {using baculovirus expression  
30 vectors}, yeast cells or mammalian cells). Suitable host cells are discussed further in Goeddel, *supra*. Alternatively, the recombinant expression vector can be transcribed

and translated *in vitro*, for example using T7 promoter regulatory sequences and T7 polymerase.

Expression of proteins in prokaryotes is most often carried out in *E. coli* with vectors containing constitutive or inducible promoters directing the expression of either fusion or non-fusion proteins. Fusion vectors add a number of amino acids to a protein encoded therein, usually to the amino terminus of the recombinant protein. Such fusion vectors typically serve three purposes: 1) to increase expression of recombinant protein; 2) to increase the solubility of the recombinant protein; and 3) to aid in the purification of the recombinant protein by acting as a ligand in affinity purification. Often, in fusion expression vectors, a proteolytic cleavage site is introduced at the junction of the fusion moiety and the recombinant protein to enable separation of the recombinant protein from the fusion moiety subsequent to purification of the fusion protein. Such enzymes, and their cognate recognition sequences, include Factor Xa, thrombin and enterokinase. Typical fusion expression vectors include pGEX (Pharmacia Biotech Inc; Smith and Johnson, 1988, *Gene* 67:31-40), pMAL (New England Biolabs, Beverly, MA) and pRIT5 (Pharmacia, Piscataway, NJ) which fuse glutathione S-transferase (GST), maltose E binding protein, or protein A, respectively, to the target recombinant protein.

Examples of suitable inducible non-fusion *E. coli* expression vectors include pTrc (Amann *et al.*, 1988, *Gene* 69:301-315) and pET 11d (Studier *et al.*, p. 60-89, In *Gene Expression Technology: Methods in Enzymology* vol.185, Academic Press, San Diego, CA, 1991). Target gene expression from the pTrc vector relies on host RNA polymerase transcription from a hybrid trp-lac fusion promoter. Target gene expression from the pET 11d vector relies on transcription from a T7 gn10-lac fusion promoter mediated by a co-expressed viral RNA polymerase (T7 gn1). This viral polymerase is supplied by host strains BL21(DE3) or HMS174(DE3) from a resident prophage harboring a T7 gn1 gene under the transcriptional control of the lacUV 5 promoter.

One strategy to maximize recombinant protein expression in *E. coli* is to express the protein in a host bacteria with an impaired capacity to proteolytically cleave the recombinant protein (Gottesman, p. 119-128, In *Gene Expression Technology: Methods in Enzymology* vol. 185, Academic Press, San Diego, CA, 1990). Another strategy is to alter the nucleic acid sequence of the nucleic acid to be inserted into an expression vector so that the individual codons for each amino acid are those preferentially utilized

in *E. coli* (Wada *et al.*, 1992, *Nucleic Acids Res.* 20:2111-2118). Such alteration of nucleic acid sequences of the invention can be carried out by standard DNA synthesis techniques.

In another embodiment, the expression vector is a yeast expression vector.

- 5 Examples of vectors for expression in yeast *S. cerevisiae* include pYepSec1 (Baldari *et al.*, 1987, *EMBO J.* 6:229-234), pMFa (Kurjan and Herskowitz, 1982, *Cell* 30:933-943), pJRY88 (Schultz *et al.*, 1987, *Gene* 54:113-123), pYES2 (Invitrogen Corporation, San Diego, CA), and pPicZ (Invitrogen Corp, San Diego, CA).

Alternatively, the expression vector is a baculovirus expression vector.

- 10 Baculovirus vectors available for expression of proteins in cultured insect cells (*e.g.*, Sf 9 cells) include the pAc series (Smith *et al.*, 1983, *Mol. Cell Biol.* 3:2156-2165) and the pVL series (Lucklow and Summers, 1989, *Virology* 170:31-39).

- In yet another embodiment, a nucleic acid of the invention is expressed in mammalian cells using a mammalian expression vector. Examples of mammalian  
15 expression vectors include pCDM8 (Seed, 1987, *Nature* 329:840) and pMT2PC (Kaufman *et al.*, 1987, *EMBO J.* 6:187-195). When used in mammalian cells, the expression vector's control functions are often provided by viral regulatory elements. For example, commonly used promoters are derived from polyoma, Adenovirus 2, cytomegalovirus and Simian Virus 40. For other suitable expression systems for both  
20 prokaryotic and eukaryotic cells see chapters 16 and 17 of Sambrook *et al.*, *supra*.

- In another embodiment, the recombinant mammalian expression vector is capable of directing expression of the nucleic acid preferentially in a particular cell type (*e.g.*, tissue-specific regulatory elements are used to express the nucleic acid). Tissue-specific regulatory elements are known in the art. Non-limiting examples of suitable  
25 tissue-specific promoters include the albumin promoter (liver-specific; Pinkert *et al.*, 1987, *Genes Dev.* 1:268-277), lymphoid-specific promoters (Calame and Eaton, 1988, *Adv. Immunol.* 43:235-275), in particular promoters of T cell receptors (Winoto and Baltimore, 1989, *EMBO J.* 8:729-733) and immunoglobulins (Banerji *et al.*, 1983, *Cell* 33:729-740; Queen and Baltimore, 1983, *Cell* 33:741-748), neuron-specific promoters  
30 (*e.g.*, the neurofilament promoter; Byrne and Ruddle, 1989, *Proc. Natl. Acad. Sci. USA* 86:5473-5477), pancreas-specific promoters (Edlund *et al.*, 1985, *Science* 230:912-916), and mammary gland-specific promoters (*e.g.*, milk whey promoter; U.S. Patent No.

4,873,316 and European Application Publication No. 264,166). Developmentally-regulated promoters are also encompassed, for example the murine hox promoters (Kessel and Gruss, 1990, *Science* 249:374-379) and the  $\alpha$ -fetoprotein promoter (Camper and Tilghman, 1989, *Genes Dev.* 3:537-546).

5           The invention further provides a recombinant expression vector comprising a DNA molecule of the invention cloned into the expression vector in an antisense orientation. That is, the DNA molecule is operably linked to a regulatory sequence in a manner which allows for expression (by transcription of the DNA molecule) of an RNA molecule which is antisense to the mRNA encoding a polypeptide of the invention.

10          Regulatory sequences operably linked to a nucleic acid cloned in the antisense orientation can be chosen which direct the continuous expression of the antisense RNA molecule in a variety of cell types, for instance viral promoters and/or enhancers, or regulatory sequences can be chosen which direct constitutive, tissue-specific or cell type specific expression of antisense RNA. The antisense expression vector can be in the  
15          form of a recombinant plasmid, phagemid, or attenuated virus in which antisense nucleic acids are produced under the control of a high efficiency regulatory region, the activity of which can be determined by the cell type into which the vector is introduced. For a discussion of the regulation of gene expression using antisense genes see Weintraub *et al.*, 1986, *Trends in Genetics*, Vol. 1(1).

20          Another aspect of the invention pertains to host cells into which a recombinant expression vector of the invention has been introduced. The terms "host cell" and "recombinant host cell" are used interchangeably herein. It is understood that such terms refer not only to the particular subject cell but to the progeny or potential progeny of such a cell. Because certain modifications may occur in succeeding generations due  
25          to either mutation or environmental influences, such progeny may not, in fact, be identical to the parent cell, but are still included within the scope of the term as used herein.

          A host cell can be any prokaryotic (*e.g.*, *E. coli*) or eukaryotic cell (*e.g.*, insect cells, yeast or mammalian cells).

30          Vector DNA can be introduced into prokaryotic or eukaryotic cells via conventional transformation or transfection techniques. As used herein, the terms "transformation" and "transfection" are intended to refer to a variety of art-recognized

techniques for introducing foreign nucleic acid into a host cell, including calcium phosphate or calcium chloride co-precipitation, DEAE-dextran-mediated transfection, lipofection, or electroporation. Suitable methods for transforming or transfecting host cells can be found in Sambrook, *et al.* (*supra*), and other laboratory manuals.

5           For stable transfection of mammalian cells, it is known that, depending upon the expression vector and transfection technique used, only a small fraction of cells may integrate the foreign DNA into their genome. In order to identify and select these integrants, a gene that encodes a selectable marker (*e.g.*, for resistance to antibiotics) is generally introduced into the host cells along with the gene of interest. Preferred  
10   selectable markers include those which confer resistance to drugs, such as G418, hygromycin and methotrexate. Cells stably transfected with the introduced nucleic acid can be identified by drug selection (*e.g.*, cells that have incorporated the selectable marker gene will survive, while the other cells die).

          A host cell of the invention, such as a prokaryotic or eukaryotic host cell in  
15   culture, can be used to produce a polypeptide corresponding to a marker of the invention. Accordingly, the invention further provides methods for producing a polypeptide corresponding to a marker of the invention using the host cells of the invention. In one embodiment, the method comprises culturing the host cell of invention (into which a recombinant expression vector encoding a polypeptide of the  
20   invention has been introduced) in a suitable medium such that the marker is produced. In another embodiment, the method further comprises isolating the marker polypeptide from the medium or the host cell.

          The host cells of the invention can also be used to produce nonhuman transgenic animals. For example, in one embodiment, a host cell of the invention is a fertilized  
25   oocyte or an embryonic stem cell into which a sequences encoding a polypeptide corresponding to a marker of the invention have been introduced. Such host cells can then be used to create non-human transgenic animals in which exogenous sequences encoding a marker protein of the invention have been introduced into their genome or homologous recombinant animals in which endogenous gene(s) encoding a polypeptide  
30   corresponding to a marker of the invention sequences have been altered. Such animals are useful for studying the function and/or activity of the polypeptide corresponding to the marker and for identifying and/or evaluating modulators of polypeptide activity. As

- 67 -

used herein, a "transgenic animal" is a non-human animal, preferably a mammal, more preferably a rodent such as a rat or mouse, in which one or more of the cells of the animal includes a transgene. Other examples of transgenic animals include non-human primates, sheep, dogs, cows, goats, chickens, amphibians, etc. A transgene is exogenous DNA which is integrated into the genome of a cell from which a transgenic animal develops and which remains in the genome of the mature animal, thereby directing the expression of an encoded gene product in one or more cell types or tissues of the transgenic animal. As used herein, an "homologous recombinant animal" is a non-human animal, preferably a mammal, more preferably a mouse, in which an endogenous gene has been altered by homologous recombination between the endogenous gene and an exogenous DNA molecule introduced into a cell of the animal, *e.g.*, an embryonic cell of the animal, prior to development of the animal.

A transgenic animal of the invention can be created by introducing a nucleic acid encoding a polypeptide corresponding to a marker of the invention into the male pronuclei of a fertilized oocyte, *e.g.*, by microinjection, retroviral infection, and allowing the oocyte to develop in a pseudopregnant female foster animal. Intronic sequences and polyadenylation signals can also be included in the transgene to increase the efficiency of expression of the transgene. A tissue-specific regulatory sequence(s) can be operably linked to the transgene to direct expression of the polypeptide of the invention to particular cells. Methods for generating transgenic animals via embryo manipulation and microinjection, particularly animals such as mice, have become conventional in the art and are described, for example, in U.S. Patent Nos. 4,736,866 and 4,870,009, U.S. Patent No. 4,873,191 and in Hogan, *Manipulating the Mouse Embryo*, Cold Spring Harbor Laboratory Press, Cold Spring Harbor, N.Y., 1986. Similar methods are used for production of other transgenic animals. A transgenic founder animal can be identified based upon the presence of the transgene in its genome and/or expression of mRNA encoding the transgene in tissues or cells of the animals. A transgenic founder animal can then be used to breed additional animals carrying the transgene. Moreover, transgenic animals carrying the transgene can further be bred to other transgenic animals carrying other transgenes.



To create an homologous recombinant animal, a vector is prepared which contains at least a portion of a gene encoding a polypeptide corresponding to a marker of the invention into which a deletion, addition or substitution has been introduced to thereby alter, *e.g.*, functionally disrupt, the gene. In a preferred embodiment, the vector

5 is designed such that, upon homologous recombination, the endogenous gene is functionally disrupted (*i.e.*, no longer encodes a functional protein; also referred to as a "knock out" vector). Alternatively, the vector can be designed such that, upon homologous recombination, the endogenous gene is mutated or otherwise altered but still encodes functional protein (*e.g.*, the upstream regulatory region can be altered to

10 thereby alter the expression of the endogenous protein). In the homologous recombination vector, the altered portion of the gene is flanked at its 5' and 3' ends by additional nucleic acid of the gene to allow for homologous recombination to occur between the exogenous gene carried by the vector and an endogenous gene in an embryonic stem cell. The additional flanking nucleic acid sequences are of sufficient

15 length for successful homologous recombination with the endogenous gene. Typically, several kilobases of flanking DNA (both at the 5' and 3' ends) are included in the vector (see, *e.g.*, Thomas and Capecchi, 1987, *Cell* 51:503 for a description of homologous recombination vectors). The vector is introduced into an embryonic stem cell line (*e.g.*, by electroporation) and cells in which the introduced gene has homologously

20 recombined with the endogenous gene are selected (see, *e.g.*, Li *et al.*, 1992, *Cell* 69:915). The selected cells are then injected into a blastocyst of an animal (*e.g.*, a mouse) to form aggregation chimeras (see, *e.g.*, Bradley, *Teratocarcinomas and Embryonic Stem Cells: A Practical Approach*, Robertson, Ed., IRL, Oxford, 1987, pp. 113-152). A chimeric embryo can then be implanted into a suitable pseudopregnant

25 female foster animal and the embryo brought to term. Progeny harboring the homologously recombined DNA in their germ cells can be used to breed animals in which all cells of the animal contain the homologously recombined DNA by germline transmission of the transgene. Methods for constructing homologous recombination vectors and homologous recombinant animals are described further in Bradley (1991)

30 *Current Opinion in Bio/Technology* 2:823-829 and in PCT Publication NOS. WO 90/11354, WO 91/01140, WO 92/0968, and WO 93/04169.

In another embodiment, transgenic non-human animals can be produced which contain selected systems which allow for regulated expression of the transgene. One example of such a system is the *cre/loxP* recombinase system of bacteriophage P1. For a description of the *cre/loxP* recombinase system, see, e.g., Lakso *et al.* (1992) *Proc. Natl. Acad. Sci. USA* 89:6232-6236. Another example of a recombinase system is the FLP recombinase system of *Saccharomyces cerevisiae* (O'Gorman *et al.*, 1991, *Science* 251:1351-1355). If a *cre/loxP* recombinase system is used to regulate expression of the transgene, animals containing transgenes encoding both the *Cre* recombinase and a selected protein are required. Such animals can be provided through the construction of "double" transgenic animals, e.g., by mating two transgenic animals, one containing a transgene encoding a selected protein and the other containing a transgene encoding a recombinase.

Clones of the non-human transgenic animals described herein can also be produced according to the methods described in Wilmut *et al.* (1997) *Nature* 385:810-813 and PCT Publication NOS. WO 97/07668 and WO 97/07669.

#### IV. Pharmaceutical Compositions

The nucleic acid molecules, polypeptides, and antibodies (also referred to herein as "active compounds") corresponding to a marker of the invention can be incorporated into pharmaceutical compositions suitable for administration. Such compositions typically comprise the nucleic acid molecule, protein, or antibody and a pharmaceutically acceptable carrier. As used herein the language "pharmaceutically acceptable carrier" is intended to include any and all solvents, dispersion media, coatings, antibacterial and antifungal agents, isotonic and absorption delaying agents, and the like, compatible with pharmaceutical administration. The use of such media and agents for pharmaceutically active substances is well known in the art. Except insofar as any conventional media or agent is incompatible with the active compound, use thereof in the compositions is contemplated. Supplementary active compounds can also be incorporated into the compositions.

The invention includes methods for preparing pharmaceutical compositions for modulating the expression or activity of a polypeptide or nucleic acid corresponding to a marker of the invention. Such methods comprise formulating a pharmaceutically

acceptable carrier with an agent which modulates expression or activity of a polypeptide or nucleic acid corresponding to a marker of the invention. Such compositions can further include additional active agents. Thus, the invention further includes methods for preparing a pharmaceutical composition by formulating a pharmaceutically

5 acceptable carrier with an agent which modulates expression or activity of a polypeptide or nucleic acid corresponding to a marker of the invention and one or more additional active compounds.

The invention also provides methods (also referred to herein as "screening assays") for identifying modulators, *i.e.*, candidate or test compounds or agents (*e.g.*,  
10 peptides, peptidomimetics, peptoids, small molecules or other drugs) which (a) bind to the marker, or (b) have a modulatory (*e.g.*, stimulatory or inhibitory) effect on the activity of the marker or, more specifically, (c) have a modulatory effect on the interactions of the marker with one or more of its natural substrates (*e.g.*, peptide, protein, hormone, co-factor, or nucleic acid), or (d) have a modulatory effect on the  
15 expression of the marker. Such assays typically comprise a reaction between the marker and one or more assay components. The other components may be either the test compound itself, or a combination of test compound and a natural binding partner of the marker.

The test compounds of the present invention may be obtained from any available  
20 source, including systematic libraries of natural and/or synthetic compounds. Test compounds may also be obtained by any of the numerous approaches in combinatorial library methods known in the art, including: biological libraries; peptoid libraries (libraries of molecules having the functionalities of peptides, but with a novel, non-peptide backbone which are resistant to enzymatic degradation but which nevertheless  
25 remain bioactive; see, *e.g.*, Zuckermann *et al.*, 1994, *J. Med. Chem.* 37:2678-85); spatially addressable parallel solid phase or solution phase libraries; synthetic library methods requiring deconvolution; the 'one-bead one-compound' library method; and synthetic library methods using affinity chromatography selection. The biological library and peptoid library approaches are limited to peptide libraries, while the other  
30 four approaches are applicable to peptide, non-peptide oligomer or small molecule libraries of compounds (Lam, 1997, *Anticancer Drug Des.* 12:145).

Examples of methods for the synthesis of molecular libraries can be found in the art, for example in: DeWitt *et al.* (1993) *Proc. Natl. Acad. Sci. U.S.A.* 90:6909; Erb *et al.* (1994) *Proc. Natl. Acad. Sci. USA* 91:11422; Zuckermann *et al.* (1994). *J. Med. Chem.* 37:2678; Cho *et al.* (1993) *Science* 261:1303; Carrell *et al.* (1994) *Angew. Chem. Int. Ed. Engl.* 33:2059; Carell *et al.* (1994) *Angew. Chem. Int. Ed. Engl.* 33:2061; and in Gallop *et al.* (1994) *J. Med. Chem.* 37:1233.

Libraries of compounds may be presented in solution (*e.g.*, Houghten, 1992, *Biotechniques* 13:412-421), or on beads (Lam, 1991, *Nature* 354:82-84), chips (Fodor, 1993, *Nature* 364:555-556), bacteria and/or spores, (Ladner, USP 5,223,409), plasmids (Cull *et al.*, 1992, *Proc Natl Acad Sci USA* 89:1865-1869) or on phage (Scott and Smith, 1990, *Science* 249:386-390; Devlin, 1990, *Science* 249:404-406; Cwirla *et al.*, 1990, *Proc. Natl. Acad. Sci.* 87:6378-6382; Felici, 1991, *J. Mol. Biol.* 222:301-310; Ladner, *supra.*).

In one embodiment, the invention provides assays for screening candidate or test compounds which are substrates of a marker or biologically active portion thereof. In another embodiment, the invention provides assays for screening candidate or test compounds which bind to a marker or biologically active portion thereof. Determining the ability of the test compound to directly bind to a marker can be accomplished, for example, by coupling the compound with a radioisotope or enzymatic label such that binding of the compound to the marker can be determined by detecting the labeled marker compound in a complex. For example, compounds (*e.g.*, marker substrates) can be labeled with  $^{125}\text{I}$ ,  $^{35}\text{S}$ ,  $^{14}\text{C}$ , or  $^3\text{H}$ , either directly or indirectly, and the radioisotope detected by direct counting of radioemission or by scintillation counting. Alternatively, assay components can be enzymatically labeled with, for example, horseradish peroxidase, alkaline phosphatase, or luciferase, and the enzymatic label detected by determination of conversion of an appropriate substrate to product.

In another embodiment, the invention provides assays for screening candidate or test compounds which modulate the activity of a marker or a biologically active portion thereof. In all likelihood, the marker can, *in vivo*, interact with one or more molecules, such as but not limited to, peptides, proteins, hormones, cofactors and nucleic acids. For the purposes of this discussion, such cellular and extracellular molecules are referred to herein as "binding partners" or marker "substrate".

One necessary embodiment of the invention in order to facilitate such screening is the use of the marker to identify its natural *in vivo* binding partners. There are many ways to accomplish this which are known to one skilled in the art. One example is the use of the marker protein as "bait protein" in a two-hybrid assay or three-hybrid assay  
5 (see, *e.g.*, U.S. Patent No. 5,283,317; Zervos *et al*, 1993, *Cell* 72:223-232; Madura *et al*, 1993, *J. Biol. Chem.* 268:12046-12054; Bartel *et al*, 1993, *Biotechniques* 14:920-924; Iwabuchi *et al*, 1993 *Oncogene* 8:1693-1696; Brent WO94/10300) in order to identify other proteins which bind to or interact with the marker (binding partners) and, therefore, are possibly involved in the natural function of the marker. Such marker  
10 binding partners are also likely to be involved in the propagation of signals by the marker or downstream elements of a marker-mediated signaling pathway. Alternatively, such marker binding partners may also be found to be inhibitors of the marker.

The two-hybrid system is based on the modular nature of most transcription factors, which consist of separable DNA-binding and activation domains. Briefly, the  
15 assay utilizes two different DNA constructs. In one construct, the gene that encodes a marker protein fused to a gene encoding the DNA binding domain of a known transcription factor (*e.g.*, GAL-4). In the other construct, a DNA sequence, from a library of DNA sequences, that encodes an unidentified protein ("prey" or "sample") is fused to a gene that codes for the activation domain of the known transcription factor. If  
20 the "bait" and the "prey" proteins are able to interact, *in vivo*, forming a marker-dependent complex, the DNA-binding and activation domains of the transcription factor are brought into close proximity. This proximity allows transcription of a reporter gene (*e.g.*, LacZ) which is operably linked to a transcriptional regulatory site responsive to the transcription factor. Expression of the reporter gene can be readily detected and cell  
25 colonies containing the functional transcription factor can be isolated and used to obtain the cloned gene which encodes the protein which interacts with the marker protein.

In a further embodiment, assays may be devised through the use of the invention for the purpose of identifying compounds which modulate (*e.g.*, affect either positively or negatively) interactions between a marker and its substrates and/or binding partners.  
30 Such compounds can include, but are not limited to, molecules such as antibodies, peptides, hormones, oligonucleotides, nucleic acids, and analogs thereof. Such compounds may also be obtained from any available source, including systematic

libraries of natural and/or synthetic compounds. The preferred assay components for use in this embodiment is an ovarian cancer marker identified herein, the known binding partner and/or substrate of same, and the test compound. Test compounds can be supplied from any source.

5           The basic principle of the assay systems used to identify compounds that interfere with the interaction between the marker and its binding partner involves preparing a reaction mixture containing the marker and its binding partner under conditions and for a time sufficient to allow the two products to interact and bind, thus forming a complex. In order to test an agent for inhibitory activity, the reaction mixture  
10 is prepared in the presence and absence of the test compound. The test compound can be initially included in the reaction mixture, or can be added at a time subsequent to the addition of the marker and its binding partner. Control reaction mixtures are incubated without the test compound or with a placebo. The formation of any complexes between the marker and its binding partner is then detected. The formation of a complex in the  
15 control reaction, but less or no such formation in the reaction mixture containing the test compound, indicates that the compound interferes with the interaction of the marker and its binding partner. Conversely, the formation of more complex in the presence of compound than in the control reaction indicates that the compound may enhance interaction of the marker and its binding partner.

20           The assay for compounds that interfere with the interaction of the marker with its binding partner may be conducted in a heterogeneous or homogeneous format. Heterogeneous assays involve anchoring either the marker or its binding partner onto a solid phase and detecting complexes anchored to the solid phase at the end of the reaction. In homogeneous assays, the entire reaction is carried out in a liquid phase. In  
25 either approach, the order of addition of reactants can be varied to obtain different information about the compounds being tested. For example, test compounds that interfere with the interaction between the markers and the binding partners (*e.g.*, by competition) can be identified by conducting the reaction in the presence of the test substance, *i.e.*, by adding the test substance to the reaction mixture prior to or  
30 simultaneously with the marker and its interactive binding partner. Alternatively, test compounds that disrupt preformed complexes, *e.g.*, compounds with higher binding constants that displace one of the components from the complex, can be tested by adding

the test compound to the reaction mixture after complexes have been formed. The various formats are briefly described below.

In a heterogeneous assay system, either the marker or its binding partner is anchored onto a solid surface or matrix, while the other corresponding non-anchored component may be labeled, either directly or indirectly. In practice, microtitre plates are often utilized for this approach. The anchored species can be immobilized by a number of methods, either non-covalent or covalent, that are typically well known to one who practices the art. Non-covalent attachment can often be accomplished simply by coating the solid surface with a solution of the marker or its binding partner and drying.

10 Alternatively, an immobilized antibody specific for the assay component to be anchored can be used for this purpose. Such surfaces can often be prepared in advance and stored.

In related embodiments, a fusion protein can be provided which adds a domain that allows one or both of the assay components to be anchored to a matrix. For example, glutathione-S-transferase/marker fusion proteins or glutathione-S-transferase/binding partner can be adsorbed onto glutathione sepharose beads (Sigma Chemical, St. Louis, MO) or glutathione derivatized microtiter plates, which are then combined with the test compound or the test compound and either the non-adsorbed marker or its binding partner, and the mixture incubated under conditions conducive to complex formation (*e.g.*, physiological conditions). Following incubation, the beads or microtiter plate wells are washed to remove any unbound assay components, the immobilized complex assessed either directly or indirectly, for example, as described above. Alternatively, the complexes can be dissociated from the matrix, and the level of marker binding or activity determined using standard techniques.

15  
20

Other techniques for immobilizing proteins on matrices can also be used in the screening assays of the invention. For example, either a marker or a marker binding partner can be immobilized utilizing conjugation of biotin and streptavidin. Biotinylated marker protein or target molecules can be prepared from biotin-NHS (N-hydroxy-succinimide) using techniques known in the art (*e.g.*, biotinylation kit, Pierce Chemicals, Rockford, IL), and immobilized in the wells of streptavidin-coated 96 well plates (Pierce Chemical). In certain embodiments, the protein-immobilized surfaces can be prepared in advance and stored.

25  
30

In order to conduct the assay, the corresponding partner of the immobilized assay component is exposed to the coated surface with or without the test compound. After the reaction is complete, unreacted assay components are removed (*e.g.*, by washing) and any complexes formed will remain immobilized on the solid surface. The detection  
5 of complexes anchored on the solid surface can be accomplished in a number of ways. Where the non-immobilized component is pre-labeled, the detection of label immobilized on the surface indicates that complexes were formed. Where the non-immobilized component is not pre-labeled, an indirect label can be used to detect complexes anchored on the surface; *e.g.*, using a labeled antibody specific for the  
10 initially non-immobilized species (the antibody, in turn, can be directly labeled or indirectly labeled with, *e.g.*, a labeled anti-Ig antibody). Depending upon the order of addition of reaction components, test compounds which modulate (inhibit or enhance) complex formation or which disrupt preformed complexes can be detected.

In an alternate embodiment of the invention, a homogeneous assay may be used.  
15 This is typically a reaction, analogous to those mentioned above, which is conducted in a liquid phase in the presence or absence of the test compound. The formed complexes are then separated from unreacted components, and the amount of complex formed is determined. As mentioned for heterogeneous assay systems, the order of addition of reactants to the liquid phase can yield information about which test compounds  
20 modulate (inhibit or enhance) complex formation and which disrupt preformed complexes.

In such a homogeneous assay, the reaction products may be separated from unreacted assay components by any of a number of standard techniques, including but not limited to: differential centrifugation, chromatography, electrophoresis and  
25 immunoprecipitation. In differential centrifugation, complexes of molecules may be separated from uncomplexed molecules through a series of centrifugal steps, due to the different sedimentation equilibria of complexes based on their different sizes and densities (see, for example, Rivas, G., and Minton, A.P., *Trends Biochem Sci* 1993 Aug;18(8):284-7). Standard chromatographic techniques may also be utilized to separate  
30 complexed molecules from uncomplexed ones. For example, gel filtration chromatography separates molecules based on size, and through the utilization of an appropriate gel filtration resin in a column format, for example, the relatively larger



- 76 -

complex may be separated from the relatively smaller uncomplexed components. Similarly, the relatively different charge properties of the complex as compared to the uncomplexed molecules may be exploited to differentially separate the complex from the remaining individual reactants, for example through the use of ion-exchange chromatography resins. Such resins and chromatographic techniques are well known to one skilled in the art (see, *e.g.*, Heegaard, 1998, *J Mol. Recognit.* 11:141-148; Hage and Tweed, 1997, *J. Chromatogr. B. Biomed. Sci. Appl.*, 699:499-525). Gel electrophoresis may also be employed to separate complexed molecules from unbound species (see, *e.g.*, Ausubel *et al* (eds.), In: *Current Protocols in Molecular Biology*, J. Wiley & Sons, New York, 1999). In this technique, protein or nucleic acid complexes are separated based on size or charge, for example. In order to maintain the binding interaction during the electrophoretic process, nondenaturing gels in the absence of reducing agent are typically preferred, but conditions appropriate to the particular interactants will be well known to one skilled in the art. Immunoprecipitation is another common technique utilized for the isolation of a protein-protein complex from solution (see, *e.g.*, Ausubel *et al* (eds.), In: *Current Protocols in Molecular Biology*, J. Wiley & Sons, New York, 1999). In this technique, all proteins binding to an antibody specific to one of the binding molecules are precipitated from solution by conjugating the antibody to a polymer bead that may be readily collected by centrifugation. The bound assay components are released from the beads (through a specific proteolysis event or other technique well known in the art which will not disturb the protein-protein interaction in the complex), and a second immunoprecipitation step is performed, this time utilizing antibodies specific for the correspondingly different interacting assay component. In this manner, only formed complexes should remain attached to the beads. Variations in complex formation in both the presence and the absence of a test compound can be compared, thus offering information about the ability of the compound to modulate interactions between the marker and its binding partner.

Also within the scope of the present invention are methods for direct detection of interactions between the marker and its natural binding partner and/or a test compound in a homogeneous or heterogeneous assay system without further sample manipulation. For example, the technique of fluorescence energy transfer may be utilized (see, *e.g.*, Lakowicz *et al*, U.S. Patent No. 5,631,169; Stavrianopoulos *et al*, U.S. Patent No.

4,868,103). Generally, this technique involves the addition of a fluorophore label on a first 'donor' molecule (*e.g.*, marker or test compound) such that its emitted fluorescent energy will be absorbed by a fluorescent label on a second, 'acceptor' molecule (*e.g.*, marker or test compound), which in turn is able to fluoresce due to the absorbed energy.

5 Alternately, the 'donor' protein molecule may simply utilize the natural fluorescent energy of tryptophan residues. Labels are chosen that emit different wavelengths of light, such that the 'acceptor' molecule label may be differentiated from that of the 'donor'. Since the efficiency of energy transfer between the labels is related to the distance separating the molecules, spatial relationships between the molecules can be

10 assessed. In a situation in which binding occurs between the molecules, the fluorescent emission of the 'acceptor' molecule label in the assay should be maximal. An FET binding event can be conveniently measured through standard fluorometric detection means well known in the art (*e.g.*, using a fluorimeter). A test substance which either enhances or hinders participation of one of the species in the preformed complex will

15 result in the generation of a signal variant to that of background. In this way, test substances that modulate interactions between a marker and its binding partner can be identified in controlled assays.

In another embodiment, modulators of marker expression are identified in a method wherein a cell is contacted with a candidate compound and the expression of

20 mRNA or protein, corresponding to a marker in the cell, is determined. The level of expression of mRNA or protein in the presence of the candidate compound is compared to the level of expression of mRNA or protein in the absence of the candidate compound. The candidate compound can then be identified as a modulator of marker expression based on this comparison. For example, when expression of marker mRNA

25 or protein is greater (statistically significantly greater) in the presence of the candidate compound than in its absence, the candidate compound is identified as a stimulator of marker mRNA or protein expression. Conversely, when expression of marker mRNA or protein is less (statistically significantly less) in the presence of the candidate compound than in its absence, the candidate compound is identified as an inhibitor of

30 marker mRNA or protein expression. The level of marker mRNA or protein expression in the cells can be determined by methods described herein for detecting marker mRNA or protein.

In another aspect, the invention pertains to a combination of two or more of the assays described herein. For example, a modulating agent can be identified using a cell-based or a cell free assay, and the ability of the agent to modulate the activity of a marker protein can be further confirmed *in vivo*, *e.g.*, in a whole animal model for  
5 cellular transformation and/or tumorigenesis.

This invention further pertains to novel agents identified by the above-described screening assays. Accordingly, it is within the scope of this invention to further use an agent identified as described herein in an appropriate animal model. For example, an agent identified as described herein (*e.g.*, an marker modulating agent, an antisense  
10 marker nucleic acid molecule, an marker-specific antibody, or an marker-binding partner) can be used in an animal model to determine the efficacy, toxicity, or side effects of treatment with such an agent. Alternatively, an agent identified as described herein can be used in an animal model to determine the mechanism of action of such an agent. Furthermore, this invention pertains to uses of novel agents identified by the  
15 above-described screening assays for treatments as described herein.

It is understood that appropriate doses of small molecule agents and protein or polypeptide agents depends upon a number of factors within the knowledge of the ordinarily skilled physician, veterinarian, or researcher. The dose(s) of these agents will vary, for example, depending upon the identity, size, and condition of the subject or  
20 sample being treated, further depending upon the route by which the composition is to be administered, if applicable, and the effect which the practitioner desires the agent to have upon the nucleic acid or polypeptide of the invention. Exemplary doses of a small molecule include milligram or microgram amounts per kilogram of subject or sample weight (*e.g.* about 1 microgram per kilogram to about 500 milligrams per kilogram,  
25 about 100 micrograms per kilogram to about 5 milligrams per kilogram, or about 1 microgram per kilogram to about 50 micrograms per kilogram). Exemplary doses of a protein or polypeptide include gram, milligram or microgram amounts per kilogram of subject or sample weight (*e.g.* about 1 microgram per kilogram to about 5 grams per kilogram, about 100 micrograms per kilogram to about 500 milligrams per kilogram, or  
30 about 1 milligram per kilogram to about 50 milligrams per kilogram). It is furthermore understood that appropriate doses of one of these agents depend upon the potency of the agent with respect to the expression or activity to be modulated. Such appropriate doses

can be determined using the assays described herein. When one or more of these agents is to be administered to an animal (*e.g.* a human) in order to modulate expression or activity of a polypeptide or nucleic acid of the invention, a physician, veterinarian, or researcher can, for example, prescribe a relatively low dose at first, subsequently  
5 increasing the dose until an appropriate response is obtained. In addition, it is understood that the specific dose level for any particular animal subject will depend upon a variety of factors including the activity of the specific agent employed, the age, body weight, general health, gender, and diet of the subject, the time of administration, the route of administration, the rate of excretion, any drug combination, and the degree  
10 of expression or activity to be modulated.

A pharmaceutical composition of the invention is formulated to be compatible with its intended route of administration. Examples of routes of administration include parenteral, *e.g.*, intravenous, intradermal, subcutaneous, oral (*e.g.*, inhalation), transdermal (topical), transmucosal, and rectal administration. Solutions or suspensions  
15 used for parenteral, intradermal, or subcutaneous application can include the following components: a sterile diluent such as water for injection, saline solution, fixed oils, polyethylene glycols, glycerine, propylene glycol or other synthetic solvents; antibacterial agents such as benzyl alcohol or methyl parabens; antioxidants such as ascorbic acid or sodium bisulfite; chelating agents such as ethylenediamine-tetraacetic  
20 acid; buffers such as acetates, citrates or phosphates and agents for the adjustment of tonicity such as sodium chloride or dextrose. pH can be adjusted with acids or bases, such as hydrochloric acid or sodium hydroxide. The parenteral preparation can be enclosed in ampules, disposable syringes or multiple dose vials made of glass or plastic.

Pharmaceutical compositions suitable for injectable use include sterile aqueous  
25 solutions (where water soluble) or dispersions and sterile powders for the extemporaneous preparation of sterile injectable solutions or dispersions. For intravenous administration, suitable carriers include physiological saline, bacteriostatic water, Cremophor EL (BASF; Parsippany, NJ) or phosphate buffered saline (PBS). In all cases, the composition must be sterile and should be fluid to the extent that easy  
30 syringability exists. It must be stable under the conditions of manufacture and storage and must be preserved against the contaminating action of microorganisms such as bacteria and fungi. The carrier can be a solvent or dispersion medium containing, for

example, water, ethanol, polyol (for example, glycerol, propylene glycol, and liquid polyethylene glycol, and the like), and suitable mixtures thereof. The proper fluidity can be maintained, for example, by the use of a coating such as lecithin, by the maintenance of the required particle size in the case of dispersion and by the use of surfactants.

- 5 Prevention of the action of microorganisms can be achieved by various antibacterial and antifungal agents, for example, parabens, chlorobutanol, phenol, ascorbic acid, thimerosal, and the like. In many cases, it will be preferable to include isotonic agents, for example, sugars, polyalcohols such as mannitol, sorbitol, or sodium chloride in the composition. Prolonged absorption of the injectable compositions can be brought about  
10 by including in the composition an agent which delays absorption, for example, aluminum monostearate and gelatin.

- Sterile injectable solutions can be prepared by incorporating the active compound (*e.g.*, a polypeptide or antibody) in the required amount in an appropriate solvent with one or a combination of ingredients enumerated above, as required,  
15 followed by filtered sterilization. Generally, dispersions are prepared by incorporating the active compound into a sterile vehicle which contains a basic dispersion medium, and then incorporating the required other ingredients from those enumerated above. In the case of sterile powders for the preparation of sterile injectable solutions, the preferred methods of preparation are vacuum drying and freeze-drying which yields a  
20 powder of the active ingredient plus any additional desired ingredient from a previously sterile-filtered solution thereof.

- Oral compositions generally include an inert diluent or an edible carrier. They can be enclosed in gelatin capsules or compressed into tablets. For the purpose of oral therapeutic administration, the active compound can be incorporated with excipients and  
25 used in the form of tablets, troches, or capsules. Oral compositions can also be prepared using a fluid carrier for use as a mouthwash, wherein the compound in the fluid carrier is applied orally and swished and expectorated or swallowed.

- Pharmaceutically compatible binding agents, and/or adjuvant materials can be included as part of the composition. The tablets, pills, capsules, troches, and the like can  
30 contain any of the following ingredients, or compounds of a similar nature: a binder such as microcrystalline cellulose, gum tragacanth or gelatin; an excipient such as starch or lactose, a disintegrating agent such as alginic acid, Primogel, or corn starch; a

lubricant such as magnesium stearate or Sterotes; a glidant such as colloidal silicon dioxide; a sweetening agent such as sucrose or saccharin; or a flavoring agent such as peppermint, methyl salicylate, or orange flavoring.

For administration by inhalation, the compounds are delivered in the form of an aerosol spray from a pressurized container or dispenser which contains a suitable propellant, *e.g.*, a gas such as carbon dioxide, or a nebulizer.

Systemic administration can also be by transmucosal or transdermal means. For transmucosal or transdermal administration, penetrants appropriate to the barrier to be permeated are used in the formulation. Such penetrants are generally known in the art, and include, for example, for transmucosal administration, detergents, bile salts, and fusidic acid derivatives. Transmucosal administration can be accomplished through the use of nasal sprays or suppositories. For transdermal administration, the active compounds are formulated into ointments, salves, gels, or creams as generally known in the art.

The compounds can also be prepared in the form of suppositories (*e.g.*, with conventional suppository bases such as cocoa butter and other glycerides) or retention enemas for rectal delivery.

In one embodiment, the active compounds are prepared with carriers that will protect the compound against rapid elimination from the body, such as a controlled release formulation, including implants and microencapsulated delivery systems. Biodegradable, biocompatible polymers can be used, such as ethylene vinyl acetate, polyanhydrides, polyglycolic acid, collagen, polyorthoesters, and polylactic acid. Methods for preparation of such formulations will be apparent to those skilled in the art. The materials can also be obtained commercially from Alza Corporation and Nova Pharmaceuticals, Inc. Liposomal suspensions (including liposomes having monoclonal antibodies incorporated therein or thereon) can also be used as pharmaceutically acceptable carriers. These can be prepared according to methods known to those skilled in the art, for example, as described in U.S. Patent No. 4,522,811.

It is especially advantageous to formulate oral or parenteral compositions in dosage unit form for ease of administration and uniformity of dosage. Dosage unit form as used herein refers to physically discrete units suited as unitary dosages for the subject to be treated; each unit containing a predetermined quantity of active compound

calculated to produce the desired therapeutic effect in association with the required pharmaceutical carrier. The specification for the dosage unit forms of the invention are dictated by and directly dependent on the unique characteristics of the active compound and the particular therapeutic effect to be achieved, and the limitations inherent in the art of compounding such an active compound for the treatment of individuals.

For antibodies, the preferred dosage is 0.1 mg/kg to 100 mg/kg of body weight (generally 10 mg/kg to 20 mg/kg). If the antibody is to act in the brain, a dosage of 50 mg/kg to 100 mg/kg is usually appropriate. Generally, partially human antibodies and fully human antibodies have a longer half-life within the human body than other antibodies. Accordingly, lower dosages and less frequent administration is often possible. Modifications such as lipidation can be used to stabilize antibodies and to enhance uptake and tissue penetration (*e.g.*, into the ovarian epithelium). A method for lipidation of antibodies is described by Cruikshank *et al.* (1997) *J. Acquired Immune Deficiency Syndromes and Human Retrovirology* 14:193.

The nucleic acid molecules corresponding to a marker of the invention can be inserted into vectors and used as gene therapy vectors. Gene therapy vectors can be delivered to a subject by, for example, intravenous injection, local administration (U.S. Patent 5,328,470), or by stereotactic injection (see, *e.g.*, Chen *et al.*, 1994, *Proc. Natl. Acad. Sci. USA* 91:3054-3057). The pharmaceutical preparation of the gene therapy vector can include the gene therapy vector in an acceptable diluent, or can comprise a slow release matrix in which the gene delivery vehicle is imbedded. Alternatively, where the complete gene delivery vector can be produced intact from recombinant cells, *e.g.* retroviral vectors, the pharmaceutical preparation can include one or more cells which produce the gene delivery system.

The pharmaceutical compositions can be included in a container, pack, or dispenser together with instructions for administration.

#### V. Predictive Medicine

The present invention pertains to the field of predictive medicine in which diagnostic assays, prognostic assays, pharmacogenomics, and monitoring clinical trails are used for prognostic (predictive) purposes to thereby treat an individual prophylactically. Accordingly, one aspect of the present invention relates to diagnostic

assays for determining the level of expression of polypeptides or nucleic acids corresponding to one or more markers of the invention, in order to determine whether an individual is at risk of developing ovarian cancer. Such assays can be used for prognostic or predictive purposes to thereby prophylactically treat an individual prior to  
5 the onset of the cancer.

Yet another aspect of the invention pertains to monitoring the influence of agents (*e.g.*, drugs or other compounds administered either to inhibit ovarian cancer or to treat or prevent any other disorder {*i.e.* in order to understand any ovarian carcinogenic effects that such treatment may have} ) on the expression or activity of a marker of the  
10 invention in clinical trials. These and other agents are described in further detail in the following sections.

#### A. Diagnostic Assays

An exemplary method for detecting the presence or absence of a polypeptide or  
15 nucleic acid corresponding to a marker of the invention in a biological sample involves obtaining a biological sample (*e.g.* an ovary-associated body fluid) from a test subject and contacting the biological sample with a compound or an agent capable of detecting the polypeptide or nucleic acid (*e.g.*, mRNA, genomic DNA, or cDNA). The detection methods of the invention can thus be used to detect mRNA, protein, cDNA, or genomic  
20 DNA, for example, in a biological sample *in vitro* as well as *in vivo*. For example, *in vitro* techniques for detection of mRNA include Northern hybridizations and *in situ* hybridizations. *In vitro* techniques for detection of a polypeptide corresponding to a marker of the invention include enzyme linked immunosorbent assays (ELISAs), Western blots, immunoprecipitations and immunofluorescence. *In vitro* techniques for  
25 detection of genomic DNA include Southern hybridizations. Furthermore, *in vivo* techniques for detection of a polypeptide corresponding to a marker of the invention include introducing into a subject a labeled antibody directed against the polypeptide. For example, the antibody can be labeled with a radioactive marker whose presence and location in a subject can be detected by standard imaging techniques.

30 A general principle of such diagnostic and prognostic assays involves preparing a sample or reaction mixture that may contain a marker, and a probe, under appropriate conditions and for a time sufficient to allow the marker and probe to interact and bind,



thus forming a complex that can be removed and/or detected in the reaction mixture. These assays can be conducted in a variety of ways.

For example, one method to conduct such an assay would involve anchoring the marker or probe onto a solid phase support, also referred to as a substrate, and detecting  
5 target marker/probe complexes anchored on the solid phase at the end of the reaction. In one embodiment of such a method, a sample from a subject, which is to be assayed for presence and/or concentration of marker, can be anchored onto a carrier or solid phase support. In another embodiment, the reverse situation is possible, in which the probe can be anchored to a solid phase and a sample from a subject can be allowed to  
10 react as an unanchored component of the assay.

There are many established methods for anchoring assay components to a solid phase. These include, without limitation, marker or probe molecules which are immobilized through conjugation of biotin and streptavidin. Such biotinylated assay components can be prepared from biotin-NHS (N-hydroxy-succinimide) using  
15 techniques known in the art (*e.g.*, biotinylation kit, Pierce Chemicals, Rockford, IL), and immobilized in the wells of streptavidin-coated 96 well plates (Pierce Chemical). In certain embodiments, the surfaces with immobilized assay components can be prepared in advance and stored.

Other suitable carriers or solid phase supports for such assays include any  
20 material capable of binding the class of molecule to which the marker or probe belongs. Well-known supports or carriers include, but are not limited to, glass, polystyrene, nylon, polypropylene, nylon, polyethylene, dextran, amylases, natural and modified celluloses, polyacrylamides, gabbros, and magnetite.

In order to conduct assays with the above mentioned approaches, the non-  
25 immobilized component is added to the solid phase upon which the second component is anchored. After the reaction is complete, uncomplexed components may be removed (*e.g.*, by washing) under conditions such that any complexes formed will remain immobilized upon the solid phase. The detection of marker/probe complexes anchored to the solid phase can be accomplished in a number of methods outlined herein.

In a preferred embodiment, the probe, when it is the unanchored assay component, can be labeled for the purpose of detection and readout of the assay, either directly or indirectly, with detectable labels discussed herein and which are well-known to one skilled in the art.

5           It is also possible to directly detect marker/probe complex formation without further manipulation or labeling of either component (marker or probe), for example by utilizing the technique of fluorescence energy transfer (see, for example, Lakowicz *et al.*, U.S. Patent No. 5,631,169; Stavrianopoulos, *et al.*, U.S. Patent No. 4,868,103). A fluorophore label on the first, 'donor' molecule is selected such that, upon excitation  
10 with incident light of appropriate wavelength, its emitted fluorescent energy will be absorbed by a fluorescent label on a second 'acceptor' molecule, which in turn is able to fluoresce due to the absorbed energy. Alternately, the 'donor' protein molecule may simply utilize the natural fluorescent energy of tryptophan residues. Labels are chosen that emit different wavelengths of light, such that the 'acceptor' molecule label may be  
15 differentiated from that of the 'donor'. Since the efficiency of energy transfer between the labels is related to the distance separating the molecules, spatial relationships between the molecules can be assessed. In a situation in which binding occurs between the molecules, the fluorescent emission of the 'acceptor' molecule label in the assay should be maximal. An FET binding event can be conveniently measured through  
20 standard fluorometric detection means well known in the art (*e.g.*, using a fluorimeter).

          In another embodiment, determination of the ability of a probe to recognize a marker can be accomplished without labeling either assay component (probe or marker) by utilizing a technology such as real-time Biomolecular Interaction Analysis (BIA) (see, *e.g.*, Sjolander, S. and Urbaniczky, C., 1991, *Anal. Chem.* 63:2338-2345 and  
25 Szabo *et al.*, 1995, *Curr. Opin. Struct. Biol.* 5:699-705). As used herein, "BIA" or "surface plasmon resonance" is a technology for studying biospecific interactions in real time, without labeling any of the interactants (*e.g.*, BIAcore). Changes in the mass at the binding surface (indicative of a binding event) result in alterations of the refractive index of light near the surface (the optical phenomenon of surface plasmon resonance (SPR)),  
30 resulting in a detectable signal which can be used as an indication of real-time reactions between biological molecules.

Alternatively, in another embodiment, analogous diagnostic and prognostic assays can be conducted with marker and probe as solutes in a liquid phase. In such an assay, the complexed marker and probe are separated from uncomplexed components by any of a number of standard techniques, including but not limited to: differential

5 centrifugation, chromatography, electrophoresis and immunoprecipitation. In differential centrifugation, marker/probe complexes may be separated from uncomplexed assay components through a series of centrifugal steps, due to the different sedimentation equilibria of complexes based on their different sizes and densities (see, for example, Rivas, G., and Minton, A.P., 1993, *Trends Biochem Sci.* 18(8):284-7).

10 Standard chromatographic techniques may also be utilized to separate complexed molecules from uncomplexed ones. For example, gel filtration chromatography separates molecules based on size, and through the utilization of an appropriate gel filtration resin in a column format, for example, the relatively larger complex may be separated from the relatively smaller uncomplexed components. Similarly, the

15 relatively different charge properties of the marker/probe complex as compared to the uncomplexed components may be exploited to differentiate the complex from uncomplexed components, for example through the utilization of ion-exchange chromatography resins. Such resins and chromatographic techniques are well known to one skilled in the art (see, *e.g.*, Heegaard, N.H., 1998, *J. Mol. Recognit.* Winter 11(1-6):141-8; Hage, D.S., and Tweed, S.A. *J Chromatogr B Biomed Sci Appl* 1997 Oct 20 10;699(1-2):499-525). Gel electrophoresis may also be employed to separate complexed assay components from unbound components (see, *e.g.*, Ausubel *et al.*, ed., *Current Protocols in Molecular Biology*, John Wiley & Sons, New York, 1987-1999). In this technique, protein or nucleic acid complexes are separated based on size or

25 charge, for example. In order to maintain the binding interaction during the electrophoretic process, non-denaturing gel matrix materials and conditions in the absence of reducing agent are typically preferred. Appropriate conditions to the particular assay and components thereof will be well known to one skilled in the art.

In a particular embodiment, the level of mRNA corresponding to the marker can

30 be determined both by *in situ* and by *in vitro* formats in a biological sample using methods known in the art. The term "biological sample" is intended to include tissues, cells, biological fluids and isolates thereof, isolated from a subject, as well as tissues,

cells and fluids present within a subject. Many expression detection methods use isolated RNA. For *in vitro* methods, any RNA isolation technique that does not select against the isolation of mRNA can be utilized for the purification of RNA from ovarian cells (see, *e.g.*, Ausubel *et al.*, ed., *Current Protocols in Molecular Biology*, John Wiley & Sons, New York 1987-1999). Additionally, large numbers of tissue samples can readily be processed using techniques well known to those of skill in the art, such as, for example, the single-step RNA isolation process of Chomczynski (1989, U.S. Patent No. 4,843,155).

The isolated mRNA can be used in hybridization or amplification assays that include, but are not limited to, Southern or Northern analyses, polymerase chain reaction analyses and probe arrays. One preferred diagnostic method for the detection of mRNA levels involves contacting the isolated mRNA with a nucleic acid molecule (probe) that can hybridize to the mRNA encoded by the gene being detected. The nucleic acid probe can be, for example, a full-length cDNA, or a portion thereof, such as an oligonucleotide of at least 7, 15, 30, 50, 100, 250 or 500 nucleotides in length and sufficient to specifically hybridize under stringent conditions to a mRNA or genomic DNA encoding a marker of the present invention. Other suitable probes for use in the diagnostic assays of the invention are described herein. Hybridization of an mRNA with the probe indicates that the marker in question is being expressed.

In one format, the mRNA is immobilized on a solid surface and contacted with a probe, for example by running the isolated mRNA on an agarose gel and transferring the mRNA from the gel to a membrane, such as nitrocellulose. In an alternative format, the probe(s) are immobilized on a solid surface and the mRNA is contacted with the probe(s), for example, in an Affymetrix gene chip array. A skilled artisan can readily adapt known mRNA detection methods for use in detecting the level of mRNA encoded by the markers of the present invention.

An alternative method for determining the level of mRNA corresponding to a marker of the present invention in a sample involves the process of nucleic acid amplification, *e.g.*, by rtPCR (the experimental embodiment set forth in Mullis, 1987, U.S. Patent No. 4,683,202), ligase chain reaction (Barany, 1991, *Proc. Natl. Acad. Sci. USA*, 88:189-193), self sustained sequence replication (Guatelli *et al.*, 1990, *Proc. Natl. Acad. Sci. USA* 87:1874-1878), transcriptional amplification system (Kwoh *et al.*, 1989,

*Proc. Natl. Acad. Sci. USA* 86:1173-1177), Q-Beta Replicase (Lizardi *et al.*, 1988, *Bio/Technology* 6:1197), rolling circle replication (Lizardi *et al.*, U.S. Patent No. 5,854,033) or any other nucleic acid amplification method, followed by the detection of the amplified molecules using techniques well known to those of skill in the art. These  
5 detection schemes are especially useful for the detection of nucleic acid molecules if such molecules are present in very low numbers. As used herein, amplification primers are defined as being a pair of nucleic acid molecules that can anneal to 5' or 3' regions of a gene (plus and minus strands, respectively, or vice-versa) and contain a short region in between. In general, amplification primers are from about 10 to 30 nucleotides in  
10 length and flank a region from about 50 to 200 nucleotides in length. Under appropriate conditions and with appropriate reagents, such primers permit the amplification of a nucleic acid molecule comprising the nucleotide sequence flanked by the primers.

For *in situ* methods, mRNA does not need to be isolated from the ovarian cells prior to detection. In such methods, a cell or tissue sample is prepared/processed using  
15 known histological methods. The sample is then immobilized on a support, typically a glass slide, and then contacted with a probe that can hybridize to mRNA that encodes the marker.

As an alternative to making determinations based on the absolute expression level of the marker, determinations may be based on the normalized expression level of  
20 the marker. Expression levels are normalized by correcting the absolute expression level of a marker by comparing its expression to the expression of a gene that is not a marker, *e.g.*, a housekeeping gene that is constitutively expressed. Suitable genes for normalization include housekeeping genes such as the actin gene, or epithelial cell-specific genes. This normalization allows the comparison of the expression level in one  
25 sample, *e.g.*, a patient sample, to another sample, *e.g.*, a non-ovarian cancer sample, or between samples from different sources.

Alternatively, the expression level can be provided as a relative expression level. To determine a relative expression level of a marker, the level of expression of the marker is determined for 10 or more samples of normal versus cancer cell isolates,  
30 preferably 50 or more samples, prior to the determination of the expression level for the sample in question. The mean expression level of each of the genes assayed in the larger number of samples is determined and this is used as a baseline expression level

for the marker. The expression level of the marker determined for the test sample (absolute level of expression) is then divided by the mean expression value obtained for that marker. This provides a relative expression level.

Preferably, the samples used in the baseline determination will be from ovarian  
5 cancer or from non-ovarian cancer cells of ovarian tissue. The choice of the cell source is dependent on the use of the relative expression level. Using expression found in normal tissues as a mean expression score aids in validating whether the marker assayed is ovarian specific (versus normal cells). In addition, as more data is accumulated, the mean expression value can be revised, providing improved relative expression values  
10 based on accumulated data. Expression data from ovarian cells provides a means for grading the severity of the ovarian cancer state.

In another embodiment of the present invention, a polypeptide corresponding to a marker is detected. A preferred agent for detecting a polypeptide of the invention is an antibody capable of binding to a polypeptide corresponding to a marker of the invention,  
15 preferably an antibody with a detectable label. Antibodies can be polyclonal, or more preferably, monoclonal. An intact antibody, or a fragment thereof (*e.g.*, Fab or F(ab')<sub>2</sub>) can be used. The term "labeled", with regard to the probe or antibody, is intended to encompass direct labeling of the probe or antibody by coupling (*i.e.*, physically linking) a detectable substance to the probe or antibody, as well as indirect labeling of the probe  
20 or antibody by reactivity with another reagent that is directly labeled. Examples of indirect labeling include detection of a primary antibody using a fluorescently labeled secondary antibody and end-labeling of a DNA probe with biotin such that it can be detected with fluorescently labeled streptavidin.

Proteins from ovarian cells can be isolated using techniques that are well known  
25 to those of skill in the art. The protein isolation methods employed can, for example, be such as those described in Harlow and Lane (Harlow and Lane, 1988, *Antibodies: A Laboratory Manual*, Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York).

A variety of formats can be employed to determine whether a sample contains a  
30 protein that binds to a given antibody. Examples of such formats include, but are not limited to, enzyme immunoassay (EIA), radioimmunoassay (RIA), Western blot analysis and enzyme linked immunoabsorbant assay (ELISA). A skilled artisan can

readily adapt known protein/antibody detection methods for use in determining whether ovarian cells express a marker of the present invention.

In one format, antibodies, or antibody fragments, can be used in methods such as Western blots or immunofluorescence techniques to detect the expressed proteins. In such uses, it is generally preferable to immobilize either the antibody or proteins on a solid support. Suitable solid phase supports or carriers include any support capable of binding an antigen or an antibody. Well-known supports or carriers include glass, polystyrene, polypropylene, polyethylene, dextran, nylon, amyloses, natural and modified celluloses, polyacrylamides, gabbros, and magnetite.

One skilled in the art will know many other suitable carriers for binding antibody or antigen, and will be able to adapt such support for use with the present invention. For example, protein isolated from ovarian cells can be run on a polyacrylamide gel electrophoresis and immobilized onto a solid phase support such as nitrocellulose. The support can then be washed with suitable buffers followed by treatment with the detectably labeled antibody. The solid phase support can then be washed with the buffer a second time to remove unbound antibody. The amount of bound label on the solid support can then be detected by conventional means.

The invention also encompasses kits for detecting the presence of a polypeptide or nucleic acid corresponding to a marker of the invention in a biological sample (*e.g.* an ovary-associated body fluid such as a urine sample). Such kits can be used to determine if a subject is suffering from or is at increased risk of developing ovarian cancer. For example, the kit can comprise a labeled compound or agent capable of detecting a polypeptide or an mRNA encoding a polypeptide corresponding to a marker of the invention in a biological sample and means for determining the amount of the polypeptide or mRNA in the sample (*e.g.*, an antibody which binds the polypeptide or an oligonucleotide probe which binds to DNA or mRNA encoding the polypeptide). Kits can also include instructions for interpreting the results obtained using the kit.

For antibody-based kits, the kit can comprise, for example: (1) a first antibody (*e.g.*, attached to a solid support) which binds to a polypeptide corresponding to a marker of the invention; and, optionally, (2) a second, different antibody which binds to either the polypeptide or the first antibody and is conjugated to a detectable label.

For oligonucleotide-based kits, the kit can comprise, for example: (1) an oligonucleotide, *e.g.*, a detectably labeled oligonucleotide, which hybridizes to a nucleic acid sequence encoding a polypeptide corresponding to a marker of the invention or (2) a pair of primers useful for amplifying a nucleic acid molecule corresponding to a marker of the invention. The kit can also comprise, *e.g.*, a buffering agent, a preservative, or a protein stabilizing agent. The kit can further comprise components necessary for detecting the detectable label (*e.g.*, an enzyme or a substrate). The kit can also contain a control sample or a series of control samples which can be assayed and compared to the test sample. Each component of the kit can be enclosed within an individual container and all of the various containers can be within a single package, along with instructions for interpreting the results of the assays performed using the kit.

#### B. Pharmacogenomics

Agents or modulators which have a stimulatory or inhibitory effect on expression of a marker of the invention can be administered to individuals to treat (prophylactically or therapeutically) ovarian cancer in the patient. In conjunction with such treatment, the pharmacogenomics (*i.e.*, the study of the relationship between an individual's genotype and that individual's response to a foreign compound or drug) of the individual may be considered. Differences in metabolism of therapeutics can lead to severe toxicity or therapeutic failure by altering the relation between dose and blood concentration of the pharmacologically active drug. Thus, the pharmacogenomics of the individual permits the selection of effective agents (*e.g.*, drugs) for prophylactic or therapeutic treatments based on a consideration of the individual's genotype. Such pharmacogenomics can further be used to determine appropriate dosages and therapeutic regimens. Accordingly, the level of expression of a marker of the invention in an individual can be determined to thereby select appropriate agent(s) for therapeutic or prophylactic treatment of the individual.

Pharmacogenomics deals with clinically significant variations in the response to drugs due to altered drug disposition and abnormal action in affected persons. See, *e.g.*, Linder (1997) *Clin. Chem.* 43(2):254-266. In general, two types of pharmacogenetic conditions can be differentiated. Genetic conditions transmitted as a single factor altering the way drugs act on the body are referred to as "altered drug action." Genetic



conditions transmitted as single factors altering the way the body acts on drugs are referred to as "altered drug metabolism". These pharmacogenetic conditions can occur either as rare defects or as polymorphisms. For example, glucose-6-phosphate dehydrogenase (G6PD) deficiency is a common inherited enzymopathy in which the  
5 main clinical complication is hemolysis after ingestion of oxidant drugs (anti-malarials, sulfonamides, analgesics, nitrofurans) and consumption of fava beans.

As an illustrative embodiment, the activity of drug metabolizing enzymes is a major determinant of both the intensity and duration of drug action. The discovery of genetic polymorphisms of drug metabolizing enzymes (*e.g.*, N-acetyltransferase 2 (NAT  
10 2) and cytochrome P450 enzymes CYP2D6 and CYP2C19) has provided an explanation as to why some patients do not obtain the expected drug effects or show exaggerated drug response and serious toxicity after taking the standard and safe dose of a drug. These polymorphisms are expressed in two phenotypes in the population, the extensive metabolizer (EM) and poor metabolizer (PM). The prevalence of PM is different among  
15 different populations. For example, the gene coding for CYP2D6 is highly polymorphic and several mutations have been identified in PM, which all lead to the absence of functional CYP2D6. Poor metabolizers of CYP2D6 and CYP2C19 quite frequently experience exaggerated drug response and side effects when they receive standard doses. If a metabolite is the active therapeutic moiety, a PM will show no therapeutic  
20 response, as demonstrated for the analgesic effect of codeine mediated by its CYP2D6-formed metabolite morphine. The other extreme are the so called ultra-rapid metabolizers who do not respond to standard doses. Recently, the molecular basis of ultra-rapid metabolism has been identified to be due to CYP2D6 gene amplification.

Thus, the level of expression of a marker of the invention in an individual can be  
25 determined to thereby select appropriate agent(s) for therapeutic or prophylactic treatment of the individual. In addition, pharmacogenetic studies can be used to apply genotyping of polymorphic alleles encoding drug-metabolizing enzymes to the identification of an individual's drug responsiveness phenotype. This knowledge, when applied to dosing or drug selection, can avoid adverse reactions or therapeutic failure  
30 and thus enhance therapeutic or prophylactic efficiency when treating a subject with a modulator of expression of a marker of the invention.

### C. Monitoring Clinical Trials

Monitoring the influence of agents (*e.g.*, drug compounds) on the level of expression of a marker of the invention can be applied not only in basic drug screening, but also in clinical trials. For example, the effectiveness of an agent to affect marker expression can be monitored in clinical trials of subjects receiving treatment for ovarian cancer. In a preferred embodiment, the present invention provides a method for monitoring the effectiveness of treatment of a subject with an agent (*e.g.*, an agonist, antagonist, peptidomimetic, protein, peptide, nucleic acid, small molecule, or other drug candidate) comprising the steps of (i) obtaining a pre-administration sample from a subject prior to administration of the agent; (ii) detecting the level of expression of one or more selected markers of the invention in the pre-administration sample; (iii) obtaining one or more post-administration samples from the subject; (iv) detecting the level of expression of the marker(s) in the post-administration samples; (v) comparing the level of expression of the marker(s) in the pre-administration sample with the level of expression of the marker(s) in the post-administration sample or samples; and (vi) altering the administration of the agent to the subject accordingly. For example, increased administration of the agent can be desirable to increase expression of the marker(s) to higher levels than detected, *i.e.*, to increase the effectiveness of the agent. Alternatively, decreased administration of the agent can be desirable to decrease expression of the marker(s) to lower levels than detected, *i.e.*, to decrease the effectiveness of the agent.

### D. Surrogate Markers

The markers of the invention may serve as surrogate markers for one or more disorders or disease states or for conditions leading up to disease states, and in particular, ovarian cancer. As used herein, a "surrogate marker" is an objective biochemical marker which correlates with the absence or presence of a disease or disorder, or with the progression of a disease or disorder (*e.g.*, with the presence or absence of a tumor). The presence or quantity of such markers is independent of the disease. Therefore, these markers may serve to indicate whether a particular course of treatment is effective in lessening a disease state or disorder. Surrogate markers are of particular use when the presence or extent of a disease state or disorder is difficult to

assess through standard methodologies (*e.g.*, early stage tumors), or when an assessment of disease progression is desired before a potentially dangerous clinical endpoint is reached (*e.g.*, an assessment of cardiovascular disease may be made using cholesterol levels as a surrogate marker, and an analysis of HIV infection may be made using HIV RNA levels as a surrogate marker, well in advance of the undesirable clinical outcomes of myocardial infarction or fully-developed AIDS). Examples of the use of surrogate markers in the art include: Koomen *et al.* (2000) *J. Mass. Spectrom.* 35: 258-264; and James (1994) *AIDS Treatment News Archive* 209.

The markers of the invention are also useful as pharmacodynamic markers. As used herein, a “pharmacodynamic marker” is an objective biochemical marker which correlates specifically with drug effects. The presence or quantity of a pharmacodynamic marker is not related to the disease state or disorder for which the drug is being administered; therefore, the presence or quantity of the marker is indicative of the presence or activity of the drug in a subject. For example, a pharmacodynamic marker may be indicative of the concentration of the drug in a biological tissue, in that the marker is either expressed or transcribed or not expressed or transcribed in that tissue in relationship to the level of the drug. In this fashion, the distribution or uptake of the drug may be monitored by the pharmacodynamic marker. Similarly, the presence or quantity of the pharmacodynamic marker may be related to the presence or quantity of the metabolic product of a drug, such that the presence or quantity of the marker is indicative of the relative breakdown rate of the drug *in vivo*. Pharmacodynamic markers are of particular use in increasing the sensitivity of detection of drug effects, particularly when the drug is administered in low doses. Since even a small amount of a drug may be sufficient to activate multiple rounds of marker transcription or expression, the amplified marker may be in a quantity which is more readily detectable than the drug itself. Also, the marker may be more easily detected due to the nature of the marker itself; for example, using the methods described herein, antibodies may be employed in an immune-based detection system for a protein marker, or marker-specific radiolabeled probes may be used to detect a mRNA marker. Furthermore, the use of a pharmacodynamic marker may offer mechanism-based prediction of risk due to drug treatment beyond the range of possible direct observations. Examples of the use of pharmacodynamic markers in the art include:

Matsuda *et al.* US 6,033,862; Hattis *et al.* (1991) *Env. Health Perspect.* 90: 229-238; Schentag (1999) *Am. J. Health-Syst. Pharm.* 56 Suppl. 3: S21-S24; and Nicolau (1999) *Am. J. Health-Syst. Pharm.* 56 Suppl. 3: S16-S20.

The markers of the invention are also useful as pharmacogenomic markers. As  
5 used herein, a “pharmacogenomic marker” is an objective biochemical marker which  
correlates with a specific clinical drug response or susceptibility in a subject (see, e.g.,  
McLeod *et al.* (1999) *Eur. J. Cancer* 35(12): 1650-1652). The presence or quantity of  
the pharmacogenomic marker is related to the predicted response of the subject to a  
specific drug or class of drugs prior to administration of the drug. By assessing the  
10 presence or quantity of one or more pharmacogenomic markers in a subject, a drug  
therapy which is most appropriate for the subject, or which is predicted to have a greater  
degree of success, may be selected. For example, based on the presence or quantity of  
RNA or protein for specific tumor markers in a subject, a drug or course of treatment  
may be selected that is optimized for the treatment of the specific tumor likely to be  
15 present in the subject. Similarly, the presence or absence of a specific sequence  
mutation in marker DNA may correlate with drug response. The use of  
pharmacogenomic markers therefore permits the application of the most appropriate  
treatment for each subject without having to administer the therapy.

#### 20 E. Computer Readable Means and Arrays

Computer readable media comprising a marker of the present invention is also  
provided. As used herein, “computer readable media” refers to any medium that can be  
read and accessed directly by a computer. Such media include, but are not limited to:  
magnetic storage media, such as floppy discs, hard disc storage medium, and magnetic  
25 tape; optical storage media such as CD-ROM; electrical storage media such as RAM and  
ROM; and hybrids of these categories such as magnetic/optical storage media. The  
skilled artisan will readily appreciate how any of the presently known computer readable  
mediums can be used to create a manufacture comprising computer readable medium  
having recorded thereon a marker of the present invention.

As used herein, "recorded" refers to a process for storing information on computer readable medium. Those skilled in the art can readily adopt any of the presently known methods for recording information on computer readable medium to generate manufactures comprising the markers of the present invention.

5           A variety of data processor programs and formats can be used to store the marker information of the present invention on computer readable medium. For example, the nucleic acid sequence corresponding to the markers can be represented in a word processing text file, formatted in commercially-available software such as WordPerfect and MicroSoft Word, or represented in the form of an ASCII file, stored in a database  
10 application, such as DB2, Sybase, Oracle, or the like. Any number of dataprocessor structuring formats (*e.g.*, text file or database) may be adapted in order to obtain computer readable medium having recorded thereon the markers of the present invention.

By providing the markers of the invention in computer readable form, one can  
15 routinely access the marker sequence information for a variety of purposes. For example, one skilled in the art can use the nucleotide or amino acid sequences of the present invention in computer readable form to compare a target sequence or target structural motif with the sequence information stored within the data storage means. Search means are used to identify fragments or regions of the sequences of the invention  
20 which match a particular target sequence or target motif.

The invention also includes an array comprising a marker of the present invention. The array can be used to assay expression of one or more genes in the array. In one embodiment, the array can be used to assay gene expression in a tissue to ascertain tissue specificity of genes in the array. In this manner, up to about 7600 genes  
25 can be simultaneously assayed for expression. This allows a profile to be developed showing a battery of genes specifically expressed in one or more tissues.

In addition to such qualitative determination, the invention allows the quantitation of gene expression. Thus, not only tissue specificity, but also the level of expression of a battery of genes in the tissue is ascertainable. Thus, genes can be  
30 grouped on the basis of their tissue expression *per se* and level of expression in that tissue. This is useful, for example, in ascertaining the relationship of gene expression between or among tissues. Thus, one tissue can be perturbed and the effect on gene

expression in a second tissue can be determined. In this context, the effect of one cell type on another cell type in response to a biological stimulus can be determined. Such a determination is useful, for example, to know the effect of cell-cell interaction at the level of gene expression. If an agent is administered therapeutically to treat one cell type but has an undesirable effect on another cell type, the invention provides an assay to determine the molecular basis of the undesirable effect and thus provides the opportunity to co-administer a counteracting agent or otherwise treat the undesired effect. Similarly, even within a single cell type, undesirable biological effects can be determined at the molecular level. Thus, the effects of an agent on expression of other than the target gene can be ascertained and counteracted.

In another embodiment, the array can be used to monitor the time course of expression of one or more genes in the array. This can occur in various biological contexts, as disclosed herein, for example development and differentiation, tumor progression, progression of other diseases, *in vitro* processes, such a cellular transformation and senescence, autonomic neural and neurological processes, such as, for example, pain and appetite, and cognitive functions, such as learning or memory.

The array is also useful for ascertaining the effect of the expression of a gene on the expression of other genes in the same cell or in different cells. This provides, for example, for a selection of alternate molecular targets for therapeutic intervention if the ultimate or downstream target cannot be regulated.

The array is also useful for ascertaining differential expression patterns of one or more genes in normal and abnormal cells. This provides a battery of genes that could serve as a molecular target for diagnosis or therapeutic intervention.

## VI. Experimental Protocol

### A. Subtracted Libraries

Subtracted libraries are generated using a PCR based method that allows the isolation of clones expressed at higher levels in one population of mRNA (tester) compared to another population (driver). Both tester and driver mRNA populations are converted into cDNA by reverse transcription, and then PCR amplified using the SMART PCR kit from Clontech. Tester and driver cDNAs are then hybridized using

the PCR-Select cDNA subtraction kit from Clontech. This technique results in both subtraction and normalization, which is an equalization of copy number of low-abundance and high-abundance sequences. After generation of the subtractive libraries, a group of 96 or more clones from each library is tested to confirm differential  
5 expression by reverse Southern hybridization.

A first group of regular cDNA libraries was constructed. Library johOa was constructed from a pool of 5 normal ovarian epithelial cell cultures. Library johOb was constructed from a pool of 5 ascites short cultured samples from ovarian cancer patients. Library johOc was constructed from a pool of 6 serous late stage (III/IV) tumor samples.  
10 Three subtracted libraries were generated from these libraries: johOd, johOe and johOf. The johOd library was a subtracted ascites library, where the tester was johOb, and the driver was johOa. The johOe and johOf libraries were both subtracted stage III/IV serous tumor libraries. The tester for both of these libraries was johOc, and the driver was a pooled RNA from normal tissues. The tissues used for this driver pool were:  
15 kidney, small intestine, prostate, lung, heart, muscle, spleen, pancreas, liver, and lymphocyte. Library cMhOg was the same as the johOc and johOf libraries, with the exception that normal ovary was added to the driver. cMhOh, i, j, and k are all stage I/II subtracted libraries made from pooled tumor RNAs of different histological types (h=serous, I=endometrioid, j=clear cell, k=mucinous). The driver was the same for these  
20 4 libraries. It consisted of normal ovarian epithelial RNA and PBML RNA.

SEQ ID NOS: 1-2795 (Tables 1 and 1A) were identified through the above-described subtractive library hybridization techniques. In Tables 1 and 1A, SEQ ID NOS: 1-773 were from Library johOd; SEQ ID NOS: 774-1331 were from Library johOe; SEQ ID NOS: 1332-2795 were from Libraries johOf.  
25 SEQ ID NOS: 2796-10795 (Table 6) and 10796-10808 (Table 6A) were also identified through the above-described subtractive library hybridization techniques. In Table 6, SEQ ID NOS: 2796-3789 were from Library cMhOg; SEQ ID NOS: 3790-6301 were from Library cMhoh; SEQ ID NOS: 6302-8108 were from Libraries cMhoi; SEQ ID NOS: 8109-9981 were from Library cMhoj; SEQ ID NOS: 9982-10795 were  
30 from Libraries cMhok.

### B. Transcript Profiling

Nylon arrays were prepared by spotting purified PCR product onto a nylon membrane using a robotic gridding system linked to a sample database. Several thousand clones were spotted on each nylon filter.

5 RNA or DNA from clinical samples (tumor and normal), and cell lines as well as from subtracted libraries, were used for hybridization against the nylon arrays. The RNA or DNA is labeled utilizing an *in vitro* reverse transcription reaction that contains a radiolabeled nucleotide that is incorporated during the reaction. Hybridization experiments were carried out by combining labeled RNA or DNA samples with nylon  
10 filters in a hybridization chamber. Duplicate, independent hybridization experiments were performed to generate transcriptional profiling data. See, *Nature Genetics*, 21 (1999).

### C. Proteomics

15 Proteins that are secreted by normal and transformed cells in culture are analyzed to identify those proteins that are likely to be secreted by cancerous cells into body fluids. Supernatants are isolated and MWT-CO filters are used to simplify the mixture of proteins. The proteins are then digested with trypsin. The tryptic peptides are loaded onto a microcapillary HPLC column where they are separated, and eluted  
20 directly into an ion trap mass spectrometer, through a custom-made electrospray ionization source. Throughout the gradient, sequence data is acquired through fragmentation of the four most intense ions (peptides) that elute off the column, while dynamically excluding those that have already been fragmented. In this way, approximately 2000 scans worth of sequence data are obtained, corresponding to  
25 approximately 50 to 200 different proteins in the sample. These data are searched against databases using correlation analysis tools, such as MS-Tag, to identify the proteins in the supernatants.

The markers of Tables 7A-7E were identified through the above-described proteomics protocol. In particular, the proteins set forth in Tables 7A-7E were identified  
30 and their expression was analyzed in seven short term cultures of ovarian cancer cells (jov891N, jov915N, jov915p6N, jov928N, jov860N, jov908N and jov926N) and six



ovarian cancer cell lines (ov17TotN, ov167TotN, ov177TotN, ov202TotN, ov207TotN and ov266TotN).

#### D. Identification of Novel Genes

5 Sequences which displayed an increase in expression in [any one of] twelve late stage ovarian tumor samples over the corresponding average expression of four non-tumor samples were blasted against both public and proprietary sequence databases in order to identify other sequences with significant overlap. Contiguous sequences were then assembled into full length genes (cDNAs). Those cDNAs in which the potential  
10 open reading frame was still open at the 5' end were experimentally extended by either 5' RACE PCR or extracted from full length cDNA libraries by a PCR reaction between the vector and 5' end of the assembled electronic sequence. To predict whether an assembled gene encodes a potential integral membrane protein or not, hydropathy predictions of the predicted open reading frame was performed. If the open reading  
15 frame contained a predicted signal peptide in the N-terminal portion and a single membrane spanning domain, it was labeled as being a potential type I transmembrane protein. If the predicted amino acid sequence contained a transmembrane domain in the N-terminal portion of the protein, it was labeled as being a potential type II transmembrane protein. If the predicted amino acid sequence was a short hydrophobic  
20 protein (<50 amino acids), such as CD52 (CAMPATH), it was labeled as a potential integral membrane protein. If the predicted amino acid sequence contained multiple membrane spanning regions it was labeled as a type III transmembrane protein.

The novel genes of Table 8 were identified through the above-described procedure.

25

#### E. Northern Blot Analysis

Northern blots were performed for several of the genes of Table 8 to analyze for expression in normal human tissues. A clone was picked and served as a template for generation of probes for Northern blots. The probes were radiolabeled using <sup>32</sup>PdCTP  
30 using standard procedures and hybridized to Clontech (Palo Alto, California) human multiple tissue northern. Clontech Human MTN blot (catalog # 7760-1) contains heart, brain, placenta, lung, liver, skeletal muscle, kidney and pancreas. Human 12-Lane

MTN blot (catalog # 7780-1) contains brain, heart, skeletal muscle, colon, thymus, spleen, kidney, liver, small intestine, placenta, lung, peripheral blood leukocytes. Human MTN blot II (#7759-1) contains spleen, thymus, prostate, testis, ovary, small intestine, colon, and peripheral blood leukocytes. The hybridization and wash

5 conditions used were as described in the Clontech Multiple Tissue Northern (MTN) Blot User Manual (Catalogue number PT1200-1). Kodak biomax film was exposed to the Northern blot membrane for 10-72 hours, which were then developed.

Tables 10A-10N summarize the Northern blot analysis performed for several of the novel genes of Table 8.

10

#### F. Gene Expression Analysis

Total RNA from normal human tissue was obtained from commercial sources. The integrity of the RNA was verified by agarose gel electrophoresis and ethidium bromide staining. Cell lines were purchased from ATCC and grown under the

15 conditions recommended by ATCC. Total RNA from a number of various breast, ovarian and prostate adenocarcinoma cell lines was prepared using commercial kits (Qiagen). First strand cDNA was prepared using oligo-dT primer and standard conditions. Each RNA preparation was treated with DNase I (Ambion) at 37°C for 1 hour.

20 Novel gene expression was measured by TaqMan<sup>®</sup> quantitative PCR (Perkin Elmer Applied Biosystems) in cDNA prepared from the following normal human tissues: prostate, cerebellum, breast, ovary, kidney, trachea, adipose, small intestine, thyroid, testis, placenta, spinal cord, cervix, esophagus, spleen, thymus, brain, lung, skeletal muscle, heart, mammary gland, salivary gland, stomach, uterus, adrenal gland,

25 bladder, medulla hippocampus, and liver from one or two adult donors. Furthermore, novel gene expression was analyzed in the following cell lines: ZR-75-30, CAMA-1, MDA-MB-157, MDA-MB-175VII, MDA-MB-231, MDA-MB-361, SK-BR-3, BT-483, BT-549, DU4475, Hs578Bst, Hs578T, MDA-MB-453, T-47D, ES-2, Caov-3, SK-OV-3, NIH:OVCA-3, HTB-78, CRL-1572, CRL-10303, CA-HPV-10, CA-HPV-7, DU145,

30 MCF-7 and MDA-MB-468.

- 102 -

PCR Probes were designed by PrimerExpress software (PE Biosystems) based on the disclosed sequences of each novel human kinase gene. The primers and probes for expression analysis of the novel genes in Table 8 are given below:

5 Marker 10

Forward primer: F GATGACTTGAGAGAAGGTGCACAGT  
Reverse primer: R AAGGACAAGTGTGTTTGGCTTCA  
TaqMan probe: P TTTGATGCAGGCTGCTGGTCTTGG

10 Marker 15

Forward primer: F TGCAGCAGCCTGTGTATGC  
Reverse primer: R AAACAGCGACACGACAGTGAA  
TaqMan probe: P TTGGCTCCGGTATCGTCAACACGG

15 Marker 19

Forward primer: F AGTTCATCACGATATCAGGGAAGAT  
Reverse primer: R TGAATGATTACTGCCGATGTAGCT  
TaqMan probe: P CAAAGAGCCGTACGTCCACTGCCAGA

20 Marker 5

Forward primer: F GGCTGCTTTGCTGCAACTG  
Reverse primer: R CAGAGCGGGCAGCAGAATA  
TaqMan Probe P ACCCCGCACAGACAAGCCTTACTCC

25 Marker 8

Forward primer F TGTGTGCTGAAGGCTACATGTTG  
Reverse primer R TCTCCATGGCTGGTTTCCA  
TaqMan Probe P TTCTTACACGTCAGGTATTTGTAATCGCCCT

30 Marker 25

Forward Primer F CTCCCACCCCTTCTTCAATG  
Reverse primer R AGCTGTACTCTGCCGTTTCTC  
TaqMan Probe P ACCTTCGACTATGACATCGCGCTGCT

Marker 39

Forward primer	F	CCCGGAATGTGGTTTATGGTATT
Reverse primer	R	GACCGTCTTGTTGTGGAGTGAAG
TaqMan Probe	P	CCTTTCCTTGACCTCTATCGCAACCCGAA

5

An internal reference gene 18S rRNA was used. Primers and probe were purchased pretested from PE Applied Biosystems. Each gene probe was labeled using FAM (6-carboxyfluorescein), and the  $\beta$ 2-microglobulin reference probe was labeled with a different fluorescent dye, VIC. The differential labeling of the target gene and internal reference gene thus enabled measurement in same well. Forward and reverse primers and the probes for both 18S rRNA and target gene were added to the TaqMan® Universal PCR Master Mix (PE Applied Biosystems). Although the final concentration of primer and probe could vary, each was internally consistent within a given experiment. A typical experiment contained 900 nM of forward and reverse primers plus 250nM probe for the target gene whereas primers and probe for 18S rRNA were used according to manufacturer's recommendations. TaqMan matrix experiments were carried out on an ABI PRISM 7700 Sequence Detection System (PE Applied Biosystems). The thermal cycler conditions were as follows: hold for 2 min at 50°C and 10 min at 95°C, followed by two-step PCR for 40 cycles of 95°C for 15 sec followed by 60°C for 1 min.

The following method was used to quantitatively calculate gene expression in the various tissues relative to 18S rRNA expression in the same tissue. The threshold cycle (Ct) value is defined as the cycle at which a statistically significant increase in fluorescence is detected. A lower Ct value is indicative of a higher mRNA concentration. The Ct value of a given gene (Ct<sub>marker</sub>) is normalized by subtracting the Ct value of the 18S rRNA gene to obtain a  $\Delta$ Ct value using the following formula:  $\Delta$ Ct = Ct<sub>marker</sub> - Ct<sub>18S rRNA</sub>. Expression is then calibrated against a no template control sample. The  $\Delta$ Ct value for the calibrator sample is then subtracted from  $\Delta$ Ct for each tissue sample according to the following formula:  $\Delta\Delta$ Ct =  $\Delta$ Ct<sub>sample</sub> -  $\Delta$ Ct<sub>calibrator</sub>. Relative expression is then calculated using the arithmetic formula given by  $2^{-\Delta\Delta$ Ct}. Table 9 graphically represents the results of the TaqMan® expression study.

### G. LightCycler

The LightCycler Instrument from Boehringer Mannheim GmbH, Mannheim is a thermocycler for the rapid analysis of PCR applications. Fluorimetric analysis of the PCR products formed is performed as "real time" measurement either continuously or at a specifically defined time during each PCR cycle. The three detection channels of the LightCycler are fitted with filter combinations which allow analysis at the given emission wavelengths, thereby enabling exact sample measurement to be carried out in parallel with the fluorophores. SYBR Green I is a dye specific for double-stranded DNA. Its inherent fluorescence is enhanced by binding to the minor groove to ds DNA. The addition of SYBR Green I to PCR reactions allows the detection of PCR products formed by the binding of this fluorophore during each phase of DNA synthesis. The point of time of fluorimetric measurement is determined at the end of the elongation phase. The LightCycler- FastStart DNA Master SYBR Green I kit manufactured by Roche was used in order to quantify the copy number of a specific target. A panel of tumors and normal tissues were used to detect the expression levels of specific markers of the present invention in ovarian tumor samples compared to normal. The results are set forth in Table 11.

### H. RT-PCR

The Gibco BRL Superscript first strand synthesis system was used for RT-PCR to synthesize first strand cDNA from total RNA of ovarian tumors as well as normal ovary. Gene specific primers were designed for clones of the present invention using software program Oligo5.1. Finished sequence for these clones was available by in house sequencing efforts. Following the use of this system, target cDNA was amplified with the gene specific primers. Presence of a band in a sample indicates that the gene is upregulated in that particular tissue or tumor. Table 11 summarizes the RT-PCR data.

### VII. Summary of the Data Provided in the Tables

The level of expression of numerous potential markers (*i.e.* "the markers of the invention") in cells obtained from seven patients afflicted with ovarian cancer, and in cells of six ovarian cancer cell lines (*i.e.* a total of thirteen sample sources) were

compared with levels of expression of the same markers in non-cancerous ovarian cell samples. Markers for which significant differences in the levels of expression in cancer-related samples and non-cancerous samples were observed are listed in the Tables.

Tables 1 and 1A list markers that were identified in subtractive libraries and  
5 which are preferentially expressed in ovarian cancer cells over normal (*i.e.* non-cancerous) ovarian cells.

Table 2A lists markers, expression of which was increased by at least 5-fold in at least one of twenty-three ovarian cancer samples tested, relative to its expression in normal (*i.e.* non-cancerous) ovarian samples. Table 2B lists markers, expression of  
10 which was increased by at least 2-fold in all twenty-three ovarian cancer samples tested, relative to its expression in normal ovarian samples. Table 2C lists markers, expression of which was increased by at least 5-fold in at least 6 of the 23 ovarian cancer samples tested, relative to its expression in normal ovarian cells. Table 2D lists markers,  
15 expression of which was increased by at least 5-fold in at least 6 of the 23 ovarian cancer samples, relative to expression in normal ovarian samples, and which can serve as antigens for embodiments of the invention based upon proteomic studies, sequence analysis and/or literature references.

Table 3A lists markers, expression of which was decreased by at least 5-fold in at least one of twenty-three ovarian cancer samples tested, relative to its expression in  
20 normal (*i.e.* non-cancerous) ovarian cells. Table 3B lists markers, expression of which was decreased by at least 2-fold in all twenty-three ovarian cancer samples tested, relative to its expression in normal ovarian cells. Table 3C lists markers, expression of which was decreased by at least 5-fold in at least 6 of the 23 ovarian cancer samples tested, relative to its expression in normal (*i.e.* non-cancerous) ovarian cells.

25 Tables 4 and 5 list markers, expression of which was either increased (Table 4) or decreased (Table 5) in ovarian cancer samples, relative to expression in normal (*i.e.*, non-cancerous) ovarian samples. In particular, expression of the markers in 37 tumors (7 endometroid tumors, 5 clear cell tumors and 25 serous tumors) was evaluated. A ranking system based on the sum of the number of tumors multiplied by the fold  
30 regulation (for 2-fold, 3-fold, 5-fold and 10-fold regulation), divided by the total number of tumors, was employed. A rank score was generated for four categories, endometroid tumors, clear cell tumors, serous tumors and overall.

The markers of Table 4 had a score of greater than 1.5 for endometroid tumors, greater than 1.5 for clear cell tumors, greater than 1 for serous tumors, or greater than 0.8 overall. Table 4A shows the markers of Table 4 with a score of greater than 3 in any of the four categories.

- 5        The markers of Table 5 had a score of greater than 2.5 for endometroid tumors, greater than 2.5 for clear cell tumors, greater than 2 for serous tumors, or greater than 1.75 overall. Table 5A shows the markers of Table 5 with a score of greater than 3 in any of the four categories.

- 10       Tables 6 and 6A list markers that were identified in subtractive libraries and which are preferentially expressed in ovarian cancer cells over normal (i.e. non-cancerous ovarian cells).

Table 7A-7E show markers of the present invention obtained through proteomic analysis as described in Section VI., subsection C., above.

- 15       Table 8 lists the nucleotide sequences of 24 novel genes identified as described in Section VI., subsection D, above.

Table 9 depicts the results of the TaqMan® expression analysis obtained as described in Section VI., subsection F, above.

Tables 10A-10N contain Northern blot analysis data obtained as described in Section VI., subsection E, above.

- 20       Table 10A shows Marker 5 expression in normal human tissue samples. The highest level of expression is seen in placenta, followed by trachea, prostate, mammary gland, and lung, with lower levels in kidney, salivary gland, small intestine, and bladder, and an even lower level of expression in normal ovary tissue.

- 25       Table 10B shows that Marker 5 is expressed in several cancer cell lines. The highest level of expression is seen in SK-BR-3, followed by T-47D, BT-483, and ZR-75-30.

- 30       Table 10C shows Marker 8 expression in wide range of normal human tissue samples. The highest level of expression was seen in cerebellum, followed by placenta, prostate, and lung. Lower levels of expression were seen in kidney, spleen, testis, whole brain, and trachea, followed by mammary gland, small intestine, and thymus, which were higher than the level of expression in normal ovary tissue.

Table 10D shows that Marker 8 is expressed in all the cancer cell lines tested, with the highest levels of expression in DU4475, followed by MDA-MB-361line.

Table 10E shows Marker 10 expression was detected in all tissue samples tested. The highest level of expression was seen in trachea, followed by testis and prostate.  
5 Lower levels of expression were seen in whole brain, salivary gland, cerebellum, and small intestine. Expression in normal ovary tissue was among the lowest levels observed.

Table 10F shows that Marker 10 expression was detected in all cancer cell lines tested, with the highest levels of expression in MDA-MB-361, followed by MDA-MB-  
10 468, and HTB-78.

Table 10G shows a limited distribution of expression of Marker 15 in the panel of normal tissues tested, with significant expression only in placenta, and much lower levels of expression in whole brain, cerebellum, and prostate. No detectable levels of expression were seen in normal ovarian tissue.

15 Table 10H shows that Marker 15 expression was detected in all cancer cell lines tested, with the highest levels of expression seen in HTB-78, followed by MDA-MB-361, SK-BR-3, Caov-3, and MDA-MB-231.

Table 10I shows that expression of Marker 19 in the panel of normal human tissues tested was much higher in testis than in prostate and whole brain. Lower, but  
20 detectable, levels of expression were seen in a number of other tissues, with ovary among the lowest.

Table 10J shows that expression of Marker 19 was seen in 22 of the 26 cancer cell lines tested. The highest levels of expression were seen in BT549 and DU145, followed by NIH-Ovar-3 and HTB-78. Lower levels of expression were seen in MDA-  
25 MB-453, MDA-MB-361, and T-470.

Table 10K shows that high levels of expression of Marker 25 in the panel of normal human tested were seen in placenta, prostate, and trachea, followed by kidney, lung, and small intestine. Lower levels of expression were seen in salivary gland, spleen, thymus, and bladder. Expression in normal ovarian tissue was just above  
30 background.



Table 10L shows that expression of Marker 25 was detected in 20 of the 26 cancer cell lines tested. The highest level of expression was seen in T-470, followed by S-BR-3. Lower levels of expression were seen in Caov-3, MDA-MB-468, and HTB-78, followed by MDA-MB-453, MDA-MB-361, BT-483, DU4475, and NIH-Ovcar-3.

5 Table 10M shows that the highest level of expression of Marker 039 was seen in whole brain, followed by cerebellum, with a lower level in prostate. Even lower levels were seen in a number of tissues, including kidney, liver, spleen, testis, thymus, trachea, and lung. Expression in normal ovarian tissue was among the lowest.

Table 10N shows that Marker 39 expression was detected in most of the cancer  
10 cell lines tested, with the highest level seen in SK-BR-3, followed by MDA-MB-361 and T470. Lower levels of expression were seen in all other cell lines tested, except for MDA-MB-157, Hs578Bst, Hs578T, and ES-2, in which no expression was detected.

Table 11 depicts the results of LightCycler data and RT-PCR data obtained as described in Section VI., subsections G. and H., respectively, above.

15 Tables 1-1, 2A-1, 2D-1, 3A-1, 4-1, 5-1 and 6-1 depict the accession number ("ACC Num") and database ("DATABASE") of the markers of the present invention with the corresponding GenBank GI number ("GI NBR"). One skilled in the art may thus obtain from the Tables of the invention, both GenBank accession number as well as the GenBank GI number for a marker of the present invention, thereby identifying the  
20 nucleotide and/or polypeptide sequence of that marker. For example, the markers of Tables 1 and 1A are referenced in Table 1-1 by both GenBank accession number and GenBank GI number.

Those skilled in the art will readily understand the data set forth in the Tables of the present invention. In particular, the following definitions will be understood to  
25 mean:

- 1) "ID #" or "#" is an arbitrary designation assigned to the marker.
- 2) "Image Clone ID" is the identification number assigned to the marker by the IMAGE Consortium (Lennon *et al.*, 1996, *Genomics* 33:151-152; see, *e.g.*, "<http://www-bio.llnl.gov/bbrp/image/image.html>" for further information). All referenced Image  
30 Clone sequences are expressly incorporated by reference.
- 3) "GenBank Accession Number" or "Accession No." or "acc" or "Accession #" or "Acc Num" is the identification number assigned to the marker in the relevant database

(see, e.g. "[http://www.ncbi.nlm.nih.gov/genbank/query\\_form.html](http://www.ncbi.nlm.nih.gov/genbank/query_form.html)" and "[www.derwent.com](http://www.derwent.com)" for further information). "GenBank Gi" or "GI NB" is the GI identification number assigned to the marker in the GenBank database (see *supra*). All referenced database sequences are expressly incorporated herein by reference.

- 5 4) "Secreted?" or "Secreted" indicates whether the protein corresponding to the marker has been demonstrated to be secreted in protein profiling experiments.
- 5) "Secretion Predicted?" indicates whether the protein corresponding to the marker is predicted, using the SIGNALP computer software described herein, to have at least one portion which is exposed to the extracellular medium upon expression of the protein.
- 10 6) "Ave-Normal-Exp" indicates the average marker expression in the non-cancerous samples.
- 7) "Max expression" and "Min-expression" indicates the highest (or lowest) marker expression value of all samples.
- 8) "Max fold up" and "Max fold-up down" indicates the highest fold positive (or
- 15 negative) induction of regulation of the marker of all samples.
- 9) "Count-up tumors" and "Count-down tumors" indicate the total number of the twenty-three tumor samples that the marker was up (or down) regulated.
- 10) "Count-up cell lines" and "Count-down cell lines" indicated the total number of the six cell lines where the marker was up (or down) regulated.
- 20 11) "Chromosome" indicates the chromosome on which the genomic sequence corresponding to the marker is located, where this location is known.
- 12) "Location" indicates the location on the chromosome at which the genomic sequence corresponding to the marker is located, where this location is known. The genes were mapped using radiation hybrid panel data that can be found in the art, for
- 25 example at "<http://www.sanger.ac.uk/HGP/Rhmap/>" and at "<http://www.ncbi.nlm.nih.gov/genemap99/>".
- 13) "Tissue Prominence" indicates up to three tissues in which expression of the marker is predicted, based on expression in the predicted tissues, of expressed sequence tags located in close proximity to the marker. The marker may also, or instead, be expressed
- 30 in tissues that are not listed in this section (*i.e.* this list is not exhaustive).

- 110 -

14) "Database" or "dbase" refers to the relevant database where the nucleotide sequence may be found according to its accession number. These public databases include GenBank, dbEST (a division of GenBank), and NUCPATENT (a GENESEQ database, available through Derwent). For examples, see

- 5 <http://www.ncbi.nlm.nih.gov/Entrez/nucleotide.html> for GenBank and  
[www.derwent.com](http://www.derwent.com) for GENESEQ. All referenced database sequences are expressly  
incorporated herein by reference.

- The contents of all references, patents, published patent applications, and  
database records including, GenBank, IMAGE consortium and GENESEQ database  
10 records, cited throughout this application are hereby incorporated by reference.

#### Other Embodiments

- Those skilled in the art will recognize, or be able to ascertain using no more than  
routine experimentation, many equivalents to the specific embodiments of the invention  
15 described herein. Such equivalents are intended to be encompassed by the following  
claims.

What is claimed is:

- 111 -

**Claims**

1. A method of assessing whether a patient is afflicted with ovarian cancer, the method comprising comparing:
  - 5 a) the level of expression of a marker in a patient sample, wherein the marker is selected from the group consisting of the markers listed in Tables 1-11, and
  - b) the normal level of expression of the marker in a control non-ovarian cancer sample,wherein a significant difference between the level of expression of the marker in  
10 the patient sample and the normal level is an indication that the patient is afflicted with ovarian cancer.
2. The method of claim 1, wherein the marker is selected from the group consisting of the markers listed in Table 2C.  
15
3. The method of claim 1, wherein the marker is selected from the group consisting of the markers listed in Table 2D.
4. The method of claim 1, wherein the marker is selected from the group  
20 consisting of the markers listed in Table 3C.
5. The method of claim 1, wherein the marker is selected from the group consisting of the markers listed in Table 4A.
- 25 6. The method of claim 1, wherein the marker is selected from the group consisting of the markers listed in Table 5A.
7. The method of claim 1, wherein the marker is selected from the group consisting of the markers listed in Tables 6 and 6A.  
30
8. The method of claim 1, wherein the marker is selected from the group consisting of the markers listed in Table 8.

- 112 -

9. The method of claim 1, wherein the marker is selected from the group consisting of the markers listed in Tables 7A-7E.

10. The method of claim 1, wherein the marker corresponds to a secreted  
5 protein.

11. The method of claim 10, wherein the marker is selected from the group consisting of the markers listed in Tables 7A-7E.

10 12. The method of claim 1, wherein the marker corresponds to a transcribed polynucleotide or portion thereof, wherein the polynucleotide comprises the marker.

13. The method of claim 1, wherein at least one tissue corresponding to the marker in the Tables is an epithelial tissue.

15

14. The method of claim 13, wherein at least one tissue corresponding to the marker in the Tables is an ovarian tissue.

15. The method of claim 1, wherein the marker is over- or under-expressed by at  
20 least two-fold in at least about 20% of ovarian cancer patients.

16. The method of claim 1, wherein the marker is not significantly expressed in non-ovarian tissues.

25 17. The method of claim 1, wherein the patient sample is an ovary-associated body fluid.

18. The method of claim 13, wherein the ovary-associated body fluid is selected from the group consisting of blood fluid, lymph, ascitic fluid, gynecological fluid, cystic  
30 fluid, urine, and a fluid collected by peritoneal rinsing.

- 113 -

19. The method of claim 1, wherein the sample comprises cells obtained from the patient.

20. The method of claim 19, wherein the cells are in a fluid selected from the group consisting of a fluid collected by peritoneal rinsing, a fluid collected by uterine rinsing, a uterine fluid, a uterine exudate, a pleural fluid, a cystic fluid, and an ovarian exudate.

21. The method of claim 1, wherein the level of expression of the marker in the sample is assessed by detecting the presence in the sample of a protein corresponding to the marker.

22. The method of claim 21, wherein the marker is selected from the group consisting of the markers listed in Tables 7A-7E and 8.

15

23. The method of claim 21, wherein the presence of the protein is detected using a reagent which specifically binds with the protein.

24. The method of claim 23, wherein the reagent is selected from the group consisting of an antibody, an antibody derivative, and an antibody fragment.

20

25. The method of claim 1, wherein the level of expression of the marker in the sample is assessed by detecting the presence in the sample of a transcribed polynucleotide or portion thereof, wherein the transcribed polynucleotide comprises the marker.

25

26. The method of claim 25, wherein the transcribed polynucleotide is an mRNA.

27. The method of claim 25, wherein the transcribed polynucleotide is a cDNA.

30

- 114 -

28. The method of claim 25, wherein the step of detecting further comprises amplifying the transcribed polynucleotide.

29. The method of claim 1, wherein the level of expression of the marker in the  
5 sample is assessed by detecting the presence in the sample of a transcribed polynucleotide which anneals with the marker or anneals with a portion of a polynucleotide wherein the polynucleotide comprises the marker, under stringent hybridization conditions.

10 30. The method of claim 1, wherein the level of expression of the marker in the sample differs from the normal level of expression of the marker in a patient not afflicted with ovarian cancer by a factor of at least about 2.

15 31. The method of claim 1, wherein the level of expression of the marker in the sample differs from the normal level of expression of the marker in a patient not afflicted with ovarian cancer by a factor of at least about 5.

32. The method of claim 1, comprising comparing:  
a) the level of expression in the sample of each of a plurality of markers  
20 independently selected from the markers listed in Tables 1-11, and  
b) the normal level of expression of each of the plurality of markers in samples of the same type obtained from control humans not afflicted with ovarian cancer, wherein the level of expression of more than one of the markers is significantly altered, relative to the corresponding normal levels of expression of the markers, is an  
25 indication that the patient is afflicted with ovarian cancer.

33. The method of claim 32, wherein the plurality comprises at least three of the markers.

30 34. The method of claim 32, wherein the plurality comprises at least five of the markers.

35. A method of assessing whether a patient is afflicted with ovarian cancer, the method comprising comparing:

5 a) the level of expression of a marker in a sample obtained from the patient, wherein the marker is selected from the group consisting of the markers listed in Tables 1-11 and

b) the normal level of expression of the marker in samples of the same type obtained from control humans not afflicted with ovarian cancer, wherein a significantly different level of expression of the marker in the sample, relative to the normal level, is an indication that the patient is afflicted with ovarian  
10 cancer.

36. A method for monitoring the progression of ovarian cancer in a patient, the method comprising:

15 a) detecting in a patient sample at a first point in time, the expression of a marker, wherein the marker is selected from the group consisting of the markers listed in Tables 1-11;

b) repeating step a) at a subsequent point in time; and

c) comparing the level of expression detected in steps a) and b), and therefrom monitoring the progression of ovarian cancer in the patient.

20

37. The method of claim 36, wherein the marker is selected from the group consisting of the markers listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8.

38. The method of claim 36, wherein the marker is selected from the group  
25 consisting of the markers listed in Tables 3A, 5, 7C and 7E.

39. The method of claim 36, wherein the marker corresponds to a secreted protein.

30 40. The method of claim 36, wherein marker corresponds to a transcribed polynucleotide or portion thereof, wherein the polynucleotide comprises the marker.



- 116 -

41. The method of claim 36, wherein the patient sample is an ovary-associated body fluid.

42. The method of claim 36, wherein the sample comprises cells obtained from  
5 the patient.

43. The method of claim 36, wherein between the first point in time and the subsequent point in time, the patient has undergone surgery to remove a tumor.

10 44. A method of assessing the efficacy of a test compound for inhibiting an ovarian cancer in a patient, the method comprising comparing:

a) expression of a marker in a first sample obtained from the patient and maintained in the presence of the test compound, wherein the marker is selected from the group consisting of the markers listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and  
15 8, and

b) expression of the marker in a second sample obtained from the patient and maintained in the absence of the test compound,  
wherein a significantly lower level of expression of the marker in the first sample, relative to the second sample, is an indication that the test compound is  
20 efficacious for inhibiting ovarian cancer in the patient.

45. The method of claim 44, wherein the first and second samples are portions of a single sample obtained from the patient.

25 46. The method of claim 44, wherein the first and second samples are portions of pooled samples obtained from the patient.

47. A method of assessing the efficacy of a test compound for inhibiting ovarian cancer in a patient, the method comprising comparing:

a) expression of a marker in a first sample obtained from the patient and maintained in the presence of the test compound, wherein the marker is selected from  
5 the group consisting of the markers listed in Tables 3A, 5, 7C and 7E, and

b) expression of the marker in a second sample obtained from the patient and maintained in the absence of the test compound,

wherein a significantly enhanced level of expression of the marker in the first sample, relative to the second sample, is an indication that the test compound is  
10 efficacious for inhibiting ovarian cancer in the patient.

48. A method of assessing the efficacy of a therapy for inhibiting ovarian cancer in a patient, the method comprising comparing:

a) expression of a marker in the first sample obtained from the patient prior to  
15 providing at least a portion of the therapy to the patient, wherein the marker is selected from the group consisting of the markers listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8, and

b) expression of the marker in a second sample obtained from the patient following provision of the portion of the therapy,

20 wherein a significantly lower level of expression of the marker in the second sample, relative to the first sample, is an indication that the therapy is efficacious for inhibiting ovarian cancer in the patient.

49. A method of assessing the efficacy of a therapy for inhibiting ovarian cancer  
25 in a patient, the method comprising comparing:

a) expression of a marker in the first sample obtained from the patient prior to providing at least a portion of the therapy to the patient, wherein the marker is selected from the group consisting of the markers listed in Tables 3A, 5, 7C and 7E, and

- 118 -

b) expression of the marker in a second sample obtained from the patient following provision of the portion of the therapy,

wherein a significantly enhanced level of expression of the marker in the second sample, relative to the first sample, is an indication that the therapy is efficacious for  
5 inhibiting ovarian cancer in the patient.

50. A method of selecting a composition for inhibiting ovarian cancer in a patient, the method comprising:

- a) obtaining a sample comprising cancer cells from the patient;
- 10 b) separately maintaining aliquots of the sample in the presence of a plurality of test compositions;
- c) comparing expression of a marker in each of the aliquots, wherein the marker is selected from the group consisting of the markers listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8; and
- 15 d) selecting one of the test compositions which induces a lower level of expression of the marker in the aliquot containing that test composition, relative to other test compositions.

51. A method of selecting a composition for inhibiting ovarian cancer in a  
20 patient, the method comprising:

- a) obtaining a sample comprising cancer cells from the patient;
- b) separately maintaining aliquots of the sample in the presence of a plurality of test compositions;
- c) comparing expression of a marker in each of the aliquots, wherein the marker  
25 is selected from the group consisting of the markers listed in Tables 3A, 5, 7C and 7E; and
- d) selecting one of the test compositions which induces an enhanced level of expression of the marker in the aliquot containing that test composition, relative to other test compositions.

30

52. A method of inhibiting ovarian cancer in a patient, the method comprising:

- a) obtaining a sample comprising cancer cells from the patient;
- b) separately maintaining aliquots of the sample in the presence of a plurality of test compositions;
- 5 c) comparing expression of a marker in each of the aliquots, wherein the marker is selected from the group consisting of the markers listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8, and
- d) administering to the patient at least one of the test compositions which induces a lower level of expression of the marker in the aliquot containing that test composition,
- 10 relative to other test compositions.

53. A method of selecting a composition for inhibiting ovarian cancer in a patient, the method comprising:

- a) obtaining a sample comprising cancer cells from the patient;
- 15 b) separately maintaining aliquots of the sample in the presence of a plurality of test compositions;
- c) comparing expression of a marker in each of the aliquots, wherein the marker is selected from the group consisting of the markers listed in Tables 3A, 5, 7C and 7E, and
- 20 d) administering to the patient at least one of the test compositions which induces an enhanced level of expression of the marker in the aliquot containing that test composition, relative to other test compositions.

54. A kit for assessing the suitability of each of a plurality of compounds for inhibiting ovarian cancer in a patient, the kit comprising:

- a) the plurality of compounds; and
- b) a reagent for assessing expression of a marker selected from the group consisting of the markers listed in Tables 1-11.

55. A kit for assessing whether a patient is afflicted with ovarian cancer, the kit comprising reagents for assessing expression of a marker selected from the group consisting of the markers listed in Tables 1-11.

56. A method of making an isolated hybridoma which produces an antibody useful for assessing whether a patient is afflicted with ovarian cancer, the method comprising:

- 5 isolating a protein corresponding to a marker selected from the group consisting of the markers listed in Tables 1-11;  
immunizing a mammal using the isolated protein;  
isolating splenocytes from the immunized mammal;  
fusing the isolated splenocytes with an immortalized cell line to form  
10 hybridomas; and  
screening individual hybridomas for production of an antibody which specifically binds with the protein to isolate the hybridoma.

57. The method of claim 56, wherein the marker is selected from the group  
15 consisting of the members listed in Tables 7A-7E and 8.

58. An antibody produced by a hybridoma made by the method of claim 56.

59. A kit for assessing the presence of human ovarian cancer cells, the kit  
20 comprising an antibody, wherein the antibody specifically binds with a protein corresponding to a marker selected from the group consisting of the markers listed in Tables 1-11.

60. A kit for assessing the presence of ovarian cancer cells, the kit comprising a  
25 nucleic acid probe wherein the probe specifically binds with a transcribed polynucleotide corresponding to a marker selected from the group consisting of the markers listed in Tables 1-11.

- 121 -

61. A method of assessing the ovarian cell carcinogenic potential of a test compound, the method comprising:

a) maintaining separate aliquots of ovarian cells in the presence and absence of the test compound; and

5        b) comparing expression of a marker in each of the aliquots, wherein the marker is selected from the group consisting of the markers listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8, and

      wherein a significantly enhanced level of expression of the marker in the aliquot maintained in the presence of the test compound, relative to the aliquot maintained in the  
10       absence of the test compound, is an indication that the test compound possesses human ovarian cell carcinogenic potential.

62. A method of assessing the ovarian cell carcinogenic potential of a test compound, the method comprising:

15       a) maintaining separate aliquots of ovarian cells in the presence and absence of the test compound; and

      b) comparing expression of a marker in each of the aliquots, wherein the marker is selected from the group consisting of the markers listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8, and

20       wherein a significantly lower level of expression of the marker in the aliquot maintained in the presence of the test compound, relative to the aliquot maintained in the absence of the test compound, is an indication that the test compound possesses ovarian cell carcinogenic potential.

25       63. A kit for assessing the ovarian cell carcinogenic potential of a test compound, the kit comprising ovarian cells and a reagent for assessing expression of a marker, wherein the marker is selected from the group consisting of the markers listed in Tables 1-11.

30       64. A method of treating a patient afflicted with ovarian cancer, the method comprising providing to cells of the cancer a protein corresponding to a marker selected from the markers listed in Tables 3A, 5, 7C and 7E.

- 122 -

65. The method of claim 62, wherein the protein is provided to the cells by providing a vector comprising a polynucleotide encoding the protein to the cells.

66. A method of treating a patient afflicted with ovarian cancer, the method comprising providing to cells of the patient an antisense oligonucleotide complementary to a polynucleotide corresponding to a marker selected from the markers listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8.

67. A method of inhibiting ovarian cancer in a patient at risk for developing ovarian cancer, the method comprising inhibiting expression of a gene corresponding to a marker selected from the markers listed in Tables 1, 1A, 2A, 4, 6, 6A, 7A, 7B, 7D and 8.

68. A method of inhibiting ovarian cancer in a patient at risk for developing ovarian cancer, the method comprising enhancing expression of a gene corresponding to a marker selected from the markers listed in Tables 3A, 5, 7C and 7E.

69. An isolated nucleic acid molecule selected from the group consisting of:  
a) a nucleic acid molecule comprising a nucleotide sequence which is at least 90% homologous to a nucleotide sequence of Table 8, or a complement thereof;  
b) a nucleic acid molecule comprising a fragment of a nucleic acid molecule comprising the nucleotide sequence of Table 8, or a complement thereof; and  
c) a nucleic acid molecule comprising the nucleotide sequence of Table 8, or a complement thereof.

25

70. A vector which contains a nucleic acid molecule of claim 69.

71. A host cell which contains a nucleic acid molecule of claim 69.

- 123 -

72. An isolated polypeptide which is encoded by a nucleic acid molecule comprising a nucleotide sequence which is at least 90% homologous to a nucleic acid comprising the nucleotide sequence of Table 8.

5        73. An antibody which selectively binds to a polypeptide of claim 72.



Table 1

Patent Sequence Number	Accession Number		
Sequence 1	AA001066	Sequence 50	AA111907
Sequence 2	AA007157	Sequence 51	AA112043
Sequence 3	AA010954	Sequence 52	AA112308
Sequence 4	AA015792	Sequence 53	AA112375
Sequence 5	AA019769	Sequence 54	AA113860
Sequence 6	AA019948	Sequence 55	AA114120
Sequence 7	AA022925	Sequence 56	AA115118
Sequence 8	AA022937	Sequence 57	AA115368
Sequence 9	AA024405	Sequence 58	AA122286
Sequence 10	AA029750	Sequence 59	AA122348
Sequence 11	AA031509	Sequence 60	AA126109
Sequence 12	AA033876	Sequence 61	AA127105
Sequence 13	AA034237	Sequence 62	AA127132
Sequence 14	AA039967	Sequence 63	AA127418
Sequence 15	AA040073	Sequence 64	AA128305
Sequence 16	AA040122	Sequence 65	AA129461
Sequence 17	AA045732	Sequence 66	AA130252
Sequence 18	AA045861	Sequence 67	AA130547
Sequence 19	AA046835	Sequence 68	AA130786
Sequence 20	AA047026	Sequence 69	AA131041
Sequence 21	AA047417	Sequence 70	AA131065
Sequence 22	AA053486	Sequence 71	AA131104
Sequence 23	AA054658	Sequence 72	AA131155
Sequence 24	AA055606	Sequence 73	AA131160
Sequence 25	AA056113	Sequence 74	AA132182
Sequence 26	AA056176	Sequence 75	AA132568
Sequence 27	AA056363	Sequence 76	AA132598
Sequence 28	AA056431	Sequence 77	AA133351
Sequence 29	AA065336	Sequence 78	AA133927
Sequence 30	AA069781	Sequence 79	AA134105
Sequence 31	AA069784	Sequence 80	AA134210
Sequence 32	AA069839	Sequence 81	AA135032
Sequence 33	AA069983	Sequence 82	AA135919
Sequence 34	AA071255	Sequence 83	AA136383
Sequence 35	AA075135	Sequence 84	AA136789
Sequence 36	AA081655	Sequence 85	AA143609
Sequence 37	AA082245	Sequence 86	AA146773
Sequence 38	AA083471	Sequence 87	AA147806
Sequence 39	AA083510	Sequence 88	AA148160
Sequence 40	AA085862	Sequence 89	AA148268
Sequence 41	AA085872	Sequence 90	AA148771
Sequence 42*	AA085947	Sequence 91	AA149056
Sequence 43	AA088770	Sequence 92	AA150307
Sequence 44	AA100333	Sequence 93	AA151310
Sequence 45	AA100719	Sequence 94	AA151775
Sequence 46	AA100793	Sequence 95	AA152037
Sequence 47	AA100852	Sequence 96	AA152416
Sequence 48	AA101270	Sequence 97	AA155853
Sequence 49	AA101561	Sequence 98	AA155926
		Sequence 99	AA157405
		Sequence 100	AA157725
		Sequence 101	AA157788

**Table 1**

Sequence 102	AA158165	Sequence 154	AA301631
Sequence 103	AA158171	Sequence 155	AA304669
Sequence 104	AA159272	Sequence 156	AA304961
Sequence 105	AA160114	Sequence 157	AA305193
Sequence 106	AA160685	Sequence 158	AA305438
Sequence 107	AA161410	Sequence 159	AA306542
Sequence 108	AA164405	Sequence 160	AA306708
Sequence 109	AA164465	Sequence 161	AA306945
Sequence 110	AA165083	Sequence 162	AA307239
Sequence 111	AA165629	Sequence 163	AA307477
Sequence 112	AA166973	Sequence 164	AA307504
Sequence 113	AA171510	Sequence 165	AA307697
Sequence 114	AA173031	Sequence 166	AA307779
Sequence 115	AA173470	Sequence 167	AA308062
Sequence 116	AA173630	Sequence 168	AA308801
Sequence 117	AA179462	Sequence 169	AA309028
Sequence 118	AA187003	Sequence 170	AA309988
Sequence 119	AA187958	Sequence 171	AA311006
Sequence 120	AA188591	Sequence 172	AA311481
Sequence 121	AA192108	Sequence 173	AA312012
Sequence 122	AA199710	Sequence 174	AA313684
Sequence 123	AA203224	Sequence 175	AA314146
Sequence 124	AA203284	Sequence 176	AA315049
Sequence 125	AA205851	Sequence 177	AA315308
Sequence 126	AA209431	Sequence 178	AA315426
Sequence 127	AA209531	Sequence 179	AA316682
Sequence 128	AA214075	Sequence 180	AA319958
Sequence 129	AA216612	Sequence 181	AA320346
Sequence 130	AA224230	Sequence 182	AA320991
Sequence 131	AA224985	Sequence 183	AA328544
Sequence 132	AA226502	Sequence 184	AA330457
Sequence 133	AA229225	Sequence 185	AA338793
Sequence 134	AA232626	Sequence 186	AA340069
Sequence 135	AA233843	Sequence 187	AA341170
Sequence 136	AA242891	Sequence 188	AA342394
Sequence 137	AA250725	Sequence 189	AA348250
Sequence 138	AA250982	Sequence 190	AA349148
Sequence 139	AA256959	Sequence 191	AA351443
Sequence 140	AA259077	Sequence 192	AA351880
Sequence 141	AA262440	Sequence 193	AA356158
Sequence 142	AA263110	Sequence 194	AA356187
Sequence 143	AA283165	Sequence 195	AA356195
Sequence 144	AA285260	Sequence 196	AA357374
Sequence 145	AA287112	Sequence 197	AA367446
Sequence 146	AA292191	Sequence 198	AA375236
Sequence 147	AA292334	Sequence 199	AA377718
Sequence 148	AA292385	Sequence 200	AA380997
Sequence 149	AA292771	Sequence 201	AA383917
Sequence 150	AA293273	Sequence 202	AA385147
Sequence 151	AA293572	Sequence 203	AA389641
Sequence 152	AA295348	Sequence 204	AA393164
Sequence 153	AA295485	Sequence 205	AA393236

Table 1

Sequence 206	AA394242	Sequence 258	AA568217
Sequence 207	AA398732	Sequence 259	AA573742
Sequence 208	AA401864	Sequence 260	AA573893
Sequence 209	AA410508	Sequence 261	AA574237
Sequence 210	AA410580	Sequence 262	AA576866
Sequence 211	AA410942	Sequence 263	AA579034
Sequence 212	AA411334	Sequence 264	AA579816
Sequence 213	AA411599	Sequence 265	AA581220
Sequence 214	AA418061	Sequence 266	AA581264
Sequence 215	AA418473	Sequence 267	AA582093
Sequence 216	AA418970	Sequence 268	AA583091
Sequence 217	AA420789	Sequence 269	AA584411
Sequence 218	AA421682	Sequence 270	AA586776
Sequence 219	AA421850	Sequence 271	AA587110
Sequence 220	AA424529	Sequence 272	AA587233
Sequence 221	AA428421	Sequence 273	AA587700
Sequence 222	AA429754	Sequence 274	AA609259
Sequence 223	AA441787	Sequence 275	AA609837
Sequence 224	AA451633	Sequence 276	AA613907
Sequence 225	AA453309	Sequence 277	AA614529
Sequence 226	AA453559	Sequence 278	AA618033
Sequence 227	AA453570	Sequence 279	AA628487
Sequence 228	AA454871	Sequence 280	AA631204
Sequence 229	AA454913	Sequence 281	AA631811
Sequence 230	AA456892	Sequence 282	AA640901
Sequence 231	AA457048	Sequence 283	AA641841
Sequence 232	AA463426	Sequence 284	AA642215
Sequence 233	AA465039	Sequence 285	AA643602
Sequence 234	AA477173	Sequence 286	AA651720
Sequence 235	AA480921	Sequence 287	AA664996
Sequence 236	AA484050	Sequence 288	AA668297
Sequence 237	AA484756	Sequence 289	AA668836
Sequence 238	AA487483	Sequence 290	AA675923
Sequence 239	AA489640	Sequence 291	AA687833
Sequence 240	AA493886	Sequence 292	AA704992
Sequence 241	AA494493	Sequence 293	AA732702
Sequence 242	AA496518	Sequence 294	AA745241
Sequence 243	AA501749	Sequence 295	AA746481
Sequence 244	AA501822	Sequence 296	AA758889
Sequence 245	AA501945	Sequence 297	AA772570
Sequence 246	AA504490	Sequence 298	AA772790
Sequence 247	AA507234	Sequence 299	AA776709
Sequence 248	AA513640	Sequence 300	AA776811
Sequence 249	AA526227	Sequence 301	AA777384
Sequence 250	AA526889	Sequence 302	AA778116
Sequence 251	AA527139	Sequence 303	AA779868
Sequence 252	AA527188	Sequence 304	AA781343
Sequence 253	AA531428	Sequence 305	AA809984
Sequence 254	AA532633	Sequence 306	AA810945
Sequence 255	AA535471	Sequence 307	AA811200
Sequence 256	AA554757	Sequence 308	AA825768
Sequence 257	AA565996	Sequence 309	AA828073

Table 1

Sequence 310	AA828722	Sequence 362	AF035316
Sequence 311	AA843176	Sequence 363	AF037204
Sequence 312	AA843661	Sequence 364	AF038451
Sequence 313	AA876526	Sequence 365	AF038662
Sequence 314	AA883255	Sequence 366	AF038963
Sequence 315	AA906652	Sequence 367	AF043431
Sequence 316	AA917638	Sequence 368	AF044956
Sequence 317	AA927734	Sequence 369	AF045941
Sequence 318	AA954939	Sequence 370	AF046997
Sequence 319	AA962622	Sequence 371	AF051894
Sequence 320	AA991285	Sequence 372	AF052124
Sequence 321	AB000115	Sequence 373	AF052578
Sequence 322	AB004047	Sequence 374	AF053233
Sequence 323	AB006746	Sequence 375	AF054838
Sequence 324	AB007619	Sequence 376	AF055012
Sequence 325	AB007860	Sequence 377	AF061736
Sequence 326	AB007965	Sequence 378	AF061738
Sequence 327	AB011101	Sequence 379	AF064603
Sequence 328	AB011169	Sequence 380	AF064854
Sequence 329	AB012701	Sequence 381	AF065388
Sequence 330	AB014536	Sequence 382	AF067168
Sequence 331	AB014565	Sequence 383	AF067174
Sequence 332	AB019568	Sequence 384	AF067817
Sequence 333	AB020623	Sequence 385	AF070523
Sequence 334	AB020629	Sequence 386	AF070561
Sequence 335	AB020693	Sequence 387	AF070562
Sequence 336	AB021288	Sequence 388	AF070596
Sequence 337	AB022663	Sequence 389	AF070664
Sequence 338	AB023214	Sequence 390	AF070674
Sequence 339	AB023230	Sequence 391	AF077048
Sequence 340	AC02059	Sequence 392	AF077051
Sequence 341	AC03653	Sequence 393	AF077200
Sequence 342	AC13415	Sequence 394	AF077671
Sequence 343	AF000982	Sequence 395	AF080246
Sequence 344	AF002985	Sequence 396	AF081484
Sequence 345	AF005654	Sequence 397	AF083470
Sequence 346	AF006086	Sequence 398	AF084523
Sequence 347	AF007791	Sequence 399	AF085355
Sequence 348	AF013758	Sequence 400	AF086003
Sequence 349	AF013988	Sequence 401	AF086080
Sequence 350	AF021232	Sequence 402	AF086183
Sequence 351	AF026939	Sequence 403	AF086545
Sequence 352	AF026941	Sequence 404	AF091263
Sequence 353	AF026942	Sequence 405	AF111713
Sequence 354	AF026943	Sequence 406	AF118023
Sequence 355	AF026944	Sequence 407	AF124438
Sequence 356	AF028832	Sequence 408	AF124439
Sequence 357	AF030455	Sequence 409	AF131808
Sequence 358	AF030514	Sequence 410	AF131820
Sequence 359	AF031469	Sequence 411	AF131848
Sequence 360	AF033095	Sequence 412	AF132966
Sequence 361	AF035286	Sequence 413	AF132968

Table 1

Sequence 414	AF146277	Sequence 466	AI382020
Sequence 415	AF147331	Sequence 467	AI400372
Sequence 416	AF150100	Sequence 468	AI417973
Sequence 417	AF150266	Sequence 469	AI431963
Sequence 418	AF151873	Sequence 470	AI453405
Sequence 419	AF151877	Sequence 471	AI457157
Sequence 420	AF151978	Sequence 472	AI457624
Sequence 421	AF167160	Sequence 473	AI459679
Sequence 422	AI023413	Sequence 474	AI460010
Sequence 423	AI027888	Sequence 475	AI469095
Sequence 424	AI031811	Sequence 476	AI469715
Sequence 425	AI033687	Sequence 477	AI471539
Sequence 426	AI042140	Sequence 478	AI476335
Sequence 427	AI075324	Sequence 479	AI479289
Sequence 428	AI075876	Sequence 480	AI499285
Sequence 429	AI126802	Sequence 481	AI521180
Sequence 430	AI127556	Sequence 482	AI538061
Sequence 431	AI129360	Sequence 483	AI567204
Sequence 432	AI139456	Sequence 484	AI587104
Sequence 433	AI140291	Sequence 485	AI587328
Sequence 434	AI144215	Sequence 486	AI609624
Sequence 435	AI161378	Sequence 487	AI610607
Sequence 436	AI188638	Sequence 488	AI612873
Sequence 437	AI215617	Sequence 489	AI627444
Sequence 438	AI216969	Sequence 490	AI632869
Sequence 439	AI241578	Sequence 491	AI633164
Sequence 440	AI250167	Sequence 492	AI636014
Sequence 441	AI253330	Sequence 493	AI637620
Sequence 442	AI253335	Sequence 494	AI676218
Sequence 443	AI253369	Sequence 495	AI683871
Sequence 444	AI253436	Sequence 496	AI684170
Sequence 445	AI261671	Sequence 497	AI693877
Sequence 446	AI262264	Sequence 498	AI694088
Sequence 447	AI267162	Sequence 499	AI732534
Sequence 448	AI267379	Sequence 500	AI743595
Sequence 449	AI267502	Sequence 501	AI744489
Sequence 450	AI267622	Sequence 502	AI745058
Sequence 451	AI279131	Sequence 503	AI753108
Sequence 452	AI285943	Sequence 504	AI791322
Sequence 453	AI289173	Sequence 505	AI798474
Sequence 454	AI290876	Sequence 506	AI803838
Sequence 455	AI292104	Sequence 507	AI811960
Sequence 456	AI300033	Sequence 508	AI813617
Sequence 457	AI300074	Sequence 509	AI815829
Sequence 458	AI312113	Sequence 510	AI826957
Sequence 459	AI336032	Sequence 511	AI831002
Sequence 460	AI337069	Sequence 512	AI863041
Sequence 461	AI340262	Sequence 513	AI867294
Sequence 462	AI346975	Sequence 514	AI912076
Sequence 463	AI354639	Sequence 515	AI915553
Sequence 464	AI366381	Sequence 516	AJ001381
Sequence 465	AI369024	Sequence 517	AJ003401

**Table 1**

Sequence 518	AJ010071	Sequence 570	D90311
Sequence 519	AJ132502	Sequence 571	D90453
Sequence 520	AL044356	Sequence 572	E01197
Sequence 521	AL044825	Sequence 573	E01198
Sequence 522	AL047024	Sequence 574	E01630
Sequence 523	AL048393	Sequence 575	E01954
Sequence 524	AL049313	Sequence 576	E01971
Sequence 525	AL049923	Sequence 577	E01972
Sequence 526	AL049954	Sequence 578	E02628
Sequence 527	AL050024	Sequence 579	E03569
Sequence 528	AL050272	Sequence 580	E03879
Sequence 529	AL050395	Sequence 581	E08663
Sequence 530	AL096714	Sequence 582	F06593
Sequence 531	AL096748	Sequence 583	F28779
Sequence 532	AL096842	Sequence 584	H25806
Sequence 533	AL110124	Sequence 585	H47546
Sequence 534	C17346	Sequence 586	H48873
Sequence 535	D00017	Sequence 587	H66467
Sequence 536	D00068	Sequence 588	H88415
Sequence 537	D11960	Sequence 589	J00196
Sequence 538	D12502	Sequence 590	J03575
Sequence 539	D12763	Sequence 591	J03858
Sequence 540	D13380	Sequence 592	J03909
Sequence 541	D13645	Sequence 593	J04164
Sequence 542	D13866	Sequence 594	K00422
Sequence 543	D14697	Sequence 595	K01763
Sequence 544	D21260	Sequence 596	L00693
Sequence 545	D23660	Sequence 597	L02426
Sequence 546	D26155	Sequence 598	L06328
Sequence 547	D26599	Sequence 599	L09159
Sequence 548	D28759	Sequence 600	L10413
Sequence 549	D29640	Sequence 601	L11066
Sequence 550	D31763	Sequence 602	L20688
Sequence 551	D31767	Sequence 603	L20941
Sequence 552	D31883	Sequence 604	L28997
Sequence 553	D38524	Sequence 605	L38995
Sequence 554	D42040	Sequence 606	L41490
Sequence 555	D45248	Sequence 607	M10119
Sequence 556	D49396	Sequence 608	M13536
Sequence 557	D50372	Sequence 609	M14328
Sequence 558	D50420	Sequence 610	M14764
Sequence 559	D55653	Sequence 611	M15329
Sequence 560	D81522	Sequence 612	M16660
Sequence 561	D83077	Sequence 613	M17017
Sequence 562	D83767	Sequence 614	M18216
Sequence 563	D86958	Sequence 615	M19723
Sequence 564	D86979	Sequence 616	M22918
Sequence 565	D87666	Sequence 617	M23613
Sequence 566	D87667	Sequence 618	M24194
Sequence 567	D87735	Sequence 619	M24594
Sequence 568	D88532	Sequence 620	M26152
Sequence 569	D89053	Sequence 621	M29540

Table 1

Sequence 622	M29541	Sequence 674	U03886
Sequence 623	M29551	Sequence 675	U04313
Sequence 624	M33146	Sequence 676	U07550
Sequence 625	M34064	Sequence 677	U07857
Sequence 626	M34455	Sequence 678	U08815
Sequence 627	M35198	Sequence 679	U09559
Sequence 628	M36693	Sequence 680	U09847
Sequence 629	M37716	Sequence 681	U10439
Sequence 630	M55268	Sequence 682	U14966
Sequence 631	M55542	Sequence 683	U18321
Sequence 632	M55543	Sequence 684	U19878
Sequence 633	M57567	Sequence 685	U23942
Sequence 634	M60333	Sequence 686	U25789
Sequence 635	M61715	Sequence 687	U28249
Sequence 636	M62831	Sequence 688	U28964
Sequence 637	M63121	Sequence 689	U32500
Sequence 638	M63838	Sequence 690	U32944
Sequence 639	M68520	Sequence 691	U33760
Sequence 640	M77945	Sequence 692	U37230
Sequence 641	M80563	Sequence 693	U37518
Sequence 642	M81757	Sequence 694	U38292
Sequence 643	M83248	Sequence 695	U38784
Sequence 644	M83654	Sequence 696	U41371
Sequence 645	M86553	Sequence 697	U41515
Sequence 646	M87284	Sequence 698	U52513
Sequence 647	M87434	Sequence 699	U56255
Sequence 648	M87503	Sequence 700	U57847
Sequence 649	M92357	Sequence 701	U61083
Sequence 650	M96982	Sequence 702	U68758
Sequence 651	M97501	Sequence 703	U73524
Sequence 652	M97935	Sequence 704	U77085
Sequence 653	N36346	Sequence 705	U78722
Sequence 654	N51262	Sequence 706	U79751
Sequence 655	N57413	Sequence 707	U94586
Sequence 656	N78477	Sequence 708	V00572
Sequence 657	N92060	Sequence 709	V00594
Sequence 658	Q21065	Sequence 710	V04202
Sequence 659	Q94780	Sequence 711	V17906
Sequence 660	R13925	Sequence 712	V36078
Sequence 661	R51732	Sequence 713	V68140
Sequence 662	R56461	Sequence 714	V86134
Sequence 663	R66489	Sequence 715	W02908
Sequence 664	R75621	Sequence 716	W05711
Sequence 665	S45630	Sequence 717	W07308
Sequence 666	S70290	Sequence 718	W25547
Sequence 667	S75295	Sequence 719	W28837
Sequence 668	S76638	Sequence 720	W37272
Sequence 669	T34641	Sequence 721	W38644
Sequence 670	T50925	Sequence 722	W39262
Sequence 671	T52715	Sequence 723	W39498
Sequence 672	T54951	Sequence 724	W52254
Sequence 673	T70793	Sequence 725	W74319

Table 1

Sequence 726	W77987	Sequence 778	AA022965
Sequence 727	W80480	Sequence 779	AA024522
Sequence 728	X00637	Sequence 780	AA028164
Sequence 729	X01742	Sequence 781	AA035775
Sequence 730	X02530	Sequence 782	AA037294
Sequence 731	X02661	Sequence 783	AA039967
Sequence 732	X04316	Sequence 784	AA045637
Sequence 733	X04371	Sequence 785	AA046815
Sequence 734	X04470	Sequence 786	AA046853
Sequence 735	X05908	Sequence 787	AA047052
Sequence 736	X07819	Sequence 788	AA047213
Sequence 737	X13238	Sequence 789	AA057071
Sequence 738	X15674	Sequence 790	AA058933
Sequence 739	X15729	Sequence 791	AA064952
Sequence 740	X16354	Sequence 792	AA075089
Sequence 741	X16356	Sequence 793	AA076291
Sequence 742	X16455	Sequence 794	AA078508
Sequence 743	X17025	Sequence 795	AA080864
Sequence 744	X20432	Sequence 796	AA083345
Sequence 745	X30167	Sequence 797	AA083693
Sequence 746	X33937	Sequence 798	AA085497
Sequence 747	X35726	Sequence 799	AA086463
Sequence 748	X41105	Sequence 800	AA093935
Sequence 749	X51841	Sequence 801	AA100291
Sequence 750	X54941	Sequence 802	AA101207
Sequence 751	X56932	Sequence 803	AA102403
Sequence 752	X57351	Sequence 804	AA111856
Sequence 753	X59710	Sequence 805	AA115174
Sequence 754	X65614	Sequence 806	AA122134
Sequence 755	X67951	Sequence 807	AA122291
Sequence 756	X68060	Sequence 808	AA125780
Sequence 757	X68277	Sequence 809	AA127322
Sequence 758	X72790	Sequence 810	AA130432
Sequence 759	X76488	Sequence 811	AA131801
Sequence 760	X83544	Sequence 812	AA132445
Sequence 761	X85134	Sequence 813	AA134109
Sequence 762	X87949	Sequence 814	AA135924
Sequence 763	X93036	Sequence 815	AA136322
Sequence 764	X99699	Sequence 816	AA143034
Sequence 765	X99920	Sequence 817	AA150057
Sequence 766	Y09267	Sequence 818	AA151651
Sequence 767	Y13323	Sequence 819	AA156335
Sequence 768	Y17392	Sequence 820	AA157333
Sequence 769	Z12830	Sequence 821	AA158987
Sequence 770	Z36815	Sequence 822	AA165439
Sequence 771	Z47087	Sequence 823	AA165632
Sequence 772	Z48570	Sequence 824	AA166618
Sequence 773	Z71389	Sequence 825	AA172067
Sequence 774	AA002223	Sequence 826	AA173031
Sequence 775	AA018843	Sequence 827	AA178870
Sequence 776	AA021647	Sequence 828	AA181874
Sequence 777	AA022842	Sequence 829	AA195194



Table 1

Sequence 830	AA203206	Sequence 882	AA446099
Sequence 831	AA203289	Sequence 883	AA446403
Sequence 832	AA204768	Sequence 884	AA447735
Sequence 833	AA206621	Sequence 885	AA449054
Sequence 834	AA213914	Sequence 886	AA449205
Sequence 835	AA218919	Sequence 887	AA449520
Sequence 836	AA224050	Sequence 888	AA452273
Sequence 837	AA224244	Sequence 889	AA455007
Sequence 838	AA227596	Sequence 890	AA455104
Sequence 839	AA229018	Sequence 891	AA459527
Sequence 840	AA229161	Sequence 892	AA460226
Sequence 841	AA236445	Sequence 893	AA461287
Sequence 842	AA236680	Sequence 894	AA464526
Sequence 843	AA243537	Sequence 895	AA468398
Sequence 844	AA252436	Sequence 896	AA469135
Sequence 845	AA252869	Sequence 897	AA469453
Sequence 846	AA256330	Sequence 898	AA470690
Sequence 847	AA262700	Sequence 899	AA479427
Sequence 848	AA278358	Sequence 900	AA480336
Sequence 849	AA287076	Sequence 901	AA483454
Sequence 850	AA291551	Sequence 902	AA487669
Sequence 851	AA293273	Sequence 903	AA488423
Sequence 852	AA295982	Sequence 904	AA488635
Sequence 853	AA301675	Sequence 905	AA488843
Sequence 854	AA301722	Sequence 906	AA489772
Sequence 855	AA302964	Sequence 907	AA503972
Sequence 856	AA303199	Sequence 908	AA508506
Sequence 857	AA304927	Sequence 909	AA513550
Sequence 858	AA305042	Sequence 910	AA513783
Sequence 859	AA305635	Sequence 911	AA514989
Sequence 860	AA315030	Sequence 912	AA516400
Sequence 861	AA315943	Sequence 913	AA520993
Sequence 862	AA317144	Sequence 914	AA521110
Sequence 863	AA326060	Sequence 915	AA523639
Sequence 864	AA327358	Sequence 916	AA523697
Sequence 865	AA336387	Sequence 917	AA528106
Sequence 866	AA346413	Sequence 918	AA528190
Sequence 867	AA352580	Sequence 919	AA528226
Sequence 868	AA363162	Sequence 920	AA534830
Sequence 869	AA375754	Sequence 921	AA548722
Sequence 870	AA399230	Sequence 922	AA551236
Sequence 871	AA400249	Sequence 923	AA551243
Sequence 872	AA401629	Sequence 924	AA558778
Sequence 873	AA402885	Sequence 925	AA563834
Sequence 874	AA406401	Sequence 926	AA576432
Sequence 875	AA421682	Sequence 927	AA580069
Sequence 876	AA422057	Sequence 928	AA580294
Sequence 877	AA424445	Sequence 929	AA582588
Sequence 878	AA424901	Sequence 930	AA584304
Sequence 879	AA424984	Sequence 931	AA588772
Sequence 880	AA425182	Sequence 932	AA593075
Sequence 881	AA428607	Sequence 933	AA595585

**Table 1**

Sequence 934	AA601895	Sequence 986	AA897461
Sequence 935	AA628700	Sequence 987	AA902582
Sequence 936	AA630326	Sequence 988	AA902644
Sequence 937	AA630642	Sequence 989	AA909144
Sequence 938	AA631178	Sequence 990	AA913281
Sequence 939	AA631218	Sequence 991	AA916756
Sequence 940	AA633550	Sequence 992	AA922420
Sequence 941	AA634808	Sequence 993	AA927283
Sequence 942	AA639199	Sequence 994	AA933075
Sequence 943	AA639791	Sequence 995	AA935979
Sequence 944	AA644273	Sequence 996	AA937947
Sequence 945	AA648897	Sequence 997	AA948295
Sequence 946	AA664732	Sequence 998	AA969131
Sequence 947	AA677550	Sequence 999	AA971881
Sequence 948	AA687308	Sequence 1000	AA973019
Sequence 949	AA705002	Sequence 1001	AA988923
Sequence 950	AA706685	Sequence 1002	AA989465
Sequence 951	AA708266	Sequence 1003	AA994023
Sequence 952	AA713687	Sequence 1004	AB002310
Sequence 953	AA719618	Sequence 1005	AB002330
Sequence 954	AA719674	Sequence 1006	AB007944
Sequence 955	AA720572	Sequence 1007	AB012911
Sequence 956	AA721752	Sequence 1008	AB017019
Sequence 957	AA723612	Sequence 1009	AB018266
Sequence 958	AA730571	Sequence 1010	AB018305
Sequence 959	AA742282	Sequence 1011	AB018347
Sequence 960	AA748437	Sequence 1012	AB019568
Sequence 961	AA749187	Sequence 1013	AB023158
Sequence 962	AA761602	Sequence 1014	AB028976
Sequence 963	AA768355	Sequence 1015	AB029005
Sequence 964	AA769127	Sequence 1016	AC28164
Sequence 965	AA774030	Sequence 1017	AD001528
Sequence 966	AA774247	Sequence 1018	AF000231
Sequence 967	AA779631	Sequence 1019	AF006088
Sequence 968	AA808747	Sequence 1020	AF006516
Sequence 969	AA809854	Sequence 1021	AF012072
Sequence 970	AA810859	Sequence 1022	AF026947
Sequence 971	AA825673	Sequence 1023	AF028832
Sequence 972	AA825768	Sequence 1024	AF030424
Sequence 973	AA826517	Sequence 1025	AF031379
Sequence 974	AA827331	Sequence 1026	AF035287
Sequence 975	AA827764	Sequence 1027	AF035309
Sequence 976	AA829511	Sequence 1028	AF038197
Sequence 977	AA831603	Sequence 1029	AF038404
Sequence 978	AA836991	Sequence 1030	AF043431
Sequence 979	AA837254	Sequence 1031	AF044670
Sequence 980	AA846480	Sequence 1032	AF044958
Sequence 981	AA846840	Sequence 1033	AF047184
Sequence 982	AA853515	Sequence 1034	AF052164
Sequence 983	AA883212	Sequence 1035	AF052496
Sequence 984	AA886885	Sequence 1036	AF052578
Sequence 985	AA889485	Sequence 1037	AF054990

Table 1

Sequence 1038	AF059524	Sequence 1090	AI338977
Sequence 1039	AF070561	Sequence 1091	AI339946
Sequence 1040	AF070626	Sequence 1092	AI373032
Sequence 1041	AF070655	Sequence 1093	AI374954
Sequence 1042	AF070674	Sequence 1094	AI380539
Sequence 1043	AF075040	Sequence 1095	AI417583
Sequence 1044	AF077030	Sequence 1096	AI432644
Sequence 1045	AF078847	Sequence 1097	AI433157
Sequence 1046	AF080246	Sequence 1098	AI457792
Sequence 1047	AF081282	Sequence 1099	AI469112
Sequence 1048	AF081484	Sequence 1100	AI471114
Sequence 1049	AF084523	Sequence 1101	AI471534
Sequence 1050	AF086163	Sequence 1102	AI473927
Sequence 1051	AF095791	Sequence 1103	AI479305
Sequence 1052	AF100756	Sequence 1104	AI499243
Sequence 1053	AF107406	Sequence 1105	AI525796
Sequence 1054	AF119297	Sequence 1106	AI525843
Sequence 1055	AF131858	Sequence 1107	AI537677
Sequence 1056	AF132940	Sequence 1108	AI541029
Sequence 1057	AF151857	Sequence 1109	AI560129
Sequence 1058	AI028733	Sequence 1110	AI583108
Sequence 1059	AI031901	Sequence 1111	AI584068
Sequence 1060	AI033739	Sequence 1112	AI587208
Sequence 1061	AI040324	Sequence 1113	AI589867
Sequence 1062	AI051172	Sequence 1114	AI610676
Sequence 1063	AI076805	Sequence 1115	AI630362
Sequence 1064	AI087005	Sequence 1116	AI633006
Sequence 1065	AI089913	Sequence 1117	AI634443
Sequence 1066	AI092007	Sequence 1118	AI635096
Sequence 1067	AI127326	Sequence 1119	AI682105
Sequence 1068	AI147251	Sequence 1120	AI683338
Sequence 1069	AI148933	Sequence 1121	AI684800
Sequence 1070	AI149846	Sequence 1122	AI684991
Sequence 1071	AI167855	Sequence 1123	AI689369
Sequence 1072	AI183965	Sequence 1124	AI689617
Sequence 1073	AI189258	Sequence 1125	AI689883
Sequence 1074	AI220148	Sequence 1126	AI693745
Sequence 1075	AI224374	Sequence 1127	AI701001
Sequence 1076	AI240095	Sequence 1128	AI733038
Sequence 1077	AI246677	Sequence 1129	AI735638
Sequence 1078	AI248538	Sequence 1130	AI741506
Sequence 1079	AI266582	Sequence 1131	AI742722
Sequence 1080	AI268864	Sequence 1132	AI742738
Sequence 1081	AI270183	Sequence 1133	AI743552
Sequence 1082	AI271795	Sequence 1134	AI753784
Sequence 1083	AI273008	Sequence 1135	AI754296
Sequence 1084	AI273841	Sequence 1136	AI754652
Sequence 1085	AI274756	Sequence 1137	AI754732
Sequence 1086	AI275528	Sequence 1138	AI765975
Sequence 1087	AI283096	Sequence 1139	AI769970
Sequence 1088	AI298059	Sequence 1140	AI819225
Sequence 1089	AI335653	Sequence 1141	AI820563

Table 1

Sequence 1142	AI827818	Sequence 1194	D30655
Sequence 1143	AI828682	Sequence 1195	D50310
Sequence 1144	AI830067	Sequence 1196	D51497
Sequence 1145	AI861989	Sequence 1197	D53031
Sequence 1146	AI887129	Sequence 1198	D62116
Sequence 1147	AI887632	Sequence 1199	D63878
Sequence 1148	AI890281	Sequence 1200	D78611
Sequence 1149	AI924046	Sequence 1201	D82348
Sequence 1150	AI924096	Sequence 1202	D83032
Sequence 1151	AI924823	Sequence 1203	D85433
Sequence 1152	AI963471	Sequence 1204	D87437
Sequence 1153	AI963604	Sequence 1205	D87667
Sequence 1154	AI972556	Sequence 1206	D89092
Sequence 1155	AI979048	Sequence 1207	D90041
Sequence 1156	AI984656	Sequence 1208	E02628
Sequence 1157	AJ010442	Sequence 1209	E05732
Sequence 1158	AJ132694	Sequence 1210	F00551
Sequence 1159	AJ224442	Sequence 1211	H08920
Sequence 1160	AL036299	Sequence 1212	H25080
Sequence 1161	AL042979	Sequence 1213	H30306
Sequence 1162	AL047305	Sequence 1214	H44647
Sequence 1163	AL049247	Sequence 1215	H81376
Sequence 1164	AL049313	Sequence 1216	H93521
Sequence 1165	AL049381	Sequence 1217	H94496
Sequence 1166	AL049932	Sequence 1218	J03464
Sequence 1167	AL050041	Sequence 1219	J03799
Sequence 1168	AL050161	Sequence 1220	J04027
Sequence 1169	AL050265	Sequence 1221	J04177
Sequence 1170	AL050268	Sequence 1222	K01228
Sequence 1171	AL050367	Sequence 1223	K01566
Sequence 1172	AL079286	Sequence 1224	L07395
Sequence 1173	AL079312	Sequence 1225	L09159
Sequence 1174	AL079314	Sequence 1226	L11315
Sequence 1175	AL080113	Sequence 1227	L13806
Sequence 1176	AL110164	Sequence 1228	L15702
Sequence 1177	AL117412	Sequence 1229	L16510
Sequence 1178	AL117612	Sequence 1230	L24804
Sequence 1179	AL119009	Sequence 1231	L25931
Sequence 1180	AW014693	Sequence 1232	L28809
Sequence 1181	AW014985	Sequence 1233	M10036
Sequence 1182	AW021794	Sequence 1234	M10905
Sequence 1183	C01521	Sequence 1235	M11353
Sequence 1184	D01096	Sequence 1236	M12267
Sequence 1185	D13119	Sequence 1237	M13536
Sequence 1186	D13627	Sequence 1238	M14483
Sequence 1187	D13630	Sequence 1239	M14630
Sequence 1188	D13639	Sequence 1240	M17885
Sequence 1189	D13665	Sequence 1241	M18366
Sequence 1190	D14530	Sequence 1242	M21575
Sequence 1191	D21260	Sequence 1243	M23254
Sequence 1192	D25278	Sequence 1244	M24194
Sequence 1193	D26361	Sequence 1245	M24486

Table 1

Sequence 1246	M26512	Sequence 1298	V24305
Sequence 1247	M28372	Sequence 1299	V81394
Sequence 1248	M31159	Sequence 1300	V84510
Sequence 1249	M32220	Sequence 1301	W19427
Sequence 1250	M36341	Sequence 1302	W65357
Sequence 1251	M36693	Sequence 1303	W75963
Sequence 1252	M38690	Sequence 1304	W80525
Sequence 1253	M58485	Sequence 1305	X01630
Sequence 1254	M59849	Sequence 1306	X04098
Sequence 1255	M62831	Sequence 1307	X04408
Sequence 1256	M64241	Sequence 1308	X06700
Sequence 1257	M69043	Sequence 1309	X14420
Sequence 1258	M77142	Sequence 1310	X51742
Sequence 1259	M77830	Sequence 1311	X60111
Sequence 1260	M86667	Sequence 1312	X69398
Sequence 1261	M88108	Sequence 1313	X72755
Sequence 1262	M93651	Sequence 1314	X74979
Sequence 1263	M95542	Sequence 1315	X76180
Sequence 1264	N43970	Sequence 1316	X78627
Sequence 1265	Q12759	Sequence 1317	X79067
Sequence 1266	Q14635	Sequence 1318	X80910
Sequence 1267	R11045	Sequence 1319	X87949
Sequence 1268	R76376	Sequence 1320	Y00052
Sequence 1269	R84450	Sequence 1321	Y00062
Sequence 1270	S74728	Sequence 1322	Y00282
Sequence 1271	S82081	Sequence 1323	Y00503
Sequence 1272	T07459	Sequence 1324	Y15286
Sequence 1273	T19883	Sequence 1325	Y17171
Sequence 1274	T21168	Sequence 1326	Z13009
Sequence 1275	T22605	Sequence 1327	Z24724
Sequence 1276	T37405	Sequence 1328	Z29083
Sequence 1277	T67129	Sequence 1329	Z29331
Sequence 1278	T69703	Sequence 1330	Z46606
Sequence 1279	T78615	Sequence 1331	Z48501
Sequence 1280	T89937	Sequence 1332	AA001460
Sequence 1281	U03851	Sequence 1333	AA001543
Sequence 1282	U12404	Sequence 1334	AA001792
Sequence 1283	U14967	Sequence 1335	AA004925
Sequence 1284	U14971	Sequence 1336	AA010897
Sequence 1285	U20659	Sequence 1337	AA017162
Sequence 1286	U25789	Sequence 1338	AA019019
Sequence 1287	U30825	Sequence 1339	AA022980
Sequence 1288	U47077	Sequence 1340	AA024595
Sequence 1289	U49844	Sequence 1341	AA024940
Sequence 1290	U63846	Sequence 1342	AA024996
Sequence 1291	U65928	Sequence 1343	AA025750
Sequence 1292	U72516	Sequence 1344	AA026598
Sequence 1293	U79282	Sequence 1345	AA029271
Sequence 1294	U90716	Sequence 1346	AA029725
Sequence 1295	U90904	Sequence 1347	AA029930
Sequence 1296	U94364	Sequence 1348	AA033832
Sequence 1297	V20437	Sequence 1349	AA035471

Table 1

Sequence 1350	AA035616	Sequence 1402	AA088344
Sequence 1351	AA036752	Sequence 1403	AA088351
Sequence 1352	AA037377	Sequence 1404	AA088693
Sequence 1353	AA039778	Sequence 1405	AA088783
Sequence 1354	AA039948	Sequence 1406	AA088829
Sequence 1355	AA040688	Sequence 1407	AA090106
Sequence 1356	AA040820	Sequence 1408	AA096032
Sequence 1357	AA041259	Sequence 1409	AA099819
Sequence 1358	AA043477	Sequence 1410	AA099923
Sequence 1359	AA044209	Sequence 1411	AA099976
Sequence 1360	AA044233	Sequence 1412	AA100764
Sequence 1361	AA044791	Sequence 1413	AA101010
Sequence 1362	AA045054	Sequence 1414	AA102013
Sequence 1363	AA045147	Sequence 1415	AA102564
Sequence 1364	AA045768	Sequence 1416	AA102830
Sequence 1365	AA046848	Sequence 1417	AA112186
Sequence 1366	AA053021	Sequence 1418	AA112645
Sequence 1367	AA053316	Sequence 1419	AA113305
Sequence 1368	AA053919	Sequence 1420	AA115218
Sequence 1369	AA054069	Sequence 1421	AA115315
Sequence 1370	AA055479	Sequence 1422	AA121656
Sequence 1371	AA055591	Sequence 1423	AA121718
Sequence 1372	AA055637	Sequence 1424	AA125809
Sequence 1373	AA057243	Sequence 1425	AA125939
Sequence 1374	AA058712	Sequence 1426	AA126452
Sequence 1375	AA059128	Sequence 1427	AA126718
Sequence 1376	AA065169	Sequence 1428	AA127436
Sequence 1377	AA069850	Sequence 1429	AA127666
Sequence 1378	AA071167	Sequence 1430	AA128063
Sequence 1379	AA075158	Sequence 1431	AA128636
Sequence 1380	AA075515	Sequence 1432	AA128641
Sequence 1381	AA075663	Sequence 1433	AA130778
Sequence 1382	AA076397	Sequence 1434	AA130982
Sequence 1383	AA076421	Sequence 1435	AA131827
Sequence 1384	AA078387	Sequence 1436	AA132056
Sequence 1385	AA078570	Sequence 1437	AA132163
Sequence 1386	AA078872	Sequence 1438	AA132574
Sequence 1387	AA079480	Sequence 1439	AA132992
Sequence 1388	AA080889	Sequence 1440	AA133351
Sequence 1389	AA081073	Sequence 1441	AA133474
Sequence 1390	AA081608	Sequence 1442	AA134460
Sequence 1391	AA081834	Sequence 1443	AA134527
Sequence 1392	AA081917	Sequence 1444	AA134589
Sequence 1393	AA082258	Sequence 1445	AA135696
Sequence 1394	AA082441	Sequence 1446	AA137017
Sequence 1395	AA083270	Sequence 1447	AA142941
Sequence 1396	AA083345	Sequence 1448	AA143001
Sequence 1397	AA083522	Sequence 1449	AA143074
Sequence 1398	AA083573	Sequence 1450	AA143746
Sequence 1399	AA083638	Sequence 1451	AA146900
Sequence 1400	AA083774	Sequence 1452	AA147200
Sequence 1401	AA088318	Sequence 1453	AA147247

Table 1

Sequence 1454	AA147781	Sequence 1506	AA188826
Sequence 1455	AA148027	Sequence 1507	AA190873
Sequence 1456	AA148136	Sequence 1508	AA191422
Sequence 1457	AA149810	Sequence 1509	AA192094
Sequence 1458	AA150377	Sequence 1510	AA193308
Sequence 1459	AA150837	Sequence 1511	AA194577
Sequence 1460	AA150928	Sequence 1512	AA195246
Sequence 1461	AA151274	Sequence 1513	AA195865
Sequence 1462	AA151594	Sequence 1514	AA196424
Sequence 1463	AA151755	Sequence 1515	AA196982
Sequence 1464	AA152476	Sequence 1516	AA203691
Sequence 1465	AA155754	Sequence 1517	AA204867
Sequence 1466	AA156066	Sequence 1518	AA206578
Sequence 1467	AA157163	Sequence 1519	AA206991
Sequence 1468	AA157993	Sequence 1520	AA209508
Sequence 1469	AA158738	Sequence 1521	AA216753
Sequence 1470	AA159110	Sequence 1522	AA219665
Sequence 1471	AA159576	Sequence 1523	AA223121
Sequence 1472	AA161003	Sequence 1524	AA223820
Sequence 1473	AA161076	Sequence 1525	AA224109
Sequence 1474	AA161467	Sequence 1526	AA224407
Sequence 1475	AA164193	Sequence 1527	AA227118
Sequence 1476	AA164473	Sequence 1528	AA229325
Sequence 1477	AA164729	Sequence 1529	AA229611
Sequence 1478	AA164873	Sequence 1530	AA232959
Sequence 1479	AA165027	Sequence 1531	AA233835
Sequence 1480	AA165068	Sequence 1532	AA233843
Sequence 1481	AA165087	Sequence 1533	AA234092
Sequence 1482	AA165174	Sequence 1534	AA234307
Sequence 1483	AA165282	Sequence 1535	AA236776
Sequence 1484	AA165293	Sequence 1536	AA242985
Sequence 1485	AA165638	Sequence 1537	AA243338
Sequence 1486	AA166618	Sequence 1538	AA244342
Sequence 1487	AA167041	Sequence 1539	AA249154
Sequence 1488	AA167750	Sequence 1540	AA255502
Sequence 1489	AA171630	Sequence 1541	AA256591
Sequence 1490	AA173506	Sequence 1542	AA261990
Sequence 1491	AA174097	Sequence 1543	AA262939
Sequence 1492	AA179187	Sequence 1544	AA278445
Sequence 1493	AA180137	Sequence 1545	AA278482
Sequence 1494	AA180224	Sequence 1546	AA278642
Sequence 1495	AA180383	Sequence 1547	AA278956
Sequence 1496	AA181075	Sequence 1548	AA279048
Sequence 1497	AA181258	Sequence 1549	AA280099
Sequence 1498	AA181684	Sequence 1550	AA280221
Sequence 1499	AA182415	Sequence 1551	AA280828
Sequence 1500	AA182540	Sequence 1552	AA282915
Sequence 1501	AA186577	Sequence 1553	AA284334
Sequence 1502	AA187817	Sequence 1554	AA284555
Sequence 1503	AA188045	Sequence 1555	AA284670
Sequence 1504	AA188140	Sequence 1556	AA284671
Sequence 1505	AA188384	Sequence 1557	AA284870

Table 1

Sequence 1558	AA284906	Sequence 1610	AA314872
Sequence 1559	AA285290	Sequence 1611	AA315363
Sequence 1560	AA286699	Sequence 1612	AA315379
Sequence 1561	AA286872	Sequence 1613	AA317243
Sequence 1562	AA287219	Sequence 1614	AA317393
Sequence 1563	AA287642	Sequence 1615	AA318969
Sequence 1564	AA287815	Sequence 1616	AA327201
Sequence 1565	AA291438	Sequence 1617	AA331991
Sequence 1566	AA291485	Sequence 1618	AA332672
Sequence 1567	AA291971	Sequence 1619	AA333358
Sequence 1568	AA292334	Sequence 1620	AA335273
Sequence 1569	AA293127	Sequence 1621	AA336666
Sequence 1570	AA293133	Sequence 1622	AA337192
Sequence 1571	AA293273	Sequence 1623	AA337489
Sequence 1572	AA293286	Sequence 1624	AA338793
Sequence 1573	AA293353	Sequence 1625	AA339957
Sequence 1574	AA293572	Sequence 1626	AA340341
Sequence 1575	AA293629	Sequence 1627	AA341446
Sequence 1576	AA293759	Sequence 1628	AA341465
Sequence 1577	AA293804	Sequence 1629	AA342969
Sequence 1578	AA296780	Sequence 1630	AA343629
Sequence 1579	AA297402	Sequence 1631	AA344084
Sequence 1580	AA298505	Sequence 1632	AA345329
Sequence 1581	AA299640	Sequence 1633	AA346393
Sequence 1582	AA301062	Sequence 1634	AA346698
Sequence 1583	AA301800	Sequence 1635	AA347887
Sequence 1584	AA303461	Sequence 1636	AA350059
Sequence 1585	AA303568	Sequence 1637	AA351507
Sequence 1586	AA306718	Sequence 1638	AA355003
Sequence 1587	AA306862	Sequence 1639	AA356682
Sequence 1588	AA306876	Sequence 1640	AA357574
Sequence 1589	AA307198	Sequence 1641	AA358887
Sequence 1590	AA307325	Sequence 1642	AA359705
Sequence 1591	AA308065	Sequence 1643	AA364352
Sequence 1592	AA308274	Sequence 1644	AA367451
Sequence 1593	AA308744	Sequence 1645	AA367773
Sequence 1594	AA310739	Sequence 1646	AA368542
Sequence 1595	AA310771	Sequence 1647	AA369400
Sequence 1596	AA311228	Sequence 1648	AA373230
Sequence 1597	AA311460	Sequence 1649	AA374754
Sequence 1598	AA311571	Sequence 1650	AA375312
Sequence 1599	AA311801	Sequence 1651	AA375815
Sequence 1600	AA311848	Sequence 1652	AA393525
Sequence 1601	AA311905	Sequence 1653	AA394115
Sequence 1602	AA312218	Sequence 1654	AA398443
Sequence 1603	AA312240	Sequence 1655	AA398585
Sequence 1604	AA312435	Sequence 1656	AA398739
Sequence 1605	AA313108	Sequence 1657	AA399165
Sequence 1606	AA313223	Sequence 1658	AA399628
Sequence 1607	AA313653	Sequence 1659	AA401329
Sequence 1608	AA313994	Sequence 1660	AA401334
Sequence 1609	AA314431	Sequence 1661	AA402191



Table 1

Sequence 1662	AA402289	Sequence 1714	AA476522
Sequence 1663	AA402775	Sequence 1715	AA477018
Sequence 1664	AA403319	Sequence 1716	AA477567
Sequence 1665	AA404613	Sequence 1717	AA477973
Sequence 1666	AA405124	Sequence 1718	AA478230
Sequence 1667	AA406239	Sequence 1719	AA479646
Sequence 1668	AA410580	Sequence 1720	AA479648
Sequence 1669	AA410982	Sequence 1721	AA479848
Sequence 1670	AA411021	Sequence 1722	AA481078
Sequence 1671	AA411252	Sequence 1723	AA481710
Sequence 1672	AA411764	Sequence 1724	AA482430
Sequence 1673	AA417794	Sequence 1725	AA482432
Sequence 1674	AA419263	Sequence 1726	AA482779
Sequence 1675	AA419284	Sequence 1727	AA483258
Sequence 1676	AA420751	Sequence 1728	AA483726
Sequence 1677	AA420758	Sequence 1729	AA483858
Sequence 1678	AA421248	Sequence 1730	AA484181
Sequence 1679	AA421682	Sequence 1731	AA486047
Sequence 1680	AA422060	Sequence 1732	AA486859
Sequence 1681	AA422143	Sequence 1733	AA488141
Sequence 1682	AA425004	Sequence 1734	AA488385
Sequence 1683	AA425468	Sequence 1735	AA488517
Sequence 1684	AA425737	Sequence 1736	AA489323
Sequence 1685	AA429794	Sequence 1737	AA489380
Sequence 1686	AA430400	Sequence 1738	AA489382
Sequence 1687	AA430436	Sequence 1739	AA491204
Sequence 1688	AA431428	Sequence 1740	AA492143
Sequence 1689	AA433988	Sequence 1741	AA493371
Sequence 1690	AA436315	Sequence 1742	AA494321
Sequence 1691	AA436411	Sequence 1743	AA494552
Sequence 1692	AA443024	Sequence 1744	AA501657
Sequence 1693	AA449394	Sequence 1745	AA502136
Sequence 1694	AA451779	Sequence 1746	AA505780
Sequence 1695	AA453878	Sequence 1747	AA512933
Sequence 1696	AA454668	Sequence 1748	AA514395
Sequence 1697	AA454953	Sequence 1749	AA514974
Sequence 1698	AA454962	Sequence 1750	AA515143
Sequence 1699	AA455245	Sequence 1751	AA516376
Sequence 1700	AA455785	Sequence 1752	AA521006
Sequence 1701	AA456454	Sequence 1753	AA523522
Sequence 1702	AA456557	Sequence 1754	AA524748
Sequence 1703	AA457255	Sequence 1755	AA524950
Sequence 1704	AA457579	Sequence 1756	AA525141
Sequence 1705	AA459167	Sequence 1757	AA526028
Sequence 1706	AA459210	Sequence 1758	AA527275
Sequence 1707	AA459527	Sequence 1759	AA527557
Sequence 1708	AA460570	Sequence 1760	AA533506
Sequence 1709	AA460816	Sequence 1761	AA534349
Sequence 1710	AA461005	Sequence 1762	AA534586
Sequence 1711	AA468657	Sequence 1763	AA534608
Sequence 1712	AA469447	Sequence 1764	AA535496
Sequence 1713	AA469453	Sequence 1765	AA541651

Table 1

Sequence 1766	AA548056	Sequence 1818	AA628536
Sequence 1767	AA548600	Sequence 1819	AA628547
Sequence 1768	AA550854	Sequence 1820	AA630611
Sequence 1769	AA550855	Sequence 1821	AA631326
Sequence 1770	AA551351	Sequence 1822	AA633909
Sequence 1771	AA551391	Sequence 1823	AA634260
Sequence 1772	AA554437	Sequence 1824	AA634298
Sequence 1773	AA554735	Sequence 1825	AA640505
Sequence 1774	AA555102	Sequence 1826	AA641289
Sequence 1775	AA564272	Sequence 1827	AA644625
Sequence 1776	AA564870	Sequence 1828	AA648944
Sequence 1777	AA565420	Sequence 1829	AA651720
Sequence 1778	AA568936	Sequence 1830	AA652478
Sequence 1779	AA569816	Sequence 1831	AA652505
Sequence 1780	AA569851	Sequence 1832	AA653775
Sequence 1781	AA569916	Sequence 1833	AA658374
Sequence 1782	AA573761	Sequence 1834	AA663005
Sequence 1783	AA573787	Sequence 1835	AA669154
Sequence 1784	AA577537	Sequence 1836	AA677560
Sequence 1785	AA578881	Sequence 1837	AA677750
Sequence 1786	AA579591	Sequence 1838	AA678185
Sequence 1787	AA579890	Sequence 1839	AA678251
Sequence 1788	AA580835	Sequence 1840	AA687495
Sequence 1789	AA582093	Sequence 1841	AA703208
Sequence 1790	AA582866	Sequence 1842	AA703667
Sequence 1791	AA583055	Sequence 1843	AA703907
Sequence 1792	AA583498	Sequence 1844	AA704208
Sequence 1793	AA583567	Sequence 1845	AA706347
Sequence 1794	AA583773	Sequence 1846	AA714010
Sequence 1795	AA584921	Sequence 1847	AA715984
Sequence 1796	AA586755	Sequence 1848	AA716651
Sequence 1797	AA587140	Sequence 1849	AA719530
Sequence 1798	AA587315	Sequence 1850	AA721642
Sequence 1799	AA587873	Sequence 1851	AA729381
Sequence 1800	AA593983	Sequence 1852	AA731946
Sequence 1801	AA594366	Sequence 1853	AA736817
Sequence 1802	AA595624	Sequence 1854	AA742713
Sequence 1803	AA595771	Sequence 1855	AA743278
Sequence 1804	AA599454	Sequence 1856	AA744681
Sequence 1805	AA600227	Sequence 1857	AA745953
Sequence 1806	AA600771	Sequence 1858	AA759195
Sequence 1807	AA601172	Sequence 1859	AA767779
Sequence 1808	AA602395	Sequence 1860	AA769697
Sequence 1809	AA602871	Sequence 1861	AA773998
Sequence 1810	AA603125	Sequence 1862	AA775058
Sequence 1811	AA603177	Sequence 1863	AA776593
Sequence 1812	AA604324	Sequence 1864	AA777384
Sequence 1813	AA604853	Sequence 1865	AA778672
Sequence 1814	AA610279	Sequence 1866	AA779949
Sequence 1815	AA610476	Sequence 1867	AA781487
Sequence 1816	AA610734	Sequence 1868	AA788907
Sequence 1817	AA614482	Sequence 1869	AA806278

Table 1

Sequence 1870	AA806735	Sequence 1922	AB007867
Sequence 1871	AA808769	Sequence 1923	AB007900
Sequence 1872	AA810149	Sequence 1924	AB007916
Sequence 1873	AA811609	Sequence 1925	AB007923
Sequence 1874	AA813604	Sequence 1926	AB007957
Sequence 1875	AA826307	Sequence 1927	AB011103
Sequence 1876	AA833766	Sequence 1928	AB011143
Sequence 1877	AA833900	Sequence 1929	AB011151
Sequence 1878	AA837457	Sequence 1930	AB011166
Sequence 1879	AA843531	Sequence 1931	AB014533
Sequence 1880	AA845737	Sequence 1932	AB014542
Sequence 1881	AA846698	Sequence 1933	AB014560
Sequence 1882	AA846856	Sequence 1934	AB015630
Sequence 1883	AA852896	Sequence 1935	AB015856
Sequence 1884	AA856902	Sequence 1936	AB018281
Sequence 1885	AA857824	Sequence 1937	AB018284
Sequence 1886	AA857882	Sequence 1938	AB018285
Sequence 1887	AA861665	Sequence 1939	AB018289
Sequence 1888	AA865960	Sequence 1940	AB018305
Sequence 1889	AA868529	Sequence 1941	AB018327
Sequence 1890	AA873271	Sequence 1942	AB018331
Sequence 1891	AA877189	Sequence 1943	AB018337
Sequence 1892	AA884922	Sequence 1944	AB019409
Sequence 1893	AA886453	Sequence 1945	AB019563
Sequence 1894	AA906652	Sequence 1946	AB019568
Sequence 1895	AA906865	Sequence 1947	AB019691
Sequence 1896	AA918993	Sequence 1948	AB020682
Sequence 1897	AA926926	Sequence 1949	AB020718
Sequence 1898	AA928934	Sequence 1950	AB021288
Sequence 1899	AA932501	Sequence 1951	AB023154
Sequence 1900	AA933987	Sequence 1952	AB023219
Sequence 1901	AA935947	Sequence 1953	AB024704
Sequence 1902	AA937302	Sequence 1954	AB027467
Sequence 1903	AA937773	Sequence 1955	AB028069
Sequence 1904	AA947835	Sequence 1956	AB028624
Sequence 1905	AA954939	Sequence 1957	AB028969
Sequence 1906	AA962587	Sequence 1958	AB028986
Sequence 1907	AA962632	Sequence 1959	AB029000
Sequence 1908	AA972525	Sequence 1960	AB029004
Sequence 1909	AA976489	Sequence 1961	AB029028
Sequence 1910	AA983380	Sequence 1962	AC03044
Sequence 1911	AA984586	Sequence 1963	AC31479
Sequence 1912	AA992596	Sequence 1964	AF000670
Sequence 1913	AB002305	Sequence 1965	AF000974
Sequence 1914	AB002330	Sequence 1966	AF001893
Sequence 1915	AB002357	Sequence 1967	AF004562
Sequence 1916	AB002806	Sequence 1968	AF006043
Sequence 1917	AB003476	Sequence 1969	AF007135
Sequence 1918	AB004066	Sequence 1970	AF007151
Sequence 1919	AB006077	Sequence 1971	AF007170
Sequence 1920	AB006534	Sequence 1972	AF009615
Sequence 1921	AB006755	Sequence 1973	AF013759

**Table 1**

Sequence 1974	AF013988	Sequence 2026	AF062318
Sequence 1975	AF015283	Sequence 2027	AF063611
Sequence 1976	AF015767	Sequence 2028	AF064019
Sequence 1977	AF016507	Sequence 2029	AF068235
Sequence 1978	AF016582	Sequence 2030	AF068846
Sequence 1979	AF017790	Sequence 2031	AF070523
Sequence 1980	AF019767	Sequence 2032	AF070537
Sequence 1981	AF021351	Sequence 2033	AF070555
Sequence 1982	AF021819	Sequence 2034	AF070561
Sequence 1983	AF022229	Sequence 2035	AF070596
Sequence 1984	AF023266	Sequence 2036	AF070600
Sequence 1985	AF025439	Sequence 2037	AF070626
Sequence 1986	AF026166	Sequence 2038	AF070649
Sequence 1987	AF026939	Sequence 2039	AF070662
Sequence 1988	AF027205	Sequence 2040	AF070672
Sequence 1989	AF031385	Sequence 2041	AF071202
Sequence 1990	AF034607	Sequence 2042	AF071219
Sequence 1991	AF035286	Sequence 2043	AF071593
Sequence 1992	AF035309	Sequence 2044	AF073298
Sequence 1993	AF035313	Sequence 2045	AF075587
Sequence 1994	AF037204	Sequence 2046	AF077030
Sequence 1995	AF038661	Sequence 2047	AF077045
Sequence 1996	AF039019	Sequence 2048	AF077200
Sequence 1997	AF039291	Sequence 2049	AF077202
Sequence 1998	AF039843	Sequence 2050	AF077207
Sequence 1999	AF040990	Sequence 2051	AF081192
Sequence 2000	AF041483	Sequence 2052	AF081484
Sequence 2001	AF042385	Sequence 2053	AF083190
Sequence 2002	AF042729	Sequence 2054	AF085355
Sequence 2003	AF044588	Sequence 2055	AF086003
Sequence 2004	AF045184	Sequence 2056	AF086116
Sequence 2005	AF047438	Sequence 2057	AF086178
Sequence 2006	AF047472	Sequence 2058	AF086205
Sequence 2007	AF048977	Sequence 2059	AF086207
Sequence 2008	AF050171	Sequence 2060	AF086336
Sequence 2009	AF050199	Sequence 2061	AF086517
Sequence 2010	AF050639	Sequence 2062	AF087135
Sequence 2011	AF052124	Sequence 2063	AF087990
Sequence 2012	AF052135	Sequence 2064	AF088036
Sequence 2013	AF052149	Sequence 2065	AF091076
Sequence 2014	AF052164	Sequence 2066	AF092563
Sequence 2015	AF052169	Sequence 2067	AF095287
Sequence 2016	AF052180	Sequence 2068	AF095791
Sequence 2017	AF052514	Sequence 2069	AF097709
Sequence 2018	AF054183	Sequence 2070	AF100741
Sequence 2019	AF054187	Sequence 2071	AF100756
Sequence 2020	AF054840	Sequence 2072	AF100928
Sequence 2021	AF055012	Sequence 2073	AF104222
Sequence 2022	AF055033	Sequence 2074	AF104913
Sequence 2023	AF057299	Sequence 2075	AF104923
Sequence 2024	AF059252	Sequence 2076	AF107405
Sequence 2025	AF061258	Sequence 2077	AF120334

Table 1

Sequence 2078	AF124438	Sequence 2130	AI090524
Sequence 2079	AF124439	Sequence 2131	AI090623
Sequence 2080	AF125525	Sequence 2132	AI091425
Sequence 2081	AF131799	Sequence 2133	AI092971
Sequence 2082	AF131814	Sequence 2134	AI095477
Sequence 2083	AF139461	Sequence 2135	AI123229
Sequence 2084	AF139658	Sequence 2136	AI125642
Sequence 2085	AF144755	Sequence 2137	AI125874
Sequence 2086	AF147331	Sequence 2138	AI127013
Sequence 2087	AF150962	Sequence 2139	AI127556
Sequence 2088	AF151832	Sequence 2140	AI140291
Sequence 2089	AF151868	Sequence 2141	AI141130
Sequence 2090	AF151898	Sequence 2142	AI141847
Sequence 2091	AF151907	Sequence 2143	AI143899
Sequence 2092	AF152097	Sequence 2144	AI144100
Sequence 2093	AF159295	Sequence 2145	AI148251
Sequence 2094	AF176702	Sequence 2146	AI149429
Sequence 2095	AF190744	Sequence 2147	AI149592
Sequence 2096	AI004664	Sequence 2148	AI186028
Sequence 2097	AI004915	Sequence 2149	AI186042
Sequence 2098	AI016073	Sequence 2150	AI190341
Sequence 2099	AI016323	Sequence 2151	AI192367
Sequence 2100	AI016791	Sequence 2152	AI192629
Sequence 2101	AI018451	Sequence 2153	AI198930
Sequence 2102	AI018625	Sequence 2154	AI216969
Sequence 2103	AI022779	Sequence 2155	AI217003
Sequence 2104	AI023799	Sequence 2156	AI223292
Sequence 2105	AI026164	Sequence 2157	AI241706
Sequence 2106	AI027516	Sequence 2158	AI251743
Sequence 2107	AI031636	Sequence 2159	AI252466
Sequence 2108	AI033037	Sequence 2160	AI253330
Sequence 2109	AI034115	Sequence 2161	AI253335
Sequence 2110	AI037859	Sequence 2162	AI253338
Sequence 2111	AI041670	Sequence 2163	AI253375
Sequence 2112	AI042034	Sequence 2164	AI253379
Sequence 2113	AI042290	Sequence 2165	AI253436
Sequence 2114	AI051971	Sequence 2166	AI262380
Sequence 2115	AI056917	Sequence 2167	AI263674
Sequence 2116	AI057124	Sequence 2168	AI267162
Sequence 2117	AI066419	Sequence 2169	AI267185
Sequence 2118	AI078041	Sequence 2170	AI267209
Sequence 2119	AI081116	Sequence 2171	AI267289
Sequence 2120	AI081472	Sequence 2172	AI267307
Sequence 2121	AI081913	Sequence 2173	AI267321
Sequence 2122	AI082244	Sequence 2174	AI267454
Sequence 2123	AI082648	Sequence 2175	AI267502
Sequence 2124	AI084731	Sequence 2176	AI268293
Sequence 2125	AI085381	Sequence 2177	AI269060
Sequence 2126	AI087291	Sequence 2178	AI269369
Sequence 2127	AI087819	Sequence 2179	AI270183
Sequence 2128	AI088178	Sequence 2180	AI270472
Sequence 2129	AI089981	Sequence 2181	AI271786

Table 1

Sequence 2182	AI272827	Sequence 2234	AI589301
Sequence 2183	AI274047	Sequence 2235	AI597938
Sequence 2184	AI276341	Sequence 2236	AI608591
Sequence 2185	AI276839	Sequence 2237	AI608787
Sequence 2186	AI278611	Sequence 2238	AI608968
Sequence 2187	AI280022	Sequence 2239	AI609193
Sequence 2188	AI283548	Sequence 2240	AI609281
Sequence 2189	AI288965	Sequence 2241	AI623804
Sequence 2190	AI290565	Sequence 2242	AI628689
Sequence 2191	AI291683	Sequence 2243	AI636635
Sequence 2192	AI292286	Sequence 2244	AI650837
Sequence 2193	AI298472	Sequence 2245	AI654096
Sequence 2194	AI298941	Sequence 2246	AI660245
Sequence 2195	AI304857	Sequence 2247	AI669253
Sequence 2196	AI308959	Sequence 2248	AI670084
Sequence 2197	AI312552	Sequence 2249	AI674313
Sequence 2198	AI333055	Sequence 2250	AI678152
Sequence 2199	AI333116	Sequence 2251	AI678703
Sequence 2200	AI335249	Sequence 2252	AI679044
Sequence 2201	AI336326	Sequence 2253	AI679321
Sequence 2202	AI345325	Sequence 2254	AI683140
Sequence 2203	AI366549	Sequence 2255	AI683338
Sequence 2204	AI367850	Sequence 2256	AI683793
Sequence 2205	AI375624	Sequence 2257	AI688798
Sequence 2206	AI376561	Sequence 2258	AI692866
Sequence 2207	AI399636	Sequence 2259	AI694087
Sequence 2208	AI417384	Sequence 2260	AI696819
Sequence 2209	AI421720	Sequence 2261	AI697501
Sequence 2210	AI424841	Sequence 2262	AI734922
Sequence 2211	AI431507	Sequence 2263	AI735069
Sequence 2212	AI433180	Sequence 2264	AI739337
Sequence 2213	AI434084	Sequence 2265	AI739377
Sequence 2214	AI434401	Sequence 2266	AI743595
Sequence 2215	AI436016	Sequence 2267	AI743691
Sequence 2216	AI436448	Sequence 2268	AI750198
Sequence 2217	AI446503	Sequence 2269	AI750909
Sequence 2218	AI453199	Sequence 2270	AI751119
Sequence 2219	AI459028	Sequence 2271	AI751364
Sequence 2220	AI469237	Sequence 2272	AI751565
Sequence 2221	AI492520	Sequence 2273	AI752319
Sequence 2222	AI492769	Sequence 2274	AI752553
Sequence 2223	AI494344	Sequence 2275	AI752929
Sequence 2224	AI523940	Sequence 2276	AI753108
Sequence 2225	AI524677	Sequence 2277	AI753671
Sequence 2226	AI538682	Sequence 2278	AI754437
Sequence 2227	AI557059	Sequence 2279	AI755181
Sequence 2228	AI561260	Sequence 2280	AI758869
Sequence 2229	AI567988	Sequence 2281	AI761927
Sequence 2230	AI569715	Sequence 2282	AI763126
Sequence 2231	AI581291	Sequence 2283	AI791906
Sequence 2232	AI583211	Sequence 2284	AI793120
Sequence 2233	AI583570	Sequence 2285	AI799521

**Table 1**

Sequence 2286	AI804346	Sequence 2338	AL049959
Sequence 2287	AI808109	Sequence 2339	AL049987
Sequence 2288	AI811021	Sequence 2340	AL049999
Sequence 2289	AI811845	Sequence 2341	AL050011
Sequence 2290	AI814139	Sequence 2342	AL050089
Sequence 2291	AI814674	Sequence 2343	AL050141
Sequence 2292	AI815868	Sequence 2344	AL050171
Sequence 2293	AI822030	Sequence 2345	AL050187
Sequence 2294	AI827641	Sequence 2346	AL050198
Sequence 2295	AI859619	Sequence 2347	AL050217
Sequence 2296	AI864580	Sequence 2348	AL050392
Sequence 2297	AI878968	Sequence 2349	AL080062
Sequence 2298	AI879179	Sequence 2350	AL080186
Sequence 2299	AI879367	Sequence 2351	AL080235
Sequence 2300	AI879992	Sequence 2352	AL096857
Sequence 2301	AI888377	Sequence 2353	AL096858
Sequence 2302	AI911704	Sequence 2354	AL110197
Sequence 2303	AI911997	Sequence 2355	AL110235
Sequence 2304	AI912084	Sequence 2356	AL117237
Sequence 2305	AI916284	Sequence 2357	AL117499
Sequence 2306	AI916584	Sequence 2358	AL117534
Sequence 2307	AI923224	Sequence 2359	AL118999
Sequence 2308	AI924096	Sequence 2360	AL119085
Sequence 2309	AI928185	Sequence 2361	AL119157
Sequence 2310	AI929819	Sequence 2362	AW020479
Sequence 2311	AI936748	Sequence 2363	AW044114
Sequence 2312	AI950087	Sequence 2364	AW102841
Sequence 2313	AI955808	Sequence 2365	C02094
Sequence 2314	AJ001258	Sequence 2366	C16886
Sequence 2315	AJ002030	Sequence 2367	C18886
Sequence 2316	AJ006026	Sequence 2368	D00017
Sequence 2317	AJ011001	Sequence 2369	D00022
Sequence 2318	AJ011915	Sequence 2370	D00068
Sequence 2319	AJ012499	Sequence 2371	D00099
Sequence 2320	AJ223183	Sequence 2372	D00422
Sequence 2321	AL035802	Sequence 2373	D10495
Sequence 2322	AL035987	Sequence 2374	D13119
Sequence 2323	AL036801	Sequence 2375	D13287
Sequence 2324	AL037646	Sequence 2376	D13665
Sequence 2325	AL038985	Sequence 2377	D13866
Sequence 2326	AL039150	Sequence 2378	D14662
Sequence 2327	AL041780	Sequence 2379	D14697
Sequence 2328	AL044019	Sequence 2380	D14710
Sequence 2329	AL046804	Sequence 2381	D14812
Sequence 2330	AL049055	Sequence 2382	D15049
Sequence 2331	AL049227	Sequence 2383	D16431
Sequence 2332	AL049229	Sequence 2384	D16937
Sequence 2333	AL049296	Sequence 2385	D17188
Sequence 2334	AL049464	Sequence 2386	D17268
Sequence 2335	AL049953	Sequence 2387	D17409
Sequence 2336	AL049954	Sequence 2388	D17793
Sequence 2337	AL049955	Sequence 2389	D21063

Table 1

Sequence 2390	D23660	Sequence 2442	E01650
Sequence 2391	D25542	Sequence 2443	E01797
Sequence 2392	D28759	Sequence 2444	E01813
Sequence 2393	D29677	Sequence 2445	E01827
Sequence 2394	D31767	Sequence 2446	E01979
Sequence 2395	D31784	Sequence 2447	E02628
Sequence 2396	D31883	Sequence 2448	E02651
Sequence 2397	D31890	Sequence 2449	E03569
Sequence 2398	D37991	Sequence 2450	E06721
Sequence 2399	D38491	Sequence 2451	E07218
Sequence 2400	D38583	Sequence 2452	F28779
Sequence 2401	D43948	Sequence 2453	F30276
Sequence 2402	D43950	Sequence 2454	F31082
Sequence 2403	D45248	Sequence 2455	H03854
Sequence 2404	D45887	Sequence 2456	H05412
Sequence 2405	D45915	Sequence 2457	H08994
Sequence 2406	D49489	Sequence 2458	H13339
Sequence 2407	D49547	Sequence 2459	H16426
Sequence 2408	D50310	Sequence 2460	H39960
Sequence 2409	D50371	Sequence 2461	H48742
Sequence 2410	D55192	Sequence 2462	H59372
Sequence 2411	D55649	Sequence 2463	H60722
Sequence 2412	D56120	Sequence 2464	H69238
Sequence 2413	D59253	Sequence 2465	H72481
Sequence 2414	D78586	Sequence 2466	H75695
Sequence 2415	D79826	Sequence 2467	H78517
Sequence 2416	D79983	Sequence 2468	H79084
Sequence 2417	D79986	Sequence 2469	H84729
Sequence 2418	D79997	Sequence 2470	H85709
Sequence 2419	D80006	Sequence 2471	H89654
Sequence 2420	D80012	Sequence 2472	J00269
Sequence 2421	D80087	Sequence 2473	J02621
Sequence 2422	D80253	Sequence 2474	J03005
Sequence 2423	D81635	Sequence 2475	J03040
Sequence 2424	D82128	Sequence 2476	J03171
Sequence 2425	D82348	Sequence 2477	J03191
Sequence 2426	D83197	Sequence 2478	J03210
Sequence 2427	D83327	Sequence 2479	J03464
Sequence 2428	D83784	Sequence 2480	J03473
Sequence 2429	D86227	Sequence 2481	J03799
Sequence 2430	D87437	Sequence 2482	J04080
Sequence 2431	D87442	Sequence 2483	J04164
Sequence 2432	D87470	Sequence 2484	J04177
Sequence 2433	D87666	Sequence 2485	J04765
Sequence 2434	D87667	Sequence 2486	J05013
Sequence 2435	D87682	Sequence 2487	J05021
Sequence 2436	D87735	Sequence 2488	J05192
Sequence 2437	D87969	Sequence 2489	J05633
Sequence 2438	D89052	Sequence 2490	K00558
Sequence 2439	D90226	Sequence 2491	K01566
Sequence 2440	D90373	Sequence 2492	K02765
Sequence 2441	E00882	Sequence 2493	L00160



Table 1

Sequence 2494	L02547	Sequence 2546	M24194
Sequence 2495	L05092	Sequence 2547	M25246
Sequence 2496	L05186	Sequence 2548	M26041
Sequence 2497	L07633	Sequence 2549	M26152
Sequence 2498	L11066	Sequence 2550	M26325
Sequence 2499	L11932	Sequence 2551	M27913
Sequence 2500	L12711	Sequence 2552	M27971
Sequence 2501	L13848	Sequence 2553	M28373
Sequence 2502	L14599	Sequence 2554	M31159
Sequence 2503	L19161	Sequence 2555	M31212
Sequence 2504	L19184	Sequence 2556	M31899
Sequence 2505	L19597	Sequence 2557	M32110
Sequence 2506	L20941	Sequence 2558	M32790
Sequence 2507	L23959	Sequence 2559	M32798
Sequence 2508	L26081	Sequence 2560	M33308
Sequence 2509	L27560	Sequence 2561	M34064
Sequence 2510	L28010	Sequence 2562	M37583
Sequence 2511	L28809	Sequence 2563	M38106
Sequence 2512	L33404	Sequence 2564	M55409
Sequence 2513	L33930	Sequence 2565	M55542
Sequence 2514	L34155	Sequence 2566	M58485
Sequence 2515	L34839	Sequence 2567	M60457
Sequence 2516	L38486	Sequence 2568	M60854
Sequence 2517	L42024	Sequence 2569	M62403
Sequence 2518	L43575	Sequence 2570	M62810
Sequence 2519	L44349	Sequence 2571	M64241
Sequence 2520	L54057	Sequence 2572	M67468
Sequence 2521	M10036	Sequence 2573	M69181
Sequence 2522	M10119	Sequence 2574	M74002
Sequence 2523	M10905	Sequence 2575	M75126
Sequence 2524	M11146	Sequence 2576	M76729
Sequence 2525	M13573	Sequence 2577	M78113
Sequence 2526	M13955	Sequence 2578	M81757
Sequence 2527	M14083	Sequence 2579	M83248
Sequence 2528	M14483	Sequence 2580	M84739
Sequence 2529	M14630	Sequence 2581	M87503
Sequence 2530	M14631	Sequence 2582	M88279
Sequence 2531	M15182	Sequence 2583	M92357
Sequence 2532	M15800	Sequence 2584	N20576
Sequence 2533	M16247	Sequence 2585	N34255
Sequence 2534	M16553	Sequence 2586	N35187
Sequence 2535	M16660	Sequence 2587	N35421
Sequence 2536	M16937	Sequence 2588	N39717
Sequence 2537	M17597	Sequence 2589	N40823
Sequence 2538	M17885	Sequence 2590	N40852
Sequence 2539	M20372	Sequence 2591	N67927
Sequence 2540	M22146	Sequence 2592	N76180
Sequence 2541	M22382	Sequence 2593	N76677
Sequence 2542	M22590	Sequence 2594	N77080
Sequence 2543	M22918	Sequence 2595	N84497
Sequence 2544	M22920	Sequence 2596	N86776
Sequence 2545	M23613	Sequence 2597	N91638

Table 1

Sequence 2598	N92086	Sequence 2650	U20896
Sequence 2599	N99205	Sequence 2651	U22431
Sequence 2600	Q37741	Sequence 2652	U22815
Sequence 2601	Q48043	Sequence 2653	U24105
Sequence 2602	Q65676	Sequence 2654	U24153
Sequence 2603	Q90526	Sequence 2655	U27768
Sequence 2604	R06046	Sequence 2656	U33760
Sequence 2605	R17092	Sequence 2657	U33833
Sequence 2606	R47228	Sequence 2658	U34877
Sequence 2607	R55150	Sequence 2659	U39361
Sequence 2608	R55398	Sequence 2660	U41515
Sequence 2609	R68132	Sequence 2661	U46570
Sequence 2610	R72676	Sequence 2662	U50733
Sequence 2611	R73306	Sequence 2663	U51586
Sequence 2612	R78333	Sequence 2664	U56255
Sequence 2613	R92367	Sequence 2665	U59305
Sequence 2614	R93637	Sequence 2666	U60975
Sequence 2615	R99649	Sequence 2667	U61083
Sequence 2616	S41458	Sequence 2668	U61397
Sequence 2617	S42303	Sequence 2669	U63846
Sequence 2618	S54005	Sequence 2670	U67784
Sequence 2619	S66431	Sequence 2671	U68723
Sequence 2620	S70154	Sequence 2672	U68727
Sequence 2621	S70290	Sequence 2673	U68758
Sequence 2622	S79895	Sequence 2674	U70735
Sequence 2623	S82076	Sequence 2675	U77085
Sequence 2624	T02792	Sequence 2676	U79258
Sequence 2625	T24119	Sequence 2677	U79274
Sequence 2626	T49314	Sequence 2678	U79278
Sequence 2627	T53479	Sequence 2679	U80213
Sequence 2628	T58797	Sequence 2680	U81234
Sequence 2629	T64560	Sequence 2681	U82130
Sequence 2630	T66112	Sequence 2682	U86602
Sequence 2631	T92160	Sequence 2683	U87309
Sequence 2632	T92396	Sequence 2684	U90028
Sequence 2633	U00947	Sequence 2685	U90441
Sequence 2634	U04815	Sequence 2686	U90902
Sequence 2635	U07151	Sequence 2687	U90917
Sequence 2636	U07857	Sequence 2688	U94831
Sequence 2637	U08470	Sequence 2689	V00478
Sequence 2638	U10323	Sequence 2690	V00503
Sequence 2639	U10439	Sequence 2691	V05728
Sequence 2640	U12465	Sequence 2692	V11636
Sequence 2641	U13665	Sequence 2693	V57903
Sequence 2642	U13877	Sequence 2694	V59662
Sequence 2643	U14550	Sequence 2695	V59746
Sequence 2644	U14966	Sequence 2696	V84428
Sequence 2645	U15008	Sequence 2697	V86232
Sequence 2646	U16306	Sequence 2698	V87930
Sequence 2647	U17104	Sequence 2699	W07215
Sequence 2648	U17496	Sequence 2700	W19127
Sequence 2649	U19769	Sequence 2701	W19407

Table 1

Sequence 2702	W19441	Sequence 2754	X71087
Sequence 2703	W25547	Sequence 2755	X73608
Sequence 2704	W26197	Sequence 2756	X73902
Sequence 2705	W38952	Sequence 2757	X74039
Sequence 2706	W56388	Sequence 2758	X74801
Sequence 2707	W68015	Sequence 2759	X74979
Sequence 2708	W73140	Sequence 2760	X76013
Sequence 2709	W73168	Sequence 2761	X76180
Sequence 2710	W76204	Sequence 2762	X78627
Sequence 2711	W87522	Sequence 2763	X81109
Sequence 2712	W87891	Sequence 2764	X82676
Sequence 2713	X00351	Sequence 2765	X84939
Sequence 2714	X00497	Sequence 2766	X85373
Sequence 2715	X01742	Sequence 2767	X93036
Sequence 2716	X01924	Sequence 2768	X93207
Sequence 2717	X03084	Sequence 2769	X94323
Sequence 2718	X04098	Sequence 2770	X94754
Sequence 2719	X04408	Sequence 2771	X97324
Sequence 2720	X04470	Sequence 2772	X99920
Sequence 2721	X05276	Sequence 2773	Y00503
Sequence 2722	X05908	Sequence 2774	Y00757
Sequence 2723	X06700	Sequence 2775	Y00815
Sequence 2724	X07819	Sequence 2776	Y09188
Sequence 2725	X13425	Sequence 2777	Y11435
Sequence 2726	X14420	Sequence 2778	Y12065
Sequence 2727	X15729	Sequence 2779	Y13247
Sequence 2728	X15880	Sequence 2780	Y13286
Sequence 2729	X16869	Sequence 2781	Y15286
Sequence 2730	X17206	Sequence 2782	Y17114
Sequence 2731	X24068	Sequence 2783	Z18538
Sequence 2732	X37385	Sequence 2784	Z18954
Sequence 2733	X37509	Sequence 2785	Z19054
Sequence 2734	X40178	Sequence 2786	Z21507
Sequence 2735	X51466	Sequence 2787	Z26317
Sequence 2736	X53505	Sequence 2788	Z29093
Sequence 2737	X54304	Sequence 2789	Z31696
Sequence 2738	X54941	Sequence 2790	Z32564
Sequence 2739	X55110	Sequence 2791	Z36531
Sequence 2740	X55885	Sequence 2792	Z37986
Sequence 2741	X56932	Sequence 2793	Z46629
Sequence 2742	X56998	Sequence 2794	Z47087
Sequence 2743	X56999	Sequence 2795	Z74615
Sequence 2744	X57766		
Sequence 2745	X62744		
Sequence 2746	X63432		
Sequence 2747	X66360		
Sequence 2748	X67698		
Sequence 2749	X68277		
Sequence 2750	X68880		
Sequence 2751	X69398		
Sequence 2752	X69838		
Sequence 2753	X70340		

Table 1

Sequence 340: found in patent publication W098/39446

AGGCGTNCCTCTGACTGCCCACTCAGTGGCNCACCGGGAGCTGNTTTGGNGCTTTGGG  
GANCCTNAACANTTNCNTCTTTCAAACCTNACTGGC

Sequence 1962: found in patent publication W098/42738

AGGTACCCGCTCTCCTGCTTCAGTAAATCTCCACTCGATCTCAGTGGGTTTCCTGTCCAT  
AGGATCCACAAGTTTGACCTGGCGGTGGAGCAAGGGGGCTTCACTAGGGATCATGGTTCC  
CCGGTAATCCATGGTCTTGCCAATGTAGCCGGTAATGTGTGTTTCCAGCCCTCCACGACCA  
CCCAGTTTCGCTGCCGGATAACTTGAACCACTTTGCCCTGCTTCCCGGCATCCTTGCCTT  
CTAGGATCTCCACCGTGTCCCCACAGAACAGATACCAAGTCTTCATCAGAGATGGGTTCC  
ACAACCACTGGGCGCCGNCCTGATCCATGGGGGGTTCCTCCTTGTCTGCAACAGAGCC  
TGGGGGGCTCATCCATAACGGTATGAGGGGGGAGAGTGACCTTGGATGCCAAGGCCAGC  
AGGGGCAAGAAAGACCCATGCCTGGAGGTTGNAAGAAAATCCCTTTGCCAGCAAAAACGC  
TTCGAAACCCTTNCCTTGTCAAGCTTTTTCACTTTTTCCGNGGCACCTTTGGGATTTTA  
GCACATTGGGGGCCCTTAAGNGTTCCTTCCCC

Sequence 341: found in patent publication W099/039941

CCCTTAGCGNGGTCGCGGCCGAGGCACAATTCGATTATTCACANGAAAGGGCAAACCTGTT  
NNTGTTNGCTGGCAGGAGNAGGTGCATATATACCAGCACTTCAAGTNNGGTATTTCCATT  
CAGGACATTTTATCTCTGTGCAAAGACCGGAGTAGAAGCTGATGAGTGGATCAAGATATT  
ACGCTGGAAATTGTCAAAATAAGAAAACAGCTCAACCAAGGGGAAGGCACCGATCCGAT  
CTCGGTCTGTTTATCTTTAATAGATCTTTCTTGCCAAGGAATGCTCTGGCCCAGGAGCAA  
GGTGAATGCTTCCCTGACGCTGCGATCTGCAGCAGACTNCAATGAAAACCGACTAAGG  
ATTTTCTTTCAAAAACAAATCAGAAGCAGATGCTGATTGGGACCCATATACCAGTTGCT  
GACTCACCGTTGCTGCCCTTNCATGGATGTTGCCATCTGCTTGAGAACACTGAAGCAATC  
ACCATTCTNGATANGAAAGTGCTTAAACCCCACTCTTAGGGCTGCTCACTTCTTAGAAC  
ACACAAAGGGAAGAGGAAAGGGGT

Sequence 342: found in patent publication W099/18126

CCGCGGTGGCGGCCCGCCGGCAGGTACCTACAGTGACACAGATCCCCTCCCGCCATCCT  
GGTCACACTGAATAACAAAGGGAAGAGAGGAGTAAGAACTGTAGTATCTAGAAATTCTCA  
GCACAGTGAAGGAAAGTGATCTTCTACTTTGTATTAGGCCTAAAAAGGAGGGGACGGG  
CCCGGCACAGTGGCTCACACCTGTAATCCAGCACTCTGGGAGGCCAAGGAGGGCAGATC  
ACCTGAGGTTGGGAGTTGAGACCAGCCTGACCAACATGGAGAAACCCTGTCTCTACTAA  
AAACACAAAATTAGCCAGGCATGGTGGCATGCGCCTATAAACCAAGCTACTCAAGAGGCT  
GAGGCAGGAGAATTGCTTGAACCCAGGAGGCAGAGGTTGTGGTGAGCCAAGATCGAGCCA  
TCACACTCCAGCCTGGGCAACAAGAGCAAACTCTGTCTCAAAAAAAAAAACAGGAGAGG  
AGGGAG

Sequence 1016: found in patent publication W099/38881

CTACTTAGGGCGAATTGGAGCTCCCCGCGGTGGCGGCCGAGGTCAAGCTTCGACCCCGCG  
TCCGTGATAAACTACTTTTGGGTTTTATTTTATTGAGGCATTTTTTTATTGTTTGAATG  
ATTCCGGCTTGTAATATATCAGCCTCTACAATGAAATGCAGAAGAGTTCATTTTTCTAG  
ATCTGTTTTTCATTAGAAATATTGACAAATAACACATTGTCAACCTGGATCCTTTGACA  
TTTACTTAACTCTGGCATGTTTCAAAAAAGTAGAACTCTAAGAGACCATTACCATTTT  
TCACAGATGTATAGGGGATGTATTCTAAAACTGACAGAAAAGAGAAATNTGATAGTCAAC  
ACTGTAACTTTTACTGNGTAATTGCCAAATACACTTTTCCAAATTTGTCCCAACAGCC  
TNTAAGCCAGCTTTCTTCTATATTTATAA

Sequence 1963: found in patent publication W099/46289

AACTGGACAGAGTAAGGGAATTCAGCATCCTCTTCTGCTTGTCTGTTACCCACAG

**Table 1**

ATCAAACCCTCAATTCTAGTTGGGGATGCTGTCTAGCCCCACACCATGACTGAAGCCTTA  
AGCACTGTTGCGCCTCATGTGCTTTGGATCAGCAACCCCAGTGGTATTCTACCAGAGCAT  
TGTGGGAAAGCAGATGTATAGTCAGGTCCCAACAGCAAAATTGTTGGGTGTGAGAGTTCTA  
AAGTATAGGGGTGAGGGAAGAGAAGGATATGAACTCCTCTGACCTTAAGCCAGCATTTCAT  
TTAACTTTTATGTCTACTTAACAAGAGAACCTGNAGAAAACCTACCGTATTCAAGAGATA  
ATCAAAATCAGTGTTTTAGCCAGGCGATGACAGAGAAGCACCATTCCCTCACCTCCATTC  
TTGTAATGTCTGTAATAAATTTAGTGCGTCAGGATGGATGAACCCAAGATCCAGTGAAT  
GATTCAGCTGTTCCAAGCCTTACATTTTCCATCATTTCATCATCCATTCTCATTGAGTA  
ACCTCTTGCACTATTGTGGTTAATTTTATGTAAAACCAGTTTATGTTTTTTTTTTAATAT  
GTGCCTATGTAATAAAAGTCTACACACTGGCAAAAAAAAAAAAAAAAAAAAAAGTCCTN

TABLE 1A

patent seq name	acc	dbase	Sequence 52	AA112308	dbEST
Sequence 1	AA001066	dbEST	Sequence 53	AA112375	dbEST
Sequence 2	AA007157	dbEST	Sequence 54	AA113860	dbEST
Sequence 3	AA010954	dbEST	Sequence 55	AA114120	dbEST
Sequence 4	AA015792	dbEST	Sequence 56	AA115118	dbEST
Sequence 5	AA019769	dbEST	Sequence 57	AA115368	dbEST
Sequence 6	AA019948	dbEST	Sequence 58	AA122286	dbEST
Sequence 7	AA022925	dbEST	Sequence 59	AA122348	dbEST
Sequence 8	AA022937	dbEST	Sequence 60	AA126109	dbEST
Sequence 9	AA024405	dbEST	Sequence 61	AA127105	dbEST
Sequence 10	AA029750	dbEST	Sequence 62	AA127132	dbEST
Sequence 11	AA031509	dbEST	Sequence 63	AA127418	dbEST
Sequence 12	AA033876	dbEST	Sequence 64	AA128305	dbEST
Sequence 13	AA034237	dbEST	Sequence 65	AA129461	dbEST
Sequence 14	AA039967	dbEST	Sequence 66	AA130252	dbEST
Sequence 15	AA040073	dbEST	Sequence 67	AA130547	dbEST
Sequence 16	AA040122	dbEST	Sequence 68	AA130786	dbEST
Sequence 17	AA045732	dbEST	Sequence 69	AA131041	dbEST
Sequence 18	AA045861	dbEST	Sequence 70	AA131065	dbEST
Sequence 19	AA046835	dbEST	Sequence 71	AA131104	dbEST
Sequence 20	AA047026	dbEST	Sequence 72	AA131155	dbEST
Sequence 21	AA047417	dbEST	Sequence 73	AA131160	dbEST
Sequence 22	AA053486	dbEST	Sequence 74	AA132182	dbEST
Sequence 23	AA054658	dbEST	Sequence 75	AA132568	dbEST
Sequence 24	AA055606	dbEST	Sequence 76	AA132598	dbEST
Sequence 25	AA056113	dbEST	Sequence 77	AA133351	dbEST
Sequence 26	AA056176	dbEST	Sequence 78	AA133927	dbEST
Sequence 27	AA056363	dbEST	Sequence 79	AA134105	dbEST
Sequence 28	AA056431	dbEST	Sequence 80	AA134210	dbEST
Sequence 29	AA065336	dbEST	Sequence 81	AA135032	dbEST
Sequence 30	AA069781	dbEST	Sequence 82	AA135919	dbEST
Sequence 31	AA069784	dbEST	Sequence 83	AA136383	dbEST
Sequence 32	AA069839	dbEST	Sequence 84	AA136789	dbEST
Sequence 33	AA069983	dbEST	Sequence 85	AA143609	dbEST
Sequence 34	AA071255	dbEST	Sequence 86	AA146773	dbEST
Sequence 35	AA075135	dbEST	Sequence 87	AA147806	dbEST
Sequence 36	AA081655	dbEST	Sequence 88	AA148160	dbEST
Sequence 37	AA082245	dbEST	Sequence 89	AA148268	dbEST
Sequence 38	AA083471	dbEST	Sequence 90	AA148771	dbEST
Sequence 39	AA083510	dbEST	Sequence 91	AA149056	dbEST
Sequence 40	AA085862	dbEST	Sequence 92	AA150307	dbEST
Sequence 41	AA085872	dbEST	Sequence 93	AA151310	dbEST
Sequence 42	AA085947	dbEST	Sequence 94	AA151775	dbEST
Sequence 43	AA088770	dbEST	Sequence 95	AA152037	dbEST
Sequence 44	AA100333	dbEST	Sequence 96	AA152416	dbEST
Sequence 45	AA100719	dbEST	Sequence 97	AA155853	dbEST
Sequence 46	AA100793	dbEST	Sequence 98	AA155926	dbEST
Sequence 47	AA100852	dbEST	Sequence 99	AA157405	dbEST
Sequence 48	AA101270	dbEST	Sequence 100	AA157725	dbEST
Sequence 49	AA101561	dbEST	Sequence 101	AA157788	dbEST
Sequence 50	AA111907	dbEST	Sequence 102	AA158165	dbEST
Sequence 51	AA112043	dbEST	Sequence 103	AA158171	dbEST

TABLE 1A

Sequence 104	AA159272	dbEST	Sequence 156	AA304961	dbEST
Sequence 105	AA160114	dbEST	Sequence 157	AA305193	dbEST
Sequence 106	AA160685	dbEST	Sequence 158	AA305438	dbEST
Sequence 107	AA161410	dbEST	Sequence 159	AA306542	dbEST
Sequence 108	AA164405	dbEST	Sequence 160	AA306708	dbEST
Sequence 109	AA164465	dbEST	Sequence 161	AA306945	dbEST
Sequence 110	AA165083	dbEST	Sequence 162	AA307239	dbEST
Sequence 111	AA165629	dbEST	Sequence 163	AA307477	dbEST
Sequence 112	AA166973	dbEST	Sequence 164	AA307504	dbEST
Sequence 113	AA171510	dbEST	Sequence 165	AA307697	dbEST
Sequence 114	AA173031	dbEST	Sequence 166	AA307779	dbEST
Sequence 115	AA173470	dbEST	Sequence 167	AA308062	dbEST
Sequence 116	AA173630	dbEST	Sequence 168	AA308801	dbEST
Sequence 117	AA179462	dbEST	Sequence 169	AA309028	dbEST
Sequence 118	AA187003	dbEST	Sequence 170	AA309988	dbEST
Sequence 119	AA187958	dbEST	Sequence 171	AA311006	dbEST
Sequence 120	AA188591	dbEST	Sequence 172	AA311481	dbEST
Sequence 121	AA192108	dbEST	Sequence 173	AA312012	dbEST
Sequence 122	AA199710	dbEST	Sequence 174	AA313684	dbEST
Sequence 123	AA203224	dbEST	Sequence 175	AA314146	dbEST
Sequence 124	AA203284	dbEST	Sequence 176	AA315049	dbEST
Sequence 125	AA205851	dbEST	Sequence 177	AA315308	dbEST
Sequence 126	AA209431	dbEST	Sequence 178	AA315426	dbEST
Sequence 127	AA209531	dbEST	Sequence 179	AA316682	dbEST
Sequence 128	AA214075	dbEST	Sequence 180	AA319958	dbEST
Sequence 129	AA216612	dbEST	Sequence 181	AA320346	dbEST
Sequence 130	AA224230	dbEST	Sequence 182	AA320991	dbEST
Sequence 131	AA224985	dbEST	Sequence 183	AA328544	dbEST
Sequence 132	AA226502	dbEST	Sequence 184	AA330457	dbEST
Sequence 133	AA229225	dbEST	Sequence 185	AA338793	dbEST
Sequence 134	AA232626	dbEST	Sequence 186	AA340069	dbEST
Sequence 135	AA233843	dbEST	Sequence 187	AA341170	dbEST
Sequence 136	AA242891	dbEST	Sequence 188	AA342394	dbEST
Sequence 137	AA250725	dbEST	Sequence 189	AA348250	dbEST
Sequence 138	AA250982	dbEST	Sequence 190	AA349148	dbEST
Sequence 139	AA256959	dbEST	Sequence 191	AA351443	dbEST
Sequence 140	AA259077	dbEST	Sequence 192	AA351880	dbEST
Sequence 141	AA262440	dbEST	Sequence 193	AA356158	dbEST
Sequence 142	AA263110	dbEST	Sequence 194	AA356187	dbEST
Sequence 143	AA283165	dbEST	Sequence 195	AA356195	dbEST
Sequence 144	AA285260	dbEST	Sequence 196	AA357374	dbEST
Sequence 145	AA287112	dbEST	Sequence 197	AA367446	dbEST
Sequence 146	AA292191	dbEST	Sequence 198	AA375236	dbEST
Sequence 147	AA292334	dbEST	Sequence 199	AA377718	dbEST
Sequence 148	AA292385	dbEST	Sequence 200	AA380997	dbEST
Sequence 149	AA292771	dbEST	Sequence 201	AA383917	dbEST
Sequence 150	AA293273	dbEST	Sequence 202	AA385147	dbEST
Sequence 151	AA293572	dbEST	Sequence 203	AA389641	dbEST
Sequence 152	AA295348	dbEST	Sequence 204	AA393164	dbEST
Sequence 153	AA295485	dbEST	Sequence 205	AA393236	dbEST
Sequence 154	AA301631	dbEST	Sequence 206	AA394242	dbEST
Sequence 155	AA304669	dbEST	Sequence 207	AA398732	dbEST

TABLE 1A

Sequence 208	AA401864	dbEST	Sequence 260	AA573893	dbEST
Sequence 209	AA410508	dbEST	Sequence 261	AA574237	dbEST
Sequence 210	AA410580	dbEST	Sequence 262	AA576866	dbEST
Sequence 211	AA410942	dbEST	Sequence 263	AA579034	dbEST
Sequence 212	AA411334	dbEST	Sequence 264	AA579816	dbEST
Sequence 213	AA411599	dbEST	Sequence 265	AA581220	dbEST
Sequence 214	AA418061	dbEST	Sequence 266	AA581264	dbEST
Sequence 215	AA418473	dbEST	Sequence 267	AA582093	dbEST
Sequence 216	AA418970	dbEST	Sequence 268	AA583091	dbEST
Sequence 217	AA420789	dbEST	Sequence 269	AA584411	dbEST
Sequence 218	AA421682	dbEST	Sequence 270	AA586776	dbEST
Sequence 219	AA421850	dbEST	Sequence 271	AA587110	dbEST
Sequence 220	AA424529	dbEST	Sequence 272	AA587233	dbEST
Sequence 221	AA428421	dbEST	Sequence 273	AA587700	dbEST
Sequence 222	AA429754	dbEST	Sequence 274	AA609259	dbEST
Sequence 223	AA441787	dbEST	Sequence 275	AA609837	dbEST
Sequence 224	AA451633	dbEST	Sequence 276	AA613907	dbEST
Sequence 225	AA453309	dbEST	Sequence 277	AA614529	dbEST
Sequence 226	AA453559	dbEST	Sequence 278	AA618033	dbEST
Sequence 227	AA453570	dbEST	Sequence 279	AA628487	dbEST
Sequence 228	AA454871	dbEST	Sequence 280	AA631204	dbEST
Sequence 229	AA454913	dbEST	Sequence 281	AA631811	dbEST
Sequence 230	AA456892	dbEST	Sequence 282	AA640901	dbEST
Sequence 231	AA457048	dbEST	Sequence 283	AA641841	dbEST
Sequence 232	AA463426	dbEST	Sequence 284	AA642215	dbEST
Sequence 233	AA465039	dbEST	Sequence 285	AA643602	dbEST
Sequence 234	AA477173	dbEST	Sequence 286	AA651720	dbEST
Sequence 235	AA480921	dbEST	Sequence 287	AA664996	dbEST
Sequence 236	AA484050	dbEST	Sequence 288	AA668297	dbEST
Sequence 237	AA484756	dbEST	Sequence 289	AA668836	dbEST
Sequence 238	AA487483	dbEST	Sequence 290	AA675923	dbEST
Sequence 239	AA489640	dbEST	Sequence 291	AA687833	dbEST
Sequence 240	AA493886	dbEST	Sequence 292	AA704992	dbEST
Sequence 241	AA494493	dbEST	Sequence 293	AA732702	dbEST
Sequence 242	AA496518	dbEST	Sequence 294	AA745241	dbEST
Sequence 243	AA501749	dbEST	Sequence 295	AA746481	dbEST
Sequence 244	AA501822	dbEST	Sequence 296	AA758889	dbEST
Sequence 245	AA501945	dbEST	Sequence 297	AA772570	dbEST
Sequence 246	AA504490	dbEST	Sequence 298	AA772790	dbEST
Sequence 247	AA507234	dbEST	Sequence 299	AA776709	dbEST
Sequence 248	AA513640	dbEST	Sequence 300	AA776811	dbEST
Sequence 249	AA526227	dbEST	Sequence 301	AA777384	dbEST
Sequence 250	AA526889	dbEST	Sequence 302	AA778116	dbEST
Sequence 251	AA527139	dbEST	Sequence 303	AA779868	dbEST
Sequence 252	AA527188	dbEST	Sequence 304	AA781343	dbEST
Sequence 253	AA531428	dbEST	Sequence 305	AA809984	dbEST
Sequence 254	AA532633	dbEST	Sequence 306	AA810945	dbEST
Sequence 255	AA535471	dbEST	Sequence 307	AA811200	dbEST
Sequence 256	AA554757	dbEST	Sequence 308	AA825768	dbEST
Sequence 257	AA565996	dbEST	Sequence 309	AA828073	dbEST
Sequence 258	AA568217	dbEST	Sequence 310	AA828722	dbEST
Sequence 259	AA573742	dbEST	Sequence 311	AA843176	dbEST



TABLE 1A

Sequence 312	AA843661	dbEST	Sequence 364	AF038451	ANUC
Sequence 313	AA876526	dbEST	Sequence 365	AF038662	ANUC
Sequence 314	AA883255	dbEST	Sequence 366	AF038963	ANUC
Sequence 315	AA906652	dbEST	Sequence 367	AF043431	ANUC
Sequence 316	AA917638	dbEST	Sequence 368	AF044956	ANUC
Sequence 317	AA927734	dbEST	Sequence 369	AF045941	ANUC
Sequence 318	AA954939	dbEST	Sequence 370	AF046997	ANUC
Sequence 319	AA962622	dbEST	Sequence 371	AF051894	ANUC
Sequence 320	AA991285	dbEST	Sequence 372	AF052124	ANUC
Sequence 321	AB000115	ANUC	Sequence 373	AF052578	ANUC
Sequence 322	AB004047	ANUC	Sequence 374	AF053233	ANUC
Sequence 323	AB006746	ANUC	Sequence 375	AF054838	ANUC
Sequence 324	AB007619	ANUC	Sequence 376	AF055012	ANUC
Sequence 325	AB007860	ANUC	Sequence 377	AF061736	ANUC
Sequence 326	AB007965	ANUC	Sequence 378	AF061738	ANUC
Sequence 327	AB011101	ANUC	Sequence 379	AF064603	ANUC
Sequence 328	AB011169	ANUC	Sequence 380	AF064854	ANUC
Sequence 329	AB012701	ANUC	Sequence 381	AF065388	ANUC
Sequence 330	AB014536	ANUC	Sequence 382	AF067168	ANUC
Sequence 331	AB014565	ANUC	Sequence 383	AF067174	ANUC
Sequence 332	AB019568	ANUC	Sequence 384	AF067817	ANUC
Sequence 333	AB020623	ANUC	Sequence 385	AF070523	ANUC
Sequence 334	AB020629	ANUC	Sequence 386	AF070561	ANUC
Sequence 335	AB020693	ANUC	Sequence 387	AF070562	ANUC
Sequence 336	AB021288	ANUC	Sequence 388	AF070596	ANUC
Sequence 337	AB022663	ANUC	Sequence 389	AF070664	ANUC
Sequence 338	AB023214	ANUC	Sequence 390	AF070674	ANUC
Sequence 339	AB023230	ANUC	Sequence 391	AF077048	ANUC
Sequence 340	AC02059	PREPATNUC	Sequence 392	AF077051	ANUC
Sequence 341	AC03653	PREPATNUC	Sequence 393	AF077200	ANUC
Sequence 342	AC13415	PREPATNUC	Sequence 394	AF077671	ANUC
Sequence 343	AF000982	ANUC	Sequence 395	AF080246	ANUC
Sequence 344	AF002985	ANUC	Sequence 396	AF081484	ANUC
Sequence 345	AF005654	ANUC	Sequence 397	AF083470	ANUC
Sequence 346	AF006086	ANUC	Sequence 398	AF084523	ANUC
Sequence 347	AF007791	ANUC	Sequence 399	AF085355	ANUC
Sequence 348	AF013758	ANUC	Sequence 400	AF086003	ANUC
Sequence 349	AF013988	ANUC	Sequence 401	AF086080	ANUC
Sequence 350	AF021232	ANUC	Sequence 402	AF086183	ANUC
Sequence 351	AF026939	ANUC	Sequence 403	AF086545	ANUC
Sequence 352	AF026941	ANUC	Sequence 404	AF091263	ANUC
Sequence 353	AF026942	ANUC	Sequence 405	AF111713	ANUC
Sequence 354	AF026943	ANUC	Sequence 406	AF118023	ANUC
Sequence 355	AF026944	ANUC	Sequence 407	AF124438	ANUC
Sequence 356	AF028832	ANUC	Sequence 408	AF124439	ANUC
Sequence 357	AF030455	ANUC	Sequence 409	AF131808	ANUC
Sequence 358	AF030514	ANUC	Sequence 410	AF131820	ANUC
Sequence 359	AF031469	ANUC	Sequence 411	AF131848	ANUC
Sequence 360	AF033095	ANUC	Sequence 412	AF132966	ANUC
Sequence 361	AF035286	ANUC	Sequence 413	AF132968	ANUC
Sequence 362	AF035316	ANUC	Sequence 414	AF146277	ANUC
Sequence 363	AF037204	ANUC	Sequence 415	AF147331	ANUC

TABLE 1A

Sequence 416	AF150100	ANUC	Sequence 468	AI417973	dbEST
Sequence 417	AF150266	dbEST	Sequence 469	AI431963	dbEST
Sequence 418	AF151873	ANUC	Sequence 470	AI453405	dbEST
Sequence 419	AF151877	ANUC	Sequence 471	AI457157	dbEST
Sequence 420	AF151978	ANUC	Sequence 472	AI457624	dbEST
Sequence 421	AF167160	ANUC	Sequence 473	AI459679	dbEST
Sequence 422	AI023413	dbEST	Sequence 474	AI460010	dbEST
Sequence 423	AI027888	dbEST	Sequence 475	AI469095	dbEST
Sequence 424	AI031811	dbEST	Sequence 476	AI469715	dbEST
Sequence 425	AI033687	dbEST	Sequence 477	AI471539	dbEST
Sequence 426	AI042140	dbEST	Sequence 478	AI476335	dbEST
Sequence 427	AI075324	dbEST	Sequence 479	AI479289	dbEST
Sequence 428	AI075876	dbEST	Sequence 480	AI499285	dbEST
Sequence 429	AI126802	dbEST	Sequence 481	AI521180	dbEST
Sequence 430	AI127556	dbEST	Sequence 482	AI538061	dbEST
Sequence 431	AI129360	dbEST	Sequence 483	AI567204	dbEST
Sequence 432	AI139456	dbEST	Sequence 484	AI587104	dbEST
Sequence 433	AI140291	dbEST	Sequence 485	AI587328	dbEST
Sequence 434	AI144215	dbEST	Sequence 486	AI609624	dbEST
Sequence 435	AI161378	dbEST	Sequence 487	AI610607	dbEST
Sequence 436	AI188638	dbEST	Sequence 488	AI612873	dbEST
Sequence 437	AI215617	dbEST	Sequence 489	AI627444	dbEST
Sequence 438	AI216969	dbEST	Sequence 490	AI632869	dbEST
Sequence 439	AI241578	dbEST	Sequence 491	AI633164	dbEST
Sequence 440	AI250167	dbEST	Sequence 492	AI636014	dbEST
Sequence 441	AI253330	dbEST	Sequence 493	AI637620	dbEST
Sequence 442	AI253335	dbEST	Sequence 494	AI676218	dbEST
Sequence 443	AI253369	dbEST	Sequence 495	AI683871	dbEST
Sequence 444	AI253436	dbEST	Sequence 496	AI684170	dbEST
Sequence 445	AI261671	dbEST	Sequence 497	AI693877	dbEST
Sequence 446	AI262264	dbEST	Sequence 498	AI694088	dbEST
Sequence 447	AI267162	dbEST	Sequence 499	AI732534	dbEST
Sequence 448	AI267379	dbEST	Sequence 500	AI743595	dbEST
Sequence 449	AI267502	dbEST	Sequence 501	AI744489	dbEST
Sequence 450	AI267622	dbEST	Sequence 502	AI745058	dbEST
Sequence 451	AI279131	dbEST	Sequence 503	AI753108	dbEST
Sequence 452	AI285943	dbEST	Sequence 504	AI791322	dbEST
Sequence 453	AI289173	dbEST	Sequence 505	AI798474	dbEST
Sequence 454	AI290876	dbEST	Sequence 506	AI803838	dbEST
Sequence 455	AI292104	dbEST	Sequence 507	AI811960	dbEST
Sequence 456	AI300033	dbEST	Sequence 508	AI813617	dbEST
Sequence 457	AI300074	dbEST	Sequence 509	AI815829	dbEST
Sequence 458	AI312113	dbEST	Sequence 510	AI826957	dbEST
Sequence 459	AI336032	dbEST	Sequence 511	AI831002	dbEST
Sequence 460	AI337069	dbEST	Sequence 512	AI863041	dbEST
Sequence 461	AI340262	dbEST	Sequence 513	AI867294	dbEST
Sequence 462	AI346975	dbEST	Sequence 514	AI912076	dbEST
Sequence 463	AI354639	dbEST	Sequence 515	AI915553	dbEST
Sequence 464	AI366381	dbEST	Sequence 516	AJ001381	ANUC
Sequence 465	AI369024	dbEST	Sequence 517	AJ003401	dbEST
Sequence 466	AI382020	dbEST	Sequence 518	AJ010071	ANUC
Sequence 467	AI400372	dbEST	Sequence 519	AJ132502	ANUC

TABLE 1A

Sequence 520	AL044356	dbEST	Sequence 572	E01197	ANUC
Sequence 521	AL044825	dbEST	Sequence 573	E01198	ANUC
Sequence 522	AL047024	dbEST	Sequence 574	E01630	ANUC
Sequence 523	AL048393	dbEST	Sequence 575	E01954	ANUC
Sequence 524	AL049313	ANUC	Sequence 576	E01971	ANUC
Sequence 525	AL049923	ANUC	Sequence 577	E01972	ANUC
Sequence 526	AL049954	ANUC	Sequence 578	E02628	ANUC
Sequence 527	AL050024	ANUC	Sequence 579	E03569	ANUC
Sequence 528	AL050272	ANUC	Sequence 580	E03879	ANUC
Sequence 529	AL050395	ANUC	Sequence 581	E08663	ANUC
Sequence 530	AL096714	ANUC	Sequence 582	F06593	dbEST
Sequence 531	AL096748	ANUC	Sequence 583	F28779	dbEST
Sequence 532	AL096842	ANUC	Sequence 584	H25806	dbEST
Sequence 533	AL110124	ANUC	Sequence 585	H47546	dbEST
Sequence 534	C17346	dbEST	Sequence 586	H48873	dbEST
Sequence 535	D00017	ANUC	Sequence 587	H66467	dbEST
Sequence 536	D00068	ANUC	Sequence 588	H88415	dbEST
Sequence 537	D11960	dbEST	Sequence 589	J00196	ANUC
Sequence 538	D12502	ANUC	Sequence 590	J03575	ANUC
Sequence 539	D12763	ANUC	Sequence 591	J03858	ANUC
Sequence 540	D13380	ANUC	Sequence 592	J03909	ANUC
Sequence 541	D13645	ANUC	Sequence 593	J04164	ANUC
Sequence 542	D13866	ANUC	Sequence 594	K00422	ANUC
Sequence 543	D14697	ANUC	Sequence 595	K01763	ANUC
Sequence 544	D21260	ANUC	Sequence 596	L00693	ANUC
Sequence 545	D23660	ANUC	Sequence 597	L02426	ANUC
Sequence 546	D26155	ANUC	Sequence 598	L06328	ANUC
Sequence 547	D26599	ANUC	Sequence 599	L09159	ANUC
Sequence 548	D28759	ANUC	Sequence 600	L10413	ANUC
Sequence 549	D29640	ANUC	Sequence 601	L11066	ANUC
Sequence 550	D31763	ANUC	Sequence 602	L20688	ANUC
Sequence 551	D31767	ANUC	Sequence 603	L20941	ANUC
Sequence 552	D31883	ANUC	Sequence 604	L28997	ANUC
Sequence 553	D38524	ANUC	Sequence 605	L38995	ANUC
Sequence 554	D42040	ANUC	Sequence 606	L41490	ANUC
Sequence 555	D45248	ANUC	Sequence 607	M10119	ANUC
Sequence 556	D49396	ANUC	Sequence 608	M13536	ANUC
Sequence 557	D50372	ANUC	Sequence 609	M14328	ANUC
Sequence 558	D50420	ANUC	Sequence 610	M14764	ANUC
Sequence 559	D55653	ANUC	Sequence 611	M15329	ANUC
Sequence 560	D81522	dbEST	Sequence 612	M16660	ANUC
Sequence 561	D83077	ANUC	Sequence 613	M17017	ANUC
Sequence 562	D83767	ANUC	Sequence 614	M18216	ANUC
Sequence 563	D86958	ANUC	Sequence 615	M19723	ANUC
Sequence 564	D86979	ANUC	Sequence 616	M22918	ANUC
Sequence 565	D87666	ANUC	Sequence 617	M23613	ANUC
Sequence 566	D87667	ANUC	Sequence 618	M24194	ANUC
Sequence 567	D87735	ANUC	Sequence 619	M24594	ANUC
Sequence 568	D88532	ANUC	Sequence 620	M26152	ANUC
Sequence 569	D89053	ANUC	Sequence 621	M29540	ANUC
Sequence 570	D90311	ANUC	Sequence 622	M29541	ANUC
Sequence 571	D90453	ANUC	Sequence 623	M29551	ANUC

TABLE 1A

Sequence 624	M33146	ANUC	Sequence 676	U07550	ANUC
Sequence 625	M34064	ANUC	Sequence 677	U07857	ANUC
Sequence 626	M34455	ANUC	Sequence 678	U08815	ANUC
Sequence 627	M35198	ANUC	Sequence 679	U09559	ANUC
Sequence 628	M36693	ANUC	Sequence 680	U09847	ANUC
Sequence 629	M37716	ANUC	Sequence 681	U10439	ANUC
Sequence 630	M55268	ANUC	Sequence 682	U14966	ANUC
Sequence 631	M55542	ANUC	Sequence 683	U18321	ANUC
Sequence 632	M55543	ANUC	Sequence 684	U19878	ANUC
Sequence 633	M57567	ANUC	Sequence 685	U23942	ANUC
Sequence 634	M60333	ANUC	Sequence 686	U25789	ANUC
Sequence 635	M61715	ANUC	Sequence 687	U28249	ANUC
Sequence 636	M62831	ANUC	Sequence 688	U28964	ANUC
Sequence 637	M63121	ANUC	Sequence 689	U32500	ANUC
Sequence 638	M63838	ANUC	Sequence 690	U32944	ANUC
Sequence 639	M68520	ANUC	Sequence 691	U33760	ANUC
Sequence 640	M77945	dbEST	Sequence 692	U37230	ANUC
Sequence 641	M80563	ANUC	Sequence 693	U37518	ANUC
Sequence 642	M81757	ANUC	Sequence 694	U38292	ANUC
Sequence 643	M83248	ANUC	Sequence 695	U38784	ANUC
Sequence 644	M83654	ANUC	Sequence 696	U41371	ANUC
Sequence 645	M86553	ANUC	Sequence 697	U41515	ANUC
Sequence 646	M87284	ANUC	Sequence 698	U52513	ANUC
Sequence 647	M87434	ANUC	Sequence 699	U56255	ANUC
Sequence 648	M87503	ANUC	Sequence 700	U57847	ANUC
Sequence 649	M92357	ANUC	Sequence 701	U61083	ANUC
Sequence 650	M96982	ANUC	Sequence 702	U68758	ANUC
Sequence 651	M97501	ANUC	Sequence 703	U73524	ANUC
Sequence 652	M97935	ANUC	Sequence 704	U77085	ANUC
Sequence 653	N36346	dbEST	Sequence 705	U78722	ANUC
Sequence 654	N51262	dbEST	Sequence 706	U79751	ANUC
Sequence 655	N57413	dbEST	Sequence 707	U94586	ANUC
Sequence 656	N78477	dbEST	Sequence 708	V00572	ANUC
Sequence 657	N92060	dbEST	Sequence 709	V00594	ANUC
Sequence 658	Q21065	NUCPATENT	Sequence 710	V04202	NUCPATENT
Sequence 659	Q94780	NUCPATENT	Sequence 711	V17906	NUCPATENT
Sequence 660	R13925	dbEST	Sequence 712	V36078	NUCPATENT
Sequence 661	R51732	dbEST	Sequence 713	V68140	NUCPATENT
Sequence 662	R56461	dbEST	Sequence 714	V86134	NUCPATENT
Sequence 663	R66489	dbEST	Sequence 715	W02908	dbEST
Sequence 664	R75621	dbEST	Sequence 716	W05711	dbEST
Sequence 665	S45630	ANUC	Sequence 717	W07308	dbEST
Sequence 666	S70290	ANUC	Sequence 718	W25547	dbEST
Sequence 667	S75295	ANUC	Sequence 719	W28837	dbEST
Sequence 668	S76638	ANUC	Sequence 720	W37272	dbEST
Sequence 669	T34641	dbEST	Sequence 721	W38644	dbEST
Sequence 670	T50925	NUCPATENT	Sequence 722	W39262	dbEST
Sequence 671	T52715	dbEST	Sequence 723	W39498	dbEST
Sequence 672	T54951	dbEST	Sequence 724	W52254	dbEST
Sequence 673	T70793	dbEST	Sequence 725	W74319	dbEST
Sequence 674	U03886	ANUC	Sequence 726	W77987	dbEST
Sequence 675	U04313	ANUC	Sequence 727	W80480	dbEST

TABLE 1A

Sequence 728	X00637	ANUC	Sequence 780	AA028164	dbEST
Sequence 729	X01742	ANUC	Sequence 781	AA035775	dbEST
Sequence 730	X02530	ANUC	Sequence 782	AA037294	dbEST
Sequence 731	X02661	ANUC	Sequence 783	AA039967	dbEST
Sequence 732	X04316	NUCPATENT	Sequence 784	AA045637	dbEST
Sequence 733	X04371	ANUC	Sequence 785	AA046815	dbEST
Sequence 734	X04470	ANUC	Sequence 786	AA046853	dbEST
Sequence 735	X05908	ANUC	Sequence 787	AA047052	dbEST
Sequence 736	X07819	ANUC	Sequence 788	AA047213	dbEST
Sequence 737	X13238	ANUC	Sequence 789	AA057071	dbEST
Sequence 738	X15674	ANUC	Sequence 790	AA058933	dbEST
Sequence 739	X15729	ANUC	Sequence 791	AA064952	dbEST
Sequence 740	X16354	ANUC	Sequence 792	AA075089	dbEST
Sequence 741	X16356	ANUC	Sequence 793	AA076291	dbEST
Sequence 742	X16455	ANUC	Sequence 794	AA078508	dbEST
Sequence 743	X17025	ANUC	Sequence 795	AA080864	dbEST
Sequence 744	X20432	NUCPATENT	Sequence 796	AA083345	dbEST
Sequence 745	X30167	NUCPATENT	Sequence 797	AA083693	dbEST
Sequence 746	X33937	NUCPATENT	Sequence 798	AA085497	dbEST
Sequence 747	X35726	NUCPATENT	Sequence 799	AA086463	dbEST
Sequence 748	X41105	NUCPATENT	Sequence 800	AA093935	dbEST
Sequence 749	X51841	ANUC	Sequence 801	AA100291	dbEST
Sequence 750	X54941	ANUC	Sequence 802	AA101207	dbEST
Sequence 751	X56932	ANUC	Sequence 803	AA102403	dbEST
Sequence 752	X57351	ANUC	Sequence 804	AA111856	dbEST
Sequence 753	X59710	ANUC	Sequence 805	AA115174	dbEST
Sequence 754	X65614	ANUC	Sequence 806	AA122134	dbEST
Sequence 755	X67951	ANUC	Sequence 807	AA122291	dbEST
Sequence 756	X68060	ANUC	Sequence 808	AA125780	dbEST
Sequence 757	X68277	ANUC	Sequence 809	AA127322	dbEST
Sequence 758	X72790	ANUC	Sequence 810	AA130432	dbEST
Sequence 759	X76488	ANUC	Sequence 811	AA131801	dbEST
Sequence 760	X83544	ANUC	Sequence 812	AA132445	dbEST
Sequence 761	X85134	ANUC	Sequence 813	AA134109	dbEST
Sequence 762	X87949	ANUC	Sequence 814	AA135924	dbEST
Sequence 763	X93036	ANUC	Sequence 815	AA136322	dbEST
Sequence 764	X99699	ANUC	Sequence 816	AA143034	dbEST
Sequence 765	X99920	ANUC	Sequence 817	AA150057	dbEST
Sequence 766	Y09267	ANUC	Sequence 818	AA151651	dbEST
Sequence 767	Y13323	ANUC	Sequence 819	AA156335	dbEST
Sequence 768	Y17392	ANUC	Sequence 820	AA157333	dbEST
Sequence 769	Z12830	ANUC	Sequence 821	AA158987	dbEST
Sequence 770	Z36815	ANUC	Sequence 822	AA165439	dbEST
Sequence 771	Z47087	ANUC	Sequence 823	AA165632	dbEST
Sequence 772	Z48570	ANUC	Sequence 824	AA166618	dbEST
Sequence 773	Z71389	ANUC	Sequence 825	AA172067	dbEST
Sequence 774	AA002223	dbEST	Sequence 826	AA173031	dbEST
Sequence 775	AA018843	dbEST	Sequence 827	AA178870	dbEST
Sequence 776	AA021647	dbEST	Sequence 828	AA181874	dbEST
Sequence 777	AA022842	dbEST	Sequence 829	AA195194	dbEST
Sequence 778	AA022965	dbEST	Sequence 830	AA203206	dbEST
Sequence 779	AA024522	dbEST	Sequence 831	AA203289	dbEST

TABLE 1A

Sequence 832	AA204768	dbEST	Sequence 884	AA447735	dbEST
Sequence 833	AA206621	dbEST	Sequence 885	AA449054	dbEST
Sequence 834	AA213914	dbEST	Sequence 886	AA449205	dbEST
Sequence 835	AA218919	dbEST	Sequence 887	AA449520	dbEST
Sequence 836	AA224050	dbEST	Sequence 888	AA452273	dbEST
Sequence 837	AA224244	dbEST	Sequence 889	AA455007	dbEST
Sequence 838	AA227596	dbEST	Sequence 890	AA455104	dbEST
Sequence 839	AA229018	dbEST	Sequence 891	AA459527	dbEST
Sequence 840	AA229161	dbEST	Sequence 892	AA460226	dbEST
Sequence 841	AA236445	dbEST	Sequence 893	AA461287	dbEST
Sequence 842	AA236680	dbEST	Sequence 894	AA464526	dbEST
Sequence 843	AA243537	dbEST	Sequence 895	AA468398	dbEST
Sequence 844	AA252436	dbEST	Sequence 896	AA469135	dbEST
Sequence 845	AA252869	dbEST	Sequence 897	AA469453	dbEST
Sequence 846	AA256330	dbEST	Sequence 898	AA470690	dbEST
Sequence 847	AA262700	dbEST	Sequence 899	AA479427	dbEST
Sequence 848	AA278358	dbEST	Sequence 900	AA480336	dbEST
Sequence 849	AA287076	dbEST	Sequence 901	AA483454	dbEST
Sequence 850	AA291551	dbEST	Sequence 902	AA487669	dbEST
Sequence 851	AA293273	dbEST	Sequence 903	AA488423	dbEST
Sequence 852	AA295982	dbEST	Sequence 904	AA488635	dbEST
Sequence 853	AA301675	dbEST	Sequence 905	AA488843	dbEST
Sequence 854	AA301722	dbEST	Sequence 906	AA489772	dbEST
Sequence 855	AA302964	dbEST	Sequence 907	AA503972	dbEST
Sequence 856	AA303199	dbEST	Sequence 908	AA508506	dbEST
Sequence 857	AA304927	dbEST	Sequence 909	AA513550	dbEST
Sequence 858	AA305042	dbEST	Sequence 910	AA513783	dbEST
Sequence 859	AA305635	dbEST	Sequence 911	AA514989	dbEST
Sequence 860	AA315030	dbEST	Sequence 912	AA516400	dbEST
Sequence 861	AA315943	dbEST	Sequence 913	AA520993	dbEST
Sequence 862	AA317144	dbEST	Sequence 914	AA521110	dbEST
Sequence 863	AA326060	dbEST	Sequence 915	AA523639	dbEST
Sequence 864	AA327358	dbEST	Sequence 916	AA523697	dbEST
Sequence 865	AA336387	dbEST	Sequence 917	AA528106	dbEST
Sequence 866	AA346413	dbEST	Sequence 918	AA528190	dbEST
Sequence 867	AA352580	dbEST	Sequence 919	AA528226	dbEST
Sequence 868	AA363162	dbEST	Sequence 920	AA534830	dbEST
Sequence 869	AA375754	dbEST	Sequence 921	AA548722	dbEST
Sequence 870	AA399230	dbEST	Sequence 922	AA551236	dbEST
Sequence 871	AA400249	dbEST	Sequence 923	AA551243	dbEST
Sequence 872	AA401629	dbEST	Sequence 924	AA558778	dbEST
Sequence 873	AA402885	dbEST	Sequence 925	AA563834	dbEST
Sequence 874	AA406401	dbEST	Sequence 926	AA576432	dbEST
Sequence 875	AA421682	dbEST	Sequence 927	AA580069	dbEST
Sequence 876	AA422057	dbEST	Sequence 928	AA580294	dbEST
Sequence 877	AA424445	dbEST	Sequence 929	AA582588	dbEST
Sequence 878	AA424901	dbEST	Sequence 930	AA584304	dbEST
Sequence 879	AA424984	dbEST	Sequence 931	AA588772	dbEST
Sequence 880	AA425182	dbEST	Sequence 932	AA593075	dbEST
Sequence 881	AA428607	dbEST	Sequence 933	AA595585	dbEST
Sequence 882	AA446099	dbEST	Sequence 934	AA601895	dbEST
Sequence 883	AA446403	dbEST	Sequence 935	AA628700	dbEST

TABLE 1A

Sequence 936	AA630326	dbEST	Sequence 988	AA902644	dbEST
Sequence 937	AA630642	dbEST	Sequence 989	AA909144	dbEST
Sequence 938	AA631178	dbEST	Sequence 990	AA913281	dbEST
Sequence 939	AA631218	dbEST	Sequence 991	AA916756	dbEST
Sequence 940	AA633550	dbEST	Sequence 992	AA922420	dbEST
Sequence 941	AA634808	dbEST	Sequence 993	AA927283	dbEST
Sequence 942	AA639199	dbEST	Sequence 994	AA933075	dbEST
Sequence 943	AA639791	dbEST	Sequence 995	AA935979	dbEST
Sequence 944	AA644273	dbEST	Sequence 996	AA937947	dbEST
Sequence 945	AA648897	dbEST	Sequence 997	AA948295	dbEST
Sequence 946	AA664732	dbEST	Sequence 998	AA969131	dbEST
Sequence 947	AA677550	dbEST	Sequence 999	AA971881	dbEST
Sequence 948	AA687308	dbEST	Sequence 1000	AA973019	dbEST
Sequence 949	AA705002	dbEST	Sequence 1001	AA988923	dbEST
Sequence 950	AA706685	dbEST	Sequence 1002	AA989465	dbEST
Sequence 951	AA708266	dbEST	Sequence 1003	AA994023	dbEST
Sequence 952	AA713687	dbEST	Sequence 1004	AB002310	ANUC
Sequence 953	AA719618	dbEST	Sequence 1005	AB002330	ANUC
Sequence 954	AA719674	dbEST	Sequence 1006	AB007944	ANUC
Sequence 955	AA720572	dbEST	Sequence 1007	AB012911	ANUC
Sequence 956	AA721752	dbEST	Sequence 1008	AB017019	ANUC
Sequence 957	AA723612	dbEST	Sequence 1009	AB018266	ANUC
Sequence 958	AA730571	dbEST	Sequence 1010	AB018305	ANUC
Sequence 959	AA742282	dbEST	Sequence 1011	AB018347	ANUC
Sequence 960	AA748437	dbEST	Sequence 1012	AB019568	ANUC
Sequence 961	AA749187	dbEST	Sequence 1013	AB023158	ANUC
Sequence 962	AA761602	dbEST	Sequence 1014	AB028976	ANUC
Sequence 963	AA768355	dbEST	Sequence 1015	AB029005	ANUC
Sequence 964	AA769127	dbEST	Sequence 1016	AC28164	PREPATNUC
Sequence 965	AA774030	dbEST	Sequence 1017	AD001528	ANUC
Sequence 966	AA774247	dbEST	Sequence 1018	AF000231	ANUC
Sequence 967	AA779631	dbEST	Sequence 1019	AF006088	ANUC
Sequence 968	AA808747	dbEST	Sequence 1020	AF006516	ANUC
Sequence 969	AA809854	dbEST	Sequence 1021	AF012072	ANUC
Sequence 970	AA810859	dbEST	Sequence 1022	AF026947	ANUC
Sequence 971	AA825673	dbEST	Sequence 1023	AF028832	ANUC
Sequence 972	AA825768	dbEST	Sequence 1024	AF030424	ANUC
Sequence 973	AA826517	dbEST	Sequence 1025	AF031379	ANUC
Sequence 974	AA827331	dbEST	Sequence 1026	AF035287	ANUC
Sequence 975	AA827764	dbEST	Sequence 1027	AF035309	ANUC
Sequence 976	AA829511	dbEST	Sequence 1028	AF038197	ANUC
Sequence 977	AA831603	dbEST	Sequence 1029	AF038404	ANUC
Sequence 978	AA836991	dbEST	Sequence 1030	AF043431	ANUC
Sequence 979	AA837254	dbEST	Sequence 1031	AF044670	ANUC
Sequence 980	AA846480	dbEST	Sequence 1032	AF044958	ANUC
Sequence 981	AA846840	dbEST	Sequence 1033	AF047184	ANUC
Sequence 982	AA853515	dbEST	Sequence 1034	AF052164	ANUC
Sequence 983	AA883212	dbEST	Sequence 1035	AF052496	dbEST
Sequence 984	AA886885	dbEST	Sequence 1036	AF052578	ANUC
Sequence 985	AA889485	dbEST	Sequence 1037	AF054990	ANUC
Sequence 986	AA897461	dbEST	Sequence 1038	AF059524	ANUC
Sequence 987	AA902582	dbEST	Sequence 1039	AF070561	ANUC



TABLE 1A

Sequence 1040	AF070626	ANUC	Sequence 1092	AI373032	dbEST
Sequence 1041	AF070655	ANUC	Sequence 1093	AI374954	dbEST
Sequence 1042	AF070674	ANUC	Sequence 1094	AI380539	dbEST
Sequence 1043	AF075040	ANUC	Sequence 1095	AI417583	dbEST
Sequence 1044	AF077030	ANUC	Sequence 1096	AI432644	dbEST
Sequence 1045	AF078847	ANUC	Sequence 1097	AI433157	dbEST
Sequence 1046	AF080246	ANUC	Sequence 1098	AI457792	dbEST
Sequence 1047	AF081282	ANUC	Sequence 1099	AI469112	dbEST
Sequence 1048	AF081484	ANUC	Sequence 1100	AI471114	dbEST
Sequence 1049	AF084523	ANUC	Sequence 1101	AI471534	dbEST
Sequence 1050	AF086163	ANUC	Sequence 1102	AI473927	dbEST
Sequence 1051	AF095791	ANUC	Sequence 1103	AI479305	dbEST
Sequence 1052	AF100756	ANUC	Sequence 1104	AI499243	dbEST
Sequence 1053	AF107406	ANUC	Sequence 1105	AI525796	dbEST
Sequence 1054	AF119297	ANUC	Sequence 1106	AI525843	dbEST
Sequence 1055	AF131858	ANUC	Sequence 1107	AI537677	dbEST
Sequence 1056	AF132940	ANUC	Sequence 1108	AI541029	dbEST
Sequence 1057	AF151857	ANUC	Sequence 1109	AI560129	dbEST
Sequence 1058	AI028733	dbEST	Sequence 1110	AI583108	dbEST
Sequence 1059	AI031901	dbEST	Sequence 1111	AI584068	dbEST
Sequence 1060	AI033739	dbEST	Sequence 1112	AI587208	dbEST
Sequence 1061	AI040324	dbEST	Sequence 1113	AI589867	dbEST
Sequence 1062	AI051172	dbEST	Sequence 1114	AI610676	dbEST
Sequence 1063	AI076805	dbEST	Sequence 1115	AI630362	dbEST
Sequence 1064	AI087005	dbEST	Sequence 1116	AI633006	dbEST
Sequence 1065	AI089913	dbEST	Sequence 1117	AI634443	dbEST
Sequence 1066	AI092007	dbEST	Sequence 1118	AI635096	dbEST
Sequence 1067	AI127326	dbEST	Sequence 1119	AI682105	dbEST
Sequence 1068	AI147251	dbEST	Sequence 1120	AI683338	dbEST
Sequence 1069	AI148933	dbEST	Sequence 1121	AI684800	dbEST
Sequence 1070	AI149846	dbEST	Sequence 1122	AI684991	dbEST
Sequence 1071	AI167855	dbEST	Sequence 1123	AI689369	dbEST
Sequence 1072	AI183965	dbEST	Sequence 1124	AI689617	dbEST
Sequence 1073	AI189258	dbEST	Sequence 1125	AI689883	dbEST
Sequence 1074	AI220148	dbEST	Sequence 1126	AI693745	dbEST
Sequence 1075	AI224374	dbEST	Sequence 1127	AI701001	dbEST
Sequence 1076	AI240095	dbEST	Sequence 1128	AI733038	dbEST
Sequence 1077	AI246677	dbEST	Sequence 1129	AI735638	dbEST
Sequence 1078	AI248538	dbEST	Sequence 1130	AI741506	dbEST
Sequence 1079	AI266582	dbEST	Sequence 1131	AI742722	dbEST
Sequence 1080	AI268864	dbEST	Sequence 1132	AI742738	dbEST
Sequence 1081	AI270183	dbEST	Sequence 1133	AI743552	dbEST
Sequence 1082	AI271795	dbEST	Sequence 1134	AI753784	dbEST
Sequence 1083	AI273008	dbEST	Sequence 1135	AI754296	dbEST
Sequence 1084	AI273841	dbEST	Sequence 1136	AI754652	dbEST
Sequence 1085	AI274756	dbEST	Sequence 1137	AI754732	dbEST
Sequence 1086	AI275528	dbEST	Sequence 1138	AI765975	dbEST
Sequence 1087	AI283096	dbEST	Sequence 1139	AI769970	dbEST
Sequence 1088	AI298059	dbEST	Sequence 1140	AI819225	dbEST
Sequence 1089	AI335653	dbEST	Sequence 1141	AI820563	dbEST
Sequence 1090	AI338977	dbEST	Sequence 1142	AI827818	dbEST
Sequence 1091	AI339946	dbEST	Sequence 1143	AI828682	dbEST



TABLE 1A

Sequence 1144	AI830067	dbEST	Sequence 1196	D51497	dbEST
Sequence 1145	AI861989	dbEST	Sequence 1197	D53031	dbEST
Sequence 1146	AI887129	dbEST	Sequence 1198	D62116	dbEST
Sequence 1147	AI887632	dbEST	Sequence 1199	D63878	ANUC
Sequence 1148	AI890281	dbEST	Sequence 1200	D78611	ANUC
Sequence 1149	AI924046	dbEST	Sequence 1201	D82348	ANUC
Sequence 1150	AI924096	dbEST	Sequence 1202	D83032	ANUC
Sequence 1151	AI924823	dbEST	Sequence 1203	D85433	ANUC
Sequence 1152	AI963471	dbEST	Sequence 1204	D87437	ANUC
Sequence 1153	AI963604	dbEST	Sequence 1205	D87667	ANUC
Sequence 1154	AI972556	dbEST	Sequence 1206	D89092	ANUC
Sequence 1155	AI979048	dbEST	Sequence 1207	D90041	ANUC
Sequence 1156	AI984656	dbEST	Sequence 1208	E02628	ANUC
Sequence 1157	AJ010442	ANUC	Sequence 1209	E05732	ANUC
Sequence 1158	AJ132694	ANUC	Sequence 1210	F00551	dbEST
Sequence 1159	AJ224442	ANUC	Sequence 1211	H08920	dbEST
Sequence 1160	AL036299	dbEST	Sequence 1212	H25080	dbEST
Sequence 1161	AL042979	dbEST	Sequence 1213	H30306	dbEST
Sequence 1162	AL047305	dbEST	Sequence 1214	H44647	dbEST
Sequence 1163	AL049247	ANUC	Sequence 1215	H81376	dbEST
Sequence 1164	AL049313	ANUC	Sequence 1216	H93521	dbEST
Sequence 1165	AL049381	ANUC	Sequence 1217	H94496	dbEST
Sequence 1166	AL049932	ANUC	Sequence 1218	J03464	ANUC
Sequence 1167	AL050041	ANUC	Sequence 1219	J03799	ANUC
Sequence 1168	AL050161	ANUC	Sequence 1220	J04027	ANUC
Sequence 1169	AL050265	ANUC	Sequence 1221	J04177	ANUC
Sequence 1170	AL050268	ANUC	Sequence 1222	K01228	ANUC
Sequence 1171	AL050367	ANUC	Sequence 1223	K01566	ANUC
Sequence 1172	AL079286	ANUC	Sequence 1224	L07395	ANUC
Sequence 1173	AL079312	ANUC	Sequence 1225	L09159	ANUC
Sequence 1174	AL079314	ANUC	Sequence 1226	L11315	ANUC
Sequence 1175	AL080113	ANUC	Sequence 1227	L13806	ANUC
Sequence 1176	AL110164	ANUC	Sequence 1228	L15702	ANUC
Sequence 1177	AL117412	ANUC	Sequence 1229	L16510	ANUC
Sequence 1178	AL117612	ANUC	Sequence 1230	L24804	ANUC
Sequence 1179	AL119009	dbEST	Sequence 1231	L25931	ANUC
Sequence 1180	AW014693	dbEST	Sequence 1232	L28809	ANUC
Sequence 1181	AW014985	dbEST	Sequence 1233	M10036	ANUC
Sequence 1182	AW021794	dbEST	Sequence 1234	M10905	ANUC
Sequence 1183	C01521	dbEST	Sequence 1235	M11353	ANUC
Sequence 1184	D01096	ANUC	Sequence 1236	M12267	ANUC
Sequence 1185	D13119	ANUC	Sequence 1237	M13536	ANUC
Sequence 1186	D13627	ANUC	Sequence 1238	M14483	ANUC
Sequence 1187	D13630	ANUC	Sequence 1239	M14630	ANUC
Sequence 1188	D13639	ANUC	Sequence 1240	M17885	ANUC
Sequence 1189	D13665	ANUC	Sequence 1241	M18366	ANUC
Sequence 1190	D14530	ANUC	Sequence 1242	M21575	ANUC
Sequence 1191	D21260	ANUC	Sequence 1243	M23254	ANUC
Sequence 1192	D25278	ANUC	Sequence 1244	M24194	ANUC
Sequence 1193	D26361	ANUC	Sequence 1245	M24486	ANUC
Sequence 1194	D30655	ANUC	Sequence 1246	M26512	ANUC
Sequence 1195	D50310	ANUC	Sequence 1247	M28372	ANUC

TABLE 1A

Sequence 1248	M31159	ANUC	Sequence 1300	V84510	NUCPATENT
Sequence 1249	M32220	ANUC	Sequence 1301	W19427	dbEST
Sequence 1250	M36341	ANUC	Sequence 1302	W65357	dbEST
Sequence 1251	M36693	ANUC	Sequence 1303	W75963	dbEST
Sequence 1252	M38690	ANUC	Sequence 1304	W80525	dbEST
Sequence 1253	M58485	ANUC	Sequence 1305	X01630	ANUC
Sequence 1254	M59849	ANUC	Sequence 1306	X04098	ANUC
Sequence 1255	M62831	ANUC	Sequence 1307	X04408	ANUC
Sequence 1256	M64241	ANUC	Sequence 1308	X06700	ANUC
Sequence 1257	M69043	ANUC	Sequence 1309	X14420	ANUC
Sequence 1258	M77142	ANUC	Sequence 1310	X51742	NUCPATENT
Sequence 1259	M77830	ANUC	Sequence 1311	X60111	ANUC
Sequence 1260	M86667	ANUC	Sequence 1312	X69398	ANUC
Sequence 1261	M88108	ANUC	Sequence 1313	X72755	ANUC
Sequence 1262	M93651	ANUC	Sequence 1314	X74979	ANUC
Sequence 1263	M95542	ANUC	Sequence 1315	X76180	ANUC
Sequence 1264	N43970	dbEST	Sequence 1316	X78627	ANUC
Sequence 1265	Q12759	NUCPATENT	Sequence 1317	X79067	ANUC
Sequence 1266	Q14635	NUCPATENT	Sequence 1318	X80910	ANUC
Sequence 1267	R11045	dbEST	Sequence 1319	X87949	ANUC
Sequence 1268	R76376	dbEST	Sequence 1320	Y00052	ANUC
Sequence 1269	R84450	dbEST	Sequence 1321	Y00062	ANUC
Sequence 1270	S74728	ANUC	Sequence 1322	Y00282	ANUC
Sequence 1271	S82081	ANUC	Sequence 1323	Y00503	ANUC
Sequence 1272	T07459	dbEST	Sequence 1324	Y15286	ANUC
Sequence 1273	T19883	NUCPATENT	Sequence 1325	Y17171	ANUC
Sequence 1274	T21168	NUCPATENT	Sequence 1326	Z13009	ANUC
Sequence 1275	T22605	NUCPATENT	Sequence 1327	Z24724	ANUC
Sequence 1276	T37405	NUCPATENT	Sequence 1328	Z29083	ANUC
Sequence 1277	T67129	dbEST	Sequence 1329	Z29331	ANUC
Sequence 1278	T69703	dbEST	Sequence 1330	Z46606	ANUC
Sequence 1279	T78615	dbEST	Sequence 1331	Z48501	ANUC
Sequence 1280	T89937	dbEST	Sequence 1332	AA001460	dbEST
Sequence 1281	U03851	ANUC	Sequence 1333	AA001543	dbEST
Sequence 1282	U12404	ANUC	Sequence 1334	AA001792	dbEST
Sequence 1283	U14967	ANUC	Sequence 1335	AA004925	dbEST
Sequence 1284	U14971	ANUC	Sequence 1336	AA010897	dbEST
Sequence 1285	U20659	ANUC	Sequence 1337	AA017162	dbEST
Sequence 1286	U25789	ANUC	Sequence 1338	AA019019	dbEST
Sequence 1287	U30825	ANUC	Sequence 1339	AA022980	dbEST
Sequence 1288	U47077	ANUC	Sequence 1340	AA024595	dbEST
Sequence 1289	U49844	ANUC	Sequence 1341	AA024940	dbEST
Sequence 1290	U63846	ANUC	Sequence 1342	AA024996	dbEST
Sequence 1291	U65928	ANUC	Sequence 1343	AA025750	dbEST
Sequence 1292	U72516	ANUC	Sequence 1344	AA026598	dbEST
Sequence 1293	U79282	ANUC	Sequence 1345	AA029271	dbEST
Sequence 1294	U90716	ANUC	Sequence 1346	AA029725	dbEST
Sequence 1295	U90904	ANUC	Sequence 1347	AA029930	dbEST
Sequence 1296	U94364	ANUC	Sequence 1348	AA033832	dbEST
Sequence 1297	V20437	NUCPATENT	Sequence 1349	AA035471	dbEST
Sequence 1298	V24305	NUCPATENT	Sequence 1350	AA035616	dbEST
Sequence 1299	V81394	NUCPATENT	Sequence 1351	AA036752	dbEST

TABLE 1A

Sequence 1352	AA037377	dbEST	Sequence 1404	AA088693	dbEST
Sequence 1353	AA039778	dbEST	Sequence 1405	AA088783	dbEST
Sequence 1354	AA039948	dbEST	Sequence 1406	AA088829	dbEST
Sequence 1355	AA040688	dbEST	Sequence 1407	AA090106	dbEST
Sequence 1356	AA040820	dbEST	Sequence 1408	AA096032	dbEST
Sequence 1357	AA041259	dbEST	Sequence 1409	AA099819	dbEST
Sequence 1358	AA043477	dbEST	Sequence 1410	AA099923	dbEST
Sequence 1359	AA044209	dbEST	Sequence 1411	AA099976	dbEST
Sequence 1360	AA044233	dbEST	Sequence 1412	AA100764	dbEST
Sequence 1361	AA044791	dbEST	Sequence 1413	AA101010	dbEST
Sequence 1362	AA045054	dbEST	Sequence 1414	AA102013	dbEST
Sequence 1363	AA045147	dbEST	Sequence 1415	AA102564	dbEST
Sequence 1364	AA045768	dbEST	Sequence 1416	AA102830	dbEST
Sequence 1365	AA046848	dbEST	Sequence 1417	AA112186	dbEST
Sequence 1366	AA053021	dbEST	Sequence 1418	AA112645	dbEST
Sequence 1367	AA053316	dbEST	Sequence 1419	AA113305	dbEST
Sequence 1368	AA053919	dbEST	Sequence 1420	AA115218	dbEST
Sequence 1369	AA054069	dbEST	Sequence 1421	AA115315	dbEST
Sequence 1370	AA055479	dbEST	Sequence 1422	AA121656	dbEST
Sequence 1371	AA055591	dbEST	Sequence 1423	AA121718	dbEST
Sequence 1372	AA055637	dbEST	Sequence 1424	AA125809	dbEST
Sequence 1373	AA057243	dbEST	Sequence 1425	AA125939	dbEST
Sequence 1374	AA058712	dbEST	Sequence 1426	AA126452	dbEST
Sequence 1375	AA059128	dbEST	Sequence 1427	AA126718	dbEST
Sequence 1376	AA065169	dbEST	Sequence 1428	AA127436	dbEST
Sequence 1377	AA069850	dbEST	Sequence 1429	AA127666	dbEST
Sequence 1378	AA071167	dbEST	Sequence 1430	AA128063	dbEST
Sequence 1379	AA075158	dbEST	Sequence 1431	AA128636	dbEST
Sequence 1380	AA075515	dbEST	Sequence 1432	AA128641	dbEST
Sequence 1381	AA075663	dbEST	Sequence 1433	AA130778	dbEST
Sequence 1382	AA076397	dbEST	Sequence 1434	AA130982	dbEST
Sequence 1383	AA076421	dbEST	Sequence 1435	AA131827	dbEST
Sequence 1384	AA078387	dbEST	Sequence 1436	AA132056	dbEST
Sequence 1385	AA078570	dbEST	Sequence 1437	AA132163	dbEST
Sequence 1386	AA078872	dbEST	Sequence 1438	AA132574	dbEST
Sequence 1387	AA079480	dbEST	Sequence 1439	AA132992	dbEST
Sequence 1388	AA080889	dbEST	Sequence 1440	AA133351	dbEST
Sequence 1389	AA081073	dbEST	Sequence 1441	AA133474	dbEST
Sequence 1390	AA081608	dbEST	Sequence 1442	AA134460	dbEST
Sequence 1391	AA081834	dbEST	Sequence 1443	AA134527	dbEST
Sequence 1392	AA081917	dbEST	Sequence 1444	AA134589	dbEST
Sequence 1393	AA082258	dbEST	Sequence 1445	AA135696	dbEST
Sequence 1394	AA082441	dbEST	Sequence 1446	AA137017	dbEST
Sequence 1395	AA083270	dbEST	Sequence 1447	AA142941	dbEST
Sequence 1396	AA083345	dbEST	Sequence 1448	AA143001	dbEST
Sequence 1397	AA083522	dbEST	Sequence 1449	AA143074	dbEST
Sequence 1398	AA083573	dbEST	Sequence 1450	AA143746	dbEST
Sequence 1399	AA083638	dbEST	Sequence 1451	AA146900	dbEST
Sequence 1400	AA083774	dbEST	Sequence 1452	AA147200	dbEST
Sequence 1401	AA088318	dbEST	Sequence 1453	AA147247	dbEST
Sequence 1402	AA088344	dbEST	Sequence 1454	AA147781	dbEST
Sequence 1403	AA088351	dbEST	Sequence 1455	AA148027	dbEST

TABLE 1A

Sequence 1456	AA148136	dbEST	Sequence 1508	AA191422	dbEST
Sequence 1457	AA149810	dbEST	Sequence 1509	AA192094	dbEST
Sequence 1458	AA150377	dbEST	Sequence 1510	AA193308	dbEST
Sequence 1459	AA150837	dbEST	Sequence 1511	AA194577	dbEST
Sequence 1460	AA150928	dbEST	Sequence 1512	AA195246	dbEST
Sequence 1461	AA151274	dbEST	Sequence 1513	AA195865	dbEST
Sequence 1462	AA151594	dbEST	Sequence 1514	AA196424	dbEST
Sequence 1463	AA151755	dbEST	Sequence 1515	AA196982	dbEST
Sequence 1464	AA152476	dbEST	Sequence 1516	AA203691	dbEST
Sequence 1465	AA155754	dbEST	Sequence 1517	AA204867	dbEST
Sequence 1466	AA156066	dbEST	Sequence 1518	AA206578	dbEST
Sequence 1467	AA157163	dbEST	Sequence 1519	AA206991	dbEST
Sequence 1468	AA157993	dbEST	Sequence 1520	AA209508	dbEST
Sequence 1469	AA158738	dbEST	Sequence 1521	AA216753	dbEST
Sequence 1470	AA159110	dbEST	Sequence 1522	AA219665	dbEST
Sequence 1471	AA159576	dbEST	Sequence 1523	AA223121	dbEST
Sequence 1472	AA161003	dbEST	Sequence 1524	AA223820	dbEST
Sequence 1473	AA161076	dbEST	Sequence 1525	AA224109	dbEST
Sequence 1474	AA161467	dbEST	Sequence 1526	AA224407	dbEST
Sequence 1475	AA164193	dbEST	Sequence 1527	AA227118	dbEST
Sequence 1476	AA164473	dbEST	Sequence 1528	AA229325	dbEST
Sequence 1477	AA164729	dbEST	Sequence 1529	AA229611	dbEST
Sequence 1478	AA164873	dbEST	Sequence 1530	AA232959	dbEST
Sequence 1479	AA165027	dbEST	Sequence 1531	AA233835	dbEST
Sequence 1480	AA165068	dbEST	Sequence 1532	AA233843	dbEST
Sequence 1481	AA165087	dbEST	Sequence 1533	AA234092	dbEST
Sequence 1482	AA165174	dbEST	Sequence 1534	AA234307	dbEST
Sequence 1483	AA165282	dbEST	Sequence 1535	AA236776	dbEST
Sequence 1484	AA165293	dbEST	Sequence 1536	AA242985	dbEST
Sequence 1485	AA165638	dbEST	Sequence 1537	AA243338	dbEST
Sequence 1486	AA166618	dbEST	Sequence 1538	AA244342	dbEST
Sequence 1487	AA167041	dbEST	Sequence 1539	AA249154	dbEST
Sequence 1488	AA167750	dbEST	Sequence 1540	AA255502	dbEST
Sequence 1489	AA171630	dbEST	Sequence 1541	AA256591	dbEST
Sequence 1490	AA173506	dbEST	Sequence 1542	AA261990	dbEST
Sequence 1491	AA174097	dbEST	Sequence 1543	AA262939	dbEST
Sequence 1492	AA179187	dbEST	Sequence 1544	AA278445	dbEST
Sequence 1493	AA180137	dbEST	Sequence 1545	AA278482	dbEST
Sequence 1494	AA180224	dbEST	Sequence 1546	AA278642	dbEST
Sequence 1495	AA180383	dbEST	Sequence 1547	AA278956	dbEST
Sequence 1496	AA181075	dbEST	Sequence 1548	AA279048	dbEST
Sequence 1497	AA181258	dbEST	Sequence 1549	AA280099	dbEST
Sequence 1498	AA181684	dbEST	Sequence 1550	AA280221	dbEST
Sequence 1499	AA182415	dbEST	Sequence 1551	AA280828	dbEST
Sequence 1500	AA182540	dbEST	Sequence 1552	AA282915	dbEST
Sequence 1501	AA186577	dbEST	Sequence 1553	AA284334	dbEST
Sequence 1502	AA187817	dbEST	Sequence 1554	AA284555	dbEST
Sequence 1503	AA188045	dbEST	Sequence 1555	AA284670	dbEST
Sequence 1504	AA188140	dbEST	Sequence 1556	AA284671	dbEST
Sequence 1505	AA188384	dbEST	Sequence 1557	AA284870	dbEST
Sequence 1506	AA188826	dbEST	Sequence 1558	AA284906	dbEST
Sequence 1507	AA190873	dbEST	Sequence 1559	AA285290	dbEST

TABLE 1A

Sequence 1560	AA286699	dbEST	Sequence 1612	AA315379	dbEST
Sequence 1561	AA286872	dbEST	Sequence 1613	AA317243	dbEST
Sequence 1562	AA287219	dbEST	Sequence 1614	AA317393	dbEST
Sequence 1563	AA287642	dbEST	Sequence 1615	AA318969	dbEST
Sequence 1564	AA287815	dbEST	Sequence 1616	AA327201	dbEST
Sequence 1565	AA291438	dbEST	Sequence 1617	AA331991	dbEST
Sequence 1566	AA291485	dbEST	Sequence 1618	AA332672	dbEST
Sequence 1567	AA291971	dbEST	Sequence 1619	AA333358	dbEST
Sequence 1568	AA292334	dbEST	Sequence 1620	AA335273	dbEST
Sequence 1569	AA293127	dbEST	Sequence 1621	AA336666	dbEST
Sequence 1570	AA293133	dbEST	Sequence 1622	AA337192	dbEST
Sequence 1571	AA293273	dbEST	Sequence 1623	AA337489	dbEST
Sequence 1572	AA293286	dbEST	Sequence 1624	AA338793	dbEST
Sequence 1573	AA293353	dbEST	Sequence 1625	AA339957	dbEST
Sequence 1574	AA293572	dbEST	Sequence 1626	AA340341	dbEST
Sequence 1575	AA293629	dbEST	Sequence 1627	AA341446	dbEST
Sequence 1576	AA293759	dbEST	Sequence 1628	AA341465	dbEST
Sequence 1577	AA293804	dbEST	Sequence 1629	AA342969	dbEST
Sequence 1578	AA296780	dbEST	Sequence 1630	AA343629	dbEST
Sequence 1579	AA297402	dbEST	Sequence 1631	AA344084	dbEST
Sequence 1580	AA298505	dbEST	Sequence 1632	AA345329	dbEST
Sequence 1581	AA299640	dbEST	Sequence 1633	AA346393	dbEST
Sequence 1582	AA301062	dbEST	Sequence 1634	AA346698	dbEST
Sequence 1583	AA301800	dbEST	Sequence 1635	AA347887	dbEST
Sequence 1584	AA303461	dbEST	Sequence 1636	AA350059	dbEST
Sequence 1585	AA303568	dbEST	Sequence 1637	AA351507	dbEST
Sequence 1586	AA306718	dbEST	Sequence 1638	AA355003	dbEST
Sequence 1587	AA306862	dbEST	Sequence 1639	AA356682	dbEST
Sequence 1588	AA306876	dbEST	Sequence 1640	AA357574	dbEST
Sequence 1589	AA307198	dbEST	Sequence 1641	AA358887	dbEST
Sequence 1590	AA307325	dbEST	Sequence 1642	AA359705	dbEST
Sequence 1591	AA308065	dbEST	Sequence 1643	AA364352	dbEST
Sequence 1592	AA308274	dbEST	Sequence 1644	AA367451	dbEST
Sequence 1593	AA308744	dbEST	Sequence 1645	AA367773	dbEST
Sequence 1594	AA310739	dbEST	Sequence 1646	AA368542	dbEST
Sequence 1595	AA310771	dbEST	Sequence 1647	AA369400	dbEST
Sequence 1596	AA311228	dbEST	Sequence 1648	AA373230	dbEST
Sequence 1597	AA311460	dbEST	Sequence 1649	AA374754	dbEST
Sequence 1598	AA311571	dbEST	Sequence 1650	AA375312	dbEST
Sequence 1599	AA311801	dbEST	Sequence 1651	AA375815	dbEST
Sequence 1600	AA311848	dbEST	Sequence 1652	AA393525	dbEST
Sequence 1601	AA311905	dbEST	Sequence 1653	AA394115	dbEST
Sequence 1602	AA312218	dbEST	Sequence 1654	AA398443	dbEST
Sequence 1603	AA312240	dbEST	Sequence 1655	AA398585	dbEST
Sequence 1604	AA312435	dbEST	Sequence 1656	AA398739	dbEST
Sequence 1605	AA313108	dbEST	Sequence 1657	AA399165	dbEST
Sequence 1606	AA313223	dbEST	Sequence 1658	AA399628	dbEST
Sequence 1607	AA313653	dbEST	Sequence 1659	AA401329	dbEST
Sequence 1608	AA313994	dbEST	Sequence 1660	AA401334	dbEST
Sequence 1609	AA314431	dbEST	Sequence 1661	AA402191	dbEST
Sequence 1610	AA314872	dbEST	Sequence 1662	AA402289	dbEST
Sequence 1611	AA315363	dbEST	Sequence 1663	AA402775	dbEST

TABLE 1A

Sequence 1664	AA403319	dbEST	Sequence 1716	AA477567	dbEST
Sequence 1665	AA404613	dbEST	Sequence 1717	AA477973	dbEST
Sequence 1666	AA405124	dbEST	Sequence 1718	AA478230	dbEST
Sequence 1667	AA406239	dbEST	Sequence 1719	AA479646	dbEST
Sequence 1668	AA410580	dbEST	Sequence 1720	AA479648	dbEST
Sequence 1669	AA410982	dbEST	Sequence 1721	AA479848	dbEST
Sequence 1670	AA411021	dbEST	Sequence 1722	AA481078	dbEST
Sequence 1671	AA411252	dbEST	Sequence 1723	AA481710	dbEST
Sequence 1672	AA411764	dbEST	Sequence 1724	AA482430	dbEST
Sequence 1673	AA417794	dbEST	Sequence 1725	AA482432	dbEST
Sequence 1674	AA419263	dbEST	Sequence 1726	AA482779	dbEST
Sequence 1675	AA419284	dbEST	Sequence 1727	AA483258	dbEST
Sequence 1676	AA420751	dbEST	Sequence 1728	AA483726	dbEST
Sequence 1677	AA420758	dbEST	Sequence 1729	AA483858	dbEST
Sequence 1678	AA421248	dbEST	Sequence 1730	AA484181	dbEST
Sequence 1679	AA421682	dbEST	Sequence 1731	AA486047	dbEST
Sequence 1680	AA422060	dbEST	Sequence 1732	AA486859	dbEST
Sequence 1681	AA422143	dbEST	Sequence 1733	AA488141	dbEST
Sequence 1682	AA425004	dbEST	Sequence 1734	AA488385	dbEST
Sequence 1683	AA425468	dbEST	Sequence 1735	AA488517	dbEST
Sequence 1684	AA425737	dbEST	Sequence 1736	AA489323	dbEST
Sequence 1685	AA429794	dbEST	Sequence 1737	AA489380	dbEST
Sequence 1686	AA430400	dbEST	Sequence 1738	AA489382	dbEST
Sequence 1687	AA430436	dbEST	Sequence 1739	AA491204	dbEST
Sequence 1688	AA431428	dbEST	Sequence 1740	AA492143	dbEST
Sequence 1689	AA433988	dbEST	Sequence 1741	AA493371	dbEST
Sequence 1690	AA436315	dbEST	Sequence 1742	AA494321	dbEST
Sequence 1691	AA436411	dbEST	Sequence 1743	AA494552	dbEST
Sequence 1692	AA443024	dbEST	Sequence 1744	AA501657	dbEST
Sequence 1693	AA449394	dbEST	Sequence 1745	AA502136	dbEST
Sequence 1694	AA451779	dbEST	Sequence 1746	AA505780	dbEST
Sequence 1695	AA453878	dbEST	Sequence 1747	AA512933	dbEST
Sequence 1696	AA454668	dbEST	Sequence 1748	AA514395	dbEST
Sequence 1697	AA454953	dbEST	Sequence 1749	AA514974	dbEST
Sequence 1698	AA454962	dbEST	Sequence 1750	AA515143	dbEST
Sequence 1699	AA455245	dbEST	Sequence 1751	AA516376	dbEST
Sequence 1700	AA455785	dbEST	Sequence 1752	AA521006	dbEST
Sequence 1701	AA456454	dbEST	Sequence 1753	AA523522	dbEST
Sequence 1702	AA456557	dbEST	Sequence 1754	AA524748	dbEST
Sequence 1703	AA457255	dbEST	Sequence 1755	AA524950	dbEST
Sequence 1704	AA457579	dbEST	Sequence 1756	AA525141	dbEST
Sequence 1705	AA459167	dbEST	Sequence 1757	AA526028	dbEST
Sequence 1706	AA459210	dbEST	Sequence 1758	AA527275	dbEST
Sequence 1707	AA459527	dbEST	Sequence 1759	AA527557	dbEST
Sequence 1708	AA460570	dbEST	Sequence 1760	AA533506	dbEST
Sequence 1709	AA460816	dbEST	Sequence 1761	AA534349	dbEST
Sequence 1710	AA461005	dbEST	Sequence 1762	AA534586	dbEST
Sequence 1711	AA468657	dbEST	Sequence 1763	AA534608	dbEST
Sequence 1712	AA469447	dbEST	Sequence 1764	AA535496	dbEST
Sequence 1713	AA469453	dbEST	Sequence 1765	AA541651	dbEST
Sequence 1714	AA476522	dbEST	Sequence 1766	AA548056	dbEST
Sequence 1715	AA477018	dbEST	Sequence 1767	AA548600	dbEST

TABLE 1A

Sequence 1768	AA550854	dbEST	Sequence 1820	AA630611	dbEST
Sequence 1769	AA550855	dbEST	Sequence 1821	AA631326	dbEST
Sequence 1770	AA551351	dbEST	Sequence 1822	AA633909	dbEST
Sequence 1771	AA551391	dbEST	Sequence 1823	AA634260	dbEST
Sequence 1772	AA554437	dbEST	Sequence 1824	AA634298	dbEST
Sequence 1773	AA554735	dbEST	Sequence 1825	AA640505	dbEST
Sequence 1774	AA555102	dbEST	Sequence 1826	AA641289	dbEST
Sequence 1775	AA564272	dbEST	Sequence 1827	AA644625	dbEST
Sequence 1776	AA564870	dbEST	Sequence 1828	AA648944	dbEST
Sequence 1777	AA565420	dbEST	Sequence 1829	AA651720	dbEST
Sequence 1778	AA568936	dbEST	Sequence 1830	AA652478	dbEST
Sequence 1779	AA569816	dbEST	Sequence 1831	AA652505	dbEST
Sequence 1780	AA569851	dbEST	Sequence 1832	AA653775	dbEST
Sequence 1781	AA569916	dbEST	Sequence 1833	AA658374	dbEST
Sequence 1782	AA573761	dbEST	Sequence 1834	AA663005	dbEST
Sequence 1783	AA573787	dbEST	Sequence 1835	AA669154	dbEST
Sequence 1784	AA577537	dbEST	Sequence 1836	AA677560	dbEST
Sequence 1785	AA578881	dbEST	Sequence 1837	AA677750	dbEST
Sequence 1786	AA579591	dbEST	Sequence 1838	AA678185	dbEST
Sequence 1787	AA579890	dbEST	Sequence 1839	AA678251	dbEST
Sequence 1788	AA580835	dbEST	Sequence 1840	AA687495	dbEST
Sequence 1789	AA582093	dbEST	Sequence 1841	AA703208	dbEST
Sequence 1790	AA582866	dbEST	Sequence 1842	AA703667	dbEST
Sequence 1791	AA583055	dbEST	Sequence 1843	AA703907	dbEST
Sequence 1792	AA583498	dbEST	Sequence 1844	AA704208	dbEST
Sequence 1793	AA583567	dbEST	Sequence 1845	AA706347	dbEST
Sequence 1794	AA583773	dbEST	Sequence 1846	AA714010	dbEST
Sequence 1795	AA584921	dbEST	Sequence 1847	AA715984	dbEST
Sequence 1796	AA586755	dbEST	Sequence 1848	AA716651	dbEST
Sequence 1797	AA587140	dbEST	Sequence 1849	AA719530	dbEST
Sequence 1798	AA587315	dbEST	Sequence 1850	AA721642	dbEST
Sequence 1799	AA587873	dbEST	Sequence 1851	AA729381	dbEST
Sequence 1800	AA593983	dbEST	Sequence 1852	AA731946	dbEST
Sequence 1801	AA594366	dbEST	Sequence 1853	AA736817	dbEST
Sequence 1802	AA595624	dbEST	Sequence 1854	AA742713	dbEST
Sequence 1803	AA595771	dbEST	Sequence 1855	AA743278	dbEST
Sequence 1804	AA599454	dbEST	Sequence 1856	AA744681	dbEST
Sequence 1805	AA600227	dbEST	Sequence 1857	AA745953	dbEST
Sequence 1806	AA600771	dbEST	Sequence 1858	AA759195	dbEST
Sequence 1807	AA601172	dbEST	Sequence 1859	AA767779	dbEST
Sequence 1808	AA602395	dbEST	Sequence 1860	AA769697	dbEST
Sequence 1809	AA602871	dbEST	Sequence 1861	AA773998	dbEST
Sequence 1810	AA603125	dbEST	Sequence 1862	AA775058	dbEST
Sequence 1811	AA603177	dbEST	Sequence 1863	AA776593	dbEST
Sequence 1812	AA604324	dbEST	Sequence 1864	AA777384	dbEST
Sequence 1813	AA604853	dbEST	Sequence 1865	AA778672	dbEST
Sequence 1814	AA610279	dbEST	Sequence 1866	AA779949	dbEST
Sequence 1815	AA610476	dbEST	Sequence 1867	AA781487	dbEST
Sequence 1816	AA610734	dbEST	Sequence 1868	AA788907	dbEST
Sequence 1817	AA614482	dbEST	Sequence 1869	AA806278	dbEST
Sequence 1818	AA628536	dbEST	Sequence 1870	AA806735	dbEST
Sequence 1819	AA628547	dbEST	Sequence 1871	AA808769	dbEST



TABLE 1A

Sequence 1872	AA810149	dbEST	Sequence 1924	AB007916	ANUC
Sequence 1873	AA811609	dbEST	Sequence 1925	AB007923	ANUC
Sequence 1874	AA813604	dbEST	Sequence 1926	AB007957	ANUC
Sequence 1875	AA826307	dbEST	Sequence 1927	AB011103	ANUC
Sequence 1876	AA833766	dbEST	Sequence 1928	AB011143	ANUC
Sequence 1877	AA833900	dbEST	Sequence 1929	AB011151	ANUC
Sequence 1878	AA837457	dbEST	Sequence 1930	AB011166	ANUC
Sequence 1879	AA843531	dbEST	Sequence 1931	AB014533	ANUC
Sequence 1880	AA845737	dbEST	Sequence 1932	AB014542	ANUC
Sequence 1881	AA846698	dbEST	Sequence 1933	AB014560	ANUC
Sequence 1882	AA846856	dbEST	Sequence 1934	AB015630	ANUC
Sequence 1883	AA852896	dbEST	Sequence 1935	AB015856	ANUC
Sequence 1884	AA856902	dbEST	Sequence 1936	AB018281	ANUC
Sequence 1885	AA857824	dbEST	Sequence 1937	AB018284	ANUC
Sequence 1886	AA857882	dbEST	Sequence 1938	AB018285	ANUC
Sequence 1887	AA861665	dbEST	Sequence 1939	AB018289	ANUC
Sequence 1888	AA865960	dbEST	Sequence 1940	AB018305	ANUC
Sequence 1889	AA868529	dbEST	Sequence 1941	AB018327	ANUC
Sequence 1890	AA873271	dbEST	Sequence 1942	AB018331	ANUC
Sequence 1891	AA877189	dbEST	Sequence 1943	AB018337	ANUC
Sequence 1892	AA884922	dbEST	Sequence 1944	AB019409	ANUC
Sequence 1893	AA886453	dbEST	Sequence 1945	AB019563	ANUC
Sequence 1894	AA906652	dbEST	Sequence 1946	AB019568	ANUC
Sequence 1895	AA906865	dbEST	Sequence 1947	AB019691	ANUC
Sequence 1896	AA918993	dbEST	Sequence 1948	AB020682	ANUC
Sequence 1897	AA926926	dbEST	Sequence 1949	AB020718	ANUC
Sequence 1898	AA928934	dbEST	Sequence 1950	AB021288	ANUC
Sequence 1899	AA932501	dbEST	Sequence 1951	AB023154	ANUC
Sequence 1900	AA933987	dbEST	Sequence 1952	AB023219	ANUC
Sequence 1901	AA935947	dbEST	Sequence 1953	AB024704	ANUC
Sequence 1902	AA937302	dbEST	Sequence 1954	AB027467	ANUC
Sequence 1903	AA937773	dbEST	Sequence 1955	AB028069	ANUC
Sequence 1904	AA947835	dbEST	Sequence 1956	AB028624	ANUC
Sequence 1905	AA954939	dbEST	Sequence 1957	AB028969	ANUC
Sequence 1906	AA962587	dbEST	Sequence 1958	AB028986	ANUC
Sequence 1907	AA962632	dbEST	Sequence 1959	AB029000	ANUC
Sequence 1908	AA972525	dbEST	Sequence 1960	AB029004	ANUC
Sequence 1909	AA976489	dbEST	Sequence 1961	AB029028	ANUC
Sequence 1910	AA983380	dbEST	Sequence 1962	AC03044	PREPATNUC
Sequence 1911	AA984586	dbEST	Sequence 1963	AC31479	PREPATNUC
Sequence 1912	AA992596	dbEST	Sequence 1964	AF000670	ANUC
Sequence 1913	AB002305	ANUC	Sequence 1965	AF000974	ANUC
Sequence 1914	AB002330	ANUC	Sequence 1966	AF001893	ANUC
Sequence 1915	AB002357	ANUC	Sequence 1967	AF004562	ANUC
Sequence 1916	AB002806	ANUC	Sequence 1968	AF006043	ANUC
Sequence 1917	AB003476	ANUC	Sequence 1969	AF007135	ANUC
Sequence 1918	AB004066	ANUC	Sequence 1970	AF007151	ANUC
Sequence 1919	AB006077	ANUC	Sequence 1971	AF007170	ANUC
Sequence 1920	AB006534	ANUC	Sequence 1972	AF009615	ANUC
Sequence 1921	AB006755	ANUC	Sequence 1973	AF013759	ANUC
Sequence 1922	AB007867	ANUC	Sequence 1974	AF013988	ANUC
Sequence 1923	AB007900	ANUC	Sequence 1975	AF015283	ANUC



TABLE 1A

Sequence 1976	AF015767	ANUC	Sequence 2028	AF064019	ANUC
Sequence 1977	AF016507	ANUC	Sequence 2029	AF068235	ANUC
Sequence 1978	AF016582	ANUC	Sequence 2030	AF068846	ANUC
Sequence 1979	AF017790	ANUC	Sequence 2031	AF070523	ANUC
Sequence 1980	AF019767	ANUC	Sequence 2032	AF070537	ANUC
Sequence 1981	AF021351	ANUC	Sequence 2033	AF070555	ANUC
Sequence 1982	AF021819	ANUC	Sequence 2034	AF070561	ANUC
Sequence 1983	AF022229	ANUC	Sequence 2035	AF070596	ANUC
Sequence 1984	AF023266	ANUC	Sequence 2036	AF070600	ANUC
Sequence 1985	AF025439	ANUC	Sequence 2037	AF070626	ANUC
Sequence 1986	AF026166	ANUC	Sequence 2038	AF070649	ANUC
Sequence 1987	AF026939	ANUC	Sequence 2039	AF070662	ANUC
Sequence 1988	AF027205	ANUC	Sequence 2040	AF070672	ANUC
Sequence 1989	AF031385	ANUC	Sequence 2041	AF071202	ANUC
Sequence 1990	AF034607	ANUC	Sequence 2042	AF071219	ANUC
Sequence 1991	AF035286	ANUC	Sequence 2043	AF071593	ANUC
Sequence 1992	AF035309	ANUC	Sequence 2044	AF073298	ANUC
Sequence 1993	AF035313	ANUC	Sequence 2045	AF075587	ANUC
Sequence 1994	AF037204	ANUC	Sequence 2046	AF077030	ANUC
Sequence 1995	AF038661	ANUC	Sequence 2047	AF077045	ANUC
Sequence 1996	AF039019	ANUC	Sequence 2048	AF077200	ANUC
Sequence 1997	AF039291	ANUC	Sequence 2049	AF077202	ANUC
Sequence 1998	AF039843	ANUC	Sequence 2050	AF077207	ANUC
Sequence 1999	AF040990	ANUC	Sequence 2051	AF081192	ANUC
Sequence 2000	AF041483	ANUC	Sequence 2052	AF081484	ANUC
Sequence 2001	AF042385	ANUC	Sequence 2053	AF083190	ANUC
Sequence 2002	AF042729	ANUC	Sequence 2054	AF085355	ANUC
Sequence 2003	AF044588	ANUC	Sequence 2055	AF086003	ANUC
Sequence 2004	AF045184	ANUC	Sequence 2056	AF086116	ANUC
Sequence 2005	AF047438	ANUC	Sequence 2057	AF086178	ANUC
Sequence 2006	AF047472	ANUC	Sequence 2058	AF086205	ANUC
Sequence 2007	AF048977	ANUC	Sequence 2059	AF086207	ANUC
Sequence 2008	AF050171	ANUC	Sequence 2060	AF086336	ANUC
Sequence 2009	AF050199	ANUC	Sequence 2061	AF086517	ANUC
Sequence 2010	AF050639	ANUC	Sequence 2062	AF087135	ANUC
Sequence 2011	AF052124	ANUC	Sequence 2063	AF087990	ANUC
Sequence 2012	AF052135	ANUC	Sequence 2064	AF088036	ANUC
Sequence 2013	AF052149	ANUC	Sequence 2065	AF091076	ANUC
Sequence 2014	AF052164	ANUC	Sequence 2066	AF092563	ANUC
Sequence 2015	AF052169	ANUC	Sequence 2067	AF095287	ANUC
Sequence 2016	AF052180	ANUC	Sequence 2068	AF095791	ANUC
Sequence 2017	AF052514	ANUC	Sequence 2069	AF097709	ANUC
Sequence 2018	AF054183	ANUC	Sequence 2070	AF100741	ANUC
Sequence 2019	AF054187	ANUC	Sequence 2071	AF100756	ANUC
Sequence 2020	AF054840	ANUC	Sequence 2072	AF100928	ANUC
Sequence 2021	AF055012	ANUC	Sequence 2073	AF104222	ANUC
Sequence 2022	AF055033	ANUC	Sequence 2074	AF104913	ANUC
Sequence 2023	AF057299	ANUC	Sequence 2075	AF104923	ANUC
Sequence 2024	AF059252	ANUC	Sequence 2076	AF107405	ANUC
Sequence 2025	AF061258	ANUC	Sequence 2077	AF120334	ANUC
Sequence 2026	AF062318	ANUC	Sequence 2078	AF124438	ANUC
Sequence 2027	AF063611	ANUC	Sequence 2079	AF124439	ANUC

TABLE 1A

Sequence 2080	AF125525	ANUC	Sequence 2132	AI091425	dbEST
Sequence 2081	AF131799	ANUC	Sequence 2133	AI092971	dbEST
Sequence 2082	AF131814	ANUC	Sequence 2134	AI095477	dbEST
Sequence 2083	AF139461	ANUC	Sequence 2135	AI123229	dbEST
Sequence 2084	AF139658	ANUC	Sequence 2136	AI125642	dbEST
Sequence 2085	AF144755	ANUC	Sequence 2137	AI125874	dbEST
Sequence 2086	AF147331	ANUC	Sequence 2138	AI127013	dbEST
Sequence 2087	AF150962	ANUC	Sequence 2139	AI127556	dbEST
Sequence 2088	AF151832	ANUC	Sequence 2140	AI140291	dbEST
Sequence 2089	AF151868	ANUC	Sequence 2141	AI141130	dbEST
Sequence 2090	AF151898	ANUC	Sequence 2142	AI141847	dbEST
Sequence 2091	AF151907	ANUC	Sequence 2143	AI143899	dbEST
Sequence 2092	AF152097	ANUC	Sequence 2144	AI144100	dbEST
Sequence 2093	AF159295	ANUC	Sequence 2145	AI148251	dbEST
Sequence 2094	AF176702	ANUC	Sequence 2146	AI149429	dbEST
Sequence 2095	AF190744	ANUC	Sequence 2147	AI149592	dbEST
Sequence 2096	AI004664	dbEST	Sequence 2148	AI186028	dbEST
Sequence 2097	AI004915	dbEST	Sequence 2149	AI186042	dbEST
Sequence 2098	AI016073	dbEST	Sequence 2150	AI190341	dbEST
Sequence 2099	AI016323	dbEST	Sequence 2151	AI192367	dbEST
Sequence 2100	AI016791	dbEST	Sequence 2152	AI192629	dbEST
Sequence 2101	AI018451	dbEST	Sequence 2153	AI198930	dbEST
Sequence 2102	AI018625	dbEST	Sequence 2154	AI216969	dbEST
Sequence 2103	AI022779	dbEST	Sequence 2155	AI217003	dbEST
Sequence 2104	AI023799	dbEST	Sequence 2156	AI223292	dbEST
Sequence 2105	AI026164	dbEST	Sequence 2157	AI241706	dbEST
Sequence 2106	AI027516	dbEST	Sequence 2158	AI251743	dbEST
Sequence 2107	AI031636	dbEST	Sequence 2159	AI252466	dbEST
Sequence 2108	AI033037	dbEST	Sequence 2160	AI253330	dbEST
Sequence 2109	AI034115	dbEST	Sequence 2161	AI253335	dbEST
Sequence 2110	AI037859	dbEST	Sequence 2162	AI253338	dbEST
Sequence 2111	AI041670	dbEST	Sequence 2163	AI253375	dbEST
Sequence 2112	AI042034	dbEST	Sequence 2164	AI253379	dbEST
Sequence 2113	AI042290	dbEST	Sequence 2165	AI253436	dbEST
Sequence 2114	AI051971	dbEST	Sequence 2166	AI262380	dbEST
Sequence 2115	AI056917	dbEST	Sequence 2167	AI263674	dbEST
Sequence 2116	AI057124	dbEST	Sequence 2168	AI267162	dbEST
Sequence 2117	AI066419	dbEST	Sequence 2169	AI267185	dbEST
Sequence 2118	AI078041	dbEST	Sequence 2170	AI267209	dbEST
Sequence 2119	AI081116	dbEST	Sequence 2171	AI267289	dbEST
Sequence 2120	AI081472	dbEST	Sequence 2172	AI267307	dbEST
Sequence 2121	AI081913	dbEST	Sequence 2173	AI267321	dbEST
Sequence 2122	AI082244	dbEST	Sequence 2174	AI267454	dbEST
Sequence 2123	AI082648	dbEST	Sequence 2175	AI267502	dbEST
Sequence 2124	AI084731	dbEST	Sequence 2176	AI268293	dbEST
Sequence 2125	AI085381	dbEST	Sequence 2177	AI269060	dbEST
Sequence 2126	AI087291	dbEST	Sequence 2178	AI269369	dbEST
Sequence 2127	AI087819	dbEST	Sequence 2179	AI270183	dbEST
Sequence 2128	AI088178	dbEST	Sequence 2180	AI270472	dbEST
Sequence 2129	AI089981	dbEST	Sequence 2181	AI271786	dbEST
Sequence 2130	AI090524	dbEST	Sequence 2182	AI272827	dbEST
Sequence 2131	AI090623	dbEST	Sequence 2183	AI274047	dbEST

TABLE 1A

Sequence 2184	AI276341	dbEST	Sequence 2236	AI608591	dbEST
Sequence 2185	AI276839	dbEST	Sequence 2237	AI608787	dbEST
Sequence 2186	AI278611	dbEST	Sequence 2238	AI608968	dbEST
Sequence 2187	AI280022	dbEST	Sequence 2239	AI609193	dbEST
Sequence 2188	AI283548	dbEST	Sequence 2240	AI609281	dbEST
Sequence 2189	AI288965	dbEST	Sequence 2241	AI623804	dbEST
Sequence 2190	AI290565	dbEST	Sequence 2242	AI628689	dbEST
Sequence 2191	AI291683	dbEST	Sequence 2243	AI636635	dbEST
Sequence 2192	AI292286	dbEST	Sequence 2244	AI650837	dbEST
Sequence 2193	AI298472	dbEST	Sequence 2245	AI654096	dbEST
Sequence 2194	AI298941	dbEST	Sequence 2246	AI660245	dbEST
Sequence 2195	AI304857	dbEST	Sequence 2247	AI669253	dbEST
Sequence 2196	AI308959	dbEST	Sequence 2248	AI670084	dbEST
Sequence 2197	AI312552	dbEST	Sequence 2249	AI674313	dbEST
Sequence 2198	AI333055	dbEST	Sequence 2250	AI678152	dbEST
Sequence 2199	AI333116	dbEST	Sequence 2251	AI678703	dbEST
Sequence 2200	AI335249	dbEST	Sequence 2252	AI679044	dbEST
Sequence 2201	AI336326	dbEST	Sequence 2253	AI679321	dbEST
Sequence 2202	AI345325	dbEST	Sequence 2254	AI683140	dbEST
Sequence 2203	AI366549	dbEST	Sequence 2255	AI683338	dbEST
Sequence 2204	AI367850	dbEST	Sequence 2256	AI683793	dbEST
Sequence 2205	AI375624	dbEST	Sequence 2257	AI688798	dbEST
Sequence 2206	AI376561	dbEST	Sequence 2258	AI692866	dbEST
Sequence 2207	AI399636	dbEST	Sequence 2259	AI694087	dbEST
Sequence 2208	AI417384	dbEST	Sequence 2260	AI696819	dbEST
Sequence 2209	AI421720	dbEST	Sequence 2261	AI697501	dbEST
Sequence 2210	AI424841	dbEST	Sequence 2262	AI734922	dbEST
Sequence 2211	AI431507	dbEST	Sequence 2263	AI735069	dbEST
Sequence 2212	AI433180	dbEST	Sequence 2264	AI739337	dbEST
Sequence 2213	AI434084	dbEST	Sequence 2265	AI739377	dbEST
Sequence 2214	AI434401	dbEST	Sequence 2266	AI743595	dbEST
Sequence 2215	AI436016	dbEST	Sequence 2267	AI743691	dbEST
Sequence 2216	AI436448	dbEST	Sequence 2268	AI750198	dbEST
Sequence 2217	AI446503	dbEST	Sequence 2269	AI750909	dbEST
Sequence 2218	AI453199	dbEST	Sequence 2270	AI751119	dbEST
Sequence 2219	AI459028	dbEST	Sequence 2271	AI751364	dbEST
Sequence 2220	AI469237	dbEST	Sequence 2272	AI751565	dbEST
Sequence 2221	AI492520	dbEST	Sequence 2273	AI752319	dbEST
Sequence 2222	AI492769	dbEST	Sequence 2274	AI752553	dbEST
Sequence 2223	AI494344	dbEST	Sequence 2275	AI752929	dbEST
Sequence 2224	AI523940	dbEST	Sequence 2276	AI753108	dbEST
Sequence 2225	AI524677	dbEST	Sequence 2277	AI753671	dbEST
Sequence 2226	AI538682	dbEST	Sequence 2278	AI754437	dbEST
Sequence 2227	AI557059	dbEST	Sequence 2279	AI755181	dbEST
Sequence 2228	AI561260	dbEST	Sequence 2280	AI758869	dbEST
Sequence 2229	AI567988	dbEST	Sequence 2281	AI761927	dbEST
Sequence 2230	AI569715	dbEST	Sequence 2282	AI763126	dbEST
Sequence 2231	AI581291	dbEST	Sequence 2283	AI791906	dbEST
Sequence 2232	AI583211	dbEST	Sequence 2284	AI793120	dbEST
Sequence 2233	AI583570	dbEST	Sequence 2285	AI799521	dbEST
Sequence 2234	AI589301	dbEST	Sequence 2286	AI804346	dbEST
Sequence 2235	AI597938	dbEST	Sequence 2287	AI808109	dbEST

TABLE 1A

Sequence 2288	AI811021	dbEST	Sequence 2340	AL049999	ANUC
Sequence 2289	AI811845	dbEST	Sequence 2341	AL050011	ANUC
Sequence 2290	AI814139	dbEST	Sequence 2342	AL050089	ANUC
Sequence 2291	AI814674	dbEST	Sequence 2343	AL050141	ANUC
Sequence 2292	AI815868	dbEST	Sequence 2344	AL050171	ANUC
Sequence 2293	AI822030	dbEST	Sequence 2345	AL050187	ANUC
Sequence 2294	AI827641	dbEST	Sequence 2346	AL050198	ANUC
Sequence 2295	AI859619	dbEST	Sequence 2347	AL050217	ANUC
Sequence 2296	AI864580	dbEST	Sequence 2348	AL050392	ANUC
Sequence 2297	AI878968	dbEST	Sequence 2349	AL080062	ANUC
Sequence 2298	AI879179	dbEST	Sequence 2350	AL080186	ANUC
Sequence 2299	AI879367	dbEST	Sequence 2351	AL080235	ANUC
Sequence 2300	AI879992	dbEST	Sequence 2352	AL096857	ANUC
Sequence 2301	AI888377	dbEST	Sequence 2353	AL096858	ANUC
Sequence 2302	AI911704	dbEST	Sequence 2354	AL110197	ANUC
Sequence 2303	AI911997	dbEST	Sequence 2355	AL110235	ANUC
Sequence 2304	AI912084	dbEST	Sequence 2356	AL117237	ANUC
Sequence 2305	AI916284	dbEST	Sequence 2357	AL117499	ANUC
Sequence 2306	AI916584	dbEST	Sequence 2358	AL117534	ANUC
Sequence 2307	AI923224	dbEST	Sequence 2359	AL118999	dbEST
Sequence 2308	AI924096	dbEST	Sequence 2360	AL119085	dbEST
Sequence 2309	AI928185	dbEST	Sequence 2361	AL119157	dbEST
Sequence 2310	AI929819	dbEST	Sequence 2362	AW020479	dbEST
Sequence 2311	AI936748	dbEST	Sequence 2363	AW044114	dbEST
Sequence 2312	AI950087	dbEST	Sequence 2364	AW102841	dbEST
Sequence 2313	AI955808	dbEST	Sequence 2365	C02094	dbEST
Sequence 2314	AJ001258	ANUC	Sequence 2366	C16886	dbEST
Sequence 2315	AJ002030	ANUC	Sequence 2367	C18886	dbEST
Sequence 2316	AJ006026	ANUC	Sequence 2368	D00017	ANUC
Sequence 2317	AJ011001	ANUC	Sequence 2369	D00022	ANUC
Sequence 2318	AJ011915	ANUC	Sequence 2370	D00068	ANUC
Sequence 2319	AJ012499	ANUC	Sequence 2371	D00099	ANUC
Sequence 2320	AJ223183	ANUC	Sequence 2372	D00422	ANUC
Sequence 2321	AL035802	dbEST	Sequence 2373	D10495	ANUC
Sequence 2322	AL035987	dbEST	Sequence 2374	D13119	ANUC
Sequence 2323	AL036801	dbEST	Sequence 2375	D13287	ANUC
Sequence 2324	AL037646	dbEST	Sequence 2376	D13665	ANUC
Sequence 2325	AL038985	dbEST	Sequence 2377	D13866	ANUC
Sequence 2326	AL039150	dbEST	Sequence 2378	D14662	ANUC
Sequence 2327	AL041780	dbEST	Sequence 2379	D14697	ANUC
Sequence 2328	AL044019	dbEST	Sequence 2380	D14710	ANUC
Sequence 2329	AL046804	dbEST	Sequence 2381	D14812	ANUC
Sequence 2330	AL049055	dbEST	Sequence 2382	D15049	ANUC
Sequence 2331	AL049227	ANUC	Sequence 2383	D16431	ANUC
Sequence 2332	AL049229	ANUC	Sequence 2384	D16937	ANUC
Sequence 2333	AL049296	ANUC	Sequence 2385	D17188	ANUC
Sequence 2334	AL049464	ANUC	Sequence 2386	D17268	ANUC
Sequence 2335	AL049953	ANUC	Sequence 2387	D17409	ANUC
Sequence 2336	AL049954	ANUC	Sequence 2388	D17793	ANUC
Sequence 2337	AL049955	ANUC	Sequence 2389	D21063	ANUC
Sequence 2338	AL049959	ANUC	Sequence 2390	D23660	ANUC
Sequence 2339	AL049987	ANUC	Sequence 2391	D25542	ANUC

TABLE 1A

Sequence 2392	D28759	ANUC	Sequence 2444	E01813	ANUC
Sequence 2393	D29677	ANUC	Sequence 2445	E01827	ANUC
Sequence 2394	D31767	ANUC	Sequence 2446	E01979	ANUC
Sequence 2395	D31784	ANUC	Sequence 2447	E02628	ANUC
Sequence 2396	D31883	ANUC	Sequence 2448	E02651	ANUC
Sequence 2397	D31890	ANUC	Sequence 2449	E03569	ANUC
Sequence 2398	D37991	ANUC	Sequence 2450	E06721	ANUC
Sequence 2399	D38491	ANUC	Sequence 2451	E07218	ANUC
Sequence 2400	D38583	ANUC	Sequence 2452	F28779	dbEST
Sequence 2401	D43948	ANUC	Sequence 2453	F30276	dbEST
Sequence 2402	D43950	ANUC	Sequence 2454	F31082	dbEST
Sequence 2403	D45248	ANUC	Sequence 2455	H03854	dbEST
Sequence 2404	D45887	ANUC	Sequence 2456	H05412	dbEST
Sequence 2405	D45915	ANUC	Sequence 2457	H08994	dbEST
Sequence 2406	D49489	ANUC	Sequence 2458	H13339	dbEST
Sequence 2407	D49547	ANUC	Sequence 2459	H16426	dbEST
Sequence 2408	D50310	ANUC	Sequence 2460	H39960	dbEST
Sequence 2409	D50371	ANUC	Sequence 2461	H48742	dbEST
Sequence 2410	D55192	dbEST	Sequence 2462	H59372	dbEST
Sequence 2411	D55649	ANUC	Sequence 2463	H60722	dbEST
Sequence 2412	D56120	dbEST	Sequence 2464	H69238	dbEST
Sequence 2413	D59253	ANUC	Sequence 2465	H72481	dbEST
Sequence 2414	D78586	ANUC	Sequence 2466	H75695	dbEST
Sequence 2415	D79826	dbEST	Sequence 2467	H78517	dbEST
Sequence 2416	D79983	ANUC	Sequence 2468	H79084	dbEST
Sequence 2417	D79986	ANUC	Sequence 2469	H84729	dbEST
Sequence 2418	D79997	ANUC	Sequence 2470	H85709	dbEST
Sequence 2419	D80006	ANUC	Sequence 2471	H89654	dbEST
Sequence 2420	D80012	ANUC	Sequence 2472	J00269	ANUC
Sequence 2421	D80087	dbEST	Sequence 2473	J02621	ANUC
Sequence 2422	D80253	dbEST	Sequence 2474	J03005	ANUC
Sequence 2423	D81635	dbEST	Sequence 2475	J03040	ANUC
Sequence 2424	D82128	dbEST	Sequence 2476	J03171	ANUC
Sequence 2425	D82348	ANUC	Sequence 2477	J03191	ANUC
Sequence 2426	D83197	ANUC	Sequence 2478	J03210	ANUC
Sequence 2427	D83327	ANUC	Sequence 2479	J03464	ANUC
Sequence 2428	D83784	ANUC	Sequence 2480	J03473	ANUC
Sequence 2429	D86227	ANUC	Sequence 2481	J03799	ANUC
Sequence 2430	D87437	ANUC	Sequence 2482	J04080	ANUC
Sequence 2431	D87442	ANUC	Sequence 2483	J04164	ANUC
Sequence 2432	D87470	ANUC	Sequence 2484	J04177	ANUC
Sequence 2433	D87666	ANUC	Sequence 2485	J04765	ANUC
Sequence 2434	D87667	ANUC	Sequence 2486	J05013	ANUC
Sequence 2435	D87682	ANUC	Sequence 2487	J05021	ANUC
Sequence 2436	D87735	ANUC	Sequence 2488	J05192	ANUC
Sequence 2437	D87969	ANUC	Sequence 2489	J05633	ANUC
Sequence 2438	D89052	ANUC	Sequence 2490	K00558	ANUC
Sequence 2439	D90226	ANUC	Sequence 2491	K01566	ANUC
Sequence 2440	D90373	ANUC	Sequence 2492	K02765	ANUC
Sequence 2441	E00882	ANUC	Sequence 2493	L00160	ANUC
Sequence 2442	E01650	ANUC	Sequence 2494	L02547	ANUC
Sequence 2443	E01797	ANUC	Sequence 2495	L05092	ANUC

TABLE 1A

Sequence 2496	L05186	ANUC	Sequence 2548	M26041	ANUC
Sequence 2497	L07633	ANUC	Sequence 2549	M26152	ANUC
Sequence 2498	L11066	ANUC	Sequence 2550	M26325	ANUC
Sequence 2499	L11932	ANUC	Sequence 2551	M27913	ANUC
Sequence 2500	L12711	ANUC	Sequence 2552	M27971	ANUC
Sequence 2501	L13848	ANUC	Sequence 2553	M28373	ANUC
Sequence 2502	L14599	ANUC	Sequence 2554	M31159	ANUC
Sequence 2503	L19161	ANUC	Sequence 2555	M31212	ANUC
Sequence 2504	L19184	ANUC	Sequence 2556	M31899	ANUC
Sequence 2505	L19597	ANUC	Sequence 2557	M32110	ANUC
Sequence 2506	L20941	ANUC	Sequence 2558	M32790	ANUC
Sequence 2507	L23959	ANUC	Sequence 2559	M32798	ANUC
Sequence 2508	L26081	ANUC	Sequence 2560	M33308	ANUC
Sequence 2509	L27560	ANUC	Sequence 2561	M34064	ANUC
Sequence 2510	L28010	ANUC	Sequence 2562	M37583	ANUC
Sequence 2511	L28809	ANUC	Sequence 2563	M38106	ANUC
Sequence 2512	L33404	ANUC	Sequence 2564	M55409	ANUC
Sequence 2513	L33930	ANUC	Sequence 2565	M55542	ANUC
Sequence 2514	L34155	ANUC	Sequence 2566	M58485	ANUC
Sequence 2515	L34839	ANUC	Sequence 2567	M60457	ANUC
Sequence 2516	L38486	ANUC	Sequence 2568	M60854	ANUC
Sequence 2517	L42024	ANUC	Sequence 2569	M62403	ANUC
Sequence 2518	L43575	ANUC	Sequence 2570	M62810	ANUC
Sequence 2519	L44349	dbEST	Sequence 2571	M64241	ANUC
Sequence 2520	L54057	ANUC	Sequence 2572	M67468	ANUC
Sequence 2521	M10036	ANUC	Sequence 2573	M69181	ANUC
Sequence 2522	M10119	ANUC	Sequence 2574	M74002	ANUC
Sequence 2523	M10905	ANUC	Sequence 2575	M75126	ANUC
Sequence 2524	M11146	ANUC	Sequence 2576	M76729	ANUC
Sequence 2525	M13573	ANUC	Sequence 2577	M78113	dbEST
Sequence 2526	M13955	ANUC	Sequence 2578	M81757	ANUC
Sequence 2527	M14083	ANUC	Sequence 2579	M83248	ANUC
Sequence 2528	M14483	ANUC	Sequence 2580	M84739	ANUC
Sequence 2529	M14630	ANUC	Sequence 2581	M87503	ANUC
Sequence 2530	M14631	ANUC	Sequence 2582	M88279	ANUC
Sequence 2531	M15182	ANUC	Sequence 2583	M92357	ANUC
Sequence 2532	M15800	ANUC	Sequence 2584	N20576	dbEST
Sequence 2533	M16247	ANUC	Sequence 2585	N34255	dbEST
Sequence 2534	M16553	ANUC	Sequence 2586	N35187	dbEST
Sequence 2535	M16660	ANUC	Sequence 2587	N35421	dbEST
Sequence 2536	M16937	ANUC	Sequence 2588	N39717	dbEST
Sequence 2537	M17597	ANUC	Sequence 2589	N40823	dbEST
Sequence 2538	M17885	ANUC	Sequence 2590	N40852	dbEST
Sequence 2539	M20372	ANUC	Sequence 2591	N67927	dbEST
Sequence 2540	M22146	ANUC	Sequence 2592	N76180	dbEST
Sequence 2541	M22382	ANUC	Sequence 2593	N76677	dbEST
Sequence 2542	M22590	ANUC	Sequence 2594	N77080	dbEST
Sequence 2543	M22918	ANUC	Sequence 2595	N84497	dbEST
Sequence 2544	M22920	ANUC	Sequence 2596	N86776	dbEST
Sequence 2545	M23613	ANUC	Sequence 2597	N91638	dbEST
Sequence 2546	M24194	ANUC	Sequence 2598	N92086	dbEST
Sequence 2547	M25246	ANUC	Sequence 2599	N99205	dbEST

TABLE 1A

Sequence 2600	Q37741	NUCPATENT	Sequence 2652	U22815	ANUC
Sequence 2601	Q48043	NUCPATENT	Sequence 2653	U24105	ANUC
Sequence 2602	Q65676	NUCPATENT	Sequence 2654	U24153	ANUC
Sequence 2603	Q90526	NUCPATENT	Sequence 2655	U27768	ANUC
Sequence 2604	R06046	dbEST	Sequence 2656	U33760	ANUC
Sequence 2605	R17092	dbEST	Sequence 2657	U33833	ANUC
Sequence 2606	R47228	dbEST	Sequence 2658	U34877	ANUC
Sequence 2607	R55150	dbEST	Sequence 2659	U39361	ANUC
Sequence 2608	R55398	dbEST	Sequence 2660	U41515	ANUC
Sequence 2609	R68132	dbEST	Sequence 2661	U46570	ANUC
Sequence 2610	R72676	dbEST	Sequence 2662	U50733	ANUC
Sequence 2611	R73306	dbEST	Sequence 2663	U51586	ANUC
Sequence 2612	R78333	dbEST	Sequence 2664	U56255	ANUC
Sequence 2613	R92367	dbEST	Sequence 2665	U59305	ANUC
Sequence 2614	R93637	dbEST	Sequence 2666	U60975	ANUC
Sequence 2615	R99649	dbEST	Sequence 2667	U61083	ANUC
Sequence 2616	S41458	ANUC	Sequence 2668	U61397	ANUC
Sequence 2617	S42303	ANUC	Sequence 2669	U63846	ANUC
Sequence 2618	S54005	ANUC	Sequence 2670	U67784	ANUC
Sequence 2619	S66431	ANUC	Sequence 2671	U68723	ANUC
Sequence 2620	S70154	ANUC	Sequence 2672	U68727	ANUC
Sequence 2621	S70290	ANUC	Sequence 2673	U68758	ANUC
Sequence 2622	S79895	ANUC	Sequence 2674	U70735	ANUC
Sequence 2623	S82076	ANUC	Sequence 2675	U77085	ANUC
Sequence 2624	T02792	NUCPATENT	Sequence 2676	U79258	ANUC
Sequence 2625	T24119	dbEST	Sequence 2677	U79274	ANUC
Sequence 2626	T49314	dbEST	Sequence 2678	U79278	ANUC
Sequence 2627	T53479	dbEST	Sequence 2679	U80213	ANUC
Sequence 2628	T58797	dbEST	Sequence 2680	U81234	ANUC
Sequence 2629	T64560	dbEST	Sequence 2681	U82130	ANUC
Sequence 2630	T66112	NUCPATENT	Sequence 2682	U86602	ANUC
Sequence 2631	T92160	NUCPATENT	Sequence 2683	U87309	ANUC
Sequence 2632	T92396	dbEST	Sequence 2684	U90028	ANUC
Sequence 2633	U00947	ANUC	Sequence 2685	U90441	ANUC
Sequence 2634	U04815	ANUC	Sequence 2686	U90902	ANUC
Sequence 2635	U07151	ANUC	Sequence 2687	U90917	ANUC
Sequence 2636	U07857	ANUC	Sequence 2688	U94831	ANUC
Sequence 2637	U08470	ANUC	Sequence 2689	V00478	ANUC
Sequence 2638	U10323	ANUC	Sequence 2690	V00503	ANUC
Sequence 2639	U10439	ANUC	Sequence 2691	V05728	NUCPATENT
Sequence 2640	U12465	ANUC	Sequence 2692	V11636	NUCPATENT
Sequence 2641	U13665	ANUC	Sequence 2693	V57903	NUCPATENT
Sequence 2642	U13877	ANUC	Sequence 2694	V59662	NUCPATENT
Sequence 2643	U14550	ANUC	Sequence 2695	V59746	NUCPATENT
Sequence 2644	U14966	ANUC	Sequence 2696	V84428	NUCPATENT
Sequence 2645	U15008	ANUC	Sequence 2697	V86232	NUCPATENT
Sequence 2646	U16306	ANUC	Sequence 2698	V87930	NUCPATENT
Sequence 2647	U17104	ANUC	Sequence 2699	W07215	dbEST
Sequence 2648	U17496	ANUC	Sequence 2700	W19127	dbEST
Sequence 2649	U19769	ANUC	Sequence 2701	W19407	dbEST
Sequence 2650	U20896	ANUC	Sequence 2702	W19441	dbEST
Sequence 2651	U22431	ANUC	Sequence 2703	W25547	dbEST



TABLE 1A

Sequence 2704	W26197	dbEST	Sequence 2756	X73902	ANUC
Sequence 2705	W38952	dbEST	Sequence 2757	X74039	ANUC
Sequence 2706	W56388	dbEST	Sequence 2758	X74801	ANUC
Sequence 2707	W68015	dbEST	Sequence 2759	X74979	ANUC
Sequence 2708	W73140	dbEST	Sequence 2760	X76013	ANUC
Sequence 2709	W73168	dbEST	Sequence 2761	X76180	ANUC
Sequence 2710	W76204	dbEST	Sequence 2762	X78627	ANUC
Sequence 2711	W87522	dbEST	Sequence 2763	X81109	ANUC
Sequence 2712	W87891	dbEST	Sequence 2764	X82676	ANUC
Sequence 2713	X00351	ANUC	Sequence 2765	X84939	NUCPATENT
Sequence 2714	X00497	ANUC	Sequence 2766	X85373	ANUC
Sequence 2715	X01742	ANUC	Sequence 2767	X93036	ANUC
Sequence 2716	X01924	NUCPATENT	Sequence 2768	X93207	ANUC
Sequence 2717	X03084	ANUC	Sequence 2769	X94323	ANUC
Sequence 2718	X04098	ANUC	Sequence 2770	X94754	ANUC
Sequence 2719	X04408	ANUC	Sequence 2771	X97324	ANUC
Sequence 2720	X04470	ANUC	Sequence 2772	X99920	ANUC
Sequence 2721	X05276	ANUC	Sequence 2773	Y00503	ANUC
Sequence 2722	X05908	ANUC	Sequence 2774	Y00757	ANUC
Sequence 2723	X06700	ANUC	Sequence 2775	Y00815	ANUC
Sequence 2724	X07819	ANUC	Sequence 2776	Y09188	ANUC
Sequence 2725	X13425	ANUC	Sequence 2777	Y11435	ANUC
Sequence 2726	X14420	ANUC	Sequence 2778	Y12065	ANUC
Sequence 2727	X15729	ANUC	Sequence 2779	Y13247	ANUC
Sequence 2728	X15880	ANUC	Sequence 2780	Y13286	ANUC
Sequence 2729	X16869	ANUC	Sequence 2781	Y15286	ANUC
Sequence 2730	X17206	ANUC	Sequence 2782	Y17114	ANUC
Sequence 2731	X24068	NUCPATENT	Sequence 2783	Z18538	ANUC
Sequence 2732	X37385	NUCPATENT	Sequence 2784	Z18954	ANUC
Sequence 2733	X37509	NUCPATENT	Sequence 2785	Z19054	ANUC
Sequence 2734	X40178	NUCPATENT	Sequence 2786	Z21507	ANUC
Sequence 2735	X51466	ANUC	Sequence 2787	Z26317	ANUC
Sequence 2736	X53505	ANUC	Sequence 2788	Z29093	ANUC
Sequence 2737	X54304	ANUC	Sequence 2789	Z31696	ANUC
Sequence 2738	X54941	ANUC	Sequence 2790	Z32564	ANUC
Sequence 2739	X55110	ANUC	Sequence 2791	Z36531	ANUC
Sequence 2740	X55885	ANUC	Sequence 2792	Z37986	ANUC
Sequence 2741	X56932	ANUC	Sequence 2793	Z46629	ANUC
Sequence 2742	X56998	ANUC	Sequence 2794	Z47087	ANUC
Sequence 2743	X56999	ANUC	Sequence 2795	Z74615	ANUC
Sequence 2744	X57766	ANUC			
Sequence 2745	X62744	ANUC			
Sequence 2746	X63432	ANUC			
Sequence 2747	X66360	ANUC			
Sequence 2748	X67698	ANUC			
Sequence 2749	X68277	ANUC			
Sequence 2750	X68880	ANUC			
Sequence 2751	X69398	ANUC			
Sequence 2752	X69838	ANUC			
Sequence 2753	X70340	ANUC			
Sequence 2754	X71087	ANUC			
Sequence 2755	X73608	ANUC			



TABLE 1A

Sequence 340: found in patent publication W098/39446

AGGCGTNCCTCTGACTGCCCACTCAGTGGCNNCACNNGGAGCTGNTTTGGNGCTTTGGG  
GANCTNAACANTTNCNTCTTTCAAACCTNACTGGC

Sequence 1962: found in patent publication W098/42738

AGGTACCCGCTCTCCTGCTTCAGTAAATCTCCACTCGATCTCAGTGGGTTTCCTGTCCAT  
AGGATCCACAAGTTTGACCTGGCGGTGGAGCAAGGGGGCTTCACTAGGGATCATGGTTCC  
CCGGTAATCCATGGTCTTGCCAATGTAGCCGGTAATGTGTGTTTCAGCCCTCCACGACCA  
CCCAGTTTCGCTGCCGGATAACTTGAACCACTTTGCCCTGCTTCCCGGCATCCTTGCCTT  
CTAGGATCTCCACCGTGTCCCCACAGAACAGATACCAAGTCTTCATCAGAGATGGGTTCC  
ACAACCACTGGGCCGCCGCTGATCCATGGGGGGTCTTCTCTTGTCTGCAACAGAGCC  
TGGGGGGCTCATCCATAACGGTNATGGGGGGGCAGAGTGACCTTGGATGCCAAGGCCAGC  
AGGGGCAAGAAAGACCCATGCCTGGAGGTTGNAAGAAAATCCCTTTGCCAGCAAAAACGC  
TTCGAAACCCTTNCCTTGTCAAGCTTTTCACTTTTCCGNGGCACCTTTGGGATTTTA  
GCACATTGGGGCCCTTAAGNGTTCCTTCCCC

Sequence 341: found in patent publication W099/039941

CCCTTAGCGNGGTCGCGGCCGAGGCACAATTCGATTATTCACANGAAAGGGCAAACCTGTT  
NNTGTTNGCTGGCAGGAGNAGGTGCATATATACCAGCACTTCAAGTNNGGTATTTCCATT  
CAGGACATTTTATCTCTGTGCAAAGACCGGAGTAGAAGCTGATGAGTGGATCAAGATATT  
ACGCTGGAAATTGTCACAAATAAGAAAACAGCTCAACCAAGGGGAAGGCACCGATCCGAT  
CTCGGTGCTTCATCTTTAAATAGATCTTTCTTGCCAAGGAATGCTCTGGCCCAGGAGCAA  
GGTGGAATGCTTCCCTGACGCTGCGATCTGCAGCAGACTNCAAATGAAAACCGACTAAGG  
ATTTCTTTCAAAAACAAATCAGAAGCAGATGCTGATTGGGACCCATATACCACGTTGCT  
GACTCACCGTTGCTGCCCTTNCATGGATGTTGCCATCTGCTTGAGAACACTGAAGCAATC  
ACCATTCTNGATANGAAAGTGCTTAACCCCCACTCTTAGGGCTGCTCACTTCTTAGAAC  
ACACAAAGGGGAAGAGGAAAGGGGT

Sequence 342: found in patent publication W099/18126

CCGCGGTGGCGGCCGCCCGGGCAGGTACCTACAGTGACACAGATCCCCCTCCCGCCATCCT  
GGTCACACTGAATAACAAAGGGAAGAGAGGAGTAAGAACTGTAGTATCTAGAAATTCTCA  
GCACAGTGAAGGAAAGTGATCTTCTACTTTGTATTTCAGGCCTAAAAAAGGAGGGGACGGG  
CCCGGCACAGTGGCTCACACCTGTAATCCCAGCACTCTGGGAGGCCAAGGAGGGCAGATC  
ACCTGAGGTTGGGAGTTTGAGACCAGCCTGACCAACATGGAGAAACCCTGTCTCTACTAA  
AAACACAAAATTAGCCAGGCATGGTGGCATGCGCCTATAAACCAAGCTACTCAAGAGGCT  
GAGGCAGGAGAATTGCTTGAACCCAGGAGGCAGAGGTTGTGGTGAGCCAAGATCGAGCCA  
TCACACTCCAGCCTGGGCAACAAGAGCAAACTCTGTCTCAAAAAAAAAAACAGGAGAGG  
AGGGAG

Sequence 1016: found in patent publication W099/38881

CTACTTAGGGCGAATTGGAGCTCCCCGCGGTGGCGGCCGAGGTCAAGCTTCGACCCCGCG  
TCCGTGATAAACTACTTTTGGGTTTTATTTTCATTGAGGCACTTTTTTTATTGTTTGAATG  
ATTCCGGCTTGTAATATATCAGCCTCTACAATGAAATGCAGAAGAGTTCATTTTTCTAG  
ATCTGTTTTTCATTAGAAATATTGACAAATAACACATTGTCAACCTGGATCCTTTGACA  
TTTACTTAACTCTGGCATGTTACAAAAAGTAGAACTCTAAGAGACCATTACCATTTT  
TCACAGATGTATAGGGGATGTATTCTAAAACTGACAGAAAAGAGAATNTGATAGTCAAC  
ACTGTAACTTTTACTGNGTAATTGCCAAATACACTTTTCCAAATTTGTCCCAACAGCC  
TNTAAGCCAGCTTTCTTCTATATTTATAA

Sequence 1963: found in patent publication W099/46289

AACTGGACAGAGTAAGGGAATTCCAGCATCCTCTTCCTGCTTGCTCGTGTTACCCACAG  
ATCAAACCTCAATTCTAGTTGGGGATGCTGTCTAGCCCCACACCATGACTGAAGCCTTA

TABLE 1A

AGCACTGTTGCGCCTCATGTGCTTTGGATCAGCAACCCCAGTGGTATTCTACCAGAGCAT  
TGTGGGAAAGCAGATGTATAGTCAGGTCCCAACAGCAAATTGTTGGGTGTGAGAGTTCTA  
AAGTATAGGGGTGAGGGAAGAGAAGGATATGAACTCCTCTGACCTTAAGCCAGCATTTCAT  
TTAACTTTTATGTCTACTTAACAAGAGAACCTGNAGAAAACTACCGTATTCAAGAGATA  
ATCAAAATCAGTGTTTTAGCCAGGCGATGACAGAGAAGCACCATTCTCACCCTCCATTC  
TTGTAATGTCTGTAATAAATTTTCAGTGCGTCAGGATGGATGAACCCAAGATCCAGTGAAT  
GATTCAGCTGTTCCAAGCCTTACATTTTCCATCATTATCATCCATTCTCATTTCAGTGA  
ACCTCTTGCACTATTGTGGTTAATTTTATGTAAAACCAGTTTATGTTTTTTTTTTAATAT  
GTGCCTATGTAATAAAAGTCTACACACTGGCAAAAAAAAAAAAAAAAAAAAAAGTCCTN

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AA001066	DBEst	1437265
AA007157	DBEst	1463141
AA010954	DBEst	1472001
AA015792	DBEst	1476959
AA019769	DBEst	1483105
AA019948	DBEst	1483755
AA022925	DBEst	1487005
AA022937	DBEst	1487036
AA024405	DBEst	1489413
AA029750	DBEst	1496068
AA031509	DBEst	1501463
AA033876	DBEst	1505694
AA034237	DBEst	1506265
AA039967	DBEst	1516280
AA040073	DBEst	1516350
AA040122	DBEst	1516400
AA045732	DBEst	1525626
AA045861	DBEst	1525757
AA046835	DBEst	1524734
AA047026	DBEst	1525061
AA047417	DBEst	1525463
AA053486	DBEst	1544124
AA054658	DBEst	1545600
AA055606	DBEst	1547963
AA056113	DBEst	1548469
AA056176	DBEst	1548514
AA056363	DBEst	1548703
AA056431	DBEst	1548771
AA065336	DBEst	1929216
AA069781	DBEst	1577149
AA069784	DBEst	1577152
AA069839	DBEst	1577199
AA069983	DBEst	1577343
AA071255	DBEst	1578610
AA075135	DBEst	1615139
AA081655	DBEst	1623857
AA082245	DBEst	1624304
AA083471	DBEst	1625557
AA083510	DBEst	1625570
AA085862	DBEst	1629449
AA085872	DBEst	1629244
AA085947	DBEst	1629482
AA088770	DBEst	1634335
AA100333	DBEst	1646685
AA100719	DBEst	1647074
AA100793	DBEst	1647210
AA100852	DBEst	1647269
AA101270	DBEst	1647951
AA101561	DBEst	1648449
AA111907	DBEst	1663978
AA112043	DBEst	1664189
AA112308	DBEst	1664577
AA112375	DBEst	1664785
AA113860	DBEst	1667753
AA114120	DBEst	1667996
AA115118	DBEst	1669966
AA115368	DBEst	1670548
AA122286	DBEst	1678525
AA122348	DBEst	1678587
AA126109	DBEst	1685775
AA127105	DBEst	1686466
AA127132	DBEst	1686492
AA127418	DBEst	1686707

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AA128305	DBEst	1688355
AA129461	DBEst	1689332
AA130252	DBEst	1691414
AA130547	DBEst	1692100
AA130786	DBEst	1692452
AA131041	DBEst	1692593
AA131065	DBEst	1692555
AA131104	DBEst	1692612
AA131155	DBEst	1692646
AA131160	DBEst	1692668
AA132182	DBEst	1693860
AA132568	DBEst	1694075
AA132598	DBEst	1694087
AA133351	DBEst	1690319
AA133927	DBEst	1690994
AA134105	DBEst	1691317
AA134210	DBEst	1691566
AA135032	DBEst	1696143
AA135919	DBEst	1697105
AA136383	DBEst	1697611
AA136789	DBEst	1697998
AA143609	DBEst	1713177
AA146773	DBEst	1716163
AA147806	DBEst	1717195
AA148160	DBEst	1717542
AA148268	DBEst	1717666
AA148771	DBEst	1721626
AA149056	DBEst	1719347
AA150307	DBEst	1721837
AA151310	DBEst	1719502
AA151775	DBEst	1720675
AA152037	DBEst	1720875
AA152416	DBEst	1718626
AA155853	DBEst	1727471
AA155926	DBEst	1727544
AA157405	DBEst	1729013
AA157725	DBEst	1729350
AA157788	DBEst	1732599
AA158165	DBEst	1732959
AA158171	DBEst	1732965
AA159272	DBEst	1734074
AA160114	DBEst	1734680
AA160685	DBEst	1736087
AA161410	DBEst	1735771
AA164405	DBEst	1740715
AA164465	DBEst	1740624
AA165083	DBEst	1740311
AA165629	DBEst	1741662
AA166973	DBEst	1745366
AA171510	DBEst	1750569
AA173031	DBEst	1754310
AA173470	DBEst	1753798
AA173630	DBEst	1753763
AA179462	DBEst	1760830
AA187003	DBEst	1775129
AA187958	DBEst	1774167
AA188591	DBEst	1775616
AA192108	DBEst	1781932
AA199710	DBEst	1795418
AA203224	DBEst	1798950
AA203284	DBEst	1799010
AA205851	DBEst	1801222
AA209431	DBEst	1807445

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AA209531	DBEst	1807492
AA214075	DBEst	1812697
AA216612	DBEst	1817292
AA224230	DBEst	1844755
AA224985	DBEst	1846276
AA226502	DBEst	1847857
AA229225	DBEst	1851057
AA232626	DBEst	1855763
AA233843	DBEst	1856920
AA242891	DBEst	1873684
AA250725	DBEst	1885707
AA250982	DBEst	1885943
AA256959	DBEst	1891227
AA259077	DBEst	1894348
AA262440	DBEst	1897800
AA263110	DBEst	1898920
AA283165	DBEst	1926099
AA285260	DBEst	1929570
AA287112	DBEst	1934119
AA292191	DBEst	1940291
AA292334	DBEst	1940314
AA292385	DBEst	1940380
AA292771	DBEst	1941593
AA293273	DBEst	1941423
AA293572	DBEst	1941239
AA295348	DBEst	1947743
AA295485	DBEst	1947839
AA301631	DBEst	1954115
AA304669	DBEst	1957001
AA304961	DBEst	1957288
AA305193	DBEst	1957520
AA305438	DBEst	1957763
AA306542	DBEst	1958871
AA306708	DBEst	1959036
AA306945	DBEst	1959275
AA307239	DBEst	1959567
AA307477	DBEst	1960025
AA307504	DBEst	1959872
AA307697	DBEst	1960187
AA307779	DBEst	1960177
AA308062	DBEst	1960391
AA308801	DBEst	1961131
AA309028	DBEst	1961354
AA309988	DBEst	1962337
AA311006	DBEst	1963405
AA311481	DBEst	1963975
AA312012	DBEst	1964341
AA313684	DBEst	1966083
AA314146	DBEst	1966495
AA315049	DBEst	1967529
AA315308	DBEst	1967637
AA315426	DBEst	1967755
AA316682	DBEst	1969010
AA319958	DBEst	1972449
AA320346	DBEst	1972675
AA320991	DBEst	1973319
AA328544	DBEst	1980860
AA330457	DBEst	1982700
AA338793	DBEst	1991103
AA340069	DBEst	1992307
AA341170	DBEst	1993406
AA342394	DBEst	1994715
AA348250	DBEst	2000486

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AA349148	DBEst	2001385
AA351443	DBEst	2003763
AA351880	DBEst	2004198
AA356158	DBEst	2008784
AA356187	DBEst	2008516
AA356195	DBEst	2008524
AA357374	DBEst	2009714
AA367446	DBEst	2019764
AA375236	DBEst	2027555
AA377718	DBEst	2030037
AA380997	DBEst	2033336
AA383917	DBEst	2036256
AA385147	DBEst	2037466
AA389641	DBEst	2042627
AA393164	DBEst	2046134
AA393236	DBEst	2046205
AA394242	DBEst	2047227
AA398732	DBEst	2051854
AA401864	DBEst	2055883
AA410508	DBEst	2069614
AA410580	DBEst	2069686
AA410942	DBEst	2070196
AA411334	DBEst	2068883
AA411599	DBEst	2069132
AA418061	DBEst	2079935
AA418473	DBEst	2080273
AA418970	DBEst	2080798
AA420789	DBEst	2094677
AA421682	DBEst	2100499
AA421850	DBEst	2100809
AA424529	DBEst	2103499
AA428421	DBEst	2112235
AA429754	DBEst	2112972
AA441787	DBEst	2153671
AA451633	DBEst	2165302
AA453309	DBEst	2166978
AA453559	DBEst	2167228
AA453570	DBEst	2167239
AA454871	DBEst	2177647
AA454913	DBEst	2177689
AA456892	DBEst	2179612
AA457048	DBEst	2179768
AA463426	DBEst	2188310
AA465039	DBEst	2189923
AA477173	DBEst	2205857
AA480921	DBEst	2210473
AA484050	DBEst	2212863
AA484756	DBEst	2214141
AA487483	DBEst	2217647
AA489640	DBEst	2219242
AA493886	DBEst	2223727
AA494493	DBEst	2224280
AA496518	DBEst	2229839
AA501749	DBEst	2236716
AA501822	DBEst	2236789
AA501945	DBEst	2236912
AA504490	DBEst	2240650
AA507234	DBEst	2243673
AA513640	DBEst	2252052
AA526227	DBEst	2268296
AA526889	DBEst	2268958
AA527139	DBEst	2269208
AA527188	DBEst	2269257

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AA531428	DBEst	2274134
AA532633	DBEst	2276887
AA535471	DBEst	2279724
AA554757	DBEst	2325296
AA565996	DBEst	2337635
AA568217	DBEst	2341271
AA573742	DBEst	2348257
AA573893	DBEst	2348408
AA574237	DBEst	2348752
AA576866	DBEst	2354340
AA579034	DBEst	2357218
AA579816	DBEst	2358000
AA581220	DBEst	2358992
AA581264	DBEst	2359036
AA582093	DBEst	2360771
AA583091	DBEst	2360451
AA584411	DBEst	2369020
AA586776	DBEst	2397590
AA587110	DBEst	2397924
AA587233	DBEst	2398047
AA587700	DBEst	2401875
AA609259	DBEst	2457687
AA609837	DBEst	2458265
AA613907	DBEst	2466041
AA614529	DBEst	2466725
AA618033	DBEst	2505238
AA628487	DBEst	2540874
AA631204	DBEst	2553815
AA631811	DBEst	2554422
AA640901	DBEst	2566151
AA641841	DBEst	2567059
AA642215	DBEst	2567433
AA643602	DBEst	2568820
AA651720	DBEst	2583372
AA664996	DBEst	2619609
AA668297	DBEst	2629796
AA668836	DBEst	2630335
AA675923	DBEst	2775270
AA687833	DBEst	2674739
AA704992	DBEst	2714910
AA732702	DBEst	2753309
AA745241	DBEst	2785227
AA746481	DBEst	2786467
AA758889	DBEst	2806752
AA772570	DBEst	2824353
AA772790	DBEst	2825632
AA776709	DBEst	2836043
AA776811	DBEst	2836142
AA777384	DBEst	2836715
AA778116	DBEst	2837517
AA779868	DBEst	2839199
AA781343	DBEst	2840674
AA809984	DBEst	2879390
AA810945	DBEst	2880556
AA811200	DBEst	2880811
AA825768	DBEst	2899080
AA828073	DBEst	2900436
AA828722	DBEst	2901821
AA843176	DBEst	2929694
AA843661	DBEst	2930179
AA876526	DBEst	2985603
AA883255	DBEst	2992785
AA906652	DBEst	3042238

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AA917638	DBEst	3057528
AA927734	DBEst	3076554
AA954939	DBEst	3118634
AA962622	DBEst	3134786
AA991285	DBEst	3177774
AB000115	GenBank	2564034
AB004047	GenBank	2116654
AB006746	GenBank	3510296
AB007619	GenBank	2465179
AB007860	GenBank	2662080
AB007965	GenBank	3413940
AB011101	GenBank	3043581
AB011169	GenBank	3043717
AB012701	GenBank	6714554
AB014536	GenBank	3327085
AB014565	GenBank	3327143
AB019568	GenBank	3885371
AB020623	GenBank	3985929
AB020629	GenBank	4240129
AB020693	GenBank	4240260
AB021288	GenBank	4038732
AB022663	GenBank	5019617
AB023214	GenBank	4589637
AB023230	GenBank	4589675
AC02059	N/A	N/A
AC03653	N/A	N/A
AC13415	N/A	N/A
AF000982	GenBank	2580549
AF002985	GenBank	2580585
AF005654	GenBank	2337951
AF006086	GenBank	2282037
AF007791	GenBank	3779196
AF013758	GenBank	3046899
AF013988	GenBank	2318114
AF021232	GenBank	3452182
AF026939	GenBank	2612967
AF026941	GenBank	2612970
AF026942	GenBank	2612971
AF026943	GenBank	2612972
AF026944	GenBank	2612973
AF028832	GenBank	3287488
AF030455	GenBank	3169829
AF030514	GenBank	3219692
AF031469	GenBank	4104081
AF033095	GenBank	2645728
AF035286	GenBank	2661038
AF035316	GenBank	2661078
AF037204	GenBank	2906012
AF038451	GenBank	3779225
AF038662	GenBank	3132899
AF038963	GenBank	4405794
AF043431	GenBank	3452280
AF044956	GenBank	5326827
AF045941	GenBank	3893854
AF046997	GenBank	3170363
AF051894	GenBank	3095110
AF052124	GenBank	3360431
AF052578	GenBank	2967847
AF053233	GenBank	2996191
AF054838	GenBank	2997740
AF055012	GenBank	3005735
AF061736	GenBank	4335936
AF061738	GenBank	4335940



TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AF064603	GenBank	3152659
AF064854	GenBank	4206066
AF065388	GenBank	3152700
AF067168	GenBank	4894369
AF067174	GenBank	4894381
AF067817	GenBank	3928846
AF070523	GenBank	3764088
AF070561	GenBank	3387928
AF070562	GenBank	3387930
AF070596	GenBank	3387973
AF070664	GenBank	4454703
AF070674	GenBank	3978243
AF077048	GenBank	4689143
AF077051	GenBank	4689149
AF077200	GenBank	4679013
AF077671	GenBank	3386485
AF080246	GenBank	3406799
AF081484	GenBank	3420928
AF083470	GenBank	3719293
AF084523	GenBank	3550342
AF085355	GenBank	5114044
AF086003	GenBank	3483348
AF086080	GenBank	3483425
AF086183	GenBank	3483528
AF086545	GenBank	3483890
AF091263	GenBank	4140646
AF111713	GenBank	5326796
AF118023	GenBank	4836400
AF124438	GenBank	4838431
AF124439	GenBank	4838433
AF131808	GenBank	4406640
AF131820	GenBank	4406655
AF131848	GenBank	4406690
AF132966	GenBank	4680702
AF132968	GenBank	4680706
AF146277	GenBank	4960046
AF147331	GenBank	4761682
AF150100	GenBank	5107187
AF150266	DBEst	5133702
AF151873	GenBank	4929698
AF151877	GenBank	4929706
AF151978	GenBank	5732679
AF167160	GenBank	5733691
AI023413	DBEst	3239819
AI027888	DBEst	3246587
AI031811	DBEst	3250023
AI033687	DBEst	3254640
AI042140	DBEst	3281334
AI075324	DBEst	3399895
AI075876	DBEst	3405054
AI126802	DBEst	3595316
AI127556	DBEst	3596070
AI129360	DBEst	3597874
AI139456	DBEst	3645428
AI140291	DBEst	3647748
AI144215	DBEst	3666024
AI161378	DBEst	3693062
AI188638	DBEst	3739847
AI215617	DBEst	3784658
AI216969	DBEst	3789623
AI241578	DBEst	3836975
AI250167	DBEst	3846696
AI253330	DBEst	3850451

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AI253335	DBEst	3850456
AI253369	DBEst	3850490
AI253436	DBEst	3850391
AI261671	DBEst	3869874
AI262264	DBEst	3870467
AI267162	DBEst	3886329
AI267379	DBEst	3886546
AI267502	DBEst	3886669
AI267622	DBEst	3886789
AI279131	DBEst	3917365
AI285943	DBEst	3924176
AI289173	DBEst	3932437
AI290876	DBEst	3933650
AI292104	DBEst	3934878
AI300033	DBEst	3959379
AI300074	DBEst	3959420
AI312113	DBEst	4017718
AI336032	DBEst	4072959
AI337069	DBEst	4073996
AI340262	DBEst	4077189
AI346975	DBEst	4084181
AI354639	DBEst	4094792
AI366381	DBEst	4126070
AI369024	DBEst	4147777
AI382020	DBEst	4194801
AI400372	DBEst	4243459
AI417973	DBEst	4261477
AI431963	DBEst	4306858
AI453405	DBEst	4281647
AI457157	DBEst	4310026
AI457624	DBEst	4310493
AI459679	DBEst	4312560
AI460010	DBEst	4312891
AI469095	DBEst	4331185
AI469715	DBEst	4331805
AI471539	DBEst	4333629
AI476335	DBEst	4329380
AI479289	DBEst	4372457
AI499285	DBEst	4391267
AI521180	DBEst	4435315
AI538061	DBEst	4452196
AI567204	DBEst	4525656
AI587104	DBEst	4573545
AI587328	DBEst	4573769
AI609624	DBEst	4618791
AI610607	DBEst	4619774
AI612873	DBEst	4622040
AI627444	DBEst	4664244
AI632869	DBEst	4684199
AI633164	DBEst	4684494
AI636014	DBEst	4687344
AI637620	DBEst	4689854
AI676218	DBEst	4876698
AI683871	DBEst	4894053
AI684170	DBEst	4895464
AI693877	DBEst	4971217
AI694088	DBEst	4971428
AI732534	DBEst	5053647
AI743595	DBEst	5111883
AI744489	DBEst	5112777
AI745058	DBEst	5113346
AI753108	DBEst	5131372
AI791322	DBEst	5339038

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AI798474	DBEst	5363946
AI803838	DBEst	5369310
AI811960	DBEst	5398526
AI813617	DBEst	5424832
AI815829	DBEst	5431375
AI826957	DBEst	5447628
AI831002	DBEst	5451673
AI863041	DBEst	5527148
AI867294	DBEst	5540310
AI912076	DBEst	5631931
AI915553	DBEst	5635408
AJ001381	GenBank	2764616
AJ003401	DBEst	2769433
AJ010071	GenBank	3483016
AJ132502	GenBank	5629914
AL044356	DBEst	5432578
AL044825	DBEst	5433037
AL047024	DBEst	5435080
AL048393	DBEst	5936479
AL049313	GenBank	4500086
AL049923	GenBank	4884169
AL049954	GenBank	4884203
AL050024	GenBank	4884093
AL050272	GenBank	4886498
AL050395	GenBank	4914616
AL096714	GenBank	5419847
AL096748	GenBank	5419879
AL096842	GenBank	5524930
AL110124	GenBank	5817017
C17346	DBEst	1572053
D00017	GenBank	219909
D00068	GenBank	220080
D11960	DBEst	2148277
D12502	GenBank	219494
D12763	GenBank	220076
D13380	GenBank	220033
D13645	GenBank	286008
D13866	GenBank	433410
D14697	GenBank	285964
D21260	GenBank	434760
D23660	GenBank	432358
D26155	GenBank	505086
D26599	GenBank	565648
D28759	GenBank	633074
D29640	GenBank	473930
D31763	GenBank	498151
D31767	GenBank	505091
D31883	GenBank	505093
D38524	GenBank	633070
D42040	GenBank	577292
D45248	GenBank	1008914
D49396	GenBank	682747
D50372	GenBank	2605593
D50420	GenBank	2618577
D55653	GenBank	871882
D81522	DBEst	1179399
D83077	GenBank	1304131
D83767	GenBank	1913784
D86958	GenBank	1503989
D86979	GenBank	6634000
D87666	GenBank	1620016
D87667	GenBank	1620019
D87735	GenBank	1620021

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
D88532	GenBank	1661000
D89053	GenBank	4165017
D90311	GenBank	219496
D90453	GenBank	219897
E01197	GenBank	2169456
E01198	GenBank	2169457
E01630	GenBank	2169883
E01954	GenBank	2170202
E01971	GenBank	2170219
E01972	GenBank	2170220
E02628	GenBank	2170856
E03569	GenBank	2171785
E03879	GenBank	2172093
E08663	GenBank	2176776
F06593	DBEst	672186
F28779	DBEst	4814405
H25806	DBEst	894929
H47546	DBEst	923598
H48873	DBEst	988713
H66467	DBEst	1025207
H88415	DBEst	1070675
J00196	GenBank	188242
J03575	GenBank	189737
J03858	GenBank	179439
J03909	GenBank	186264
J04164	GenBank	177801
K00422	GenBank	184322
K01763	GenBank	184316
L00693	GenBank	180228
L02426	GenBank	403455
L06328	GenBank	340200
L09159	GenBank	307374
L10413	GenBank	388755
L11066	GenBank	307322
L20688	GenBank	404044
L20941	GenBank	507251
L28997	GenBank	607027
L38995	GenBank	704415
L41490	GenBank	927064
M10119	GenBank	182517
M13536	GenBank	180248
M14328	GenBank	182113
M14764	GenBank	189204
M15329	GenBank	186277
M16660	GenBank	184420
M17017	GenBank	179579
M18216	GenBank	178690
M19723	GenBank	186726
M22918	GenBank	189019
M23613	GenBank	189271
M24194	GenBank	187701
M24594	GenBank	186262
M26152	GenBank	1160968
M29540	GenBank	180222
M29541	GenBank	189103
M29551	GenBank	180708
M33146	GenBank	181070
M34064	GenBank	416292
M34455	GenBank	185790
M35198	GenBank	9446401
M36693	GenBank	338285
M37716	GenBank	338266
M55268	GenBank	177837

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
M55542	GenBank	183001
M55543	GenBank	829176
M57567	GenBank	178986
M60333	GenBank	188268
M61715	GenBank	340367
M62831	GenBank	182260
M63121	GenBank	339755
M63838	GenBank	184568
M68520	GenBank	180177
M77945	DBEst	273682
M80563	GenBank	179916
M81757	GenBank	337732
M83248	GenBank	189150
M83654	GenBank	179660
M86553	GenBank	179958
M87284	GenBank	338651
M87434	GenBank	338653
M87503	GenBank	184652
M92357	GenBank	306463
M96982	GenBank	338262
M97501	GenBank	180621
M97935	GenBank	2281070
N36346	DBEst	1157488
N51262	DBEst	1192428
N57413	DBEst	1201303
N78477	DBEst	1241178
N92060	DBEst	1264369
Q21065	N/A	N/A
Q94780	N/A	N/A
R13925	DBEst	767001
R51732	DBEst	813634
R56461	DBEst	826567
R66489	DBEst	839127
R75621	DBEst	850303
S45630	GenBank	256398
S70290	GenBank	546602
S75295	GenBank	913392
S76638	GenBank	243420
T34641	DBEst	616739
T50925	DBEst	652785
T52715	DBEst	654575
T54951	DBEst	656812
T70793	DBEst	685314
U03886	GenBank	458225
U04313	GenBank	453368
U07550	GenBank	469170
U07857	GenBank	469048
U08815	GenBank	508722
U09559	GenBank	791184
U09847	GenBank	495565
U10439	GenBank	577169
U14966	GenBank	550012
U18321	GenBank	603763
U19878	GenBank	755465
U23942	GenBank	1698395
U25789	GenBank	808089
U28249	GenBank	897916
U28964	GenBank	899458
U32500	GenBank	1000750
U32944	GenBank	1209060
U33760	GenBank	995823
U37230	GenBank	1574941
U37518	GenBank	1149557

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
U38292	GenBank	1790879
U38784	GenBank	1574947
U41371	GenBank	1173904
U41515	GenBank	1209723
U52513	GenBank	1777781
U56255	GenBank	1399688
U57847	GenBank	1373420
U61083	GenBank	4097430
U68758	GenBank	4097815
U73524	GenBank	1644401
U77085	GenBank	1684789
U78722	GenBank	1699000
U79751	GenBank	2257753
U94586	GenBank	1946691
V00572	GenBank	35434
V00594	GenBank	37120
V04202	N/A	N/A
V17906	N/A	N/A
V36078	N/A	N/A
V68140	N/A	N/A
V86134	N/A	N/A
W02908	DBEst	1274885
W05711	DBEst	1278502
W07308	DBEst	1281506
W25547	DBEst	1303421
W28837	DBEst	1308785
W37272	DBEst	1318866
W38644	DBEst	1320349
W39262	DBEst	1320979
W39498	DBEst	1321206
W52254	DBEst	1349394
W74319	DBEst	1384468
W77987	DBEst	1388521
W80480	DBEst	1391538
X00637	GenBank	32429
X01742	GenBank	35324
X02530	GenBank	33917
X02661	GenBank	23795
X04316	N/A	N/A
X04371	GenBank	23792
X04470	GenBank	28638
X05908	GenBank	34387
X07819	GenBank	35798
X13238	GenBank	1200056
X15674	GenBank	35995
X15729	GenBank	38317
X16354	GenBank	37197
X16356	GenBank	37203
X16455	GenBank	29854
X17025	GenBank	488749
X20432	N/A	N/A
X30167	N/A	N/A
X33937	N/A	N/A
X35726	N/A	N/A
X41105	N/A	N/A
X51841	GenBank	33910
X54941	GenBank	29976
X56932	GenBank	23690
X57351	GenBank	311373
X59710	GenBank	35049
X65614	GenBank	36177
X67951	GenBank	287640
X68060	GenBank	37230

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
X68277	GenBank	29980
X72790	GenBank	311401
X76488	GenBank	434305
X83544	GenBank	1089849
X85134	GenBank	755749
X87949	GenBank	1143491
X93036	GenBank	1085025
X99699	GenBank	1869900
X99920	GenBank	1694827
Y09267	GenBank	1834492
Y13323	GenBank	5042231
Y17392	GenBank	3212109
Z12830	GenBank	551637
Z36815	GenBank	533929
Z47087	GenBank	860989
Z48570	GenBank	695580
Z71389	GenBank	2239127
AA002223	DBEst	1445158
AA018843	DBEst	1482235
AA021647	DBEst	1485308
AA022842	DBEst	1487015
AA022965	DBEst	1487064
AA024522	DBEst	1489238
AA028164	DBEst	1494289
AA035775	DBEst	1507603
AA037294	DBEst	1512438
AA039967	DBEst	1516280
AA045637	DBEst	1525513
AA046815	DBEst	1524920
AA046853	DBEst	1524752
AA047052	DBEst	1524950
AA047213	DBEst	1525113
AA057071	DBEst	1549810
AA058933	DBEst	1551788
AA064952	DBEst	1559216
AA075089	DBEst	1615078
AA076291	DBEst	1616160
AA078508	DBEst	1837982
AA080864	DBEst	1623371
AA083345	DBEst	1625405
AA083693	DBEst	1625753
AA085497	DBEst	1628765
AA086463	DBEst	1629080
AA093935	DBEst	1639528
AA100291	DBEst	1646582
AA101207	DBEst	1647860
AA102403	DBEst	1647188
AA111856	DBEst	1663943
AA115174	DBEst	1670371
AA122134	DBEst	1678255
AA122291	DBEst	1678547
AA125780	DBEst	1685521
AA127322	DBEst	1686638
AA130432	DBEst	1691715
AA131801	DBEst	1693290
AA132445	DBEst	1694012
AA134109	DBEst	1691321
AA135924	DBEst	1697110
AA136322	DBEst	1697597
AA143034	DBEst	1712411
AA150057	DBEst	1721279
AA151651	DBEst	1720206
AA156335	DBEst	1727969

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AA157333	DBEst	1728942
AA158987	DBEst	1733823
AA165439	DBEst	1741455
AA165632	DBEst	1741665
AA166618	DBEst	1745207
AA172067	DBEst	1751125
AA173031	DBEst	1754310
AA178870	DBEst	1760393
AA181874	DBEst	1765359
AA195194	DBEst	1784884
AA203206	DBEst	1798916
AA203289	DBEst	1799038
AA204768	DBEst	1802618
AA206621	DBEst	1802009
AA213914	DBEst	1812716
AA218919	DBEst	1832993
AA224050	DBEst	1844591
AA224244	DBEst	1844769
AA227596	DBEst	1849140
AA229018	DBEst	1851983
AA229161	DBEst	1851090
AA236445	DBEst	1858734
AA236680	DBEst	1860973
AA243537	DBEst	1874328
AA252436	DBEst	1887407
AA252869	DBEst	1885537
AA256330	DBEst	1891867
AA262700	DBEst	1898112
AA278358	DBEst	1921666
AA287076	DBEst	1934137
AA291551	DBEst	1939545
AA293273	DBEst	1941423
AA295982	DBEst	1948378
AA301675	DBEst	1954018
AA301722	DBEst	1954065
AA302964	DBEst	1955294
AA303199	DBEst	1955604
AA304927	DBEst	1957254
AA305042	DBEst	1957368
AA305635	DBEst	1957960
AA315030	DBEst	1967520
AA315943	DBEst	1968272
AA317144	DBEst	1969699
AA326060	DBEst	1978315
AA327358	DBEst	1979623
AA336387	DBEst	1988636
AA346413	DBEst	1998651
AA352580	DBEst	2004900
AA363162	DBEst	2015480
AA375754	DBEst	2028074
AA399230	DBEst	2053028
AA400249	DBEst	2054315
AA401629	DBEst	2055827
AA402885	DBEst	2056782
AA406401	DBEst	2064402
AA421682	DBEst	2100499
AA422057	DBEst	2100890
AA424445	DBEst	2103415
AA424901	DBEst	2107006
AA424984	DBEst	2107137
AA425182	DBEst	2105974
AA428607	DBEst	2112800
AA446099	DBEst	2158764



TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AA446403	DBEst	2159068
AA447735	DBEst	2161405
AA449054	DBEst	2163074
AA449205	DBEst	2162668
AA449520	DBEst	2163270
AA452273	DBEst	2165942
AA455007	DBEst	2177783
AA455104	DBEst	2177880
AA459527	DBEst	2184434
AA460226	DBEst	2185042
AA461287	DBEst	2186407
AA464526	DBEst	2189410
AA468398	DBEst	2194932
AA469135	DBEst	2195669
AA469453	DBEst	2194248
AA470690	DBEst	2197999
AA479427	DBEst	2207983
AA480336	DBEst	2208487
AA483454	DBEst	2212267
AA487669	DBEst	2217833
AA488423	DBEst	2215854
AA488635	DBEst	2216066
AA488843	DBEst	2218445
AA489772	DBEst	2220656
AA503972	DBEst	2238939
AA508506	DBEst	2246009
AA513550	DBEst	2251962
AA513783	DBEst	2252204
AA514989	DBEst	2254589
AA516400	DBEst	2253762
AA520993	DBEst	2261536
AA521110	DBEst	2261653
AA523639	DBEst	2264567
AA523697	DBEst	2264625
AA528106	DBEst	2270175
AA528190	DBEst	2270259
AA528226	DBEst	2270295
AA534830	DBEst	2279083
AA548722	DBEst	2319004
AA551236	DBEst	2321488
AA551243	DBEst	2321495
AA558778	DBEst	2329255
AA563834	DBEst	2335473
AA576432	DBEst	2353932
AA580069	DBEst	2355396
AA580294	DBEst	2355621
AA582588	DBEst	2359948
AA584304	DBEst	2368913
AA588772	DBEst	2402503
AA593075	DBEst	2408837
AA595585	DBEst	2410935
AA601895	DBEst	2436048
AA628700	DBEst	2541087
AA630326	DBEst	2552937
AA630642	DBEst	2553253
AA631178	DBEst	2553789
AA631218	DBEst	2553829
AA633550	DBEst	2556764
AA634808	DBEst	2558022
AA639199	DBEst	2562978
AA639791	DBEst	2563570
AA644273	DBEst	2569491
AA648897	DBEst	2575326

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AA664732	DBEst	2619345
AA677550	DBEst	2658072
AA687308	DBEst	2675499
AA705002	DBEst	2714920
AA706685	DBEst	2716603
AA708266	DBEst	2718184
AA713687	DBEst	2725961
AA719618	DBEst	2732717
AA719674	DBEst	2732773
AA720572	DBEst	2736707
AA721752	DBEst	2737814
AA723612	DBEst	2741319
AA730571	DBEst	2751775
AA742282	DBEst	2784282
AA748437	DBEst	2788395
AA749187	DBEst	2789145
AA761602	DBEst	2810532
AA768355	DBEst	2819370
AA769127	DBEst	2820365
AA774030	DBEst	2825919
AA774247	DBEst	2825545
AA779631	DBEst	2838962
AA808747	DBEst	2878153
AA809854	DBEst	2879260
AA810859	DBEst	2880470
AA825673	DBEst	2898985
AA825768	DBEst	2899080
AA826517	DBEst	2898339
AA827331	DBEst	2899772
AA827764	DBEst	2901323
AA829511	DBEst	2902610
AA831603	DBEst	2904702
AA836991	DBEst	2912190
AA837254	DBEst	2912453
AA846480	DBEst	2932620
AA846840	DBEst	2932980
AA853515	DBEst	2940254
AA883212	DBEst	2992742
AA886885	DBEst	3001993
AA889485	DBEst	3016364
AA897461	DBEst	3034081
AA902582	DBEst	3037705
AA902644	DBEst	3037767
AA909144	DBEst	3048549
AA913281	DBEst	3052673
AA916756	DBEst	3056148
AA922420	DBEst	3069729
AA927283	DBEst	3076180
AA933075	DBEst	3087008
AA935979	DBEst	3093136
AA937947	DBEst	3096058
AA948295	DBEst	3109548
AA969131	DBEst	3144311
AA971881	DBEst	3147171
AA973019	DBEst	3148199
AA988923	DBEst	3174494
AA989465	DBEst	3174829
AA994023	DBEst	3180568
AB002310	GenBank	2224564
AB002330	GenBank	2224604
AB007944	GenBank	3413911
AB012911	GenBank	3062802
AB017019	GenBank	4512256

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AB018266	GenBank	3882166
AB018305	GenBank	3882244
AB018347	GenBank	3882328
AB019568	GenBank	3885371
AB023158	GenBank	4589525
AB028976	GenBank	5689442
AB029005	GenBank	5689500
AC28164	N/A	N/A
AD001528	GenBank	2198556
AF000231	GenBank	2149974
AF006088	GenBank	2282041
AF006516	GenBank	2245670
AF012072	GenBank	2895096
AF026947	GenBank	2736255
AF028832	GenBank	3287488
AF030424	GenBank	2623155
AF031379	GenBank	4894208
AF035287	GenBank	2661040
AF035309	GenBank	2661070
AF038197	GenBank	2795918
AF038404	GenBank	2707904
AF043431	GenBank	3452280
AF044670	GenBank	4191318
AF044958	GenBank	4164447
AF047184	GenBank	2909859
AF052164	GenBank	3360475
AF052496	DBEst	3090893
AF052578	GenBank	2967847
AF054990	GenBank	3005703
AF059524	GenBank	4091867
AF070561	GenBank	3387928
AF070626	GenBank	3283892
AF070655	GenBank	4454685
AF070674	GenBank	3978243
AF075040	GenBank	3377580
AF077030	GenBank	4689107
AF078847	GenBank	5531808
AF080246	GenBank	3406799
AF081282	GenBank	4336324
AF081484	GenBank	3420928
AF084523	GenBank	3550342
AF086163	GenBank	3483508
AF095791	GenBank	3777595
AF100756	GenBank	5410297
AF107406	GenBank	5531905
AF119297	GenBank	4633508
AF131858	GenBank	4406705
AF132940	GenBank	4680650
AF151857	GenBank	4929666
AI028733	DBEst	3246042
AI031901	DBEst	3250113
AI033739	DBEst	3254692
AI040324	DBEst	3279518
AI051172	DBEst	3306706
AI076805	DBEst	3404634
AI087005	DBEst	3425428
AI089913	DBEst	3428972
AI092007	DBEst	3427205
AI127326	DBEst	3595840
AI147251	DBEst	3674933
AI148933	DBEst	3677402
AI149846	DBEst	3678315
AI167855	DBEst	3701025

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AI183965	DBEst	3734603
AI189258	DBEst	3740467
AI220148	DBEst	3802351
AI224374	DBEst	3807087
AI240095	DBEst	3835492
AI246677	DBEst	3842074
AI248538	DBEst	3843935
AI266582	DBEst	3884740
AI268864	DBEst	3888031
AI270183	DBEst	3889350
AI271795	DBEst	3890962
AI273008	DBEst	3895276
AI273841	DBEst	3896109
AI274756	DBEst	3897030
AI275528	DBEst	3897802
AI283096	DBEst	3921329
AI298059	DBEst	3957795
AI335653	DBEst	4072580
AI338977	DBEst	4075904
AI339946	DBEst	4076873
AI373032	DBEst	4152898
AI374954	DBEst	4174944
AI374954	DBEst	4174944
AI380539	DBEst	4190392
AI417583	DBEst	4261087
AI432644	DBEst	4283347
AI433157	DBEst	4287209
AI457792	DBEst	4310661
AI469112	DBEst	4331202
AI471114	DBEst	4333204
AI471534	DBEst	4333624
AI473927	DBEst	4326972
AI479305	DBEst	4372473
AI499243	DBEst	4391225
AI525796	DBEst	4439931
AI525843	DBEst	4439978
AI537677	DBEst	4451812
AI541029	DBEst	4458402
AI560129	DBEst	4510470
AI583108	DBEst	4569005
AI584068	DBEst	4569965
AI587208	DBEst	4573649
AI589867	DBEst	4598915
AI610676	DBEst	4619843
AI630362	DBEst	4681692
AI633006	DBEst	4684336
AI634443	DBEst	4685773
AI635096	DBEst	4686426
AI682105	DBEst	4892287
AI683338	DBEst	4893520
AI684800	DBEst	4896094
AI684991	DBEst	4896285
AI689369	DBEst	4900663
AI689617	DBEst	4900911
AI689883	DBEst	4901177
AI693745	DBEst	4971085
AI701001	DBEst	4988901
AI733038	DBEst	5054151
AI735638	DBEst	5057162
AI741506	DBEst	5109794
AI742722	DBEst	5111010
AI742738	DBEst	5111026
AI743552	DBEst	5111840

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AI753784	DBEst	5132136
AI754296	DBEst	5132560
AI754652	DBEst	5132916
AI754732	DBEst	5132996
AI765975	DBEst	5232484
AI769970	DBEst	5236479
AI819225	DBEst	5438304
AI820563	DBEst	5439642
AI827818	DBEst	5448489
AI828682	DBEst	5449353
AI830067	DBEst	5450738
AI861989	DBEst	5526096
AI887129	DBEst	5592293
AI887632	DBEst	5592796
AI890281	DBEst	5595445
AI924046	DBEst	5660010
AI924096	DBEst	5660060
AI924823	DBEst	5660787
AI963471	DBEst	5756184
AI963604	DBEst	5756382
AI972556	DBEst	5769302
AI979048	DBEst	5804078
AI984656	DBEst	5811933
AJ010442	GenBank	3954884
AJ132694	GenBank	4454210
AJ224442	GenBank	2911586
AL036299	DBEst	5405889
AL042979	DBEst	5422409
AL047305	DBEst	4727252
AL049247	GenBank	4499985
AL049313	GenBank	4500086
AL049381	GenBank	4500168
AL049932	GenBank	4884176
AL050041	GenBank	4884283
AL050161	GenBank	4884375
AL050265	GenBank	4886440
AL050268	GenBank	4886442
AL050367	GenBank	4914600
AL079286	GenBank	5102746
AL079312	GenBank	5102890
AL079314	GenBank	5102893
AL080113	GenBank	5262540
AL110164	GenBank	5817069
AL117412	GenBank	5912102
AL117612	GenBank	5912188
AL119009	DBEst	5924908
AW014693	DBEst	5863450
AW014985	DBEst	5863742
AW021794	DBEst	5875324
C01521	DBEst	1433751
D01096	GenBank	220128
D13119	GenBank	285909
D13627	GenBank	286010
D13630	GenBank	286000
D13639	GenBank	285990
D13665	GenBank	393318
D14530	GenBank	414348
D21260	GenBank	434760
D25278	GenBank	434780
D26361	GenBank	452516
D30655	GenBank	485387
D50310	GenBank	1183161
D51497	DBEst	951733

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
D53031	DBEst	954928
D62116	DBEst	965892
D63878	GenBank	961447
D78611	GenBank	1655421
D82348	GenBank	1311461
D83032	GenBank	1374697
D85433	GenBank	1841371
D87437	GenBank	1665768
D87667	GenBank	1620019
D89092	GenBank	2780747
D90041	GenBank	219413
E02628	GenBank	2170856
E05732	GenBank	2173919
F00551	DBEst	707254
H08920	DBEst	873742
H25080	DBEst	893979
H30306	DBEst	901216
H44647	DBEst	920699
H81376	DBEst	1059465
H93521	DBEst	1099849
H94496	DBEst	1102129
J03464	GenBank	179595
J03799	GenBank	186840
J04027	GenBank	950413
J04177	GenBank	179729
K01228	GenBank	180391
K01566	GenBank	187721
L07395	GenBank	190218
L09159	GenBank	307374
L11315	GenBank	403386
L13806	GenBank	306554
L15702	GenBank	291921
L16510	GenBank	291887
L24804	GenBank	438651
L25931	GenBank	438638
L28809	GenBank	454151
M10036	GenBank	339840
M10905	GenBank	182696
M11353	GenBank	184092
M12267	GenBank	189328
M13536	GenBank	180248
M14483	GenBank	339692
M14630	GenBank	339690
M17885	GenBank	190231
M18366	GenBank	179131
M21575	GenBank	1311702
M23254	GenBank	511636
M24194	GenBank	187701
M24486	GenBank	190785
M26512	GenBank	177796
M28372	GenBank	643575
M31159	GenBank	183115
M32220	GenBank	186619
M36341	GenBank	178984
M36693	GenBank	338285
M38690	GenBank	1048988
M58485	GenBank	180154
M59849	GenBank	182591
M62831	GenBank	182260
M64241	GenBank	190813
M69043	GenBank	187290
M77142	GenBank	339700
M77830	GenBank	4689438

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
M86667	GenBank	189066
M88108	GenBank	189499
M93651	GenBank	338038
M95542	GenBank	184271
N43970	DBEst	1182498
Q12759	N/A	N/A
Q14635	N/A	N/A
R11045	DBEst	763780
R76376	DBEst	851058
R84450	DBEst	942856
S74728	GenBank	797409
S82081	GenBank	1488412
T07459	DBEst	318608
T19883	DBEst	597628
T21168	DBEst	2596291
T22605	DBEst	2597187
T37405	DBEst	621222
T67129	DBEst	676569
T69703	DBEst	680851
T78615	DBEst	697124
T89937	DBEst	718450
U03851	GenBank	433307
U12404	GenBank	531170
U14967	GenBank	550014
U14971	GenBank	550022
U20659	GenBank	929920
U25789	GenBank	808089
U30825	GenBank	1049077
U47077	GenBank	9027566
U49844	GenBank	1235901
U63846	GenBank	1480921
U65928	GenBank	1549382
U72516	GenBank	1673521
U79282	GenBank	1710254
U90716	GenBank	1946350
U90904	GenBank	1913882
U94364	GenBank	2618769
V20437	N/A	N/A
V24305	N/A	N/A
V81394	N/A	N/A
V84510	N/A	N/A
W19427	DBEst	1295328
W65357	DBEst	1373499
W75963	DBEst	1386337
W80525	DBEst	1391689
X01630	GenBank	28871
X04098	GenBank	28338
X04408	GenBank	31914
X06700	GenBank	30053
X14420	GenBank	30057
X51742	GenBank	52697
X60111	GenBank	34768
X69398	GenBank	396175
X72755	GenBank	311375
X74979	GenBank	400462
X76180	GenBank	452649
X78627	GenBank	607129
X79067	GenBank	483524
X80910	GenBank	531475
X87949	GenBank	1143491
Y00052	GenBank	30308
Y00062	GenBank	34275
Y00282	GenBank	36048

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
Y00503	GenBank	34038
Y15286	GenBank	2584788
Y17171	GenBank	3093333
Z13009	GenBank	31072
Z24724	GenBank	505034
Z29083	GenBank	435654
Z29331	GenBank	483539
Z46606	GenBank	575250
Z48501	GenBank	693936
AA001460	DBEst	1436925
AA001543	DBEst	1437008
AA001792	DBEst	1445606
AA004925	DBEst	1448503
AA010897	DBEst	1471994
AA017162	DBEst	1479361
AA019019	DBEst	1482429
AA022980	DBEst	1487079
AA024595	DBEst	1489500
AA024940	DBEst	1489864
AA024996	DBEst	1489901
AA025750	DBEst	1491134
AA026598	DBEst	1492433
AA029271	DBEst	1496712
AA029725	DBEst	1497138
AA029930	DBEst	1496355
AA033832	DBEst	1505650
AA035471	DBEst	1507128
AA035616	DBEst	1507426
AA036752	DBEst	1509790
AA037377	DBEst	1512540
AA039778	DBEst	1516057
AA039948	DBEst	1516243
AA040688	DBEst	1517002
AA040820	DBEst	1517098
AA041259	DBEst	1517683
AA043477	DBEst	1521333
AA044209	DBEst	1522066
AA044233	DBEst	1522109
AA044791	DBEst	1522994
AA045054	DBEst	1523256
AA045147	DBEst	1523487
AA045768	DBEst	1525870
AA046848	DBEst	1524747
AA053021	DBEst	1544277
AA053316	DBEst	1545775
AA053919	DBEst	1544863
AA054069	DBEst	1545012
AA055479	DBEst	1547884
AA055591	DBEst	1547956
AA055637	DBEst	1547976
AA057243	DBEst	1550096
AA058712	DBEst	1551520
AA059128	DBEst	1552146
AA065169	DBEst	1559064
AA069850	DBEst	1577210
AA071167	DBEst	1578528
AA075158	DBEst	1615214
AA075515	DBEst	1615385
AA075663	DBEst	1615533
AA076397	DBEst	1616448
AA076421	DBEst	1616290
AA078387	DBEst	1837861
AA078570	DBEst	1838051



TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AA078872	DBEst	1617825
AA079480	DBEst	1618390
AA080889	DBEst	1623378
AA081073	DBEst	1622991
AA081608	DBEst	1623666
AA081834	DBEst	1623893
AA081917	DBEst	1623975
AA082258	DBEst	1624341
AA082441	DBEst	1624498
AA083270	DBEst	1625391
AA083345	DBEst	1625405
AA083522	DBEst	1625582
AA083573	DBEst	1625633
AA083638	DBEst	1625697
AA083774	DBEst	1625832
AA088318	DBEst	1633822
AA088344	DBEst	1633856
AA088351	DBEst	1633889
AA088693	DBEst	1634214
AA088783	DBEst	1634295
AA088829	DBEst	1634323
AA090106	DBEst	1636590
AA096032	DBEst	1641617
AA099819	DBEst	1645918
AA099923	DBEst	1646071
AA099976	DBEst	1646109
AA100764	DBEst	1647117
AA101010	DBEst	1647531
AA102013	DBEst	1645759
AA102564	DBEst	1647756
AA102830	DBEst	1648675
AA112186	DBEst	1664473
AA112645	DBEst	1665346
AA113305	DBEst	1665010
AA115218	DBEst	1670047
AA115315	DBEst	1670632
AA121656	DBEst	1679269
AA121718	DBEst	1679447
AA125809	DBEst	1688020
AA125939	DBEst	1687931
AA126452	DBEst	1686119
AA126718	DBEst	1686236
AA127436	DBEst	1686832
AA127666	DBEst	1686935
AA128063	DBEst	1687342
AA128636	DBEst	1688579
AA128641	DBEst	1688584
AA130778	DBEst	1692444
AA130982	DBEst	1692473
AA131827	DBEst	1693380
AA132056	DBEst	1693545
AA132163	DBEst	1693672
AA132574	DBEst	1694081
AA132992	DBEst	1694561
AA133351	DBEst	1690319
AA133474	DBEst	1690442
AA134460	DBEst	1692042
AA134527	DBEst	1692092
AA134589	DBEst	1695586
AA135696	DBEst	1696707
AA137017	DBEst	1698226
AA142941	DBEst	1712319
AA143001	DBEst	1712506

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AA143074	DBEst	1712578
AA143746	DBEst	1713134
AA146900	DBEst	1716290
AA147200	DBEst	1716573
AA147247	DBEst	1716883
AA147781	DBEst	1717161
AA148027	DBEst	1717451
AA148136	DBEst	1717510
AA149810	DBEst	1720917
AA150377	DBEst	1721908
AA150837	DBEst	1722412
AA150928	DBEst	1722439
AA151274	DBEst	1719600
AA151594	DBEst	1720081
AA151755	DBEst	1720310
AA152476	DBEst	1718704
AA155754	DBEst	1727371
AA156066	DBEst	1727700
AA157163	DBEst	1728787
AA157993	DBEst	1732804
AA158738	DBEst	1733549
AA159110	DBEst	1733921
AA159576	DBEst	1735127
AA161003	DBEst	1735290
AA161076	DBEst	1735364
AA161467	DBEst	1735906
AA164193	DBEst	1741344
AA164473	DBEst	1740650
AA164729	DBEst	1740889
AA164873	DBEst	1741032
AA165027	DBEst	1740273
AA165068	DBEst	1740296
AA165087	DBEst	1740315
AA165174	DBEst	1740402
AA165282	DBEst	1740510
AA165293	DBEst	1740521
AA165638	DBEst	1741671
AA166618	DBEst	1745207
AA167041	DBEst	1745434
AA167750	DBEst	1744900
AA171630	DBEst	1750889
AA173506	DBEst	1753638
AA174097	DBEst	1754302
AA179187	DBEst	1760556
AA180137	DBEst	1761403
AA180224	DBEst	1761553
AA180383	DBEst	1761692
AA181075	DBEst	1764592
AA181258	DBEst	1764785
AA181684	DBEst	1765213
AA182415	DBEst	1766238
AA182540	DBEst	1766256
AA186577	DBEst	1774676
AA187817	DBEst	1774011
AA188045	DBEst	1774295
AA188140	DBEst	1774332
AA188384	DBEst	1775418
AA188826	DBEst	1775853
AA190873	DBEst	1779393
AA191422	DBEst	1780101
AA192094	DBEst	1782111
AA193308	DBEst	1782719
AA194577	DBEst	1784338

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AA195246	DBEst	1784946
AA195865	DBEst	1791530
AA196424	DBEst	1791999
AA196982	DBEst	1792573
AA203691	DBEst	1799458
AA204867	DBEst	1802927
AA206578	DBEst	1801958
AA206991	DBEst	1801246
AA209508	DBEst	1807460
AA216753	DBEst	1817452
AA219665	DBEst	1833722
AA223121	DBEst	1843680
AA223820	DBEst	1844362
AA224109	DBEst	1844668
AA224407	DBEst	1845029
AA227118	DBEst	1848672
AA229325	DBEst	1851167
AA229611	DBEst	1851608
AA232959	DBEst	1855951
AA233835	DBEst	1856856
AA233843	DBEst	1856920
AA234092	DBEst	1858897
AA234307	DBEst	1858618
AA236776	DBEst	1860841
AA242985	DBEst	1873780
AA243338	DBEst	1874149
AA244342	DBEst	1875177
AA249154	DBEst	1879783
AA255502	DBEst	1892406
AA256591	DBEst	1892130
AA261990	DBEst	1897971
AA262939	DBEst	1898659
AA278445	DBEst	1919782
AA278482	DBEst	1919801
AA278642	DBEst	1919962
AA278956	DBEst	1920495
AA279048	DBEst	1920577
AA280099	DBEst	1921573
AA280221	DBEst	1921759
AA280828	DBEst	1923508
AA282915	DBEst	1925910
AA284334	DBEst	1928614
AA284555	DBEst	1927484
AA284670	DBEst	1927581
AA284671	DBEst	1927582
AA284870	DBEst	1927464
AA284906	DBEst	1927448
AA285290	DBEst	1929600
AA286699	DBEst	1933581
AA286872	DBEst	1933932
AA287219	DBEst	1934280
AA287642	DBEst	1933325
AA287815	DBEst	1933514
AA291438	DBEst	1939417
AA291485	DBEst	1939506
AA291971	DBEst	1940027
AA292334	DBEst	1940314
AA293127	DBEst	1941167
AA293133	DBEst	1941173
AA293273	DBEst	1941423
AA293286	DBEst	1941377
AA293353	DBEst	1940750
AA293572	DBEst	1941239

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AA293629	DBEst	1941280
AA293759	DBEst	1941542
AA293804	DBEst	1941727
AA296780	DBEst	1949277
AA297402	DBEst	1949735
AA298505	DBEst	1950908
AA299640	DBEst	1951971
AA301062	DBEst	1953433
AA301800	DBEst	1954133
AA303461	DBEst	1955795
AA303568	DBEst	1955901
AA306718	DBEst	1959046
AA306862	DBEst	1959190
AA306876	DBEst	1959204
AA307198	DBEst	1959526
AA307325	DBEst	1959653
AA308065	DBEst	1960394
AA308274	DBEst	1960673
AA308744	DBEst	1961143
AA310739	DBEst	1963088
AA310771	DBEst	1963242
AA311228	DBEst	1963628
AA311460	DBEst	1963786
AA311571	DBEst	1964055
AA311801	DBEst	1964150
AA311848	DBEst	1964177
AA311905	DBEst	1964306
AA312218	DBEst	1964618
AA312240	DBEst	1964578
AA312435	DBEst	1964763
AA313108	DBEst	1965456
AA313223	DBEst	1965552
AA313653	DBEst	1965983
AA313994	DBEst	1966555
AA314431	DBEst	1966760
AA314872	DBEst	1967221
AA315363	DBEst	1967742
AA315379	DBEst	1967707
AA317243	DBEst	1969570
AA317393	DBEst	1969772
AA318969	DBEst	1971371
AA327201	DBEst	1979467
AA331991	DBEst	1984254
AA332672	DBEst	1984936
AA333358	DBEst	1985601
AA335273	DBEst	1987516
AA336666	DBEst	1988905
AA337192	DBEst	1989429
AA337489	DBEst	1989954
AA338793	DBEst	1991103
AA339957	DBEst	1992224
AA340341	DBEst	1992579
AA341446	DBEst	1993684
AA341465	DBEst	1993733
AA342969	DBEst	1995205
AA343629	DBEst	1995868
AA344084	DBEst	1996321
AA345329	DBEst	1997564
AA346393	DBEst	1998631
AA346698	DBEst	1999168
AA347887	DBEst	2000122
AA350059	DBEst	2002398
AA351507	DBEst	2003827

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AA355003	DBEst	2007559
AA356682	DBEst	2009000
AA357574	DBEst	2009894
AA358887	DBEst	2011207
AA359705	DBEst	2012096
AA364352	DBEst	2016692
AA367451	DBEst	2019769
AA367773	DBEst	2020090
AA368542	DBEst	2021085
AA369400	DBEst	2021719
AA373230	DBEst	2025550
AA374754	DBEst	2027074
AA375312	DBEst	2027642
AA375815	DBEst	2028133
AA393525	DBEst	2046493
AA394115	DBEst	2047129
AA398443	DBEst	2051755
AA398585	DBEst	2051827
AA398739	DBEst	2051861
AA399165	DBEst	2052960
AA399628	DBEst	2052642
AA401329	DBEst	2053554
AA401334	DBEst	2053559
AA402191	DBEst	2056138
AA402289	DBEst	2056202
AA402775	DBEst	2056528
AA403319	DBEst	2056820
AA404613	DBEst	2058825
AA405124	DBEst	2063536
AA406239	DBEst	2064220
AA410580	DBEst	2069686
AA410982	DBEst	2070088
AA411021	DBEst	2070171
AA411252	DBEst	2068793
AA411764	DBEst	2070352
AA417794	DBEst	2079604
AA419263	DBEst	2078976
AA419284	DBEst	2079014
AA420751	DBEst	2094630
AA420758	DBEst	2094637
AA421248	DBEst	2100135
AA421682	DBEst	2100499
AA422060	DBEst	2100893
AA422143	DBEst	2101011
AA425004	DBEst	2107073
AA425468	DBEst	2106322
AA425737	DBEst	2107249
AA429794	DBEst	2113001
AA430400	DBEst	2110974
AA430436	DBEst	2110992
AA431428	DBEst	2115136
AA433988	DBEst	2138902
AA436315	DBEst	2141229
AA436411	DBEst	2141325
AA443024	DBEst	2155699
AA449394	DBEst	2162785
AA451779	DBEst	2165448
AA453878	DBEst	2167547
AA454668	DBEst	2177444
AA454953	DBEst	2177729
AA454962	DBEst	2177738
AA455245	DBEst	2178021
AA455785	DBEst	2178561

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AA456454	DBEst	2179030
AA456557	DBEst	2179133
AA457255	DBEst	2179975
AA457579	DBEst	2180299
AA459167	DBEst	2184074
AA459210	DBEst	2184117
AA459527	DBEst	2184434
AA460570	DBEst	2185690
AA460816	DBEst	2185936
AA461005	DBEst	2186125
AA468657	DBEst	2195191
AA469447	DBEst	2194242
AA469453	DBEst	2194248
AA476522	DBEst	2204733
AA477018	DBEst	2205229
AA477567	DBEst	2206201
AA477973	DBEst	2206607
AA478230	DBEst	2206864
AA479646	DBEst	2205532
AA479648	DBEst	2205534
AA479848	DBEst	2205734
AA481078	DBEst	2210630
AA481710	DBEst	2211262
AA482430	DBEst	2210108
AA482432	DBEst	2210110
AA482779	DBEst	2211624
AA483258	DBEst	2212071
AA483726	DBEst	2212539
AA483858	DBEst	2212671
AA484181	DBEst	2212994
AA486047	DBEst	2216263
AA486859	DBEst	2217023
AA488141	DBEst	2215572
AA488385	DBEst	2215816
AA488517	DBEst	2215948
AA489323	DBEst	2218925
AA489380	DBEst	2218982
AA489382	DBEst	2218984
AA491204	DBEst	2220377
AA492143	DBEst	2221705
AA493371	DBEst	2223212
AA494321	DBEst	2224108
AA494552	DBEst	2224339
AA501657	DBEst	2236624
AA502136	DBEst	2237103
AA505780	DBEst	2241917
AA512933	DBEst	2251356
AA514395	DBEst	2253995
AA514974	DBEst	2254574
AA515143	DBEst	2254743
AA516376	DBEst	2253738
AA521006	DBEst	2261549
AA523522	DBEst	2264234
AA524748	DBEst	2265676
AA524950	DBEst	2265878
AA525141	DBEst	2266069
AA526028	DBEst	2268097
AA527275	DBEst	2269344
AA527557	DBEst	2269626
AA533506	DBEst	2277602
AA534349	DBEst	2278602
AA534586	DBEst	2278839
AA534608	DBEst	2278861

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AA535496	DBEst	2279749
AA541651	DBEst	2288085
AA548056	DBEst	2318338
AA548600	DBEst	2318882
AA550854	DBEst	2321106
AA550855	DBEst	2321107
AA551351	DBEst	2321603
AA551391	DBEst	2321643
AA554437	DBEst	2324976
AA554735	DBEst	2325274
AA555102	DBEst	2325641
AA564272	DBEst	2335911
AA564870	DBEst	2336509
AA565420	DBEst	2337059
AA568936	DBEst	2341990
AA569816	DBEst	2343796
AA569851	DBEst	2343831
AA569916	DBEst	2343896
AA573761	DBEst	2348276
AA573787	DBEst	2348302
AA577537	DBEst	2355011
AA578881	DBEst	2357065
AA579591	DBEst	2357775
AA579890	DBEst	2355217
AA580835	DBEst	2358607
AA582093	DBEst	2360771
AA582866	DBEst	2360226
AA583055	DBEst	2360415
AA583498	DBEst	2368107
AA583567	DBEst	2368176
AA583773	DBEst	2368382
AA584921	DBEst	2367701
AA586755	DBEst	2397569
AA587140	DBEst	2397954
AA587315	DBEst	2398129
AA587873	DBEst	2402048
AA593983	DBEst	2409333
AA594366	DBEst	2409716
AA595624	DBEst	2410974
AA595771	DBEst	2411121
AA599454	DBEst	2433079
AA600227	DBEst	2433852
AA600771	DBEst	2434396
AA601172	DBEst	2434797
AA602395	DBEst	2436373
AA602871	DBEst	2436805
AA603125	DBEst	2436986
AA603177	DBEst	2437038
AA604324	DBEst	2445233
AA604853	DBEst	2445717
AA610279	DBEst	2458707
AA610476	DBEst	2458904
AA610734	DBEst	2459162
AA614482	DBEst	2466678
AA628536	DBEst	2540923
AA628547	DBEst	2540934
AA630611	DBEst	2553222
AA631326	DBEst	2553937
AA633909	DBEst	2557123
AA634260	DBEst	2557474
AA634298	DBEst	2557512
AA640505	DBEst	2565755
AA641289	DBEst	2566539

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AA644625	DBEst	2569843
AA648944	DBEst	2575373
AA651720	DBEst	2583372
AA652478	DBEst	2584130
AA652505	DBEst	2584157
AA653775	DBEst	2589929
AA658374	DBEst	2594528
AA663005	DBEst	2616996
AA669154	DBEst	2630653
AA677560	DBEst	2658082
AA677750	DBEst	2658272
AA678185	DBEst	2658707
AA678251	DBEst	2658773
AA687495	DBEst	2675686
AA703208	DBEst	2706321
AA703667	DBEst	2713585
AA703907	DBEst	2713825
AA704208	DBEst	2714126
AA706347	DBEst	2716265
AA714010	DBEst	2726284
AA715984	DBEst	2728258
AA716651	DBEst	2728925
AA719530	DBEst	2732629
AA721642	DBEst	2736625
AA729381	DBEst	2750740
AA731946	DBEst	2753897
AA736817	DBEst	2768051
AA742713	DBEst	2782219
AA743278	DBEst	2782784
AA744681	DBEst	2783445
AA745953	DBEst	2785939
AA759195	DBEst	2807058
AA767779	DBEst	2818794
AA769697	DBEst	2820935
AA773998	DBEst	2825887
AA775058	DBEst	2834392
AA776593	DBEst	2835927
AA777384	DBEst	2836715
AA778672	DBEst	2838003
AA779949	DBEst	2839280
AA781487	DBEst	2840818
AA788907	DBEst	2849027
AA806278	DBEst	2875028
AA806735	DBEst	2875485
AA808769	DBEst	2878175
AA810149	DBEst	2879555
AA811609	DBEst	2881220
AA813604	DBEst	2882289
AA826307	DBEst	2899619
AA833766	DBEst	2908534
AA833900	DBEst	2907499
AA837457	DBEst	2912656
AA843531	DBEst	2930049
AA845737	DBEst	2931877
AA846698	DBEst	2932838
AA846856	DBEst	2932996
AA852896	DBEst	2939635
AA856902	DBEst	2945204
AA857824	DBEst	2946126
AA857882	DBEst	2946184
AA861665	DBEst	2953805
AA865960	DBEst	2958236
AA868529	DBEst	2963974



TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AA873271	DBEst	2969393
AA877189	DBEst	2986266
AA884922	DBEst	2994903
AA886453	DBEst	3001561
AA906652	DBEst	3042238
AA906865	DBEst	3042109
AA918993	DBEst	3058883
AA926926	DBEst	3075823
AA928934	DBEst	3078291
AA932501	DBEst	3087282
AA933987	DBEst	3090255
AA935947	DBEst	3093104
AA937302	DBEst	3095413
AA937773	DBEst	3095884
AA947835	DBEst	3109088
AA954939	DBEst	3118634
AA962587	DBEst	3134751
AA962632	DBEst	3134796
AA972525	DBEst	3145289
AA976489	DBEst	3152281
AA983380	DBEst	3161905
AA984586	DBEst	3163111
AA992596	DBEst	3179352
AB002305	GenBank	2224554
AB002330	GenBank	2224604
AB002357	GenBank	2224658
AB002806	GenBank	2780782
AB003476	GenBank	2081606
AB004066	GenBank	2308996
AB006077	GenBank	2564010
AB006534	GenBank	2924619
AB006755	GenBank	2979417
AB007867	GenBank	2662094
AB007900	GenBank	2662160
AB007916	GenBank	6683704
AB007923	GenBank	3413869
AB007957	GenBank	3413931
AB011103	GenBank	3043585
AB011143	GenBank	3043665
AB011151	GenBank	3043681
AB011166	GenBank	3043711
AB014533	GenBank	3327079
AB014542	GenBank	3327097
AB014560	GenBank	3327133
AB015630	GenBank	4586837
AB015856	GenBank	3953530
AB018281	GenBank	3882196
AB018284	GenBank	3882202
AB018285	GenBank	3882204
AB018289	GenBank	3882212
AB018305	GenBank	3882244
AB018327	GenBank	3882288
AB018331	GenBank	3882296
AB018337	GenBank	3882308
AB019409	GenBank	4587128
AB019563	GenBank	3885366
AB019568	GenBank	3885371
AB019691	GenBank	5051742
AB020682	GenBank	4240238
AB020718	GenBank	4240310
AB021288	GenBank	4038732
AB023154	GenBank	4589517
AB023219	GenBank	4589647

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AB024704	GenBank	4589928
AB027467	GenBank	6172222
AB028069	GenBank	4996095
AB028624	GenBank	5103045
AB028969	GenBank	5689428
AB028986	GenBank	5689462
AB029000	GenBank	5689490
AB029004	GenBank	5689498
AB029028	GenBank	5689546
AC03044	N/A	N/A
AC31479	N/A	N/A
AF000670	GenBank	3153911
AF000974	GenBank	2232135
AF001893	GenBank	2529723
AF004562	GenBank	3041872
AF006043	GenBank	2674061
AF007135	GenBank	2852610
AF007151	GenBank	2852629
AF007170	GenBank	2865251
AF009615	GenBank	2393946
AF013759	GenBank	3153208
AF013988	GenBank	2318114
AF015283	GenBank	2384720
AF015767	GenBank	2353176
AF016507	GenBank	2909776
AF016582	GenBank	2367668
AF017790	GenBank	2501872
AF019767	GenBank	3510461
AF021351	GenBank	2460207
AF021819	GenBank	2460317
AF022229	GenBank	2809382
AF023266	GenBank	4103447
AF025439	GenBank	2815605
AF026166	GenBank	4090928
AF026939	GenBank	2612967
AF027205	GenBank	2598967
AF031385	GenBank	2606093
AF034607	GenBank	4426566
AF035286	GenBank	2661038
AF035309	GenBank	2661070
AF035313	GenBank	2661075
AF037204	GenBank	2906012
AF038661	GenBank	3132897
AF039019	GenBank	2828109
AF039291	GenBank	4104738
AF039843	GenBank	2809399
AF040990	GenBank	2804783
AF041483	GenBank	3493528
AF042385	GenBank	2828148
AF042729	GenBank	7770717
AF044588	GenBank	2865520
AF045184	GenBank	3417598
AF047438	GenBank	3335131
AF047472	GenBank	2921872
AF048977	GenBank	3005586
AF050171	GenBank	5668577
AF050199	GenBank	2961556
AF050639	GenBank	4164453
AF052124	GenBank	3360431
AF052135	GenBank	3360444
AF052149	GenBank	3360459
AF052164	GenBank	3360475
AF052169	GenBank	3360480

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AF052180	GenBank	3360492
AF052514	GenBank	3510662
AF054183	GenBank	4092053
AF054187	GenBank	4092059
AF054840	GenBank	2997744
AF055012	GenBank	3005735
AF055033	GenBank	3005763
AF057299	GenBank	5739040
AF059252	GenBank	3372629
AF061258	GenBank	3108092
AF062318	GenBank	3152814
AF063611	GenBank	4731856
AF064019	GenBank	3347856
AF068235	GenBank	4321975
AF068846	GenBank	3201999
AF070523	GenBank	3764088
AF070537	GenBank	3387894
AF070555	GenBank	3387920
AF070561	GenBank	3387928
AF070596	GenBank	3387973
AF070600	GenBank	3387979
AF070626	GenBank	3283892
AF070649	GenBank	3283923
AF070662	GenBank	4454699
AF070672	GenBank	3978239
AF071202	GenBank	3335172
AF071219	GenBank	3288867
AF071593	GenBank	3249712
AF073298	GenBank	3641537
AF075587	GenBank	3319325
AF077030	GenBank	4689107
AF077045	GenBank	4689137
AF077200	GenBank	4679013
AF077202	GenBank	4679017
AF077207	GenBank	4679027
AF081192	GenBank	3420798
AF081484	GenBank	3420928
AF083190	GenBank	3599414
AF085355	GenBank	5114044
AF086003	GenBank	3483348
AF086116	GenBank	3483461
AF086178	GenBank	3483523
AF086205	GenBank	3483550
AF086207	GenBank	3483552
AF086336	GenBank	3483681
AF086517	GenBank	3483862
AF087135	GenBank	3641297
AF087990	GenBank	3523196
AF088036	GenBank	3523242
AF091076	GenBank	3859989
AF092563	GenBank	3851583
AF095287	GenBank	3766235
AF095791	GenBank	3777595
AF097709	GenBank	3777616
AF100741	GenBank	5138992
AF100756	GenBank	5410297
AF100928	GenBank	4323586
AF104222	GenBank	3983426
AF104913	GenBank	3941723
AF104923	GenBank	4680483
AF107405	GenBank	5531903
AF120334	GenBank	4191615
AF124438	GenBank	4838431

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AF124439	GenBank	4838433
AF125525	GenBank	4689281
AF131799	GenBank	4406628
AF131814	GenBank	4406648
AF139461	GenBank	4894945
AF139658	GenBank	4894940
AF144755	GenBank	5006628
AF147331	GenBank	4761682
AF150962	GenBank	5020252
AF151832	GenBank	4929616
AF151868	GenBank	4929688
AF151898	GenBank	4929748
AF151907	GenBank	4929766
AF152097	GenBank	4929772
AF159295	GenBank	5714635
AF176702	GenBank	6103642
AF190744	GenBank	6176531
AI004664	DBEst	3214174
AI004915	DBEst	3214425
AI016073	DBEst	3230409
AI016323	DBEst	3230659
AI016791	DBEst	3231127
AI018451	DBEst	3232970
AI018625	DBEst	3233144
AI022779	DBEst	3238020
AI023799	DBEst	3238843
AI026164	DBEst	3241777
AI027516	DBEst	3246446
AI031636	DBEst	3249848
AI033037	DBEst	3253990
AI034115	DBEst	3255068
AI037859	DBEst	3277053
AI041670	DBEst	3280864
AI042034	DBEst	3281228
AI042290	DBEst	3281484
AI051971	DBEst	3307962
AI056917	DBEst	3330706
AI057124	DBEst	3331000
AI066419	DBEst	3367121
AI078041	DBEst	3412449
AI081116	DBEst	3417908
AI081472	DBEst	3418264
AI081913	DBEst	3418705
AI082244	DBEst	3419036
AI082648	DBEst	3419440
AI084731	DBEst	3423154
AI085381	DBEst	3423804
AI087291	DBEst	3425714
AI087819	DBEst	3426852
AI088178	DBEst	3427256
AI089981	DBEst	3429040
AI090524	DBEst	3429583
AI090623	DBEst	3429682
AI091425	DBEst	3430484
AI092971	DBEst	3431947
AI095477	DBEst	3434453
AI123229	DBEst	3538995
AI125642	DBEst	3594156
AI125874	DBEst	3594388
AI127013	DBEst	3595527
AI127556	DBEst	3596070
AI140291	DBEst	3647748
AI141130	DBEst	3648587

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AI141847	DBEst	3649304
AI143899	DBEst	3665708
AI144100	DBEst	3665909
AI148251	DBEst	3675933
AI149429	DBEst	3677898
AI149592	DBEst	3678061
AI186028	DBEst	3736666
AI186042	DBEst	3736680
AI190341	DBEst	3741550
AI192367	DBEst	3743576
AI192629	DBEst	3743838
AI198930	DBEst	3751536
AI216969	DBEst	3789623
AI217003	DBEst	3789657
AI223292	DBEst	3805495
AI241706	DBEst	3837103
AI251743	DBEst	3848272
AI252466	DBEst	3848995
AI253330	DBEst	3850451
AI253335	DBEst	3850456
AI253338	DBEst	3850459
AI253375	DBEst	3850496
AI253379	DBEst	3850500
AI253436	DBEst	3850391
AI262380	DBEst	3870583
AI263674	DBEst	3871877
AI267162	DBEst	3886329
AI267185	DBEst	3886352
AI267209	DBEst	3886376
AI267289	DBEst	3886456
AI267307	DBEst	3886474
AI267321	DBEst	3886488
AI267454	DBEst	3886621
AI267502	DBEst	3886669
AI268293	DBEst	3887460
AI269060	DBEst	3888227
AI269369	DBEst	3888536
AI270183	DBEst	3889350
AI270472	DBEst	3889639
AI271786	DBEst	3890953
AI272827	DBEst	3895095
AI274047	DBEst	3896315
AI276341	DBEst	3898615
AI276839	DBEst	3899113
AI278611	DBEst	3916845
AI280022	DBEst	3918255
AI283548	DBEst	3921781
AI288965	DBEst	3931274
AI290565	DBEst	3933339
AI291683	DBEst	3934457
AI292286	DBEst	3935060
AI298472	DBEst	3958208
AI298941	DBEst	3958595
AI304857	DBEst	3988546
AI308959	DBEst	4003830
AI312552	DBEst	4018157
AI333055	DBEst	4069614
AI333116	DBEst	4069675
AI335249	DBEst	4072176
AI336326	DBEst	4073253
AI345325	DBEst	4082531
AI366549	DBEst	4126238
AI366549	DBEst	4126238

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AI367850	DBEst	4137595
AI375624	DBEst	4175614
AI375624	DBEst	4175614
AI376561	DBEst	4186410
AI399636	DBEst	4242723
AI417384	DBEst	4260888
AI421720	DBEst	4267651
AI424841	DBEst	4270772
AI431507	DBEst	4303669
AI433180	DBEst	4287371
AI434084	DBEst	4293703
AI434401	DBEst	4295922
AI436016	DBEst	4307232
AI436448	DBEst	4281781
AI446503	DBEst	4295666
AI453199	DBEst	4308687
AI459028	DBEst	4311607
AI469237	DBEst	4331327
AI492520	DBEst	4393523
AI492769	DBEst	4393772
AI494344	DBEst	4395347
AI523940	DBEst	4438075
AI524677	DBEst	4438812
AI538682	DBEst	4452817
AI557059	DBEst	4489422
AI561260	DBEst	4511601
AI567988	DBEst	4526440
AI569715	DBEst	4533089
AI581291	DBEst	4565667
AI583211	DBEst	4569108
AI583570	DBEst	4569467
AI589301	DBEst	4598349
AI597938	DBEst	4606986
AI608591	DBEst	4617758
AI608787	DBEst	4617954
AI608968	DBEst	4618135
AI609193	DBEst	4618360
AI609281	DBEst	4618448
AI623804	DBEst	4648735
AI628689	DBEst	4665489
AI636635	DBEst	4687965
AI650837	DBEst	4734816
AI654096	DBEst	4738075
AI660245	DBEst	4763815
AI669253	DBEst	4834027
AI670084	DBEst	4834858
AI674313	DBEst	4874793
AI678152	DBEst	4888334
AI678703	DBEst	4888885
AI679044	DBEst	4889226
AI679321	DBEst	4889503
AI683140	DBEst	4893322
AI683338	DBEst	4893520
AI683793	DBEst	4893975
AI688798	DBEst	4900092
AI692866	DBEst	4970206
AI694087	DBEst	4971427
AI696819	DBEst	4984719
AI697501	DBEst	4985401
AI734922	DBEst	5056446
AI735069	DBEst	5056668
AI739337	DBEst	5101318
AI739377	DBEst	5101358

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AI743595	DBEst	5111883
AI743691	DBEst	5111979
AI750198	DBEst	5128462
AI750909	DBEst	5129173
AI751119	DBEst	5129306
AI751364	DBEst	5129628
AI751565	DBEst	5129829
AI752319	DBEst	5130583
AI752553	DBEst	5130817
AI752929	DBEst	5131193
AI753108	DBEst	5131372
AI753671	DBEst	5131935
AI754437	DBEst	5132701
AI755181	DBEst	5133445
AI758869	DBEst	5152594
AI761927	DBEst	5177594
AI763126	DBEst	5178793
AI791906	DBEst	5339622
AI793120	DBEst	5340836
AI799521	DBEst	5364993
AI804346	DBEst	5369818
AI808109	DBEst	5394597
AI811021	DBEst	5397587
AI811845	DBEst	5398411
AI814139	DBEst	5425354
AI814674	DBEst	5425889
AI815868	DBEst	5431414
AI822030	DBEst	5441109
AI827641	DBEst	5448312
AI859619	DBEst	5513235
AI864580	DBEst	5528687
AI878968	DBEst	5553017
AI879179	DBEst	5553228
AI879367	DBEst	5553416
AI879992	DBEst	5554041
AI888377	DBEst	5593464
AI911704	DBEst	5631559
AI911997	DBEst	5631852
AI912084	DBEst	5631939
AI916284	DBEst	5636229
AI916584	DBEst	5636439
AI923224	DBEst	5659188
AI924096	DBEst	5660060
AI928185	DBEst	5664149
AI929819	DBEst	5665783
AI936748	DBEst	5675618
AI950087	DBEst	5742397
AI955808	DBEst	5748118
AJ001258	GenBank	2769648
AJ002030	GenBank	2570006
AJ006026	GenBank	3127893
AJ011001	GenBank	4456466
AJ011915	GenBank	3757675
AJ012499	GenBank	5441359
AJ223183	GenBank	3925598
AL035802	DBEst	5927582
AL035987	DBEst	5405617
AL036801	DBEst	5927917
AL037646	DBEst	5928237
AL038985	DBEst	5408101
AL039150	DBEst	5408232
AL041780	DBEst	5421127
AL044019	DBEst	5432247

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
AL046804	DBEst	5434866
AL049055	DBEst	4728364
AL049227	GenBank	4499957
AL049229	GenBank	4499961
AL049296	GenBank	4500057
AL049464	GenBank	4500256
AL049953	GenBank	4884201
AL049954	GenBank	4884203
AL049955	GenBank	4884205
AL049959	GenBank	4884211
AL049987	GenBank	4884238
AL049999	GenBank	4884252
AL050011	GenBank	4884080
AL050089	GenBank	4884107
AL050141	GenBank	4884352
AL050171	GenBank	4884383
AL050187	GenBank	4884402
AL050198	GenBank	4884436
AL050217	GenBank	4884458
AL050392	GenBank	4914613
AL080062	GenBank	5262466
AL080186	GenBank	5262664
AL080235	GenBank	5262728
AL096857	GenBank	5541862
AL096858	GenBank	5541864
AL110197	GenBank	5817115
AL110235	GenBank	5817176
AL117237	GenBank	5834563
AL117499	GenBank	5912003
AL117534	GenBank	5912062
AL118999	DBEst	5924898
AL119085	DBEst	5924984
AL119157	DBEst	5925056
AW020479	DBEst	5874009
AW044114	DBEst	5904643
AW102841	DBEst	6073454
C02094	DBEst	1434324
C16886	DBEst	1571593
C18886	DBEst	1580488
D00017	GenBank	219909
D00022	GenBank	219653
D00068	GenBank	220080
D00099	GenBank	219941
D00422	GenBank	220063
D10495	GenBank	520586
D13119	GenBank	285909
D13287	GenBank	496370
D13665	GenBank	393318
D13866	GenBank	433410
D14662	GenBank	285948
D14697	GenBank	285964
D14710	GenBank	559324
D14812	GenBank	285968
D15049	GenBank	475003
D16431	GenBank	598955
D16937	GenBank	598856
D17188	GenBank	598702
D17268	GenBank	598899
D17409	GenBank	2335046
D17793	GenBank	457407
D21063	GenBank	434752
D23660	GenBank	432358
D25542	GenBank	662389



TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
D28759	GenBank	633074
D29677	GenBank	473938
D31767	GenBank	505091
D31784	GenBank	974184
D31883	GenBank	505093
D31890	GenBank	505107
D37991	GenBank	1019367
D38491	GenBank	559327
D38583	GenBank	560790
D43948	GenBank	603950
D43950	GenBank	603954
D45248	GenBank	1008914
D45887	GenBank	665587
D45915	GenBank	1483130
D49489	GenBank	1136742
D49547	GenBank	710654
D50310	GenBank	1183161
D50371	GenBank	2605591
D55192	DBEst	957089
D55649	GenBank	1132478
D56120	DBEst	970603
D59253	GenBank	1060898
D78586	GenBank	1228048
D79826	DBEst	1180177
D79983	GenBank	1136383
D79986	GenBank	1136389
D79997	GenBank	1136409
D80006	GenBank	1136427
D80012	GenBank	1136437
D80087	DBEst	1177964
D80253	DBEst	1178130
D81635	DBEst	1179512
D82128	DBEst	1183520
D82348	GenBank	1311461
D83197	GenBank	3893154
D83327	GenBank	2687860
D83784	GenBank	1663695
D86227	GenBank	2081619
D87437	GenBank	1665768
D87442	GenBank	1665772
D87470	GenBank	1665822
D87666	GenBank	1620016
D87667	GenBank	1620019
D87682	GenBank	1663699
D87735	GenBank	1620021
D87969	GenBank	1694636
D89052	GenBank	1694672
D90226	GenBank	219946
D90373	GenBank	219477
E00882	GenBank	2169143
E01650	GenBank	2169903
E01797	GenBank	2170049
E01813	GenBank	2170065
E01827	GenBank	2170079
E01979	GenBank	2170227
E02628	GenBank	2170856
E02651	GenBank	2170879
E03569	GenBank	2171785
E06721	GenBank	2174903
E07218	GenBank	2175359
F28779	DBEst	4814405
F30276	DBEst	4815902
F31082	DBEst	4816708

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
H03854	DBEst	866787
H05412	DBEst	868964
H08994	DBEst	873816
H13339	DBEst	878159
H16426	DBEst	881246
H39960	DBEst	916012
H48742	DBEst	988582
H59372	DBEst	1012204
H60722	DBEst	1013554
H69238	DBEst	1030614
H72481	DBEst	1044297
H75695	DBEst	1049638
H78517	DBEst	1056606
H79084	DBEst	1057173
H84729	DBEst	1063923
H85709	DBEst	1067288
H89654	DBEst	1080084
J00269	GenBank	186699
J02621	GenBank	184229
J03005	GenBank	183183
J03040	GenBank	338312
J03171	GenBank	184645
J03191	GenBank	190385
J03210	GenBank	180670
J03464	GenBank	179595
J03473	GenBank	337423
J03799	GenBank	186840
J04080	GenBank	179645
J04164	GenBank	177801
J04177	GenBank	179729
J04765	GenBank	189404
J05013	GenBank	182417
J05021	GenBank	340216
J05192	GenBank	178026
J05633	GenBank	186504
K00558	GenBank	340020
K01566	GenBank	187721
K02765	GenBank	179664
L00160	GenBank	189904
L02547	GenBank	180598
L05092	GenBank	388031
L05186	GenBank	182394
L07633	GenBank	186512
L11066	GenBank	307322
L11932	GenBank	307423
L12711	GenBank	388890
L13848	GenBank	307382
L14599	GenBank	348238
L19161	GenBank	306899
L19184	GenBank	440305
L19597	GenBank	306467
L20941	GenBank	507251
L23959	GenBank	414316
L26081	GenBank	799328
L27560	GenBank	452059
L28010	GenBank	452047
L28809	GenBank	454151
L33404	GenBank	521214
L33930	GenBank	500848
L34155	GenBank	551596
L34839	GenBank	1220373
L38486	GenBank	790816
L42024	GenBank	804748

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
L43575	GenBank	899064
L44349	DBEst	1048859
L54057	GenBank	1196416
M10036	GenBank	339840
M10119	GenBank	182517
M10905	GenBank	182696
M11146	GenBank	182504
M13573	GenBank	189663
M13955	GenBank	186729
M14083	GenBank	189566
M14483	GenBank	339692
M14630	GenBank	339690
M14631	GenBank	183416
M15182	GenBank	183232
M15800	GenBank	187297
M16247	GenBank	178044
M16553	GenBank	339503
M16660	GenBank	184420
M16937	GenBank	184300
M17597	GenBank	340057
M17885	GenBank	190231
M20372	GenBank	189372
M22146	GenBank	337929
M22382	GenBank	190126
M22590	GenBank	179418
M22918	GenBank	189019
M22920	GenBank	189021
M23613	GenBank	189271
M24194	GenBank	187701
M25246	GenBank	340233
M26041	GenBank	188134
M26152	GenBank	1160968
M26325	GenBank	186688
M27913	GenBank	339807
M27971	GenBank	187621
M28373	GenBank	609448
M31159	GenBank	183115
M31212	GenBank	188589
M31899	GenBank	182178
M32110	GenBank	189421
M32790	GenBank	180804
M32798	GenBank	180856
M33308	GenBank	340236
M34064	GenBank	416292
M37583	GenBank	184059
M38106	GenBank	189169
M55409	GenBank	189596
M55542	GenBank	183001
M58485	GenBank	180154
M60457	GenBank	181249
M60854	GenBank	338446
M62403	GenBank	184815
M62810	GenBank	188563
M64241	GenBank	190813
M67468	GenBank	182672
M69181	GenBank	641957
M74002	GenBank	178996
M75126	GenBank	184020
M76729	GenBank	189519
M78113	DBEst	273850
M81757	GenBank	337732
M83248	GenBank	189150
M84739	GenBank	179881

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
M87503	GenBank	184652
M88279	GenBank	186389
M92357	GenBank	306463
N20576	DBEst	1125531
N34255	DBEst	1155397
N35187	DBEst	1156329
N35421	DBEst	1156563
N39717	DBEst	1163262
N40823	DBEst	1164420
N40852	DBEst	1164449
N67927	DBEst	1220052
N76180	DBEst	1238758
N76677	DBEst	1239255
N77080	DBEst	1239658
N84497	DBEst	1260122
N86776	DBEst	1439978
N91638	DBEst	1263947
N92086	DBEst	1264395
N99205	DBEst	1270661
Q37741	N/A	N/A
Q48043	N/A	N/A
Q65676	N/A	N/A
Q90526	N/A	N/A
R06046	DBEst	756666
R17092	DBEst	770702
R47228	DBEst	808115
R55150	DBEst	824379
R55398	DBEst	824693
R68132	DBEst	841649
R72676	DBEst	846708
R73306	DBEst	847338
R78333	DBEst	853443
R92367	DBEst	959907
R93637	DBEst	967803
R99649	DBEst	986250
S41458	GenBank	252252
S42303	GenBank	253482
S54005	GenBank	264772
S66431	GenBank	435777
S70154	GenBank	546900
S70290	GenBank	546602
S79895	GenBank	1195555
S82076	GenBank	1488423
T02792	DBEst	319308
T24119	DBEst	523315
T49314	DBEst	651174
T53479	DBEst	655339
T58797	DBEst	660634
T64560	DBEst	673605
T66112	DBEst	675157
T92160	DBEst	724073
T92396	DBEst	724309
U00947	GenBank	405049
U04815	GenBank	507157
U07151	GenBank	460624
U07857	GenBank	469048
U08470	GenBank	478884
U10323	GenBank	532312
U10439	GenBank	577169
U12465	GenBank	562073
U13665	GenBank	606922
U13877	GenBank	606943
U14550	GenBank	565079

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
U14966	GenBank	550012
U15008	GenBank	600747
U16306	GenBank	608514
U17104	GenBank	609307
U17496	GenBank	596139
U19769	GenBank	924600
U20896	GenBank	1046220
U22431	GenBank	881345
U22815	GenBank	930340
U24105	GenBank	1638873
U24153	GenBank	780807
U27768	GenBank	1216372
U33760	GenBank	995823
U33833	GenBank	1517815
U34877	GenBank	1143231
U39361	GenBank	1066081
U41515	GenBank	1209723
U46570	GenBank	1688073
U50733	GenBank	1255187
U51586	GenBank	1809247
U56255	GenBank	1399688
U59305	GenBank	1695872
U60975	GenBank	5030423
U61083	GenBank	4097430
U61397	GenBank	1518693
U63846	GenBank	1480921
U67784	GenBank	1617516
U68723	GenBank	2114391
U68727	GenBank	2052384
U68758	GenBank	4097815
U70735	GenBank	2360944
U77085	GenBank	1684789
U79258	GenBank	1710211
U79274	GenBank	1710240
U79278	GenBank	1710247
U80213	GenBank	1857418
U81234	GenBank	4098960
U82130	GenBank	1772663
U86602	GenBank	1835785
U87309	GenBank	1842092
U90028	GenBank	2745975
U90441	GenBank	2439984
U90902	GenBank	1913880
U90917	GenBank	1913898
U94831	GenBank	2276459
V00478	GenBank	28244
V00503	GenBank	30123
V05728	N/A	N/A
V11636	N/A	N/A
V57903	N/A	N/A
V59662	N/A	N/A
V59746	N/A	N/A
V84428	N/A	N/A
V86232	N/A	N/A
V87930	N/A	N/A
W07215	DBEst	1281217
W19127	DBEst	1294870
W19407	DBEst	1295308
W19441	DBEst	1295361
W25547	DBEst	1303421
W26197	DBEst	1306608
W38952	DBEst	1320872
W56388	DBEst	1358278

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
W68015	DBEst	1376884
W73140	DBEst	1383275
W73168	DBEst	1383322
W76204	DBEst	1386429
W87522	DBEst	1401728
W87891	DBEst	1401976
X00351	GenBank	28251
X00497	GenBank	32130
X01742	GenBank	35324
X01924	N/A	N/A
X03084	GenBank	29537
X04098	GenBank	28338
X04408	GenBank	31914
X04470	GenBank	28638
X05276	GenBank	37201
X05908	GenBank	34387
X06700	GenBank	30053
X07819	GenBank	35798
X13425	GenBank	31590
X14420	GenBank	30057
X15729	GenBank	38317
X15880	GenBank	30029
X16869	GenBank	31091
X17206	GenBank	34391
X24068	N/A	N/A
X37385	N/A	N/A
X37509	N/A	N/A
X40178	N/A	N/A
X51466	GenBank	31105
X53505	GenBank	36145
X54304	GenBank	34755
X54941	GenBank	29976
X55110	GenBank	35086
X55885	GenBank	34030
X56932	GenBank	23690
X56998	GenBank	37564
X56999	GenBank	37568
X57766	GenBank	456256
X62744	GenBank	36062
X63432	GenBank	28335
X66360	GenBank	36616
X67698	GenBank	37476
X68277	GenBank	29980
X68880	GenBank	31141
X69398	GenBank	396175
X69838	GenBank	287864
X70340	GenBank	37089
X71087	GenBank	288396
X73608	GenBank	793844
X73902	GenBank	452754
X74039	GenBank	456192
X74801	GenBank	671526
X74979	GenBank	400462
X76013	GenBank	531595
X76180	GenBank	452649
X78627	GenBank	607129
X81109	GenBank	535057
X82676	GenBank	3929753
X84939	GenBank	695548
X85373	GenBank	806565
X93036	GenBank	1085025
X93207	GenBank	2462486
X94323	GenBank	1213612

TABLE 1-1

<u>ACC NUM</u>	<u>DATABASE</u>	<u>GI NBR</u>
X94754	GenBank	1702931
X97324	GenBank	1806039
X99920	GenBank	1694827
Y00503	GenBank	34038
Y00757	GenBank	23910
Y00815	GenBank	34266
Y09188	GenBank	2230868
Y11435	GenBank	2910996
Y12065	GenBank	2230877
Y13247	GenBank	2117158
Y13286	GenBank	2853173
Y15286	GenBank	2584788
Y17114	GenBank	4160551
Z18538	GenBank	28711
Z18954	GenBank	396706
Z19054	GenBank	38519
Z21507	GenBank	38521
Z26317	GenBank	416177
Z29093	GenBank	732799
Z31696	GenBank	479156
Z32564	GenBank	473235
Z36531	GenBank	535184
Z37986	GenBank	780262
Z46629	GenBank	758102
Z47087	GenBank	860989
Z74615	GenBank	1418927

Table 2A

#	IMAGE_ID	Gen Bank Accession Number	Ave-Normal- expression	Max- expression-of	Count-up tumors	Count-up lines	Chromosome	Location	Tissue 1	Tissue 2	Tissue 3
26	138917	R62862	3.40	104.38	5.00	0.00	1	119.16	Heart	Placenta	Testis
35	141562	R73003	284.38	1936.47	2.00	0.00			Ear	Placenta	Breast
36	245330	N54596	17.81	1247.32	5.00	0.00	11	16.42			Aorta
38	665656	AA186638	5.96	41.56	1.00	0.00	5	640.93	Parathyroid	Placenta	Foreskin
43	41659	R52852	172.84	1148.01	1.00	1.00	8	370.3	Lymph	Heart	Bone
44	307553	N95249	29.89	401.10	1.00	0.00	12	93.93	Neural	Tonsil	Pooled
46	770192	AA434102	8.17	94.05	2.00	0.00			Ovary	Blood	Lymph
54	82871	T89346	2.44	13.61	1.00	0.00			Small intestine	Liver	Testis
72	123730	R01281	4.04	21.40	2.00	0.00	17	339.25	Spleen	Germ Cell	Pool
80	206994	R99436	0.87	4.57	1.00	0.00	1	165.59	Blood	Parathyroid	Pool
90	210575	H55066	12.82	75.35	3.00	0.00			Brain	CNS	
94	753184	AA400739	19.83	101.81	1.00	0.00			Ear	Thyroid	Stomach
96	299679	W05826	16.34	170.18	2.00	0.00			Foreskin	Whole embryo	Cervix
97	124597	R02373	17.70	142.15	3.00	0.00	3	697.77	Liver	Pool	Kidney
98	78628	T62481	10.11	128.38	6.00	0.00	2	473.05	Lymph	Synovial mem	Thyroid
105	366834	AA028418	1.91	18.11	5.00	0.00			Pancreas	Ovary	Breast
115	785595	AA449459	32.04	420.31	5.00	0.00	4	420.51	Whole embryo	Pool	Heart
116	245330	N54596	18.00	2255.66	5.00	0.00	11	16.42			
120	193913	R83836	12.68	64.55	5.00	0.00	8	316.17	Head and nec	Thyroid	Adipose
121	66532	T67005	5.59	37.42	2.00	0.00	20	336.98	Neural	Pancreas	Heart
125	183337	H42879	13.99	179.31	5.00	0.00	6	118.59	Lymph	Lung	Breast
142	789389	AA464856	55.22	1034.25	2.00	0.00	6	78.3	Thyroid	Bone	Ovary
148	810787	AA481768	52.79	372.83	1.00	0.00	19	80.38	Omentum	Nose	Eye
150	214441	H73550	3.21	285.97	7.00	0.00	14	278.45	Lymph node	Adipose	Small intestine
157	727251	AA412053	100.19	1801.71	6.00	0.00	12	39.87	Esophagus	Synovial mem	Parathyroid
159	564603	AA129552	18.52	118.34	1.00	1.00	12	17.99	Larynx	Thyroid	CNS
161	813678	AA453823	6.53	92.09	3.00	0.00	5	124.15	Eye	CNS	Brain
170	82991	T70503	9.44	51.73	5.48	0.00	6	528.66	Liver	Pooled	Whole embryo
176	785845	AA449118	11.47	64.12	5.59	0.00	10	335.81	Whole embryo	Pooled	Umbilical cord
177	241412	H81220	7.55	44.36	5.88	1.00	13	130.71	Small intestine	Tonsil	CNS
181	125134	R05416	1.08	5.80	5.39	1.00	3	395.72	Smooth musc	Germ Cell	Blood
189	813552	AA455448	25.65	337.42	13.15	3.00	7	427.43	Synovial mem	Stomach	Kidney
193	363377	AA019591	11.01	64.88	5.69	1.00	17	372.99	Placenta	Lung	CNS
201	300412	W20275	132.74	960.13	7.23	1.00	5	572.03	Aorta	Heart	Brain
202	66333	T66828	8.98	35.59	5.10	0.00	22	78.48	Pool	LID not found	Uterus
203	245006	N76278	22.03	125.65	5.70	0.00	1	70.87	Stomach	Prostate	Other
207	121994	T98244	8.31	44.02	5.29	1.00	4	0	Placenta	Tonsil	Parathyroid
215	122019	T98320	10.70	97.20	5.35	0.00					
217	241392	H91281	33.82	209.27	6.17	2.00	5	322.42	Tonsil	Pool	Testis
219	207858	H60588	5.62	35.13	6.25	1.00	6	457.79	Thymus	Blood	Other
222	185081	R91137	9.71	48.79	5.02	1.00			Adrenal gland	Pancreas	Placenta
223	240937	H00997	8.74	60.87	6.96	2.00	5	322.42	Tonsil	Pool	Testis
227	127185	R08220	41.66	295.18	7.09	0.00	6	650.09	Pool	LID not found	LID not found
228	132358	R27329	7.27	48.63	6.69	1.00			Placenta	LID not found	Other
231	357870	W92594	25.09	125.47	5.00	0.00	6	457.79	Thymus	Blood	Ear
232	136919	R38459	6.31	35.24	5.58	1.00			Breast	Ovary	Placenta
234	66377	T66907	541.40	3318.74	6.13	2.00			Pool	LID not found	Other
236	136737	R63623	8.74	100.61	11.51	1.00	12	277.88	Smooth musc	Aorta	Muscle
239	296444	W01048	8.03	102.30	16.95	8.00	1	550.58	Stomach	Germ Cell	Colon
243	120318	T97119	6.40	46.50	7.28	2.00	5	121.96	Pool	LID not found	Other
244	131130	R28397	10.74	60.32	5.62	2.00	6	475.49	Head and nec	Pooled	Pancreas
246	144924	R78513	21.71	285.65	12.24	2.00					



Table 2A

247	122126	T98484	28.60	185.68	6.98	5.00	3.00	1	75.41	Pool	LID not found	Other
248	136933	R39745	9.01	48.82	5.40	1.00	0.00			Placenta	Pool	LID not found
249	186918	H43317	1.43	22.38	15.68	4.00	0.00			Pooled	Breast	Ovary
254	132569	R26798	72.02	500.58	6.95	5.00	0.00	14	253.8	Pooled	Whole embryo	Placenta
255	122170	T98511	1.05	78.19	7.08	2.00	0.00	7	183.65	Pool	LID not found	Other
256	137647	R37684	2.90	14.73	5.07	1.00	0.00			Pooled	Placenta	LID not found
259	126234	R06362	3.19	25.58	8.01	3.00	0.00	8	374.98	Tonsil	Pool	Ovary
261	128833	R10185	7.10	38.50	5.43	2.00	0.00	12	14.49	Pancreas	Pool	LID not found
263	123065	T98529	25.64	237.06	9.25	6.00	0.00	21	225.9		Pool	LID not found
265	68753	T67652	25.31	165.76	6.55	4.00	3.00			Aorta	Pool	LID not found
266	66550	T67022	15.68	126.20	8.05	4.00	0.00	11	54.46	Breast	CNS	Pool
268	132323	R25464	180.89	1094.00	6.05	3.00	1.00	10	475.18	Placenta	LID not found	Other
275	127173	R08153	49.58	348.37	7.03	2.00	0.00	19	208.56		Thyroid	Muscle
278	416833	V06653	18.24	101.94	5.30	1.00	0.00	6	127.46	Gall bladder	LID not found	Other
281	122161	T98458	179.55	1082.64	6.03	2.00	0.00			Pool	Thyroid	Muscle
283	246074	N78944	15.36	91.97	5.99	0.00	1.00	20	212.02	Pancreas	Muscle	Eye
287	294255	N98839	78.47	702.83	8.96	6.00	2.00	6	539.64	Kidney	Pool	LID not found
289	200031	R87154	3.67	24.17	6.39	1.00	0.00			Parathyroid	Pool	LID not found
290	245489	N72510	3.23	32.05	8.93	2.00	0.00			Smooth muscle	Testis	Colon
292	795856	AA461521	3.12	41.15	13.17	5.00	0.00	1	707.84	Ear	Ear	Cervix
298	430153	AA010158	12.49	114.03	8.13	1.00	1.00	22	136.36	Pool	LID not found	Other
297	294444	V01484	3.19	16.00	5.02	1.00	0.00			Pool	Lung	LID not found
302	296741	N74055	19.71	119.09	6.04	0.00	2.00	11	240.08		LID not found	Other
303	110503	T89986	22.10	307.62	13.92	0.00	2.00			Pool	LID not found	Other
305	292986	N63753	7.54	63.74	8.45	0.00	3.00			Pool	LID not found	Other
309	201030	H48360	25.61	135.80	5.30	0.00	1.00	1	174.05	Pool	LID not found	Other
311	141768	R70462	13.88	122.97	8.81	2.00	1.00	17	307.47		LID not found	Other
316	417251	V07752	8.69	105.10	12.09	0.00	1.00	1	590.63	Testis	Prostate	
317	286149	N74360	93.73	682.47	7.28	2.00	3.00			Pool	LID not found	Other
321	198837	R92865	38.26	306.84	8.02	6.00	0.00			Pool	LID not found	Other
333	247194	N54038	3.08	18.86	6.12	2.00	0.00			Pool	LID not found	Other
339	111054	T81574	17.34	131.78	7.60	4.00	1.00			Pool	LID not found	Other
342	298141	V02824	2.08	14.68	7.07	1.00	4.00	7	148.31	Pool	LID not found	Other
343	148373	H04382	1.89	154.54	81.80	1.00	0.00			Pool	Whole embryo	
345	198222	R92962	37.51	298.59	7.96	4.00	0.00	7	547.88	Pooled	Pool	LID not found
348	238681	H79568	3.55	40.18	11.33	2.00	0.00	12	222.35	Heart	Blood	Breast
350	298168	N74365	17.45	132.42	7.59	5.00	0.00	1	85.25	Pool	LID not found	Other
355	111264	T84084	3.19	22.61	7.08	3.00	0.00	11	67.01	CNS	Foreskin	Pool
357	232227	H54822	10.15	59.26	5.84	1.00	0.00	12	248.58	Pool	LID not found	Other
362	244350	N75735	6.71	53.36	7.95	2.00	0.00			Ovary	Pool	LID not found
363	126438	R05642	38.58	354.88	8.20	1.00	2.00			Brain	LID not found	Other
369	196636	R93007	18.40	376.08	20.44	6.00	5.00	16	361.71	Pool	LID not found	Other
371	111634	T60991	40.82	266.49	6.53	1.00	1.00	15	227.19	Pool	LID not found	Other
373	194307	H50747	9.12	116.99	12.83	2.00	5.00	19	184.27	Pool	LID not found	Other
378	201517	R97031	27.36	216.95	7.93	5.00	2.00	5	540.55	Neural	Kidney	Lung
381	109123	T60978	5.98	60.96	10.19	9.00	0.00			Ovary	Lung	Adrenal gland
384	262053	H98812	68.75	374.82	5.62	0.00	1.00			Synovial mem	Thyroid	Parathyroid
392	770868	AA434487	6.01	32.14	5.35	1.00	0.00	12	247.33	Smooth muscle	Thymus	Parathyroid
403	211206	H67988	19.59	103.50	5.28	1.00	0.00	8	151.95	Spleen	Thymus	Neural
412	233457	H80129	5.31	36.08	6.80	1.00	0.00	X	317.31	Eye	Parathyroid	Muscle
414	50182	H17882	5.08	78.05	15.44	2.00	0.00			Esophagus	Uterus	Whole embryo
419	234907	H73080	48.24	723.83	15.00	2.00	0.00	4	490.87		Stomach	Larynx
423	123790	R01428	12.73	106.43	8.38	1.00	4.00			Esophagus	Colon	Adrenal gland
427	753313	AA410265	7.78	86.78	11.15	3.00	0.00	1	98.87		Colon	
446	823696	AA489743	20.09	441.68	21.98	2.00	0.00	10	421.81			

Table 2A

450	510032	AA053051	9.68	52.00	5.37	1.00	0.00	17	126.91	Adipose	Colon	Kidney
452	724112	AA411244	5.50	43.21	7.86	9.00	0.00	2	578.49	Ovary	Testis	Pool
460	153411	R47878	103.88	1436.31	13.83	6.00	0.00	6	117.99	Lymph node	Lymph	Small intestine
472	144834	R77251	8.95	48.81	5.57	1.00	0.00	6	544.54	CNS	Thyroid	Cervix
475	753682	AA086601	51.39	322.04	6.27	1.00	0.00	10	510.68	Nose	Pooled	Stomach
485	840364	AA485626	76.01	443.38	5.83	1.00	0.00	X	87.96	Ear	Done	Foreskin
490	713690	AA284328	4.15	52.21	12.57	4.00	0.00	5	353.54	Bone	Ear	Adrenal gland
498	51463	H21071	74.83	488.06	6.50	2.00	0.00	15	340.51	Colon	Pool	LID not found
604	198694	R95132	125.76	1236.62	9.85	1.00	0.00	14	271.74	Tonsil	Tonsil	Blood
513	824074	AA481227	12.79	94.86	7.41	3.00	0.00	11	264.42	Pooled	Placenta	Thymus
520	144740	R76228	4.18	40.44	9.87	2.00	0.00	2	578.48	Pooled	CNS	Foreskin
527	342593	W66537	4.62	76.81	18.61	5.00	0.00	11	40.11	Pooled	Breast	Parathyroid
528	66336	T66833	7.00	35.19	5.03	1.00	0.00	10	42.59	Nose	Aorta	Ear
538	358885	W84714	6.25	125.76	20.11	1.00	0.00	3	459.05			
557	32684	R43544	691.72	4078.06	5.89	1.00	0.00	4	644.84	Gall bladder	Eye	Bone
564	838359	AA457178	198.25	1429.49	7.21	5.00	2.00	14	15.32	Synovial mem	Nose	Skin
565	897635	AA598582	92.09	537.51	5.84	1.00	0.00	7	510.47	Peripheral ner	Nose	Pancreas
571	626502	AA188155	37.77	303.69	8.04	3.00	0.00	9	410.77	Breast	Tonsil	Aorta
575	814460	AA459247	5.28	28.23	5.36	1.00	0.00	20	107.35	Heart	CNS	Brain
581	290039	N76471	4.86	41.62	8.56	4.00	0.00	10	360.96		CNS	Germ Cell
583	241003	H81010	13.87	336.66	24.28	2.00	0.00	10		Adrenal gland	Pool	LID not found
584	194704	R89904	139.34	975.12	7.00	4.00	0.00	14	217.02	Thyroid	LID not found	Other
587	126237	R06370	3.59	30.07	8.37	2.00	0.00	21	146.87	Placenta	LID not found	Other
592	194395	R83180	4.86	43.50	8.77	2.00	0.00	12	22.39	Lymph	Kidney	Testis
594	139766	R62242	4.03	23.97	5.95	1.00	0.00	X	26.02	Pool	LID not found	Other
596	141106	R66219	30.01	315.07	10.50	6.00	5.00	6	18.54	Larynx	Thyroid	Pancreas
599	282171	N79167	37.25	254.39	6.83	1.00	2.00	13	130.76	Pool	LID not found	Other
600	243638	N49883	21.84	147.82	8.83	3.00	0.00	10	529.82		LID not found	Other
605	154462	R54664	19.87	113.40	5.77	1.00	0.00	3	198.24	Pool	Muscle	Spleen
608	109523	T81972	14.03	280.23	19.98	3.00	0.00	17	278.88	Tonsil	CNS	Pancreas
611	126401	R06568	20.42	148.44	7.27	3.00	0.00	1	111.83	Pooled	Brain	Pool
624	285973	N73551	21.73	297.33	13.69	5.00	0.00	11	45.26	Prostate	LID not found	Other
627	126465	R06668	6.32	33.06	5.23	0.00	3.00	10	540.8		Prostate	Kidney
630	120309	T97215	86.73	548.95	6.33	2.00	0.00	18	185.45	Parathyroid	Foreskin	Thyroid
639	196860	R93087	7.30	38.47	5.27	1.00	0.00	13	135.33			
640	242778	H93603	19.78	104.65	5.29	2.00	0.00	10	107.85	Heart	Placenta	Kidney
646	111136	T83556	48.05	391.53	8.15	4.00	0.00	17	331.17	Epididymis	Breast	Tonsil
648	169841	R96526	18.22	258.47	14.24	6.00	4.00	8	278.83	CNS	Pool	LID not found
649	325102	W49715	55.50	308.70	5.53	1.00	0.00	0	588.46	CNS	Uterus	Eye
652	140921	R66585	11.48	113.13	9.85	3.00	0.00	0	382.37	Pool	LID not found	Other
655	321560	W32884	22.61	144.53	6.39	2.00	0.00	1	623.7	Lymph	Pool	LID not found
658	194600	R84242	16.46	171.30	10.41	3.00	0.00	14	158.61	Placenta	LID not found	Other
660	141230	R66652	5.33	66.65	12.51	5.00	4.00	9	370.33	Ear	Lymph	Germ Cell
662	110507	T82819	83.42	442.60	5.31	1.00	0.00	3	570.8	Pool	Testis	Placenta
664	194587	R84376	56.83	430.85	7.56	5.00	0.00					
665	135673	R31591	36.62	249.43	6.48	4.00	2.00					
667	126509	R06745	4.10	32.22	7.85	2.00	0.00					
673	292522	N91307	28.39	259.36	9.13	1.00	0.00					
677	240050	H82330	105.95	601.66	5.68	1.00	2.00					
688	785382	AA453498	7.42	60.23	8.12	3.00	0.00					
689	207655	H82287	44.34	286.04	8.45	2.00	3.00					
690	136855	R36212	34.89	197.48	5.64	1.00	1.00					
693	282308	N80950	5.08	60.21	11.85	1.00	0.00					
694	212180	H66943	17.47	105.55	6.04	0.00	2.00					
696	810096	AA464987	3.63	28.47	7.85	2.00	0.00					

Table 2A

696	128260	R09873	3.07	17.78	5.80	1.00	0.00	4	451.92	CNS	Germ Cell	Uterus
699	810104	AA464970	1.35	10.20	7.54	2.00	1.00	15	129.27	Peripheral ner	Blood	Spleen
700	810866	AA459869	3.01	18.96	5.84	1.00	0.00	3	473.39	Pooled	Brain	Pool
702	139376	R65573	4.27	28.10	6.11	1.00	0.00	17	534.21	Cervix	Colon	.
703	1030928	AA620346	1.42	7.13	5.02	0.00	1.00			Tonsil	Testis	LID not found
704	810643	AA459859	5.57	39.10	7.02	5.00	0.00			Thyroid	Spleen	Aorta
706	292812	N90470	102.29	1074.94	10.51	2.00	0.00	X	245.06	CNS	Germ Cell	Colon
707	824659	AA491302	5.45	46.58	8.55	6.00	0.00			Pool	LID not found	Other
709	240318	H89795	10.25	78.85	7.50	6.00	0.00			Bone	Muscle	Lymph
710	247835	N58163	39.23	549.46	14.01	8.00	5.00	11	235.26	Synovial mem	Germ Cell	Liver
711	324225	WA7350	18.74	143.45	7.27	3.00	0.00			Ignore	Colon	Blood
719	321271	AA037410	2.82	33.33	11.83	9.00	0.00			Ignore	Colon	Blood
723	773220	AA426655	15.17	98.17	6.54	2.00	0.00	X	238.33	Small intestine	Gall bladder	Parathyroid
729	293847	N86001	12.26	103.51	8.44	2.00	5.00	7	530.17	Pool	LID not found	Other
739	327247	AA264292	1.12	6.46	5.74	1.00	0.00	9	119.58	Eye	Pooled	Placenta
740	810083	AA464962	9.45	105.15	11.12	2.00	0.00			Ignore	Synovial mem	Skin
741	240674	H90946	41.50	206.27	5.01	1.00	1.00	1	745.7	Prostate	Pool	LID not found
746	785788	AA459853	26.14	177.95	6.81	1.00	0.00			Adrenal gland	Parathyroid	Germ Cell
764	364555	AA022601	3.67	39.80	10.84	5.00	0.00			Aorta	CNS	Brain
765	241355	H80355	46.78	457.50	9.78	3.00	2.00	11	67.01	Cell bladder	Adipose	Liver
798	161458	H25546	5.22	244.53	48.84	1.00	1.00			Brain	Germ Cell	Breast
800	178534	H51481	5.80	68.54	11.83	2.00	0.00			Brain	Colon	Pool
806	292613	N68585	31.75	237.99	7.50	1.00	1.00	4	423.35	Ovary	Colon	Pool
816	129146	R10896	80.28	411.33	5.12	0.00	1.00	2	130.74	Thymus	Adipose	Parathyroid
836	593598	AA102670	3.06	23.49	7.69	2.00	0.00	5	631.73	Pancreas	Stomach	Breast
843	754031	AA479981	4.32	28.07	6.74	1.00	1.00	12	43.68	Eye	Brain	CNS
844	209137	H83934	5.84	39.73	8.80	2.00	0.00	X	345.45	Bone marrow	Pooled	Placenta
846	135221	R32652	2.34	78.39	33.57	2.00	0.00			Cervix	Placenta	Colon
851	809784	AA454743	3.78	396.88	104.82	16.00	1.00	19	274.67	Ovary	CNS	Colon
854	66731	T64905	1.81	33.71	16.64	1.00	0.00			Whole embryo	Placenta	Germ Cell
856	734479	AA410567	10.24	488.78	47.75	5.00	0.00	1	181.81	Adipose	Skin	Stomach
871	144777	R78263	48.30	252.69	5.23	0.00	1.00	1	539.01	Synovial mem	Adrenal gland	Gall bladder
872	305608	N80246	2.99	58.07	24.30	3.00	1.00	7	655.1	Parathyroid	Liver	Pancreas
878	73923	T47229	32.41	194.29	6.00	1.00	0.00	10	270.14	Parathyroid	Breast	CNS
893	243816	N39161	6.10	64.47	10.56	1.00	0.00	7	449.88	Adipose	Spleen	Placenta
900	838373	AA458801	16.75	188.15	11.83	0.00	2.00	13	220.79	Thymus	Ear	Cervix
905	770957	AA428170	13.13	72.40	5.51	0.00	1.00			Aorta	CNS	Bone
908	66564	T67056	7.48	39.15	5.24	1.00	0.00			Larynx	Pool	Colon
923	50690	H18438	7.30	36.82	5.05	1.00	0.00			Whole embryo	Brain	Testis
924	753862	AA410517	1172.13	6187.71	5.28	1.00	0.00	6	17.79	Cervix	Placenta	Skin
925	770212	AA434115	10.84	360.36	33.25	10.00	0.00	1	671.26	Adipose	Breast	Ovary
927	841332	AA487634	34.60	370.01	10.89	0.00	2.00	12	47.11	Lymph	Bone marrow	Thymus
929	186992	R93124	40.11	375.69	9.37	0.00	1.00			Liver	Esophagus	Synovial membrane
934	86220	T72422	5.52	31.38	5.68	1.00	0.00	4	644.55	Gall bladder	Liver	Pool
940	813530	AA447774	40.44	311.58	7.70	2.00	0.00	8	563.64	Head and neck	Muscle	Ovary
948	40017	R52654	185.04	936.24	5.08	1.00	0.00	7	111.22	Adrenal gland	Thyroid	CNS
954	784876	AA448015	3.03	23.48	7.74	0.00	2.00	10	471.64	CNS	Whole embryo	Testis
956	83231	T68351	3.89	18.90	5.12	0.00	1.00			Liver	Breast	Colon
981	111204	T64382	50.91	402.17	7.90	2.00	0.00	9	357.99	Eye	Pool	LID not found
982	66562	T67093	126.61	695.78	5.50	1.00	0.00	14	217.02	Thyroid	Breast	Pool
984	132848	R25641	5.08	28.44	5.79	2.00	0.00	14	7.79	Placenta	Heart	LID not found
989	244205	N52860	30.49	200.44	6.57	0.00	3.00			Eye	Colon	Whole embryo
973	711857	AA281189	111.11	671.93	6.05	1.00	0.00	8	186.21	.	Lung	Breast
974	88400	T66930	88.51	616.50	6.98	3.00	0.00	19	291.14	Breast	Pool	LID not found
985	292416	N91188	14.95	83.13	5.58	1.00	0.00	1	32.1			

Table 2A

989	195458	R01710	27.12	175.16	6.46	4.00	0.00	4	480.05	Prostate	Eye	Fore skin
991	122384	T99145	47.05	333.08	7.08	1.00	2.00	17	372.88	Placenta	LID not found	Other
996	132623	R26813	104.02	678.64	6.52	3.00	0.00	22	71.14	Pool	LID not found	Other
998	234133	N99789	88.44	641.84	6.52	2.00	0.00	17	140.84	Thyroid	Spleen	Adrenal gland
999	122359	T99150	5.35	32.07	6.00	1.00	0.00	8	123.4	Testis	Pool	LID not found
1005	186148	R92352	182.88	1375.47	7.52	3.00	0.00	7	521.82	Pool	LID not found	Other
1007	122702	T88972	4.36	41.19	8.44	1.00	1.00	6	150.81	Pool	LID not found	Other
1011	233347	H77533	2.82	15.59	5.52	2.00	0.00	13	157.6	Placenta	Pool	LID not found
1013	201207	R99288	37.43	282.46	7.55	3.00	0.00	3	197.02	Larynx	Brain	LID not found
1015	122694	T69011	3.54	24.49	6.92	1.00	0.00	19	82.52	Pool	Adipose	Skin
1016	154312	R53024	3.94	33.50	8.50	0.00	1.00	19	430.49	Germ Cell	Tail	Breast
1017	123354	T69817	5.79	40.72	7.03	1.00	0.00	X	245.08	Placenta	LID not found	Other
1020	133225	R26931	10.47	76.10	7.27	1.00	0.00	19	82.52	Blood	Lung	Placenta
1021	244329	N75729	35.08	219.10	6.25	3.00	0.00	4	40.26	Umbilical cord	Uterus	Placenta
1023	212034	H69048	59.04	437.88	7.42	4.00	2.00	3	61.75	Pool	LID not found	Other
1024	136188	R53860	18.37	169.86	9.25	8.00	0.00	1	111.69	Breast	Cervix	Blood
1029	110930	T90360	43.18	325.87	7.55	3.00	0.00	17	429.02	Pool	LID not found	Other
1033	293437	R92085	6.95	41.05	5.90	0.00	3.00	12	47.11	Small intestine	Smooth muscle	Nose
1036	133333	R26855	32.01	214.24	6.69	4.00	1.00	19	82.52	Blood	Lung	Placenta
1038	50214	H16746	1.78	13.27	7.46	1.00	0.00	4	40.26	Umbilical cord	Uterus	Placenta
1040	136165	R53900	61.67	406.70	6.64	2.00	0.00	3	61.75	Pool	LID not found	Other
1041	127076	R07998	17.84	209.29	11.73	5.00	5.00	1	111.69	Breast	Cervix	Blood
1045	127243	R08297	6.75	43.43	8.25	5.00	0.00	17	429.02	Pool	LID not found	Other
1048	136189	R53910	24.19	195.56	8.25	5.00	0.00	12	47.11	Small intestine	Smooth muscle	Nose
1049	126638	R08938	6.16	48.63	8.05	3.00	0.00	19	48.71	Blood	Colon	Testis
1051	123720	R01277	38.35	234.45	6.11	3.00	0.00	4	40.26	Umbilical cord	Uterus	Placenta
1054	306808	W24055	4.33	40.15	9.27	3.00	1.00	3	61.75	Pool	LID not found	Other
1057	245413	N77203	27.69	626.95	11.28	6.00	5.00	1	111.69	Breast	Cervix	Blood
1060	251351	H96213	5.86	66.15	11.28	1.00	1.00	17	429.02	Pool	LID not found	Other
1062	296199	W02639	18.99	112.20	6.81	4.00	3.00	12	47.11	Small intestine	Smooth muscle	Nose
1063	60394	T65770	1.86	20.07	10.25	2.00	0.00	19	48.71	Blood	Colon	Testis
1065	197051	R63153	57.86	500.64	8.66	4.00	0.00	4	40.26	Umbilical cord	Uterus	Placenta
1067	124271	R02036	1.63	10.32	6.34	1.00	0.00	3	61.75	Pool	LID not found	Other
1068	291459	N72852	7.30	51.41	7.05	3.00	1.00	1	111.69	Breast	Cervix	Blood
1069	242644	H94978	2.87	18.35	8.40	2.00	0.00	17	429.02	Pool	LID not found	Other
1070	296330	W03050	17.71	94.24	5.32	1.00	1.00	12	473.2	Pool	LID not found	Other
1073	197093	R93412	12.23	84.11	6.88	3.00	0.00	4	669.93	Pool	LID not found	Other
1074	245556	N77223	55.16	407.65	7.39	1.00	2.00	14	135.79	Pool	LID not found	Other
1075	128795	R16769	12.74	178.01	13.98	2.00	2.00	12	473.2	Pool	LID not found	Other
1076	296334	W03052	23.32	129.45	5.55	1.00	0.00	4	669.93	Pool	LID not found	Other
1081	230637	H75490	18.86	175.55	9.26	6.00	0.00	14	135.79	Pool	LID not found	Other
1083	295527	N74942	60.37	512.31	8.48	3.00	3.00	12	473.2	Pool	LID not found	Other
1106	245836	N76803	2.43	23.11	9.51	1.00	0.00	4	350.76	Pool	LID not found	Other
1108	287843	N62328	3.17	17.94	5.66	1.00	0.00	7	468.83	Aorta	CNS	Brain
1111	187618	R63758	2.89	56.67	18.95	6.00	1.00	1	119.16	Heart	Placenta	Breast
1113	186593	R94683	3.26	16.73	5.13	1.00	0.00	7	551.66	Pool	LID not found	Other
1115	136874	R39705	43.22	287.22	6.65	2.00	2.00	16	482.73	Small intestine	Adipose	Nose
1120	785604	AA460003	1.45	12.57	8.66	1.00	0.00	8	334.17	Pool	LID not found	Other
1121	292569	N91330	80.83	465.59	5.76	1.00	0.00	5	578.78	Pool	LID not found	Other
1126	265321	W04369	18.15	93.56	5.15	1.00	0.00	11	41.44	Lymph	Heart	Breast
1127	186282	R02609	6.44	37.80	5.87	1.00	0.00	4	450.16	Testis	Placenta	Pool
1129	190339	R94212	9.24	113.39	12.27	3.00	5.00	1	119.16	Heart	Placenta	Breast
1130	246144	N55492	24.28	185.77	7.69	5.00	0.00	4	450.16	Testis	Placenta	Pool
1137	186026	R94601	18.08	128.78	7.18	5.00	0.00	1	711.38	Testis	Placenta	Pool
1139	194155	H51056	7.92	43.68	5.54	1.00	0.00					

Table 2A

1140	782217	AA431988	3.50	28.76	8.21	4.00	0.00	1	143.55	Colon	Testis	Prostate
1142	287411	W03672	20.92	140.00	6.69	5.00	4.00	6	458.69	Kidney	Pool	LID not found
1149	202337	H53156	27.14	235.66	8.68	1.00	3.00			Brain	Pool	LID not found
1151	230370	H75531	357.49	1974.00	5.52	1.00	0.00	17	24.51	Adipose	Bone	Blood
1152	236515	H95141	1.83	10.29	5.61	2.00	0.00			Bone	Whole embryo	Pancreas
1164	767183	AA424575	5.80	37.86	6.42	2.00	1.00	3	422.37	Lymph	Tonsil	Blood
1165	266988	N70349	12.50	68.30	5.47	2.00	0.00			Pool	Heart	LID not found
1174	154016	AA48796	123.44	713.75	5.78	2.00	0.00	16	480.32	Cervix	Stomach	Lung
1184	563444	AA112660	224.35	1858.11	8.28	2.00	0.00	19	235.13	Nose	Synovial mem	Thyroid
1191	130153	R21614	28.66	211.73	7.39	1.00	1.00	16	423.94			
1195	411541	R52789	44.80	285.97	6.37	1.00	4.00	X	351.05	Spleen	Lymph	Tonsil
1196	212640	H68620	2.04	22.60	11.07	1.00	0.00			CNS	Prostate	Colon
1202	244189	N52474	10.02	73.50	7.33	4.00	0.00	X	320.23	Ear	Gall bladder	Blood
1203	687990	AA236617	61.49	332.30	5.40	1.00	0.00	2	668.45	Liver	Ovary	Spleen
1206	233721	H79047	115.33	793.08	6.88	1.00	0.00			Ear	Prostate	Eye
1222	613179	AA456321	13.08	211.38	18.18	3.00	0.00	12	397.89			
1227	245970	N52263	107.40	782.42	7.28	2.00	0.00	4	71.55	Gall bladder	Pooled	Fore skin
1228	285728	W02285	15.92	93.15	5.85	1.00	0.00			Gall bladder	Adipose	Aorta
1231	265680	N31467	26.17	145.89	5.56	1.00	1.00			Ignore	Head and nec	Germ Cell
1232	308041	W24429	2.84	18.79	5.71	1.00	0.00			Brain	Ovary	Lung
1235	740027	AA477514	87.23	582.63	6.68	1.00	0.00			Eye	Lymph	Testis
1240	240399	N48137	3.07	17.85	5.82	1.00	0.00	4	618.04	Pool	Lung	LID not found
1249	840788	AA486138	19.56	141.71	7.24	2.00	0.00	16	61.77	Larynx	Synovial mem	Pancreas
1250	363088	AA019482	4.30	82.85	19.25	7.00	0.00	15	145.79	Eye	Pooled	Brain
1261	840483	AA487797	37.99	330.90	8.71	2.00	0.00			Bone marrow	Adipose	Breast
1264	287421	W03677	132.17	704.71	5.33	1.00	0.00	2	607	Adipose	Muscle	Liver
1266	713762	AA292876	12.27	96.61	8.04	3.00	0.00	1	538.34	Stomach	Pancreas	Pooled
1271	828478	AA504351	46.80	284.97	6.08	1.00	0.00	19	216	Stomach	Tonsil	Pooled
1273	815526	AA456878	18.85	104.91	5.56	0.00	1.00	20	236.87	Larynx	Skin	Germ Cell
1279	825577	AA504710	4.57	34.79	7.61	1.00	0.00			Adrenal gland	Placenta	Tonsil
1284	129392	R11236	6.48	43.89	6.78	1.00	0.00	11	253.28	Aorta	Colon	Brain
1289	783729	AA443351	20.42	128.31	6.28	0.00	1.00	17	307.47			
1292	199251	R65780	18.05	136.35	7.55	3.00	0.00	3	512.91	Pool	LID not found	Other
1299	897687	AA588884	14.14	173.00	12.23	1.00	0.00	12	25.02	Peripheral ner	Cervix	Colon
1304	208531	H61878	22.28	292.14	13.11	2.00	0.00	3	157.87	Thyroid	Muscle	
1311	788334	AA453015	6.57	38.03	5.78	1.00	0.00			Germ Cell	Uterus	Cervix
1313	714106	AA284668	102.68	621.94	6.06	0.00	1.00	4	436.69	Larynx	Pancreas	Gall bladder
1320	812266	AA455062	37.95	191.67	5.05	1.00	0.00			Placenta	Uterus	Kidney
1330	263200	H99344	17.89	106.27	5.88	0.00	1.00	3	340.31	Eye	Tonsil	Ovary
1331	782811	AA448261	130.33	687.71	5.28	0.00	2.00	2	87.96		Fore skin	Pool
1335	786675	AA451804	17.56	4927.42	280.58	21.00	1.00	20	255.21	Epididymis	Ovary	Thyroid
1336	296072	W02558	8.04	59.98	7.46	1.00	0.00	3	732.12	Thymus	Tonsil	Breast
1345	244147	N51018	13.92	78.52	5.84	1.00	0.00	X	356.28	Nose	Bone	Placenta
1349	178232	H60229	3.08	23.17	7.52	2.00	0.00			Bone	Germ Cell	Colon
1350	128503	R10159	7.33	43.58	5.95	1.00	0.00			Pool	LID not found	Other
1352	194856	R84407	109.88	704.54	6.41	3.00	0.00	12	395.98	Pool	Brain	LID not found
1356	140301	R86924	17.17	145.62	8.48	6.00	4.00			Stomach	Bone	Placenta
1360	214858	H74032	142.25	964.40	6.78	3.00	0.00			Pool	LID not found	Other
1365	125788	R07694	56.28	358.42	6.37	2.00	0.00	18	299.63	Spleen	Germ Cell	Pancreas
1367	128617	R16479	12.84	66.58	5.30	0.00	1.00	2	76.15	Pool	LID not found	Other
1369	202339	H53038	10.03	56.20	5.80	2.00	0.00	X	245.06	Pool	LID not found	Other
1371	128568	R08682	37.81	219.55	5.79	0.00	1.00			Pool	LID not found	Other
1372	138239	R67891	6.52	80.26	12.32	5.00	4.00	9	252.87	Tonsil	Prostate	Placenta
1375	128627	R16484	33.68	245.63	7.29	4.00	0.00			Lung	Pool	LID not found

Table 2A

1382	200863	R98877	150.20	803.04	5.35	1.00	0.00	15	142.19	Fore skin	Placenta	Stomach
1384	208949	R88736	5.57	31.49	5.66	1.00	0.00	19	73.81	Pool	LID not found	Other
1388	123815	R01451	3.92	19.87	5.07	0.00	1.00			Breast	Pool	LID not found
1393	320630	V01784	5.51	36.06	6.54	1.00	0.00	11	292.28	Parathyroid	LID not found	Other
1398	137760	R68514	74.76	865.11	8.90	3.00	0.00					
1400	242010	H83819	27.74	600.51	21.65	6.00	5.00	19	278.5	Parathyroid	Uterus	Pool
1401	294916	N71457	68.02	407.57	5.99	1.00	2.00	17	268.7	Kidney	Colon	Testis
1402	274932	R84636	3.51	23.16	6.80	1.00	0.00			Placenta	LID not found	Other
1404	137767	R68245	72.54	805.29	8.34	5.00	0.00	5	528.2	Pancreas	Uterus	Germ Cell
1408	139354	R63735	16.98	114.43	6.74	3.00	0.00			Ovary	Kidney	Tonsil
1416	195720	R89157	11.00	64.99	5.91	3.00	2.00	1	712.78	Pool	LID not found	Other
1418	110282	T81988	21.87	176.10	8.05	4.00	0.00	15	288.07	Parathyroid	Aorta	Whole embryo
1419	128984	R07313	4.93	25.54	5.18	1.00	0.00	21	225.9	Pool	LID not found	Other
1422	110987	T90369	111.54	967.64	8.68	5.00	0.00	18	81.8	Pool	LID not found	Other
1424	193481	H47297	19.67	185.11	9.41	4.00	0.00	13	65.51	Placenta	Pool	LID not found
1428	137885	R68381	15.71	182.89	11.64	6.00	5.00	4	637.68	Placenta	Fore skin	Pool
1430	203918	H58655	5.13	32.70	6.38	0.00	1.00	18	442.88	Pool	LID not found	Other
1432	183200	H47335	63.58	650.81	10.24	3.00	4.00	19	71.08	Adipose	Thyroid	Bone
1438	123079	R02528	32.03	199.57	6.23	4.00	0.00	2	513.74	Blood	Testis	Pool
1441	210494	H65481	8.02	39.02	6.46	2.00	0.00	3	222.62	Pool	Prostate	LID not found
1442	129800	R16800	1.87	12.41	6.64	1.00	0.00	7	615.41	Fore skin	Kidney	Pool
1445	214175	H90477	36.55	289.65	7.83	4.00	0.00	5	500.38	Pool	LID not found	Other
1448	210431	H84938	90.85	459.67	5.06	1.00	0.00	13	268.41	Esophagus	Ovary	CNS
1452	322537	W16283	12.82	645.92	42.58	4.00	0.00	20	70.87	Muscle	Prostate	Pool
1453	241497	H50490	38.28	293.84	7.68	3.00	4.00	19	77.18	Placenta	LID not found	Other
1454	292542	H91317	34.54	332.54	9.83	3.00	0.00	1	109.76	Pool	CNS	Stomach
1456	39874	R54560	5.31	116.80	21.81	6.00	0.00	2	631.89	Smooth muscle	CNS	Stomach
1457	210301	H64972	8.10	54.06	8.80	1.00	0.00	1	711.82	Stomach	Germ Cell	Heart
1468	123546	R01566	34.96	232.53	6.65	2.00	0.00	X	245.08	Aorta	CNS	Placenta
1472	770935	AA434382	122.18	1320.42	10.81	3.00	1.00	11	247.68	Tonsil	Pool	LID not found
1473	210525	H65044	4.06	36.83	9.06	1.00	0.00	X	245.08	Bone	Placenta	LID not found
1474	109049	T80942	82.55	417.41	6.67	3.00	0.00	13	239.18	Peripheral ner	Lung	Breast
1479	365865	AA028030	1.38	8.22	5.95	1.00	0.00	12	513.08	Pool	LID not found	Other
1484	257823	N30839	108.28	670.95	6.31	3.00	0.00	9	148.92	Stomach	Testis	CNS
1485	240768	H91353	14.27	98.73	6.92	1.00	0.00	17	416.26	Pool	Brain	LID not found
1495	810864	AA403972	6.22	34.07	5.48	1.00	0.00			Thymus	CNS	Tonsil
1501	241274	H91216	65.83	309.64	5.55	2.00	0.00	20	193	Prostate	Smooth muscle	Thyroid
1502	198258	R85827	8.14	41.66	5.12	1.00	0.00	7	638.71	Whole embryo	Heart	Pool
1503	809603	AA458483	276.70	1508.42	5.44	1.00	0.00	6	829.01	Pool	LID not found	Other
1505	210565	H65839	102.74	591.69	5.76	3.00	1.00	7	682.13	Ear	CNS	Liver
1507	505414	AA156251	1.38	7.07	5.12	1.00	0.00	2	356.2	Umbilical cord	Muscle	Marrow
1508	278399	N48708	4.05	28.75	6.60	1.00	0.00	19	35.68	Smooth muscle	Tonsil	Blood
1510	203510	N68574	40.34	261.52	6.48	3.00	0.00	22	136.3	Synovial mem	Umbilical cord	Lymph node
1526	143756	R78489	8.32	42.26	5.08	1.00	0.00	7	528.17	Umbilical cord	Brain	Adrenal gland
1528	811028	AA465373	69.42	382.14	5.50	1.00	0.00	5	421.53	Liver	Colon	
1531	244323	N54803	122.58	782.85	6.47	4.00	2.00					
1533	246276	N77098	69.29	489.47	7.38	3.00	5.00					
1535	325380	AA284285	3.58	28.41	7.93	2.00	0.00					
1548	247281	N57984	176.76	1010.53	5.72	1.00	0.00					
1552	207288	H58620	8.53	45.76	5.36	1.00	0.00					
1556	322561	W15277	848.75	7077.92	8.34	2.00	0.00					
1564	123400	T99639	19.25	119.66	6.21	2.00	0.00					
1568	365945	AA083831	5.83	32.20	5.52	1.00	0.00					
1587	33271	R25521	12.26	85.06	6.84	1.00	0.00					
1570	78284	T50786	11.05	98.99	6.74	6.00	0.00					

Table 2A

1571	283086	N69872	21.38	111.69	5.22	1.00	0.00	6	98.75	Spleen	Blood	Liver
1590	123117	T96559	18.52	128.10	6.97	1.00	0.00	7	40.71	Smooth muscle	Nose	Forebrain
1600	33478	R44864	1.07	10.86	10.24	1.00	0.00	8	385.82	Germ Cell	Testis	Ovary
1608	223098	H84113	48.32	399.66	8.27	3.00	0.00	11	278.98	Eye	Bone	Breast
1610	124261	R02346	96.33	617.76	8.41	2.00	0.00	19	269.17	Umbilical cord	Nose	Esophagus
1616	770910	AA433951	2.57	284.36	110.73	22.00	3.00	1	671.44	Esophagus	Adipose	Colon
1638	739801	AA477893	281.52	1440.73	5.12	1.00	0.00	13	86.55	Synovial mem	Blood	Colon
1639	739126	AA421687	11.46	95.99	8.35	2.00	0.00			Cervix	Synovial mem	Marrow
1646	897982	AA598863	34.56	427.57	12.37	5.00	0.00	20	123.04	Peripheral nervous system	Adipose	
1652	949938	AA599177	222.88	2198.58	8.86	1.00	0.00	17	317.13	Adrenal gland	Muscle	Cervix
1653	241474	H80415	68.14	405.18	6.00	1.00	1.00	19	78.13	Neural	Shin	Blood
1655	815265	AA461554	24.00	218.20	9.09	1.00	0.00	11	259.81	Lymph node	Pancreas	
1668	841641	AA487700	134.93	754.06	5.59	1.00	0.00	6	118.59	Trachea	Ignore	Aorta
1672	801099	T63324	13.00	200.07	15.38	8.00	0.00	18	78.05	Trachea	Ignore	Aorta
1674	810859	AA458965	63.93	359.71	5.63	1.00	0.00	7	384.36	Gall bladder	Uterus	Whole embryo
1679	814465	AA459286	6.59	34.82	5.29	1.00	0.00	11	252.9	Pancreas	Pooled	Uterus
1687	591907	AA143436	4.84	67.25	13.60	1.00	0.00	11	373.41	Omentum	Testis	Heart
1691	358162	W95346	6.38	56.90	8.92	1.00	0.00	11	471.03	Ovary	Stomach	Ear
1692	839736	AA504943	33.81	320.07	9.47	2.00	0.00	18	217.43	Acorta	Forebrain	Pancreas
1696	782513	AA432030	48.69	416.62	8.56	5.00	0.00	21	317.36	Cervix	Uterus	Brain
1698	815542	AA456886	21.58	508.40	23.61	7.00	0.00	13	387.63	Muscle	Adrenal gland	Pooled
1699	563465	AA113331	2.24	18.67	7.43	1.00	0.00	5	39.19	Umbilical cord	Cervix	Esophagus
1700	785985	AA460480	11.42	72.35	6.33	1.00	0.00	6	114.01			
1701	788168	AA424833	10.20	57.27	5.62	1.00	0.00	22	345.25	Small intestine	Gall bladder	CNS
1703	843312	AA489555	31.61	254.18	8.04	1.00	0.00	X	528.46	Bone	Ear	Umbilical cord
1714	824070	AA491225	11.60	63.56	5.46	1.00	0.00	6	218.19	Stomach	Pooled	Tonsil
1716	899092	AA598784	177.33	992.01	5.59	0.00	1.00	12	87.5	Tonsil	Pancreas	Brain
1720	242578	H94949	6.68	45.86	6.87	3.00	0.00	19	104.97	Pool	LID not found	Other
1721	276547	N43830	4.34	28.92	8.68	1.00	0.00	20	73.14			
1729	197474	H52098	5.66	30.08	5.31	1.00	0.00	7	535.67	Whole embryo	Placenta	Heart
1733	246705	N59717	14.97	87.18	5.83	2.00	0.00			Lung	Pool	LID not found
1737	241824	H93217	10.37	56.53	5.45	1.00	0.00			Cervix	Pool	LID not found
1739	120834	T85234	108.50	1068.76	9.78	4.00	0.00	12	471.75	Testis	Tonsil	Brain
1741	124719	R02168	55.90	282.47	5.05	1.00	0.00	22	17.69	Adipose	Pool	LID not found
1747	120631	T85238	64.89	478.95	7.37	3.00	1.00	2	508.62	Pool	LID not found	Other
1751	240199	H69637	4.39	43.60	9.93	1.00	0.00	X	245.06	Brain	Lymph	Uterus
1753	297155	W03972	33.18	211.85	6.38	3.00	1.00	11	271.39	Pool	LID not found	Other
1756	206816	R88295	205.87	1376.11	6.68	2.00	0.00	5	596.86	Forebrain	Germ Cell	Brain
1761	201322	R88685	15.09	115.32	7.64	4.00	0.00	7	580.91	Prostate	Pool	Brain
1763	120162	T85274	9.20	172.76	18.79	5.00	0.00	12	130.31	Testis	Pool	LID not found
1768	154789	R55406	60.85	339.65	5.58	2.00	0.00	4	134.94	Pool	LID not found	Other
1772	134719	R28287	7.83	112.75	14.78	3.00	0.00	9	292.64	Pool	LID not found	Other
1773	245319	N76675	5.16	31.53	6.11	1.00	0.00	10	434.43	Stomach	Pool	LID not found
1781	285283	W04411	16.57	103.73	6.26	3.00	0.00					
1785	292207	N80622	75.72	781.78	10.06	5.00	0.00					
1786	67070	T70429	46.44	278.87	5.96	1.00	0.00					
1789	120124	T95160	12.95	112.46	8.69	5.00	0.00					
1791	122913	T98894	48.43	358.85	7.41	3.00	2.00					
1792	240638	H80890	3.44	23.82	6.86	2.00	0.00					
1793	113488	T78084	13.81	163.85	11.85	8.00	0.00					
1795	120823	T95462	3.18	23.26	7.32	5.00	0.00					
1797	124090	R02710	15.11	102.18	6.76	3.00	0.00					
1801	195641	R92841	14.11	88.52	6.27	3.00	0.00					
1802	68423	R15715	48.59	398.38	7.43	4.00	0.00					
1805	206094	H61608	28.07	258.86	9.22	4.00	0.00					

Table 2A

809	208769	H61037	18.90	161.93	8.57	5.00	0.00	1	727.12	Heart	Pool	LID not found	Other
811	120173	T65693	3.91	24.88	6.36	1.00	0.00	2	339.39			LID not found	Pool
814	236129	H33732	86.56	556.21	8.45	2.00	0.00					Lymph	Breast
816	160628	H25019	3.56	42.17	11.84	3.00	0.00					LID not found	Other
818	68430	R15709	4.24	22.12	5.22	1.00	0.00	3	697.77	Pool	Pool	LID not found	Other
825	198578	R64808	8.56	96.17	11.21	3.00	5.00					LID not found	Other
828	210710	H66883	64.09	424.11	6.62	2.00	1.00					LID not found	Other
833	204833	H77643	12.32	61.62	5.00	0.00	1.00	21	254.53	Lung	Pool	LID not found	Other
837	202482	H53224	26.35	154.92	5.96	1.00	0.00	3	726.64	Pool	Pool	LID not found	Other
838	112409	T65990	26.87	177.11	6.59	1.00	0.00					LID not found	Other
841	198582	N52811	99.41	622.47	6.26	1.00	0.00	3	124.02	Pool		LID not found	Other
844	198582	R94810	5.12	33.31	6.51	2.00	0.00	19	247.58	Pool	Fore skin	LID not found	Other
843	194872	R91033	81.01	515.27	6.36	3.00	0.00	10	114.69	Pool		LID not found	Other
845	202553	H35262	49.02	304.89	6.22	1.00	0.00	1	27.3	Pool		LID not found	Other
846	18846	R00220	58.94	402.53	6.93	3.00	0.00	3	49.43	Pool		LID not found	Other
849	275634	R94840	32.69	512.12	15.67	6.00	5.00	5	353.54	Bone		Ear	Adrenal gland
855	1046522	A6521150	3.41	25.32	7.41	2.00	0.00					Colon	Tonsil
861	202704	H53878	265.96	1460.54	5.49	1.00	0.00	X	92.32			LID not found	Other
862	123448	R00688	190.44	1272.50	6.68	2.00	0.00					LID not found	Other
868	183553	H53553	13.42	113.02	8.42	1.00	5.00	6	495.48	Pool		LID not found	Other
874	246824	N56494	5.08	50.15	9.87	4.00	0.00					LID not found	Other
875	205417	H59938	29.74	266.28	8.95	4.00	2.00	13	155.48	Pool		LID not found	Other
877	202602	H53920	28.65	219.15	7.60	6.00	0.00	X	93.95	Muscle		Ovary	Pool
878	126290	R09960	37.79	288.23	7.93	2.00	0.00	1	340.75	Pool		LID not found	Other
883	232568	H73321	47.88	358.75	7.49	3.00	2.00	22	84.53	Pool		LID not found	Other
889	199228	R95819	7.40	48.35	6.53	5.00	0.00	1	581.29	Thyroid		CNS	Breast
890	246688	R76301	36.63	443.91	12.12	6.00	3.00					Small intestine	Gall bladder
892	143790	R76762	12.03	77.18	8.42	3.00	1.00	18	-13.12	Thyroid		Adipose	Parathyroid
895	273435	N36882	187.81	854.31	5.09	1.00	0.00					LID not found	Gall bladder
1805	1695	R95851	17.87	160.66	8.99	5.00	0.00	9	82.35	Pool		LID not found	Other
1907	233579	H78482	65.07	488.83	7.51	4.00	0.00	6	118.49	Small intestine	Thymus	Lung	
1908	417711	W08967	111.37	651.16	5.85	1.00	0.00					LID not found	Other
1913	199243	R95869	11.23	103.65	9.25	3.00	3.00	5	501.86	Pancreas		Placenta	
1915	229330	H78363	183.58	1340.02	7.30	3.00	0.00					LID not found	Other
1917	202785	H53964	54.55	320.43	5.87	2.00	0.00	19	250.8	Breast		LID not found	Other
1918	1918	R70361	89.10	682.72	7.66	5.00	0.00	11	248.09	Tonsil		Gall bladder	Pancreas
1920	323917	A4284180	1.46	12.52	8.57	1.00	0.00					Colon	
1934	1834	AA457156	9.65	52.33	5.42	1.00	0.00	8	168.41	Ovary		Kidney	LID not found
1842	740925	AA479278	7.36	93.08	12.65	5.00	1.00	2	54.22	Fore skin		Gall bladder	Eye
1951	813426	AA450653	20.62	176.15	8.54	1.00	0.00	9	27.23	Gall bladder		Spleen	Kidney
1952	248281	N78093	2.50	37.43	14.96	8.00	2.00	2	357.93	Uterus		Placenta	Kidney
1958	235353	H78484	4.35	28.78	6.62	2.00	0.00	2	203.26	CNS		Gall bladder	Eye
1963	361943	AA001444	44.28	299.05	6.75	2.00	0.00	2	336.85	Aorta		Spleen	Parathyroid
1966	10213	AA464525	22.60	233.09	10.31	3.00	0.00	2	-13.12	Thyroid		Adipose	Gall bladder
1967	210522	H65034	5.93	30.58	5.15	0.00	0.00	18	130.31	Thyroid		Esophagus	
1975	133178	R28423	212.84	1548.07	7.27	1.00	0.00	12	37.19	Aorta		Placenta	Heart
1978	128159	R11490	5.44	41.98	7.72	1.00	0.00	22	40.65	Cervix		Blood	
1980	234011	H66168	19.37	258.81	13.41	4.00	0.00	10	111.21	NO OBSERVED TISSUES		Spleen	
2000	361639	W98268	29.87	243.29	8.15	4.00	0.00	7	215.58	Uterus		Placenta	
2006	128413	R06834	4.92	44.36	9.01	1.00	0.00						
2008	381612	AA058857	1.53	7.95	5.20	1.00	0.00						
2014	689845	AA233078	5.50	51.37	9.35	1.00	0.00						



Table 2A

2034	470379	AA031284	10.20	63.47	6.22	3.00	0.00	3	121.8	Blood	Germ Cell	Tonsil
2036	135608	R31395	4.41	32.04	7.27	0.00	1.00	3	117.28	Blood	Placenta	Brain
2043	41650	R52797	8.61	457.65	53.16	0.00	1.00	7	455.24	CNS	Skin	Muscle
2046	840333	AA485401	18.36	108.93	5.93	1.00	0.00	3	113.12	Larynx	Spleen	Lymph
2053	786888	AA460756	47.07	266.65	5.56	1.00	0.00	12	-5.34	Ignore	Adrenal gland	Blood
2056	340712	W55997	5.05	34.15	6.76	2.00	0.00	11	253.28	Adrenal gland	Forebrain	
2051	525566	AA064715	21.72	177.40	6.17	2.00	0.00	18	87.54			
2063	841008	AA466849	13.71	101.15	7.38	2.00	0.00	1	252.77	Umbilical cord	Thyroid	Pancreas
2070	815329	AA457047	30.58	154.37	5.05	1.00	0.00	3	69.22	Tonsil	Ear	Ovary
2073	24415	R39356	85.67	501.66	5.86	0.00	2.00	17	53.69	Umbilical cord	Pool	Pool
2081	841340	AA467637	14.27	109.00	7.64	2.00	0.00	6	117.98	Omentum	Larynx	Ear
2084	366341	AA025778	4.85	35.01	7.53	2.00	0.00	7	84.81	Thymus	Placenta	Skin
2085	511816	AA068745	38.62	204.15	5.29	1.00	0.00	11	271.39	Colon	Pool	LID not found
2092	140354	R65622	10.46	101.22	8.66	1.00	0.00	19	17.75			
2120	193586	H47475	24.52	179.74	7.33	2.00	0.00	5	529.34	Pool	LID not found	Other
2127	129853	R17054	3.67	22.41	6.11	2.00	0.00	5	110.31	Parathyroid	Uterus	Placenta
2132	141453	R68897	7.43	80.27	10.80	4.00	5.00	4	102.82	Pool	LID not found	Other
2136	193533	H47542	33.47	294.33	8.79	5.00	0.00	4		Umbilical cord	Aorta	Pool
2138	108651	T72691	6.16	44.92	7.29	0.00	1.00	19	277.06	Thyroid	Testis	Whole embryo
2142	111004	T90374	80.87	511.71	6.33	2.00	0.00	20	336.98	Aorta	Umbilical cord	Muscle
2145	294310	N64431	5.63	46.63	8.28	1.00	0.00	9	356.18	Placenta	LID not found	Other
2148	141765	R69788	5.88	34.36	5.77	1.00	0.00	16	96.57	CNS	Breast	Pool
2152	233289	H78655	32.08	285.95	8.29	5.00	0.00	X	245.06	Kidney	LID not found	Other
2156	142387	R69934	144.43	1282.55	8.88	3.00	0.00	3	726.84	Pool	LID not found	Other
2160	193817	H47928	3.69	21.37	5.80	1.00	0.00	12	241.3	Ovary	Nose	Lymph
2161	276266	R94591	84.38	700.57	8.30	5.00	0.00	15	215.11	Adipose	Ear	Thyroid
2168	295866	N67006	80.55	770.55	9.57	3.00	0.00	12	55.14	Liver	Cervix	Parathyroid
2169	137367	R38133	5.63	83.92	14.91	3.00	3.00	6	43.96	Blood	Breast	Placenta
2171	125789	R07695	46.22	359.21	7.77	4.00	0.00	18	427	Aorta	Adrenal gland	Forebrain
2178	294973	N89539	15.66	110.38	8.96	1.00	0.00	6	117.99	Pool	LID not found	Other
2180	142532	R70140	12.35	112.24	9.09	1.00	0.00	3	141.89	Breast	-	Prostate
2183	341805	W60845	12.68	67.31	5.31	2.00	0.00	14	112.35	-	Uterus	Kidney
2184	193713	H48115	5.72	43.03	7.52	2.00	0.00	6	137.73	Pool	LID not found	Other
2189	121275	T96731	47.39	333.67	7.04	2.00	0.00	5	367.75	Pool	LID not found	Other
2198	155128	R70318	61.78	415.93	6.73	2.00	0.00	14	192.57	Pool	Thymus	Pool
2198	115110	T80784	8.21	53.74	6.55	2.00	0.00	17	404.02	Germ Cell	LID not found	Other
2199	198190	R82412	46.38	340.79	7.51	4.00	0.00	8	474.66	Pool	LID not found	Other
2204	240748	H81337	32.24	396.25	12.28	6.00	5.00	14	241.84	Spleen	Pancreas	Ovary
2208	193724	H47883	17.13	279.24	16.30	6.00	0.00	14		Adrenal gland	-	Lung
2213	242084	H93842	25.24	176.76	7.00	5.00	0.00	14				
2214	294127	N71365	133.97	1144.80	8.54	3.00	0.00	14				
2216	771133	AA427762	2.17	13.96	6.43	2.00	0.00	14				
2217	228851	H68442	54.66	367.61	6.72	3.00	0.00	14				
2223	342522	W68559	1.98	10.16	5.12	1.00	0.00	17				
2226	206906	H48369	5.02	26.44	5.67	1.00	1.00	8				
2232	809532	AA456598	2.67	30.26	11.32	8.00	0.00	14				
2234	324342	W47576	6.63	43.63	5.06	1.00	0.00	9				
2235	1049291	AA620759	4.23	45.28	10.70	4.00	0.00	9				
2245	242780	H93604	6.56	62.08	9.45	1.00	0.00	5				
2247	290693	N72009	12.40	103.45	8.34	6.00	0.00	5				
2251	267241	N24561	69.63	337.04	5.06	1.00	0.00	8				
2252	805788	AA454745	108.79	628.23	5.72	2.00	0.00	8				
2254	322723	W15465	31.84	265.18	8.38	2.00	3.00	8				
2255	247582	N54244	163.85	1183.85	7.23	4.00	0.00	1				
2257	292424	N81202	17.17	108.48	6.38	3.00	0.00	1				

Table 2A

2261	242797	H94043	114.39	616.96	5.38	1.00	0.00	Umbilical cord Aorta	Whole embryo
2262	241179	H91121	41.59	383.46	9.22	5.00	0.00	Gall bladder	Spleen
2263	415946	W060376	109.92	985.51	8.78	3.00	0.00	Pool	LID not found Other
2266	213496	H70120	21.82	135.46	6.21	3.00	0.00	Pool	LID not found
2269	242698	H94282	7.03	38.98	5.28	2.00	0.00	Colon	Pool
2270	210744	H66856	4.23	37.73	8.92	6.00	0.00	Pool	LID not found Other
2271	268950	N24845	21.41	165.74	7.74	5.00	2.00	Fore skin	LID not found
2276	211859	H68719	105.30	931.91	8.65	3.00	4.00	Pool	LID not found Other
2278	325150	AA284258	15.39	93.18	6.05	3.00	0.00	Heart	Parathyroid
2280	66878	T67281	11.86	73.97	6.18	3.00	0.00	Pool	LID not found
2292	503083	AA148641	19.98	118.63	5.94	0.00	1.00	Uterus	Whole embryo Lung
2299	368578	AA026688	26.38	145.08	5.50	3.00	0.00	Small intestine	Gall bladder
2300	610813	AA464741	7.48	48.81	6.53	1.00	0.00	Synovial mem.	Muscle
2301	243199	H94578	3.49	25.42	7.29	1.00	0.00	Pool	LID not found Other
2302	235173	H73013	35.34	283.21	8.30	5.00	0.00	Pool	LID not found Other
2303	377587	AA055829	8.81	72.71	8.44	4.00	0.00	Placenta	Heart
2308	190491	H37774	42.49	288.01	7.01	3.00	0.00	Colon	Muscle
2312	43743	H04985	3.12	36.65	11.73	1.00	0.00	CNS	Eye
2316	841357	AA487452	35.60	181.59	5.10	0.00	1.00	Synovial mem	Bone
2327	768031	AA418846	3.56	30.39	8.54	3.00	0.00	Adrenal gland	Pancreas
2335	184038	H30688	1.95	28.40	14.59	7.00	0.00	Pancreas	Lymph
2348	51737	H23021	65.76	382.27	5.81	1.00	1.00	Fore skin	Brain
2358	810521	AA464644	6.25	152.12	24.34	2.00	0.00	Synovial mem	Tonsil
2367	175103	H39187	6.25	67.88	10.86	8.00	2.00	Parathyroid	Parathyroid
2378	243343	N38069	181.11	888.28	5.88	2.00	0.00	Head and nec	Umbilical cord Lymph
2404	85497	T71879	11.90	237.13	19.93	4.00	0.00	Smooth musc	Synovial mem Brain
2413	810444	AA457114	21.63	148.28	6.86	4.00	0.00	Liver	Thymus
2421	201727	R69749	10.81	83.80	8.67	10.00	0.00	Head and nec	Pooled
2423	741841	AA402879	13.62	98.50	7.23	2.00	0.00	CNS	Uterus
2424	66560	T67053	15.83	529.35	33.43	6.00	0.00	Ignore	Skin
2444	85128	T71284	19.00	230.89	12.16	3.00	0.00	Liver	Thymus
2446	417508	W80899	3.33	78.08	23.47	4.00	1.00	Stomach	Kidney
2448	85840	T72235	152.31	788.55	5.18	1.00	0.00	Liver	Gall bladder
2449	193712	R89492	11.24	66.46	5.91	2.00	2.00	Liver	Pool
2450	840687	AA488073	9.89	778.31	78.69	22.00	0.00	Stomach	Nose
2451	118914	T04283	51.83	670.76	12.94	7.00	0.00	Stomach	Pancreas
2452	704020	AA279147	4.64	27.22	5.88	2.00	0.00	Placenta	Tonsil
2453	295137	N71853	17.50	88.21	5.61	1.00	0.00	Fore skin	Kidney
2477	178890	H51574	4.22	31.31	7.43	2.00	0.00	Smooth musc	Tonsil
2478	39593	R51912	1.23	135.62	110.61	5.00	0.00	Pooled	Brain
2479	380737	AA054358	4.28	23.81	5.59	1.00	0.00	Eye	Heart
2493	159608	H15842	28.37	320.35	12.15	1.00	0.00	Peripheral ner	Breast
2497	285839	N67034	13.84	359.35	26.69	3.00	0.00	Adipose	Breast
2499	120681	T95657	40.85	307.34	7.52	5.00	0.00	Adipose	Breast
2501	128775	R10007	18.80	125.68	7.57	2.00	0.00	Pool	Uterus
2502	29063	R40970	7.26	40.07	5.52	0.00	1.00	Synovial mem	Umbilical cord Bone
2505	141854	R70598	54.14	1070.61	19.77	5.00	0.00	Esophagus	Skin
2507	186109	R92236	44.41	311.44	7.01	3.00	0.00	Pool	Brain
2509	127119	R08083	18.59	134.53	7.24	3.00	0.00	Pool	Tonsil
2514	31251	R42852	3.49	80.80	23.15	5.00	0.00	Ear	Bone
2515	120685	T85804	1.83	11.57	5.98	1.00	0.00	Synovial mem	Lymph
2516	234378	N28268	7.27	169.22	23.29	7.00	0.00	Adipose	Pooled
2518	120413	T95953	8.44	64.51	7.66	2.00	0.00	Spleen	Pool
2524	134255	R31107	2.47	14.88	6.02	3.00	0.00	Placenta	Germ Cell

Table 2A

2525	232899	H75578	135.71	1157.00	8.53	4.00	0.00	11	104.49	Spleen	Adrenal gland	Aorta
2526	126741	R07663	4.04	22.65	5.61	1.00	0.00			Pool	Pool	LID not found
2527	208264	H65569	4.01	20.77	5.18	1.00	0.00	13	74	Tonsil	Pool	LID not found
2528	282392	N88390	4.69	23.56	5.03	1.00	0.00	9	-276.55	Whole embryo	Germ Cell	Pancreas
2529	211876	H68724	12.97	78.72	6.07	2.00	2.00	13	85.4	Pool	LID not found	Heart
2530	184634	R55184	4.38	68.52	15.27	9.00	0.00	16	362.18	Pool	Breast	Heart
2531	128355	R06544	30.98	208.48	6.67	2.00	0.00			Pool	LID not found	Other
2543	122889	R00151	20.06	158.02	7.88	3.00	1.00	1	56.78			
2544	138601	R63407	4.74	36.34	7.67	2.00	0.00			Small intestine	Tonsil	Pool
2545	243260	H65823	94.16	824.13	8.75	5.00	0.00	13	107.1	Pool	Lung	LID not found
2553	128735	R14802	28.54	258.66	8.78	5.00	1.00	5	378.73	Blood	Pool	LID not found
2555	120383	T96035	3.82	54.07	14.15	3.00	0.00	12	487.41	Pool	LID not found	Other
2561	110565	T90201	15.24	158.45	10.27	7.00	0.00	12	-1.31	Pool	LID not found	Other
2563	120863	T86077	5.29	61.32	11.60	5.00	0.00	12	246.56	Placenta	LID not found	Other
2564	135094	R31426	137.52	1087.72	7.91	3.00	0.00	10	426.58	Thymus	Adipose	Tonsil
2566	366966	AA026562	51.89	513.37	9.89	4.00	4.00	17	363.86	Placenta	LID not found	Other
2569	137417	R38239	127.80	877.81	6.87	3.00	0.00	14	281.08			
2573	282879	N62885	400.94	2154.88	5.37	3.00	0.00	7	448.84			
2578	120873	T86215	26.88	185.98	7.19	3.00	0.00			Thymus	Skin	Aorta
2585	245585	N72540	13.52	83.07	6.14	0.00	1.00	1	130.77	Heart	Lung	LID not found
2590	204735	H57242	36.75	344.96	9.39	5.00	0.00	9	387.64	Ovary	Breast	Pool
2591	123459	R06808	4.31	47.67	11.05	1.00	0.00	4	495.88	Pool	Brain	LID not found
2592	214205	H77797	40.13	537.54	13.40	6.00	5.00	8	315.52	Pool	LID not found	Other
2597	297084	W03793	2.97	15.00	5.04	1.00	0.00	10	310.17	Pool	Kidney	Pool
2598	208904	H63760	3.89	22.89	5.89	1.00	0.00	11	250.51	Pool	LID not found	Other
2601	187856	R66208	25.79	220.20	8.54	3.00	0.00	8	24.47	Muscle	Heart	Pool
2606	191572	H37880	102.78	630.02	6.13	2.00	0.00	4	102.82			
2611	243414	N48139	20.60	142.27	6.91	5.00	0.00	17	486.39	Pool	Prostate	LID not found
2613	203038	H54188	1.57	13.53	8.59	1.00	0.00			Small intestine	Fore skin	Germ Cell
2614	244343	N72321	28.63	201.79	6.81	3.00	0.00			Pool	LID not found	Other
2622	185132	R91215	78.02	697.42	8.94	3.00	0.00	7	593.58	Pool	Prostate	Pool
2625	213969	H72388	3.49	18.26	5.52	2.00	0.00			Pool	LID not found	Other
2627	113300	T63928	7.44	39.45	5.30	1.00	1.00			Pool	LID not found	Other
2629	203122	H54423	60.77	453.86	7.47	2.00	3.00			Pool	LID not found	Other
2630	292471	N91231	41.37	242.75	5.87	2.00	0.00	7	593.58	Pool	Prostate	Pool
2633	180111	R98436	3.81	22.55	5.92	4.00	0.00			Pool	LID not found	Other
2634	293421	N82136	62.21	373.00	8.00	1.00	2.00	13	226.98	Fore skin	Ear	Placenta
2636	271050	N29914	3.14	56.77	18.07	1.00	0.00	8	478.59	Pool	LID not found	Other
2637	203287	H54609	21.94	135.43	8.17	4.00	2.00	10	510.78	Pool	LID not found	Other
2638	248652	N57713	172.05	1108.66	6.78	1.00	2.00	17	68.81	Pool	LID not found	Other
2641	245299	N53453	28.31	241.70	8.25	7.00	0.00	7	521.82	Testis	Pool	LID not found
2646	201651	R98262	14.85	100.28	8.70	4.00	0.00	7	655.1	Parathyroid	Liver	LID not found
2647	305608	N90248	2.26	30.08	13.30	10.00	0.00			Aorta	Pancreas	Pancreas
2648	232837	H73329	5.33	28.12	5.28	1.00	0.00			Whole embryo	Pool	LID not found
2650	284311	W02016	2.84	15.44	5.45	1.00	0.00	8	436.84	Thyroid	Ear	Bone
2651	613410	AA458840	57.11	348.14	6.10	1.00	0.00	6	137.43	Pool	LID not found	Other
2653	203388	H54811	3.70	34.01	9.18	2.00	0.00	18	14.7	Pool	LID not found	Other
2654	201203	R98287	158.79	1072.08	8.75	2.00	0.00			Pool	LID not found	Other
2657	199623	R98561	1.99	18.66	9.37	2.00	0.00			Pool	LID not found	Other
2662	202414	H52623	78.26	525.41	6.71	2.00	1.00	1	81.13	Thymus	Colon	Breast
2664	244044	N38801	16.16	258.25	15.88	5.00	0.00			Eye	Pool	LID not found
2665	199602	R98586	13.16	92.70	7.05	4.00	0.00			Pool	LID not found	Other
2666	293306	N91731	41.14	280.22	6.81	1.00	1.00			Pool	LID not found	Other
2670	249478	N58638	169.54	1088.24	6.30	2.00	0.00			Pool	LID not found	Other
2675	288787	N74086	20.84	167.00	7.88	4.00	0.00					

Table 2A

2676	243426	N49439	27.64	193.16	6.99	5.00	3.00	Pool	LID not found	Other
2681	199709	R06094	8.65	45.30	5.24	1.00	0.00	96.76 Ear	Pool	LID not found
2683	305227	W19461	172.09	829.50	5.40	1.00	0.00	631.66 Ignore	CNS	.
2685	207448	H60119	13.74	151.89	11.06	4.00	5.00	192.45 Marrow	Ovary	Pool
2686	294953	N89533	191.71	1146.72	5.98	1.00	0.00	371.28	Foreskin	Blood
2687	324815	W49563	9.08	74.88	8.25	1.00	1.00	57.43 Thyroid	Heart	Spleen
2688	115443	T67341	18.22	117.14	6.43	0.00	2.00	701.1 Germ Cell	Larynx	Pancreas
2689	204335	H59915	81.58	1116.13	12.19	10.00	0.00	347.38 Thyroid	Muscle	.
2690	753775	AA406242	13.19	149.08	11.31	1.00	0.00	58.61 Foreskin	Esophagus	Blood
2702	38684	R52542	8.30	44.79	5.40	1.00	0.00	585.14	Pool	LID not found
2712	122636	T89886	48.48	306.49	6.32	2.00	1.00	308.38 Kidney	Head and nec Bone marrow	Neural
2714	741497	AA401137	6.39	446.78	69.69	16.00	2.00	390	Gall bladder	Stomach
2716	810891	AA459519	87.81	1099.85	12.52	3.00	0.00	138.07 Eye	Tonsil	Breast
2727	24642	T60232	8.68	63.58	7.34	2.00	0.00	Thymus	Skin	Stomach
2731	160793	H24707	21.08	148.31	7.03	0.00	2.00	80.45 Blood	Eye	Pool
2743	470178	AA028963	12.86	76.40	5.94	3.00	0.00	553.7 Placenta	Aorta	Germ Cell
2755	813149	AA456895	11.37	98.21	8.46	3.00	0.00	107.37 Ignore	Parathyroid	Whole embryo
2756	243741	N49629	4.84	255.58	55.06	8.00	0.00	-10.31 Eye	Colon	Whole embryo
2759	50413	H17975	2.89	18.54	6.21	1.00	0.00	41.94 Ear	Bone	Lung
2760	823859	AA490688	8.33	52.06	6.25	1.00	0.00	58.38 Germ Cell	Whole embryo	.
2762	41565	R66447	3.95	45.07	11.41	2.00	0.00	Foreskin	Brain	Breast
2763	261204	H98218	13.30	620.78	46.88	0.00	2.00	250.6 Liver	Gall bladder	.
2779	155718	R72075	4.68	49.41	10.58	0.00	0.00	280.97 Tonsil	Bone	Pool
2783	85509	T71888	5.75	30.52	5.31	1.00	0.00	414.87	Gall bladder	Adipose
2780	815503	AA450869	26.57	134.44	5.06	0.00	1.00	24.9 Ovary	Cervix	Pancreas
2799	700527	AA291163	42.21	288.07	6.83	1.00	0.00	Placenta	LID not found	Other
2802	843028	AA488406	9.55	337.36	35.33	3.00	0.00	239.18 Periphereal ner	Thymus	Thyroid
2812	136802	R36175	19.30	156.71	8.12	2.00	0.00	458.37 Small intestine	Skin	Muscle
2813	207029	H46420	83.05	4603.26	7.05	0.00	0.00	144.44 Whole embryo	Brain	Pool
2815	51702	H22856	47.23	253.57	5.37	0.00	1.00	249.3 Pool	LID not found	Other
2819	40580	R55130	4.88	25.65	5.26	1.00	0.00	Periphereal ner	CNS	.
2820	195200	R91266	7.97	40.67	5.10	1.00	0.00	Placenta	LID not found	Other
2821	22731	T75041	14.95	465.42	31.13	2.00	0.00	362.25 Bone marrow	Nose	Aorta
2824	131365	R23089	16.85	96.08	5.70	3.00	2.00	198.83 Aorta	Spleen	.
2829	703581	AA278759	81.01	710.88	6.77	0.00	1.00	147.41 CNS	Uterus	Tonsil
2830	84820	T74808	145.69	849.43	5.83	0.00	0.00	819.8 Umbilical cord	Gall bladder	Thymus
2832	212772	H70099	4.82	38.61	8.01	2.00	0.00	Thymus	Skin	Whole embryo
2836	814054	AA465479	5.26	43.37	6.25	5.00	0.00	245.8 Aorta	Heart	Aorta
2850	796757	AA460727	96.74	632.23	6.54	0.00	1.00	Esophagus	Breast	Tonsil
2852	248543	N77515	104.96	620.12	5.91	3.00	0.00	432.9 Adipose	Pool	Tonsil
2860	471198	AA034213	73.56	426.01	5.78	1.00	0.00	Germ Cell	LID not found	Other
2866	712023	AA281816	3.09	16.09	5.21	0.00	1.00	Pool	LID not found	Other
2882	123408	R00395	20.55	119.71	5.83	3.00	0.00	180.04 Small intestine	Head and nec Placenta	.
2883	127098	R08109	10.81	66.23	6.13	1.00	0.00	218.82 Pool	LID not found	Other
2888	242687	H93550	4.17	25.20	6.05	1.00	0.00	284.32 Tonsil	Uterus	Testis
2892	142733	R71414	4.91	32.70	6.86	1.00	1.00	378.6 Prostate	Pool	LID not found
2893	327506	W32731	47.97	281.72	5.87	3.00	0.00	362.64 Testis	Pool	Brain
2896	202621	H53791	4.76	29.69	6.24	3.00	0.00	Breast	Whole embryo	Tonsil
2897	292612	N90368	142.16	1082.50	7.81	2.00	0.00	635.65 Pool	LID not found	Other
2898	108836	T77812	3.22	20.05	6.22	3.00	0.00	Pool	LID not found	Other
2908	156270	R72681	26.97	157.44	5.84	5.00	0.00	Placenta	LID not found	Other
2920	207860	H60503	9.74	101.79	10.45	4.00	4.00	Pool	LID not found	Other
2923	127216	R08275	14.93	110.80	7.42	4.00	2.00	Pool	LID not found	Other
2924	141522	R73075	9.98	67.30	6.76	1.00	0.00	Pool	LID not found	Other
2928	197637	R87184	137.97	907.45	6.58	4.00	0.00	Pool	LID not found	Other

Table 2A

2932	160730	H24816	7.81	44.90	5.90	1.00	0.00	Breast	LID not found	Other
2936	241171	H80423	8.70	62.63	7.20	3.00	0.00	Parathyroid	Breast	Germ Cell
2938	128331	R12287	19.07	150.41	7.89	3.00	0.00	Pool	LID not found	Other
2944	194906	R90957	12.00	143.03	11.92	6.00	0.00	Pool	LID not found	Other
2947	127768	R08680	35.47	233.77	6.59	3.00	0.00	Pool	LID not found	Other
2949	251260	H86334	6.87	70.96	12.52	3.00	0.00	Pool	Forebrain	Muscle
2954	241539	H90603	69.25	587.51	8.48	5.00	0.00	Pool	LID not found	Other
2960	194985	R91004	24.84	325.60	13.11	6.00	5.00	Pool	LID not found	Other
2961	211202	H67686	193.76	1244.58	6.42	3.00	0.00	Pool	LID not found	Other
2962	246116	N52394	11.89	74.15	6.24	4.00	0.00	Pool	LID not found	Other
2968	195037	R91031	12.01	155.50	12.94	5.00	5.00	Pool	Testis	LID not found
2970	292492	N80381	26.69	146.11	5.47	2.00	0.00	Pool	Brain	Brain
2971	127409	R08761	22.48	216.48	9.63	5.00	0.00	Pool	LID not found	Other
2986	135226	H68721	76.49	478.38	6.09	3.00	0.00	Pool	Smooth muscle	Esophagus
2980	214614	H73661	35.76	224.01	6.26	3.00	0.00	Pool	Umbilical cord	Liver
2991	415178	W95104	2.17	14.15	6.53	4.00	1.00	Pool	LID not found	Other
2994	295389	W05000	74.27	492.35	6.63	1.00	0.00	Pool	Breast	Pool
2997	245401	N77198	6.96	50.95	7.32	4.00	3.00	Pool	LID not found	Other
3001	212098	H68932	5.76	37.09	6.44	1.00	1.00	Pool	LID not found	Other
3005	243245	H95044	4.85	38.31	7.90	3.00	1.00	Pool	LID not found	Other
3009	295497	W23541	3.31	19.24	5.81	2.00	0.00	Pool	Head and neck	Breast
3011	808587	AA456611	27.81	173.02	6.22	1.00	0.00	Pool	Muscle	Tonsil
3014	235419	H77707	15.20	125.29	8.25	4.00	3.00	Pool	Whole embryo	Tonsil
3015	756332	AA404276	66.84	455.79	8.02	6.00	0.00	Pool	Uterus	Pool
3017	234080	H89004	31.02	264.66	8.53	1.00	3.00	Pool	Eye	LID not found
3021	243317	H95086	161.20	1617.25	10.03	3.00	2.00	Pool	Kidney	LID not found
3022	239711	H79850	10.95	136.03	12.42	4.00	5.00	Pool	Heart	Pool
3023	214503	H71224	23.34	228.49	9.83	5.00	0.00	Pool	LID not found	Other
3024	810754	AA457728	21.29	197.78	9.29	3.00	1.00	Pool	Ovary	Pooled
3025	245784	N76873	46.42	299.06	6.44	1.00	0.00	Pool	LID not found	Other
3030	245534	N7321	8.35	43.04	6.78	1.00	1.00	Pool	Forebrain	LID not found
3031	271378	N34751	9.18	70.64	7.71	6.00	0.00	Pool	Placenta	LID not found
3033	233399	H77697	54.37	446.37	8.21	2.00	3.00	Pool	Prostate	LID not found
3038	134537	R27733	3.80	263.90	69.52	1.00	0.00	Pool	Prostate	LID not found
3046	245774	N76878	27.27	155.30	5.70	3.00	1.00	Pool	Prostate	LID not found
3048	212712	H69653	67.41	474.67	7.04	4.00	3.00	Pool	Prostate	LID not found
3055	269951	N26072	59.86	384.91	6.62	2.00	1.00	Pool	Thyroid	Pool
3061	362894	AA018134	113.95	687.63	6.03	1.00	0.00	Pool	Adrenal gland	Pool
3063	128277	R06284	12.84	98.69	7.88	4.00	0.00	Pool	Whole embryo	Liver
3069	245409	N45012	104.92	582.10	5.55	2.00	0.00	Pool	Small intestine	Lymph node
3071	195034	R88764	17.51	212.01	12.10	7.00	0.00	Pool	Testis	Thymus
3084	548857	AA115919	123.77	648.70	5.24	2.00	0.00	Pool	Adipose	Tonsil
3090	823864	AA480680	10.94	58.73	5.37	1.00	0.00	Pool	Adipose	Adipose
3100	213136	H69582	6.00	33.11	5.52	1.00	0.00	Pool	Adipose	Adipose
3102	245296	N72452	5.30	37.09	7.00	1.00	0.00	Pool	Adipose	Adipose
3106	240706	H60215	17.79	102.85	5.78	0.00	1.00	Pool	Adipose	Adipose
3150	229365	H74265	4.25	24.78	5.83	1.00	0.00	Pool	Adipose	Adipose
3152	207794	H59000	3.40	39.81	11.69	1.00	0.00	Pool	Adipose	Adipose
3158	240009	H82419	56.27	305.55	5.43	1.00	0.00	Pool	Adipose	Adipose
3169	128795	R07167	7.20	100.24	13.92	1.00	0.00	Pool	Adipose	Adipose
3170	840942	AA486627	30.39	438.19	14.42	5.00	0.00	Pool	Adipose	Adipose
3173	246765	N53169	55.02	308.72	5.81	1.00	0.00	Pool	Adipose	Adipose
3183	416833	W86653	18.89	95.85	5.07	1.00	0.00	Pool	Adipose	Adipose
3185	249688	H84153	7.63	58.41	7.66	1.00	0.00	Pool	Adipose	Adipose
3201	220096	H82535	3.03	49.04	18.18	2.00	1.00	Pool	Adipose	Adipose

Table 2A

3203	789232	AA450227	85.13	486.26	5.48	0.00	1.00	1	550.58	Marrow	Synovial mem	Cervix	
3207	814306	AA459318	24.32	217.93	8.98	3.00	0.00	8	400.57	Small intestine	Stomach	Pancreas	
3209	28475	R13434	28.18	158.75	5.63	1.00	0.00	1	188.67	Small intestine	Smooth muscle	Gall bladder	
3214	823590	AA487051	8.08	138.70	16.97	9.00	1.00	18	181.31	Prostate	CNS	Lung	
3217	R67147	R67147	2.00	19.14	9.59	2.00	0.00						
3220	897770	AA598508	16.06	632.20	38.37	13.00	1.00						
3222	592801	AA160852	4.84	94.31	20.34	4.00	0.00	14	207.23	Larynx	Nose	Pancreas	
3231	897594	AA498337	41.59	230.89	5.55	1.00	0.00	8	118.71	Synovial mem	Skin	Esophagus	
3234	809598	AA442864	31.28	208.43	8.60	2.00	0.00	9	373.87	Thymus	CNS	Adipose	
3241	247816	N73030	53.87	284.92	5.28	1.00	1.00	6	118.59	Bone	CNS	Spleen	
3242	80109	T63324	18.89	276.39	16.74	7.00	0.00	2	219.22	Stomach	Nose	Gall bladder	
3245	755506	AA419108	108.93	567.02	5.21	0.00	1.00	16	98.33	CNS	Thyroid	Cervix	
3247	767345	AA418564	7.17	44.04	6.14	2.00	0.00	7	548.17	Testis	Lymph	Breast	
3250	824547	AA480920	19.55	103.99	5.32	1.00	0.00	20	11.08	Adipose	Thyroid	Bone	
3252	841684	AA487560	32.66	257.47	7.88	0.00	1.00			Mouth	Adrenal gland	Whole embryo	
3258	786087	AA448755	13.25	115.35	8.71	4.00	0.00	14	111.85	Whole embryo	Ovary	Muscle	
3264	544839	AA074677	3.03	34.92	11.51	4.00	2.00	18	380.27	Testis	Pancreas	Placenta	
3265	202158	H52361	10.54	72.86	8.92	3.00	0.00	12	248.58	Adipose	Colon	Heart	
3266	138021	R63197	6.92	36.13	5.22	1.00	0.00			Placenta	LID not found	Other	
3267	121239	T98708	12.55	157.22	12.53	2.00	4.00			Placenta	LID not found	Other	
3271	188214	R92577	72.28	552.40	7.84	4.00	0.00			Placenta	Placenta	Pool	
3272	207098	H48502	307.87	2229.97	7.24	5.00	0.00			Placenta	Prostate	LID not found	
3273	245073	N55583	72.28	552.40	7.84	4.00	0.00			Placenta	Prostate	LID not found	
3276	134229	R31865	3.44	22.15	8.44	1.00	0.00			Placenta	Prostate	LID not found	
3279	123508	R00828	4.61	44.43	9.84	3.00	0.00			Placenta	Prostate	LID not found	
3283	121355	T98780	18.60	126.34	7.61	4.00	1.00			Placenta	Prostate	LID not found	
3285	282452	N88424	78.98	454.30	5.80	2.00	0.00			Placenta	Prostate	LID not found	
3287	123817	R00848	188.14	1165.32	6.19	3.00	0.00			Placenta	Prostate	LID not found	
3291	121341	T98870	6.48	37.82	5.85	2.00	0.00			Placenta	Prostate	LID not found	
3295	211024	H65775	28.50	208.43	7.31	4.00	0.00			Placenta	Prostate	LID not found	
3296	139764	R62241	106.10	691.48	6.52	1.00	0.00			Placenta	Prostate	LID not found	
3301	120544	T95503	2.25	11.35	5.04	1.00	0.00			Placenta	Prostate	LID not found	
3304	139835	R62288	10.56	63.44	7.90	1.00	0.00			Placenta	Prostate	LID not found	
3307	121412	T96809	43.61	311.39	7.14	4.00	0.00			Placenta	Prostate	LID not found	
3315	121415	T96819	12.00	98.20	8.18	0.00	2.00			Placenta	Prostate	LID not found	
3317	358182	W95346	5.20	57.48	11.04	1.00	0.00			Placenta	Prostate	LID not found	
3318	307255	W21482	9.09	65.35	7.18	5.00	0.00			Placenta	Prostate	LID not found	
3320	142884	R71190	4.16	60.30	14.50	5.00	0.00			Placenta	Prostate	LID not found	
3322	261587	N24268	28.01	163.84	5.64	2.00	0.00			Placenta	Prostate	LID not found	
3323	207379	H58884	70.32	458.48	6.52	3.00	0.00			Placenta	Prostate	LID not found	
3329	470378	AA031388	8.23	71.20	8.65	3.00	1.00			Placenta	Prostate	LID not found	
3332	135450	R32751	97.59	883.38	9.05	4.00	0.00			Placenta	Prostate	LID not found	
3336	138999	R62853	5.57	32.68	5.87	1.00	2.00			Placenta	Prostate	LID not found	
3340	135454	R32754	6.77	39.47	5.83	1.00	0.00			Placenta	Prostate	LID not found	
3343	244058	N45440	3.36	22.13	6.58	0.00	2.00			Placenta	Prostate	LID not found	
3347	121206	T97078	15.26	105.16	6.89	4.00	2.00			Placenta	Prostate	LID not found	
3348	135203	R32839	13.55	74.38	5.48	0.00	3.00			Placenta	Prostate	LID not found	
3353	197833	R66358	81.26	418.18	5.15	1.00	0.00			Placenta	Prostate	LID not found	
3381	295600	N88843	6.53	54.17	8.29	7.00	0.00			Placenta	Prostate	LID not found	
3365	204038	H58997	55.60	1040.08	18.71	6.00	5.00			Placenta	Prostate	LID not found	
3368	295324	W04231	78.49	541.43	8.81	2.00	0.00			Placenta	Prostate	LID not found	
3368	487819	AA043484	3.41	17.55	5.15	1.00	0.00			Placenta	Prostate	LID not found	
3374	233214	H75898	116.84	721.17	6.18	1.00	0.00			Placenta	Prostate	LID not found	
3377	201519	R97269	7.31	112.00	15.33	4.00	2.00			Placenta	Prostate	LID not found	
3381	203781	H56207	50.48	288.06	5.80	1.00	0.00			Placenta	Prostate	LID not found	

Table 2A

3384	309515	N94385	2.84	45.06	17.09	2.00	0.00	19	101.92	Ignore	Bone	Pancreas
3386	288010	N67041	78.50	428.36	5.46	1.00	0.00			Forebrain	Eye	Heart
3387	344133	W73782	112.51	704.79	8.26	1.00	0.00			Pancreas	Colon	Heart
3390	230341	H80858	14.08	90.95	6.46	5.00	1.00	17	347.35	Gall bladder	LID not found	Other
3391	343646	W68471	1.09	8.98	8.27	1.00	0.00	3	461.06	Pool	Heart	Uterus
3393	200418	R87234	23.48	228.54	8.74	0.00	3.00			Adipose	LID not found	Other
3395	108330	T70850	1.49	15.29	10.25	3.00	0.00	17	272.44	Pool	Testis	Breast
3397	203405	H58424	14.15	98.25	6.94	6.00	1.00			LID not found	Other	
3398	243460	N48213	22.41	136.79	6.10	2.00	0.00	15	283.15	Pool	Ovary	Blood
3401	245366	N54893	12.82	285.33	20.69	8.00	5.00			LID not found	Other	
3404	190340	H38088	28.42	178.78	6.29	1.00	0.00	10	178.79	Pool	Stomach	Forebrain
3405	203850	H58438	28.60	200.93	7.03	4.00	0.00			LID not found	Other	
3409	283128	N81997	5.30	31.43	5.94	2.00	1.00	4	444.57	Head and neck	Thymus	Eye
3412	758356	AA404268	9.57	59.37	8.20	1.00	0.00	20	193	Pool	LID not found	Other
3414	283358	N92034	98.84	593.71	6.13	3.00	0.00	14	108.8	Pool	LID not found	Other
3421	203956	H56876	2.51	17.91	7.14	1.00	0.00	2	466.93	Peripheral nervous system	Ovary	
3426	284167	W01893	18.46	111.09	5.71	2.00	0.00			Pool	LID not found	Other
3427	135220	R32944	28.97	148.09	5.49	1.00	0.00			Pool	LID not found	Other
3429	234469	H65358	6.83	36.19	5.30	1.00	0.00			Pool	LID not found	Other
3433	243194	H94571	3.34	17.83	5.27	0.00	1.00			Pool	LID not found	Other
3437	204558	H56981	51.27	369.39	7.21	3.00	0.00	12	500.87	Blood	Tonsil	Whole embryo
3441	206781	R88074	12.65	158.14	12.35	6.00	5.00			Pool	LID not found	Other
3442	287437	W03668	43.11	357.35	8.29	3.00	0.00	4	134.94	Pool	LID not found	Other
3444	344568	W73140	2.78	102.48	38.91	6.00	1.00			Larynx	Heart	Testis
3446	310034	W24161	6.82	47.02	6.90	4.00	0.00			Umbilical cord	Forebrain	Prostate
3450	283358	N92035	26.95	163.02	6.26	1.00	0.00			Heart	LID not found	Other
3451	228580	H66011	39.52	232.23	5.88	1.00	0.00			Uterus	Tonsil	Whole embryo
3453	204624	H57017	42.99	263.42	8.13	3.00	0.00			Blood	LID not found	Other
3463	183476	H45617	35.02	384.82	10.99	3.00	0.00			Bone marrow	Skin	Breast
3467	134270	R31168	24.16	430.88	17.83	1.00	1.00	1	695.02	Small intestine	Thyroid	Parathyroid
3470	767638	AA418251	3.03	25.64	8.47	1.00	0.00	7	152.22	Thyroid	Umbilical cord	Pool
3472	756490	AA439406	7.11	48.64	6.84	2.00	0.00	17	325.76	Whole embryo	Liver	Tonsil
3475	810083	AA465021	2.84	22.95	8.10	1.00	0.00	16	28.14	Cervix	Blood	Lymph
3480	134748	R28284	40.70	207.56	5.10	1.00	0.00			Blood	Kidney	Muscle
3482	136188	R33154	9.84	100.01	102.48	4.00	1.00	4	34.95	Uterus	Heart	Placenta
3492	49518	H15707	198.87	1000.80	5.08	1.00	0.00	8	387.23	Smooth muscle	Umbilical cord	Pancreas
3494	288696	N62820	24.35	288.90	12.19	2.00	1.00	1	738.14	CNS	Eye	Colon
3499	234191	H64324	18.88	143.44	7.60	1.00	1.00	2	120.76	Adipose	Tonsil	Bone
3503	132122	R26070	15.52	93.92	6.05	1.00	0.00	1	132.58	Forebrain	Bone	Brain
3512	753467	AA406551	19.44	195.97	10.08	0.00	1.00	1	142.99	Adipose	Nose	
3514	261971	N27227	0.63	7.82	9.54	1.00	0.00	16	490.28	Larynx	Adrenal gland	Germ Cell
3516	360556	AA026631	13.90	71.52	5.15	0.00	1.00	5	524.57	Smooth muscle	Pancreas	
3522	754509	AA410591	6.56	48.64	7.41	1.00	1.00	7	554.6	Forebrain	Pool	Placenta
3526	199945	R67068	44.83	572.27	12.77	1.00	2.00	20	212.78	Placenta	Lymph node	Aorta
3530	296880	W01240	20.92	105.51	5.04	1.00	0.00	X	354.25	Thyroid	CNS	Liver
3538	209814	R68851	9.24	107.94	11.68	1.00	2.00	3	571.11	Small intestine	Smooth muscle	Kidney
3544	587927	AA085597	39.89	231.92	5.81	1.00	0.00	3	880.68	Thymus	Forebrain	Gall bladder
3556	841278	AA486836	13.35	148.85	11.00	2.00	0.00			Parathyroid	Adrenal gland	Lung
3561	148421	H12312	21.72	197.39	9.09	3.00	3.00	4	249.05	Placenta	Tonsil	LID not found
3576	282731	N91428	25.94	157.35	6.07	1.00	1.00			Cervix	Ovary	Kidney
3580	122822	T86688	28.51	168.38	5.91	0.00	2.00	17	96.5	Gall bladder	Stomach	Adrenal gland
3588	112885	T87139	24.97	140.54	5.93	1.00	2.00			Blood	Tonsil	Pool
3594	842863	AA469261	52.96	307.36	5.80	3.00	0.00	8	499.22	Larynx	Peripheral nerve	Nose
3603	530814	AA070226	187.57	1242.16	6.82	1.00	0.00	5	155.62	Small intestine	Liver	
3605	897689	AA498810	30.52	178.88	5.79	2.00	0.00			Bone	Whole embryo	Eye

Table 2A

3626	76869	T51162	9.45	54.67	5.78	1.00	0.00	Lymph node	Nose	Ovary
3629	824568	AA490881	14.10	103.31	7.33	2.00	0.00	Prostate	Blood	Breast
3639	814326	AA49588	12.01	123.11	10.25	4.00	1.00	Thymus	Ovary	Blood
3649	122345	T89191	2.44	27.91	11.42	3.00	0.00	Eye	Pool	LID not found
3652	246335	N59766	44.57	265.36	5.95	3.00	0.00	Tonsil	Uterus	Testis
3653	469762	AA027864	2.11	19.23	9.11	0.00	2.00	Uterus	Brain	LID not found
3654	111684	T84998	2.86	14.64	5.00	1.00	0.00	Breast	Muscle	-
3657	275612	R83354	3.97	31.30	7.86	2.00	0.00	Whole embryo-		Kidney
3664	185358	R89539	23.61	239.05	10.13	5.00	0.00	Pool	LID not found	Other
3668	207813	H59056	4.73	140.02	29.63	2.00	0.00	313.88		
3670	121072	T96523	15.40	85.85	5.59	1.00	1.00	Pool	LID not found	Other
3672	195381	R89999	29.51	287.96	8.78	5.00	0.00	354.68	Pool	Testis
3674	234318	H93238	51.27	376.88	7.35	4.00	0.00	569.13	CNS	
3688	185553	R91821	11.20	188.27	16.81	6.00	5.00	283.63	Pool	LID not found
3691	127514	R08868	46.71	325.72	6.97	3.00	0.00	378.75	Pool	LID not found
3696	195346	R89216	8.03	57.99	7.22	2.00	4.00	327.26	Pool	LID not found
3699	127542	R08683	35.08	273.00	7.78	4.00	0.00	Pool	Brain	LID not found
3700	244289	N75715	20.92	170.05	8.13	3.00	0.00	8cns	Pool	Lung
3704	242011	H93319	63.22	845.51	10.21	6.00	4.00	249.15	Pool	LID not found
3708	113431	T76571	16.94	120.58	7.12	2.00	0.00	Pool	LID not found	Other
3710	111200	T84381	22.65	165.90	7.33	5.00	0.00	Pool	Pool	Kidney
3711	130758	R22088	4.18	28.71	6.87	4.00	0.00	Placenta	Pool	LID not found
3712	185784	R89285	25.92	255.51	8.85	5.00	0.00	Pool	Pool	LID not found
3713	308989	W25368	4.33	150.87	34.82	17.00	1.00	Esophagus	Breast	Prostate
3715	127943	R09153	38.48	235.67	6.13	2.00	2.00	Thyroid	Parathyroid	Spleen
3718	130791	R22113	13.46	113.53	8.43	5.00	0.00	Forekin	Pool	Testis
3722	113538	T79129	14.65	110.89	7.56	3.00	0.00	Breast	Lung	Placenta
3727	130901	R22085	29.02	180.51	6.22	2.00	2.00	Pool	LID not found	Other
3730	247710	N68198	248.83	1359.00	5.46	3.00	0.00	382.83	Pool	LID not found
3732	262498	N81290	7.16	44.85	6.23	3.00	0.00	251.63	Pool	LID not found
3736	185568	R89471	4.65	28.41	6.46	3.00	0.00	Pool	Colon	LID not found
3738	127710	R08498	98.35	647.47	6.58	3.00	0.00	550.79	Skin	-
3741	470061	AA028041	27.86	155.53	5.59	1.00	0.00	Placenta	Heart	Breast
3743	130824	R22239	54.71	366.58	6.70	1.00	1.00	Tonsil	Pool	LID not found
3745	234419	H95342	1.76	16.53	8.41	1.00	0.00	Parathyroid	Blood	Eye
3747	257391	N30708	60.87	536.61	8.82	4.00	0.00	Adrenal gland	Pooled	CNS
3750	132871	R27505	13.70	143.05	10.44	3.00	2.00	Whole embryo	Testis	Forekin
3755	194524	R86333	87.50	522.40	5.97	2.00	0.00	Blood	Germ Cell	Placenta
3758	246194	N77006	55.38	330.01	5.85	3.00	0.00	Thymus	Pooled	Aorta
3760	356635	W84612	7.10	50.26	7.07	2.00	0.00	Pool	LID not found	Other
3766	247466	N64285	5.49	57.99	10.57	1.00	0.00	Pool	Parathyroid	-
3768	385177	AA024666	13.77	76.55	5.56	1.00	0.00	Pool	LID not found	Other
3769	226997	H71314	4.12	27.39	6.64	2.00	0.00	Pool	Parathyroid	-
3772	480755	AA133167	2.70	16.21	6.01	1.00	0.00	Pool	Peripheral ner	Uterus
3774	246698	N78306	15.47	126.08	8.09	5.00	0.00	Pool	LID not found	Other
3778	428788	AA004671	0.94	5.58	5.92	1.00	0.00	Pool	Parathyroid	-
3783	296095	N73811	13.48	69.57	5.16	0.00	1.00	Pool	Parathyroid	-
3785	213509	H72247	5.17	52.72	10.20	3.00	0.00	Pool	LID not found	Other
3788	811066	AA485443	15.02	81.86	5.45	1.00	0.00	Pool	Adipose	Umbilical cord
3790	292770	N63846	3.42	25.89	7.57	1.00	0.00	Pool	Stomach	Pooled
3792	277003	N34987	3.11	20.89	6.76	2.00	0.00	Pool	Parathyroid	Lung
3793	213698	H72280	42.19	265.06	8.28	3.00	0.00	Pool	Parathyroid	Lung
3798	293584	N64143	11.82	82.55	7.10	3.00	0.00	Pool	Adipose	LID not found
3800	428789	AA009773	10.37	61.56	5.93	1.00	0.00	Pool	Germ Cell	Colon
3801	213535	H72259	14.24	101.45	7.13	6.00	0.00	Pool	LID not found	Other



Table 2A

3806	293835	N95107	18.39	138.39	7.53	5.00	0.00	12	485.73	Adipose	Pool	LID not found
3809	214823	H71854	36.33	234.35	6.45	3.00	2.00	5	283.38	Stomach	Umbilical cord	Spleen
3819	198256	R94456	25.79	186.78	7.63	5.00	0.00					
3822	294740	N69252	10.84	96.38	9.06	1.00	0.00					
3826	292559	N80384	9.28	64.38	6.94	1.00	4.00					
3829	243385	N48130	133.02	756.41	5.69	1.00	2.00	22	175.49	Pool	LID not found	Other
3830	297919	N70072	35.59	178.79	5.02	1.00	0.00	19	66.07	Epididymis	Brain	Heart
3834	212820	H70554	2.53	26.50	10.49	0.00	1.00					
3840	320712	V31675	4.40	36.67	9.01	4.00	0.00	15	344.93	CNS	Lymph	Parathyroid
3848	243675	N50014	91.43	897.64	9.82	4.00	0.00	20	293.98	Eye	Breast	Parathyroid
3852	813707	AA453774	4.32	169.71	39.25	0.00	1.00	1	630.42	Stomach	Kidney	Tonsil
3876	721792	AA393408	5.84	63.28	10.83	1.00	0.00	1	247.44	Thyroid	Stomach	Testis
3887	289502	N76669	168.72	877.18	5.20	1.00	0.00	1	15.07	Epididymis	Gall bladder	CNS
3902	187280	R90744	3.16	22.64	7.16	1.00	0.00	14	278.45	Blood	Placenta	Pool
3908	289337	N92646	18.04	2000.27	105.08	11.00	0.00	2	827.13			
3912	813714	AA453850	54.89	354.12	8.45	1.00	0.00	16	170.16	Synovial mem	Nose	Blood
3919	119914	T94293	56.45	466.48	8.26	4.00	0.00	15	283.4	Small intestine	Stomach	CNS
3940	841470	AA487346	22.82	169.89	7.51	2.00	0.00	10	426.08	Small intestine	Stomach	CNS
3942	310105	W24248	6.79	37.19	5.48	1.00	0.00	21	247.7	Stomach	Germ cell	Kidney
3944	130843	R22308	2.82	16.94	8.00	1.00	0.00	8	31.44	Thyroid	Small intestine	Skin
3945	301081	W07798	9.88	51.16	5.18	1.00	0.00	1	294.09	Ear	Thymus	Whole embryo
3948	890335	AA598950	76.84	550.84	7.17	2.00	0.00	7	875.52	Placenta	Uterus	Kidney
3961	134783	R31701	12.88	99.79	7.78	2.00	0.00	1	339.21	Pancreas	Nose	Thyroid
3977	381204	AA017528	6.46	33.28	5.15	1.00	0.00	1	634.12	Umbilical cord	CNS	Bone
3981	70827	T46924	3.54	1338.88	378.45	6.00	0.00	6	643.74	Head and nec	Adrenal gland	Colon
3987	841679	AA486899	21.70	276.81	10.51	2.00	0.00	3	148.77	Thymus	Spleen	Thyroid
3988	823775	AA489256	16.88	177.19	27.84	2.00	0.00	9	368.81	Smooth musc	Larynx	Esophagus
4003	825085	AA489246	7.02	195.43	22.22	9.00	0.00	5	401.58	Smooth musc	Stomach	Breast
4010	815501	AA456888	7.77	47.44	6.11	1.00	0.00	16	588.76	Thyroid	Ear	Parathyroid
4015	712341	AA405000	12.15	483.45	27.84	2.00	0.00	2	650.84	Smooth musc	Placenta	Umbilical cord
4018	203132	H54629	11.87	447.40	37.70	6.00	0.00	1	594.85	Stomach	Blood	Placenta
4025	594050	AA114226	5.52	33.51	6.07	1.00	0.00	4	355.71	Placenta	LID not found	Other
4027	321388	W02272	5.04	41.69	8.27	1.00	0.00	1	594.85	Placenta	Whole embryo	Lung
4028	69672	T53626	6.48	64.26	9.92	0.00	1.00	1	355.71	Tonsil	Lung	Kidney
4032	235155	H79353	7.53	62.04	8.24	5.00	0.00	4	355.71	Placenta	Eye	Pool
4035	121214	T97080	1.88	10.76	5.43	1.00	0.00	3	560.21	Brain	Placenta	LID not found
4040	136051	R62826	1.72	16.20	9.42	2.00	0.00	19	180.72	Adrenal gland	Placenta	Tonsil
4043	137793	R68537	4.64	45.25	8.10	2.00	0.00	4	301.37	Parathyroid	Placenta	Kidney
4047	124128	R01348	4.60	38.45	8.36	1.00	0.00	4	301.37	Neural	Synovial mem	Skin
4049	247901	N77671	125.22	844.82	7.55	5.00	0.00	6	457.37	Breast	Germ cell	Pancreas
4056	137981	R63134	6.50	51.40	7.81	1.00	0.00	8	457.37	Brain	Germ cell	Lung
4060	136114	R33285	4.56	28.58	6.27	2.00	0.00	19	84.49	Adrenal gland	Placenta	Heart
4064	137989	R63137	5.87	47.56	8.10	2.00	0.00	8	42.15	Head and nec	Pooled	Other
4067	247216	N59057	25.34	257.17	10.15	2.00	0.00	4	301.37	Placenta	LID not found	Other
4069	347887	W81562	75.85	403.17	5.26	1.00	0.00	6	457.37	Placenta	Germ cell	Pancreas
4070	308092	N95361	27.95	172.23	6.25	1.00	0.00	19	84.49	Head and nec	Pooled	Heart
4078	201301	R99573	10.68	74.59	6.69	3.00	0.00	8	42.15	Head and nec	Pooled	Other
4085	140197	R68101	3.32	52.65	15.86	3.00	0.00	4	301.37	Placenta	LID not found	Other
4088	138579	R63285	6.40	63.99	10.00	2.00	0.00	6	457.37	Placenta	Germ cell	Pancreas
4092	135999	R33570	14.14	101.19	7.16	3.00	0.00	4	301.37	Bone	Colon	Whole embryo
4098	141209	R66533	25.56	253.22	9.91	2.00	0.00	1	638.07	Pool	LID not found	Other
4099	121501	T97309	12.11	87.60	7.23	3.00	2.00	4	422.24	Placenta	LID not found	Other
4107	121808	T97427	12.66	124.60	9.84	6.00	2.00	1	638.07	Pool	LID not found	Other
4116	135853	R33699	7.53	61.45	8.16	1.00	0.00	4	422.24	Placenta	LID not found	Other
4124	136909	R39730	150.48	1600.38	10.63	1.00	0.00					

Table 2A

4125	468019	AA053285	2.22	15.91	7.63	1.00	0.00	10	45.1 Uterus	Liver	Brain
4129	205849	R88107	21.08	131.91	6.26	1.00	0.00	11	159.67 Pool	LID not found	Other
4130	309665	W30810	4.80	29.05	6.05	0.00	0.00	3	121.68 Cdon	-	Pancreas
4134	120701	T95668	7.57	40.56	5.35	0.00	1.00	15	258.87 Foreskin	Placenta	Uterus
4137	200780	R88191	12.58	145.94	11.60	4.00	5.00	14	276.5 Pool	Whole embryo	Testis
4139	840683	AA488072	9.23	150.41	16.29	0.00	1.00	10	432.98 Heart	Lung	Pool
4140	261936	H88856	12.97	81.78	6.31	0.00	2.00		Small intestine	Cervix	Gall bladder
4144	762635	AA447569	16.88	116.75	6.93	3.00	0.00		Bone marrow	Larynx	Esophagus
4149	204814	H57111	9.53	53.97	5.67	3.00	0.00	10	439.93 Heart	Testis	Pool
4150	132789	R27412	5.16	42.80	6.30	4.00	0.00	7	494.74 Placenta	Muscle	Pool
4151	586888	AA128089	8.98	35.65	5.11	1.00	0.00	3	160.78 Cervix	Colon	Blood
4153	230191	H64849	7.65	91.50	11.96	4.00	5.00	10	336.96 Pool	LID not found	Other
4156	468422	AA044662	7.03	47.35	6.73	2.00	0.00	19	32.22 Smooth musc	Adrenal gland	Foreskin
4157	204888	H57273	9.16	287.13	31.34	7.00	2.00	11	289.86 Pool	LID not found	Other
4158	345751	W72621	7.76	46.42	5.99	1.00	0.00	3	211.28 Smooth musc	Placenta	Uterus
4164	324712	AA284234	3.74	23.64	6.32	1.00	0.00	11	255.63 Pooled	-	Blood
4165	204444	H58001	31.62	169.79	5.37	1.00	0.00	15	16.25 Pool	LID not found	Other
4171	128532	R10043	3.37	18.09	5.37	1.00	0.00	10	434.43 Kidney	Breast	Colon
4173	209468	H65231	12.17	107.52	8.83	2.00	4.00	10	198.24 Thyroid	LID not found	Other
4174	209381	H64095	4.16	32.30	7.72	1.00	0.00	3	77.59 Lung	Pool	Heart
4175	770593	AA434160	0.63	3.66	5.81	0.00	1.00	1	235.02 Gall bladder	Ovary	Uterus
4176	810558	AA464588	87.06	1082.48	12.43	1.00	1.00	19	Aorta	Germ Cell	Skin
4178	293500	N94080	5.97	32.51	5.44	1.00	0.00		Foreskin	Germ Cell	Pancreas
4180	268211	N29986	3.84	44.57	12.24	4.00	0.00		Aorta	Tongue	Lung
4189	234647	H77736	9.38	59.93	6.39	5.00	0.00		39.26 Peripheral nervous system	LID not found	Other
4191	783597	AA468839	18.58	134.35	7.23	1.00	0.00	14	615.85 Umbilical cord	Thymus	Pooled
4192	1048810	AA821342	22.57	354.24	15.69	1.00	0.00	1	Pool	LID not found	Other
4193	201334	R98591	27.16	282.10	9.64	7.00	0.00		473.63 Placenta	Cervix	Stomach
4195	130120	R21425	20.57	106.94	5.20	0.00	1.00	5	117.11 Pool	LID not found	Other
4198	136780	R36070	107.78	626.68	5.81	1.00	1.00	20	66.76 Lung	LID not found	Other
4205	204489	H58574	29.59	289.99	9.80	6.00	5.00	2	Parathyroid	Kidney	Uterus
4207	148743	H12777	8.59	63.82	7.43	4.00	0.00	8	78.97 Pool	LID not found	Other
4208	66481	R16134	5.34	67.51	12.64	2.00	0.00		Pool	LID not found	Other
4209	206882	R98905	19.32	100.72	5.21	0.00	1.00		672.08 Heart	Testis	LID not found
4210	263675	N84274	12.54	120.05	9.59	3.00	3.00	1	51.17 Breast	Germ Cell	Kidney
4212	277138	N40919	2.46	13.50	5.46	1.00	0.00		Pool	LID not found	Other
4214	345942	W76945	56.40	323.85	5.74	2.00	2.00	2	318.41 Pool	LID not found	Other
4215	155072	R71393	2.56	74.61	29.10	1.00	0.00	11	348.77 Gall bladder	Pancreas	Uterus
4217	206885	R98913	10.72	87.46	8.16	4.00	0.00	2	Head and nec	Esophagus	Prostate
4221	207370	H58834	8.67	65.81	7.59	4.00	0.00	2	597.59 Eye	Bone	Placenta
4228	470393	AA031513	6.34	1316.45	207.75	12.00	1.00	12	45.69 Gall bladder	Kidney	Testis
4228	511428	AA128115	4.84	114.71	23.89	13.00	0.00	20	333.71 Germ Cell	Ovary	Prostate
4235	341310	W58032	7.07	41.82	5.93	1.00	0.00	11	351.63 Umbilical cord	Larynx	Esophagus
4244	727292	AA401893	8.70	67.82	7.81	2.00	0.00	17	371.63 Lung	Eye	Kidney
4246	725680	AA398334	4.31	53.22	12.35	10.00	0.00	15	271.67 Brain	LID not found	Other
4250	569115	AA143201	224.97	1902.36	8.46	0.00	1.00	11	272.9 Thymus	Ovary	Skin
4255	309161	N99243	0.98	5.31	5.35	0.00	1.00	20	334.7 Ear	Bone	Ovary
4256	163200	H44858	68.69	334.70	5.02	1.00	0.00	2	282.87 Cervix	Brain	Kidney
4260	72050	T52435	107.26	1045.66	9.75	0.00	1.00	11	78.13 Pooled	Placenta	Pool
4275	741139	AA402207	3.77	243.17	64.45	17.00	0.00	13	146.12 Spleen	Pool	CNS
4282	667482	AA227594	2.85	781.29	274.41	11.00	2.00				
4310	346552	W74377	5.38	40.76	7.57	1.00	0.00				
4312	131639	R24635	5.59	571.41	102.25	15.00	0.00				
4320	341931	AA058828	2.22	13.36	6.02	1.00	0.00				
4321	71727	T51350	9.85	55.59	5.64	1.00	0.00				

Table 2A

4322	49873	H26922	0.90	6.47	7.16	1.00	0.00	13	329.7	Synovial mem	Adipose	Breast
4323	787861	AA452378	7.02	48.04	6.84	1.00	0.00	12	277.86	Smooth musc	Aorta	Muscle
4324	124052	R02800	2.92	23.38	8.01	1.00	0.00			Pool	LID not found	Other
4326	132142	R26184	7.16	68.24	9.54	1.00	0.00	7	152.54	Placenta	Aorta	Eye
4340	120343	T87139	3.59	22.14	6.17	1.00	0.00			Pool	Bone	Pool
4347	825923	AA188901	23.06	164.27	7.12	1.00	1.00	14	16.43	Adipose	Adrenal gland	
4351	808464	AA443093	6.07	108.14	17.82	3.00	0.00	10	526.52	Cervix	Ear	Skin
4380	161992	H26176	45.02	247.82	5.50	1.00	1.00	10	436.45	Head and nec	Stomach	Breast
4393	787016	AA463565	21.12	118.97	5.63	2.00	0.00	11	65.95	Gall bladder	Fore skin	Parathyroid
4396	110503	T62617	47.17	339.52	7.20	0.00	2.00		240.08			
4397	785605	AA448988	51.23	508.87	8.93	3.00	0.00	16	170.16	Neural	Smooth musc	Placenta
4405	840840	AA486626	282.86	7043.36	24.88	2.00	0.00	8	439.33	Ear	Head and neck	
4412	140515	R66057	30.54	455.14	14.90	8.00	3.00			Placenta	Pool	LID not found
4413	361974	AA001448	12.23	323.96	28.48	3.00	0.00	7	624.62	Ear	Adrenal gland	Parathyroid
4414	72391	T51680	7.94	161.34	20.31	13.00	0.00					
4418	248624	N78263	2.16	13.50	6.24	0.00	1.00	15	167.89	Ear	Muscle	Pool
4420	139189	R68706	96.86	677.35	6.99	3.00	0.00	11	272.02		LID not found	Other
4424	209593	H97748	20.35	173.21	8.51	3.00	0.00	19	250.6	Pool	Pool	Whole embryo
4427	198981	H83233	55.12	382.12	8.93	1.00	0.00	2	198.38	Testis	Prostate	Placenta
4428	139250	R68738	28.75	502.50	17.48	8.00	5.00			Brain	LID not found	Other
4432	194670	R68682	5.87	30.28	5.16	1.00	0.00	9	350.75	Thyroid	Pool	Eye
4433	210698	H68877	7.36	39.07	5.31	0.00	1.00	8	327.75	Pool	LID not found	Other
4435	128118	R08870	2.91	14.79	5.08	1.00	0.00	1	387.26	Pool	LID not found	Other
4440	211319	H66650	9.03	91.86	10.16	2.00	3.00				Pool	LID not found
4442	293417	N92134	130.10	916.66	7.05	3.00	0.00			Pool	LID not found	Other
4448	194921	R91060	4.37	32.72	7.48	2.00	4.00			Fore skin	Tonsil	Pancreas
4458	251461	H98001	9.54	50.74	5.32	1.00	0.00	1	81.13	Thymus	Colon	Breast
4461	292833	N80491	14.03	188.93	13.47	6.00	0.00	2	229.2	CNS	Pool	Whole embryo
4462	144852	R78527	121.99	615.16	5.04	1.00	0.00			Thyroid	Pool	Whole embryo
4464	195052	R91176	31.82	652.86	20.65	6.00	5.00			Pancreas	Breast	Lung
4468	470348	AA028381	53.50	897.06	16.77	4.00	4.00	1		Who's embryo	Cervix	Eye
4469	298763	W01171	41.02	270.05	6.56	5.00	0.00			Pool	LID not found	Other
4470	111722	T91088	1.95	11.52	5.90	1.00	0.00			Tonsil	Heart	Lymph
4471	130572	R22420	6.40	32.82	5.14	1.00	0.00			Ignore	Skin	Germ Cell
4477	741841	AA402879	13.06	100.83	7.72	2.00	1.00			Pool	LID not found	Other
4480	195091	R91244	22.26	407.93	18.31	6.00	5.00	3	694.2	Ear	Breast	Pool
4492	248108	N78103	19.14	157.51	8.23	5.00	0.00	9	404.02	Pool	LID not found	Other
4496	185139	R91271	26.88	380.59	14.16	6.00	3.00	6	569.46	Pool	LID not found	Other
4498	210923	H70862	11.51	63.48	5.52	1.00	0.00			Esophagus	Placenta	Prostate
4503	130977	R22926	50.11	260.63	5.20	1.00	0.00			Pool	Kidney	LID not found
4507	128993	R10311	2.99	17.48	5.85	1.00	0.00			Pool	LID not found	Other
4510	111825	T84685	32.44	226.10	6.87	1.00	3.00			Pool	LID not found	Other
4512	195974	R92648	32.93	277.38	8.42	1.00	3.00	8	565.78		Pool	LID not found
4517	243405	N49436	177.13	1055.80	5.96	3.00	0.00	4	558.68	Fore skin	Pool	LID not found
4518	298041	W02424	41.53	500.26	12.05	4.00	5.00			Heart	LID not found	Other
4519	327221	AA284305	2.43	13.32	5.47	1.00	0.00	1	705.71	Pool	LID not found	Other
4521	232723	H72700	56.50	462.62	8.19	2.00	0.00	17	322.28	Ovary	Pool	LID not found
4528	235026	H73304	55.65	308.83	5.55	1.00	0.00	11	220.06			
4531	126239	R06372	10.28	64.93	6.32	1.00	0.00	13	83.11			
4532	504300	AA140640	2.93	15.54	5.31	1.00	0.00	9	422.75	Blood	Tonsil	Pool
4534	288182	W02630	6.55	33.19	5.07	1.00	0.00	5	504.31	Placenta	Heart	Whole embryo
4535	341569	W58368	8.33	254.68	30.58	0.00	2.00			Pool	LID not found	Other
4542	295514	W23548	134.96	977.10	7.24	3.00	0.00	6	150.71	Blood	Placenta	Pool
4543	135752	R33082	4.99	43.96	8.82	2.00	0.00				Eye	Fore skin
4545	244784	N54401	12.56	101.24	8.06	1.00	1.00	6	391.02	Bone		

Table 2A

4550	796198	AA461108	4.86	97.87	20.14	6.00	5.00	13	297.84	Prostate	Pool	LID not found
4557	130027	R19406	127.45	858.30	6.73	3.00	0.00	14	192.96	Lung	-	LID not found
4558	301678	N78558	45.80	914.31	19.98	6.00	5.00	14	207.23	Pool	LID not found	Other
4561	293990	N95656	73.34	583.30	7.68	5.00	0.00	1	695.13	Pooled	Tonsil	Pleocenta
4562	139358	R62339	3.30	16.00	5.45	1.00	0.00	3	141.89	Germ Cell	Testis	Pancreas
4565	296180	W00899	51.70	346.25	6.70	3.00	0.00	9	410.91	Bone	CNS	Parathyroid
4566	121420	T97267	49.03	298.34	6.04	1.00	0.00	15	258.04	Synovial mem	Eye	Bone
4567	280286	N49224	4.76	26.00	5.46	2.00	0.00	9	592.98	Thymus	Codon	Aorta
4571	809694	AA454702	3.87	130.14	33.66	4.00	0.00	1	592.98	Adrenal gland	Pool	LID not found
4573	347036	W81128	5.22	76.40	15.02	10.00	0.00	7	20.09	Pooled	Brain	Uterus
4579	271008	N34362	14.81	91.74	6.20	1.00	0.00	5	357.64	Tonsil	Pool	Uterus
4586	212406	H69471	70.98	408.97	5.78	1.00	3.00	6	129.58	Smooth musc	Testis	Pool
4590	235008	H79130	22.78	282.62	11.54	6.00	5.00	11	221.51	Pancreas	Uterus	Parathyroid
4593	230180	H74330	28.48	153.65	5.39	3.00	0.00	14	249.31	Brain	Testis	LID not found
4605	243969	N45364	9.36	60.96	6.51	1.00	4.00	3	461.79	Lymph	Placenta	Ovary
4606	344141	W69781	35.47	636.95	17.96	6.00	5.00	17	405.59	Umbilical cord	Aorta	Pancreas
4607	810809	AA458882	9.80	83.87	8.47	1.00	0.00	13	288.41	Esophagus	Ovary	CNS
4612	286678	W02256	3.92	31.84	8.12	2.00	0.00	16	411.33	Foreskin	Pancreas	Heart
4626	592243	AA155695	7.89	47.80	6.06	1.00	0.00	1	266.38	Small intestine	Eye	Foreskin
4636	123561	R00822	39.06	387.80	9.42	2.00	3.00	6	16.54	Larynx	Thyroid	Pancreas
4642	382787	AA088598	6.72	151.19	22.31	1.00	0.00	1	740.99	Adipose	Foreskin	Pool
4674	770014	AA427667	5.32	28.30	5.32	1.00	0.00	11	225.28	Head and nec	Thyroid	Blood
4684	130541	R22412	63.99	366.28	5.72	0.00	1.00	12	240.88	Bone	Muscle	Pooled
4700	810761	AA480851	7.32	183.05	22.26	5.00	0.00	19	223.76	Pancreas	Nose	Parathyroid
4708	773301	AA425356	9.09	118.85	12.85	3.00	0.00	20	293.34	Smooth musc	Adrenal gland	Colon
4715	687397	AA235332	54.17	399.84	7.39	0.00	2.00	22	154.71	Ovary	LID not found	Other
4716	638568	AA456931	185.57	1028.72	5.34	2.00	0.00	1	127.63	Slin	Placenta	
4722	842836	AA486275	25.47	144.29	5.66	1.00	0.00	11	165.06	Larynx	Head and nec	Skin
4729	298198	N74383	15.94	91.14	5.72	1.00	0.00	12	227.72	Smooth musc	Testis	-
4734	380831	AA036148	17.76	141.22	7.95	3.00	0.00	2	388.44	Pancreas	Omentum	Nose
4742	796268	AA460827	5.71	32.71	5.73	2.00	0.00	21	215.71	Thyroid	Parathyroid	Pooled
4747	814378	AA459038	46.27	1271.96	27.49	7.00	1.00	6	42.01	Ear	Pool	Adipose
4750	814595	AA480908	31.24	178.65	5.72	2.00	0.00	22	45.84	Germ Cell	Blood	Adrenal gland
4751	724378	AA250771	203.37	1453.55	7.15	4.00	2.00	5	337.8	Gall bladder	Bone	Carvix
4752	207358	H58873	20.78	150.27	7.23	5.00	3.00	2	585.19	Ear	Eye	Stomach
4753	771196	AA443506	105.78	582.39	5.51	1.00	0.00	10	188.11	Pool	Lung	-
4754	592540	AA160507	7.94	93.24	11.74	4.00	1.00	9	359.82	Placenta	LID not found	Other
4755	795771	AA460330	28.48	146.10	5.13	1.00	0.00	7	510.43	Foreskin	Lung	Uterus
4762	810131	AA464250	40.21	384.90	9.67	3.00	0.00	2	304.52	Ear	Tonsil	Blood
4771	785566	AA452868	32.42	256.88	7.92	1.00	0.00	14	67.03	Kidney	Germ Cell	Pool
4772	785975	AA448599	8.61	49.78	5.78	2.00	0.00					
4774	897656	AA598817	12.66	207.32	16.37	5.00	3.00					
4776	199180	R95740	8.81	65.64	7.45	3.00	0.00					
4785	754436	AA410207	13.11	356.01	27.15	6.00	0.00					
4795	39808	R54050	23.99	285.25	11.89	1.00	0.00					
4796	897656	AA498804	12.33	67.84	5.49	1.00	0.00					
4798	781362	AA448400	11.81	104.95	8.81	3.00	0.00					
4803	214185	H17772	42.80	370.44	8.65	4.00	2.00					
4806	136592	R63342	60.99	389.12	6.36	3.00	0.00					
4807	293380	N82046	5.97	63.96	10.90	0.00	1.00					
4808	140057	R65963	2.68	16.25	6.80	2.00	0.00					
4813	298600	N74882	2.54	223.42	87.88	4.00	0.00					
4816	139331	R63782	19.71	205.86	10.45	4.00	4.00					
4820	136289	R33780	38.27	284.91	7.44	4.00	0.00					
4821	122892	R00332	6.77	47.00	6.94	2.00	0.00					

Table 2A

4824	141288	R04408	45.31	1092.18	24.11	4.00	5.00	1	642.27	Placenta	Fore skin	Testis
4835	285492	W23522	258.83	1641.13	6.34	1.00	0.00	19	218.53	CNS	Uterus	Eye
4839	127354	R08359	12.19	102.76	8.43	2.00	2.00	6	115.59	Placenta	Testis	LID not found
4840	141366	R04449	16.79	217.56	12.96	5.00	3.00	6	Parathyroid	Fore skin	Fore skin	Eye
4854	121521	T97794	63.38	451.13	7.12	3.00	0.00	6	381.95	Placenta	LID not found	Other
4856	139656	R63980	2.85	20.17	7.08	4.00	0.00	6	44.4	Fore skin	Spleen	Blood
4858	241066	H01404	57.15	295.74	5.17	0.00	1.00	6	Pool	LID not found	Other	
4859	295594	W02403	183.89	1448.47	7.47	2.00	0.00	5	628.75	Adrenal gland	Ovary	Prostate
4867	121543	T97809	4.93	36.42	7.39	1.00	1.00	5	Neural	Thyroid	Muscle	
4874	357091	W93510	29.88	162.71	5.48	1.00	0.00	7	500.33	Placenta	Germ Cell	Ovary
4881	136303	R34013	58.80	298.31	5.06	2.00	0.00	7	Pool	LID not found	Other	
4882	66474	R16009	9.06	87.27	9.83	5.00	0.00	20	237.59	Pool	LID not found	Other
4884	136632	R34957	23.51	119.77	5.09	1.00	0.00	10	188.13	CNS	Prostate	Colon
4886	38465	R48470	40.64	316.04	7.78	4.00	0.00	3	597.76	Pool	LID not found	Other
4891	121954	T97870	36.96	278.50	7.48	5.00	1.00	13	87.45	Bone marrow	Blood	Uterus
4892	247042	N57848	102.34	727.83	7.11	3.00	0.00	8	327.75	Pool	LID not found	Other
4901	294018	N64033	4.32	21.75	5.03	1.00	0.00	20	200.08	Pool	LID not found	Other
4904	67002	T69708	2.74	15.22	5.55	2.00	0.00	14	112.05	Pool	LID not found	Other
4905	293932	N95842	24.01	275.72	11.49	4.00	5.00	14	98.28	Pool	LID not found	Other
4908	294701	N84840	2.87	15.13	5.68	1.00	0.00	13	Pool	LID not found	Other	
4909	207421	H58866	26.50	182.32	6.12	1.00	0.00	8	Pool	LID not found	Other	
4911	196640	R03009	47.50	298.08	6.28	3.00	2.00	13	87.45	Bone marrow	Blood	Uterus
4915	134712	R28280	5.75	31.21	5.43	1.00	1.00	20	327.75	Pool	LID not found	Other
4917	248481	N78198	3.63	28.12	7.74	3.00	0.00	14	98.28	Pool	LID not found	Other
4918	275950	R63394	9.51	48.64	5.11	1.00	0.00	13	87.45	Bone marrow	Blood	Uterus
4919	202921	H54384	2.99	18.19	6.09	2.00	0.00	8	327.75	Pool	LID not found	Other
4921	200873	R98774	49.91	705.25	14.13	2.00	3.00	20	200.08	Pool	LID not found	Other
4922	286905	N70286	4.42	62.30	14.08	0.00	0.00	14	98.28	Pool	LID not found	Other
4923	135247	R31521	37.87	211.71	5.57	4.00	0.00	14	98.28	Pool	LID not found	Other
4924	416099	W85382	4.56	25.75	5.65	0.00	1.00	13	87.45	Bone marrow	Blood	Uterus
4927	240480	H80746	70.90	676.54	9.54	4.00	0.00	8	327.75	Pool	LID not found	Other
4929	233952	H65942	3.20	50.39	15.74	3.00	0.00	20	200.08	Pool	LID not found	Other
4930	296901	W04272	6.76	53.78	7.95	2.00	0.00	14	98.28	Pool	LID not found	Other
4931	139680	R63998	51.21	285.30	5.57	2.00	1.00	11	593.19	Placenta	Colon	LID not found
4933	204251	H59188	3.12	16.19	5.20	2.00	0.00	1	373.42	Testis	Adrenal gland	Germ Cell
4935	234664	H77737	1.99	10.08	5.07	1.00	0.00	5	473.67	Pool	LID not found	Other
4937	200838	R88947	3.35	37.75	11.26	4.00	0.00	5	473.67	Pool	LID not found	Other
4938	285483	W05028	124.53	3956.30	31.77	12.00	2.00	12	319.85	Brain	Pool	LID not found
4939	127400	R08563	2.94	28.70	10.11	2.00	0.00	12	319.85	Brain	Pool	LID not found
4943	295604	N68845	17.20	142.30	8.27	4.00	0.00	12	319.85	Brain	Pool	LID not found
4945	200840	R88848	2.48	16.72	6.80	1.00	0.00	12	319.85	Brain	Pool	LID not found
4947	141894	R69645	2.25	15.27	6.78	2.00	0.00	12	319.85	Brain	Pool	LID not found
4948	249887	H85454	0.85	6.99	10.76	0.00	1.00	12	319.85	Brain	Pool	LID not found
4951	295873	N73510	7.57	63.84	8.43	3.00	2.00	12	319.85	Brain	Pool	LID not found
4954	294445	W01511	19.57	135.95	6.85	1.00	1.00	12	319.85	Brain	Pool	LID not found
4961	200637	R99004	22.51	434.35	19.30	6.00	5.00	12	319.85	Brain	Pool	LID not found
4964	416587	W86431	4.90	50.28	10.25	1.00	0.00	12	319.85	Brain	Pool	LID not found
4966	486279	AA044205	7.99	316.80	39.67	12.00	1.00	12	319.85	Brain	Pool	LID not found
4967	233627	H65984	128.22	937.80	7.31	2.00	0.00	12	319.85	Brain	Pool	LID not found
4969	247482	N54161	13.46	99.18	7.37	2.00	4.00	12	319.85	Brain	Pool	LID not found
4975	240981	H90899	34.10	259.48	7.61	3.00	0.00	12	319.85	Brain	Pool	LID not found
4976	795877	AA480152	79.76	500.54	6.27	0.00	1.00	12	319.85	Brain	Pool	LID not found
4977	201393	R99627	41.46	515.90	12.44	5.00	5.00	12	319.85	Brain	Pool	LID not found
4978	294942	N71473	169.43	938.84	5.54	2.00	0.00	12	319.85	Brain	Pool	LID not found
4980	341634	W60847	196.87	1153.92	5.86	0.00	1.00	12	319.85	Brain	Pool	LID not found

Table 2A

4881	207293	H56670	4.45	34.40	7.74	1.00	0.00	15	50.93	Pool	LID not found	Other
4984	284479	N52350	2.32	11.77	5.07	1.00	0.00	3	180.11	CNS	Brain	LID not found
4885	203551	H58033	68.70	538.39	7.72	4.00	0.00	2	194			
4888	380029	AA063574	2.32	21.23	8.15	0.00	2.00	7	475.57	Kidney	LID not found	Other
4994	813258	AA455911	3.27	79.35	24.25	5.00	0.00	16	408.29	Small intestine	Esophagus	
4888	251018	H97778	25.48	586.14	23.00	8.00	1.00	3	483.34	Synovial mem	Germ Cell	
5002	442955	H06113	41.35	288.53	6.45	2.00	0.00	17	404.41	Synovial mem	Adrenal gland	
5007	741087	AA478436	145.78	2893.05	19.85	4.00	1.00	1	27	Eye	Spleen	
5011	345586	H76376	3.33	75.85	22.75	2.00	0.00	12	248.83			
5015	160838	H24688	7.91	95.86	12.12	4.00	0.00	1	671.23	Breast	Adipose	Stomach
5016	811162	AA486471	30.70	195.58	6.37	2.00	0.00	19	35.88	Eye	Breast	Muscle
5020	341246	H56858	20.53	107.67	5.25	1.00	0.00	5	153.89	Eye	Lymph	Spleen
5024	752631	AA417854	3.48	73.61	14.99	7.00	0.00	4	24.07	Larynx	Stomach	Pod
5031	296529	H01011	23.25	234.52	10.09	1.00	0.00	18	193.03	Ear	CNS	Germ Cell
5032	808464	AA456160	4.21	90.82	21.59	4.00	0.00	10	628.52	Cervix	Ear	Skin
5048	121722	H98152	13.81	330.56	23.94	0.00	2.00	5	504.31	Placenta	Heart	Whole embryo
5051	769921	AA430504	7.89	156.62	19.47	8.00	0.00	20	254.9	Larynx	Head and nec	Colon
5070	823590	AA487051	7.16	116.02	18.20	9.00	1.00					
5080	813757	AA463816	5.57	90.47	18.24	8.00	0.00					
5084	210317	H65528	8.09	44.17	5.46	0.00	1.00	9	301.16	Mamow	Whole embryo	Pool
5089	857822	AA588572	2.62	301.01	114.81	9.00	1.00					
5090	786872	AA451891	3.12	31.08	9.95	2.00	0.00	22	133.9	Parathyroid	Germ Cell	Whole embryo
5093	67634	T49539	33.32	238.64	7.10	1.00	0.00					
5094	814546	AA480859	3.82	20.52	5.66	2.00	0.00	14	151.92	Pooled	Spleen	Heart
5097	359781	AA011320	6.93	38.09	5.50	1.00	0.00	9	417.73			
5108	233365	H78888	13.96	882.94	7.22	1.00	0.00	13	141.14	Parathyroid	Foreskin	Heart
5111	80500	T64625	273.12	3218.30	11.78	5.00	0.00	8	104.03	Lymph node	Synovial mem	Cervix
5118	788091	AA453105	34.44	280.94	7.58	3.00	0.00	2	587.39	Omentum	CNS	Parathyroid
5124	201886	R99423	18.58	97.87	5.27	1.00	0.00	X	85.4	Uterus	Parathyroid	Cervix
5125	503097	AA151486	58.37	308.18	5.30	1.00	0.00	5	350.47	Cervix	Muscle	Bone
5130	843133	AA486524	151.58	1317.88	8.70	2.00	0.00	7	494.74	Placenta	Tonsil	Ovary
5131	812286	AA455082	21.38	220.91	10.34	2.00	0.00	9	392.05	Tonsil	Muscle	Pod
5138	246789	N53177	6.09	51.98	8.54	4.00	0.00	19	213.91	Smooth musc	Blood	Adipose
5139	824393	AA469714	8.29	35.89	5.71	1.00	0.00	2	304.52	Ear	Tonsil	Blood
5147	824922	AA488017	1.07	14.57	13.67	1.00	0.00	2	222.73	Gall bladder	Tonsil	Foreskin
5164	248649	N58542	9.19	55.80	6.07	1.00	1.00	3	340.31	Eye	Foreskin	Pool
5176	244684	N54296	25.33	136.01	5.37	0.00	1.00	20	284.27	Pool	LID not found	Other
5184	285494	N21309	12.17	95.74	7.87	0.00	1.00	10	441.32	Breast	LID not found	Other
5187	128301	R11529	3.41	17.25	5.08	0.00	1.00					
5188	188390	H43857	5.17	33.01	6.38	1.00	0.00	3	395.51	Placenta	Pooled	Heart
5190	111844	T64889	2.16	12.07	5.58	1.00	0.00	8	540.74	Breast	LID not found	Other
5191	131316	R23055	7.59	47.43	6.25	0.00	2.00	7	607.47	Pool	LID not found	Other
5193	138861	R62780	23.90	197.08	8.24	1.00	0.00	9	123.91	Esophagus	Larynx	Blood
5196	182365	H27590	17.53	128.29	7.20	3.00	0.00	2	216.99	Pool	Placenta	LID not found
5198	206472	H00973	4.01	21.84	5.45	2.00	0.00	4	450.11	Uterus	Placenta	Parathyroid
5200	199327	R65916	101.49	685.32	6.57	2.00	0.00					
5202	109108	T60718	2.79	20.56	7.37	1.00	0.00					
5205	196005	R91904	6.21	41.77	6.73	1.00	0.00					
5207	131388	R23087	3.30	18.58	5.93	2.00	0.00					
5208	210622	H84244	9.85	198.27	20.12	6.00	5.00					
5213	134948	R31831	7.01	54.85	7.82	5.00	0.00					
5214	111765	T84865	51.96	288.27	5.55	1.00	0.00					
5216	186522	R91557	13.24	73.07	5.52	1.00	0.00					
5220	209187	H63668	28.20	687.08	23.68	5.00	3.00					

Table 2A

5222	249143	N52408	187.72	1040.36	5.54	2.00	0.00	19	250.6	Aorta	Uterus
5227	345670	W76603	10.56	413.19	39.13	1.00	0.00	18	66.29	Pooled	Heart
5228	202740	H53803	33.88	486.27	14.45	6.00	3.00	12	484.53	Breast	Colon
5230	121681	T97816	10.69	73.24	6.85	3.00	3.00			LID not found	Other
5232	203772	H56088	53.03	455.36	6.59	4.00	0.00			LID not found	Other
5233	210873	H66312	6.32	2018.02	281.89	1.00	3.00	17	322.57	Pool	Other
5244	201662	R97050	4.34	28.98	6.90	1.00	0.00			LID not found	Other
5248	284968	N98518	3.75	16.88	5.03	0.00	2.00	13	131.78	Pool	Other
5252	181506	H37646	28.55	308.99	10.82	4.00	0.00			LID not found	Other
5255	132159	R26034	9.71	68.16	7.02	5.00	1.00			Pooled	Brain
5262	244154	N75668	70.65	581.16	8.19	1.00	0.00	12	457.31	Pancreas	Testis
5264	185314	R92032	2.34	12.02	5.15	1.00	0.00	13	85.4	Testis	LID not found
5265	138218	R33841	5.13	35.99	7.01	2.00	0.00			Testis	LID not found
5266	239692	H00558	18.28	101.49	5.56	1.00	1.00	7	588.48	Aorta	Pool
5267	128567	R14884	30.46	283.75	9.32	5.00	0.00			LID not found	Other
5272	185820	R92285	46.14	752.28	16.30	6.00	5.00			LID not found	Other
5273	130371	R21785	19.07	133.41	6.99	0.00	1.00			Testis	Placenta
5275	129922	R19183	85.35	775.83	9.08	4.00	2.00			LID not found	Other
5277	368389	AA025807	228.73	1357.48	5.93	3.00	0.00	5	626.75	Ear	Snail intestine
5280	185821	R92292	54.59	680.52	12.47	5.00	5.00	11	124.08	Pool	Other
5288	770838	AA427732	7.49	40.16	5.36	1.00	0.00	5	14.04	Cervix	Pancreas
5289	235055	H73608	14.90	92.63	6.22	0.00	1.00	6	164.21	Spleen	LID not found
5296	811010	AA485365	34.61	212.20	6.15	1.00	0.00	10	342.2	Brain	Lung
5298	133236	R28926	16.89	105.73	6.26	2.00	0.00	12	55.14	Ear	Gall bladder
5305	233446	H77714	26.02	571.38	21.96	4.00	5.00			Muscle	Breast
5308	213871	AA040269	4.08	23.40	9.34	1.00	0.00	8	540.74	Neural	Pool
5310	257011	N26802	12.75	215.41	16.89	5.00	4.00	7	640.87	Pancreas	Ovary
5317	280386	N48375	3.17	18.81	5.84	0.00	1.00			CNS	Muscle
5319	220077	H82532	10.68	56.21	5.26	1.00	0.00	2	545.08	Placenta	Pool
5327	765371	AA453273	31.49	185.46	5.89	1.00	0.00			Brain	LID not found
5328	341854	W58342	2.58	39.90	15.47	1.00	0.00	7	140.47	Mouth	Bone
5329	233719	H79046	1.58	95.64	80.38	1.00	0.00	13	107	Acrida	Whole embryo
5336	111070	T83394	28.74	184.18	6.41	1.00	2.00	17	329.98	Brain	LID not found
5341	243848	N48986	47.92	249.31	5.20	1.00	0.00			Pool	Other
5342	140635	R66845	8.53	50.86	5.97	1.00	0.00	7	413.5	Pool	Other
5344	428773	AA004664	3.90	75.90	19.45	4.00	0.00	9	53.56	Neural	Heart
5345	294150	N89803	47.71	280.39	5.88	1.00	0.00	X	85.88	Placenta	Heart
5348	243656	N48895	166.09	1077.98	5.79	2.00	2.00	4	403.14	Pool	Whole embryo
5350	130835	R22252	11.90	81.49	5.17	0.00	1.00	10	480.31	Parathyroid	LID not found
5353	298748	N74059	20.02	117.13	5.85	2.00	0.00			Uterus	Pool
5354	115223	T86603	37.37	216.81	5.80	1.00	0.00			Foreskin	Placenta
5390	502062	AA126825	4.28	25.03	5.83	0.00	1.00	20	70.67	Pool	Other
5396	125685	R07594	31.68	278.82	6.80	1.00	0.00	9	323.12	Tonsil	Other
5397	811128	AA485683	3.12	42.07	13.50	0.00	2.00	13	155.38	Whole embryo	CNS
5371	247230	N57927	3.06	35.04	11.38	3.00	0.00			Ovary	Spleen
5372	376652	AA045257	12.80	68.37	5.34	0.00	0.00			Colon	Heart
5373	195487	R02163	28.43	177.29	6.71	0.00	1.00	2	300.43	Foreskin	LID not found
5383	812975	AA484805	9.33	55.98	6.00	1.00	0.00	13	84.78	CNS	Heart
5384	150702	H02340	1.26	36.02	29.13	5.00	0.00	9	14.64	Unilateral cord	LID not found
5390	68982	T67549	5.86	36.04	6.15	4.00	0.00	17	338.92	Pooled	Whole embryo
5392	897497	AA497085	3.43	32.94	9.61	1.00	0.00	2	318.41		Blood
5395	813873	AA453831	117.19	672.25	5.74	1.00	0.00	7	115.72	Blood	Kidney
5396	48897	H09814	14.47	158.57	10.96	1.00	0.00	1	582.84	Bone marrow	Synovial membrane
5408	813841	AA453728	41.53	311.01	7.49	0.00	1.00	1	557.85	Smooth muscle	Pancreas

Table 2A

5408	212849	H70473	2.40	20.42	6.50	2.00	0.00	Liver	Pool	Colon
5418	810873	AA459197	31.54	381.89	12.11	4.00	0.00	Esophagus	Ovary	Colon
5428	293104	N91890	29.90	414.75	13.87	1.00	0.00	Gall bladder	Thyroid	Pool
5438	204897	H57180	5.40	32.79	6.07	0.00	1.00	Lymph	Kidney	Pool
5440	211780	H71868	127.85	739.00	5.78	1.00	0.00	Head and nec	Umbilical cord	Thyroid
5443	31842	R41839	6.87	52.68	7.67	3.00	0.00	Adipose	Skin	Umbilical cord
5446	51918	H22563	5.17	39.43	7.63	1.00	1.00	Thyroid	Forsklin	Blood
5454	742082	AA405769	2.11	29.83	14.11	3.00	1.00	Adipose	Liver	Kidney
5455	245199	N76581	5.58	57.99	10.21	2.00	0.00	Adipose	Stomach	Stomach
5456	150623	H02158	15.68	137.58	8.78	5.00	0.00	Trachea	Parathyroid	Thyroid
5458	153008	R48988	2.83	16.41	6.23	1.00	0.00	Adipose	Parathyroid	Brain
5460	27787	R40400	12.29	87.57	7.12	1.00	0.00	Adipose	CNS	Brain
5475	789382	AA453338	2.10	11.87	5.57	1.00	0.00	Adipose	Germ Cell	Breast
5478	49260	H18573	0.87	4.80	5.27	1.00	0.00	Adipose	Breast	Whole embryo
5480	788185	AA453410	9.66	67.53	8.99	0.00	1.00	Adipose	Eye	Testis
5485	51448	H21041	10.99	74.74	6.80	2.00	0.00	Thyroid	Umbilical cord	Spleen
5488	51865	H23187	0.54	3.11	5.77	0.00	1.00	Thyroid	Parathyroid	Parathyroid
5494	838082	AA598776	37.84	315.88	8.33	0.00	1.00	Bone	Kidney	CNS
5500	843049	AA485883	6.11	59.83	8.78	1.00	0.00	Larynx	Lymph	Colon
5509	368511	AA026609	37.73	292.78	7.76	2.00	0.00	Larynx	Muscle	Blood
5510	700302	AA283693	4.06	111.12	27.35	1.00	0.00	Omentum	Breast	Umbilical cord vein
5516	785827	AA461508	5.22	57.32	10.98	1.00	0.00	Tonsil	Brain	Muscle
5524	704459	AA278883	3.39	50.52	14.88	1.00	0.00	Blood	Tonsil	CNS
5533	48799	H14841	105.24	814.30	5.84	2.00	0.00	Eye	Brain	CNS
5538	782313	AA446478	49.15	815.43	11.79	5.00	0.00	Ovary	Stomach	Ear
5543	232873	H75547	229.90	1373.78	5.98	3.00	0.00	Small intestine	Stomach	Stomach
5548	504228	AA132090	4.64	48.76	10.51	4.00	0.00	Thymus	Lymph	Blood
5552	208413	H62162	2.75	77.76	28.27	7.00	1.00	Kidney	Heart	Prostate
5555	469281	AA026112	1.80	10.03	5.56	1.00	0.00	Pancreas	Heart	Uterus
5560	814526	AA459588	4.49	58.09	12.83	7.00	0.00	Thymus	Ovary	Blood
5561	48182	H09614	28.72	182.95	6.37	0.00	1.00	Ear	Forsklin	Breast
5562	784319	AA447098	58.78	302.27	5.32	1.00	0.00	Muscle	Heart	Heart
5563	627939	AA195959	58.12	325.27	5.60	0.00	2.00	Blood	Parathyroid	Lymph
5564	306013	N91385	3.75	27.58	7.36	1.00	0.00	Blood	Parathyroid	Stomach
5567	825583	AA504817	7.99	40.21	5.03	1.00	0.00	Head and nec	Colon	Pancreas
5568	785283	AA478543	22.82	120.51	5.28	0.00	1.00	Pooled	Aorta	Whole embryo
5570	248703	N59716	20.00	114.85	6.75	1.00	0.00	Thyroid	Small intestine	Thymus
5578	121981	T87889	6.70	49.29	7.35	0.00	1.00	Thyroid	LID not found	Other
5583	124091	R02718	2.45	13.82	5.65	1.00	0.00	Pool	Pool	LID not found
5590	201764	R99311	5.55	61.58	11.10	5.00	0.00	Liver	Whole embryo	Prostate
5598	138775	R35253	11.61	60.78	5.23	1.00	0.00	Placenta	Eye	CNS
5599	124079	R02820	13.57	72.53	5.35	0.00	1.00	Thyroid	Tonsil	Ovary
5605	108422	T77647	15.84	161.17	11.44	3.00	0.00	CNS	Parathyroid	Forsklin
5610	294040	N68497	29.72	233.68	7.88	1.00	3.00	Pancreas	Muscle	Tonsil
5615	128785	R10015	8.08	66.55	8.23	2.00	0.00	Ear	Small intestine	Thymus
5617	245428	N77205	81.49	492.14	6.04	1.00	0.00	Ear	LID not found	Other
5619	121736	T86075	18.48	215.95	11.69	8.00	2.00	Pool	Brain	LID not found
5629	22711	T74714	6.74	37.14	5.51	1.00	0.00	CNS	Umbilical cord	Whole embryo
5630	243317	N45248	36.59	245.26	6.35	4.00	0.00	CNS	Placenta	Other
5632	140334	R68994	7.95	53.74	6.76	3.00	3.00	Placenta	LID not found	Kidney
5633	122178	T98615	11.93	122.49	10.26	6.00	0.00	Placenta	Placenta	Heart
5644	138501	R38181	85.42	729.64	8.54	3.00	0.00	Placenta	Breast	Heart
5646	292749	N80491	14.30	72.90	5.10	1.00	0.00	Blood	Germ Cell	Pool
5647	124822	R05837	2.25	17.82	7.93	2.00	0.00	Ovary	Brain	Pool



Table 2A

5653	121776	T98098	3.24	17.87	5.45	1.00	0.00	3	46.9	Blood	Placenta	Prostate
5657	127636	R08301	11.75	77.84	6.82	1.00	1.00		Pool	Pool	LID not found	Other
5660	138496	R68634	8.36	62.36	7.46	1.00	1.00	19	284.08			
5665	200545	H48445	7.60	42.69	5.82	1.00	0.00		Pool		LID not found	Other
5673	200904	H48407	17.21	131.41	7.84	3.00	0.00	15	243	Pool	LID not found	Other
5674	295044	N71585	10.93	98.74	8.04	5.00	3.00		Pool		LID not found	Other
5677	207881	H80317	49.61	685.43	13.82	4.00	2.00		Pool		LID not found	Other
5678	137139	R38006	7.83	50.87	6.47	2.00	1.00	6	364.54	Tonsil	Agria	Whole embryo
5682	295108	W01645	3.52	44.23	12.55	8.00	0.00	11	236.72	Pool	LID not found	Other
5685	207068	H60491	16.17	132.68	8.74	6.00	0.00	2	544.88	Pool	LID not found	Other
5687	138444	R88272	3.68	61.83	16.80	1.00	0.00	10	345.66	Eye	CNS	Placenta
5690	295590	W02401	44.35	282.25	6.36	2.00	1.00	20	336.98	Pool	LID not found	Other
5693	207952	H60523	19.24	145.45	7.56	2.00	0.00		Pool		LID not found	Other
5696	271952	N35301	5.53	91.00	18.46	5.00	2.00	15	220.59	Pool	Gall bladder	Forebrain
5697	201228	R98386	10.11	180.28	17.83	5.00	5.00		Pancreas		LID not found	Other
5698	288094	W02591	101.25	656.99	8.49	2.00	0.00		Pool		LID not found	Other
5703	159725	H23983	5.01	56.47	11.26	1.00	1.00	2	508.5	Breast	LID not found	Other
5709	208984	H60868	74.43	471.73	8.34	1.00	0.00		Adrenal gland	Pool	LID not found	Other
5713	201314	R99682	32.16	202.66	8.30	2.00	1.00		Pool		LID not found	Other
5715	144880	R78580	1.92	12.53	6.52	2.00	0.00	19	35.68	Pool	Placenta	Pool
5719	131239	R24258	3.51	25.05	7.13	2.00	0.00	1	15.07			
5720	810403	AA464202	15.02	82.07	5.46	0.00	1.00	11	227.6	Parathyroid	Ovary	Breast
5721	201317	R99680	13.71	206.53	15.08	6.00	5.00		Pool		LID not found	Other
5722	268559	W00793	73.17	507.83	8.94	1.00	0.00		Pool		LID not found	Other
5724	281476	N15121	10.49	385.12	36.70	1.00	0.00	4	673.59			
5725	210548	H65052	29.39	174.78	6.95	6.00	4.00	11	167.19	Liver	Spleen	Blood
5726	109314	T60848	3.36	19.20	5.71	1.00	0.00	19	32.22	Synovial mem	Adrenal gland	Placenta
5728	380075	AA013240	2.11	12.37	5.88	1.00	0.00	5	402.9	Germ Cell	Eye	Brain
5731	135888	R32408	5.83	32.43	5.47	1.00	0.00	3	453.05	Placenta	Germ Cell	Uterus
5733	213118	H69576	128.72	862.04	6.70	1.00	1.00		Pool		LID not found	Other
5737	201784	R99938	11.19	64.31	5.75	1.00	0.00		Pool		LID not found	Other
5738	286802	W01028	68.30	436.73	6.59	2.00	1.00	7	64.01	Placenta	Aorta	Pool
5739	240430	H78087	14.99	99.87	6.88	1.00	0.00		Pool		LID not found	Other
5741	208840	H61684	6.56	44.63	6.80	4.00	0.00		Pool		LID not found	Other
5747	295918	W04152	51.22	527.18	10.29	6.00	2.00		Placenta	Tonsil		Pool
5748	795883	AA460168	0.78	12.86	16.97	1.00	1.00	19	287.41	Smooth musc	Pancreas	Skin
5749	211351	H75690	2.98	20.28	8.81	1.00	0.00	2	88.08	Pool	LID not found	Other
5750	122787	T99674	0.88	37.21	6.14	1.00	0.00		Pool		LID not found	Other
5753	201818	H48318	2.64	18.03	6.82	1.00	4.00		Pool		LID not found	Other
5755	284592	N76193	20.72	132.39	6.39	0.00	1.00	11	204.26	Placenta	Pool	Cervix
5757	208434	H62166	20.58	162.07	7.87	2.00	2.00		Pool		LID not found	Other
5759	233645	H78007	76.18	550.08	7.22	2.00	2.00		Pool		LID not found	Other
5783	28012	R40784	25.40	139.90	5.51	1.00	0.00	X	238.33	Small intestine	Gall bladder	Parathyroid
5787	23932	R39484	4.33	63.77	14.73	1.00	0.00	12	457.31	Pancreas	CNS	Testis
5788	810813	AA45684	10.58	303.83	28.72	8.00	4.00		Larynx		Nose	Head and neck
5790	303048	N91584	2048.08	12924.63	6.32	1.00	0.00	9	70.89	Prostate	Lung	LID not found
5791	46358	H08938	10.99	89.80	8.17	1.00	0.00		Lung		Brain	LID not found
5794	122774	T98783	34.32	203.01	5.91	1.00	1.00	2	108.91			
5796	177737	H46683	0.79	39.24	5.78	1.00	0.00		Blood		Prostate	Muscle
5804	742101	AA405891	14.67	288.69	18.68	10.00	3.00	2	412.17	Thyroid	Skin	Ovary
5812	810724	AA480815	1.15	13.72	11.95	1.00	1.00	6	116.56	Marrow	Smooth musc	Larynx
5820	511068	AA100296	1.89	15.26	8.08	3.00	0.00		Ear		Umbilical cord	Tonsil
5824	46508	H15634	13.41	105.49	7.87	3.00	0.00	19	71.09	Whole embryo	Pool	Kidney
5826	51447	H20872	7.47	119.82	16.02	6.00	0.00	1	576.3	Peripheral ner	CNS	Blood
5827	767069	AA424516	8.20	240.83	26.18	16.00	0.00	11	400.33	Ear	Brain	Colon

Table 2A

5831	769028	AA426311	5.01	151.45	30.25	0.00	1.00	17	320.28	Whole embryo	Pool	Adipose
5839	809598	AA458472	38.49	297.75	8.16	6.00	0.00	6	118.71	Thymus	Skin	Brain
5847	23173	T75436	7.87	39.39	5.14	1.00	0.00	4	457.01	Adrenal gland	Eye	Cervix
5856	137017	R35665	0.77	10.38	13.44	0.00	1.00	7	252.33	Placenta	Tonsil	Thyroid
5857	897906	AA598652	8.76	73.88	8.43	8.00	0.00	3	701.95	Kidney	Aorta	Forebrain
5859	667883	AA258398	38.68	259.23	6.70	2.00	1.00	12	311.24	Forebrain	Umbilical cord	Esophagus
5861	842714	AA460305	673.76	5862.61	8.70	2.00	0.00	12	384.44	Smooth muscle	Gall bladder	Eye
5866	300037	AA048411	60.04	303.58	5.06	1.00	0.00	17	382.93	Parathyroid	Kidney	Breast
5877	430318	AA010609	48.12	232.69	5.04	1.00	0.00	22	119.23	Parathyroid	Thyroid	Lymph
5878	815239	AA481277	3.85	28.18	7.31	1.00	0.00	19	239.99	Skin	Nose	Cervix
5883	828135	AA521346	8.15	41.31	5.07	1.00	0.00	6	132.5	Stomach	LID not found	Other
5884	209518	H65260	1.69	8.88	5.25	1.00	0.00	X	351.05	Blood	Pancreas	Bone
5891	810942	AA459380	42.16	222.37	5.27	1.00	0.00	2	223.55	Neural	Tonsil	Germ cell
5892	127408	R08755	31.42	167.01	5.32	3.00	1.00	22	51.16	Larynx	Pooled	Forebrain
5901	827144	AA521243	8.79	135.51	15.41	3.00	0.00	X	138.23	CNS	Forebrain	Lymph
5903	725503	AA292995	24.28	144.22	5.95	2.00	0.00	17	345.1	Placenta	Pool	Heart
5906	49117	H14804	7.64	40.03	5.31	1.00	0.00	2	601.35	Nose	Gall bladder	Umbilical cord
5911	198444	R82485	6.04	35.89	5.94	0.00	1.00	11	257.7	Lymph node	Tonsil	Blood
5913	840691	AA486387	54.57	459.40	8.42	1.00	0.00	1	292.92	Cervix	LID not found	Other
5921	850607	AA608548	216.08	1241.33	5.74	2.00	0.00	22	114.61	Small intestine	Pooled	Thyroid
5923	815284	AA481547	9.15	58.89	6.43	1.00	0.00			CNS	Lung	Ovary
5927	195753	R89083	132.07	999.77	7.57	2.00	0.00			Breast	Pool	LID not found
5934	205445	H57830	18.29	123.48	6.75	0.00	1.00			Pooled	Placenta	Pool
5942	811600	AA458533	3.45	71.40	20.87	5.00	0.00			Brain	Pool	LID not found
5954	109271	T81281	69.82	507.82	7.27	3.00	1.00			Testis	Uterus	Pool
5982	132630	R25980	6.16	33.31	5.41	1.00	0.00	2	545.1	Forebrain	Heart	Whole embryo
5984	191518	H36148	7.53	94.44	12.55	6.00	0.00			Pool	LID not found	Other
5965	184351	H60871	75.77	481.07	6.35	2.00	0.00	10	338.83	Placenta	LID not found	Other
5968	195853	R82310	13.98	253.22	18.12	6.00	5.00	10	508.28	Pool	LID not found	Other
5970	122899	T89881	30.71	201.82	6.37	1.00	2.00	15	243.99	Placenta	Pooled	Tonsil
5972	242700	H84183	18.35	157.40	9.62	3.00	0.00	1	582.98	Dvary	Colon	Tonsil
5975	131446	R23952	3.77	25.27	6.70	2.00	0.00			Placenta	CNS	Pancreas
5976	244781	N54407	13.56	256.46	19.06	4.00	5.00	10	536.56	Liver	Eye	Pool
5983	134235	R31154	68.98	542.83	8.10	5.00	0.00	4	470.33	Muscle	Thyroid	Stomach
5984	198125	R92347	24.08	233.28	9.89	6.00	0.00	3	430.63	Pool	LID not found	Other
6000	186303	R62435	2.85	16.82	5.94	1.00	0.00	10	536.56	Liver	Eye	Pool
6007	143322	R74357	4.83	86.93	17.99	9.00	0.00	4	470.33	Muscle	Thyroid	Stomach
6008	198350	R92545	118.58	1124.58	9.40	3.00	0.00	3	430.63	Pool	LID not found	Other
6014	202209	H52534	12.93	79.88	6.18	2.00	0.00	1	27			
6016	198345	R82455	121.53	1332.25	10.98	3.00	0.00	X	85.38	Pool	LID not found	Other
6018	243770	N39325	80.01	787.20	8.75	1.00	0.00	8	247			
6021	66815	T64956	18.90	121.28	6.42	2.00	0.00	21	144.83	Gall bladder	Eye	Parathyroid
6023	136317	R34121	36.44	325.85	8.94	5.00	0.00	6	104.03	Lymph node	Prostate	Uterus
6027	129342	R16431	16.11	138.54	10.25	4.00	0.00			Pool	Tonsil	Colon
6028	194401	R3017	104.75	1073.88	10.25	4.00	0.00	6	463.92	Pool	LID not found	Other
6031	198185	R91948	5.99	28.94	5.00	1.00	0.00	4	678.51	Pool	Whole embryo	Parathyroid
6032	243784	N33927	22.22	210.26	9.46	5.00	5.00	8	477.99	Pool	LID not found	Other
6034	120097	T95151	4.17	53.48	12.87	5.00	0.00			CNS	Brain	Eye
6035	293683	N94181	55.73	758.70	13.61	4.00	0.00				Aorta	Whole embryo
6036	204570	H63223	3.12	32.78	10.50	2.00	0.00					
6037	66697	T64881	5.14	25.73	5.01	1.00	0.00					
6038	242070	H94238	5.01	41.59	8.30	1.00	3.00					
6039	233318	H76863	3.73	34.93	9.37	1.00	0.00					
6040	290054	N64671	7.33	37.08	5.08	1.00	0.00					
6042	230247	H94934	11.15	82.07	7.36	0.00	3.00					

Table 2A

6043	203514	H55968	3.76	23.26	6.19	2.00	0.00	11	288.08	Gall bladder	Fore skin	Whole embryo
6044	295741	W02483	2.20	15.80	7.19	1.00	0.00	17	406.09	Parathyroid	Ovary	Brain
6047	131824	R25114	2.48	12.87	5.18	0.00	1.00	11	425.84	Placenta	Heart	Lung
6048	280122	N49231	59.47	1328.43	22.30	6.00	5.00					
6052	282404	N49774	8.11	58.31	6.40	2.00	0.00	12	134.37	CNS	Breast	Kidney
6055	782203	AA431972	13.54	92.37	5.95	0.00	1.00					
6057	240138	H79613	2.84	16.53	5.89	3.00	0.00					
6058	121811	T97580	78.78	494.97	6.28	1.00	0.00	16	162.07	Aorta	Stomach	Pooled
6063	307314	N85217	202.10	1072.03	5.30	2.00	0.00			Fore skin	Esophagus	Cervix
6064	241847	H83393	8.06	92.05	11.42	3.00	0.00			Ignote	Pool	LID not found
6067	416856	W06660	2.87	15.83	5.51	0.00	1.00	21	215.71	Smooth muscle	Pooled	Prostate
6069	284341	N52254	3.17	45.78	14.43	0.00	1.00	15	243.11			
6071	203474	H55784	140.53	977.85	6.86	1.00	1.00	7	94.41	Fore skin	LID not found	Other
6075	286000	N23753	18.46	106.59	5.77	2.00	0.00	2	394.55	Thymus	Breast	Whole embryo
6083	810203	H25846	60.84	644.48	10.83	5.00	0.00			Ovary	LID not found	Other
6085	244722	N52535	35.31	321.13	9.10	4.00	0.00	2	558.69	Pool	LID not found	Other
6086	470646	AA031770	77.49	420.72	5.43	0.00	1.00	X	275.71	Ovary	Uterus	Tonsil
6089	241087	H80336	16.11	118.61	7.36	2.00	3.00			Pool	LID not found	Other
6092	792537	AA448484	4.58	44.76	9.76	3.00	0.00			Prostate	Germ Cell	Testis
6095	417305	W80001	11.16	64.46	5.77	4.00	0.00	6	64.16	Germ Cell	Pool	LID not found
6097	295992	N73555	6.18	37.53	6.07	1.00	0.00			Fore skin	Uterus	Cervix
6099	416611	W88466	5.53	81.61	16.57	3.00	0.00			Pancreas	Pool	LID not found
6101	245082	N72384	6.58	33.87	5.15	0.00	1.00			Prostate	Pool	LID not found
6103	771023	AA427878	4.85	39.05	6.39	7.00	4.00	1	28.08	Lung	Ovary	LID not found
6110	233289	H77508	98.89	843.44	6.55	2.00	0.00	19	271.02	Pool	LID not found	Other
6113	241507	H80724	74.60	423.77	5.68	2.00	0.00			Adipose	Bone	Ear
6115	488202	AA046424	4.16	30.84	7.36	1.00	0.00			Larynx	Kidney	Breast
6116	770866	AA434390	12.25	218.85	17.86	1.00	1.00	9	377.24	Thyroid	Stomach	Ovary
6118	810809	AA464739	120.65	838.42	6.93	1.00	1.00			Pool	LID not found	Other
6122	129616	R16656	10.16	67.75	6.67	3.00	0.00			Pool	Ovary	Germ Cell
6124	771058	AA427521	4.82	28.32	6.09	1.00	0.00			Pool	LID not found	Other
6125	246478	N73227	27.56	247.28	8.97	3.00	0.00			Cervix	Pool	LID not found
6127	195370	R89581	6.39	43.67	6.84	6.00	0.00	15	243.19	Pool	LID not found	Other
6129	230509	H81048	151.03	1118.42	7.41	2.00	1.00	9	252.77	Liver	Ear	Skin
6130	233071	H75832	5.96	37.27	8.28	1.00	0.00			Pooled	Thyroid	Tonsil
6131	212438	H65528	58.01	551.59	8.51	5.00	0.00			Pool	LID not found	Other
6133	246786	N53167	1.87	14.96	6.02	1.00	0.00			Pool	LID not found	Other
6136	418289	W80749	3.02	36.73	12.17	1.00	0.00	14	260.52	Fore skin	Pool	Brain
6137	244201	N52978	2.38	42.67	18.17	2.00	0.00	12	53.22	Pool	LID not found	Other
6138	134312	R31218	3.40	33.60	9.87	1.00	0.00	14	260.52	Placenta	LID not found	Other
6140	66856	T67223	7.54	55.44	7.36	2.00	5.00	1	150.31		Testis	Lung
6143	132217	R26396	5.46	43.00	7.87	0.00	2.00			Placenta	Eye	Placenta
6144	140267	R87803	151.34	1382.44	9.00	2.00	3.00			Nose	Fore skin	LID not found
6148	266161	N21592	11.28	80.12	7.98	7.00	0.00	9	357.99	Nose	Parathyroid	Pooled
6156	322223	W38022	3.13	60.77	18.44	16.00	0.00	12	415.54	Fore skin	Heart	LID not found
6160	376943	AA046112	8.40	60.79	6.05	0.00	1.00			Fore skin	LID not found	Other
6168	272548	N35892	7.41	70.52	9.52	3.00	0.00					
6171	133192	R28456	2.56	15.47	6.03	1.00	0.00	12	242.02	CNS	Heart	Brain
6191	357285	W83682	3.75	57.93	15.45	2.00	0.00			Pool	LID not found	Other
6194	121649	T97650	2.89	14.82	5.05	1.00	0.00	11	99.65	Placenta	Pool	Kidney
6203	357396	W93847	2.98	45.29	15.22	1.00	1.00	19	34.68	Marrow	Breast	Ovary
6211	782233	AA431721	8.53	48.58	5.70	2.00	0.00	11	277.15	Blood	Fore skin	Colon
6215	251712	N67839	24.87	153.46	6.17	2.00	0.00	6	503.11	Placenta	Parathyroid	Prostate
6219	428431	AA004415	1.46	14.81	10.00	2.00	0.00					

Table 2A

6220	213850	H72388	20.72	218.92	10.56	0.00	1.00	11	131.77	Whole embryo	Ovary	Kidney
6236	112456	T90971	3.27	17.94	5.48	1.00	0.00	10	400.44	Spleen	Tonsil	Fore skin
6240	301043	N81017	2.99	24.86	8.31	5.00	0.00	2	743.9	Heart	Lung	Testis
6244	278666	H68542	23.57	137.85	5.82	1.00	0.00	19	68.86	Bone marrow	Head and nec	Pituitary
6247	327038	AA264267	2.66	17.03	5.96	1.00	0.00	10	549.88	Testis	LID not found	Other
6248	297731	N65908	44.42	830.09	18.68	3.00	0.00	1	32.1	Pool	LID not found	Other
6251	809703	AA454710	4.39	33.40	7.61	2.00	0.00	18	115.51	Heart	LID not found	Other
6254	428764	AA009677	2.47	22.30	9.01	4.00	0.00	9	417.73	Ear	Heart	Brain
6259	154769	R55630	6.94	49.19	7.09	1.00	0.00	12	84.71	Testis	Pool	LID not found
6260	203114	H54419	104.49	1043.38	9.89	4.00	0.00	1	165.69	Heart	Testis	LID not found
6263	327239	AA284307	11.58	84.71	5.59	1.00	0.00	2	714.07	Prostate	Heart	LID not found
6265	365590	AA028666	27.63	139.08	5.03	1.00	0.00	5	154.93	Head and nec	Stomach	Umbilical cord
6268	203554	H56028	213.57	2457.62	11.51	4.00	0.00	5	633.32	Lung	Heart	Ovary
6275	234131	H70808	252.22	1518.18	8.02	1.00	0.00	7	94.72	Cervix	Umbilical cord	Pooled
6277	327179	W02753	4.39	40.09	8.14	1.00	0.00	7	110.66	Pooled	Heart	Kidney
6285	327195	AA284249	15.86	91.15	5.75	1.00	0.00	1	165.69	Heart	Testis	LID not found
6292	245542	N53564	25.65	140.20	5.47	1.00	1.00	2	487.75	Prostate	Heart	LID not found
6298	320495	W16059	3.54	31.61	8.93	1.00	0.00	2	714.07	Ear	Umbilical cord	Aorta
6300	208389	H64150	187.99	1297.51	6.80	2.00	3.00	2	487.75	Prostate	Heart	LID not found
6303	321470	W23303	11.36	57.93	5.10	1.00	0.00	5	154.93	Head and nec	Stomach	Umbilical cord
6304	590148	AA156030	13.06	70.29	5.38	1.00	0.00	5	633.32	Lung	Heart	Ovary
6305	384583	AA022949	4.71	164.27	34.90	7.00	3.00	5	154.93	Head and nec	Stomach	Umbilical cord
6308	505000	AA151285	3.90	28.85	7.65	3.00	0.00	7	94.72	Cervix	Umbilical cord	Pooled
6316	479494	AA011347	11.40	61.38	5.39	1.00	0.00	7	110.66	Pooled	Heart	Kidney
6318	343695	W69170	4.94	44.88	8.11	1.00	0.00	1	145.6	Thyroid	Liver	Adrenal gland
6320	232274	H68518	8.08	56.97	7.05	1.00	0.00	10	384.98	Head and nec	Uterus	Brain
6321	501868	AA128005	19.13	114.62	5.99	1.00	0.00	14	280.52	Nose	Stomach	Lymph
6324	207636	H59083	35.97	228.38	8.35	2.00	3.00	11	18.88	Nose	Skin	Adipose
6326	346360	W74254	3.42	20.51	5.99	1.00	0.00	1	145.6	Thyroid	Liver	Adrenal gland
6328	488145	AA058709	4.34	32.31	7.45	3.00	0.00	19	194.27	Cervix	Colon	Tonsil
6331	502860	AA128382	21.80	150.03	6.88	1.00	0.00	11	268.99	Smooth muscle	Spleen	Lymph
6341	744960	AA625915	2.38	17.53	7.36	1.00	0.00	1	708.84	Adipose	Cervix	Parathyroid
6346	530545	AA112978	8.24	41.48	5.04	1.00	0.00	1	93.22	Penpheral net	Bone	Germ Cell
6351	509641	AA058323	478.17	3249.58	6.80	2.00	0.00	12	427.01	Liver	Esophagus	Ovary
6352	85384	T71985	52.99	317.62	5.99	1.00	0.00	3	156.88	Umbilical cord	Brain	Testis
6353	435858	AA701545	6.04	34.49	5.71	2.00	0.00	20	229.43	Muscle	CNS	Brain
6354	1031799	AA609655	15.85	93.01	5.67	1.00	0.00	8	425.57	Ear	Kidney	Germ Cell
6355	455121	AA676804	16.59	116.35	7.01	0.00	1.00	11	276.98	Blood	Fore skin	Prostate
6359	236034	H61243	21.79	621.89	28.53	18.00	0.00	11	221.51	Spleen	Kidney	Tonsil
6360	79000	T61938	10.97	193.22	17.61	4.00	0.00	3	162.83	Stomach	Colon	Prostate
6369	744800	AA644448	7.62	75.32	9.89	5.00	0.00	12	427.01	Liver	Esophagus	Ovary
6377	877782	AA626787	8.34	47.69	5.72	1.00	0.00	12	427.01	Liver	Esophagus	Ovary
6385	725308	AA281556	2.25	12.42	5.52	2.00	0.00	3	156.88	Umbilical cord	Brain	Testis
6394	32687	R43605	2.94	40.75	13.86	1.00	0.00	2	357.89	Lymph node	Nose	Gall bladder
6397	743880	AA534464	9.24	58.62	6.35	1.00	0.00	17	475.65	Eye	Stomach	Breast
6405	740742	AA479891	23.06	118.87	5.15	0.00	1.00	4	427.9	Tonsil	Eye	Pituitary
6407	731426	AA412064	6.23	34.41	5.52	1.00	0.00	20	229.43	Muscle	CNS	Brain
6408	79726	T62552	6.67	131.55	19.71	8.00	0.00	8	425.57	Ear	Kidney	Germ Cell
6410	742763	AA400186	10.99	58.95	5.38	1.00	0.00	11	276.98	Blood	Fore skin	Prostate
6412	67330	T49238	15.28	78.34	5.12	1.00	0.00	11	221.51	Spleen	Kidney	Tonsil
6418	41647	R52784	3.33	35.58	10.70	0.00	1.00	3	162.83	Stomach	Colon	Prostate
6419	252663	H68329	8.76	210.92	21.61	2.00	1.00	12	427.01	Liver	Esophagus	Ovary
6420	67385	T49309	39.49	225.30	5.71	0.00	1.00	4	427.9	Tonsil	Eye	Pituitary
6424	78576	T62849	12.81	77.90	6.08	1.00	0.00	8	425.57	Ear	Kidney	Germ Cell
6425	809846	AA454813	28.02	164.76	6.59	0.00	1.00	11	276.98	Blood	Fore skin	Prostate

Table 2A

6426	365175	AA074222	28.46	168.58	5.92	0.00	1.00	240.06	Pancreas	Aorta	Germ Cell
6427	858535	AA635577	11.50	155.33	13.51	0.00	2.00	194.42	Cervix	Liver	
6440	34070	R44850	66.21	392.88	5.93	1.00	2.00	19.4	Brain	LID not found	Other
6443	784174	AA432108	72.21	417.09	5.78	0.00	1.00	347.25	Ear	Stomach	Umbilical cord
6445	50130	H16989	2.33	13.13	5.63	0.00	1.00	223.28	Brain	LID not found	Other
6450	246430	N33031	4.14	80.65	19.48	11.00	1.00	419.22	-	Pool	LID not found
6452	509089	AA148213	24.25	160.39	6.61	0.00	1.00	542.75	Umbilical cord	Blood	Kidney
6455	270559	N33258	5.01	37.12	7.41	1.00	0.00	96.5	Gall bladder	Stomach	Adrenal gland
6456	20064	R44955	137.13	888.45	6.48	1.00	2.00	291.03	Brain	Heart	LID not found
6464	33817	R44717	2.59	14.14	5.45	1.00	0.00	Brain	Breast	Lung	
6469	49918	H15286	202.82	1149.89	5.67	1.00	1.00	323.38	Breast	Prostate	Pool
6479	592594	AA159578	39.38	483.21	12.27	4.00	0.00	195.75	Pancreas	Uterus	Heart
6480	32598	R43721	4.87	29.81	6.12	0.00	1.00	644.77	Brain	LID not found	Other
6482	858187	AA630628	186.23	861.50	5.18	0.00	2.00	121.59	Cervix	Lymph	Adrenal gland
6490	743230	AA400234	77.96	616.12	7.90	5.00	0.00	237.93	Trachea	Parathyroid	Thyroid
6497	49704	H15886	10.03	83.28	8.30	0.00	2.00	387.03	CNS	Blood	Eye
6498	431803	AA678021	170.86	1007.34	5.89	1.00	0.00	674.5	Gall bladder	Bone	Muscle
6506	840753	AA486072	1502.29	10587.00	7.65	2.00	0.00	271.02	-	-	-
6509	45852	H08862	5.19	162.57	31.34	0.00	1.00	839.73	CNS	Adrenal gland	Eye
6521	47358	H10983	1.09	21.25	19.46	2.00	0.00	524.72	Brain	LID not found	Other
6523	284908	N68750	30.89	190.53	6.17	0.00	1.00	477.69	CNS	Thyroid	Skin
6524	291633	N73448	113.94	920.08	8.08	0.00	2.00	556.38	Forebrain	LID not found	Other
6525	47580	H11718	62.53	433.10	6.93	2.00	0.00	240.78	Brain	Bone	Brain
6528	39973	R52522	94.50	764.45	8.09	2.00	1.00	-6.83	Forebrain	LID not found	Other
6532	484271	AA088214	41.62	309.91	7.45	2.00	1.00	317.39	Breast	Uterus	Lung
6546	418651	W87714	35.68	230.40	6.46	2.00	0.00	840.65	CNS	Forebrain	Tonsil
6554	280000	N38891	8.63	50.49	7.82	3.00	0.00	-	-	Prostate	-
6562	128515	R10875	8.07	124.83	20.56	4.00	0.00	488.94	Mouth	Uterus	-
6570	414894	W83087	39.07	269.43	6.90	1.00	0.00	86.18	Ignore	Ear	Bone
6575	347320	W80688	4.39	47.45	10.81	3.00	0.00	221.61	-	Smooth muscle	Colon
6578	487151	AA043790	33.60	363.84	10.77	3.00	4.00	93.95	Breast	Placenta	Testis
6591	189186	R65841	254.66	1824.02	7.16	3.00	0.00	46.83	Tonsil	LID not found	Other
6598	429210	AA007276	4.28	27.01	6.34	1.00	0.00	274.89	-	Heart	Brain
6599	365955	AA063588	187.30	1123.88	6.72	2.00	0.00	575.4	Esophagus	Forebrain	Forebrain
6602	271830	N35156	38.38	239.91	6.09	1.00	0.00	77.44	Ovary	Adipose	Gall bladder
6604	324148	W46629	60.99	318.27	5.22	0.00	1.00	476.7	Gall bladder	Colon	Tonsil
6608	292082	N73308	160.67	1020.03	6.35	0.00	0.00	117.84	Small intestine	Bone	Muscle
6622	809585	AA458629	14.13	93.33	6.61	2.00	0.00	227.72	Adrenal gland	LID not found	Other
6629	487371	AA046700	12.41	140.16	11.29	1.00	0.00	204.36	Pool	Thymus	Forebrain
6633	429932	AA033991	2.58	14.77	5.73	1.00	0.00	118.71	Thymus	Skin	Adipose
6640	811139	AA485739	109.53	2237.65	20.43	3.00	0.00	-	-	Heart	Whole embryo
6644	755689	AA486438	10.58	53.68	5.07	1.00	0.00	650.68	Stomach	Blood	Parathyroid
6645	429927	AA034058	1.29	372.24	288.89	1.00	0.00	107.16	Omentum	CNS	Eye
6648	854444	AA669055	20.16	191.50	9.50	6.00	0.00	287.98	Small intestine	Uterus	Placenta
6672	72395	T51539	11.69	89.65	7.67	1.00	0.00	526.2	Placenta	Tonsil	Parathyroid
6673	365425	AA025248	117.03	728.09	6.22	2.00	0.00	-	-	Pool	LID not found
6674	302997	N91145	21.72	126.14	5.81	1.00	0.00	32.75	Smooth muscle	Synovial mem	Blood
6678	381323	AA017544	4.99	119.62	23.95	14.00	0.00	-	-	Heart	-
6680	611443	AA176581	11.86	111.11	9.37	4.00	0.00	-	-	Eye	-
6683	400718	AA115761	2.35	46.54	19.79	2.00	0.00	-	-	Muscle	-
6684	745343	AA625855	3.55	70.07	19.71	2.00	0.00	-	-	Pooled	-
6691	306420	N92699	3.99	32.06	8.04	2.00	0.00	-	-	Uterus	-
6693	490970	AA138666	2.74	37.03	13.49	1.00	0.00	-	-	Small intestine	-
6694	755358	AA453485	4.85	41.04	6.46	1.00	0.00	-	-	Adipose	-
6698	273546	N33274	83.64	501.27	5.99	0.00	2.00	-	-	Uterus	-

Table 2A

6698	269433	N26175	112.45	735.87	6.54	2.00	0.00	8	440.23	Muscle	Foreskin	Heart
6712	447588	AA702422	51.96	309.09	5.95	1.00	2.00	12	397.57	Uterus	Lymph	Heart
6718	771241	AA443587	186.58	1484.29	8.89	2.00	2.00	15	230.36	Liver	Pooled Lymph	Lymph
6720	788256	AA454098	27.04	167.65	6.20	0.00	1.00	4	24.02	Brain	LID not found	Other
6721	323500	W45688	295.35	1983.30	6.72	1.00	0.00	4	24.02	Small intestine	Lymph	Tonsil
6728	34098	R46930	9.85	88.35	9.16	2.00	3.00	4	24.02	Brain	Lymph	Tonsil
6730	796278	AA460838	11.29	83.60	7.41	0.00	2.00	4	-6.08	Brain	Whole embryo	LID not found
6731	841396	AA487543	6.54	48.98	7.49	4.00	0.00	4	-6.08	Spleen	Pod	Testis
6736	40009	R54034	1.88	10.99	5.85	1.00	0.00	1	570.96	Spleen	Pod	Testis
6753	757222	AA498149	6.26	449.41	71.83	2.00	0.00	1	570.96	Spleen	Pod	Testis
6759	45807	H08210	7.48	44.72	5.98	1.00	0.00	6	377.78	Stomach	Placenta	Pancreas
6764	22374	T82459	65.53	358.58	5.47	0.00	1.00	2	545.17	Colon	Breast	Pooled
6770	148225	H13888	9.30	108.85	11.68	1.00	0.00	19	234.91	Brain	LID not found	Other
6771	154172	R52030	8.68	147.36	16.97	6.00	0.00	15	238.91	Brain	Adipose	Ear
6772	22376	T82461	30.98	244.90	7.90	0.00	3.00	13	134.61	Brain	Testis	LID not found
6774	725321	AA291749	14.37	323.40	22.50	13.00	0.00	6	48.24	Umbilical cord	Umbilical cord	Skin
6776	33860	R44840	43.05	441.17	10.25	3.00	5.00	3	41.9	Aorta	Whole embryo	Brain
6782	47359	H11003	6.25	78.77	12.61	3.00	3.00	15	241.89	Prostate	Brain	LID not found
6784	41789	R39197	14.53	78.43	8.40	0.00	0.00	6	40.6	Brain	Gall bladder	Blood
6785	48567	H15114	2.70	20.39	7.55	4.00	0.00	21	177.48	Ignore	Parathyroid	Thymus
6787	491751	AA150500	23.89	209.55	8.78	1.00	0.00	10	359.28	CNS	Brain	Testis
6788	22393	T87228	11.27	74.60	6.62	0.00	0.00	8	318.21	Brain	LID not found	Other
6789	884482	AA629707	3.20	23.10	7.22	1.00	0.00	8	318.21	Brain	LID not found	Other
6792	33894	R44530	2.89	15.20	5.25	1.00	0.00	8	318.21	Brain	LID not found	Other
6796	45312	H08226	11.17	80.40	5.41	1.00	0.00	8	318.21	Brain	LID not found	Other
6800	47460	H11454	38.58	365.98	9.49	4.00	3.00	8	318.21	Brain	LID not found	Other
6804	22334	T88939	23.60	245.83	10.41	2.00	2.00	8	318.21	Brain	LID not found	Other
6808	49588	H15153	0.98	5.74	5.86	0.00	1.00	8	318.21	Brain	LID not found	Other
6813	380240	AA012839	13.32	77.41	5.61	0.00	1.00	8	318.21	Brain	LID not found	Other
6814	26568	R27815	13.58	213.43	15.71	3.00	3.00	8	318.21	Brain	LID not found	Other
6815	41495	R39111	17.06	172.52	10.11	1.00	0.00	8	318.21	Brain	LID not found	Other
6816	49810	R54073	4.76	39.60	8.32	0.00	1.00	8	318.21	Brain	LID not found	Other
6817	841314	H15288	3.11	19.28	6.18	1.00	0.00	8	318.21	Brain	LID not found	Other
6819	810801	AA487218	12.77	102.08	7.98	2.00	0.00	8	318.21	Brain	LID not found	Other
6820	52021	AA458878	9.55	245.62	25.73	8.00	0.00	8	318.21	Brain	LID not found	Other
6822	281194	H98215	20.67	135.94	6.58	1.00	0.00	8	318.21	Brain	LID not found	Other
6828	283108	H99837	6.38	49.90	7.82	1.00	0.00	8	318.21	Brain	LID not found	Other
6829	51599	H18932	2.30	40.76	17.73	3.00	0.00	8	318.21	Brain	LID not found	Other
6830	592728	AA160670	2.80	14.48	5.17	1.00	0.00	8	318.21	Brain	LID not found	Other
6834	278557	N82914	6.63	43.66	6.59	1.00	0.00	8	318.21	Brain	LID not found	Other
6835	611324	AA176819	222.27	1138.84	5.12	1.00	0.00	8	318.21	Brain	LID not found	Other
6840	41842	R52681	8.88	88.81	10.23	2.00	5.00	8	318.21	Brain	LID not found	Other
6842	272816	N38130	270.30	2000.42	7.40	0.00	0.00	8	318.21	Brain	LID not found	Other
6844	129020	R10823	13.63	105.83	7.76	0.00	1.00	8	318.21	Brain	LID not found	Other
6845	41345	R58985	19.86	110.16	5.51	1.00	0.00	8	318.21	Brain	LID not found	Other
6849	73244	T58987	32.67	187.52	5.13	1.00	0.00	8	318.21	Brain	LID not found	Other
6851	510790	AA102053	66.78	380.33	5.70	0.00	1.00	8	318.21	Brain	LID not found	Other
6858	40449	R53258	43.97	362.11	8.23	4.00	3.00	8	318.21	Brain	LID not found	Other
6859	841185	AA487070	65.89	631.29	9.58	2.00	0.00	8	318.21	Brain	LID not found	Other
6863	950355	AA600184	11.65	92.19	7.92	1.00	0.00	8	318.21	Brain	LID not found	Other
6865	781442	AA428603	4.50	32.34	7.19	1.00	0.00	8	318.21	Brain	LID not found	Other
6869	45801	H08208	8.48	63.24	7.46	1.00	4.00	8	318.21	Brain	LID not found	Other
6875	41070	R56100	9.76	51.53	5.28	1.00	0.00	8	318.21	Brain	LID not found	Other
6878	809738	AA454713	19.84	109.53	5.52	0.00	1.00	8	318.21	Brain	LID not found	Other

Table 2A

6880	50477	H17004	45.26	244.57	5.40	0.00	3.00	2	568.11	Brain	Heart
6885	50615	H17513	61.16	399.80	6.53	0.00	1.00	6	118.05	Pool	Fore skin
6890	897576	AA498871	12.94	71.13	5.50	0.00	1.00			Lung	LID not found
6896	40491	R55873	16.09	113.46	7.05	2.00	3.00	2	480.97	Eye	Fore skin
6897	71312	T47825	55.70	682.13	12.26	1.00	2.00	9	392.46	Cervix	Ovary
6901	74738	T57359	35.63	190.20	5.34	2.00	1.00			LID not found	Other
6902	564514	AA121697	20.51	105.24	5.13	0.00	1.00			Colon	Prostate
6903	809455	AA443089	225.23	1135.07	5.04	1.00	0.00	X	88.99	Parathyroid	Cervix
6909	86160	T72336	6.50	58.06	8.93	1.00	0.00			Pool	LID not found
6912	33814	R44714	3.01	23.17	7.70	2.00	0.00	9	382.37	Pool	LID not found
6914	121256	T86605	24.45	151.19	8.18	1.00	0.00			CNS	Placenta
6915	281870	N51838	1.80	9.09	5.56	3.00	0.00			Placenta	Colon
6931	133860	R27975	7.46	44.19	5.93	3.00	0.00	19	101.7	Pool	Parathyroid
6939	247089	N57856	0.94	20.21	21.42	1.00	0.00	X	143.33	Fore skin	Pool
6944	281280	N72228	8.49	84.31	9.93	1.00	0.00	8	401.87	Fore skin	Kidney
6948	271078	N29918	9.42	102.30	10.86	6.00	0.00			LID not found	Other
6963	133884	R28650	4.19	58.73	15.92	5.00	0.00	12	69.28	Cervix	Pancreas
6972	239815	H79538	55.96	335.56	6.00	0.00	1.00	3	162.18	Pool	LID not found
6982	810890	AA459278	7.92	42.08	5.32	1.00	0.00	1	88.34	Fore skin	Ovary
6984	270889	N32502	9.52	57.27	6.01	1.00	0.00	9	24.47	Pool	LID not found
6986	125118	R03293	38.74	202.03	5.21	1.00	0.00			Neural	LID not found
6990	758338	AA404278	32.23	169.82	5.27	1.00	0.00	9	407.26	Placenta	Breast
6995	134011	R31282	4.99	35.81	7.17	1.00	0.00			Brain	LID not found
7000	291827	N72876	7.25	51.54	7.11	0.00	1.00			Fore skin	LID not found
7004	268478	N25798	4.63	27.76	6.00	0.00	1.00	20	257.89	Head and nec	Adipose
7010	323988	V46433	15.18	104.59	6.89	2.00	0.00	11	227.2	Uterus	Pool
7012	294535	N71028	6.27	55.78	8.89	3.00	0.00	2	519.95	Brain	Brain
7014	428756	AA05219	8.72	72.55	7.47	1.00	0.00	2	334.11	Tonsil	Fore skin
7017	489800	AA098820	7.65	77.01	10.06	1.00	0.00	11	363.53	Lung	LID not found
7020	207275	H59618	0.97	5.72	5.90	1.00	0.00			Pool	LID not found
7034	307138	N93721	2.30	14.14	6.16	1.00	0.00	10	508.28	Fore skin	Parathyroid
7036	123811	R01448	228.40	1373.17	6.01	1.00	2.00			Prostate	Placenta
7037	344854	W72972	7.37	132.81	18.02	8.00	0.00			Colon	Lung
7043	268258	N30006	61.41	614.11	7.54	2.00	0.00			Prostate	Prostate
7046	359689	W96473	4.72	28.95	5.70	1.00	0.00			Breast	Tonsil
7047	324333	AA284109	40.55	225.22	5.55	1.00	0.00			Fore skin	LID not found
7053	321886	W37628	17.28	114.54	6.64	1.00	0.00			Prostate	Prostate
7057	299182	N70553	4.89	30.57	6.25	0.00	1.00			Fore skin	LID not found
7060	207932	H60514	131.17	903.67	6.89	2.00	2.00			Fore skin	LID not found
7081	321908	W37860	5.20	98.20	18.89	9.00	3.00	5	287.05	Pancreas	Pool
7084	769800	AA425900	2.06	24.85	12.08	7.00	0.00	3	161.52	Gall bladder	Lymph
7072	250883	N23464	26.30	198.48	7.01	0.00	1.00	12	246.58	Kidney	LID not found
7075	143145	R73661	161.71	1126.96	6.97	2.00	0.00			Fore skin	LID not found
7078	417059	W87801	12.82	85.35	6.66	2.00	4.00	11	130.57	CNS	Placenta
7078	503334	AA134111	14.22	88.89	6.25	1.00	0.00	2	174.16	Umbilical cord	Larynx
7080	877827	AA625632	1739.50	9069.75	5.23	1.00	0.00			Broad	Placenta
7087	345034	W72294	7.49	552.19	73.72	5.00	0.00			Umbilical cord	Stomach
7090	245174	N54458	53.28	298.56	5.82	0.00	1.00	8	99.74	Pancreas	LID not found
7092	240878	H80286	120.80	1407.02	11.65	4.00	3.00	1	111.21	Pancreas	LID not found
7093	347613	W81504	12.41	71.32	5.75	0.00	1.00	3	43.88	Prostate	Fore skin
7094	488276	AA65759	8.44	91.42	9.65	1.00	0.00			Fore skin	Fore skin
7096	460114	AA678840	6.39	41.43	6.48	1.00	0.00	13	268.41	Whole embryo	Testis
7097	298104	N70759	23.56	169.43	7.19	0.00	1.00	4	17.79	Ovary	Parathyroid
7106	769712	AA428958	10.64	53.49	5.03	1.00	0.00			Pool	Pool

Table 2A

7108	67625	T48530	7.62	72.88	9.56	0.00	1.00	18	347.76	Uterus	Placenta	Pool
7110	387495	AA701081	3.32	18.27	5.50	2.00	0.00	5	28.67	Brain	LID not found	Other
7111	47355	H10961	2.28	19.97	8.82	1.00	0.00	19	283.98	Peripheral ner	Cervix	Umbilical cord
7114	264646	N20338	2.28	12.36	5.42	1.00	0.00	20	333.71	CNS	Kidney	Heart
7132	68636	T49802	15.29	161.47	10.56	2.00	2.00	17	24.51	Parathyroid	CNS	Thyroid
7134	51743	H23081	179.73	1030.55	5.73	1.00	0.00	11	239.66	Blood	Eye	Ovary
7139	344430	W73473	7.12	470.67	66.11	6.00	2.00	7	424.13	CNS	Adrenal gland	Pancreas
7148	70152	T50041	44.92	235.92	5.25	1.00	0.00	2	511.27	Peripheral ner	Brain	Aorta
7152	60226	T64216	0.72	4.13	5.71	0.00	1.00	7	186.21	Synovial mem	Smooth musc	Ear
7153	146866	R60779	5.70	51.01	8.95	1.00	0.00	14	251	Liver	Spleen	Pool
7159	450453	T62577	60.57	309.75	5.11	1.00	0.00	1	619.07	Brain	LID not found	Other
7165	51070	AA682815	5.13	35.84	6.98	0.00	2.00	12	246.56	Placenta	LID not found	Other
7167	51070	H17115	6.52	37.80	5.80	1.00	0.00	1	284.5	Testis	Brain	LID not found
7172	77193	T50121	14.64	77.62	5.30	1.00	0.00	14	173.38	Brain	LID not found	Other
7187	154472	R54846	61.81	368.99	6.99	1.00	0.00	9	21.93	Brain	LID not found	Other
7184	248412	N66568	3.15	15.91	5.05	1.00	0.00	8	433.98	Small intestine	Colon	LID not found
7209	41548	R66415	14.47	147.85	10.23	5.00	5.00	1	289.5	Uterus	Placenta	Kidney
7215	60882	T40868	49.42	314.59	6.37	0.00	2.00	5	283.38	Aorta	Kidney	Kidney
7216	40108	R52633	14.47	90.66	6.27	0.00	2.00	12	277.88	Brain	LID not found	Other
7221	51939	H22956	2.78	45.22	16.25	6.00	1.00	5	580.59	Eye	Synovial mem	Kidney
7232	40036	R53442	77.36	828.19	8.12	0.00	4.00	4	558.61	Pool	Brain	Lung
7245	51397	H18417	2.07	10.70	5.16	1.00	0.00	6	328.36	Umbilical cord	Smooth musc	Thyroid
7246	511809	AA088861	85.14	892.20	10.48	5.00	5.00	5	542.07	Brain	LID not found	Other
7247	279577	N48899	41.02	257.30	6.27	0.00	2.00	11	282.56	Liver	Whole embryo	LID not found
7248	32095	R42898	84.42	625.57	7.41	4.00	4.00	1	562.43	Synovial mem	Thyroid	Blood
7254	324122	W46577	70.03	484.02	6.91	0.00	2.00	20	74.7	Ear	Placenta	Parathyroid
7277	52432	H23256	1.88	10.30	5.55	0.00	1.00	12	469.76	Blood	Germ Cell	Aorta
7278	855523	AA684180	30.82	1750.51	57.17	3.00	0.00	16	423.94	Bone	Blood	CNS
7285	52056	H24347	2.15	17.94	8.33	1.00	1.00	8	334.17	CNS	LID not found	Other
7287	744593	AA621266	45.43	237.63	5.23	0.00	1.00	3	728.84	Brain	LID not found	Other
7288	40768	R56432	118.71	1367.55	11.52	4.00	2.00	6	172.31	Nose	Spleen	Placenta
7289	39770	R64558	2.09	13.28	6.36	0.00	1.00	7	586.57	Heart	LID not found	Other
7292	61639	T41032	3.96	28.10	6.59	1.00	0.00	1	740.89	CNS	Uterus	Pool
7294	324715	W47362	6.80	107.28	15.77	8.00	0.00	10	165.42	Nose	Whole embryo	Cervix
7298	122183	T98628	12.28	66.14	5.40	1.00	0.00	6	118.49	Small intestine	Thymus	Lung
7302	810597	AA464576	32.78	171.37	5.23	0.00	1.00	8	91.28	Foreskin	Whole embryo	Germ Cell
7306	430264	AA010557	15.95	129.35	8.11	3.00	5.00	12				
7314	415250	W91885	51.41	315.39	8.14	2.00	0.00					
7322	809609	AA458486	19.20	107.04	5.57	1.00	0.00					
7324	280231	N62273	7.16	43.02	6.01	1.00	0.00					
7332	271280	N34637	211.90	1086.74	5.13	1.00	0.00					
7335	418279	W90323	30.51	351.39	11.52	1.00	0.00					
7338	234121	H70603	24.14	180.89	7.49	2.00	0.00					
7339	281162	N50982	7.04	39.13	5.56	1.00	0.00					
7340	289867	N62080	17.77	210.46	11.85	6.00	5.00					
7363	178656	H49517	427.45	2222.49	5.20	1.00	0.00					
7366	810727	AA457718	9.43	277.32	29.09	1.00	0.00					
7367	201217	R98293	98.79	571.85	5.79	3.00	0.00					
7368	505597	AA147654	3.78	28.48	7.00	1.00	0.00					
7370	428749	AA004652	8.89	71.11	8.00	0.00	4.00					
7383	201855	H48251	107.45	868.70	8.08	2.00	0.00					
7418	855547	AA684195	65.02	863.56	13.26	3.00	0.00					
7418	272038	N31848	2.92	92.77	31.75	1.00	0.00					



Table 2A

7427	376475	AA041386	104.18	533.03	5.12	1.00	0.00	17	53.69	Lung	Placenta	Ovary
7432	857284	AA688689	17.25	137.53	7.97	0.00	1.00	6	380.74	Brain	Foreskin	Codon
7434	359285	AA016234	36.46	518.68	14.23	0.00	1.00			Foreskin	Heart	Kidney
7440	293925	N63943	128.20	761.43	5.89	1.00	2.00	16	162.07	Pool	LID not found	Other
7449	418750	W86521	12.78	74.51	5.83	1.00	0.00	10	498.96	Umbilical cord	Thyroid	Ear
7455	210486	H65478	191.51	1175.82	6.14	1.00	0.00	13	136.15	Lymph	Blood	Heart
7456	344586	W73144	14.93	512.53	34.96	6.00	0.00			Testis	LID not found	Other
7462	762782	AA448167	9.71	53.61	5.54	1.00	0.00	3	18.66	Foreskin	Heart	Pool
7466	272890	N36123	65.27	366.56	5.62	2.00	0.00	9	416.03	Ovary	Whole embryo	Prostate
7468	741880	AA402883	6.99	159.54	22.82	4.00	0.00	14	123.72	Small intestine	Head and nec	Esophagus
7480	855910	AA630328	43.72	219.10	5.01	1.00	0.00			CNS	Eye	Kidney
7482	284355	N52136	5.42	47.59	6.78	1.00	1.00	19	244.11	Placenta	Pool	Parathyroid
7488	362059	AA001432	3.88	66.74	17.21	0.00	1.00	18	185.94	Pancreas	Podied	Eye
7489	24918	R39039	6.32	141.98	17.05	5.00	1.00	9	226.16	Eye	Colon	Aorta
7493	510381	AA055565	37.85	312.68	6.24	2.00	0.00	10	32.49	Eye	Skin	Esophagus
7494	808828	AA455521	13.11	127.97	9.76	1.00	0.00	8	417.63	Tonsil	Blood	Parathyroid
7496	46051	H09064	11.43	63.49	5.55	0.00	1.00	7	624.62	Brain	LID not found	Other
7500	50847	H17929	10.39	60.06	5.78	1.00	0.00	8	276.6	Brain	LID not found	Other
7504	47459	H11453	3.99	111.46	27.93	7.00	0.00	1	192.55			
7512	45531	H08541	1.57	8.90	5.85	2.00	1.00	2	487.92		Kidney	LID not found
7516	50749	H17322	4.74	56.00	11.82	4.00	5.00	8	399.21	Brain	Blood	Brain
7517	770424	AA430675	20.00	128.48	6.43	1.00	0.00	9	137.48	Ovary	Cervix	
7523	773479	AA427889	303.30	1788.94	5.90	0.00	1.00	6	116.42	Adrenal gland	Lymph	
7532	50173	H17464	23.31	224.14	9.62	2.00	5.00	2	356.74	Brain	LID not found	Other
7533	32304	R42894	6.10	66.28	10.87	3.00	0.00	3	180.69	CNS	Brain	Eye
7534	742115	AA405800	22.26	175.16	7.87	5.00	0.00	16	48.85	Pooled	Blood	Adrenal gland
7540	22541	T89084	1.02	6.25	6.10	0.00	1.00	8	447.85	Brain	LID not found	Other
7541	140574	R66139	4.20	72.19	17.20	2.00	0.00	16	365.99	Brain	Breast	Prostate
7546	855880	AA630320	18.89	106.74	5.62	1.00	0.00	8	436.63			
7551	69835	T48692	18.11	106.23	5.87	1.00	0.00	10	489.34	Ear	Placenta	CNS
7554	41648	R52798	18.22	94.88	5.21	0.00	1.00			Brain	Testis	Lung
7555	196488	R81539	17.69	110.72	6.26	2.00	0.00	13	147.98	Whole embryo	Pancreas	CNS
7556	34010	R44647	41.18	374.58	9.10	6.00	5.00	16	183.03	Lymph	Brain	LID not found
7557	135630	R31582	10.59	64.53	6.10	2.00	1.00	4	450.16	Placenta	Pancreas	Ovary
7571	743081	AA405901	78.80	579.47	7.35	0.00	1.00			Prostate	LID not found	Kidney
7576	52716	H28245	1.27	15.84	12.58	0.00	2.00	11	37.97	Skin	Pooled	Germ Cell
7579	773106	AA435316	18.63	103.70	5.57	1.00	0.00	5	578.31	Adipose	Stomach	Prostate
7582	291985	N73101	23.68	799.69	33.80	4.00	1.00	19	284.15	Pool	Brain	Lung
7583	50299	H17854	2.84	15.39	5.42	1.00	0.00			Brain	LID not found	Other
7592	51828	H22846	152.46	1025.58	6.73	1.00	3.00	4	437.92	Tonil		Kidney
7595	248027	N68372	98.78	613.02	6.21	2.00	2.00			Brain	LID not found	Other
7600	48993	H28738	15.04	81.94	5.45	0.00	1.00	3	49.58	Umbilical cord	vein	Synovial membrane
7601	41607	R64177	222.18	1286.36	5.79	1.00	1.00	2	356.03	Ovary	Blood	Pooled
7603	279085	N51705	46.41	280.13	6.04	2.00	0.00	14	192.86	Pooled	CNS	Germ Cell
7612	202912	H99659	84.43	445.79	5.28	1.00	0.00			Thymus	CNS	Pancreas
7621	67187	T52652	341.51	1997.97	5.85	1.00	0.00			Liver	Eye	Ovary
7624	84211	T72915	26.41	168.43	5.76	1.00	0.00			Ear	Muscle	Whole embryo
7627	609155	AA176867	15.19	96.55	6.36	0.00	1.00	1	35.21	CNS	Testis	Germ Cell
7631	757265	AA426113	15.72	136.39	8.87	1.00	2.00	10	360.96	Brain	Pool	LID not found
7633	50689	H17046	35.92	328.24	9.14	2.00	5.00			Aorta		Whole embryo
7646	772918	AA479913	21.46	165.22	7.89	0.00	1.00	11	259.91			
7657	26414	R20862	2.72	26.43	9.72	3.00	0.00	1	31.2	Eye	Cervix	Blood
7660	840062	AA405249	47.31	339.43	7.18	1.00	0.00	8	407.58	Codon	LID not found	Other
7662	510736	AA089748	5.97	42.15	7.08	1.00	0.00					

Table 2A

7672	46565	H09759	6.80	35.10	5.32	1.00	0.00	540.74	Brain	Whole embryo	LID not found
7673	41103	R66769	61.03	633.49	10.38	1.00	3.00	Brain	Brain	LID not found	Other
7677	45318	H08734	15.21	76.21	5.01	0.00	1.00	Parathyroid	Parathyroid	Tonsil	Whole embryo
7688	284105	N20577	79.93	400.02	5.00	1.00	0.00	Eye	Forebrain	Forebrain	Lymph
7700	281687	H69780	324.08	1849.12	5.71	1.00	0.00	271.02	Pooled	Germ Cell	Placenta
7707	135108	R33037	7.07	49.53	7.00	1.00	0.00	95.03	Forebrain	LID not found	Other
7708	282334	H99384	220.65	1603.13	7.27	3.00	0.00	353.04	Adrenal gland	Tonsil	Brain
7712	272531	N35889	23.74	130.74	5.51	2.00	0.00	18.43	Ear	CNS	Pooled
7714	762269	AA431746	10.11	52.08	5.15	1.00	0.00	672.99	Adrenal gland	Tonsil	Brain
7715	795730	AA460282	9.82	52.06	5.30	1.00	0.00	846.44	Blood	Muscle	Lung
7718	366763	AA029331	11.83	97.51	8.24	2.00	0.00	351.05	Forebrain	Pancreas	Whole embryo
7720	269787	N27145	6.78	39.34	5.80	1.00	0.00	223.68	Forebrain	Spleen	Germ Cell
7722	809487	AA443105	11.92	67.44	5.86	0.00	1.00	Placenta	Placenta	Tonsil	Heart
7731	360025	AA063573	2.14	12.19	5.89	1.00	0.00	584.18	Pool	LID not found	Other
7740	415111	W63147	6.44	43.25	6.71	3.00	0.00	732.12	Pool	LID not found	Other
7742	430336	AA010619	0.92	6.25	8.96	1.00	0.00	250.28	Prostate	Pool	LID not found
7748	426726	AA011678	10.31	67.94	6.59	0.00	1.00	Tonsil	Pool	LID not found	Other
7750	427931	AA001983	4.22	22.29	5.28	1.00	0.00	395.81	CNS	Placenta	Germ Cell
7751	288741	N59219	5.70	43.43	7.62	2.00	0.00	46.33	Forebrain	Heart	LID not found
7752	375853	AA039857	20.17	121.67	6.03	1.00	1.00	253.28	Aorta	Colon	Brain
7762	129882	R17086	6.39	65.40	7.79	1.00	0.00	226.12	CNS	Forebrain	Pancreas
7771	291059	N72116	135.89	899.23	6.62	3.00	0.00	108.41	Pool	LID not found	Other
7774	427893	AA001359	9.85	56.52	5.74	1.00	0.00	Neural	Neural	Prostate	LID not found
7760	429489	AA001480	8.29	43.83	5.30	0.00	1.00	17.4	Uterus	Pool	Brain
7791	415562	W60701	5.98	358.44	60.12	20.00	1.00	535.4	Uterus	Eye	Testis
7798	502593	AA136049	1.00	7.22	7.25	1.00	0.00	111.36	Adipose	Ovary	Pool
7804	293097	N68738	42.15	344.25	6.17	3.00	3.00	685.61	Heart	Adrenal gland	Whole embryo
7806	505183	AA151111	112.19	820.11	5.53	1.00	0.00	476.39	Germ Cell	Breast	Heart
7814	810002	AA454864	21.22	191.72	9.04	6.00	1.00	106.17	Pool	LID not found	Other
7816	298062	N70734	3.63	27.97	7.71	2.00	0.00	17.11	Tonsil	Whole embryo	Pool
7819	364865	AA035745	5.14	27.49	5.34	1.00	0.00	269.6	Muscle	CNS	Forebrain
7820	247482	N54157	17.34	100.47	5.78	0.00	2.00	504.52	Liver	Pool	Forebrain
7821	376333	AA041293	10.15	68.51	8.75	1.00	0.00	592.45	Tonsil	Kidney	Pool
7827	197285	R68970	1.90	12.09	6.38	1.00	0.00	599.98	Colon	Pool	LID not found
7833	480784	AA133194	1.33	8.69	6.53	1.00	0.00	394.49	Ovary	Prostate	Parathyroid
7836	211367	H66670	482.02	3163.12	6.54	0.00	2.00	77.53	Smooth muscle	Nose	Other
7840	289428	N63949	9.38	52.36	5.58	1.00	0.00	LID not found	Other	LID not found	Other
7842	129227	R11047	27.31	185.70	6.80	0.00	1.00	87.87	Blood	Lymph	Placenta
7843	196579	R01689	90.78	592.31	6.52	2.00	0.00	367.22	Uterus	Heart	Pool
7844	211865	H66710	7.06	120.70	17.10	3.00	0.00	726.84	Esophagus	Spleen	Prostate
7846	810283	AA464728	51.66	303.66	5.88	1.00	0.00	28.79	Pooled	Placenta	Blood
7847	417867	W60128	62.61	526.71	8.41	2.00	0.00	116.59	Liver	Gall bladder	Blood
7852	288452	N74825	7.86	104.35	13.28	4.00	3.00	36.48	Blood	Thyroid	Kidney
7854	810320	AA464140	3.47	25.03	7.21	2.00	0.00	445.78	Adipose	Spleen	Umbilical cord
7858	259017	N32811	72.73	450.82	6.20	2.00	0.00	444.53	Stomach	Ear	Ear
7860	293110	N63646	44.70	732.45	16.39	3.00	3.00				
7862	770860	AA434388	280.85	2734.72	9.74	2.00	0.00				
7863	357138	W63523	7.98	41.66	5.22	1.00	0.00				
7867	201071	R99847	123.03	1038.86	6.44	3.00	0.00				
7878	72426	T51592	135.03	1085.12	8.04	2.00	0.00				
7878	212115	H68885	179.08	1043.98	5.83	2.00	0.00				
7883	724177	AA401441	51.32	599.25	11.68	2.00	0.00				
7900	72616	T51995	5.01	41.17	6.22	2.00	0.00				
7901	624443	AA181333	1.35	18.06	11.88	1.00	0.00				
7908	73561	T55592	12.72	78.28	8.15	1.00	0.00				

Table 2A

7915	342161	W83749	16.94	317.86	18.76	1.00	0.00	18	414.52	Thyroid	Fore skin	CNS
7916	32517	R43271	6.35	38.06	5.89	2.00	0.00	18	193.03	Pancreas	Placenta	Kidney
7922	141485	R73564	3.36	28.83	8.54	1.00	0.00			Brain	Ovary	Uterus
7927	52226	H23265	3.65	81.73	16.93	7.00	1.00			Head end nec	Thymus	Cervix
7930	69184	T54144	187.62	1000.39	5.33	0.00	1.00			CNS	Lung	LID not found
7936	80912	T70032	64.91	520.69	8.02	0.00	4.00	1	151.73	Pooled	Germ Cell	Whole embryo
7937	194236	R83277	97.87	759.02	7.77	1.00	3.00			Uterus	Whole embryo	Lung
7938	41358	R59167	38.17	457.85	12.00	0.00	1.00	8	486.43	Smooth musc	Larynx	Fore skin
7940	72016	T52325	94.28	859.21	8.11	0.00	2.00			Small intestine	Smooth musc	Ear
7949	970813	AA683102	25.48	161.66	9.69	0.00	1.00	5	344.82	Thyroid	Umbilical cord	Blood
7954	433307	AA689732	4.03	39.06	9.69	0.00	1.00					
7961	243159	H94471	4.48	52.60	11.73	6.00	2.00					
7968	85224	T71578	54.33	323.37	5.95	0.00	1.00	19	61.79	Ear	Ovary	Heart
7977	33998	R44538	26.16	207.03	7.92	2.00	0.00	5	489.95	Brain	LID not found	Other
7985	34014	R44564	2.16	14.37	6.65	2.00	1.00			Pancreas	Brain	LID not found
7992	46919	H10226	12.04	84.81	5.38	2.00	2.00	1	553.84	Brain	Pool	LID not found
7993	48287	H09620	4.66	31.48	6.75	1.00	2.00	2	718.03	Stomach	Esophagus	Synovial membrane
7994	770704	AA476284	106.98	1442.59	13.48	0.00	3.00	4	672.97	Eye	CNS	Uterus
7995	839698	AA504658	30.40	421.81	13.88	1.00	0.00	17	60.19	Placenta	Brain	LID not found
8001	41850	R52766	2.93	48.37	18.84	4.00	0.00	22	67.29	Spleen	Lung	LID not found
8005	49887	H15250	2.28	22.85	10.10	2.00	0.00	3	459.29	Brain	LID not found	Other
8012	71591	T48011	28.48	162.89	5.71	2.00	0.00	3	92.81	Adipose	Skin	Tonsil
8013	50865	H18458	2.18	13.55	6.20	0.00	1.00	6	443.43	Parathyroid	Prostate	Tonsil
8021	50240	H17083	8.03	48.47	6.16	2.00	3.00	1	49.6	Thyroid	Placenta	LID not found
8023	78736	T81888	9.63	66.65	6.92	0.00	2.00	1	49.6	Brain	LID not found	Other
8025	52543	H23482	6.32	35.56	8.28	1.00	0.00	2	358.3	Brain	LID not found	Other
8028	69893	T48849	289.13	3420.38	11.83	2.00	0.00	4	74.59	Brain	LID not found	Other
8029	50860	H17981	5.06	41.26	8.15	4.00	4.00	1	841.15	Prostate	Skin	Esophagus
8032	45999	H09317	98.43	1055.20	10.94	3.00	4.00	12	228.02	Larynx	Skin	Esophagus
8039	730288	AA412509	16.30	105.82	6.49	1.00	0.00	9	59.14	Fore skin	Uterus	Brain
8045	51020	H18312	88.66	821.16	9.28	2.00	2.00	3	546.7	Small intestine	Head end nec	Synovial membrane
8055	34304	R44210	23.89	225.88	9.48	4.00	0.00	16	471.03	Thyroid	Breast	
8057	47452	H11448	48.02	482.48	10.48	3.00	3.00	X	245.06	CNS	Heart	Germ Cell
8058	592111	AA150532	5.84	76.09	13.03	5.00	0.00	2	557.51	Pool	Prostate	Colon
8059	785837	AA461511	15.79	170.63	10.81	1.00	0.00			Aorta	Whole embryo	Adrenal gland
8061	49281	H15653	211.47	1731.38	8.19	2.00	0.00	15	245.32	Eye	Brain	Heart
8063	773446	AA428049	56.11	307.25	5.48	1.00	0.00	17	460.48	Ignoie	Pool	LID not found
8070	795756	AA460304	115.08	797.22	8.93	1.00	0.00	1	251.6	CNS	Lung	LID not found
8086	811138	AA485730	21.28	287.53	13.51	1.00	0.00	22	17.11	Heart	CNS	Whole embryo
8087	429685	AA011598	16.54	142.66	8.63	6.00	0.00	11	20.3	Ovary	CNS	Blood
8098	365004	AA024832	507.22	3348.21	6.80	2.00	0.00	12	258.32	Ovary	CNS	Lymph
8098	416945	W87710	22.16	188.28	7.59	4.00	3.00					
8098	306921	N91962	37.34	212.50	5.89	1.00	0.00					
8103	194156	H51050	29.21	180.89	8.19	2.00	0.00					
8104	359684	AA011100	2.87	19.37	6.74	1.00	0.00					
8106	488885	AA043092	8.85	135.78	15.35	6.00	1.00					
8111	194023	H51271	27.11	166.67	6.89	1.00	0.00					
8112	282315	N51861	19.02	141.62	7.44	3.00	0.00					
8116	288938	N71147	21.75	165.66	7.62	0.00	3.00					
8118	358936	W82233	3.80	19.07	5.01	1.00	0.00					
8119	366763	AA028703	7.16	43.55	6.09	1.00	0.00					
8122	811038	AA485424	20.92	129.21	6.18	1.00	0.00					
8127	278392	N48700	2.36	13.17	5.58	1.00	0.00					
8128	770840	AA427733	4.27	75.03	17.59	3.00	0.00					

Table 2A

8144	306066	N81003	5.22	64.81	12.42	2.00	1.00	CNS	Adrenal gland	Parathyroid
8154	209189	M83858	58.77	580.30	10.04	3.00	0.00	474.75	LID not found	Other
8162	277134	N40917	5.84	30.71	5.28	1.00	0.00	Blood	Adrenal gland	CNS
8163	298809	N70848	82.67	502.00	6.07	2.00	0.00	628.75	Uterus	Pool
8165	491097	AA135889	6.21	36.83	5.93	1.00	0.00	Blood	Colon	Colon
8166	770969	AA430829	1200.45	8826.70	5.89	1.00	0.00	243.59	Synovial mem	Umbilical cord
8172	325641	W51985	7.29	46.71	6.41	0.00	1.00	63.51	LID not found	Other
8176	864844	AA629897	4.67	168.95	36.21	6.00	0.00	466.22	Gall bladder	Umbilical cord
8178	283968	N53360	89.68	938.12	10.47	2.00	0.00	58.88	Spleen	Blood
8182	711142	AA427778	140.87	813.47	5.78	1.00	0.00	71.41	Trachea	Bone marrow
8183	277187	N40932	19.94	118.12	5.92	1.00	0.00	576.51	Adrenal gland	Adipose
8189	487317	AA043800	3.63	22.09	6.26	1.00	0.00	553.01	Gall bladder	Thyroid
8191	811612	AA455013	8.83	111.24	12.80	2.00	0.00	272.85	Peripheral ner	Germ Cell
8192	138369	R88108	10.67	58.09	5.44	1.00	0.00	Uterus	Lung	LID not found
8193	427697	AA001884	36.11	198.63	5.50	1.00	0.00	Testis	Parathyroid	Kidney
8194	245688	N53355	29.29	263.98	9.01	3.00	0.00	227.72	Colon	Omentum
8195	347035	W81135	36.83	219.77	5.97	0.00	1.00	44.84	Pool	Ovary
8196	760947	AA429661	5.58	30.05	11.30	6.00	0.00	227.72	Colon	Lymph node
8197	501479	AA115328	1.91	21.54	11.30	6.00	0.00	549.88	Fore skin	Lung
8205	501854	AA127885	5.30	47.84	9.05	3.00	0.00	516.86	Stomach	Eye
8208	795755	AA460313	4.40	46.86	10.66	1.00	0.00	422.8	Ovary	Uterus
8208	1035889	AA628189	3.88	78.59	20.38	2.00	0.00	397.3	Thymus	Fore skin
8215	207838	H60286	4.92	68.36	13.91	1.00	0.00	Prostate	Pool	LID not found
8216	855521	AA864179	113.25	1550.11	13.69	5.00	0.00	37.18	Small intestine	Colon
8217	430073	AA010000	5.68	42.92	7.56	0.00	2.00	308.04	Small intestine	Colon
8218	322841	W15318	2.21	29.98	13.56	2.00	0.00	308.17	Whole embryo	Ovary
8220	361122	AA017379	5.73	57.70	10.07	1.00	0.00	524.94	Whole embryo	Ovary
8229	808957	AA455128	19.13	101.19	5.29	1.00	0.00	28.78	Gall bladder	Blood
8231	195357	R88982	183.22	1413.03	7.71	2.00	0.00	253.83	Head and nec	Blood
8233	417761	W88725	18.49	325.08	17.58	6.00	0.00	151.82	Brain	LID not found
8234	291594	N67810	155.42	822.84	5.94	2.00	0.00	680.08	Brain	LID not found
8239	192188	H41144	186.31	2383.91	12.69	2.00	0.00	336.84	Whole embryo	Fore skin
8240	586786	AA133469	13.32	81.91	6.15	2.00	0.00	467.91	Whole embryo	CNS
8248	869375	AA679907	9.56	94.38	9.85	6.00	0.00	411.77	Germ Cell	Brain
8250	291241	N72210	35.63	295.93	8.30	2.00	0.00	Eye	Ear	Brain
8254	770898	AA434411	16.63	63.65	5.03	1.00	0.00	Brain	LID not found	Other
8257	510002	AA052832	1.07	27.77	25.95	3.00	0.00	Brain	LID not found	Other
8258	277186	N40945	2.83	27.85	9.77	3.00	0.00	Brain	LID not found	Other
8263	48953	H29276	32.08	332.56	10.37	4.00	0.00	Brain	LID not found	Other
8264	40217	R55658	1.75	8.54	5.46	1.00	0.00	Brain	LID not found	Other
8265	263118	N20108	3.57	17.97	5.03	0.00	1.00	Brain	LID not found	Other
8272	47400	H10413	2.62	20.09	7.67	2.00	0.00	Brain	LID not found	Other
8276	24081	R37568	9.23	51.58	5.59	1.00	0.00	Brain	LID not found	Other
8278	755228	AA486334	3.49	55.93	16.04	1.00	0.00	Brain	LID not found	Other
8280	52704	H29227	2.86	64.43	21.80	10.00	1.00	Brain	LID not found	Other
8289	177772	H45976	6.41	49.22	7.88	2.00	0.00	Brain	LID not found	Other
8292	41835	R54212	2.70	53.36	19.74	4.00	0.00	Brain	LID not found	Other
8314	845663	AA670155	32.08	204.40	6.37	2.00	0.00	Brain	LID not found	Other
8317	460666	AA700322	4.23	87.39	20.88	3.00	0.00	Brain	LID not found	Other
8320	46383	H08601	50.03	357.26	7.14	0.00	5.00	Brain	LID not found	Other
8328	345525	W72437	49.22	522.50	10.61	0.00	1.00	Brain	LID not found	Other
8340	480519	AA122287	6.91	36.19	5.24	1.00	0.00	Brain	LID not found	Other
8353	798176	AA461088	88.06	498.57	5.84	1.00	0.00	Brain	LID not found	Other
8359	772377	AA404585	48.90	236.81	5.05	1.00	0.00	Brain	LID not found	Other

Table 2A

8362	784017	AA443706	0.72	7.06	9.87	1.00	1.00	Aorta	Lung	Brain
8366	782843	AA448261	4.33	40.52	9.36	4.00	0.00	Testis	Pool	LID not found
8374	487327	AA045524	9.70	286.81	29.78	6.00	0.00	Uterus	Heart	Breast
8376	49203	H15695	8.34	147.99	17.74	2.00	0.00	Prostate	Whole embryo	Kidney
8377	52618	H29303	71.11	490.54	6.90	1.00	1.00	494.82	Blood	Eye
8384	784212	AA448975	15.96	103.37	6.48	3.00	0.00	245.05	Liver	Parathyroid
8386	428640	AA446864	7.32	44.00	6.01	2.00	0.00	278.76	Lymph	Pool
8388	787038	AA011837	15.77	111.34	7.06	0.00	3.00	62.6	Cervix	Parathyroid
8392	73785	AA496247	82.25	417.75	5.08	1.00	0.00	61.75	Smooth muscle	Gall bladder
8395	798887	AA463516	6.75	56.21	8.33	1.00	0.00	127.04	Spleen	Brain
8399	850700	AA608572	23.97	221.39	9.24	1.00	0.00	184.58	Stomach	Whole embryo
8402	626088	AA211459	103.86	689.28	8.55	2.00	0.00	363.57	Bone	Colon
8405	796127	AA460965	104.93	530.01	5.05	0.00	1.00	277.72	Head and neck	Cervix
8413	71557	T47871	132.31	713.23	5.39	1.00	1.00	66.18	Placenta	Gall bladder
8416	51015	H18307	126.36	866.36	8.70	2.00	0.00	Whole embryo	Parathyroid	
8419	784229	AA468887	10.51	86.75	6.25	1.00	0.00	401.78	Forebrain	CNS
8422	773073	AA425302	342.85	3025.83	7.90	2.00	2.00	436.1	Ear	Kidney
8426	430172	AA010247	34.82	205.57	5.94	2.00	0.00	697.77	Thyroid	Germ Cell
8433	51165	H17143	18.15	124.13	6.48	0.00	1.00	189.43	Pool	LID not found
8438	58053	AA133001	163.83	1729.28	10.55	4.00	4.00	354.01	Aorta	Forebrain
8445	62263	T41173	6.15	53.05	8.62	1.00	0.00	428.27	Smooth muscle	Stomach
8446	253733	N22552	7.41	46.80	6.33	1.00	0.00	137.65	Ear	Uterus
8478	410350	N93772	5.67	36.30	6.40	2.00	0.00	13.47	Ear	Heart
8482	307157	N93740	157.82	1059.74	8.71	1.00	0.00	392.03	Pool	LID not found
8486	487861	AA054722	75.56	453.46	6.00	1.00	0.00	372.61	Lung	Pool
8495	504372	AA142842	53.29	284.48	6.34	1.00	0.00	31.46	Ear	GNS
8498	247261	N57950	46.41	284.35	6.13	2.00	0.00	53.48	Aorta	Whole embryo
8502	297899	N70059	30.25	228.02	7.47	2.00	1.00	Synovial mem	Breast	Germ Cell
8508	241648	H91615	13.30	131.27	9.87	3.00	0.00	491.75	Stomach	Umbilical cord
8510	782794	AA448182	14.21	77.12	5.43	1.00	0.00	28.83	Testis	Colon
8516	302180	N79889	7.88	68.22	8.66	1.00	0.00	510.24	Eye	Testis
8528	281558	N72882	107.53	554.40	5.16	1.00	0.00	410.63	Forebrain	Breast
8531	809731	AA455508	7.08	72.02	10.17	1.00	0.00	-8.01	Forebrain	Ovary
8532	258242	N30655	5.78	53.88	9.31	1.00	0.00	518.98	Brain	Heart
8538	811048	AA485428	10.80	67.50	6.25	1.00	0.00	460.63	Adipose	Stomach
8542	201440	R09105	7.19	87.44	12.18	8.00	0.00	232.44	Pool	LID not found
8547	203008	H54263	40.12	205.61	5.12	1.00	0.00	Germ Cell	Parathyroid	Prostate
8558	771257	AA443594	4.60	24.67	5.36	1.00	0.00	442.78	Pool	LID not found
8560	454083	AA676888	17.83	148.22	8.37	4.00	0.00	Whole embryo	Brain	Brain
8563	204478	H58808	3.94	20.02	5.08	1.00	0.00	Colon	Pool	LID not found
8564	810457	AA457137	5.39	58.93	10.94	6.00	0.00	Ovary	Pooled	Lung
8568	853368	AA653110	50.44	283.41	5.82	0.00	1.00	13.7	Spleen	Adrenal gland
8573	325344	AA284265	81.88	513.18	6.27	1.00	0.00	Umbilical cord	Bone	Bone
8576	754034	AA478058	22.60	236.44	10.42	4.00	0.00	Pool	LID not found	LID not found
8582	806585	AA458480	16.39	149.75	9.13	2.00	0.00	245.05	Pooled	Ovary
8584	345553	W73089	16.43	129.72	7.90	2.00	0.00	143.02	Pooled	Whole embryo
8588	430235	AA010223	6.40	60.32	9.43	2.00	0.00	553.7	Placenta	Lung
8591	428507	AA004525	28.87	227.26	7.87	1.00	1.00	235.13	Liver	Pool
8596	234955	H37628	48.67	495.61	10.18	3.00	3.00	Pool	LID not found	Other
8602	134997	R31783	6.65	43.28	6.51	3.00	0.00	Pooled	Eye	Placenta
8604	129320	R12879	13.13	89.28	6.80	1.00	0.00	443.85	Pool	LID not found
8605	325169	AA284261	7.09	48.78	6.88	2.00	0.00	-	-	LID not found
8606	324220	AA284184	285.30	3165.68	10.69	2.00	0.00	-	-	Testis
8610	252491	H87459	8.66	118.24	17.74	2.00	0.00	237.77	Adipose	Placenta

Table 2A

8613	490023	AA114966	8.14	54.87	8.74	0.00	1.00
8615	365317	AA008593	109.91	1055.38	9.60	2.00	0.00
8624	855788	AA684040	79.37	1511.88	19.05	0.00	1.00
8625	259066	N32832	47.84	282.39	5.90	1.00	2.00
8627	324946	AA284281	73.80	690.33	9.35	2.00	0.00
8628	233277	H77494	35.54	334.54	9.41	3.00	3.00
8630	810446	AA457115	17.82	123.40	6.93	2.00	0.00
8631	430188	AA010188	3.28	53.80	16.41	4.00	0.00
8632	454440	AA677308	5.90	46.40	7.87	0.00	1.00
8633	110903	T80448	97.42	587.87	6.04	2.00	0.00
8635	262023	H98683	109.20	1285.23	11.77	2.00	1.00
8636	289060	N63598	39.43	287.24	6.78	0.00	2.00
8639	430320	AA010408	30.29	330.73	10.92	4.00	0.00
8645	460584	AA700419	8.79	48.75	5.32	1.00	1.00
8654	433111	AA680407	13.60	91.23	6.71	1.00	0.00
8655	731433	AA412217	17.36	108.73	8.15	0.00	1.00
8656	85804	T72086	7.62	39.98	5.24	1.00	0.00
8658	770388	AA430865	5.21	162.53	35.06	22.00	1.00
8675	212188	H68848	39.11	222.54	5.69	1.00	1.00
8678	39843	R53527	7.40	45.85	6.20	1.00	0.00
8691	178825	H49511	2.67	90.53	33.88	2.00	0.00
8693	856464	AA630794	115.38	890.52	5.98	0.00	2.00
8694	68818	T53431	11.34	153.87	13.55	5.00	0.00
8698	84695	T74586	14.15	75.78	5.36	1.00	0.00
8699	454872	AA677185	22.04	160.52	7.28	6.00	0.00
8706	725977	AA282226	50.62	1877.11	37.08	9.00	0.00
8710	197520	H52110	28.13	168.15	6.44	1.00	0.00
8713	859807	AA688527	78.02	778.78	9.86	0.00	3.00
8717	884980	AA630016	3.72	19.19	5.17	0.00	1.00
8720	52186	H24352	2.22	59.44	26.81	0.00	1.00
8727	34400	R44353	16.40	120.24	7.33	3.00	0.00
8731	785321	AA454175	25.52	140.36	5.50	1.00	0.00
8743	843048	AA488413	11.07	63.20	5.71	1.00	0.00
8745	49275	H18701	3.65	23.12	6.33	0.00	1.00
8747	487165	AA045074	11.52	64.17	5.57	2.00	0.00
8750	897427	AA489470	10.87	66.03	8.19	2.00	0.00
8753	50613	H17511	1.98	17.07	8.71	3.00	0.00
8757	46471	H09776	153.24	952.89	6.22	1.00	0.00
8765	52191	H24355	11.70	70.59	6.03	3.00	0.00
8768	625623	AA186605	21.84	148.76	6.81	0.00	1.00
8769	47225	H11270	22.09	228.61	10.35	4.00	5.00
8773	46907	H10204	2.46	16.77	6.81	1.00	0.00
8776	50593	H17825	160.22	1147.64	7.16	2.00	5.00
8778	291818	N87816	155.92	928.81	5.94	1.00	0.00
8786	858422	AA686180	334.58	2292.62	6.85	1.00	3.00
8790	416711	W86608	7.44	37.62	5.08	0.00	1.00
8791	798357	AA458139	4.11	42.44	10.32	2.00	0.00
8792	51542	H20757	341.33	3089.48	8.89	3.00	0.00
8794	841070	AA488761	93.93	878.52	9.36	0.00	1.00
8796	72083	T52375	42.52	392.08	9.22	1.00	3.00
8797	49240	H15436	21.21	121.09	5.71	1.00	0.00
8800	51485	H24018	2.09	12.35	5.92	1.00	0.00
8809	47149	H10893	56.98	302.36	5.31	0.00	2.00
8812	73758	T54643	13.42	138.78	10.19	3.00	0.00
8817	46187	H09245	2.84	19.93	7.03	1.00	3.00

Table 2A

8618	80336	T65736	47.65	345.61	7.26	1.00	0.00	1	545.68	Thymus	Spleen	Adipose
8619	810239	AA464709	53.40	442.34	8.28	2.00	2.00			CNS	Ovary	Prostate
8620	73787	T54673	92.77	901.06	9.71	3.00	3.00	15	227.19	Spleen	LID not found	Other
8622	133454	R27457	59.50	389.84	6.55	0.00	1.00	12	205.02	Spleen	Whole embryo	Prostate
8623	1032048	AA610040	93.15	790.01	8.48	2.00	0.00	15	215.11			
8624	52917	H29231	183.04	1622.27	8.68	3.00	2.00			Colon	CNS	Whole embryo
8629	50508	H17508	88.19	1092.05	12.38	2.00	0.00			Pool	LID not found	Other
8642	415447	W80361	36.08	244.91	6.79	3.00	1.00			CNS	LID not found	Other
8644	276920	N39448	12.98	79.23	6.11	1.00	0.00			CNS	Heart	Brain
8663	347772	W81603	28.40	167.51	6.60	1.00	1.00			Kidney	LID not found	Other
8664	428011	AA004719	66.75	354.07	5.30	1.00	0.00	16	113.68	Pool	LID not found	Other
8671	415229	W81879	4.29	180.46	42.05	9.00	2.00	1	334.46	Kidney	Colon	Pool
8675	810459	AA457138	9.49	113.26	11.94	2.00	1.00	10	188.88	Pooled	Eye	Kidney
8679	203888	H56840	24.70	153.84	6.23	1.00	0.00					
8683	240509	H90767	85.64	499.68	5.83	3.00	0.00	11	254.48	Pool	Brain	LID not found
8691	416627	W86630	23.52	121.21	5.15	1.00	0.00	2	647.14	CNS	Liver	Eye
8690	299465	N71080	109.26	585.83	5.39	1.00	0.00	18	-10.98	Prostate	LID not found	
8903	204598	H57060	2.21	12.32	5.57	1.00	0.00			Pool	LID not found	Other
8911	204661	H57130	3.26	51.45	15.78	3.00	0.00	5	455.8	Cervix	Pool	LID not found
8915	345081	W74802	2.11	18.88	8.93	7.00	0.00	8	438.5	Stomach	Breast	Heart
8938	809893	AA454881	146.42	1591.08	10.87	2.00	0.00			Ovary	Whole embryo	Lung
8960	269608	N59270	10.72	67.63	6.31	1.00	0.00	4	615.42	Synovial mem	LID not found	
8963	771128	AA429398	4.51	33.08	7.33	1.00	0.00	8	-14.62	Ovary	Uterus	Brain
8964	361840	W82514	63.81	526.49	8.25	4.00	0.00	5	570.71	Eye	Pool	LID not found
8965	488054	AA053286	20.73	115.22	5.56	2.00	0.00	2	462.78	Uterus	Breast	Tonsil
8967	185347	R89225	157.87	1020.30	6.46	2.00	0.00	9	19.96			
8969	418435	W83024	90.21	570.88	6.33	2.00	0.00	X	137.28	Pool	Brain	LID not found
8982	782822	AA447553	7.87	40.61	5.16	1.00	0.00			Thyroid	Foreskin	Testis
8984	280035	N30372	24.74	146.03	5.90	1.00	0.00	7	591.02	Tonsil	Blood	Germ Cell
8986	342008	W80057	42.21	307.87	7.29	1.00	0.00	17	307.36	Esophagus	Synovial mem	Prostate
8988	882459	AA876404	10.36	74.61	7.20	2.00	0.00	5	491.07	Cervix	Ear	Uterus
8993	415085	W93382	11.55	67.99	5.89	1.00	0.00			Prostate	Pool	LID not found
8996	782840	AA448271	36.34	234.14	6.44	2.00	0.00			Testis	LID not found	Other
9000	145112	R77293	141.93	818.52	5.75	0.00	1.00	12	430.39	Epididymis	Liver	Stomach
9001	771301	AA443837	4.45	77.86	17.50	8.00	0.00			Ovary	Heart	Germ Cell
9003	156046	R72434	1.51	11.31	7.48	0.00	0.00	4	436.32	Gall bladder	Muscle	Breast
9006	785309	AA454180	4.85	93.46	19.28	2.00	0.00			Kidney	Testis	Lung
9012	230862	R85962	2.82	21.57	7.65	1.00	2.00	11	117.56	Eye	Brain	Testis
9013	488149	AA058711	4.88	32.59	6.68	1.00	0.00	X	289.73	Neural	Colon	LID not found
9014	809503	AA484562	7.23	687.41	95.10	10.00	0.00				Stomach	Adrenal gland
9016	754408	AA436167	17.44	115.28	6.61	1.00	0.00			Kidney	Pool	LID not found
9023	428211	AA007283	141.05	943.35	6.72	4.00	0.00	12	242.45	Muscle	CNS	
9024	377871	AA055979	3.87	78.47	20.28	1.00	0.00			Pooled	Ovary	Testis
9025	755861	AA486339	5.07	392.62	77.40	7.00	0.00			CNS	Brain	LID not found
9028	51842	H24308	1.83	8.58	5.25	1.00	0.00	1	553.7	Skin	Uterus	Pancreas
9029	293950	N84014	51.45	351.95	6.84	0.00	1.00	5	-13.58	Larynx	Blood	Brain
9047	45587	H07934	16.31	81.63	5.00	1.00	0.00	12	45.2	Blood	Uterus	Pancreas
9053	40871	R56251	6.28	41.81	6.70	0.00	1.00	1	165.69	Small intestine	Skin	Pancreas
9054	46166	H08076	2.42	23.15	8.58	2.00	0.00	3	331.21	Pancreas	Foreskin	Brain
9059	51210	H19227	3.88	31.82	8.20	3.00	0.00	1	674.5	Pooled	Tonsil	Heart
9060	49555	H15089	103.13	641.75	6.22	2.00	4.00	12	247.55	Neural	Breast	Eye
9067	731308	AA418759	82.93	473.20	5.71	0.00	1.00	21	266.24	Esophagus	Cervix	Liver
9070	769837	AA430387	26.28	280.51	10.87	1.00	0.00	1	576.82			
9071	45417	H08720	5.72	62.20	10.67	1.00	0.00	11	127.97	Skin	Aorta	Pool
9075	48410	H15885	80.75	718.01	8.89	1.00	3.00					

[illegible]



Table 2A

9369	252278	H87153	68.35	491.18	7.40	3.00	0.00	18	-13.12	Peripheral ner	Thymus	Head and neck
9374	770518	AA434187	42.92	283.01	6.13	1.00	0.00	7	424.26	Bone	Pooled	Blood
9375	375882	AA032221	21.05	124.44	5.91	0.00	1.00			Larynx	Esophagus	Placenta
9376	325928	AA074511	19.03	118.33	6.22	1.00	0.00	2	71.28	Muscle	Fore skin	Pool
9377	271926	N35250	414.27	3017.05	7.28	2.00	0.00	2	118.93			
9385	272658	N32261	280.66	1598.25	5.70	2.00	0.00	15	227.19			
9386	415089	W93378	18.03	98.16	5.33	1.00	0.00	1	695.02	Pooled	Testis	Brain
9389	340842	W56753	48.32	285.07	5.88	2.00	1.00	3	418.03	Tonsil	Heart	Brain
9391	502561	AA157017	48.04	288.22	6.00	1.00	0.00	1	56.78	Uterus	Brain	Pool
9393	201609	R88003	117.56	1009.81	8.59	2.00	2.00	17	452.42	Pool	LID not found	Other
9396	283013	N63768	0.58	10.04	17.18	1.00	0.00			Pool	LID not found	Other
9407	811142	AA485731	60.45	487.38	8.06	2.00	0.00			Spleen	-	Lung
9414	153541	R48320	473.01	3242.42	6.85	1.00	1.00	2	720.63	Placenta	Aorta	Uterus
9415	34405	R44357	13.08	69.07	8.21	0.00	1.00	X	117.87	Breast	Brain	Testis
9416	23903	R39520	5.54	28.93	5.28	1.00	0.00	7	675.52	Codon	Parathyroid	Pool
9421	811046	AA485427	18.58	145.68	8.78	7.00	0.00	14	278.45	Synovial mem	Stomach	Ear
9423	1034776	AA821535	11.78	60.40	5.12	1.00	0.00	11	36.84	Esophagus	Adrenal gland	Ovary
9424	40881	R36055	5.37	52.66	9.81	1.00	0.00	11	47.55	Stomach	Spleen	Germ Cell
9425	377682	AA056013	42.85	533.24	12.45	1.00	0.00	12	345.1	Brain	Prostate	Testis
9436	73438	T55407	77.48	580.28	7.23	0.00	3.00	17	238.8	Breast	Heart	Blood
9439	47763	H11938	5.83	33.64	5.77	0.00	1.00	17	597.27	Unbical cord	Spleen	Whole embryo
9456	22786	R38835	43.37	256.70	5.88	0.00	2.00	12	66.6	Neural	Brain	LID not found
9457	854284	AA668726	4.70	28.77	5.70	1.00	0.00	11	597.27	Thymus	Eye	Codon
9460	73608	T65704	19.33	98.88	6.44	0.00	1.00	12	238.8	Parathyroid	Lymph	Germ Cell
9461	344243	W68906	26.78	172.26	6.44	0.00	1.00	1	597.27	Skin	Spleen	Pooled
9464	33026	R44498	4.97	51.88	10.39	0.00	1.00	11	66.6	Heart	Testis	Heart
9470	1031078	AA810088	17.18	132.19	7.89	0.00	0.00	5	504.31	Blood	Codon	Whole embryo
9472	41643	R52882	30.49	303.67	8.86	3.00	5.00	X	278.32	Brain	LID not found	Other
9476	73241	T58007	1.61	10.24	6.37	1.00	0.00	6	542.05	CNS	Pooled	Spleen
9478	868332	AA834028	53.84	1075.17	19.97	9.00	0.00	1	572.03	Smooth musc	Nose	Esophagus
9480	23114	R38852	43.58	381.56	8.76	2.00	5.00	5	139.45	Brain	LID not found	Other
9483	214189	H90431	6.85	111.01	18.68	0.00	2.00	14	58.23	Pooled	Prostate	Breast
9487	33076	R44048	33.35	239.68	7.19	1.00	1.00	5	454.26	Brain	LID not found	Other
9488	51395	H18415	24.31	422.45	17.38	5.00	3.00	12	58.23	Pool	LID not found	Other
9493	415086	W93370	7.08	82.56	11.84	3.00	3.00	17	89.67	Colon	Heart	Brain
9497	433481	AA698573	2.85	89.93	13.51	0.00	1.00	18	182.99	Parathyroid	Fore skin	Lung
9499	377252	AA055350	18.48	222.83	13.51	0.00	1.00	14	278.24	Brain	LID not found	Other
9501	281580	H98684	5.81	83.37	14.85	0.00	1.00	8	355.78	Brain	Adipose	Eye
9502	253009	H88540	311.21	4461.26	14.34	7.00	1.00	1	740.89	Synovial mem	Lymph	Pancreas
9509	50018	H16781	87.14	525.84	6.03	2.00	1.00	15	186.48	Thyroid	Lymph	CNS
9512	890076	AA598781	50.55	444.59	8.80	0.00	1.00	2	694.44	Brain	LID not found	Other
9513	45845	H08753	181.35	1279.28	7.05	1.00	0.00	14	278.24	Spleen	Tonsil	Pool
9514	511091	AA080256	58.55	858.15	15.14	3.00	0.00	19	290.72	Stomach	Pancreas	Cervix
9515	32186	R42871	30.23	192.99	6.38	0.00	3.00	14	55.7	Tonsil	Whole embryo	Brain
9531	71825	T52331	34.64	208.27	6.01	1.00	0.00	17	347.4	Small intestine	LID not found	Other
9537	49944	H28215	50.41	278.75	5.49	0.00	1.00	20	192.65	CNS	Cervix	Germ Cell
9538	276237	R64175	15.97	85.23	5.96	0.00	2.00	9	32.97	Spleen	Whole embryo	Brain
9539	50883	H18428	58.66	337.10	6.76	0.00	1.00	9	419.74	Bone	Colon	Prostate
9540	80699	T57648	9.37	50.30	5.24	1.00	0.00	9	419.74	Bone	Colon	Prostate
9544	730633	AA412738	32.78	171.68	5.24	1.00	0.00	9	419.74	Bone	Colon	Prostate
9552	887595	AA498887	8.31	61.69	6.63	1.00	0.00	9	419.74	Bone	Colon	Prostate
9553	47188	H10403	5.85	40.69	6.66	0.00	1.00	9	419.74	Bone	Colon	Prostate
9555	511096	AA088458	24.54	155.28	6.33	0.00	3.00	9	419.74	Bone	Colon	Prostate
9558	69360	T58848	134.08	1459.02	10.88	2.00	4.00			Placenta	LID not found	Other

Table 2A

9573	40631	H29257	196.31	1278.87	6.51	3.00	2.00	Whole embryo/Brain	Tonsil
9575	249755	H85476	88.89	455.14	5.12	1.00	0.00	Muscle	Pancreas
9576	530237	AA111978	101.13	101.13	5.06	1.00	0.00	LID not found	Other
9579	428544	AA004868	38.91	211.86	5.44	0.00	2.00	Pool	Whole embryo
9582	593537	AA165410	5.71	110.06	19.28	2.00	1.00	Ovary	Other
9587	428333	AA007502	27.49	174.35	6.34	2.00	2.00	Pool	LID not found
9591	587398	AA130351	128.23	778.20	6.03	2.00	0.00	Colon	Other
9597	46508	H09143	11.39	120.19	10.55	4.00	2.00	LID not found	Other
9599	346889	H78834	1.87	18.22	9.74	1.00	1.00	Heart	LID not found
9599	416095	H85890	63.88	532.08	8.33	3.00	2.00	Pool	LID not found
9607	204442	H56000	17.82	101.92	5.72	4.00	0.00	LID not found	Other
9608	248232	H58473	153.16	1000.45	6.53	3.00	0.00	Pool	LID not found
9612	283268	N45301	9.75	57.26	5.87	0.00	1.00	CNS	Other
9615	428529	AA004846	3.43	17.42	5.09	1.00	0.00	Parathyroid	Prostate
9616	244305	N54793	4.23	23.65	5.59	1.00	0.00	Pool	Prostate
9630	359135	AA010128	10.84	75.56	6.97	0.00	1.00	Bone	Colon
9646	489220	AA058734	4.38	32.29	7.41	2.00	0.00	Pool	Uterus
9651	384352	AA022498	4.13	35.98	8.70	1.00	0.00	Parathyroid	Heart
9654	357298	H93688	10.22	61.48	6.02	1.00	0.00	Prostate	Heart
9655	283919	N50797	27.04	162.73	8.02	2.00	0.00	Synovial mem Cervix	Other
9659	191868	H40351	1.36	8.01	5.89	1.00	0.00	LID not found	Other
9680	277747	N46086	4.19	96.43	23.02	2.00	0.00	CNS	Pool
9688	279164	N46321	44.10	358.03	7.82	1.00	0.00	Umbilical cord Tonsil	Other
9672	247381	N58022	2.74	13.89	5.06	1.00	0.00	Uterus	Placenta
9682	124895	R06123	25.73	179.55	6.88	3.00	0.00	Pool	LID not found
9688	294281	N84428	2.54	16.64	6.56	1.00	0.00	LID not found	Other
9691	782688	AA47593	25.28	162.18	6.42	1.00	0.00	Bone	Breast
9694	283995	N53378	8.88	44.55	5.03	1.00	0.00	Testis	Liver
9695	207828	H60560	6.24	37.40	6.00	2.00	0.00	Eye	Tonsil
9696	489189	AA056580	18.08	105.60	5.54	0.00	1.00	LID not found	Other
9700	205239	H60824	4.17	27.32	6.56	3.00	0.00	Skin	CNS
9701	303180	N92764	6.00	35.21	5.87	1.00	0.00	Thyroid	Heart
9708	755301	AA496380	5.83	33.79	6.00	3.00	0.00	Muscle	Eye
9711	292223	N62464	96.86	533.08	5.50	1.00	0.00	Lymph	Breast
9712	22895	R38840	3.92	172.21	43.91	3.00	0.00	Testis	Blood
9714	427750	AA001897	2.61	14.52	5.57	1.00	0.00	Fore skin	Pool
9715	120929	T98107	4.41	28.74	6.74	1.00	0.00	Brain	LID not found
9716	434833	AA703141	6.24	151.91	24.34	0.00	5.00	Pool	LID not found
9724	280970	N50854	21.51	144.48	6.72	1.00	2.00	Germ Cell	CNS
9728	298448	N74623	16.41	238.92	145.45	3.00	0.00	Head and nec Esophagus	Skin
9733	359009	H92134	1.85	21.18	11.44	5.00	0.00	Eye	Pancreas
9737	795315	AA454172	96.40	574.13	5.96	2.00	0.00	Testis	Testis
9750	785564	AA458674	5.70	28.83	5.08	1.00	0.00	Pool	Pool
9753	428338	AA005153	5.38	48.88	9.07	7.00	0.00	Blood	Parathyroid
9758	795614	AA460006	1.08	9.91	9.34	1.00	0.00	Testis	LID not found
9759	212347	H88286	28.95	237.69	8.82	1.00	0.00	Pool	Brain
9761	428832	AA004528	43.78	246.10	5.82	2.00	0.00	Pool	LID not found
9767	243068	N39542	24.68	125.89	5.12	1.00	0.00	Lung	Synovial membrane
9778	375857	AA037810	8.39	45.30	5.40	0.00	1.00	Head and nec Cervix	Brain
9781	504431	AA151245	1.50	20.77	13.83	2.00	0.00	Lung	Uterus
9784	755589	AA419251	64.15	1420.22	22.14	7.00	0.00	Nose	Skin
9787	810700	AA457688	26.28	141.61	5.80	2.00	0.00	Ovary	Adipose
9788	155768	R72097	1.88	32.64	17.33	5.00	0.00	Esophagus	LID not found
9789	785274	AA454012	4.28	24.40	5.70	1.00	0.00	Neural	Stomach
									Breast
									Ovary

Table 2A

9790	810205	AA464518	3.45	61.13	17.72	1.00	0.00	8	404.08	Blood	Heart	Ovary
9792	289496	N63988	10.03	428.48	42.72	3.00	0.00	10	421.81	Smooth musc	Stomach	CNS
9804	51508	H18953	42.30	284.50	6.73	0.00	4.00			Brain	Pool	LID not found
9805	590774	AA157499	1.42	7.33	5.16	1.00	0.00	3	726.84	Parathyroid	Pancreas	Prostate
9809	773478	AA427891	19.92	106.31	5.34	1.00	0.00			Tonsil	Blood	Whole embryo
9813	270917	N32514	0.63	16.96	26.81	0.00	0.00	8	165.67	Synovial mem	Neural	Heart
9818	364934	AA025275	33.68	219.53	6.48	1.00	0.00	16	398.69	CNS	Pancreas	Parathyroid
9820	51218	H18472	8.64	48.69	5.75	0.00	1.00	10	363.47	Marrow	Ovary	Pooled
9824	47378	H11036	13.72	97.07	7.07	1.00	2.00	8	434.48	Brain	Whole embryo	Kidney
9828	33837	R45939	12.69	101.37	7.89	2.00	0.00	4	439.11		Adrenal gland	Lung
9830	713974	AA324954	12.28	62.28	5.07	1.00	0.00	5	578.76			
9831	52577	H29566	161.57	1246.09	7.71	3.00	5.00	9	358.01			
9832	32150	R43352	5.81	56.89	5.18	0.00	2.00	11	131.57	CNS	Brain	Pool
9837	761047	AA446482	26.06	135.05	5.18	0.00	1.00			Adrenal gland	Testis	Germ Cell
9844	22600	T87235	1.33	14.61	10.98	2.00	0.00	12	80.95	Brain	LID not found	Other
9849	868530	AA664389	256.71	1388.48	5.33	1.00	0.00	3	472.27	Brain	Blood	Prostate
9855	50250	H17800	5.59	32.93	5.89	2.00	0.00	13	131.31	Stomach	Head and neck	Ear
9856	51608	H18936	2.37	30.59	12.90	2.00	2.00	3	697.57	Brain	LID not found	Other
9860	34560	R44193	2.61	19.31	7.38	1.00	0.00	13	155.48	Pancreas	Lung	Brain
9861	303109	N90783	12.95	96.96	7.49	1.00	0.00	13	147.98	Whole embryo	Pancreas	CNS
9862	433155	AA680136	4.11	29.01	7.08	0.00	0.00	1	611.54	Gall bladder	Pool	Pancreas
9864	32496	R43486	4.01	28.80	7.17	2.00	0.00			Testis	Pool	Brain
9865	287745	N62244	21.29	162.81	7.65	2.00	0.00	1	740.99	Esophagus	Adrenal gland	Aorta
9870	432851	AA700556	25.79	140.68	5.45	0.00	1.00				LID not found	Other
9871	50132	H17055	2.88	15.79	5.90	1.00	0.00	13	204.51	Brain	Pool	LID not found
9872	49554	H15087	4.37	23.62	5.40	1.00	0.00	2	97.87	Brain	Brain	Whole embryo
9875	51774	H23212	14.91	115.07	7.72	2.00	0.00	5	-9.1	Foreskin	Brain	LID not found
9876	24915	R39066	3.46	24.72	7.10	1.00	0.00			Pool	Cervix	Germ Cell
9881	241738	H91691	35.14	224.01	6.37	0.00	1.00	12	244.17	Skin	CNS	Germ Cell
9884	72003	T52330	1.99	14.07	7.07	3.00	0.00	1	538.57	Eye		Lymph
9885	598725	AA129171	9.62	48.18	5.01	1.00	0.00	11	236.57	Brain	LID not found	Other
9888	40649	R56397	1.75	10.41	5.93	1.00	0.00	18	329.18	Brain	LID not found	Other
9892	33510	R43896	2.58	16.23	6.28	2.00	0.00	7	485.65	Neural	Foreskin	Brain
9895	951125	AA620556	115.42	950.64	8.24	0.00	1.00	6	21.28	Umbilical cord	Muscle	Uterus
9904	83156	T68113	243.52	1425.87	5.86	0.00	3.00			Liver	Uterus	LID not found
9909	45300	H04203	1.66	49.89	26.79	0.00	1.00	11	97.53	CNS	Brain	LID not found
9913	23586	R38369	8.07	73.87	9.16	1.00	3.00	2	675.85	Germ Cell	Lung	
9921	31818	R41730	24.68	136.75	5.54	1.00	0.00	2	228.05	Brain	Whole embryo	Lung
9923	743224	AA400229	13.12	71.96	5.48	1.00	0.00	5	654.42	Smooth musc	Testis	Tonsil
9925	47234	H10939	2.22	47.37	21.35	5.00	1.00	5	584.16	Larynx	Brain	Parathyroid
9928	796078	AA460353	0.84	9.76	11.68	1.00	0.00	9	383.2	Neural	Synovial mem	Blood
9935	950668	AA608367	202.83	1140.65	5.82	1.00	0.00			Lymph node	Foreskin	Bone
9956	26259	R20547	4.57	55.33	12.10	2.00	0.00			Brain	LID not found	Other
9957	51496	H18927	2.25	18.69	8.31	2.00	0.00			CNS	Brain	Aorta
9958	346671	W74646	5.54	30.82	5.56	1.00	0.00	16	194.85	Lymph node	Synovial membrane	
9960	25384	R12808	50.97	391.93	7.89	3.00	0.00					
9966	282720	N50079	55.00	288.40	5.43	0.00	1.00	13	120.56	Thymus	Ear	Adrenal gland
9970	841226	AA487115	48.62	270.05	5.79	1.00	0.00	1	639.63	Adipose	CNS	Muscle
9973	46561	H09757	13.82	197.65	14.30	1.00	0.00	8	91.28	Foreskin	Whole embryo	Germ Cell
9874	849988	AA600214	178.13	1782.73	9.95	1.00	0.00	8	442.4	Synovial mem	Skin	Cervix
9980	32683	R43543	3.31	22.66	6.85	0.00	1.00			Whole embryo	Brain	LID not found
9981	50879	H18424	6.40	52.52	8.21	0.00	1.00	17	332.34	Skin	Brain	Whole embryo
9988	795446	AA436816	1.45	28.38	18.61	5.00	0.00	5	39.72	Germ Cell	Testis	Eye
9991	284100	N53421	6.62	112.86	17.06	1.00	2.00	16	-10.98	Smooth musc	CNS	Pancreas

Table 2A

9996	320201	W15542	105.30	672.71	6.39	2.00	0.00	1	98.48	Parathyroid	Eye	Prostate
10003	810454	AA457118	39.05	293.81	7.52	4.00	0.00	21	217.43	Larynx	Cervix	Skin
10014	124229	R02329	5.72	54.41	9.52	2.00	0.00			Pool	LID not found	Other
10019	303049	N91589	4.53	26.38	5.82	1.00	0.00	16	22.69	Colon	Breast	Lung
10020	270277	N33555	47.05	279.16	5.93	1.00	0.00	3	141.89	Foreskin	LID not found	Other
10022	415806	W84774	3.40	17.48	5.15	1.00	0.00	3	287.63	Pool	Colon	LID not found
10027	503725	AA131530	160.27	1134.59	7.08	2.00	0.00	X	298.2	CNS	Bone	Breast
10030	120878	T95650	51.78	449.02	8.67	2.00	2.00	3	198.24	Breast	Brain	Pool
10034	308216	N92404	13.53	114.47	8.46	3.00	3.00	4	420.51	Whole embryo	Pool	Heart
10038	428299	AA007370	9.85	97.80	9.93	1.00	0.00	19	218.02	Stomach	CNS	Heart
10042	113193	T83864	49.05	282.84	5.77	2.00	2.00			CNS	Placenta	Pool
10047	277714	N49377	2.25	14.52	6.46	1.00	0.00	4	579.1	Colon	Heart	Uterus
10048	343387	W87292	4.07	23.21	5.71	1.00	0.00	8	404	Bone	Foreskin	Parathyroid
10052	321800	W37447	16.58	141.89	8.56	3.00	0.00	10	126.87	Pool	LID not found	Other
10058	429108	AA004786	4.59	27.01	5.89	1.00	0.00	10	45.1	Foreskin	Testis	Pool
10060	273677	N36172	5.59	38.72	6.92	1.00	0.00	10		Foreskin	Colon	Pool
10064	286250	N26714	20.94	122.51	5.85	1.00	0.00	22	118.23			
10078	120824	T95909	141.55	1082.71	7.65	3.00	0.00			Ear	Prostate	Heart
10081	327337	W02102	39.91	230.36	5.77	0.00	1.00	1	576.51	Ear	Brain	Eye
10084	345743	W72668	20.15	178.23	8.85	5.00	0.00	21	154.34	Neural	Foreskin	Kidney
10089	261163	H98201	286.70	1632.70	5.69	1.00	0.00	3	164.31	Adipose	Colon	Heart
10095	377363	AA055052	2.55	19.63	7.71	2.00	0.00	20	270.8	Epithelium	Neural	Pancreas
10096	856796	AA689545	20.31	122.45	6.03	1.00	0.00	14	14.91	CNS	Foreskin	Whole embryo
10099	277042	N39572	4.73	47.81	10.10	4.00	0.00			Thymus	LID not found	Other
10105	273652	N36994	108.66	652.50	6.00	1.00	2.00	3	69.98	Foreskin	Eye	CNS
10112	384610	AA022561	38.68	348.61	9.51	0.00	5.00	8	116.81	Foreskin	Brain	Eye
10113	291417	N72307	213.91	1366.82	8.39	2.00	0.00			CNS	Eye	Breast
10115	286748	N62487	9.73	95.69	9.83	1.00	0.00				Heart	LID not found
10117	342208	W63785	33.89	204.61	6.04	0.00	3.00					
10120	433350	AA700604	34.17	200.17	5.86	1.00	1.00	2	120.34		Pool	LID not found
10124	257248	N28806	8.42	92.76	11.01	1.00	0.00	16	354.33	Placenta	LID not found	Other
10128	139840	R62364	3.38	27.21	8.05	1.00	0.00	5	289.88	CNS	LID not found	Other
10129	280387	N49278	11.39	75.96	6.67	0.00	4.00	3	52.66	Brain	Liver	Breast
10136	177967	H46254	4.23	21.50	5.08	1.00	0.00	7	264.9	Skin	Liver	Breast
10139	278091	N46240	31.92	168.49	5.22	0.00	1.00			Ear	Pool	Pool
10140	418539	W086423	102.91	639.40	6.21	2.00	0.00			Tonsil	Breast	Germ Cell
10143	504623	AA149228	14.33	95.93	6.70	1.00	0.00	2	140.92	Foreskin	Aorta	Umbilical cord
10144	53092	R15785	45.48	390.24	8.58	0.00	1.00			Heart	LID not found	Other
10149	365575	AA009484	1.29	23.42	18.22	3.00	0.00			Heart	Kidney	Pancreas
10157	376080	AA040387	6.21	39.11	6.29	1.00	0.00			Foreskin	Eye	Heart
10162	272327	N32199	180.45	1654.38	10.28	2.00	0.00	1	210.4	Heart	Aorta	Whole embryo
10165	345051	W72803	13.67	64.70	5.20	0.00	1.00	19	250.87	Uterus	Bone	Whole embryo
10166	503096	AA151480	9.10	48.85	5.37	1.00	0.00			Ovary	Lung	Brain
10167	809674	AA454689	79.24	704.02	8.89	3.00	0.00			Heart	Whole embryo	Lung
10173	377168	AA055163	4.11	21.86	5.32	1.00	0.00	8	323.97	Ear	Tonsil	Brain
10174	758298	AA040437	8.66	43.97	5.08	1.00	0.00	17	41.12	Neural	Umbilical cord	Muscle
10177	531937	AA113881	10.03	58.80	5.86	1.00	0.00			Pancreas	Colon	Kidney
10181	121551	T97710	3.90	71.56	18.34	7.00	1.00	6	102			
10181	32632	R43088	8.27	47.33	7.55	0.00	3.00	2	576.38	Aorta	Brain	LID not found
10192	41793	R59200	3.00	32.46	10.81	1.00	0.00			Kidney	Pool	Brain
10195	363055	AA019320	6.52	36.12	5.54	0.00	1.00			Umbilical cord	Spleen	Foreskin
10196	80692	T57834	158.02	1113.27	7.14	1.00	3.00	6	650.09	Brain	LID not found	Other
10199	46196	H08243	382.88	1970.93	5.15	0.00	1.00	3	684.31	CNS	Spleen	Tonsil
10204	74713	T57349	11.63	106.89	9.20	0.00	0.00			Breast	Eye	Lung
10205	152433	R46202	22.18	149.69	6.75	0.00	2.00					

Table 2A

10207	73222	T57221	53.18	332.42	6.25	1.00	3.00	Pool	LID not found
10210	120600	T95113	8.28	230.57	27.84	2.00	0.00	Lymph	.
10216	32825	R43646	2.42	19.79	8.17	1.00	0.00	Brain	LID not found
10223	742635	AA400262	38.08	352.82	9.27	2.00	3.00	Colon	LID not found
10224	23216	R38089	3.59	21.50	5.68	2.00	0.00	Pool	LID not found
10225	428448	AA007889	21.09	112.40	5.33	1.00	0.00	Pool	.
10232	46380	H09789	29.89	318.81	10.67	5.00	5.00	Brain	LID not found
10235	868304	AA634008	546.62	4715.07	8.63	1.00	0.00	Small intestine	Bone
10240	41405	R56146	2.57	19.30	7.52	3.00	0.00	LID not found	Other
10246	505491	AA156461	82.10	325.37	5.24	0.00	1.00	Small intestine	Neural
10248	23554	R38274	4.76	28.74	6.25	0.00	1.00	Eye	CNS
10253	121406	T96888	2.15	11.44	5.31	1.00	0.00	Lymph	.
10261	894718	AA629558	64.78	737.72	11.39	0.00	1.00	Colon	Aorta
10262	135085	R33031	293.58	2277.83	7.76	2.00	2.00	Adrenal gland	Colon
10264	49949	H10079	5.56	43.90	7.90	0.00	1.00	Eye	Breast
10268	432042	AA078280	30.83	217.12	7.04	0.00	1.00	Ear	Fore skin
10269	306575	N94820	83.72	505.47	6.04	1.00	2.00	Uterus	Lung
10272	23726	R39555	50.72	585.34	11.76	0.00	2.00	Lung	Uterus
10275	376358	AA041300	9.85	51.42	5.33	0.00	1.00	Placenta	Colon
10276	40352	R54797	47.70	558.85	11.87	1.00	4.00	LID not found	Other
10282	742038	AA402812	11.87	59.91	5.05	0.00	1.00	Cervix	Adrenal gland
10288	509943	AA052960	36.44	281.87	7.73	2.00	1.00	Cervix	Stomach
10291	511808	AA088430	161.98	1094.79	8.76	0.00	2.00	Adrenal gland	Colon
10293	46438	H09864	29.38	226.19	7.70	1.00	4.00	Heart	Brain
10294	32483	R43456	6.48	32.69	5.04	0.00	0.00	LID not found	Other
10297	40150	R53980	48.00	620.10	12.92	5.00	5.00	Testis	Brain
10302	22845	R43678	4.27	84.15	19.88	1.00	0.00	Brain	LID not found
10309	46238	H10681	80.09	645.84	8.06	6.00	4.00	Testis	Brain
10312	950876	AA608555	5.64	69.52	12.32	3.00	0.00	Aorta	.
10314	627251	AA191468	266.88	1984.87	7.44	0.00	2.00	Small intestine	Whole embryo
10315	795084	AA460366	28.77	228.15	7.93	0.00	4.00	Whole embryo	Adrenal gland
10316	23063	R45284	3.98	25.58	8.43	1.00	1.00	Brain	Pool
10317	47059	H11016	61.65	402.97	6.54	4.00	2.00	Parathyroid	Fore skin
10319	840576	AA487898	10.42	54.16	5.20	0.00	1.00	Lung	Prostate
10321	40038	R53446	3.63	119.59	32.94	6.00	4.00	LID not found	Other
10323	839579	AA489813	32.82	378.45	10.01	3.00	0.00	LID not found	Other
10324	41595	R59556	4.57	33.88	7.41	3.00	0.00	Blood	.
10326	429936	AA034039	66.53	343.88	5.02	1.00	0.00	LID not found	Other
10327	321751	W05416	208.05	1389.36	6.68	1.00	2.00	Parathyroid	LID not found
10329	41391	R56123	4.33	62.27	14.40	4.00	0.00	Pool	Other
10340	23822	R36546	3.60	23.55	6.53	1.00	0.00	Arteries	Breast
10342	345262	W72881	4.07	32.38	7.96	3.00	0.00	Brain	LID not found
10344	76647	T50995	1.24	7.97	6.41	1.00	0.00	Nose	LID not found
10347	839370	AA489847	3.64	20.93	5.75	0.00	1.00	Ovary	Eye
10348	23800	R38381	5.40	54.60	10.11	5.00	5.00	Fore skin	Eye
10351	781461	AA429573	14.71	89.28	6.07	0.00	3.00	LID not found	Other
10355	809600	AA458473	30.61	237.42	7.75	0.00	2.00	Whole embryo	Pool
10356	32092	R42695	1.51	13.64	9.05	2.00	1.00	Ovary	Parathyroid
10361	41825	R54109	0.85	22.97	24.09	5.00	2.00	Colon	Brain
10362	131626	R23735	1.32	9.29	7.05	0.00	1.00	CNS	Brain
10363	839807	AA489788	199.76	1642.03	8.22	2.00	2.00	Synovial mem	Placenta
10364	24237	R38018	2.36	23.90	10.13	3.00	2.00	Eye	LID not found
10366	322652	W15487	148.26	1014.81	6.84	1.00	0.00	LID not found	Other
10367	808657	AA454682	134.69	1288.39	9.42	1.00	0.00	Parathyroid	Placenta
								Fore skin	Ear

Table 2A

10368	842878	AA486410	117.10	993.99	8.49	2.00	2.00	6	95.89	Cervix	Parathyroid	Bone
10370	126513	R06748	34.79	341.73	9.82	4.00	4.00			Lymph	Kidney	Pool
10374	782851	AA448285	26.21	230.85	8.81	3.00	3.00			Head and nec	Cervix	Cervix
10395	197374	R07122	48.08	235.75	5.12	1.00	0.00	15	262.02	Muscle	Heart	Stomach
10418	126739	R07128	11.78	9.98	9.98	2.00	0.00	18	62.25	Pool	LID not found	Other
10427	185801	R09104	73.91	449.48	6.08	2.00	0.00					
10435	194689	R09928	5.58	41.33	7.40	3.00	0.00	1	629.23	Pool	LID not found	Other
10438	130078	R21408	4.33	30.97	7.15	2.00	0.00			Forebrain	Placenta	Uterus
10444	360155	AA013260	0.68	4.65	6.88	1.00	0.00			CNS	Gall bladder	Eye
10446	360065	AA053815	6.18	37.03	6.01	1.00	0.00			CNS	Heart	Testis
10448	290182	N62213	2.84	15.02	5.30	1.00	0.00			CNS	Testis	LID not found
10450	505506	AA148958	11.25	70.67	5.28	1.00	0.00	3	730.78	Uterus	Pool	LID not found
10459	418323	W00764	29.44	202.03	6.86	1.00	0.00					
10460	289287	N73705	19.56	98.66	5.04	1.00	0.00	2	184	Thyroid	Adrenal gland	Lymph
10462	795378	AA453495	13.74	208.51	15.17	4.00	0.00			Germ Cell	Placenta	Forebrain
10465	428492	AA005428	14.87	78.31	5.27	1.00	0.00	12	316.62	Pool	LID not found	Other
10467	234623	H77728	7.08	37.61	5.31	1.00	0.00	16	22.79	Blood	Pool	Eye
10468	811015	AA485377	8.43	82.30	9.77	2.00	0.00	14	199.9	Aorta	Nose	Thymus
10471	301817	N91165	1185.06	10375.29	8.91	2.00	0.00			Liver	Lung	Pool
10472	742132	AA406020	21.64	246.82	11.40	8.00	0.00	8	440.23	Esophagus	Cervix	Germ Cell
10474	376482	AA039512	117.36	633.85	5.40	1.00	0.00			Thyroid	Blood	Tonsil
10476	39820	R53935	31.54	305.80	9.69	3.00	0.00			Pool	LID not found	Other
10480	267697	N59150	118.00	882.27	7.43	1.00	0.00	21	152.58	Bone	CNS	Whole embryo
10488	588815	AA157813	237.84	1553.70	6.54	4.00	0.00			Esophagus	Head and nec	Nose
10493	365041	AA074586	19.71	133.85	6.79	0.00	1.00			Head and nec	Adrenal gland	Pool
10496	755578	AA419177	60.74	675.41	11.17	0.00	1.00	16	484.48	Cervix	Lymph	Forebrain
10501	810089	AA464983	17.49	309.86	17.71	7.00	0.00	8	561.65	Breast	Cervix	Ovary
10504	418889	W06199	20.78	143.51	6.90	0.00	1.00			Blood	Eye	Whole embryo
10511	418328	W00760	34.54	192.82	5.58	0.00	2.00	4	467.5	Pool	LID not found	Other
10512	344134	W73750	10.78	55.44	5.15	1.00	0.00	22	63.7			
10517	810230	AA464694	11.15	72.95	6.54	2.00	0.00			Adipose	Ovary	Prostate
10522	785561	AA459681	15.19	81.51	5.38	1.00	0.00			Blood	Aorta	CNS
10525	415715	W04687	3.57	17.87	5.01	1.00	0.00	2	103.51	Brain	Pool	LID not found
10530	507817	AA135868	63.53	724.53	6.97	2.00	0.00	6	299.58	CNS	Uterus	Lung
10531	481644	AA150263	6.83	38.81	5.69	0.00	1.00	4	444.77	Stomach	Pool	Uterus
10533	489805	AA102068	2.28	37.27	16.35	2.00	0.00			Parathyroid	Colon	Uterus
10537	376043	AA040265	2.32	33.75	14.52	1.00	0.00	10	403.27	Heart	Whole embryo	Prostate
10542	321482	W032511	32.32	196.83	6.08	0.00	1.00	5	-10.43	Adrenal gland	Bone	Ear
10543	229849	H70887	3.13	50.86	16.25	1.00	0.00			Ovary	Testis	Pool
10545	430291	AA010600	3.68	31.42	8.54	3.00	0.00	19	274.87	Larynx	Ovary	Pancreas
10549	810960	AA459401	8.17	560.56	68.80	13.00	1.00	16	44.86			
10550	489208	AA045658	15.81	117.88	7.41	0.00	1.00	11	183.74	Pool	LID not found	Other
10553	193311	R04491	28.39	167.78	6.38	2.00	1.00	10	335.81	Whole embryo	Pool	Umbilical cord
10554	504826	AA150777	13.30	77.50	5.83	0.00	1.00			Aorta	Muscle	Germ Cell
10557	782737	AA448003	2.98	15.25	5.11	1.00	0.00	19	250.4	Pool	Spleen	Brain
10558	429060	AA005140	9.20	55.32	6.01	1.00	0.00	7	137.5	Forebrain	CNS	Pool
10559	298486	N70208	3.43	450.59	131.58	22.00	1.00	17	48.48	Thymus	Pool	Germ Cell
10562	884655	AA629809	138.10	1266.57	9.10	0.00	1.00	20	208.92	Brain	LID not found	Other
10565	531319	AA071486	12.10	147.97	12.23	0.00	1.00	5	591.55			
10567	40873	R56044	8.18	60.13	7.35	1.00	0.00	11	67.07	CNS	Adrenal gland	Placenta
10571	46898	H08818	7.81	54.10	6.93	1.00	2.00	2	410.75	Brain	LID not found	Other
10573	744047	AA629282	18.22	114.76	5.97	0.00	1.00	15	257.33	Brain	LID not found	Other
10575	46173	H09098	15.28	1683.19	110.14	11.00	1.00					
10576	33408	R44082	87.32	604.24	6.92	4.00	3.00					
10583	47161	H10995	9.11	180.45	20.81	0.00	3.00					

Table 2A

10584	48829	H10047	15.09	138.25	9.03	2.00	3.00	7	126.26	Brain	LID not found	Other
10587	45877	H05582	6.76	62.56	9.28	1.00	0.00	7	323.97	Brain	Kidney	Lung
10588	46367	H05859	6.28	39.97	6.37	1.00	1.00	11	259.81	Brain	Pool	LID not found
10593	46266	H05078	5.33	28.84	5.41	1.00	0.00	7	141.38	CNS	Brain	Uterus
10613	430968	AA078335	6.00	234.02	39.01	1.00	0.00	3	143.12	Pool	Testis	Prostate
10614	472008	AA036881	6.82	151.10	22.17	0.00	1.00	3	488.05	Lung	Lymph	Uterus
10616	50587	H17620	70.18	627.73	8.94	4.00	5.00	2	488.05	Brain	Testis	LID not found
10620	51991	H23228	74.03	710.43	9.80	3.00	0.00	8	426.08	Ovary	Blood	Prostate
10622	223350	H85554	3.16	372.02	117.56	21.00	2.00	8	426.08	Gall bladder	Liver	Ovary
10623	46273	H05082	8.47	66.85	7.80	2.00	1.00	1	130.87	-	CNS	Germ Cell
10626	878638	AA670430	15.43	112.35	7.28	1.00	0.00	7	468.83	Aorta	CNS	Brain
10627	47788	H11728	2.73	20.43	7.47	2.00	0.00	21	229.15	Pancreas	Bone	Colon
10632	49227	H15408	6.73	65.74	9.76	2.00	0.00	14	21.38	Gall bladder	Liver	Kidney
10635	26503	R20639	4.32	23.78	5.51	0.00	0.00	2	97.77	Fore skin	Brain	Whole embryo
10639	50004	AA419229	11.98	188.01	15.69	0.00	1.00	2	471.33	Brain	Pool	CNS
10640	33715	H16733	2.07	18.33	8.83	1.00	0.00	2	97.77	Brain	Brain	Brain
10641	884783	R44078	17.44	119.79	6.87	0.00	0.00	3	91.11	Adrenal gland	Ovary	LID not found
10642	49987	AA620603	0.73	70.39	98.44	1.00	0.00	7	601.67	Spleen	Brain	Pool
10644	25838	H28734	2.00	79.04	39.51	7.00	0.00	1	286.7	Neural	Marrow	Muscle
10648	739193	R37108	24.83	385.27	15.52	4.00	5.00	20	194.08	Gall bladder	Germ Cell	Eye
10651	41720	AA421218	5.15	100.31	19.49	8.00	0.00	15	258.04	Brain	Fore skin	Lung
10652	25520	R37696	16.19	195.37	12.07	7.00	1.00	2	97.77	Synovial mem	Eye	Bone
10653	214985	H72030	11.55	132.27	11.45	0.00	1.00	19	278.4	Brain	Pool	CNS
10650	34468	R44214	82.53	428.97	5.20	2.00	0.00	17	340.31	Whole embryo	Brain	LID not found
10655	51907	H23529	70.24	819.70	11.67	1.00	5.00	13	32.28	Neural	Parathyroid	Stomach
10677	48411	H08164	6.69	40.85	6.11	2.00	0.00	19	214.41	Uterus	Brain	Pool
10688	287468	N68081	9.73	67.74	9.02	0.00	1.00	13	278.42	Brain	Parathyroid	Kidney
10711	595109	AA173926	22.05	158.85	7.20	0.00	1.00	11	162.94	Ear	Breast	Testis
10716	51604	H18934	4.03	38.90	9.89	3.00	0.00	16	276.85	Fore skin	Prostate	Blood
10718	850428	AA599085	19.84	103.13	5.20	0.00	1.00	12	24.9	Lymph	Placenta	Colon
10721	841084	AA486780	2.44	25.74	10.55	1.00	0.00	12	188.25	CNS	Adrenal gland	Testis
10735	279308	N46354	3.78	20.36	5.37	1.00	0.00	4	56.88	CNS	Gall bladder	Parathyroid
10737	773639	AA431887	15.00	133.33	8.89	2.00	0.00	14	227.31	CNS	LID not found	Other
10740	51433	H20747	2.50	41.18	16.48	1.00	0.00	17	282.8	Head and nec	Cervix	Muscle
10741	51896	H23225	4.36	47.78	10.96	0.00	3.00	17	53.59	Prostate	Brain	LID not found
10744	66747	T64919	24.18	208.97	8.64	2.00	4.00	21	229.15	Ovary	Pool	Brain
10746	585197	AA164847	5.22	82.18	15.78	6.00	0.00	21	307.17	Eye	Tonsil	Testis
10747	190915	H38221	6.80	48.34	5.43	1.00	0.00	17	109.9	Larynx	Synovial mem	Heart
10751	897582	AA486878	68.44	347.95	5.08	1.00	0.00	22	438.17	Thyroid	Breast	Testis
10756	257096	N30792	81.64	587.57	6.95	2.00	0.00	10	109.9	Spleen	CNS	Heart
10763	417976	W80660	10.83	79.73	7.28	0.00	1.00	10	529.52	Pool	LID not found	Other
10764	259072	N28356	2.96	17.38	5.87	1.00	0.00	10	57.93	Nose	Muscle	Kidney
10766	305253	N95011	22.71	235.73	10.38	2.00	0.00	7	141.87	Ear	Adipose	Prostate
10767	309081	N82895	4.10	38.47	8.89	1.00	0.00	3	167.66	Neural	Germ Cell	Pool
10772	254028	N30308	4.39	117.55	28.75	1.00	0.00	15	167.66	Ear	LID not found	Other
10778	428652	AA004321	34.20	289.72	8.47	3.00	3.00	15	167.66	Prostate	Pool	LID not found
10782	121028	T86309	5.96	57.72	9.68	3.00	3.00	13	284.56	Aorta	Uterus	Uterus
10798	429198	AA005355	9.62	60.93	6.33	1.00	0.00	19	112.49	Pool	LID not found	Other
10799	491186	AA137073	5.08	39.49	7.76	1.00	0.00	14	123.72	Fore skin	Heart	LID not found
10802	428697	AA004353	36.44	211.40	5.80	0.00	1.00	2	592.27	Pool	LID not found	Other
10806	271471	N35025	114.02	607.20	5.33	1.00	0.00	14	123.72	Pool	Heart	LID not found
10814	418081	W80067	4.60	28.13	5.68	1.00	0.00	2	592.27	Pool	LID not found	Other
10816	280022	N32623	7.60	40.24	5.29	1.00	0.00	2	592.27	Pool	LID not found	Other

Table 2A

10819	810937	AA459389	1.72	27.95	16.26	3.00	0.00	22	74.71	Ignore	Germ Cell	Bone
10820	291690	N67832	4.10	27.99	6.63	1.00	0.00	16	489.89	Forebrain	Codon	.
10826	109952	T88816	10.50	56.46	5.38	0.00	1.00			Pool	LID not found	Other
10828	491367	AA148524	111.89	637.65	5.70	0.00	1.00			Heart	Uterus	Pool
10840	258686	N30069	3.43	34.18	9.86	0.00	1.00	8	383.61	Thyroid	Germ Cell	Lymph
10842	194515	R86242	8.82	56.94	6.46	2.00	0.00	2	387.37	Stomach	Kidney	Parathyroid
10848	345638	W72692	142.07	1010.76	7.11	5.00	0.00			Heart	Pool	LID not found
10852	241996	H92985	10.44	63.17	8.05	0.00	2.00	11	253.29	CNS	Stomach	Parathyroid
10855	134942	R32334	7.00	49.81	7.12	5.00	0.00	2	603.72	Lymph	Breast	Pool
10858	196543	R91570	3.53	317.54	89.86	3.00	0.00	9	-10.76	Ovary	Pool	LID not found
10860	242001	H92874	31.78	399.27	12.56	5.00	5.00			Larynx	CNS	Eye
10867	283190	N51357	2.78	46.17	16.58	1.00	0.00	15	135.17			
10868	241948	H93050	15.74	95.38	6.06	0.00	1.00	7	481.98	Parathyroid	CNS	Prostate
10875	244796	N52554	257.16	2560.65	9.96	2.00	0.00	6	99.75			
10876	429505	AA011383	12.30	96.29	7.83	2.00	3.00	6	111.13	Thymus	Germ Cell	Parathyroid
10889	307687	N92924	5.32	31.80	5.96	0.00	1.00	2	709.74	Forebrain	Testis	Kidney
10890	272690	N36174	4.18	56.65	14.04	0.00	1.00			Bone	Forebrain	Colon
10895	782838	AA448270	4.31	37.41	8.68	1.00	0.00			Parathyroid	Forebrain	Germ Cell
10897	304886	N92502	16.99	132.95	7.83	2.00	0.00	12	42.89	Esophagus	Smooth muscle	Stomach
10899	279999	N57577	83.84	545.96	6.51	4.00	0.00	1	682.97	Thyroid	Spleen	Parathyroid
10903	770681	AA476285	5.30	39.13	7.38	5.00	0.00	7	94.82		Spleen	Prostate
10904	236155	H61758	4.09	40.08	9.81	1.00	1.00	1	554.22	Skin	Forebrain	Germ Cell
10905	320455	W04713	42.55	280.02	6.58	0.00	1.00	16	464.97	CNS	Adrenal gland	Kidney
10908	287458	N25240	28.20	224.85	7.97	0.00	1.00	7	675.72	Smooth muscle	Skin	Umbilical cord
10915	785522	AA454218	3.72	18.70	5.03	0.00	0.00	10	126.92	Small intestine	Parathyroid	Eye
10917	810496	AA457153	23.68	172.84	7.29	1.00	0.00			Blood	Lymph	Testis
10918	418004	W90705	23.64	153.71	6.50	2.00	0.00			Heart	LID not found	Other
10923	307249	N93438	21.53	123.01	5.71	1.00	1.00			Thyroid	Ovary	Muscle
10933	346699	W78169	4.32	24.22	5.60	0.00	0.00			Lung	Uterus	Tonsil
10935	810911	AA459286	48.11	381.75	7.93	10.00	1.00	10	540.8		Lung	LID not found
10937	307304	N93455	13.15	82.48	6.27	2.00	0.00			Skin	Bone	Ear
10939	293895	N86025	60.68	384.81	6.34	1.00	2.00			Kidney	Spleen	Breast
10940	416744	W88518	18.62	256.87	13.78	4.00	3.00			Blood	Lymph	Spleen
10950	884283	AA668750	32.83	458.69	13.97	5.00	0.00	11	387.82	Larynx	Head and neck	Esophagus
10951	148469	H12336	18.03	146.94	19.27	5.00	1.00	7	525.98	Breast	Brain	LID not found
10953	487888	AA047478	4.52	48.81	10.36	4.00	0.00	19	238.07	Synovial mem.	Skin	
10954	377275	AA065486	5.04	97.13	19.27	5.00	1.00	6	104.19	CNS	Heart	Lymph
10960	52071	H24350	2.07	10.84	5.14	1.00	0.00	7	653.71	Thymus	Lymph	Lymph node
10963	49318	H15716	20.35	105.91	5.20	0.00	5.00	7	567.71	Brain	LID not found	Other
10968	23579	R38349	61.59	462.04	7.50	0.00	1.00	1	558.36	Muscle	Heart	Breast
10974	293984	N68053	231.06	1187.94	5.05	0.00	1.00	12	44.83	Spleen	Testis	LID not found
10975	308841	N91921	12.57	159.70	12.70	1.00	0.00	16	61.44	Brain	LID not found	Other
10984	23897	R38364	2.06	15.49	7.60	1.00	0.00	19	185.78	CNS	Nose	Forebrain
10985	756931	AA425924	7.75	493.54	63.69	13.00	0.00	19	299.02		Adipose	Pooled
10988	74518	T59016	12.12	67.66	5.58	0.00	2.00	16	21.87	Larynx	Esophagus	Gall bladder
10998	30175	R42433	3.55	28.46	8.02	2.00	0.00	9	416.74	Cervix	Eye	Brain
11008	34041	R44754	2.83	16.26	6.19	2.00	0.00			Brain	LID not found	Other
11011	52741	H29521	12.48	64.80	5.19	1.00	0.00	16	01.44	Ignore	Brain	Bone
11017	280697	N47443	7.45	37.66	5.05	1.00	0.00	19	185.78	CNS	Nose	Forebrain
11018	770588	AA434139	32.28	167.02	5.17	1.00	0.00	19	299.02		Adipose	Pooled
11024	21920	T72533	3.22	25.76	8.00	1.00	3.00	16	21.87	Larynx	Esophagus	Gall bladder
11030	203003	H54417	104.55	599.29	5.73	0.00	1.00	9	416.74	Cervix	Eye	Brain
11031	32229	R42813	177.53	1271.87	7.16	1.00	2.00			Brain	LID not found	Other
11032	21922	T72535	10.28	98.41	9.57	1.00	5.00					



Table 2A

11036	81316	T60061	17.14	282.19	18.46	1.00	1.00	12	118.94	Ear	Foreskin	CNS
11038	741958	AA402891	3.13	28.59	9.14	1.00	0.00	11	242.83	Ovary	CNS	Prostate
11039	34366	R43667	8.40	62.33	9.75	0.00	3.00	9	37.32	Brain	LID not found	Other
11041	41822	R54193	7.40	41.20	5.57	2.00	0.00			Breast	Prostate	Ovary
11043	123085	R00046	38.27	432.24	11.29	2.00	4.00			Pool	LID not found	Other
11046	743113	AA401397	4.81	43.84	9.13	1.00	0.00	6	118.71	Brain	LID not found	Other
11047	838732	AA457543	66.35	464.36	7.00	3.00	0.00	15	75.1	Smooth muscle	Thymus	Placenta
11053	49822	H28985	11.47	110.09	9.59	4.00	5.00	8	151.38	Brain	LID not found	Other
11056	298059	N70520	37.75	200.66	5.32	0.00	1.00	17	53.59	Adrenal gland	Pooled	Lymph
11060	25029	R37615	30.41	311.45	10.24	2.00	5.00	10	373.9	Heart	Pod	LID not found
11061	49284	H15533	4.02	34.87	8.67	0.00	1.00	7	52.28	Adrenal gland	Eye	Ovary
11067	340894	W57872	42.73	363.47	8.51	2.00	0.00	19	214.06	Brain	Pool	LID not found
11074	755581	AA419143	89.80	1386.32	15.41	2.00	0.00	1	338.35	Colon	Germ Cell	Kidney
11075	838287	AA457485	79.16	417.11	5.27	1.00	0.00	8	416.03	Kidney	Aorta	Prostate
11078	827040	AA181019	29.35	184.37	5.62	2.00	1.00	6	308.19	Thyroid	Gall bladder	Parathyroid
11079	46567	H09790	20.69	120.73	5.83	0.00	1.00	20	38.8	Cervix	Thyroid	Stomach
11080	588057	AA134885	6.05	45.92	7.59	1.00	0.00	12	455.53	Breast	Tonsil	Whole embryo
11081	32663	R43318	3.29	34.80	10.57	2.00	0.00	1	711.38	Breast	Lung	Germ Cell
11082	850981	AA689674	574.04	3656.31	6.37	1.00	0.00	4	644.94	Brain	LID not found	Other
11083	627226	AA195463	13.75	157.20	11.43	2.00	0.00	5	-4.6	Brain	Germ Cell	Ovary
11085	50582	H16783	17.87	216.14	12.23	2.00	0.00	5	28.06	Brain	LID not found	Other
11088	626531	AA186327	108.04	628.67	5.82	0.00	1.00	12	45.2	Ignore	Pooled	Pancreas
11092	25081	R38944	25.11	206.06	8.21	4.00	4.00	15	245.32	Neural	Pooled	Lymph
11093	50595	H16785	48.24	322.62	6.89	0.00	0.00	1	740.99	Pod	Brain	LID not found
11097	40384	R54822	11.46	69.53	6.07	2.00	0.00	10	508.5	Pooled	Lymph	Foreskin
11098	884500	AA629987	301.59	2058.98	6.83	1.00	0.00	11	295.61	Lymph	Breast	Placenta
11101	50786	H16832	6.88	38.79	5.64	1.00	0.00	12	287.91	Brain	Whole embryo	Spleen
11103	950497	AA588140	32.84	202.83	6.18	0.00	2.00	12	101.74	Synovial mem	Foreskin	CNS
11105	41186	R56130	3.07	31.21	10.18	2.00	0.00	21	157.88	Brain	LID not found	Other
11108	25396	R17747	6.12	35.72	5.84	1.00	0.00	14	139.45	Ear	Kidney	Brain
11109	50421	H17308	5.13	33.98	6.62	1.00	0.00	12	272.22	Pool	LID not found	Other
11114	302632	N80281	5.04	72.42	14.37	3.00	0.00	11	295.61	Lymph	Breast	Placenta
11115	268785	N24024	103.58	643.57	8.14	2.00	0.00	12	287.91	Brain	Whole embryo	Spleen
11116	25355	R38928	7.52	40.01	5.32	0.00	3.00	12	101.74	Synovial mem	Foreskin	CNS
11117	50805	H17634	97.71	650.07	6.85	3.00	3.00	21	157.88	Brain	LID not found	Other
11118	26186	R20755	13.33	265.82	19.83	7.00	4.00	14	139.45	Ear	Kidney	Brain
11120	510080	AA053411	27.02	163.84	6.06	0.00	1.00	12	272.22	Pool	LID not found	Other
11122	526945	AA113291	12.16	69.90	5.75	1.00	0.00	11	295.61	Lymph	Breast	Placenta
11123	731121	AA417307	14.17	73.75	5.20	0.00	1.00	12	287.91	Brain	Whole embryo	Spleen
11124	25983	R37410	0.88	7.66	7.78	2.00	0.00	12	101.74	Synovial mem	Foreskin	CNS
11126	608743	AA168372	9.77	143.60	14.70	5.00	1.00	21	157.88	Brain	LID not found	Other
11127	595078	AA164819	50.61	820.56	16.21	5.00	0.00	14	139.45	Ear	Kidney	Brain
11128	41214	R58898	0.70	5.31	7.60	3.00	0.00	12	272.22	Pool	LID not found	Other
11133	67440	T49355	33.98	335.05	9.86	1.00	1.00	12	272.22	Pool	LID not found	Other
11139	185274	R92011	34.01	263.80	7.76	3.00	0.00	3	403.91	Ear	Foreskin	Lung
11150	270385	N33063	8.00	1138.23	141.94	22.00	2.00	7	653.71	Pool	LID not found	Other
11154	127230	R08184	12.07	91.10	7.55	3.00	0.00	5	378.73	Pool	Lung	LID not found
11163	272709	N32285	7.35	59.71	8.13	4.00	0.00	17	339.79	Bone	Synovial mem	Eye
11170	417730	W68497	20.23	118.36	6.75	0.00	1.00			Pool	LID not found	Other
11171	198325	R92448	10.58	64.88	6.13	1.00	0.00			CNS	Cdon	Heart
11174	308238	N92415	12.69	83.89	6.82	1.00	0.00			CNS	Pancreas	Ovary
11188	415535	W60635	21.32	198.18	9.34	4.00	3.00			CNS	Pancreas	Ovary
11188	261681	N48057	155.59	1345.02	8.64	4.00	0.00			CNS	Pancreas	Ovary
11181	281045	N50904	5.72	40.59	7.09	1.00	0.00			CNS	Pancreas	Ovary
11200	278875	N63034	121.50	737.87	6.07	1.00	0.00			CNS	Pancreas	Ovary

Table 2A

11208	428828	AA005254	2.27	32.71	14.39	2.00	0.00	15	185.7	Brain	Uterus	Tonsil
11209	210554	H85934	3.17	24.37	7.70	1.00	0.00	14	123.62	CNS	LID not found	Other
11210	278759	N62946	4.15	53.38	12.85	2.00	0.00	2	508.52	Tonsil	Lung	Pool
11211	307740	N92947	8.18	887.19	14.50	6.00	5.00	6	567.03	Pool	Muscle	Pancreas
11212	416400	V92798	19.70	100.34	5.03	1.00	0.00	8	548.99	CNS	Pancreas	Tonsil
11220	502674	AA135886	8.31	84.03	10.11	1.00	0.00	5	343.84	Brain	Tonsil	Forebrain
11222	131599	R23727	2.43	20.49	8.43	3.00	0.00	1	832.4	Adrenal gland	Pooled	Parathyroid
11237	782575	AA447522	3.08	145.25	47.20	4.00	2.00	15	238.08	Lung	Blood	Parathyroid
11242	204740	H57308	23.86	220.28	9.23	3.00	0.00	17	305.8	Adipose	Muscle	Prostate
11246	811033	AA485432	10.41	68.40	6.57	0.00	1.00	3	180.78	CNS	Lymph	Testis
11254	795165	AA453474	18.18	290.55	15.98	5.00	0.00	10	522.77	Forebrain	Ovary	Blood
11258	753751	AA496830	4.20	70.39	16.76	1.00	0.00	12	252.49	Liver	Pool	Kidney
11258	811604	AA454610	264.67	1560.18	5.89	1.00	0.00	3	472.48	Testis	Heart	Heart
11262	284160	N53512	4.39	24.23	5.52	1.00	0.00	12	20.42	Germ Cell	Pancreas	LID not found
11265	785207	AA453588	3.83	31.10	8.01	2.00	0.00	7	675.72	Spleen	Heart	LID not found
11267	143450	R74478	9.82	89.02	9.07	4.00	0.00	3	688.28	Pool	LID not found	Other
11269	490060	AA136080	8.87	47.17	6.87	1.00	0.00	X	251.52	Pool	LID not found	Other
11283	488584	AA044814	11.82	88.78	7.51	6.00	0.00	5	345.44	Pool	Small intestine	Skin
11287	418154	V85927	74.78	538.79	7.20	3.00	0.00	8	31.44	Thyroid	Synovial mem	Blood
11288	382564	AA069372	41.81	524.80	12.55	0.00	2.00	8	90.95	Pooled	Spleen	Uterus
11289	385707	AA025434	4.72	115.70	24.53	1.00	0.00	17	308.88	Adipose	Pancreas	Tonsil
11290	809829	AA455528	77.17	908.87	11.79	3.00	0.00	3	683.89	Colon	Pool	LID not found
11295	292738	N69393	11.57	85.13	7.36	0.00	2.00	2	568	CNS	Pooled	Eye
11298	346888	V74602	17.31	111.75	6.45	0.00	3.00	21	18.92	Brain	LID not found	Other
11297	795276	AA454022	6.51	42.66	8.55	1.00	0.00	6	201.24	Pooled	Forebrain	Brain
11298	795171	AA453470	3.97	32.74	8.24	2.00	0.00	1	591.57	Brain	Brain	LID not found
11301	415630	V84780	20.12	168.22	8.38	1.00	0.00	18	81.46	Nose	Stomach	Blood
11303	244012	N38787	37.49	588.81	15.70	4.00	0.00	7	546.17	Brain	Pancreas	LID not found
11308	201172	R98487	15.17	88.43	5.83	1.00	0.00	1	785.21	Brain	Brain	LID not found
11317	324205	V47179	176.99	1085.53	6.02	2.00	0.00	5	49.48	Muscle	Brain	LID not found
11320	471859	AA035144	1.47	17.20	11.86	1.00	0.00	X	309.33	Lymph	Gall bladder	Kidney
11321	490329	AA127741	3.52	45.33	12.88	7.00	0.00	4	613.42	Parathyroid	Aorta	Heart
11324	147075	R60235	18.02	106.06	5.89	1.00	0.00	13	155.48	Brain	LID not found	Other
11327	201090	R98849	30.93	170.88	5.52	1.00	0.00	X	272.16	Larynx	Spleen	Skin
11330	362409	AA018457	4.24	45.89	10.82	0.00	0.00	18	113.09	Synovial mem	Skin	CNS
11333	416390	V86860	19.38	103.64	5.35	0.00	1.00	18	113.09	Synovial mem	Skin	CNS
11344	40598	R55367	2.18	20.79	9.54	2.00	0.00	18	113.09	Synovial mem	Skin	CNS
11347	52865	H26820	69.85	912.96	13.07	4.00	0.00	18	113.09	Synovial mem	Skin	CNS
11348	52847	H26771	56.91	407.74	6.92	5.00	2.00	18	113.09	Synovial mem	Skin	CNS
11349	32299	R42685	14.61	95.97	6.84	0.00	1.00	18	113.09	Synovial mem	Skin	CNS
11350	377461	AA055835	131.51	1261.57	9.59	0.00	0.00	18	113.09	Synovial mem	Skin	CNS
11351	50559	H18780	99.80	753.04	7.56	1.00	3.00	18	113.09	Synovial mem	Skin	CNS
11352	47074	H10417	10.85	55.53	5.07	0.00	1.00	18	113.09	Synovial mem	Skin	CNS
11360	51775	H23213	2.21	18.88	7.55	1.00	0.00	18	113.09	Synovial mem	Skin	CNS
11361	141852	R70505	5.30	28.20	5.32	0.00	0.00	18	113.09	Synovial mem	Skin	CNS
11362	286663	N62394	3.12	47.72	15.31	11.00	0.00	18	113.09	Synovial mem	Skin	CNS
11363	50114	H18743	9.93	108.85	10.74	4.00	0.00	18	113.09	Synovial mem	Skin	CNS
11364	25636	R36098	47.73	328.03	6.87	1.00	3.00	18	113.09	Synovial mem	Skin	CNS
11368	845355	AA644088	16.98	115.78	6.83	2.00	0.00	18	113.09	Synovial mem	Skin	CNS
11373	344432	V73474	36.81	247.87	6.38	1.00	0.00	18	113.09	Synovial mem	Skin	CNS
11376	33022	R44607	63.62	746.53	11.70	6.00	3.00	18	113.09	Synovial mem	Skin	CNS
11377	158437	R73570	8.80	228.88	28.12	1.00	0.00	18	113.09	Synovial mem	Skin	CNS
11389	811064	AA485442	37.64	479.89	12.75	0.00	1.00	18	113.09	Synovial mem	Skin	CNS
11390	854138	AA69272	1.38	14.83	10.85	2.00	0.00	18	113.09	Synovial mem	Skin	CNS
11391	50843	H17551	39.45	319.46	8.10	0.00	1.00	18	113.09	Synovial mem	Skin	CNS

Table 2A

11392	33653	R44707	70.43	923.03	13.11	2.00	4.00	5	358.9	Placenta	Brain	LID not found
11393	432072	AA679278	9.96	71.22	7.15	1.00	0.00	18	477.48	Cervix	Pooled	Tonsil
11395	47054	H11012	2.84	18.65	6.56	2.00	0.00	15	59.38	Brain	Parathyroid	Testis
11398	50030	H17051	3.00	26.16	8.73	2.00	0.00	1	629.81	Brain	LID not found	Other
11403	39800	R52641	48.55	578.34	11.91	4.00	3.00	1	364.02	Brain	LID not found	Other
11408	46385	H09774	24.85	200.07	8.02	2.00	5.00	12		Brain	CNS	Germ Cell
11418	558447	AA630800	31.86	544.63	17.04	10.00	0.00			Lymph node	Neural	Thymus
11432	46375	H09086	44.82	394.78	8.59	2.00	5.00			Brain	Pod	LID not found
11433	60738	T40840	2.97	15.05	5.07	1.00	0.00	1	716.02	Blood	Breast	Ovary
11434	592491	AA160498	142.34	1029.78	7.23	0.00	4.00	15	171.53	Adipose	Lymph	Pancreas
11437	24938	R38865	7.38	58.03	7.86	1.00	1.00	1	628.85	Ear	Brain	Testis
11441	51083	H17139	1.49	13.89	9.29	1.00	0.00	5	39.72	Brain	CNS	Eye
11444	49941	H29211	38.28	357.57	9.34	4.00	5.00			Brain	Pod	LID not found
11446	626018	AA187938	29.87	544.51	18.23	1.00	0.00			Brain	Pod	LID not found
11447	842896	AA486418	40.08	337.37	8.42	0.00	1.00	7	94.82	Umbilical cord	Parathyroid	Cervix
11448	838999	AA484935	6.92	45.33	6.55	0.00	1.00	4	635.64	Stomach	Colon	Brain
11452	773573	AA428239	6.72	124.52	18.54	4.00	0.00	4	109.02	Spleen	Whole embryo	Pool
11460	52057	H24323	2.94	19.63	6.69	1.00	0.00	8	300.54	Brain	Aorta	Pool
11473	841624	AA487462	17.64	91.09	5.16	1.00	0.00			Foreskin	Placenta	LID not found
11477	70384	T54474	28.38	188.40	6.64	3.00	0.00	2	56.46	Pooled	Brain	Codon
11480	838668	AA457235	8.00	121.81	15.23	2.00	0.00	1	649.11	Blood	Testis	Germ Cell
11482	731044	AA421273	18.11	121.44	6.71	0.00	1.00	12	247.55	Cervix	Lymph node	Marrow
11487	842973	AA488332	209.06	1432.51	6.85	1.00	1.00	21	152.58			
11489	70027	T48767	13.55	76.05	5.61	2.00	0.00	9	274.79	Brain	CNS	Whole embryo
11492	53081	R15764	1.72	28.75	16.69	3.00	0.00	17	317.6	Pooled	Germ Cell	Whole embryo
11496	758318	AA404273	2.00	17.37	8.67	1.00	0.00	7	408.94	Brain	LID not found	Other
11498	32530	R43286	4.34	23.38	5.38	1.00	1.00			Pooled	Lymph	Eye
11502	587992	AA130596	31.72	257.83	8.12	0.00	0.00	3	120.51	Lymph node	Head and nec	Foreskin
11509	840763	AA480092	17.02	261.44	15.36	9.00	0.00			Ovary	Testis	Whole embryo
11510	950667	AA608560	8.23	81.41	9.89	1.00	0.00	9	377.82	Cervix	Brain	Ovary
11516	28203	R20850	3.10	16.38	5.28	1.00	1.00	4	671.55	CNS	Testis	Ovary
11527	268023	N62737	3.88	32.62	8.19	1.00	0.00	5	421.53	Tonsil	Lung	Pool
11534	302025	N69753	6.30	66.71	10.59	2.00	0.00	1	825.71	Blood	Spleen	Pancreas
11539	810737	AA457723	17.34	117.63	6.78	1.00	0.00	14	280.52		Pool	Tonsil
11542	124742	R02173	7.18	51.57	7.18	0.00	2.00			Foreskin	LID not found	Other
11554	479353	AA007522	34.31	183.25	5.34	1.00	0.00	12	104.58	Foreskin	Testis	LID not found
11580	290748	N71792	4.26	22.73	5.33	1.00	0.00	6	86.9	Pod	LID not found	Other
11586	415904	W68195	5.77	42.00	7.28	3.00	0.00	X	272.18	Larynx	Spleen	Skin
11588	268286	N26863	15.68	102.37	6.63	1.00	0.00	21	245.37	Nose	Kidney	Colon
11590	415894	W65367	3.87	41.48	10.71	0.00	1.00	11	230.82	Peripheral	Neural	Thyroid
11604	291091	N67878	18.53	102.08	5.51	0.00	0.00	10	249.3	Tonsil	Placenta	Pod
11617	810050	AA455286	32.08	192.03	5.99	3.00	0.00	1	553.07	Bone marrow	Larynx	Smooth muscle
11619	298087	N69044	524.79	4910.39	9.36	2.00	0.00			Eye	Brain	Heart
11620	415321	W92041	3.44	24.85	7.23	2.00	0.00	1	554.63	Skin	CNS	Testis
11622	250033	N30316	6.26	34.38	5.49	0.00	1.00			Foreskin	Colon	
11623	810621	AA464744	6.05	32.34	5.35	2.00	0.00			Neural	Umbilical cord	Tonsil
11624	562729	AA086471	9.30	466.41	50.13	3.00	0.00			Eye	Brain	Heart
11630	293438	N68864	541.51	4939.44	8.12	1.00	0.00			Neural	Umbilical cord	Tonsil
11631	810082	AA484952	50.59	308.50	6.10	2.00	0.00			Eye	Brain	Heart
11632	377441	AA055242	6.52	126.13	19.33	2.00	0.00			Foreskin	Colon	Testis
11636	281274	H98255	6.87	46.50	6.77	0.00	1.00			Neural	Umbilical cord	Tonsil
11639	245137	N54395	2.84	30.60	10.48	4.00	0.00			Neural	Umbilical cord	Tonsil
11640	758595	AA444051	21.04	251.75	11.97	3.00	0.00			Neural	Umbilical cord	Tonsil
11643	302221	N77828	20.48	138.65	6.87	2.00	0.00			Neural	Umbilical cord	Tonsil
11645	357544	W84063	11.51	100.88	8.76	0.00	2.00			Neural	Umbilical cord	Tonsil





Table 2A

12046	12447	R01094	17.48	123.68	7.07	6.00	0.00	2	680.88	Ovary	Spleen	Codon
12049	481311	AA150188	2.52	15.80	6.31	2.00	0.00			Uterus	LID not found	Other
12051	324323	W47552	10.21	69.89	8.64	3.00	0.00	6	483.25			
12053	501876	AA128008	67.89	557.11	8.19	3.00	0.00	18	27.41	Uterus	LID not found	Other
12054	428165	AA005108	10.29	55.37	5.38	1.00	0.00	11	292.28	Neural	Pool	LID not found
12058	854896	AA830084	38.84	391.68	10.84	6.00	0.00	8	564.08	Larynx	Esophagus	Thymus
12059	204536	H58250	99.09	613.80	6.26	1.00	0.00	5	356.9			
12060	291890	N67487	67.91	404.61	5.88	1.00	0.00	1	62.05	Skin	Bone	Blood
12061	502618	AA136052	2.88	33.04	11.48	1.00	0.00	21	142.57	Uterus	Prostate	LID not found
12063	502333	AA158783	24.09	155.77	6.47	2.00	0.00	20	287.89	Head and nec	Adipose	Adipose
12066	795687	AA459937	4.98	36.48	7.33	3.00	0.00			Testis	Brain	LID not found
12069	502634	AA127017	35.10	445.08	12.68	7.00	0.00	1	15.89	Muscle	Uterus	Brain
12071	810711	AA457700	34.08	205.23	6.02	0.00	1.00	10	480.31	Breast	Small intestine	Ovary
12073	501890	AA128017	2.81	23.99	8.34	2.00	1.00			Uterus	Testis	LID not found
12077	770788	AA427621	5.18	37.49	7.24	2.00	0.00	12	298.79	Bone	Kidney	Ovary
12085	301867	N92478	12.26	68.39	5.58	1.00	0.00	15	147.26	Gall bladder	Bone	Ovary
12086	771060	AA427522	5.60	50.48	9.01	1.00	0.00			Ovary	Brain	LID not found
12087	202194	H52379	3.17	23.48	7.41	1.00	0.00	1	568.94	Pool	Bone	LID not found
12088	858135	AA830804	0.64	6.24	9.68	2.00	0.00	18	200.17	Small intestine	Skin	CNS
12089	481244	AA152299	2.47	15.89	6.47	2.00	0.00			Uterus	LID not found	Other
12092	590284	AA155813	130.40	677.23	5.19	1.00	0.00	12	68.19	Umbilical cord	Adipose	Aorta
12093	503749	AA131469	4.93	45.82	9.30	1.00	0.00	6	328.81	Uterus	LID not found	Other
12094	144825	R78521	536.48	8898.87	18.78	5.00	0.00	19	271.02			
12103	81044	T40725	41.18	403.08	9.79	5.00	1.00			Spleen	Lymph	Pool
12104	48931	H10030	3.69	19.99	5.41	0.00	1.00	14	247.37	Brain	Testis	Prostate
12108	51511	H18958	3.08	17.82	5.77	2.00	0.00	14	175.36	Heart	Brain	Breast
12113	45378	H07920	3.50	224.15	83.98	2.00	1.00			Pool	Muscle	Lymph
12119	81462	T40927	2.31	14.80	6.41	1.00	0.00			Liver	LID not found	Other
12120	50983	H18017	1.08	6.60	6.10	1.00	0.00	1	173.21	Brain	LID not found	Other
12122	433253	AA699427	10.02	73.68	7.36	1.00	0.00					
12123	46716	H10012	4.26	28.77	6.98	1.00	0.00	17	534.21	Blood	Whole embryo	Tonsil
12129	141615	R70885	25.04	189.33	7.66	2.00	0.00	20	73.18	Head and nec	Blood	Pool
12130	279790	N48355	45.38	280.10	6.17	1.00	0.00	X	337.33	CNS	Tonsil	Uterus
12131	85409	T71891	92.22	508.27	5.49	1.00	0.00	1	606.38	Larynx	Thymus	Stomach
12134	724888	AA291484	3.75	119.61	31.92	10.00	0.00	1	151.1	Ovary	Eye	Heart
12136	34321	R44949	1.63	10.11	6.20	1.00	0.00	10	90.1	Brain	LID not found	Other
12138	502367	AA134871	46.83	437.30	9.34	1.00	0.00	22	154.77	Ignore	Placenta	Esophagus
12139	51700	H22834	7.76	118.71	15.29	1.00	0.00			Parathyroid	CNS	Brain
12140	25984	R37411	1.82	73.60	40.50	4.00	3.00	9	237.9	Small intestine	Brain	LID not found
12143	77811	T61289	11.40	188.77	16.56	1.00	0.00	8	447.53	Liver	Foreskin	Pool
12145	855745	AA663981	29.67	3034.57	102.27	9.00	0.00	14	278.45			
12148	45391	H08194	4.44	47.46	10.68	1.00	0.00	7	30.48	Brain	Heart	Pool
12150	854401	AA668959	15.45	87.81	5.68	1.00	1.00	10	313.32	Pool	Lung	LID not found
12151	78946	T61792	49.64	590.37	11.89	1.00	0.00	7	492.68	Small intestine	Thymus	Synovial membrane
12152	34442	R44986	17.53	125.84	7.17	1.00	0.00	3	143.12	CNS	Whole embryo	Brain
12153	77539	T58775	41.58	351.19	8.45	3.00	0.00			Liver	LID not found	Other
12154	855755	AA663986	227.17	1789.14	7.88	1.00	1.00	18	235.13	Umbilical cord	Skin	Stomach
12160	46328	H10641	2.32	26.81	11.57	2.00	0.00	15	206.38	Brain	LID not found	Other
12162	868380	AA634109	13.50	98.24	7.27	3.00	0.00	1	578.51	Placenta	Aorta	Gall bladder
12163	83345	T68440	110.82	1275.71	11.51	3.00	3.00	9	140.54	Pool	Liver	Foreskin
12164	32331	R42922	81.54	770.47	9.45	3.00	5.00	1	258.83	Brain	LID not found	Other
12168	307471	N93505	28.77	160.20	5.38	1.00	0.00	X	115.87	Adrenal gland	Ear	Eye
12169	745019	AA628028	8.25	43.41	6.95	2.00	0.00	11	248.09	Tonsil	Gall bladder	Pancreas
12170	307660	N92801	21.63	766.95	35.45	1.00	0.00			Adipose	Pool	Heart
12171	47428	H11088	12.11	90.48	7.47	0.00	1.00	12	349.48	Stomach	Spleen	Placenta

Table 2A

12172	40100	R54590	3.74	35.78	9.56	2.00	0.00	1	674.22	CNS	Colon	Lung
12173	345626	W72051	3.24	19.91	6.15	1.00	0.00			Bone	CNS	Brain
12174	530185	AA111069	12.83	93.26	7.27	1.00	0.00	6	43.37	Thyroid	Lymph	Tonsil
12176	49443	H15427	1.53	9.47	6.17	1.00	1.00	1	597.06	Brain	LID not found	Other
12178	79565	T82854	11.13	56.82	5.11	2.00	0.00	12	51.04	Colon	Bone	Lung
12180	50443	H17325	1.49	8.81	5.90	1.00	0.00	17	55.98	Brain	LID not found	Other
12181	739155	AA421819	25.69	382.06	14.76	0.00	1.00	5	116.08	Pooled	Ovary	Whole embryo
12185	740604	AA479785	16.77	148.77	8.87	2.00	0.00			Esophagus	Gall bladder	Blood
12190	431655	AA676453	10.87	55.57	5.21	1.00	0.00	19	258.89	Lymph	Blood	Breast
12192	32587	R43595	36.21	546.94	15.10	2.00	4.00	18	62.76	Germ Cell	Foreskin	Heart
12193	83358	T68445	16.69	134.75	8.12	4.00	0.00	X	238.33	Smooth musc	Prostate	
12201	84713	T74257	8.01	47.51	5.93	2.00	0.00	4	534.21	Gall bladder	Liver	Spleen
12202	25132	R37620	28.43	155.76	5.29	2.00	0.00			Brain	LID not found	Other
12210	838829	AA481769	43.33	220.08	5.08	0.00	1.00	1	667.01	Thyroid	Testis	Heart
12213	46977	H10372	166.82	2024.13	12.15	0.00	4.00	19	71.09			
12217	67759	T49652	11.14	75.15	6.74	2.00	0.00	13	88.36	Gall bladder	Blood	Lymph
12218	773556	AA428182	185.52	1230.88	8.20	0.00	2.00					
12219	843008	AA488391	117.75	1108.01	9.41	0.00	1.00	3	137.52	Thymus	Lymph	Cervix
12227	842946	AA489324	79.39	482.56	6.08	1.00	2.00	17	319.26	Testis	Omentum	Larynx
12235	610464	AA056656	63.42	451.49	7.12	2.00	0.00	17	341.88		Pooled	Brain
12237	49303	H15877	3.04	118.26	38.92	3.00	0.00	5	552.2	Ear	Placenta	Aorta
12240	786760	AA460722	60.77	633.37	10.42	1.00	2.00			Stomach		
12244	23774	R38196	6.61	55.60	8.41	5.00	0.00	1	630.22		LID not found	Other
12245	46553	H10679	7.60	50.71	8.67	4.00	2.00	9	382.31	Brain	Pool	Parathyroid
12247	234977	H73640	9.43	60.55	6.42	0.00	2.00	8	89.33	Foreskin	Foreskin	Tonsil
12250	786285	AA460848	26.70	137.22	5.14	1.00	0.00	1	737.93	Thyroid	Uterus	Colon
12254	602277	AA156597	18.64	96.67	5.19	1.00	0.00	16	63.62	Testis	Foreskin	Thyroid
12259	897761	AA598468	120.52	741.07	6.15	0.00	1.00			Larynx	LID not found	Other
12260	22762	R38543	3.49	33.77	9.87	2.00	0.00	8	483.92	Brain		
12261	52755	H28783	9.47	219.72	23.20	9.00	0.00					
12262	612809	AA179600	5.25	30.81	5.87	0.00	1.00	10	349.63	Cervix	LID not found	
12267	838774	AA457566	2.16	27.87	12.89	2.00	0.00	1	720.56	Thymus	Cervix	Foreskin
12268	22773	R38613	187.84	2384.83	12.75	2.00	0.00	X	245.06	Whole embryo	Heart	Brain
12271	731270	AA416664	11.64	61.68	5.30	0.00	1.00			Parathyroid	Heart	Whole embryo
12272	136399	R34287	72.17	498.81	6.91	0.00	1.00			Thymus	Placenta	Cervix
12276	231116	R39179	95.08	674.88	7.10	4.00	3.00	16	502.56			
12281	41188	R56870	31.76	160.84	5.06	1.00	0.00	5	810.47	Thyroid	Ear	Lung
12282	48843	H10072	5.18	27.22	5.28	1.00	0.00	19	20.17	Heart	Brain	Breast
12283	731023	AA421266	8.09	40.62	5.02	0.00	1.00			Stomach	Tonsil	Kidney
12285	50266	H17463	4.81	67.50	13.74	1.00	0.00			Kidney	Brain	Brain
12286	286503	N87386	86.20	563.46	6.54	0.00	1.00			Ear	Lymph	Umbilical cord
12280	320209	WD4509	3.54	27.56	7.79	0.00	1.00			Germ Cell	Parathyroid	Testis
12300	298091	N70756	65.05	501.31	7.71	2.00	2.00			Lung	LID not found	Other
12304	309368	N83967	171.56	1303.92	7.60	2.00	3.00	4	450.37	Lung	LID not found	Other
12306	784214	AA446865	40.13	280.86	7.00	1.00	0.00	15	307.86	Adipose	Umbilical cord	Lymph
12307	251147	H97366	18.98	104.47	5.51	0.00	1.00	4	843.02	Foreskin	Parathyroid	Whole embryo
12316	525478	AA05042	9.73	110.92	11.40	4.00	0.00	21	170.13	Tonsil	Prostate	Parathyroid
12319	260170	N32072	24.77	136.35	12.77	2.00	0.00			Trachea	Parathyroid	Thyroid
12322	510273	AA053165	167.30	2378.38	14.22	5.00	0.00	11	249.52		LID not found	Other
12328	309895	N94488	8.91	114.36	12.84	4.00	2.00	1	81.33	Pool	LID not found	Other
12343	259275	N32847	600.52	6799.12	11.32	2.00	0.00	14	75.08	Parathyroid	Whole embryo	Testis
12344	321805	W37663	8.39	161.74	19.27	1.00	1.00			Tonsil	Eye	Aorta
12351	594323	AA169202	15.71	83.20	5.29	1.00	0.00			Uterus	Tonsil	Colon
12358	127368	R08548	26.48	133.72	5.05	0.00	1.00			Colon	Pool	LID not found
12382	510380	AA055404	78.28	666.72	8.41	0.00	5.00					

Table 2A

12363	281408	H98967	678.82	4011.17	5.91	1.00	0.00	1.00	0.00	3	139.86	Blood	Foreskin	Heart
12364	291082	N72113	593.39	5.78	0.00	0.00	1.00	0.00	1.00	17	385.71	Muscle	Tonsil	Foreskin
12370	299723	N75055	32.59	276.70	8.49	0.00	1.00	0.00	1.00	11	276.96	Cervix	Aorta	Breast
12374	950603	AA608346	34.98	181.50	5.19	0.00	1.00	0.00	1.00			Thyroid	Ear	Adrenal gland
12382	773283	AA425214	3.52	76.41	21.70	1.00	0.00	0.00	0.00			Pooled	Muscle	Breast
12384	308446	N95490	93.80	671.66	7.16	0.00	2.00	0.00	2.00			Aorta	Lung	Pool
12387	298236	N70837	3.44	21.88	6.37	2.00	0.00	0.00	0.00			Ear	Ovary	Eye
12395	286568	N87305	37.55	353.93	9.43	0.00	1.00	0.00	1.00	7	84.81	Colon	LID not found	Other
12397	511302	AA086005	3.76	32.83	5.70	1.00	0.00	0.00	0.00			Thyroid	Whole embryo	Lung
12402	300815	N80764	3.14	51.61	16.44	0.00	1.00	0.00	1.00			Ear	Prostate	Whole embryo
12404	950594	AA608531	8.90	91.05	10.23	1.00	0.00	0.00	0.00	1	689.84	Testis	LID not found	Other
12406	743314	AA400482	3.81	21.93	5.75	1.00	0.00	0.00	0.00	5	509.69	Heart	LID not found	Other
12411	304868	N93197	13.85	92.16	6.65	1.00	1.00	0.00	1.00			Parathyroid	LID not found	Other
12419	305408	N95073	59.74	502.87	8.42	0.00	3.00	0.00	3.00	10	381.57	Parathyroid	LID not found	Other
12420	758343	AA404266	126.17	650.44	5.16	0.00	1.00	0.00	1.00	10	45.1	Kidney	Tonsil	Colon
12421	511776	AA088438	55.43	1150.75	20.76	0.00	0.00	0.00	0.00			Head and neck	Prostate	Prostate
12431	784016	AA443868	42.95	296.62	6.91	1.00	1.00	0.00	1.00			Cervix	Tonsil	Prostate
12441	826716	AA191548	34.14	183.71	5.38	0.00	1.00	0.00	1.00			Ear	Adrenal gland	Parathyroid
12444	743297	AA400422	22.29	165.43	7.42	0.00	1.00	0.00	1.00			Ear	Adrenal gland	Parathyroid
12449	840530	AA487834	38.21	298.35	8.24	2.00	0.00	0.00	0.00			Ear	Adrenal gland	Parathyroid
12461	488633	AA101878	20.37	319.08	15.74	0.00	1.00	0.00	1.00			Ear	Adrenal gland	Parathyroid
12463	757327	AA437094	24.84	140.70	5.66	0.00	1.00	0.00	1.00			Ear	Adrenal gland	Parathyroid
12465	375086	AA040369	3.76	36.98	8.83	0.00	2.00	0.00	2.00	4	490.66	Heart	Pool	LID not found
12466	324762	AA284112	7.77	71.23	9.16	1.00	4.00	0.00	4.00			Pool	LID not found	Other
12468	242009	H93318	7.91	88.94	11.24	1.00	0.00	0.00	0.00			Testis	LID not found	Other
12470	742830	AA406048	7.47	38.53	5.16	1.00	0.00	0.00	0.00	14	280.52	Whole embryo	Parathyroid	Kidney
12471	796166	AA461084	120.27	788.29	6.55	0.00	2.00	0.00	2.00			Pooled	Kidney	LID not found
12475	323251	W42996	0.84	4.42	5.25	1.00	0.00	0.00	0.00	2	467.65	Pool	Tonsil	-
12476	113257	T83846	12.99	188.46	14.51	0.00	2.00	0.00	2.00			Muscle	LID not found	Other
12477	582447	AA100595	94.86	843.60	8.88	1.00	2.00	0.00	2.00	11	301.58	Brain	LID not found	Other
12482	292982	N69100	33.00	217.30	6.59	4.00	0.00	0.00	0.00			Testis	LID not found	Other
12483	1031599	AA609485	8.51	77.75	9.14	0.00	3.00	0.00	3.00			Ovary	Placenta	Brain
12485	134297	R31933	110.55	750.61	6.79	0.00	2.00	0.00	2.00	1	668.75	Stomach	Brain	Pool
12488	48142	H12105	6.34	50.11	7.91	0.00	1.00	0.00	1.00			Testis	Kidney	Pool
12489	1031807	AA609744	6.81	74.79	10.98	1.00	0.00	0.00	0.00	9	377.24	Ovary	LID not found	Other
12493	77371	T55340	84.68	464.12	5.48	0.00	1.00	0.00	1.00	6	78.75	Smooth muscle	Thymus	Thymus
12495	788205	AA453420	18.45	107.33	5.82	1.00	0.00	0.00	0.00	5	380.7	Brain	Pool	Thymus
12468	43940	H04785	6.91	36.33	5.26	0.00	1.00	0.00	1.00			Pancreas	Kidney	Heart
12497	1031919	AA609749	54.00	370.29	6.86	0.00	1.00	0.00	1.00	6	463.04	Brain	LID not found	Other
12499	1031719	AA609585	60.21	508.50	8.45	0.00	2.00	0.00	2.00			Testis	LID not found	Other
12504	41999	R39960	51.47	398.04	7.73	2.00	0.00	0.00	0.00	1	722.7	Foreskin	Bone	Eye
12508	298468	N70203	116.59	864.90	7.58	0.00	2.00	0.00	2.00	X	143.33	Testis	LID not found	Other
12507	1031748	AA609589	4.41	48.96	11.10	1.00	1.00	0.00	1.00	20	194.39	Tonsil	Heart	Brain
12511	768260	AA424950	14.35	88.83	6.18	1.00	0.00	0.00	0.00	5	320.68	Brain	LID not found	Other
12512	35366	R45567	33.88	313.65	9.23	2.00	5.00	0.00	5.00	1	82.43	Spleen	LID not found	Other
12517	74283	T55236	1156.15	5802.59	5.02	1.00	0.00	0.00	0.00			Pancreas	LID not found	Other
12518	593026	AA159358	6.05	60.40	9.99	1.00	0.00	0.00	0.00			Spleen	LID not found	Other
12525	73472	T55437	59.85	560.89	9.37	0.00	3.00	0.00	3.00	14	151.82	Ovary	LID not found	Other
12526	593972	AA165116	85.83	627.56	7.31	0.00	2.00	0.00	2.00	9	377.31	Pool	LID not found	Other
12527	811927	AA454668	40.64	427.64	10.52	5.00	0.00	0.00	0.00	12	200.16	Pancreas	Blood	Whole embryo
12528	35826	R45292	23.27	312.81	13.44	0.00	2.00	0.00	2.00	3	484.87	Germ Cell	Heart	Kidney
12531	1031767	AA609626	121.20	871.57	7.19	0.00	2.00	0.00	2.00			Testis	LID not found	Other
12535	36491	R48700	1.89	10.11	5.36	1.00	0.00	0.00	0.00	15	171.5	Breast	Pancreas	Uterus
12536	294665	N71303	228.24	2498.37	10.95	0.00	2.00	0.00	2.00			Pool	LID not found	Other
12548	665405	AA185002	19.86	111.28	5.60	1.00	0.00	0.00	0.00			Pool	Pancreas	Uterus



Table 2A

12549	38816	R49144	978.67	9313.36	9.52	0.00	1.00	0.00	2	682.76	Colon	Pool
12550	594178	AA169498	74.70	595.57	7.97	0.00	4.00	LID not found	5	413.49	Ovary	Other
12555	1031785	AA609848	30.95	469.68	15.17	0.00	4.00	LID not found	20	193	Testis	Other
12557	27315	R40208	730.67	7138.07	9.77	0.00	2.00	Adipose	5	511.76	Nose	Thyroid
12558	594226	AA189535	305.00	1526.80	5.01	0.00	1.00	Brain			LID not found	
12563	1031839	AA609895	42.40	500.01	11.79	0.00	4.00	Testis			LID not found	
12567	685358	AA194983	16.91	122.78	7.26	1.00	1.00	Aorta	22	27.42	Thyroid	Brain
12568	38676	R51504	11.96	82.28	6.88	1.00	0.00	Placenta	12	311.24	Brain	Brain
12573	32775	R43521	5.65	30.93	5.47	0.00	1.00	LID not found	4	422.79	Adipose	Other
12581	787690	AA417950	10.68	95.74	8.86	1.00	0.00	Blood	2	599.98	Forebrain	CNS
12589	767706	AA417956	55.25	1348.92	24.41	0.00	0.00	Tonsil	2	317.73	Adipose	Parathyroid
12592	280507	N47312	37.55	195.86	5.22	1.00	0.00	Adrenal gland	X	339.35	Forebrain	Pooled
12597	787721	AA417882	10.09	81.46	6.09	1.00	0.00	Pool	11	419.03	Colon	LID not found
12603	38883	R51305	2.72	30.89	11.37	0.00	4.00	Testis	3	358.39	Testis	Brain
12607	43966	H04828	24.05	401.19	16.68	4.00	5.00	Brain	18	427.39	Forebrain	Pool
12608	271744	N31585	5.54	56.58	10.21	0.00	1.00	Colon	4	295.78	Ignore	Head and neck
12620	1492104	AA088146	521.50	2629.07	5.04	0.00	1.00	Cervix	X	255.55	Heart	Lung
12623	43764	H05089	29.07	321.68	11.06	0.00	4.00	Brain	19	214.37	Whole embryo	Prostate
12630	785842	AA450338	14.39	104.77	7.28	1.00	0.00	Uterus	14		Stomach	Lung
12631	43879	H05939	3.18	25.00	7.86	3.00	0.00	Testis			Whole embryo	Cell
12638	785585	AA449444	4.74	110.99	23.41	7.00	0.00	Stomach			Umbilical cord	Aorta
12641	754449	AA410296	16.42	91.22	5.56	0.00	2.00	Prostate	2	412.17	Prostate	Thyroid
12659	39147	R51836	148.86	1173.53	7.88	0.00	2.00	Umbilical cord	1	293.34	Umbilical cord	Muscle
12661	767823	AA418728	21.70	255.92	11.79	0.00	0.00	Testis	9	358.85	Brain	Lung
12663	44154	H08154	10.87	59.93	5.46	1.00	0.00	Whole embryo	9	256.34	Forebrain	Kidney
12667	39265	R51871	95.00	662.13	6.97	2.00	0.00	Pancreas	18	78.13	Adipose	Placenta
12668	1492426	AA0878576	17.62	166.98	9.48	0.00	0.00	LID not found			Other	
12669	767843	AA418743	21.49	180.87	8.89	1.00	0.00	Pool			Parathyroid	Aorta
12680	415042	W93106	5.60	29.15	5.21	1.00	0.00	Pool			CNS	LID not found
12682	742679	AA401376	3.93	47.69	12.13	3.00	0.00	Parathyroid			CNS	LID not found
12683	280281	N47861	0.46	3.19	8.86	2.00	0.00	Heart	6	49.89	Heart	LID not found
12692	343381	W67193	36.73	269.53	7.34	0.00	2.00	Uterus	11	343.85	Aorta	Breast
12694	504959	AA149051	17.89	264.78	14.80	2.00	0.00	Thyroid	5	527.16	Neural	Brain
12706	610883	AA172188	13.40	97.79	7.30	2.00	0.00	CNS	7	443.2	CNS	LID not found
12715	281970	N48181	106.13	820.25	5.84	2.00	0.00	Heart			Heart	Tonsil
12732	342497	W68268	6.59	75.47	11.46	1.00	0.00	Esophagus	5	627.78	Nose	Thymus
12734	593104	AA158600	73.48	401.27	5.46	1.00	0.00	Liver	1	131.82	Liver	Pancreas
12742	951091	AA620458	22.43	115.93	5.17	0.00	1.00	Heart			Heart	LID not found
12748	343569	W69435	6.28	43.64	6.95	0.00	2.00	Aorta	10	45	Aorta	Forebrain
12754	564771	AA136551	47.43	241.38	5.09	0.00	1.00	Lung			Heart	Prostate
12756	343930	W69774	8.59	220.63	23.01	1.00	0.00	Heart			Heart	LID not found
12760	358699	W64247	8.40	83.02	9.88	3.00	2.00	Pericardial	19	103.68	Pericardial	nerve
12762	838761	AA457670	7.49	49.22	6.57	1.00	1.00	Testis			Testis	Brain
12769	491712	AA150459	7.19	42.95	5.97	1.00	0.00	Uterus			Uterus	LID not found
12771	253241	H89293	23.89	182.95	7.66	3.00	0.00	Bone	2	515.69	Ear	Colon
12782	785284	AA454016	5.64	30.64	5.44	1.00	0.00	Colon			Colon	LID not found
12783	324154	W46832	854.32	6022.66	7.05	2.00	0.00	Skin			Skin	Parathyroid
12784	244300	N54783	180.84	1483.88	8.21	1.00	2.00	Tonsil	1	209.39	Tonsil	LID not found
12788	294915	N71463	84.99	822.10	9.67	3.00	3.00	Tonsil			Tonsil	LID not found
12789	293414	N63596	15.24	95.82	6.27	3.00	0.00	Esophagus			Esophagus	LID not found
12797	595238	AA173430	6.79	94.50	13.92	13.00	2.00	Lung			Lung	Synovial mem
12800	841471	AA487241	117.25	1862.86	15.89	0.00	5.00	Germ Cell			Germ Cell	Other
12802	758271	AA423978	5.49	212.28	38.64	1.00	0.00	Testis	X	245.06	Parathyroid	Pool
12815	342027	W60283	453.49	3690.52	8.14	2.00	0.00	Pooled			Pooled	Testis
12817	592777	AA156234	5.12	42.39	8.27	2.00	0.00	Stomach	11	154.35	Stomach	Pancreas

Table 2A

12824	365085	AA024604	20.32	123.60	6.08	1.00	0.00	8	458.89	Parathyroid	Breast	Whole embryo
12828	898070	AA598779	14.30	83.74	5.86	2.00	0.00	6	21.48	Pooled	Gall bladder	CNS
12829	613070	AA181646	10.50	52.70	5.02	0.00	1.00			Cervix	Pool	LID not found
12833	591253	AA160780	43.19	237.98	5.51	0.00	1.00			Pancreas	LID not found	Other
12836	730772	AA436008	5.11	25.75	5.04	1.00	0.00			Testis	LID not found	Other
12840	210952	H66005	20.97	203.41	9.70	0.00	2.00	6	510.58	Germ Cell	Pancreas	Testis
12841	592599	AA159962	4.03	26.45	6.57	1.00	0.00	9	59.14	Whole embryo	Brain	LID not found
12858	773168	AA425665	301.81	2179.31	7.22	0.00	2.00	1	174.53	Small intestine	Germ Cell	Stomach
12861	625563	AA187641	28.78	381.81	13.28	0.00	0.00			Heart	Pool	LID not found
12862	795230	AA453598	2.51	17.74	7.08	2.00	0.00	4	682	Testis	LID not found	Other
12863	348119	W72749	8.05	111.90	13.89	4.00	4.00	4	350.78	Germ Cell	Pool	LID not found
12868	780538	AA429804	17.60	277.80	15.78	2.00	0.00	4	728.84	Peripheral ner	Adipose	Brain
12868	122782	T99719	1.12	5.67	5.07	1.00	0.00	3	726.94	Cervix	Umbilical cord	Ovary
12872	753028	AA438456	5.50	39.79	7.23	4.00	0.00	11	59.6	Foreskin	LID not found	Other
12878	838611	AA456975	49.82	586.03	11.78	1.00	0.00	10	288.35	Testis	Pool	LID not found
12882	781482	AA432121	148.79	796.56	5.35	2.00	0.00			Uterus	Lung	Pool
12883	786037	AA448637	21.36	139.88	6.55	1.00	0.00			Nose	LID not found	Other
12889	250868	N23598	47.56	262.69	5.52	0.00	2.00	7	640.65	CNS	Heart	Bone
12900	128881	R10099	51.07	260.36	5.10	1.00	0.00			Whole embryo	Pool	LID not found
12908	196187	R91949	5.04	34.28	6.80	5.00	0.00	12	245.31	Colon	Pool	LID not found
12910	950596	AA608532	9.17	46.14	5.03	1.00	0.00	8	438.11	Brain	Eye	Breast
12913	254549	N23887	146.51	1355.51	9.25	2.00	0.00	19	305.11	Omentum	Muscle	Other
12914	762270	AA431736	4.74	35.69	7.53	1.00	0.00	11	386.48	Uterus	Foreskin	Pool
12924	128490	R06618	17.70	108.45	6.13	1.00	0.00	14	15.73	Esophagus	Prostate	Cervix
12935	43729	H06195	58.56	309.94	5.29	0.00	1.00	22	16.88	Stomach	Prostate	Foreskin
12936	753113	AA400710	3.73	26.59	7.13	1.00	3.00			Whole embryo	LID not found	Other
12939	758066	AA461486	348.68	2323.83	6.70	1.00	2.00			Uterus	Whole embryo	Pool
12955	35481	R45592	5.81	830.35	142.87	5.00	1.00			Testis	LID not found	Other
12956	743441	AA609364	24.34	273.37	11.23	4.00	5.00	19	305.11	Omentum	Muscle	Other
12964	1409509	AA468929	11.29	109.06	9.66	4.00	2.00	11	386.48	Uterus	Foreskin	Pool
12965	613284	AA455933	3.00	32.08	10.71	0.00	1.00	14	15.73	Esophagus	Prostate	Cervix
12966	30580	R42182	95.52	514.37	5.38	4.00	0.00			Whole embryo	LID not found	Other
12968	267420	N24966	2.86	16.09	5.62	1.00	0.00			Larynx	Esophagus	Blood
12970	788248	AA454085	4.29	26.53	6.18	1.00	0.00			Bone	Liver	Adrenal gland
12973	813286	AA455934	4.73	52.97	11.20	4.00	1.00	1	305.09	Thyroid	Ovary	Aorta
12976	267725	N25578	46.57	308.82	6.63	2.00	0.00	19	158.84	Thyroid	Testis	Prostate
12983	809998	AA454854	11.58	1193.03	103.05	12.00	0.00	1	143.55	Colon	Testis	Prostate
12986	788309	AA450020	45.48	250.49	5.51	0.00	1.00	1	278.61	Adipose	Gall bladder	Blood
12988	1160618	AA877618	8.75	64.58	7.38	3.00	0.00	19	88.05	Gall bladder	Stomach	Pituitary
12989	311195	R41911	3.79	21.94	5.79	2.00	0.00	22	116.78	Thyroid	Foreskin	Gall bladder
12996	1160723	AA877845	23.81	186.20	7.82	3.00	0.00	6		Lymph	Germ Cell	Tonsil
13008	290841	N71982	54.56	395.10	7.24	2.00	1.00			Ovary	Germ Cell	Lung
13010	788355	AA453028	13.17	75.98	5.77	1.00	0.00			Pooled	Blood	Germ Cell
13012	1159983	AA877255	12.77	74.48	5.83	1.00	0.00	14	182.05	Adipose	Brain	Parathyroid
13020	1161013	AA877669	7.72	50.78	6.68	0.00	1.00	19	269.06	Pool	Whole embryo	Pituitary
13028	241176	R39325	8.08	47.71	5.90	3.00	0.00	10	318.05	Stomach	Uterus	Other
13032	258973	N32804	12.88	83.05	6.45	3.00	0.00	14	260.52	Bone	Brain	LID not found
13043	869375	AA253464	53.71	588.53	10.98	0.00	1.00	2	718.37	Cervix	LID not found	Other
13045	813414	AA456648	4.92	30.96	6.29	2.00	1.00	16	194.76	Adrenal gland	CNS	Whole embryo
13046	31281	R42884	16.93	119.56	7.06	1.00	0.00	9	307.02	CNS	Uterus	Lymph
13047	626987	AA190871	6.10	34.82	5.68	1.00	0.00			Thyroid	Cervix	LID not found
13048	276982	N39229	41.72	268.78	6.39	1.00	0.00	18	477.68	Testis	Brain	LID not found
13049	811895	AA454982	9.25	89.92	7.56	0.00	1.00					
13051	669378	AA238788	12.50	64.40	5.15	1.00	0.00					
13054	31475	R42871	2.18	15.53	7.14	1.00	0.00					

Table 2A

13056	278243	N63575	18.60	108.80	5.74	0.00	3.00	15	43.88	Tonsil	CNS	Lymph
13060	291385	N72288	22.12	152.22	6.88	0.00	2.00	2	529.87	Foreskin	LID not found	Other
13068	291416	N72300	13.47	127.57	9.47	2.00	5.00	X	93.95	Cervix	LID not found	Other
13072	310501	N98513	177.72	1041.01	5.86	1.00	2.00	X	93.95	Cervix	LID not found	Other
13102	627272	AA191336	46.56	268.84	5.77	0.00	1.00	X	93.95	Cervix	LID not found	Other
13108	627428	AA190313	136.16	1099.47	8.07	2.00	0.00	X	93.95	Cervix	LID not found	Other
13118	298369	N74108	12.07	121.41	10.06	1.00	1.00	6	314.79	Whole embryo	Tonsil	Lung
13118	850451	AA599104	14.72	80.90	5.49	0.00	1.00	3	120.51	Breast	Whole embryo	Brain
13122	642762	AA486183	12.45	108.06	8.88	2.00	1.00	17	377.83	Stomach	Placenta	Whole embryo
13132	299498	N74958	102.71	549.89	5.35	1.00	2.00	17	377.83	Stomach	Placenta	Whole embryo
13151	278712	N34895	5.13	33.58	6.54	3.00	0.00	X	93.95	Cervix	LID not found	Other
13154	345847	W70342	48.21	450.94	9.78	0.00	3.00	X	93.95	Cervix	LID not found	Other
13157	627555	AA192435	5.91	33.57	5.88	1.00	0.00	3	411.15	Pool	Adipose	Germ Cell
13160	233246	H75778	223.17	1123.34	5.03	1.00	0.00	2	93.04	Testis	LID not found	Other
13164	744391	AA621201	4.51	30.27	8.72	0.00	4.00	2	93.04	Testis	LID not found	Other
13165	611209	AA176413	11.04	58.27	5.10	1.00	0.00	19	87.54	Pool	Pancreas	LID not found
13166	742867	AA406210	3.59	21.13	5.89	0.00	2.00	3	411.15	Pool	LID not found	Other
13167	838855	AA481789	11.11	55.56	5.00	0.00	1.00	2	93.04	Testis	LID not found	Other
13168	209179	H62011	31.94	270.63	8.47	2.00	5.00	19	87.54	Pool	LID not found	Other
13171	772938	AA479928	3.55	21.65	6.09	0.00	1.00	3	583.17	Whole embryo	Heart	LID not found
13178	210531	H65832	30.52	263.32	8.93	1.00	2.00	12	45.7	Pool	LID not found	Other
13183	838518	AA481729	12.87	94.65	7.36	1.00	0.00	3	332.38	Pool	LID not found	Other
13184	233844	H78999	70.37	392.13	5.57	1.00	0.00	11	219.43	Smooth muscle	LID not found	Other
13185	641386	AA487527	26.88	172.32	6.46	1.00	0.00	3	332.38	Pool	LID not found	Other
13189	418113	W90105	25.88	210.34	8.18	0.00	4.00	11	219.43	Smooth muscle	LID not found	Other
13200	241241	H81083	112.11	1388.49	12.37	0.00	3.00	11	59.08	Pool	LID not found	Other
13214	742952	AA405580	81.91	370.06	5.98	1.00	1.00	10	534.94	Testis	LID not found	Other
13216	213575	H70163	25.89	244.04	9.43	2.00	0.00	1	82.43	Pool	LID not found	Other
13223	839545	AA489791	6.08	40.57	6.87	3.00	0.00	20	117.11	Eye	LID not found	Other
13231	839837	AA489826	269.52	1733.60	5.79	2.00	0.00	20	117.11	Eye	LID not found	Other
13235	356940	W92738	20.51	137.42	6.70	0.00	4.00	20	21.8	Heart	LID not found	Other
13239	839855	AA489840	258.81	1713.78	8.87	1.00	0.00	20	21.8	Heart	LID not found	Other
13245	260412	N47208	1.01	10.41	10.26	1.00	0.00	7	501.89	Brain	LID not found	Other
13247	838894	AA480048	15.73	88.46	5.50	1.00	0.00	8	440.43	Blood	Prostate	Other
13252	43885	H05091	8.21	298.53	28.80	1.00	0.00	8	157.53	Brain	CNS	Other
13254	593929	AA189379	7.38	74.50	10.12	7.00	0.00	8	157.53	Brain	CNS	Other
13255	753248	AA408231	10.92	78.95	7.23	2.00	0.00	2	567.28	Pool	LID not found	Other
13258	291700	N73477	62.15	396.07	6.37	3.00	0.00	11	278.14	Small intestine	Heart	Brain
13266	296022	N73571	26.65	229.65	8.62	0.00	5.00	X	245.08	Aorta	Foreskin	Whole embryo
13272	53110	R15832	42.32	280.27	8.62	0.00	1.00	1	171.28	Heart	CNS	Pool
13274	269402	N73807	666.09	4171.04	6.26	1.00	0.00	2	73.17	Brain	Lung	Pool
13276	813154	AA456289	40.77	304.93	7.48	1.00	0.00	X	245.06	Testis	Pool	LID not found
13279	279388	N48698	16.18	348.36	21.55	18.00	1.00	2	578.44	Pool	LID not found	Other
13280	31869	R43017	2.40	20.01	8.33	5.00	0.00	7	152	Brain	LID not found	Other
13282	296574	N73846	8.81	58.24	6.61	1.00	0.00	8	278.48	Pool	LID not found	Other
13283	1031027	AA609881	275.74	1853.58	6.72	1.00	0.00	10	528.52	Lung	UD not found	Other
13286	260273	H98229	12.65	91.65	7.25	0.00	3.00	9	256.24	Brain	Germ Cell	CNS
13287	813311	AA447692	13.42	108.95	6.12	4.00	0.00	8	278.48	Pool	LID not found	Other
13288	37539	R49045	2.29	22.48	9.80	4.00	0.00	8	278.48	Pool	LID not found	Other
13289	298710	N74042	45.55	408.88	8.97	0.00	4.00	10	528.52	Lung	UD not found	Other
13312	44073	H06508	2.56	13.02	5.09	1.00	0.00	8	278.48	Pool	LID not found	Other
13314	299412	N78101	47.58	341.68	7.18	0.00	1.00	10	528.52	Lung	UD not found	Other
13334	629883	AA219033	1.01	5.24	5.20	0.00	1.00	9	256.24	Brain	Germ Cell	CNS
13335	728703	AA398264	8.99	53.89	7.71	5.00	0.00	9	256.24	Brain	Germ Cell	CNS
13336	42824	R60170	2.41	261.02	108.17	6.00	3.00	9	256.24	Brain	Germ Cell	CNS

Table 2A

13344	42271	R60328	0.77	4.78	6.15	2.00	0.00	17	528.04	Germ Cell	Brain	Lung
13348	1493160	AA878880	5.08	348.78	68.25	8.00	0.00	4	453.51	Neural	Smooth muscle	Lymph
13358	785993	AA449329	4.07	43.91	10.80	4.00	0.00			Pooled	Whole embryo	Testis
13381	754485	AA410180	42.71	238.54	5.59	0.00	2.00			Fore skin	Pool	LID not found
13383	39286	R62347	57.51	470.44	8.18	2.00	5.00	8	438.5	Brain	LID not found	Other
13384	1328920	AA757351	11.31	96.68	8.55	0.00	1.00	2	580.63	Synovial mem	Whole embryo	Germ Cell
13388	1389373	AA844124	3.18	19.87	6.25	1.00	1.00	14	276.5	Parathyroid	Testis	Kidney
13370	726551	AA398112	7.11	36.98	5.20	1.00	0.00			Testis	Prostate	LID not found
13371	39191	R54444	3.86	37.45	9.69	2.00	0.00			CNS	Whole embryo	Prostate
13374	785699	AA449332	4.30	23.93	5.57	0.00	1.00			Whole embryo	Kidney	Lung
13375	44387	H08525	2.50	32.65	13.07	0.00	3.00	6	46.35	Breast	Brain	Parathyroid
13376	39311	R51361	3.51	22.87	6.51	1.00	0.00	9	53.37	Cervix	Whole embryo	Brain
13383	44300	H06377	16.57	188.06	11.35	4.00	0.00			Brain	LID not found	Other
13384	1389018	AA655158	3.96	23.98	6.06	1.00	0.00	17	340.79	Germ Cell	Heart	Kidney
13391	44409	H08385	115.41	587.95	5.18	0.00	2.00	3	65.25	Muscle	Blood	Fore skin
13400	1390860	AA844447	10.43	67.84	5.53	0.00	2.00	3	628.88	Parathyroid	Pool	LID not found
13403	39336	R51186	4.17	21.68	5.20	1.00	0.00	2	481.92	Liver	Brain	LID not found
13407	47918	H11968	110.58	678.82	6.15	0.00	4.00	2	521.66	Brain	LID not found	Other
13415	48238	H11887	15.27	136.40	8.93	1.00	1.00	2	554.85	Brain	LID not found	Other
13419	39453	R51631	16.07	85.66	5.33	0.00	1.00			Pooled	Brain	Heart
13423	48033	H11780	6.50	43.58	6.70	2.00	0.00			Germ Cell	Pool	LID not found
13429	768018	AA418747	4.52	34.80	7.70	1.00	0.00	1	118.74	Placenta	Breast	Uterus
13430	785760	AA449686	0.57	6.37	11.27	0.00	1.00			Testis	LID not found	Other
13434	726595	AA397018	77.18	476.55	6.17	0.00	4.00			Pancreas	Blood	Parathyroid
13439	48060	H11631	2.30	23.55	10.25	1.00	0.00	18	185.79	Pancreas	CNS	Parathyroid
13440	1405689	AA890863	28.65	221.88	6.33	0.00	1.00	11	277.15	CNS	LID not found	Other
13443	277487	N56888	11.93	819.40	68.68	0.00	5.00			Whole embryo	Ovary	Heart
13444	344010	W70242	15.18	98.29	8.47	1.00	0.00			Eye	Heart	Pool
13448	358800	W94383	128.17	849.98	6.74	0.00	2.00			CNS	LID not found	Other
13455	278516	N82817	5.74	29.48	5.14	1.00	0.00	3	422.58	CNS	LID not found	Other
13460	344834	W70264	5.08	29.28	5.78	0.00	0.00			Germ Cell	Pool	Pool
13468	366209	AA062985	0.51	3.20	6.32	1.00	0.00			Pooled	Lung	Heart
13480	358872	W94820	48.88	327.02	6.72	0.00	2.00			Heart	LID not found	Other
13487	289742	N62969	12.58	78.05	6.21	4.00	0.00	7	564.37	CNS	Muscle	Fore skin
13490	796916	AA463208	9.48	48.02	5.06	1.00	0.00	22	78.88	Parathyroid	Pool	LID not found
13504	416182	W95106	28.11	252.46	8.98	0.00	4.00			Kidney	Testis	Pool
13512	842895	AA486427	53.22	527.31	9.91	0.00	3.00			Cervix	Pancreas	LID not found
13537	591457	AA160692	2.97	29.48	9.91	4.00	0.00	14	282.48	Stomach	Pancreas	Pool
13538	773189	AA425700	11.48	64.37	5.61	0.00	1.00			Whole embryo	Pool	LID not found
13546	773392	AA425749	6.39	358.22	55.74	0.00	1.00			Whole embryo	Pool	LID not found
13550	786458	AA453823	51.11	410.59	8.03	3.00	2.00			Testis	Pool	LID not found
13552	510908	AA102223	291.83	2073.78	7.11	0.00	2.00			Blood	Colon	Lung
13553	582523	AA180484	12.82	122.27	9.54	1.00	0.00			Blood	Pancreas	LID not found
13560	584587	AA127385	99.85	1264.57	12.66	0.00	4.00			Ear	Ovary	Lung
13563	282827	H99704	345.70	2124.21	6.14	2.00	0.00			Aorta	Uterus	Placenta
13565	625693	AA188460	33.98	235.51	6.83	1.00	0.00	7	15.39	Umbilical cord	Pooled	Germ Cell
13568	566383	AA151775	10.08	184.20	18.27	4.00	0.00	1	192.24	Adipose	Uterus	Pool
13574	811572	AA454595	6.48	72.09	11.12	1.00	0.00	7	521.82	Pool	Ovary	LID not found
13585	593174	AA159605	58.02	510.35	9.11	1.00	5.00	19	112.29	Tonsil	Ovary	Kidney
13591	328889	W45453	5.14	45.23	8.78	1.00	0.00	16	474.57	Fore skin	Parathyroid	Thyroid
13592	642766	AA466135	73.63	573.11	7.78	0.00	4.00	22	136.79	Larynx	Cervix	Prostate
13594	773548	AA428179	130.11	1092.90	8.40	0.00	3.00	12	52.19	Blood	Ovary	Eye
13596	784005	AA443290	20.58	105.81	5.14	0.00	1.00	2	101.7	Brain	Muscle	Lymph
13600	120273	T96986	27.95	217.85	7.80	0.00	3.00			Ovary	Uterus	Brain
13604	132392	R26531	5.41	30.87	5.70	1.00	0.00	13	146.89	Ovary	Placenta	Pool

Table 2A

13605	627211	AA195398	15.59	130.75	8.39	0.00	2.00	1	889.77	Ignore	Umbilical cord	Bone
13617	593690	AA166695	15.84	91.98	5.86	1.00	0.00			Aorta	Ovary	Lung
13627	297830	N69882	207.98	1066.14	5.13	1.00	0.00	10	415.29	Ignore	Cervix	Lung
13628	628851	AA191437	46.43	320.60	6.91	0.00	2.00					
13632	240223	H89505	121.48	1329.27	10.94	0.00	4.00	11	252.9	Fore skin	LID not found	Other
13633	269031	N24028	29.55	208.24	7.05	0.00	3.00			Testis	LID not found	Other
13636	743536	AA609422	17.04	97.59	5.73	1.00	0.00	13	226.98	Pool	LID not found	Other
13639	811943	AA455012	3.47	78.84	22.73	7.00	2.00			Synovial mem	Skin	Germ Cell
13640	753188	AA406346	12.53	92.49	7.38	1.00	1.00	10	104.49	Fore skin	Pool	LID not found
13641	267085	N24848	0.76	4.02	5.30	1.00	0.00			CNS	Pool	LID not found
13645	278058	N51682	7.47	216.16	28.95	7.00	0.00	7	875.52	Parathyroid	Breast	Placenta
13648	753213	AA406363	10.85	57.48	5.30	1.00	0.00	10	364.62	Fore skin	LID not found	Other
13649	268892	N25920	9.58	135.88	14.18	3.00	2.00	10	438.89	Adipose	Ear	Fore skin
13655	798406	AA459944	6.00	40.35	6.72	2.00	0.00	6	117.99	CNS	Pooled	Germ Cell
13656	753236	AA406373	14.78	77.16	5.22	1.00	0.00			Uterus	Brain	Pool
13687	767425	AA417940	5.89	61.85	10.50	1.00	0.00	3	714.15	Stomach	Pool	Heart
13688	664975	AA194633	171.04	1259.80	7.37	1.00	0.00	17	321.68	Nose	Nose	Prostate
13672	753252	AA408233	49.91	687.44	13.37	3.00	0.00	1	553.91	Nose	LID not found	Other
13681	255897	N27366	39.35	218.99	5.59	2.00	0.00	3	391.24	Adipose	Tonsil	Adrenal gland
13687	34942	R43755	10.30	78.89	7.66	0.00	1.00					
13688	753278	AA411856	39.82	224.58	5.64	1.00	0.00	11	319.19	Tonsil	Brain	Pool
13692	685082	AA119441	9.46	48.86	5.17	1.00	0.00			Whole embryo	LID not found	Other
13698	784050	AA443719	882.59	5090.11	5.18	0.00	1.00	12	105.89	Ear	Germ Cell	Testis
13699	753940	AA478106	44.48	281.64	6.33	0.00	1.00			Breast	Fore skin	Muscle
13707	44287	H08249	3.13	21.70	6.82	1.00	0.00	17	24.51	Pancreas	Fore skin	LID not found
13715	768897	AA424754	13.60	92.80	6.82	2.00	0.00	4	615.42	Ear	Cervix	Pool
13726	744438	AA621224	1.95	10.63	5.44	2.00	0.00			Testis	LID not found	Other
13730	788507	AA452572	22.85	134.85	5.90	0.00	0.00			Synovial mem	Fore skin	Lymph
13732	1161830	AA876021	6.33	57.78	9.12	4.00	0.00	9	64.24	Thymus	Whole embryo	CNS
13735	789888	AA443284	10.23	78.18	7.64	1.00	0.00	2	500.76	Adrenal gland	Brain	CNS
13736	276430	N66104	7.57	40.01	5.28	1.00	0.00	14	173.12	Pooled	Brain	Pancreas
13740	33627	R44741	0.43	5.12	11.96	2.00	0.00	14	14.91	CNS	Fore skin	Whole embryo
13746	34149	R44762	2.78	15.74	5.65	1.00	0.00	8	404.08	Muscle	Pancreas	Lung
13754	785524	AA452801	9.32	239.92	25.75	4.00	0.00			Placenta	Parathyroid	CNS
13755	681891	AA256176	2.45	15.40	6.29	1.00	0.00			Whole embryo	Blood	
13758	31759	R43008	24.45	201.18	8.23	0.00	2.00					
13762	786541	AA452818	35.18	344.69	9.40	1.00	2.00	20	328.88	Larynx	Thyroid	Cervix
13770	786554	AA452822	178.95	917.74	5.13	1.00	0.00	3	160.11	Larynx	Fore skin	CNS
13774	31972	R43020	16.56	91.50	5.32	0.00	2.00	22	-13.35	Ignore	CNS	Ovary
13775	741888	AA402915	21.30	467.98	21.97	3.00	0.00	18	287.82	Fore skin	Pooled	Pancreas
13783	773367	AA425684	3.82	38.16	10.81	1.00	0.00	10	362.35	CNS	Whole embryo	Pool
13786	788558	AA452824	18.81	145.87	8.68	0.00	1.00	11	352.44	Larynx	Bone	
13789	813513	AA455033	6.55	75.47	11.53	0.00	2.00					
13794	788575	AA452877	8.87	100.11	11.55	3.00	0.00	21	242.21	Small intestine	Stomach	Germ Cell
13800	324492	W51794	16.27	254.33	15.83	1.00	2.00			Cervix	Liver	Adipose
13801	811954	AA456635	13.44	81.81	6.07	0.00	1.00	1	750.98	Skin	Placenta	Adrenal gland
13804	844816	AA773894	78.03	864.87	8.41	2.00	0.00	20	40.68	Umbilical cord	Pool	Kidney
13808	346009	W72140	13.64	176.48	12.94	0.00	0.00	3	56.14	Parathyroid	LID not found	Other
13815	882522	AA676466	59.77	639.22	10.69	1.00	0.00					
13820	853450	AA773883	3.25	27.09	8.33	1.00	0.00					
13822	34326	R44396	12.13	63.85	5.27	1.00	0.00					
13823	843386	AA469383	4.57	30.17	6.80	2.00	0.00					
13826	589027	AA149117	45.23	334.48	7.39	2.00	0.00					
13832	320425	W04695	48.73	322.51	6.82	0.00	4.00					

Table 2A

13842	568466	AA152340	22.14	132.95	6.01	1.00	0.00	16	191.21	Stomach	Pooled	Blood
13848	322447	W16425	12.79	103.45	8.09	0.00	3.00			Parathyroid	LID not found	Other
13855	276861	N34933	52.25	323.77	8.20	0.00	2.00	16	202.63	CNS	LID not found	Other
13860	259459	N78133	10.26	59.41	5.69	0.00	1.00			Heart	Whole embryo	Lung
13860	327480	W20462	2.66	63.69	23.65	5.00	0.00			Pooled	Kidney	Heart
13862	116049	T92200	209.53	1879.75	8.97	2.00	0.00	X	245.05	Lung	Pool	LID not found
13864	300000	N78889	87.83	683.58	7.78	0.00	2.00			Prostate	Lung	LID not found
13868	627633	AA186281	14.23	98.63	6.79	0.00	1.00			Cervix	Uterus	Kidney
13894	365508	AA009738	21.60	108.48	5.02	0.00	1.00	17	41.55	Stomach	Colon	Kidney
13896	327732	W23581	22.62	221.25	9.78	0.00	3.00			Heart	LID not found	Other
13914	429050	AA005135	5.20	41.95	8.08	1.00	0.00			Cervix	Colon	Lung
13916	841280	AA487192	211.35	1651.17	7.81	0.00	2.00	6	522.76	Pool	LID not found	Other
13924	255285	N23885	8.13	70.09	8.62	1.00	0.00	12	277.57	Pooled	Placenta	Lung
13925	587430	AA132524	172.73	1402.32	8.12	2.00	0.00	3	537.68	Nose	Stomach	Breast
13926	743030	AA406081	0.96	5.95	6.17	1.00	0.00			Pool	Colon	LID not found
13927	839903	AA490056	8.09	77.47	8.53	4.00	0.00			Testis	LID not found	Other
13928	244194	N51030	23.43	153.04	6.53	3.00	0.00	1	88.45	Eye	Pool	LID not found
13941	490991	AA120881	23.65	144.84	6.12	0.00	2.00			Uterus	LID not found	Other
13945	487938	AA044741	167.46	911.37	5.44	0.00	1.00			Eye	LID not found	Other
13947	360177	AA012911	21.41	148.56	6.85	3.00	0.00			Eye	LID not found	Other
13952	245324	N53456	10.77	57.25	5.32	0.00	1.00	5	504.31	Eye	LID not found	Other
13955	360355	AA013353	3.93	24.42	6.21	3.00	0.00			Eye	Pool	LID not found
13957	490188	AA121271	8.89	54.09	6.08	0.00	1.00	15	347.96	Uterus	Pool	LID not found
13959	490434	AA122078	5.27	254.25	48.25	3.00	0.00	4	648.75	Uterus	Pool	LID not found
13963	840503	AA485886	165.98	1064.67	5.43	1.00	0.00	12	101.74	Synovial mem	Forestin	CNS
13965	489331	AA121518	102.71	588.75	5.73	0.00	3.00	X	78.07	Uterus	LID not found	Other
13968	246552	N57659	34.27	210.83	8.15	1.00	4.00	18	82.53	Pool	LID not found	Other
13969	487068	AA045300	86.72	1083.12	15.47	3.00	0.00			Lymph	Uterus	Kidney
13970	368508	AA146979	116.80	1091.03	9.16	0.00	3.00	15	164.8	Colon	Kidney	LID not found
13980	796095	AA460376	878.64	5263.63	5.99	2.00	0.00			Whole embryo	LID not found	Other
13986	592530	AA161161	16.12	84.13	5.22	1.00	0.00	10	428.78	Head and nec	Ear	Pancreas
13988	841674	AA487233	104.88	871.77	9.27	0.00	4.00	5	529.13	Spleen	Tonsil	Ovary
13996	810209	AA484522	9.79	59.00	6.03	1.00	0.00			Spleen	Tonsil	Ovary
14005	511833	AA126958	14.87	204.49	13.66	1.00	0.00	9	121.08	Whole embryo	Colon	Pool
14006	809722	AA455463	1027.59	5787.16	5.64	2.00	0.00	22	71.14	Ovary	Colon	Lung
14016	233759	H64591	5.80	98.24	17.18	6.00	0.00	15	240.37	Ovary	Heart	Lung
14018	267864	N25650	473.68	6608.43	13.85	0.00	4.00	18	162.07	Heart	Heart	Lung
14028	767206	AA424586	159.16	962.04	6.04	0.00	2.00	3	575.4	Esophagus	Adipose	Gall bladder
14040	29920	R42312	2.87	38.06	12.83	8.00	0.00			Eye	Spleen	Brain
14047	786448	AA459883	4.04	26.83	6.64	2.00	0.00	18	316.2	Eye	Umbilical cord	Ear
14051	1031182	AA608955	8.11	44.32	5.46	1.00	0.00			Testis	LID not found	Other
14058	306276	N90595	87.43	449.68	6.67	0.00	2.00	3	459.05	Parathyroid	LID not found	Other
14068	766550	AA452125	29.77	157.14	5.28	0.00	1.00			Parathyroid	LID not found	Other
14072	43532	H05535	5.44	33.82	6.22	1.00	1.00	1	126.05	Brain	LID not found	Other
14080	37310	R49597	7.33	820.13	125.57	22.00	4.00			Parathyroid	Pool	Brain
14082	308620	N82804	248.86	1489.82	5.99	0.00	1.00	2	743.9	Lung	LID not found	Other
14088	31979	R43026	67.40	515.85	7.65	0.00	2.00	12	229.98	Brain	LID not found	Other
14090	309638	N94447	13.45	88.90	6.68	0.00	1.00			Testis	LID not found	Other
14101	121154	T88935	18.25	182.86	11.87	7.00	1.00			Pool	LID not found	Other
14106	309119	N88238	122.99	755.21	6.14	0.00	2.00			Stomach	Lung	LID not found
14108	122752	T89043	76.46	740.06	9.68	0.00	4.00	6	495.58	Pool	LID not found	Other
14112	28737	R40835	3.86	45.19	11.40	1.00	0.00	3	455.08	Brain	LID not found	Other
14114	728598	AA398141	39.58	227.49	5.75	1.00	0.00	2	97.87	Adrenal gland	Foreskin	Placenta
14124	1468220	AA884897	15.84	86.89	5.56	0.00	1.00			Cervix	Tonsil	Prostate

Table 2A

14129	754582	AA406311	6.76	37.93	5.61	1.00	0.00	17	281.02	Lymph	Blood	Tonsil
14131	41869	R66438	2.46	39.09	15.90	2.00	0.00	4	168.31	Parathyroid	Pool	Prostate
14137	754588	AA406201	44.76	236.76	5.29	1.00	0.00	3	347.25	Bone	Liver	Germ Cell
14145	754591	AA406320	32.83	204.88	8.24	1.00	0.00	3	395.51	Pooled	Salivary gland	Placenta
14156	1412238	AA404416	7.88	392.52	51.11	8.00	0.00	7	592.03	Pancreas	Kidney	Ear
14164	1412245	AA4044831	4.94	164.96	33.42	2.00	0.00	2	241.71	Pancreas	Small intestine	Colon
14169	754853	AA411204	4.06	27.03	6.66	0.00	2.00	2	323.47	CNS	Thyroid	Whole embryo
14172	1412300	AA4044864	3.49	67.84	25.19	2.00	0.00	15	145.08	Brain	LID not found	Other
14174	785690	AA449362	89.12	457.43	5.13	1.00	0.00	3	-12.44	Brain	LID not found	Other
14177	754654	AA411607	0.75	15.53	20.78	3.00	0.00	2	410.63	Breast	Brain	LID not found
14190	765910	AA449481	10.64	216.84	20.38	0.00	3.00	3	357.99	Heart	Whole embryo	Lung
14191	46887	H18098	15.52	103.43	6.66	2.00	0.00	2	668.82	Bone	Heart	LID not found
14195	41913	RS9608	1.02	7.88	7.72	2.00	0.00	7	504.11	Brain	Aorta	Lung
14198	785913	AA449490	20.80	104.27	5.01	1.00	0.00	X	83.98	Ear	Pooled	Whole embryo
14199	48610	H18179	1.79	29.99	16.80	1.00	0.00	10	374.69	Eye	Aorta	CNS
14205	768264	AA424944	81.52	487.39	5.73	0.00	2.00	3	200.65	Eye	Blood	Muscle
14207	48953	H18725	218.35	1424.48	6.52	0.00	3.00	22	37.19	Heart	LID not found	Other
14208	448432	AA777551	1.24	9.58	7.71	2.00	0.00	15	190.28	Pool	LID not found	Other
14222	587333	AA132887	155.88	1414.68	9.08	2.00	0.00	4	349.65	Eye	LID not found	Other
14232	359072	W82315	178.96	1345.43	7.80	2.00	0.00	15	14.4	Pool	Testis	Breast
14235	840728	AA487848	297.87	3930.50	13.20	1.00	0.00	5	57.43	Eye	LID not found	Other
14236	344959	W72870	7.12	43.30	6.08	1.00	0.00	19	250.6	Testis	LID not found	Other
14238	564847	AA128217	6.16	71.52	11.61	1.00	0.00	20	120.04	Colon	Lung	LID not found
14242	897865	AA599640	11.91	93.68	7.87	1.00	0.00	X	121.77	Stomach	Tonsil	Brain
14244	345116	W72820	4.88	35.24	7.07	1.00	0.00	21	154.34	Neural	Brain	Eye
14251	280468	N51585	102.34	542.72	5.30	0.00	1.00	2	245.08	Testis	LID not found	Other
14254	838998	AA487287	192.85	1473.86	7.64	2.00	3.00	2	668.82	Bone	Heart	LID not found
14270	784154	AA432096	0.71	6.78	9.54	2.00	0.00	7	504.11	Brain	Aorta	Lung
14274	730408	AA468964	27.95	179.30	6.41	2.00	0.00	X	83.98	Ear	Pooled	Whole embryo
14285	278144	N63516	11.24	59.22	5.27	0.00	1.00	10	374.69	Eye	Aorta	CNS
14289	283581	N50702	5.77	131.13	22.72	1.00	0.00	3	200.65	Eye	Blood	Muscle
14300	344590	W73597	167.79	1041.92	6.21	1.00	0.00	22	37.19	Heart	LID not found	Other
14304	427754	AA002228	35.43	189.93	5.36	2.00	1.00	15	190.28	Pool	LID not found	Other
14307	291947	N73083	0.78	8.81	11.25	1.00	0.00	4	349.65	Eye	LID not found	Other
14316	122008	T83555	20.15	102.07	5.07	0.00	1.00	15	190.28	Pool	LID not found	Other
14322	780944	AA428807	6.11	44.05	7.21	0.00	1.00	4	349.65	Eye	LID not found	Other
14324	122435	T89243	7.95	48.95	6.16	2.00	0.00	15	14.4	Pool	Testis	Breast
14325	737246	AA426028	10.03	103.89	10.37	0.00	3.00	22	68	Nose	Kidney	Breast
14328	836230	AA456874	33.42	182.38	5.46	1.00	0.00	5	57.43	Eye	LID not found	Other
14328	241447	H80407	32.70	301.43	9.22	0.00	5.00	19	250.6	Testis	LID not found	Other
14331	255651	N27637	218.53	1509.71	6.88	3.00	0.00	20	120.04	Colon	Lung	LID not found
14334	810884	AA459403	272.32	1472.70	5.41	0.00	3.00	X	121.77	Stomach	Tonsil	Brain
14337	608521	AA187565	91.21	987.23	10.80	0.00	5.00	5	57.43	Eye	LID not found	Other
14342	785534	AA458648	11.43	68.92	5.85	0.00	1.00	19	250.6	Testis	LID not found	Other
14345	608620	AA187589	57.11	556.28	9.74	0.00	3.00	20	120.04	Colon	Lung	LID not found
14350	785580	AA458689	15.24	92.28	5.40	0.00	2.00	X	121.77	Stomach	Tonsil	Brain
14352	586845	AA133554	1.28	12.88	10.13	3.00	0.00	5	57.43	Eye	LID not found	Other
14353	610006	AA189173	7.02	51.88	7.38	3.00	0.00	19	250.6	Testis	LID not found	Other
14355	488553	AA047275	34.00	231.39	6.81	2.00	0.00	20	120.04	Colon	Lung	LID not found
14356	126540	R06754	136.48	920.17	6.74	0.00	2.00	X	121.77	Stomach	Tonsil	Brain
14360	347740	W81624	25.71	153.68	5.98	0.00	3.00	5	57.43	Eye	LID not found	Other
14367	429122	AA004803	6.43	41.25	6.41	1.00	0.00	19	250.6	Testis	LID not found	Other
14374	785784	AA458851	161.02	921.82	5.72	1.00	0.00	20	120.04	Colon	Lung	LID not found
14381	628701	AA218873	138.33	845.72	6.11	1.00	0.00	21	154.34	Neural	Brain	Eye
14384	744360	AA621103	9.34	61.28	6.56	0.00	1.00					

Table 2A

14389	629385	AA219047	374.87	2285.82	6.04	1.00	1.00	3	697.77	Eye	LID not found Other
14395	244050	N34042	1.92	23.21	12.07	4.00	0.00		Pool	Pool	LID not found Other
14396	197102	R93409	5.52	34.05	6.17	1.00	0.00	2	75.2	Pool	LID not found Other
14397	629907	AA219230	56.06	346.97	6.19	1.00	2.00		Eye	Eye	LID not found Other
14401	257170	N30557	278.83	1681.08	6.03	0.00	2.00	1	589.13	Nose	Gall bladder Adrenal gland
14421	283919	N52876	5.54	31.91	5.78	1.00	0.00		CNS	CNS	LID not found Other
14422	744611	AA621294	5.12	26.38	5.15	1.00	0.00		Testis	Testis	LID not found Other
14428	665316	AA185318	48.96	282.39	5.77	1.00	0.00		Colon	Colon	Heart Pool
14431	811907	AA454654	5.78	71.10	12.34	5.00	0.00		Larynx	Larynx	Umbilical cord
14434	795376	AA453494	25.62	163.18	6.37	0.00	2.00		Testis	Testis	LID not found Other
14435	753633	AA478596	22.78	183.14	6.04	2.00	0.00	6	477.99	CNS	Pancreas Kidney
14437	283682	N52935	5.72	32.63	5.74	1.00	0.00		CNS	CNS	LID not found Other
14439	812074	AA455986	4.75	96.82	20.82	5.00	1.00		Pancreas	Pancreas	Liver Kidney
14440	753376	AA411685	5.97	107.04	17.84	8.00	1.00		Brain	Brain	Ovary Tonsil
14444	665379	AA194993	8.93	53.68	6.01	1.00	0.00	6	18.39	CNS	LID not found Other
14445	283688	N52938	165.50	830.04	5.02	1.00	0.00	9	355.78	CNS	Pool Aorta
14446	198635	R96922	9.36	132.05	14.11	0.00	1.00	19	247.58	Placenta	Thyroid Brain
14463	353329	R45517	7.80	45.38	5.75	1.00	0.00	5	527.16	Neural	LID not found Other
14465	243477	N33610	4.20	21.75	5.16	1.00	0.00	10	372.38	Pool	LID not found Other
14466	811503	AA454616	110.50	665.02	6.02	1.00	0.00	15	118.58	Thyroid	Pool Pancreas
14467	752836	AA418603	104.18	731.42	7.02	2.00	0.00		Ovary	Ovary	LID not found Other
14474	809663	AA455130	6.71	82.04	12.22	2.00	0.00	3	20.28	Thymus	Pool CNS
14476	665445	AA195080	38.85	203.49	5.26	1.00	0.00	4	450.16	Pool	LID not found Other
14477	247265	N54081	32.56	190.33	5.85	0.00	3.00		Pool	Pool	LID not found Other
14480	753596	AA478717	9.78	72.82	7.45	1.00	3.00		Eye	Eye	Head and nec Brain
14486	731469	AA412417	10.58	70.67	7.53	3.00	0.00	11	271.39	Eye	LID not found Other
14487	838778	AA457576	0.91	4.90	5.40	2.00	0.00	3	472.27	Brain	Blood Prostate
14481	34869	R44447	1.53	15.24	9.96	1.00	1.00	1	628	Pool	LID not found Other
14493	247698	N54274	498.68	2822.94	5.65	1.00	0.00	4	639.11	Pool	Germ Cell
14495	37823	R58473	0.88	7.17	8.15	2.00	0.00		Uterus	Uterus	Brain Pool
14498	788841	AA448832	8.34	105.02	12.59	3.00	0.00	8	470.89	Uterus	Umbilical cord Pooled
14499	682072	AA256464	22.68	156.78	6.91	0.00	4.00	7	94.72	Cervix	LID not found Other
14500	824744	AA187143	4.57	87.39	19.10	3.00	5.00	14	82.24	Brain	LID not found Other
14502	34745	R44409	3.15	140.68	44.81	5.00	0.00	10	250.28	Whole embryo	LID not found Other
14506	811983	AA458654	35.57	186.88	5.31	0.00	1.00	3	628.88	Adipose	Umbilical cord Pancreas
14508	788687	AA448847	8.99	68.23	6.83	1.00	0.00		Colon	Colon	Brain LID not found
14510	625011	AA181023	9.30	735.73	79.13	22.00	3.00	7	675.52	Head and nec Esophagus	Pancreas
14511	897768	R44428	3.74	24.36	6.62	1.00	0.00	22	37.18	Forasquin	Pool
14512	345077	AA598507	3.33	24.03	7.22	3.00	0.00	20	328.54	Thymus	Bone Foreskin
14526	34901	R45114	28.49	356.01	12.07	2.00	0.00	9	422.75	Gall bladder	Spleen Skin
14528	366518	AA026605	26.52	173.44	6.54	0.00	3.00	X	141.57	Peripheral nec Tonsil	Colon
14529	812012	AA455882	6.67	52.07	5.87	1.00	0.00	4	562.82	Adipose	Gall bladder Breast
14536	377468	AA055440	34.43	282.70	8.21	1.00	1.00		217.43	Cervix	Forasquin Lymph
14544	378935	AA778382	227.48	1482.41	6.52	1.00	0.00		Neural	Neural	Pool Tonsil
14552	378461	AA775816	21.01	1900.75	90.45	3.00	0.00	5	626.75	Eye	Parathyroid Brain
14569	812053	AA455980	8.88	56.95	6.41	1.00	0.00		Prostate	Prostate	Ovary
14580	1155071	AA706301	6.78	43.53	6.42	0.00	2.00	2	294.33	CNS	Blood Parathyroid
14584	385003	AA709143	165.20	1387.94	8.40	0.00	2.00		194.51		
14585	812069	AA455984	20.32	132.57	6.53	0.00	1.00				
14590	35300	R43786	24.85	258.12	10.30	1.00	0.00				
14600	594693	AA185313	13.11	78.61	5.99	1.00	0.00				
14618	797057	AA463249	234.02	1486.72	6.35	0.00	2.00				
14620	843251	AA468648	104.73	726.20	6.93	5.00	0.00				



Table 2A

14624	326287	W31919	7.82	216.37	27.71	1.00	2.00	Pancreas	LID not found Other
14628	784142	AA432081	230.65	1953.21	8.47	0.00	2.00	Esophagus	Parathyroid Whole embryo
14638	950924	AA608728	473.38	3028.61	6.40	1.00	0.00	433.36	246.56
14644	565110	AA128462	37.94	398.16	10.49	0.00	2.00	Lung	LID not found
14654	490985	AA120886	19.58	426.06	21.76	8.00	0.00	Stomach	Prostate
14659	784272	AA447478	21.78	111.38	6.11	1.00	0.00	Pool	LID not found Other
14662	595209	AA173411	31.84	193.81	6.09	1.00	0.00	Small intestine	Uterus
14666	119330	T84556	10.77	55.94	5.18	1.00	0.00	Lung	LID not found Other
14675	285780	N84145	243.53	1407.07	5.78	2.00	0.00	Ear	Foreskin LID not found
14678	757337	AA437099	5.87	45.68	7.82	2.00	0.00	CNS	Muscle Parathyroid
14684	272552	N35884	308.13	1763.62	5.72	1.00	1.00	173.27	245.06
14686	568501	AA151917	27.97	140.21	5.01	0.00	1.00	Aorta	Foreskin Whole embryo
14687	277083	N38603	24.91	132.91	5.33	1.00	0.00	123.62	Colon Heart
14688	322033	W37833	18.28	93.04	5.09	1.00	1.00	546.25	LID not found Other
14698	121580	T97921	48.78	449.24	9.21	0.00	2.00	530.53	LID not found Other
14702	730742	AA435988	5.55	33.55	6.05	1.00	0.00	Testis	LID not found Other
14704	757205	AA443978	15.48	179.89	11.62	4.00	0.00	Prostate	Uterus
14715	489109	AA056484	33.07	172.98	5.23	1.00	0.00	Lung	LID not found Other
14719	840470	AA485668	6.15	45.03	7.32	2.00	0.00	Lung	LID not found Other
14727	840514	AA485969	154.78	1177.04	7.81	2.00	0.00	Colon	LID not found Other
14729	510397	AA053682	24.74	130.91	5.29	2.00	0.00	CNS	Pooled Foreskin
14731	897722	AA598883	9.99	65.30	6.54	1.00	0.00	330.19	Uterus Testis
14735	843058	AA486504	4.62	61.72	13.36	8.00	0.00	Cervix	Prostate
14736	111735	T91088	8.14	59.84	7.35	1.00	0.00	Cervix	Pool
14743	843278	AA488659	15.07	133.96	8.89	2.00	3.00	654.24	Kidney Whole embryo
14749	564698	AA126318	8.82	49.54	5.82	0.00	1.00	327.15	Blood Aorta
14755	594946	AA172039	23.27	150.86	6.48	0.00	1.00	327.49	Prostate Pancreas
14766	731445	AA412443	4.24	45.21	10.65	2.00	0.00	Testis	Pool
14770	786227	AA460669	9.45	56.30	5.96	0.00	1.00	Adipose	Whole embryo
14772	839037	AA487501	20.39	289.95	14.22	4.00	5.00	Neural	Eye Pool
14773	526567	AA128407	3.24	75.85	23.42	8.00	0.00	22.82	Pancreas Blood
14775	841016	AA486858	25.12	128.29	5.11	1.00	0.00	149.53	Colon LID not found Other
14784	122872	R00130	177.91	1531.72	8.61	1.00	2.00	719.04	Testis LID not found
14793	1031580	AA609310	42.09	239.38	5.69	1.00	0.00	Muscle	LID not found Other
14794	327748	W23441	50.53	342.75	6.78	0.00	4.00	347.4	LID not found Other
14797	126549	R08860	240.13	1258.08	5.24	1.00	0.00	680.45	Parathyroid LID not found Other
14802	321310	W32192	12.17	62.84	5.15	0.00	1.00	Adipose	Brain LID not found
14808	35285	R45579	83.13	333.15	5.28	0.00	1.00	LID not found Other	LID not found Other
14810	324111	W46575	68.82	540.10	7.76	0.00	2.00	385.62	Pool LID not found Other
14813	127192	R08260	30.57	211.52	6.92	0.00	3.00	110.93	Small intestine Gall bladder
14814	287865	N25657	171.67	1025.19	5.97	0.00	2.00	Brain	Prostate LID not found
14816	37387	R49650	4.06	26.82	6.55	1.00	0.00	Pool	LID not found Other
14817	244959	N54825	207.89	1211.40	5.83	0.00	2.00	151.92	Foreskin Heart
14818	340737	W56308	211.56	1740.75	8.23	0.00	2.00	Testis	LID not found Other
14819	1030855	AA621781	73.00	963.30	13.20	2.00	1.00	21.19	Brain Pancreas LID not found
14824	44156	H08157	2.53	22.90	9.06	2.00	0.00	Foreskin	LID not found Other
14825	265845	N25338	300.18	2126.96	7.09	2.00	0.00	5	840.93
14826	133947	R27619	7.30	61.73	11.20	1.00	0.00	96	Parathyroid Placenta Foreskin
14832	31584	R42058	10.98	354.13	32.25	2.00	0.00	4	88.95
14833	276397	N40180	6.17	45.40	5.55	1.00	0.00	CNS	LID not found Other
14837	128777	R18983	73.32	449.88	6.14	0.00	1.00	Pool	LID not found Other
14840	37814	R69355	3.76	22.23	5.91	1.00	0.00	Brain	Pool LID not found
14842	338178	W80473	7.04	70.55	10.03	6.00	0.00	Pooled	Germ Cell Pancreas
14845	130392	R21741	175.15	2089.40	11.93	0.00	2.00	Placenta	LID not found Other

Table 2A

14849	278137	N63520	15.76	89.85	5.70	0.00	1.00	CNS	LID not found Other
14851	1030959	AA620359	6.27	43.04	6.86	1.00	0.00	Testis	LID not found Other
14852	115277	W66932	8.24	101.29	12.29	0.00	1.00	Pancreas	Whole embryo
14858	343174	W67536	7.55	67.51	8.94	0.00	0.00	Uterus	Heart
14866	343235	W67368	545.65	4144.44	7.60	2.00	0.00	Testis	LID not found Other
14867	1049168	AA620669	7.76	69.08	8.90	0.00	1.00	Brain	Testis
14872	38887	R51514	2.01	16.57	8.24	3.00	0.00	Whole embryo	Germ Cell
14876	785994	AA449321	9.14	59.24	6.48	1.00	0.00	Thyroid	Stomach
14878	281859	N48050	8.44	62.82	7.44	1.00	0.00	Testis	Pool
14882	786053	AA448653	10.04	59.29	5.90	0.00	0.00	CNS	Tonsil
14884	450574	AA704255	5.74	35.34	6.18	2.00	0.00	Brain	LID not found Other
14886	726695	AA388365	28.86	179.72	6.69	1.00	0.00	Testis	Eye
14887	42302	R61700	114.02	588.02	5.23	1.00	0.00	Brain	LID not found Other
14892	451504	AA707321	55.57	386.48	6.60	2.00	0.00	Pool	Colon
14893	767113	AA424537	12.07	124.06	10.28	0.00	2.00	Parathyroid	Smooth muscle
14896	1416782	AA084557	33.94	485.07	14.59	4.00	2.00	Esophagus	Pool
14921	767128	AA424534	7.93	41.85	5.28	0.00	1.00	Esophagus	Kidney
14924	1434948	AA857131	37.40	224.00	5.69	2.00	0.00	Esophagus	Synovial mem
14928	452588	AA778919	2.59	18.03	6.19	3.00	0.00	Testis	Brain
14930	726709	AA388267	3.56	21.44	6.02	2.00	0.00	Brain	LID not found Other
14931	42330	R61187	34.27	171.71	5.01	1.00	0.00	Larynx	Blood
14932	1434905	AA857101	4.95	280.19	58.58	10.00	0.00	Testis	Brain
14935	847444	AA199666	38.81	272.37	7.02	2.00	0.00	Testis	LID not found
14947	42660	R61231	2.54	14.11	5.56	1.00	0.00	Brain	LID not found Other
14955	42452	R61297	663.79	3511.35	6.25	1.00	0.00	Adipose	Pool
14957	768417	AA495835	10.49	66.32	6.32	0.00	1.00	Muscle	Whole embryo
14958	788154	AA448855	3.25	16.48	5.08	0.00	1.00	Lung	Placenta
14968	482595	AA705112	9.99	89.27	8.94	2.00	0.00	Testis	LID not found Other
14969	767178	AA424562	11.40	105.55	8.26	4.00	0.00	Germ Cell	Prostate
14970	726731	AA398285	127.67	879.81	6.89	2.00	0.00	Pool	Heart
14982	345761	W72671	5.80	43.08	7.43	3.00	0.00	CNS	Heart
14984	427677	AA001879	200.05	1430.29	7.15	2.00	0.00	CNS	Heart
14987	283744	N50740	9.51	57.47	6.05	1.00	0.00	Cervix	Whole embryo
14996	346368	W74257	24.68	166.07	6.32	0.00	2.00	Cervix	Heart
14998	898050	AA588947	285.58	1540.53	5.80	1.00	0.00	Cervix	Heart
15004	358217	W95636	11.19	145.79	13.03	3.00	0.00	Cervix	Heart
15011	281834	N51068	11.16	118.67	10.63	1.00	0.00	CNS	Heart
15012	346523	W73984	6.08	33.18	5.46	0.00	1.00	Heart	LID not found Other
15016	427897	AA001824	15.96	118.08	7.40	4.00	0.00	Pool	LID not found Other
15038	212784	H86891	321.40	1884.75	6.18	2.00	1.00	CNS	LID not found Other
15039	277871	N84198	30.02	151.32	5.04	0.00	1.00	Smooth muscle	Germ Cell
15040	784183	AA446661	21.16	207.35	9.80	0.00	1.00	Nose	LID not found Other
15046	255295	N23717	947.84	6391.17	6.74	1.00	0.00	Nose	LID not found Other
15048	428592	AA004487	5.52	44.94	8.14	2.00	0.00	Pool	Ear
15055	839527	AA491457	3.26	19.72	8.05	1.00	0.00	Eye	Prostate
15073	594083	AA188840	36.84	185.56	5.04	0.00	1.00	Ovary	LID not found Other
15084	591814	AA143487	21.11	117.03	5.54	2.00	0.00	Pancreas	Colon
15085	829994	AA219172	31.32	237.88	7.60	0.00	2.00	Eye	LID not found Other
15087	213464	H72232	713.96	3960.65	5.13	1.00	0.00	Pool	LID not found Other
15101	838653	AA481788	88.62	726.33	6.23	2.00	0.00	Tonsil	Eye
15113	594694	AA171760	8.39	166.58	19.86	16.00	2.00	Stomach	Colon
15114	782171	AA431210	6.49	32.60	5.02	0.00	1.00	Testis	LID not found Other
15121	594758	AA172056	19.89	233.32	11.73	3.00	0.00	Cervix	Ovary
15133	748605	AA821281	10.20	81.19	6.00	1.00	0.00	Cervix	Testis
15140	213682	H72279	22.29	184.56	8.28	2.00	0.00		Breast

Table 2A

15142	785887	AA460147	9.82	51.11	5.20	1.00	0.00	Testis	LID not found	Other
15143	345711	W71894	18.21	123.62	6.79	0.00	1.00	27.1 Syriovial mem	Pooled	Breast
15147	290227	N62271	86.33	508.21	5.66	2.00	0.00	CNS	Tonsil	LID not found
15153	609950	AA174106	5.49	30.43	5.54	2.00	0.00	Eye	LID not found	Other
15154	782409	AA431771	62.94	340.33	5.41	1.00	0.00	Testis	LID not found	Other
15155	286961	N62712	27.25	139.37	5.11	1.00	0.00	424.67 CNS	LID not found	Other
15156	569440	AA148862	21.48	150.53	7.01	2.00	0.00	487.5 Sbrnach	Uterus	Muscle
15165	838478	AA457517	36.05	285.84	7.93	2.00	0.00	437.97 Salivary gland	Adipose	CNS
15187	697422	AA489463	24.33	173.86	7.15	2.00	1.00	85.75 CNS	Spleen	Lung
15171	27404	R40031	3.13	16.95	5.42	2.00	0.00	672.07 Germ Cell	Eye	
15178	753626	AA478852	11.24	57.94	5.15	1.00	0.00	88.61 Pooled	Fore skin	Pool
15182	798328	AA461317	128.49	883.48	5.32	0.00	1.00	245.06 CNS	Whole embryo	LID not found
15190	814286	AA458993	8.18	43.55	5.32	1.00	0.00	15.07 Neural	Pooled	Placenta
15191	753385	AA410345	18.00	85.86	5.37	1.00	0.00	Pooled	Pancreas	Kidney
15205	245683	N55361	140.20	742.30	5.29	1.00	0.00	337.25 Spleen	CNS	Tonsil
15208	753684	AA408589	17.20	93.81	5.45	0.00	1.00	339.35 Placenta	Parathyroid	Fore skin
15215	813169	AA458318	14.56	80.68	5.54	1.00	0.00	Muscle	Kidney	
15223	752625	AA419608	24.86	180.85	7.27	1.00	1.00	85.2 Parathyroid	Nose	CNS
15225	243347	N38980	27.02	182.04	7.11	0.00	1.00	179.2 Pod	LID not found	Other
15229	277327	N57483	1.22	7.21	5.93	1.00	0.00	164.67 CNS	LID not found	Other
15231	767082	AA424511	7.92	44.22	5.58	0.00	1.00	Pool	LID not found	Other
15241	277039	N38577	54.60	329.80	6.04	0.00	3.00	CNS	LID not found	Other
15243	685385	AA195021	55.28	280.51	5.07	0.00	1.00	403.8 Pancreas	Cervix	Umbilical cord
15247	728483	AA398289	38.53	317.18	8.23	0.00	1.00	Uterus	Pool	LID not found
15281	247177	N57606	14.34	133.20	9.29	1.00	1.00	Pool	LID not found	Other
15282	788155	AA461090	200.57	1105.18	5.51	2.00	2.00	CNS	LID not found	Other
15285	812088	AA456001	7.45	89.03	11.86	2.00	0.00	Whole embryo	LID not found	Other
15288	187814	R83757	9.04	78.31	8.44	3.00	0.00	Kidney	Germ Cell	Ovary
15292	1091543	AA599311	325.95	2098.08	6.43	2.00	0.00	823.42 Umbilical cord	Germ Cell	Fore skin
15294	35147	R45550	10.94	67.12	6.14	0.00	1.00	349.34 Skin	Cervix	Bone
15317	813719	AA453779	12.45	1484.63	117.83	11.00	0.00	554.03 Brain	LID not found	Other
15318	35812	R45627	242.87	1224.71	5.04	0.00	0.00	Muscle	Heart	Testis
15323	768421	AA495835	21.43	159.27	7.43	0.00	5.00	246.7 Brain	Lung	LID not found
15334	35728	R45692	4.13	25.72	8.23	1.00	0.00	Pool	LID not found	Other
15335	1499234	AA865728	11.33	137.74	12.15	1.00	0.00	633.69 Lung	Brain	LID not found
15337	812172	AA456036	8.01	50.35	6.29	1.00	0.00	83.98 Skin	Germ Cell	Colon
15341	813748	AA453802	0.73	8.15	11.15	1.00	0.00	Tonsil	Testis	Pool
15345	812175	AA456044	8.40	73.50	7.82	1.00	4.00	Fore skin	Pool	LID not found
15348	214008	H70775	16.72	138.19	8.28	3.00	0.00	653.71 Cervix	Placenta	Pancreas
15362	626348	AA186553	11.81	225.11	18.06	3.00	0.00	Fore skin	LID not found	Other
15363	268697	N22897	248.28	1302.54	5.25	2.00	0.00	174.53 Umbilical cord	Ear	Thymus
15372	300024	N78903	16.17	81.94	5.07	0.00	1.00	102.24 Parathyroid	LID not found	Other
15376	321958	W27733	122.09	1402.23	11.48	0.00	2.00	245.06 CNS	LID not found	Other
15386	289168	N68970	56.84	482.61	8.14	0.00	4.00	Lung	Tonsil	LID not found
15394	841670	AA487563	14.61	95.54	6.54	0.00	0.00	481.83 Lymph nodes	Head and nec	CNS
15402	627401	AA190825	6.03	115.96	19.22	5.00	0.00	Nose	Aorta	Heart
15403	254562	N23651	9.45	74.55	7.89	1.00	0.00	157.82 Fore skin	Spleen	Adrenal gland
15404	291891	N73011	28.77	302.80	10.53	0.00	1.00	474.75	Whole embryo	LID not found
15408	322926	W45025	480.32	2443.33	5.00	0.00	1.00	460.41 Placenta	CNS	Bone
15430	488499	AA047462	34.02	205.07	6.03	0.00	1.00	Larynx	Pancreas	Blood
15432	592771	AA159994	8.15	37.15	6.04	1.00	0.00	670.02 CNS	Pool	LID not found
15439	282475	N49850	31.28	235.84	7.54	0.00	4.00	CNS	Ovary	Heart
15440	323268	W42746	6.70	34.87	5.20	0.00	1.00	Parathyroid	LID not found	Other
15444	305461	N89812	9.73	54.67	5.62	0.00	0.00	281.02 Nose		
15451	254694	N25049	42.24	252.71	5.98	2.00	0.00			

Table 2A

15455	280529	N90001	141.95	761.82	5.37	0.00	1.00	1	142.28	CNS	LID not found	Other
15459	503583	AA131240	6.48	32.95	5.10	0.00	1.00	4	201.92	Adrenal gland	Whole embryo	Parathyroid
15468	950988	AA598402	109.03	774.36	7.10	2.00	0.00		Colon	Pancreas	LID not found	Other
15470	730150	AA412497	27.98	416.52	14.89	2.00	0.00	2	96.76	Testis	LID not found	Other
15471	839609	AA487301	88.86	665.29	7.49	0.00	2.00		Eye	Esophagus	Umbilical cord	
15472	135811	R33363	25.44	164.26	6.46	1.00	0.00	10	254.83	Uterus	LID not found	Other
15475	503675	AA131450	20.94	212.25	10.14	5.00	2.00	18	22.88	Lung	LID not found	Other
15477	565683	AA133385	5.88	58.67	10.14	2.00	0.00	5	501.96	Colon	LID not found	Other
15481	509463	AA058383	205.97	1587.83	7.71	0.00	2.00		Kidney	LID not found	Other	
15484	897773	AA598515	28.51	205.12	7.19	0.00	2.00		CNS	Heart		
15487	344108	W73781	17.91	131.83	7.37	0.00	1.00	11	373.42	Forebrain	Whole embryo	
15503	839048	AA487505	9.36	69.52	7.41	2.00	0.00	X	245.06	Parathyroid	Pancreas	Eye
15509	512116	AA133590	17.48	389.56	22.29	6.00	0.00	9	358.18	Pool	Brain	LID not found
15511	328813	W45285	1.66	9.97	5.95	1.00	1.00		Adipose	Pancreas	Blood	
15524	123735	R01179	136.20	1538.66	11.30	0.00	3.00		Whole embryo	Colon	LID not found	Other
15531	591116	AA158352	35.03	255.97	7.31	2.00	0.00	X	350.62	Pool	LID not found	Other
15533	784216	AA446866	0.52	2.97	5.74	1.00	0.00		Brain	Parathyroid	Ovary	
15536	206666	H60898	10.24	75.93	7.42	1.00	0.00	18	247.79	Forebrain	Brain	Breast
15547	593638	AA166743	1276.25	6957.63	6.48	1.00	0.00	6	117.83	Pool	LID not found	Other
15568	37671	R61377	32.33	209.42	6.48	1.00	1.00		Kidney	Heart	LID not found	Other
15569	292399	N68398	28.25	228.05	8.11	0.00	3.00	1	174.53	Placenta	LID not found	Other
15573	138281	R68013	47.38	255.92	5.40	0.00	0.00		Pool	Testis	Pool	
15577	201422	R98092	13.58	96.24	7.08	0.00	0.00	6	137.8	Uterus	LID not found	Other
15578	416309	W68185	23.55	125.82	5.54	1.00	0.00		Brain	Thymus	Adrenal gland	
15585	505341	AA156235	154.72	856.76	5.34	0.00	2.00	20	309.22	Brain	Lung	Breast
15592	30459	R42143	3.57	25.47	7.13	1.00	0.00	11	367.35	Pool	LID not found	Other
15593	241705	H81680	61.33	576.55	5.39	1.00	0.00	20	12.07	Nose	Adrenal gland	Brain
15600	44303	H06390	6.73	36.27	5.08	1.00	0.00	17	375.84	Pool	LID not found	Other
15613	195888	R61401	835.50	4242.62	8.06	1.00	0.00		Testis	LID not found	Other	
15617	123326	R00311	46.42	442.67	9.54	0.00	4.00	1	620.01	Lung	LID not found	Other
15623	666492	AA233070	21.53	173.65	8.62	0.00	0.00		Testis	LID not found	Other	
15625	415165	W93407	158.31	1364.45	8.82	0.00	3.00	18	281.83	Pool	LID not found	Other
15626	416557	W68992	41.08	320.75	7.81	0.00	0.00	X	245.06	Pool	LID not found	Other
15627	1049321	AA620783	1.42	8.88	6.24	1.00	0.00	10	104.78	Lung	Adrenal gland	Eye
15630	309224	N93653	47.04	300.26	6.38	0.00	1.00	1	639.73	CNS	Colon	Blood
15635	1055487	AA620784	13.77	107.83	7.63	1.00	5.00	2	645.89	Larynx	Ear	Cervix
15638	753038	AA438460	11.07	69.74	6.82	1.00	0.00		Testis	LID not found	Other	
15637	200417	R97240	389.23	2400.78	8.17	2.00	1.00		Brain	LID not found	Other	
15645	201551	R97970	1160.17	5935.08	5.12	1.00	0.00	12	42.69	Esophagus	Spleen	Parathyroid
15646	306821	N95041	93.00	592.11	6.26	0.00	3.00	X	358.17	Pool	Germ Cell	Pool
15675	42415	R60981	10.38	95.92	9.24	0.00	1.00	20	247.78	Skin	Adrenal gland	Bone
15678	1473690	AA916728	7.54	77.58	10.29	5.00	1.00	4	671.55	CNS	Testis	Eye
15689	787236	AA424675	42.85	235.28	5.49	1.00	0.00		CNS	LID not found	Other	
15695	648046	AA206914	235.50	1202.45	5.11	1.00	0.00		Cervix	Parathyroid	Lung	
15697	787239	AA418402	22.33	168.50	7.46	0.00	4.00		Testis	Brain	Pool	
15699	43065	R61886	2.15	23.93	11.11	3.00	0.00		Aorta	CNS	Germ Cell	Pool
15713	787282	AA418382	13.06	457.14	35.00	1.00	0.00	12	42.69	Esophagus	Spleen	Parathyroid
15723	42872	R61883	4.67	38.88	8.32	1.00	1.00		CNS	LID not found	Other	
15726	786308	AA451863	9.52	117.89	12.38	6.00	0.00	20	247.78	Skin	Adrenal gland	Bone
15729	787273	AA418408	8.82	102.85	11.66	0.00	1.00		Lung	Testis	LID not found	Other
15731	43080	R61289	5.09	38.68	7.59	0.00	0.00	4	671.55	CNS	Testis	Eye
15735	587055	AA131315	20.35	200.56	9.66	4.00	0.00		CNS	LID not found	Other	
15738	726821	AA388341	6.22	36.00	5.79	0.00	1.00					
15783	281545	N51601	50.58	290.69	5.75	0.00	1.00					

Table 2A

15768	429447	AA007826	53.69	301.95	5.62	1.00	0.00	Pool	LID not found	Other
15778	568421	AA148858	46.20	281.97	6.32	0.00	3.00	Colon	Colon	Placenta
15788	322028	W37424	39.95	282.31	7.07	0.00	1.00	Parathyroid	Breast	LID not found
15803	282481	N52039	133.96	802.08	5.99	0.00	2.00	CNS	LID not found	Other
15814	220059	H85434	81.15	500.26	6.16	0.00	0.00	Eye	LID not found	Other
15822	219361	H84785	15.29	101.73	8.65	2.00	1.00	Eye	Germ Cell	
15827	284457	N52337	21.06	250.96	11.91	2.00	5.00	CNS	LID not found	Other
15831	284545	N64762	7.10	36.46	5.13	1.00	0.00	CNS	LID not found	Other
15832	430313	AA010611	5.01	44.59	8.89	3.00	0.00	Pool	LID not found	Other
15843	243602	N49717	21.65	125.16	5.73	2.00	0.00	Whole embryo	Pool	LID not found
15846	797042	AA463221	10.34	105.79	10.23	1.00	0.00	Muscle	Whole embryo	LID not found
15849	611269	AA174508	8.63	147.49	17.10	0.00	1.00	Whole embryo	LID not found	
15850	784010	AA443695	153.08	890.95	5.82	1.00	2.00	Blood	Whole embryo	LID not found
15862	277509	N56892	21.67	124.24	5.73	0.00	1.00	45 Kidney	CNS	Tonsil
15865	612613	AA179392	160.45	1487.32	9.15	2.00	0.00	Heart	Heart	LID not found
15870	610112	AA464972	40.35	251.92	6.24	0.00	2.00	Parathyroid	Ovary	LID not found
15880	1031592	AA609473	16.60	86.22	5.32	1.00	0.00	Pooled	Blood	Muscle
15882	638003	AA434482	360.37	2332.06	6.47	2.00	0.00	-2.94 Eye	LID not found	Other
15884	255182	N22033	7.28	121.39	18.68	1.00	0.00	Adrenal gland	Pooled	Nose
15893	120528	T85320	50.76	667.01	13.14	0.00	5.00	Pool	LID not found	Other
15897	322588	W15284	287.42	1721.92	5.99	0.00	2.00	Parathyroid	Colon	LID not found
15901	757180	AA443658	16.59	107.12	6.46	0.00	0.00	Thymus	Synovial mem	Cervix
15903	627288	AA191483	16.80	139.16	8.28	0.00	2.00	Cervix	Eye	LID not found
15906	730732	N48593	21.89	146.53	6.70	1.00	0.00	Skin	Pooled	Kidney
15919	796159	AA461091	16.15	86.33	5.35	1.00	0.00	Adrenal gland	CNS	Eye
15920	1049185	AA620697	7.39	185.46	25.10	0.00	0.00	Whole embryo	Testis	Tonsil
15921	743568	AA609463	70.39	406.39	5.77	0.00	1.00	Bone	Muscle	Testis
15922	757365	AA437124	116.62	937.30	8.04	0.00	4.00	217.43 Thyroid	Forearm	LID not found
15924	122068	T98287	8.44	48.09	5.46	1.00	1.00	219.18 Testis	Prostate	LID not found
15926	753794	AA410383	5.29	35.15	6.65	2.00	0.00	Pool	LID not found	Other
15927	511852	AA100874	485.63	9928.18	21.31	6.00	0.00	437.61 Spleen	Germ Cell	Lymph
15928	1056172	AA620995	21.52	129.21	6.00	3.00	0.00	Colon	Brain	LID not found
15931	279938	N57535	5.02	77.38	15.43	4.00	0.00	144.49 -	Kidney	Blood
15934	843278	AA486558	301.86	1685.70	5.58	0.00	2.00	610.42 Gall bladder	Colon	CNS
15936	1048714	AA620628	89.47	578.60	6.47	0.00	3.00	Adipose	Cervix	Adrenal gland
15940	665029	AA193579	13.56	189.67	13.99	6.00	0.00	Lung	Testis	LID not found
15943	767456	AA417994	15.03	89.67	5.97	1.00	0.00	244.38 Blood	CNS	Pool
15950	1033342	AA621381	635.67	3258.35	5.13	1.00	0.00	438.03 Spleen	Placenta	Uterus
15952	753957	AA479382	42.42	485.61	11.68	3.00	0.00	133.06	CNS	
15957	247962	N58276	85.25	489.92	5.75	2.00	0.00	355.29 Pool	LID not found	Other
15977	276441	N40211	63.28	377.14	5.96	0.00	2.00	Smooth musc	CNS	LID not found
15982	1030618	AA608824	10.75	68.83	6.40	1.00	0.00	Testis	Pool	LID not found
15986	282837	N68578	14.17	94.75	6.69	0.00	4.00	Pool	LID not found	Other
15993	214179	H17614	35.90	275.50	7.67	0.00	3.00	Pool	LID not found	Other
16008	753993	AA478962	13.41	71.01	5.30	0.00	1.00	Ear	Pooled	Parathyroid
16019	786213	AA453435	47.10	252.45	5.36	0.00	2.00	Stomach	Uterus	Aorta
16031	788526	AA462802	28.41	196.88	7.45	0.00	1.00	Parathyroid	Whole embryo	Colon
16046	35681	R45970	3.85	20.22	5.26	1.00	0.00	252.96		
16055	1470048	AA865464	37.10	244.02	6.58	1.00	0.00	Ovary	Pancreas	Stomach
16058	26759	R36505	7.25	491.44	67.78	2.00	0.00	295.78 Whole embryo	Kidney	Brain
16062	35788	R45976	3.85	31.66	8.03	0.00	1.00	88.99 Colon	CNS	Testis
16065	612242	AA455041	25.90	137.42	5.37	0.00	1.00	Pool	LID not found	Other
16067	766590	AA425116	24.97	146.71	5.68	2.00	0.00	Colon	Testis	Brain
16071	1475859	AA872020	8.71	199.68	29.75	14.00	1.00	Pooled	Pancreas	Stomach

Table 2A

15072	435934	AA701844	9.42	63.64	6.76	1.00	0.00	20	197.2	Parathyroid	Placenta	Fore skin
15079	1470530	AA864524	29.34	148.49	5.06	1.00	0.00	4	499.2	Brain	Testis	LID not found
16086	35804	R48000	4.09	158.41	38.75	0.00	1.00	7	511.72	Pancreas	Brain	LID not found
16087	1475595	AA873885	6.78	135.95	20.05	8.00	0.00	2	182.89	Blood	Tonsil	Fore skin
16080	27277	R37487	1.43	7.40	5.18	1.00	1.00					
16117	823575	AA497044	25.90	134.32	5.19	0.00	1.00					
16120	447187	AA702873	9.06	48.14	5.31	1.00	0.00					
16123	53331	R15922	1.56	10.69	6.66	1.00	0.00					
16128	450152	AA703449	7.44	51.38	6.91	1.00	0.00	17	374.6	Ear	Colon	Blood
16130	595697	AA187382	9.76	58.78	6.02	1.00	0.00	16	333.64	Nose	LID not found	Diher
16131	254749	N25085	231.79	1569.52	6.77	2.00	0.00				Placenta	Lymph
16134	510688	AA101954	7.60	39.55	5.20	1.00	0.00	11	45.93	CNS	Whole embryo	Pool
16147	784278	AA447480	14.22	171.58	12.07	1.00	1.00				LID not found	Other
16148	305556	N90218	15.77	81.95	5.20	0.00	1.00				Tonsil	Heart
16152	323806	W46341	34.87	184.98	5.30	1.00	0.00	6	111.13	Thymus	Germ Cell	Parathyroid
16174	568597	AA149987	7.57	68.85	9.23	2.00	0.00				Pool	LID not found
16176	324322	W47416	7.96	59.11	7.43	0.00	1.00				Parathyroid	Other
16183	279484	N48764	168.77	1218.43	7.22	2.00	0.00	20	11.91	Pooled	Ovary	Parathyroid
16186	810027	AA455275	7.28	51.92	7.13	3.00	0.00	1	146.01	Aorta	Colon	Whole embryo
16192	324983	W49487	61.73	387.67	6.28	0.00	3.00	6	358.37	Colon	Whole embryo	Pool
16202	588139	AA132172	117.75	750.27	6.37	0.00	4.00	2	703.17	Thymus	Whole embryo	Parathyroid
16203	322005	W37418	27.10	141.67	5.23	1.00	0.00	17	341.05	Pooled	Placenta	Pool
16204	308384	N90704	8.22	181.10	22.03	1.00	0.00	2	402.09	Germ Cell	Eye	Kidney
16208	813237	AA181787	12.54	257.51	20.53	0.00	1.00	5	108.43	Gall bladder	Testis	Kidney
16211	257206	N26899	49.64	328.85	6.82	0.00	2.00				LID not found	Other
16220	303110	N80774	60.91	313.84	5.15	0.00	2.00	6	44.4	Neural	Esophagus	Fore skin
16223	279703	N49005	4.61	40.05	8.68	2.00	0.00	1	678.07	Testis	LID not found	Other
16229	743481	AA809392	5.04	39.05	7.74	1.00	0.00	X	245.08	Testis	LID not found	Other
16231	626848	AA181426	87.81	549.50	5.63	1.00	0.00	1	32.73	Pool	Brain	Heart
16235	842840	AA486277	175.91	1425.30	8.10	2.00	0.00	5	379.83	Pool	LID not found	Other
16238	730858	AA418984	78.23	492.69	6.30	2.00	1.00	17	316.02	Colon	LID not found	Other
16248	214624	H71242	7.15	80.34	11.24	4.00	0.00	8	390.78	Ear	Whole embryo	Ovary
16263	531459	AA074079	213.29	1097.38	5.14	0.00	1.00	1	167.12	Pancreas	Eye	Pool
16264	232912	H72843	68.12	668.70	8.78	2.00	3.00	1	75.41	Parathyroid	Lung	Uterus
16268	587268	AA132650	6.24	40.28	6.45	0.00	0.00	1	263.38	Tonsil	Pool	LID not found
16271	595001	AA184782	18.64	171.97	8.67	2.00	0.00	8	103.36	Testis	LID not found	Other
16279	139304	R63714	155.10	854.15	5.51	2.00	0.00				Fore skin	Whole embryo
16281	592403	AA159497	0.80	11.45	14.28	4.00	0.00	17	102.53	Gall bladder	Adrenal gland	Tonsil
16285	731198	AA417355	3.55	49.29	13.87	6.00	0.00	2	233.5	Pooled	Whole embryo	Uterus
16288	197638	R86198	582.33	3431.32	6.10	1.00	0.00				Lymph	LID not found
16298	743197	AA401438	15.16	181.33	11.96	0.00	4.00				LID not found	Other
16303	357180	W93544	471.01	2915.09	8.19	1.00	0.00	6	39.19	Testis	Pool	LID not found
16304	233627	H78411	49.18	259.02	5.27	2.00	1.00	15	227.19	Pool	LID not found	Other
16305	528116	AA084888	6.82	53.77	6.10	1.00	0.00	19	271.02	Brain	LID not found	Other
16313	796871	AA461490	5.52	31.12	5.64	1.00	0.00	14	130.93	Brain	LID not found	Other
16315	638874	AA481795	129.64	652.97	5.04	1.00	0.00	1	153.98	Pool	LID not found	Other
16316	123585	R00835	44.66	274.86	6.11	0.00	3.00	X	74.75	Breast	Colon	
16317	565647	AA138540	38.74	208.53	5.68	1.00	0.00					
16318	731202	AA417252	4.57	32.74	7.17	1.00	0.00					
16321	68734	T84868	75.71	399.75	5.28	0.00	1.00					
16325	201213	R98471	365.60	2018.86	5.52	0.00	2.00					
16328	323041	W42450	322.08	1735.00	5.39	1.00	0.00					
16336	29583	R42218	208.26	1111.72	5.34	0.00	1.00					
16341	208789	H63241	414.28	2136.82	5.16	1.00	0.00					
16351	43733	H04789	13.44	67.89	5.05	0.00	1.00					

Table 2A

16357	233942	H66122	80.90	487.59	6.03	0.00	2.00	101.02	Pool	LID not found	Other
16375	42331	R61821	87.34	376.58	5.62	2.00	0.00	449.86	Brain	LID not found	Other
16400	35620	R45832	8.48	48.45	5.71	1.00	0.00	Colon	Colon	Pool	Pool
16402	427778	AA002258	67.42	728.00	10.80	1.00	4.00	Peripheral ner	Pool	LID not found	Other
16405	239943	H81938	159.60	842.34	5.28	0.00	2.00	Pool	Pool	LID not found	Other
16408	44092	H06282	14.23	89.83	6.31	1.00	2.00	92.24	Placenta	Pool	Lung
16409	262282	H99362	4.06	25.85	6.37	0.00	1.00	714.07	Foreskin	Pool	Kidney
16415	42803	R60014	3.64	20.46	5.62	2.00	0.00	592.03	Peripheral ner	Adipose	Brain
16421	789596	AA425056	9.20	36.99	6.41	2.00	0.00	421.53	Heart	Pool	LID not found
16430	786534	AA452118	18.17	111.38	6.13	0.00	1.00	Eye	Whole embryo	LID not found	Other
16435	42807	R60135	3.29	19.33	5.88	1.00	0.00	Brain	LID not found	Other	Other
16438	786537	AA452130	149.27	1224.89	8.21	0.00	2.00	485.68	Pooled	Whole embryo	Prostate
16439	888204	AA598594	7.74	50.13	6.48	2.00	0.00	Kidney	Brain	Pool	Pool
16454	786545	AA452134	23.63	129.22	5.47	1.00	0.00	39.72	Thyroid	Foreskin	Parathyroid
16467	42816	R60044	17.05	97.62	5.73	1.00	0.00	143.02	Thyroid	Lung	Pool
16472	1367900	AA810225	2.83	17.38	8.14	1.00	0.00	57.93	Ear	Pooled	Placenta
16478	786593	AA424820	8.45	64.20	6.41	1.00	0.00	Lung	Whole embryo	Brain	Brain
16478	786590	AA452250	6.91	67.80	9.81	1.00	0.00	201.74	Foreskin	Spleen	Brain
16479	898227	AA598625	11.78	63.87	5.42	0.00	1.00	Testis	LID not found	Other	Other
16482	726858	AA398355	38.84	241.60	6.56	1.00	0.00	Tonsil	Pool	LID not found	Other
16485	786661	AA425543	5.19	31.53	6.08	1.00	0.00	282.85	Ear	CNS	Bone
16489	37880	R81372	1.88	148.03	87.95	4.00	4.00	Kidney	Lung	Pool	Pool
16500	111489	T90789	2.88	15.78	5.48	1.00	0.00	Testis	LID not found	Other	Other
16508	728889	AA398430	5.67	78.61	13.38	2.00	0.00	Cervix	Spleen	Parathyroid	Parathyroid
16511	898259	AA598679	51.55	288.40	5.58	0.00	1.00	Whole embryo	Colon	Foreskin	Foreskin
16538	786330	AA481318	81.82	564.13	6.89	0.00	2.00	Pancreas	Lung	Pool	Pool
16547	247840	NS3670	105.64	813.61	7.70	1.00	3.00	42.61	Skin	Pooled	Foreskin
16562	786117	AA460981	7.15	68.47	9.59	2.00	0.00	443.42	Pool	LID not found	Other
16588	416644	W86445	42.56	358.85	8.48	0.00	4.00	144.58	Pancreas	Testis	Whole embryo
16611	292531	N62852	8.18	70.07	8.58	1.00	0.00	Colon	Prostate	Heart	Heart
16612	510578	AA055788	6.41	257.50	47.84	5.00	0.00	Smooth musc	Testis	Testis	LID not found
16615	511080	AA100283	7.37	48.12	6.28	1.00	0.00	206.24	Pool	LID not found	Other
16618	1031698	AA609556	37.88	270.79	7.15	1.00	0.00	Testis	Tonsil	Colon	Colon
16619	292122	N62434	14.00	145.88	10.42	3.00	0.00	191.7	Kidney	Colon	Brain
16620	728956	AA412047	6.17	32.75	5.30	1.00	0.00	875.72	CNS	Heart	Pool
16621	511117	AA088231	7.72	74.17	8.60	1.00	0.00	Skin	Stomach	Blood	Pool
16634	784032	AA443712	13.38	74.39	5.58	1.00	0.00	105.6	Lung	Breast	Pool
16638	240148	H82435	8.91	58.86	6.61	3.00	0.00	Pool	Heart	LID not found	Other
16637	810217	AA464698	51.18	448.96	8.77	0.00	1.00	427.01	Eye	Thyroid	Parathyroid
16638	841621	AA487468	10.19	72.54	7.12	4.00	0.00	Blood	Eye	Ovary	Ovary
16644	415712	W84858	30.98	354.57	11.44	5.00	0.00	381.52	Pool	LID not found	Other
16648	1030849	AA608775	55.61	333.70	6.00	1.00	3.00	Marrow	Spleen	Parathyroid	Kidney
16651	840837	AA486538	8.58	44.26	5.18	1.00	0.00	15.78	Pooled	Breast	Colon
16655	770654	AA427737	15.85	106.92	6.75	2.00	0.00	49.08	Pool	LID not found	Other
16659	283058	N83777	639.84	3692.85	5.77	2.00	0.00	545.43	CNS	Eye	Lung
16669	786723	AA443140	8.19	99.35	12.13	8.00	0.00	201.71	Germ Cell	Pancreas	Blood
16678	897587	AA496894	8.88	116.50	13.12	2.00	0.00	285.87	Small intestine	Aorta	Placenta
16679	503234	AA151621	5.20	42.86	8.24	3.00	0.00	391.77	CNS	LID not found	Other
16684	124139	R01248	15.32	133.38	8.70	2.00	4.00	340.75	Thyroid	Stomach	Parathyroid
16685	289019	N82735	4.51	27.79	8.18	0.00	1.00	463.73	Thyroid	Stomach	Parathyroid
16689	131452	R23270	25.24	176.24	6.98	0.00	1.00	463.73	Thyroid	Stomach	Parathyroid
16701	757244	AA426025	12.54	336.28	26.82	3.00	0.00	463.73	Thyroid	Stomach	Parathyroid
16703	197056	R92801	17.51	224.18	12.80	0.00	1.00	463.73	Thyroid	Stomach	Parathyroid
16713	278308	N48353	240.04	1299.80	5.42	1.00	2.00	463.73	Thyroid	Stomach	Parathyroid
16723	813730	AA453783	14.01	862.63	47.31	17.00	2.00	463.73	Thyroid	Stomach	Parathyroid

Table 2A

16728	1030543	AA608923	68.21	438.99	6.45	0.00	3.00	Testis	LID not found	Other
16729	2791178	NA68845	10.34	183.72	17.76	1.00	0.00	CNS	LID not found	Other
16730	627888	AA197344	86.40	488.88	5.67	1.00	0.00	Cervix	LID not found	Other
16738	754101	AA478818	8.96	118.83	13.25	2.00	0.00	Adipose	Heart	Whole embryo
16757	276168	N63528	10.90	63.65	5.91	1.00	2.00	CNS	LID not found	Other
16768	754126	AA476623	6.38	74.13	8.84	2.00	0.00	620.52 Ear	Pancreas	Bone
16783	39189	R54443	11.88	259.87	22.23	1.00	0.00	501.89 Brain	Prostate	-
16788	836888	AA481801	51.94	533.77	10.28	2.00	5.00	Eye	LID not found	Other
16787	44081	H08266	12.14	100.01	8.24	0.00	3.00	75.83 Whole embryo	Testis	Pool
16789	293975	N64024	53.28	351.08	6.59	0.00	4.00	337.98 Pool	LID not found	Other
16804	180179	R65537	4.74	24.86	5.27	1.00	0.00	210.72 Pancreas	Lymph	Heart
16808	37205	R49592	2.20	12.08	5.51	0.00	1.00	Brain	LID not found	Other
16810	28705	R40357	420.29	2363.88	5.82	0.00	2.00	Neural	Whole embryo	LID not found
16818	29030	R40377	35.87	255.78	7.17	2.00	0.00	127.28 Thyroid	Bone	Whole embryo
16819	53031	R16946	38.42	214.31	5.88	1.00	2.00	178.58 Brain	LID not found	Other
16824	448514	AA777637	7.53	45.84	6.06	2.00	0.00	74.01 Whole embryo	Brain	Breast
16828	161172	H25223	6.30	33.03	5.25	0.00	1.00	155.35 Kidney	Heart	Pool
16831	1323328	AA872602	2.96	37.58	12.70	2.00	0.00	84.38 Placenta	Breast	Tonsil
16837	823815	AA486357	10.81	56.92	5.27	1.00	0.00	Brain	LID not found	Other
16846	450060	AA703392	81.82	486.00	5.95	0.00	1.00	Kidney	Whole embryo	LID not found
16850	28658	R40987	8.55	82.56	9.66	0.00	0.00	Brain	LID not found	Other
16858	29077	R40983	147.22	747.28	5.08	0.00	1.00	Brain	Whole embryo	Pool
16858	53158	R16148	3.83	31.36	8.18	1.00	1.00	Kidney	Foreskin	Pool
16863	1323448	AA873804	86.53	527.84	7.83	7.00	1.00	50.5 CNS	Blood	Lung
16868	172751	H19687	4.35	23.84	5.48	3.00	0.00	238.05 Germ Cell	Brain	Breast
16870	37505	R51103	6.14	65.89	10.73	1.00	0.00	117.35 Brain	LID not found	Other
16878	1323539	AA955296	4.23	28.96	6.85	2.00	0.00	Kidney	Pool	LID not found
16880	450710	AA704448	66.32	428.08	6.45	0.00	1.00	Brain	LID not found	Other
16882	29251	R41389	11.82	203.87	17.54	1.00	0.00	250.83 Brain	LID not found	Other
16883	53081	R18144	5.27	34.96	6.64	0.00	1.00	Lung	Brain	Brain
16884	221828	H92234	2.34	23.32	8.98	3.00	0.00	742.57 Brain	Germ Cell	Eye
16885	823647	AA486884	10.07	169.12	16.80	7.00	1.00	287.88 CNS	Pooled	Whole embryo
16887	1323591	AA858028	18.86	697.83	35.12	0.00	0.00	251 Adrenal gland	Germ Cell	Gall bladder
16889	812871	AA464603	21.81	122.89	5.61	1.00	0.00	259.71 Pooled	CNS	Parathyroid
16892	380392	AA013268	71.98	531.97	7.39	0.00	4.00	697.77 Eye	Kidney	Testis
16894	38381	R49033	5.92	42.44	7.17	1.00	0.00	Brain	LID not found	Other
16896	450745	AA704492	12.94	80.00	8.18	0.00	1.00	527.46 Esophagus	Pool	LID not found
16907	257826	N27026	464.82	2447.39	5.05	1.00	0.00	154.34 Neural	Brain	Eye
16910	839862	AA480044	11.03	70.28	6.37	3.00	0.00	538.58 Synovial mem	Uterus	-
16918	729924	AA386833	22.10	113.58	5.14	1.00	0.00	Heart	Testis	LID not found
16928	340904	V57787	4.41	28.83	6.54	1.00	0.00	741.96 Testis	Pancreas	Colon
16934	731348	AA421018	6.30	74.38	11.81	1.00	0.00	CNS	Testis	LID not found
16943	731275	AA420998	316.25	1982.87	6.27	2.00	0.00	599.03	Placenta	Muscle
16956	307995	N92293	208.70	1483.93	7.11	0.00	2.00	64.16 Pooled	Whole embryo	Parathyroid
16962	595161	AA173408	24.82	212.05	8.51	2.00	0.00	99.87 Synovial mem	Parathyroid	Testis
16971	772982	AA478258	110.44	584.25	5.29	1.00	0.00	92.4 Colon	Testis	-
16978	259827	N29776	9.89	52.38	5.30	1.00	0.00	191.31 CNS	Testis	Kidney
16983	280308	AA17075	6.34	275.05	43.41	0.00	0.00	Eye	LID not found	Other
16997	565949	AA136541	71.25	489.49	5.59	0.00	2.00	Pool	Heart	LID not found
17004	127711	R09504	7.00	40.11	5.73	0.00	2.00	Testis	Kidney	LID not found
17010	742586	AA401347	14.68	80.49	5.48	1.00	0.00	Pool	LID not found	Other
17020	128810	R07268	151.49	1000.79	6.81	2.00	1.00	369.24 Heart	LID not found	Other
17025	360035	AA063577	97.68	497.57	5.09	1.00	0.00	Testis	LID not found	Other
17030	731193	AA417356	4.29	26.97	6.28	1.00	0.00	31.88 Pool	LID not found	Other
17032	243024	H95669	31.81	186.84	5.90	1.00	1.00			